OPU NEEDLE

KITAZATO®

Copyright © KITAZATO CORPORATION. All Rights Reserved.

Kitazato Corporation 100-10 Yanagishima, Fuji, Shizuoka 416-0932 Japan
TEL: +81-545-65-7122 FAX: +81-545-65-7128

EC REP Dibimed - Biomedical Supply, S. L. C./ Jorge Comín, 3 Bajo 1-2, 46015 Valencia, Spain
TEL: +34-96-305-63-95 FAX: +34-96-305-63-96

CE 2797

EN OPU NEEDLE

GENERAL INFORMATION: Refer to the label on the product or package before use.

INTENDED USE: The Kitazato OPU needle is used for puncturing into human body and picking up tissue follicular fluid or oocyte from ovarian follicle.

CONTRAINDICATIONS: • No reuse. It may cause infection when reuse. • No resterilization.

INSTRUCTIONS FOR USE: (1) Apply ultrasound probe to the puncture region and check the condition of the follicle. For transvaginal puncture, make sure of the adjustment for the transvaginal ultrasound probe guide and the guage of the needle. (2) Connect the needle hub or the female connector hookup to the syringe or the suction pump and make sure the direction to the follicle looking at the monitor of ultrasound. (3) Insert the needle tip to the follicle and aspirate the follicular fluid slowly. (4) For puncturing another follicle in the same ovary, keep the needle tip in the ovary. It is recommended to do flushing with the 10ml syringe accompanying the Double Lumen of OPU needle. Conduct the flushing slowly. (5) Remove the needle after aspirating the follicular fluid. (6) Lightly press the puncture region with gauze. Apply an adhesive plaster to the puncture region if bleeding cannot be observed. Check the condition of the patient for 15 to 30 minutes. (7) Immediately do speculum or store for the aspirated and cell.

CAUTIONS FOR USE: (1) Confirm the specification (gauge and length) of the needle before opening it. (2) Confirm adaptability with the other medical equipment to use together beforehand. (3) Pay enough attention to the needle tip and do not touch it until it is used for puncture. An efficacy of the puncture may deteriorate. (4) Confirm a route of the puncture enough before puncturing the needle. (5) When the clogging in the needle or the connecting tube is confirmed on aspiration of tissues, remove the clogging from it by pressurizing with the cannula attached to the product. However, avoid the enforced handling, and, depending on the situation, change a needle newly. (6) Read an Instructions for Use before use. (7) Do not use the product after expiration date. (8) Do not use if the package is opened or damaged. (9) This product is intended to be used by medical specialist trained in fertility treatment. (10) Do not use the product except the collection of the tissues and the cells. (11) Use the product aseptically. Carry out enough precautionary measures for bacteria. (12) Confirm leakage through solution such as saline before use. Do not use the product when air bubble is confirmed by through solution such as saline. (13) Do not recap after use. (14) Pay enough attention to a false puncture. (15) Throw the product away promptly after having used it. (16) Do not refurbish. (17) Use immediately after opening. (18) Store away from high temperature, high humidity.

WARRANTY: The manufacturer is not liable for any damage resulting from misuse or incorrect operation or failure to follow the cautions or instructions contained in this instruction. For Use or the instructions manual for any compatible medical equipment.

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