MEGAGEN

Surgical Instrument

www.imegagen.com

MegaGen Implant Co., Ltd
45, Sechen-ro 7-gil, Dasa-eup, Dalseong-gun, Daegu 42921, Republic of Korea
Tel. 82-53-222-2828 / Fax. 82-53-289-3420

EC REP Europe representative : MDSS GmbH
Schiffgraben 41 30175 Hannover, Germany

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CE 2460 ● CE conformity marking

REF ● Catalogue number | LOT ● Batch code

● Date of manufacture | SN ● Serial number

NON STERILE ● Non-sterile

● Temperature limitation

● Caution, consult accompanying documents

● Do not use if package is damaged

● Manufacturer

● Consult instructions for use

EC REP ● Authorised representative in the european community.

Rx Only ● Caution: Federal law(USA) restricts this device to sale by or on the order of a licensed dentist

IFU-II-R10/ 20-08-28

English

1. Product Summary

To recover a patient's masticatory function, it is necessary to have a dental implant which is implanted into a patient's alveolar bone as a root of a tooth and abutment which is protruded through the gum in the oral cavity. These products include various Instruments to support implant surgery.

2. Intended Use

These products are surgical instruments designed for the surgical placement of dental implants or prosthesis. It typically contain various dental drills, drill extensions, depth gauges, wrenches, drivers, osteotomes, path finders, and direction indicators designed in various types and sizes for use.

※ Surgical insertion instruments for temporary use(like Drill etc.)
→ All surgically invasive devices intended for transient use
Surgically invasive devices that perform for creating hard tissue/ soft tissue shapes suitable for implant surgery, like to drilling holes on the bone cutting tissue.

Product lines are as follows:
- Lance drill
- Tissue Punch
- Marking Drill
- Shaping Drill
- Stopper Drill
 etc.

※ The Instruments that is not inserted surgically(like Connector, Driver etc.)
→ only to assist in the procedure
Surgical non-invasive devices that perform to help for implant surgery and don't modify hard tissue/ soft tissue.

Product lines are as follows:
- Handpiece Connector
- Ratchet Connector
- Implant Carrier
- Hand Driver
- Octa Driver
- Torque Driver
- Torque Wrench
 etc.

3. Pre-operative planning

For a successful treatment, the following confirmation procedures should be performed carefully to confirm the treatment plan and evaluation, and also the product with medical and dentistry's point of view.
(1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.
(2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.
(3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.
(4) The product is supplied as non-sterile devices. All product used for surgery must be sterilized before use.
(5)Choose the surgical instrument after considering right dimension at surgical site.
(6)Depth Gauge is needed to confirm that the product is the same as shown in the catalog. In particular, it is recommended to use after checking if laser marking is in the same location.
* Marking can easily be erased when the product is left for a long time after the procedure, so caution is required.
(7) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, expired date, and etc. after opening the package.
(8) Check the product information such as reference code, size, manufacturing date and expiration date before starting the surgery.
(9) Check for compatibility with the opponent component before surgery.
(10) Check if there is any problem in the shape of the instrument

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5. Adverse Effects

(1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.
(2) Complications such as malaise, pain, edema, bleeding, Inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient overall health.
(3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.

6. Contraindications

Contraindications should be observed when selecting patients. These include, but are not limited to:
(1) Low bone quality and inadequate bone quantity
(2) Bad oral hygiene and pathological infection at the alveolar bone.
(3) Internal diseases such as hematodyscrasia and diabetes, or mal- nutrition
(4) Alcoholic, drug abuse, heavy smoking and a mental patient.
(5) Other periodontal and surgical operation prohibited patient

7. Preservation & Handling

(1) The surgeon should check each product being used during surgery is in right position and whether if the product is contaminated or not.
(2) The surgical instrument should be used for dental implant surgery following procedure as guided.
(3) The edge of surgical instrument for making the bone cavity should be checked before use.
(4) The surgical instrument for making the bone cavity should be used with saline solution to prevent necrosis.
(5) The surgeon should use the surgical instrument with suitable RPM(revolutions per minute) by considering bone condition such as bone density.
(6) It should be used surgical instruments of relevant Implant System.
(7) A small devices such as gauges with measuring function can be entered through throat during surgery. Fixation with thread or wire should be needed.
(8) The surgical instruments are reusable and dental drills should be replaced after being used 50 times in condition of TiN coating and non-coating or 50 times in condition of AlTiN coating or when the cutting efficiency declines. The lifetime of dental drill can be changed by handling method.
(9)The frequency of use could be changed depending on the handling.
(10)It can not be used when the cutting force of Drill is reduced.
(11) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
(12) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.
(13) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
(14) If the package that is not opened by dentist is defective, request a replacement from manufacture.
(15) This Product should be stored at 15℃~30℃ condition.

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8.1. Cleaning

❑Clean surgical instruments separately and go through a special cleaning procedure before sterilization as mentioned below.
1) Remove all remaining pollutants (tissues, blood, bone residuals, etc.) on the surface by immersing used surgical instruments in water. Keep surgical instruments in wet environment until the next step indicated below is initiated.
2-1) Wash surgical instruments manually (conventional method).
a. Immerse surgical instruments in the water mixed with enzyme detergent (e.g. enzyme detergent with a pH level between 6 and 9) or its similar detergents (available to use commercial products).
b. Scrub the surface of surgical instruments until all visible pollutants will be removed with a suitable soft bristled nylon brush (more than 4minutes). Clean minute parts of surgical instruments with cleansing needle.
c. Rinse the surface of surgical instruments with tap water to remove all cleaning solution.
2-2) Automatic washing using ultrasonic cleaner (recommendation method)
a. After putting surgical instruments in a washing container (e.g. plastic basket) inside of an ultrasonic cleaner, fill water completely. [You'll get better effects when using enzyme detergent (e.g. enzyme detergent with a pH level between 6 and 9) or its similar detergents (available to use commercial products.)]
b. Clean them for 5~15 minutes and drain the water.
c. Rinse the surface of the surgical instruments with tap water to remove all cleaning solution.
❑Repeat step a, b, c until surgical instruments are completely clean and perform a visual inspection.
❑Dry the surface of surgical instruments by natural drying or according to a drying cycle of cleaner. If needed, manual drying using nonabrasive tower with no lint is possible.
❑For more information about available enzyme detergents and cleaner, refer to instructions from the manufacturer and have to be followed.
❑Wash surgical devices separately and go through a special washing procedure as below before sterilization.

8.2. Sterilization

❑Surgical Instruments should be sterilized at temperature of 132℃(270°F) for 15 minutes by high pressure steam sterilization and dry for 30 minutes.
❑In case of more than 149℃(300.2°F), a kit box may be deformed so it should be careful.

9. Caution ⚠

■ Product should not be used if packaging is opened or damaged.
■ Do not use the contaminated product by mistake of surgeon during the surgery.
■ The product should not be reused products between patients without sterilization.
■ If the instrument is contaminated, even a single patient should not use it in another part.
■ Immediate action should be taken in the event of an accident.
■ Quality and accuracy cannot be guaranteed if the customer (user) uses similar products not manufactured or sold by us.
■ If the direction of fixation between the fixture and the connector is not correct, it may be deformed.
■ If the recommended number of uses is exceeded, product deformation or fractures may be indicated.
■ User shall state that the product should not be used if the packaging is not completely sealed or torn.
■ Patients aged 18 and over who have completed the growth of the jaw bone (Partial or complete edentate patient)
■Manufacturer is not responsible for any damages to the product from sterilization or improper handling other than the high pressure steam sterilization method suggested by.
* Problems can occur if not used in the recommended way by manufacture.
* Guide should be used when using Removal Tool, if don't use guide, tool could be fracture.
* Changing the instrument's path during use may cause deformation of the product.

10. Warning

The following is strongly to be recommended for safe and effective use of MegaGen Implant System.
(1) The surgeon should be a medical specialist who practiced high technique for implant surgery.
(2)The surgeon should be fully aware of the operating method and cautions for the product before the use.
(3) The surgeon should wear an aseptic operating gown, and perform surgery in the operating room remained the aseptic condition during surgery.
(4) Inappropriate selection of patient and technique can cause a failure of the implant surgery and loss of supportive bone.
(5) Patients poor oral hygiene management can lead to side effects.
(6) Careless animal rearing in the dentist can contaminate the instrument.
(7)The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
(8) The product must be used in clean and pleasant environment area.

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CE
2460

www.imegagen.com

Cleaning and Sterilization of Surgical and Prosthetic Instrument Kits

MEGAGEN

1. Product Summary

Care and maintenance of reusable instruments are crucial for a successful treatment. Well-maintained instruments are not only safeguard your patients and staff against infection, but also essential for the outcome of the total treatment.

● Please follow the legal regulations in your country as well as hygienic guidelines, to prevent from contaminations and infections by handling the used surgical kits.

● Please remember that you are responsible for the maintenance and sterility of your medical / dental products and device. It is essential to follow proper cleaning, disinfection and sterilization procedures all the time strictly.

● It is also important to follow the manufacturer's recommendation on the usage of instruments. Please keep a log of how many times the instruments have been used.

● Drill's life time is determined by various factors of the surgical cases, not the numbers of usage only. It's highly depending on the bone density and usage practice.

● Replacement of drills should be considered usually after 50 times of usage, or when the cutting efficiency declines. Check the conditions of drill in every time of usage.

The reusable drills must be sterilized before use and require specific cleaning procedures prior to sterilization.

2. Intended Use

These products are surgical instrument kit designed for the surgical placement of dental implants or prosthesis.

3. Pre-operative planning

For a successful treatment, the following confirmation procedures should be performed carefully to confirm the treatment plan and evaluation, and also the product with medical and dentistry's point of view.

(1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.

(2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.

(3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.

(4) The product is supplied as non-sterile devices. All product used for surgery must be sterilized before use.

(5)Choose the surgical instrument after considering right dimension at surgical site.

(6)Depth Gauge is needed to confirm that the product is the same as shown in the catalog. In particular, it is recommended to use after checking if laser marking is in the same location.

* Marking can easily be erased when the product is left for a long time after the procedure, so caution is required.

(7) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, expired date, and etc. after opening the package.

(8) Check the product information such as reference code, size, manufacturing date and expiration date before starting the surgery.

(9) Check for compatibility with the opponent component before surgery.

(10) Check if there is any problem in the shape of the instrument.

4. Instruction for Use & Procedure

For a detailed procedure and the features of our dental implant system components used with the products, please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.

5 Adverse Effects

(1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.

(2) Complications such as malaise, pain, edema, bleeding, Inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient overall health.

(3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.

6. Contraindications

Contraindications should be observed when selecting patients. These include, but are not limited to:

(1) Low bone quality and inadequate bone quantity

(2) Bad oral hygiene and pathological infection at the alveolar bone.

(3) Internal diseases such as hematodyscrasia and diabetes, or mal-nutrition

(4) Alcoholic, drug abuse, heavy smoking and a mental patient.

(5) Other periodontal and surgical operation prohibited patient

7. Preservation & Handling

(1) The surgeon should check each product being used during surgery is in right position and whether if the product is contaminated or not.

(2) The surgical instrument should be used only for dental implant surgery following procedure as guided.

(3) The edge of surgical instrument for making the bone cavity should be checked before use.

(4) The surgical instrument for making the bone cavity should be used with saline solution to prevent necrosis.

(5) The surgeon should use the surgical instrument with suitable RPM(revolutions per minute) by considering bone condition such as bone density.

(6) It should be used surgical instruments of relevant Implant System.

(7) A small devices such as gauges with measuring function can be entered through throat during surgery. Fixation with thread or wire should be needed.

(8) The surgical instruments are reusable and dental drills should be replaced after being used 50 times or when the cutting efficiency declines. The lifetime of dental drill can be changed by handling method.

(9)The frequency of use could be changed depending on the handling.

(10)It can not be used when the cutting force of Drill is reduced.

(11) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.

(12) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.

(13) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

(14) If the package that is not opened by dentist is defective, request a replacement from manufacture.

(15) This Product should be stored at 15℃~30℃ condition.

8. Reusable method and Storage

The used surgical instruments should be store after cleaning and sterilization.

8.1. Cleaning

❑Clean surgical instruments separately and go through a special cleaning procedure before sterilization as mentioned below.

1) Remove all remaining pollutants (tissues, blood, bone residuals, etc.) on the surface by immersing used surgical instruments in water. Keep surgical instruments in wet environment until the next step indicated below is initiated.

2-1) Wash surgical instruments manually (conventional method).

a. Immerse surgical instruments in the water mixed with enzyme detergent (e.g. enzyme detergent with a pH level between 6 and 9) or its similar detergents (available to use commercial products).

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CE 2460

www.imegagen.com

Cleaning and Sterilization of Surgical and Prosthetic Instrument Kits

MEGAGEN

b. Scrub the surface of surgical instruments until all visible pollutants will be removed with a suitable soft bristled nylon brush (more than 4minutes).
Clean minute parts of surgical instruments with cleansing needle.
c. Rinse the surface of surgical instruments with tap water to remove all
2-2) Automatic washing using ultrasonic cleaner (recommendation method)
a. After putting surgical instruments in a washing container (e.g. plastic basket) inside of an ultrasonic cleaner, fill water completely. [You'll get better effects when using enzyme detergent (e.g. enzyme detergent with a pH level between 6 and 9) or its similar detergents (available to use commercial products.)]
b. Clean them for 5~15 minutes and drain the water.
c. Rinse the surface of the surgical instruments with tap water to remove all cleaning solution.
❏Repeat step a, b, c until surgical instruments are completely clean and perform a visual inspection.
❏Dry the surface of surgical instruments by natural drying or according to a drying cycle of cleaner. If needed, manual drying using nonabrasive tower with no lint is possible.
❏For more information about available enzyme detergents and cleaner, refer to instructions from the manufacturer and have to be followed.
❏Wash surgical devices separately and go through a special washing procedure as below before sterilization.

8.2. Sterilization

❏Surgical Instruments should be sterilized at temperature of 132℃ (270℉) for 15 minutes by high pressure steam sterilization and dry for 30 minutes.
❏In case of more than 149℃(300.2℉), a kit box may be deformed so it should be careful.

9. Caution ⚠

- Product should not be used if packaging is opened or damaged.
- Do not use the contaminated product by mistake of surgeon during the surgery.
- The product should not be reused products between patients without sterilization.
- If the instrument is contaminated, even a single patient should not use it in another part.
- Immediate action should be taken in the event of an accident.
- Quality and accuracy cannot be guaranteed if the customer (user) uses similar products not manufactured or sold by us.
- If the direction of fixation between the fixture and the connector is not correct, it may be deformed.
- If the recommended number of uses is exceeded, product deformation or fractures may be indicated.
- User shall state that the product should not be used if the packaging is not completely sealed or torn.
- Patients aged 18 and over who have completed the growth of the jaw bone (Partial or complete edentate patient)
- Manufacturer is not responsible for any damages to the product from sterilization or improper handling other than the high pressure steam sterilization method suggested by.

* Problems can occur if not used in the recommended way by manufacture.
* Guide should be used when using Removal Tool, if don't use guide, tool could be fracture.
* Changing the instrument's path during use may cause deformation of the product.

10. Warning

The following is strongly to be recommended for safe and effective use of MegaGen Implant System.
(1) The surgeon should be a medical specialist who practiced high technique for implant surgery.
(2) The surgeon should be fully aware of the operating method and cautions for the product before the use.
(3) The surgeon should wear an aseptic operating gown, and perform surgery in the operating room remained the aseptic condition during surgery.
(4) Inappropriate selection of patient and technique can cause a failure of the implant surgery and loss of supportive bone.
(5) Patients poor oral hygiene management can lead to side effects.
(6) Careless animal rearing in the dentist can contaminate the instrument.
(7) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
(8) The product must be used in clean and pleasant environment area.

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	Manufacturer		Caution, consult accompanying documents
CE 2460	CE conformity marking		Consult instructions for use
REF	Catalogue number		Temperature limitation
SN	Serial number		Do not use if package is damaged
LOT	Batch code		Date of manufacture
NON STERILE	Non-sterile		Use by date
Rx Only	For USA Only (Federal law(USA) restricts this device to sale by or on the order of a licensed dentist)		
EC REP	Authorized representative in the european community		

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MEGAGEN IMPLANT CO., LTD