



Boston Scientific
Advancing science for life™


Boston Scientific Neutral Electrodes

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


Instructions for Use Neutral Electrodes

This document applies to Boston Scientific Neutral Electrodes REF DGP-PMC2-5, REF DGP-PMC2-25, REF DGP-PM2-5, and REF DGP-PM2-25.

Intended Use

Neutral electrodes are an accessory for monopolar HF-surgery and represent the large area and low impedance contact with the patient's skin required for returning the electric current to the HF-generator. Neutral electrodes are self-adhesive, ready-to-use disposable products.

The Boston Scientific Neutral Electrodes REF DGP-PMC2-5, REF DGP-PMC2-25, REF DGP-PM2-5, and REF DGP-PM2-25 are intended to be used in combination with the devices "Boston Scientific G4™ and 1A/1B Radiofrequency Generators".

Indications

Indication is any medical condition that is treated by monopolar biterminal electrosurgery.

Patient groups are children (5 kg to 15 kg) and adults (>15 kg). Electrodes are specified for one or more of these patient groups.

Contraindications

Neutral electrodes must not be applied to injured or damaged skin. Therefore, a general contraindication is applicable to patients with large-area skin injuries (e.g. burns) if there is not enough intact skin application area available.

General Warnings Safety and Precautions

Neutral electrodes may be used by medical professionals only who must have undergone training regarding monopolar HF-surgery. Ignoring these instrucions for use may cause severe damage to the patient's skin!

- A safe use of the neutral electrodes REF DGP-PMC2-5, REF DGP-PMC2-25, REF DGP-PM2-5, and REF DGP-PM2-25 has been tested against a maximum allowed continuous RF-current activation on the "Boston Scientific G4™ and 1A/1B Radiofrequency Generators" of 3 minutes. A longer continuous activation of the current is not allowed, because this could lead to undesired severe tissue damage under the neutral electrode.
- The neutral electrodes are ready-for-use. Do not use additional contact gel.
- Prior to application, make sure that the neutral electrode is intact and that the expiration date has not been exceeded.
- Any modification of the neutral electrode (e.g. cutting) voids any liability of the manufacturer.
- Neutral electrodes are single use products. Do not reapply an improperly applied electrode.
- Review the application scopes for the different types of electrodes according to the weight and power specifications of the HF-generator as indicated on the packaging (max. 250W for children, max. 150W for newborns).
- With non-split neutral electrodes, regularly check that the entire surface area of the neutral electrode has contact with the skin (especially after moving the patient) because this is not monitored by the HF-generator.
- Use neutral electrodes that are compatible with your HF-generator only. Therefore, review the Declaration of Compatibility.
- Always observe the instructions for use of the HF-generator.

1. **Positioning** - Select an application site for the neutral electrode as close as possible to the operating site, but not closer than 20 cm. Avoid proximity to ECG electrodes and ensure that the heart is never in the current path. Preferred application sites are the thigh and upper arm. Always ensure the patient is not lying on top of the neutral electrode.
2. **Avoid** - The neutral electrode must not be put on scars, inflamed or injured skin, bony parts of the body, above metal implants or above massive subcutaneous fat tissue.
3. **Prep** - Carefully remove the patient's hair at the application site. If necessary, clean it. The application site must be dry and free of skin creams. Allow any used cleaning products to dry out completely.
4. **Open** - Open the pouch and remove the neutral electrode. Once the pouch has been opened, use within 7 days. Always keep unused neutral electrodes in the pouch and close it by folding over.
5. **Peel** - Check that the electrode and connected cable (if applicable) are free from any damages. Then, peel off the neutral electrode from the protective foil. Do not use damaged electrodes.
6. **Caution** - Rectangular neutral electrodes must be applied in such a way that the long edge faces the site of surgery. Electrodes with an oval shape can be oriented in any direction due to their design.
7. **Apply** - Attach the neutral electrode with its full area and without creases. Make sure that there are no air bubbles below the neutral electrode and that the protective edge (foam edge) is applied to the patient completely and without creases.
8. **Connect-** Link the neutral electrode with the generator according to its instructions for use. Avoid contact of the cable with the patient and other cables.
9. **Remove** - Remove the neutral electrode carefully from one corner after use to avoid skin irritations. Please note that disposal must be in line with the applicable laws and requirements of the respective country.

Explanation of Packaging Symbols

A. Electrode suited for adults (weight > 15 kg)	**I.** Temperature limitation	**Q.** Does not contain PVC
B. Electrode suited for children (weight 5 kg to 15 kg)	**J.** Keep dry	**R.** Do not use if pacakge is damaged.
C. Electrode not suited for newborns (weight < 5 kg)	**K.** Keep Away from sunlight	**S.** Quantity
D. Electrode with preattached cable and 4 mm connector	**L.** Green Dot scheme	**T.** Lot Number
E. For professional users only	**M.** CE Marking of Conformity and number of the involved Notified Body	**U.** Catalog Number
F. Consult instructions for use	**N.** Medical Device	**V.** Manufacturer
G. Single use. Do not re-use	**O.** Date of Manufacture	**W.** Authorized representative in the European Community
H. Not made with natural rubber latex	**P.** Use by date	

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All rights reserved.
Nissha Medical Technologies Ltd
Nissha Medical Technologies SAS
Boston Scientific (Australia) Pty Ltd
©2020 Boston Scientific Corporation or its affiliates.