TERUMO

CE 2797

2022-11 Rev.5

Runthrough™ NS

Extension Wire

PTCA Guide Wire
Guide PTCA
PTCA Führungsdraht
Guía ACTP
Fio-Guia para PTCA
Guida per PTCA
PTCA Voerdraad

PTCA Guide Wire
PTCA Guidewire
PTCA Guidewire
PTCA-ohjainlanka
Οδηγό Σύρμα PTCA
Проводник ЧТКА
Prowadnik do PTCA


PTCA vezetődrót
Vodící drát PTCA
Vodiaci drôt pre PTCA
PTCA Kılavuz Tel
PTCA Juhtetraat
PTCA Vadītājstīga
PTCA Kreipiamoji viela

Vodilna žica PTCA
PTCA žica vodič
Ghid pentru PTCA
Водач за PTCA
Провідник ЧТКА
Žica vodilica PTCA

Fig. / Fig. / Abb. / Fig. / Fig. / Fig / Fig. / Bild / Fig. / Fig. / Kuva. / Εικ. / Рис. / Ryc / Ábra / Obr. / Obr . / Şekil. / Joonis / Att. / Pav. / Sl. / Sl. / Fig. / Фигура / Мал. / Sl.

1

Grasp as closely to the joint as possible.
Saisir au plus près du joint.
Fassen Sie so dicht wie möglich an der Verbindungsstelle an
Sosténgalas lo más cerca posible de la articulación
Segurar o mais perto possível da junta
Serrare il più vicino possibile al giunto.
Zo dicht mogelijk bij de verbinding vasthouden.
Ta tag så nära fogen som möjligt.
Tag fat så tæt på leddet som muligt.
Grip så nær koplingen som mulig.
Tartu kiinni mahdollisimman läheltä liitosta.
Πιάστε όσο γίνεται πιο κοντά στην ένωση
Захватите как можно ближе к соединению.
Należy uchwycić jak najbliżej złącza.
Olyan közel fogja meg a csatlakozáshoz, amennyire csak lehetséges.
Uchopte dle možnosti co nejblíže spoje.
Uchopte čo najbližšie k miestu spojenia.
Mümkün olduğunca ek yerine yakın tutun.
Haarake jätkukohast võimalikult lähedalt.
Satveriet pēc iespējas tuvāk pie savienojuma.
Suimkite kuo arčiau sandūros.
Primite čim bliže spoju.
Uhvatite što je moguće bliže spoju.
Apucați cât mai aproape de îmbinare posibil.
Захванете колкото е възможно по-близо до съединението.
Захоплюйте якомога ближче до з'єднання.
Uhvatite što je moguće bliže spoju.

TW34M0600-04

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ENGLISH

INSTRUCTIONS FOR USE

Read the following warnings, precautions, and directions for use carefully.

DESCRIPTION

The RUNTHROUGH NS Extension Wire is a PTFE-coated stainless steel guide wire attachment with a diameter of 0.014"(0.36mm) and length of 150cm. The shaft coating provides lubricity to reduce frictional resistance with interventional catheters. It is exclusively compatible with RUNTHROUGH NS Series ("guide wire") which have a modified proximal end that allow for connection to the RUNTHROUGH NS Extension Wire.

INTENDED PURPOSE

The RUNTHROUGH NS Extension Wire is a wire intended to extend the RUNTHROUGH NS for facilitating the replacement of interventional devices. Do not use for any other purpose.

INDICATION

RUNTHROUGH NS Extension Wire is indicated for coronary artery disease (CAD), characterized by narrowing of the coronary arteries that carry blood to the heart.

CONTRAINDICATIONS

Carefully read the list of CONTRAINDICATIONS in the instructions for RUNTHROUGH NS which is used with the RUNTHROUGH NS Extension Wire.

PATIENT TARGET GROUP

Patients who need coronary angioplasty.

INTENDED USERS

The RUNTHROUGH NS Extension Wire shall be used by interventional cardiologists who are well trained in angiography and percutaneous transluminal coronary angioplasty (PTCA).

INDIRECT CLINICAL BENEFIT

It enables coronary intervention by extending a RUNTHROUGH NS and exchanging interventional devices while the insertion route is maintained.

COMPLICATIONS

Possible complications associated with the use of the RUNTHROUGH NS Extension Wire include, but are not limited to, the following:

- Arterial perforation
- Arterial rupture
- Coronary artery injury
- Infectious disease and puncture site complication

Carefully read the COMPLICATIONS in the instructions for RUNTHROUGH NS which is used with the RUNTHROUGH NS Extension Wire.

WARNINGS AND PRECAUTIONS

WARNINGS

Prior to the procedure, carefully examine all equipment to verify its proper function and integrity. Failure to abide by the following warnings might result in damage to the vessel, release of plastic fragments from the RUNTHROUGH NS Extension Wire and/or damage to or breakage/separation of the wire, that may necessitate retrieval.

- **Do not manipulate and/or withdraw the RUNTHROUGH NS Extension Wire through a metal entry needle, a metal dilator, or a metal guide wire inserter. Manipulation and/or withdrawal through a metal entry needle, a metal dilator, or a metal guide wire inserter may result in destruction and/or separation of the RUNTHROUGH NS Extension Wire.**
- **Do not use the RUNTHROUGH NS Extension Wire with devices which contain metal parts such as Atherectomy catheters.**
- **Attachment and detachment to and from the guide wire should be done while confirming the position of guide wire tip and interventional device under high resolution fluoroscopy. The unintentional advancement of the wire into the terminus of a small lateral branch may result in perforation or damage to the inner vessel membrane.**
- **If any resistance is felt or if the tip's behaviour and/or location seems improper, stop manipulating the RUNTHROUGH NS Extension Wire and/or the interventional device and determine the cause by high resolution fluoroscopy. Failure to exercise proper caution may result in bending, kinking, separation of the guide wire's tip, damage to the interventional device, or damage to the vessel.**

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- **Do not bend the RUNTHROUGH NS Extension Wire repeatedly at the same point. It may result in deformation, breakage, or separation of the RUNTHROUGH NS Extension Wire.**
- **If an interventional device is used with a stopcock, do not operate this stopcock while the RUNTHROUGH NS Extension Wire is in the device. It may result in breakage or separation of the RUNTHROUGH NS Extension Wire.**
- **Do not use under MRI.**
- **After replacement of interventional device using the RUNTHROUGH NS Extension Wire, make sure that the extension wire is removed before next intervention. Separation of the joint may result from manipulation such as rotation while the RUNTHROUGH NS and the RUNTHROUGH NS Extension Wire are connected.**
- **The RUNTHROUGH NS Extension Wire is exclusively compatible with the RUNTHROUGH NS which have a modified proximal end that allow for connection to the RUNTHROUGH NS Extension Wire. Check the product label on the package of the RUNTHROUGH NS product to confirm that the guide wire is extendable.**
- **Do not torque or manipulate the RUNTHROUGH NS Extension Wire once it is attached. This may cause the RUNTHROUGH NS Extension Wire to detach from the guide wire.**
- **Do not insert the RUNTHROUGH NS Extension Wire in the living body. It may result in disconnection of the joint, damage to the vessel or breakage of the RUNTHROUGH NS Extension Wire.**
- **Cobalt (CAS:7440-48-4) classified as CMR* 1B is present in the shaft in a concentration above 0.1% weight by weight.**
 ***CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008).**
 Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

PRECAUTIONS

- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Carefully inspect the unit package to make sure that no damage to the sterile barrier has occurred.
- Sterile in an unopened and undamaged unit package. Do not use if the unit package or the RUNTHROUGH NS Extension Wire is damaged or soiled.
- The product should be used immediately after opening the package, and be disposed of safely and properly after use following the local regulations. The product is biohazardous as it is contaminated with blood and is physically hazardous due to its sharp edge.
- The entire procedure should be carried out aseptically.
- When using a drug or a device concurrently with the RUNTHROUGH NS Extension Wire, the operator should have a full understanding of the properties/characteristics of the drug or device so as to avoid damage to the RUNTHROUGH NS Extension Wire.
- Do not wipe the RUNTHROUGH NS Extension Wire with alcohol, or gluconic acid chlorhexidine solution because they may adversely affect the surface of the RUNTHROUGH NS Extension Wire.

DIRECTIONS FOR USE

1.Preparation of the RUNTHROUGH NS Extension Wire

1-1 Remove the RUNTHROUGH NS Extension Wire and the holder from the package.

1-2 After removing the wire stopper, remove the RUNTHROUGH NS Extension Wire from the holder.

PRECAUTIONS

- Make sure to detach the wire stopper before removing the RUNTHROUGH NS Extension Wire from the holder to avoid breakage of the product.
- The proximal end of the RUNTHROUGH NS Extension Wire exposes from the holder. Make sure of the direction of the RUNTHROUGH NS Extension Wire when jointing the guide wire.
- When reinserting the RUNTHROUGH NS Extension Wire back into the holder, take care not to damage the surface of the RUNTHROUGH NS Extension Wire with the edge of the holder.

1-3 Connect the proximal end of the guide wire with the distal end of the RUNTHROUGH NS Extension Wire for connection.

PRECAUTION

When connecting the guide wire and the RUNTHROUGH NS Extension Wire, grasp as closely to the joint as possible. Failure to do so can result in breakage of the guide wire or the RUNTHROUGH NS Extension Wire (Fig1.).

1-4 Check for secure attachment by pulling gently both the guide wire and the RUNTHROUGH NS Extension Wire on each side of the joint.

1-5 Withdraw the inserted interventional device over the extended guide wire.

1-6 Reconfirm the connection by pulling gently both side of the joint.

1-7 Before advancement of another interventional device, confirm the location of distal end of the guide wire under fluoroscopy.

1-8 Insert another interventional device from the proximal end of the RUNTHROUGH NS Extension Wire.

1-9 Advance the interventional device while maintaining position of the guide wire and guiding catheter.

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1-10 After replacement of interventional device, grip the joint of the guide wire and the RUNTHROUGH NS Extension Wire, and release the connection at the piping while maintaining the position of the guide wire.

PRECAUTIONS

- Wipe the blood residues on the proximal end of the guide wire before connecting the RUNTHROUGH NS Extension Wire and the guide wire. The thrombus may prevent smooth detachment of the joint.
- Do not advance the joint of the RUNTHROUGH NS Extension Wire and the guide wire into a Y connector. It may cause detachment of the joint.
- Do not insert the proximal end of the RUNTHROUGH NS Extension Wire into the interventional device when the device is curved. Interventional device may be damaged at by the proximal end of the RUNTHROUGH NS Extension Wire.
- If any resistance is felt at insertion of the RUNTHROUGH NS Extension Wire, stop manipulating the RUNTHROUGH NS Extension Wire and determine the cause. Successive manipulation under resistance may result in damage to vessel by the RUNTHROUGH NS Extension Wire and the interventional device.
- When wiping blood and other residues while the guide wire and the RUNTHROUGH NS Extension Wire connected, grip the joint firmly so that any load would not be applied to the joint. It may cause detachment of the connection of the guide wire and the RUNTHROUGH NS Extension Wire.
- When connection is degraded from repetition of attachment and detachment of the RUNTHROUGH NS Extension Wire and the guide wire, stop manipulation and exchange it for a new one. It may result in separation at the joint and cause damage to the RUNTHROUGH NS Extension Wire.
- Use an appropriate interventional device corresponding to the outer diameter (O.D.) of the RUNTHROUGH NS Extension Wire. Due to the variation of certain interventional device tip inner diameters, detachment of the connection of the guide wire and the RUNTHROUGH NS Extension Wire may occur during manipulation. If any resistance is felt during introduction of the RUNTHROUGH NS Extension Wire into the interventional device, it is advisable to stop using such devices.

2.For temporary storage during the procedure

2-1 Remove the blood residues on the RUNTHROUGH NS Extension Wire by rinsing in a bowl of heparinized physiological saline solution. Wipe once with sterile gauze moistened with heparinized physiological saline solution if any blood residues still remains on the RUNTHROUGH NS Extension Wire.

2-2 Reinsert the RUNTHROUGH NS Extension Wire into the holder.

PRECAUTIONS FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperature, or high humidity during storage.

REPORT OF INCIDENT

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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© TERUMO CORPORATION 2022-11
All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.