Filename Date Time
23.0 x 10.6 inches (584 mm x 269 mm)
Folded: 2.6 x 5.3 inches (66 mm x 135 mm)

Fold

Medtronic Confidential
IFUGlobalManualLarge.xsl - IFUGlobalManualLarge.fm
Template version 6.5: 09-30-2020

Medtronic

SenSight™ **B31010**
Depth Stop and Cranial Tunneler Kit

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© Medtronic 2020
All Rights Reserved.
M991029A002 Rev A

M991029A002

2021-03-01

M991029A002 Rev A

Printing instructions: *Refer to doc# A88898* for Printing Instructions

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Description

The Medtronic Model B31010 SenSight™ depth stop and cranial tunneler kit are components used during a lead implant procedure for a deep brain stimulation (DBS) system.

Intended purpose

Depth stop: The depth stop is an external device intended to retain the DBS lead position during implant.

Cranial Tunneler: The cranial tunneler is a device intended to create a subcutaneous path for tunneling the DBS lead to the intended location during implant.

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Leads: The SenSight cranial tunneler can be used with all Medtronic DBS leads.

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Package contents

The product kit contains the following sterile components:

- 1 Depth stop
- 1 Cranial tunneler
 - Sharp tip
 - Passing straw
 - Tunneling rod

Specifications


Figure 1. Depth stop and cranial tunneler.

① Sharp tip (cranial tunneler part)
② Passing tube (cranial tunneler part)
③ Tunneling rod (cranial tunneler part)
④ Depth stop

All measurements are approximate.

Depth stop

Inner diameter = 1.52 mm

The material and substance to which the patients can be exposed is polyoxymethylene.

Cranial tunneler

Passing straw diameter = 3.40 mm

The materials and substances to which the patients can be exposed are stainless steel and polytetrafluoroethylene (PTFE). The materials listed are ordered from the most to least amount.

Discuss allergies or other intolerances related to the materials with the patient before the procedure. Stainless steel contains nickel. Nickel is used in mixtures of metals, and is a known allergen.

Note: These materials were tested for category 1A or 1B carcinogenic, mutagenic, or toxic for reproduction (CMR) substances, and endocrine disrupting chemicals (EDC). No known CMRs or EDCs were found in the materials used for this component.

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Refer to the appropriate lead implant manual for procedural information related to the use of these components.

Refer to the information for prescribers manual for information on component disposal.

Using the depth stop

1. Attach the depth stop to the lead, ensuring that the lead is aligned to the proper target depth for insertion.
2. Tighten the screw on the depth stop to secure it.
3. To remove the depth stop, loosen the screw and remove the depth stop from the lead.

 ⚠ **Caution:** Ensure that the depth stop is fully loosened before removing it from the lead. If the depth stop is not loosened completely before it is removed, lead damage and inadequate stimulation may occur, requiring revision surgery.

Using the cranial tunneler

1. Remove the sharp tip from the plastic protector.
2. Attach the sharp tip to the tunneling rod.
3. Create an incision above and behind the ear.
4. Tunnel between the cranial incision and incision created above and behind the ear.
5. When tunneling is complete, carefully remove the sharp tip and tunneling rod, leaving the passing tube in place.
6. Pass the lead end to the desired location, using the passing tube provided.
7. Carefully remove the passing tube, leaving the lead in place.

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English

Spare part kits

The components contained in these Medtronic kits include spare parts. For information and instructions on the use of these parts and tools, which include indications, contraindications, warnings, precautions, and adverse events, refer to the appropriate product manual.

Important: When replacing a system component with a spare part from this kit, make sure the part in the kit matches the original component.

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Medtronic

SenSight™ Burr Hole Device Kit B32000
SenSight™ Burr Hole Device B31000

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Medtronic and Medtronic logo are trademarks of Medtronic. Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

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Table of contents

Description 7
Intended purpose 7
Compatibility 7
Before the procedure 7
Package contents 7
Device Specifications 9
Instructions for use 9
Preparing the burr hole 9
Securing the base of the burr hole device 9
Inserting the support clip into the base 11
Closing the support clip and anchoring the lead 12
Opening the support clip 12
Snapping the cover into the base 13
Removing the cover from the base 13

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Information available for the Deep Brain Stimulation Therapies:

The information for prescribers manual provides information about contraindications, warnings, precautions, adverse events, sterilization, and component disposal. The information for prescribers addendums contain indication-specific information. For customers in Japan, the appropriate package insert provides information about safety, contraindications, warnings, precautions, and adverse events.

The indications sheet provides information about indications and related information. For customers in Japan, the appropriate package insert provides information about indications.

The system eligibility and battery longevity manual describes programming considerations and provides battery longevity information to aid in the appropriate neurostimulator selection.

MRI guidelines provide information about any MRI conditions and MRI-specific contraindications, warnings, and precautions for MRI scans with the neurostimulation system.

Product manuals, such as programming guides, recharging guides, and implant manuals provide device descriptions, package contents, device specifications, product-specific warnings and precautions, and instructions for use.

! USA **The clinical summary provides information about the clinical study results for the neurostimulation system.**

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Description

The Medtronic Models B32000 and B31000 SenSight™ burr hole devices are used with a neurostimulation system for deep brain stimulation (DBS).

Note: The burr hole device alone is MR Conditional. Refer to the MRI Guidelines.

Intended purpose

The burr hole device is an implanted device intended to cover a 14-mm burr hole on the skull, and anchor a single compatible lead following cranial surgery.

Compatibility

⚠ **Warning:** Only use SenSight components with other SenSight components, unless specific compatibility is stated in this manual. Do not use incompatible components. Using incompatible components may result in component damage, intermittent or loss of stimulation, and require surgical replacement or revision.

Before the procedure

⚠ **Warning:** Do not reuse, reprocess, or resterilize single use products. Reusing, reprocessing, or resterilizing may compromise the functional integrity of the products and may create a risk of contamination, which could result in patient injury, illness, or death.

⚠ **Warning:** Medtronic has sterilized the contents of the package according to the process indicated on the package label. Do not use a component if the following circumstances have occurred:

- The storage package or sterile seal has been pierced or altered and component sterility cannot be guaranteed and infection may occur.
- The component shows signs of damage and the component may not function properly.
- The use-by date has expired and component sterility cannot be guaranteed and infection may occur.

Package contents

- 1 Burr hole device (Model B32000 or B31000)
 - 1 Base
 - 1 Centering tool (preattached to the base)
 - 2 Bone screws (preinserted in the legs of the centering tool)
 - 1 Support clip
 - 1 Cover
- 3 Tools (Model B32000 kits only)
 - 1 Insertion tool (gray and blue; use with Nexframe™ Stereotactic Systems)
 - 1 Insertion tool (clear; use with stereotactic frames)
 - 1 Disposable screwdriver
- Product literature
- ! USA Warranty card

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
Figure 1. *Model B32000 kit components.*

① Centering tool (preattached)
② Bone screw (preinserted)
③ Base
④ Support clip
⑤ Screwdriver
⑥ Insertion tool (gray and blue)
⑦ Insertion tool (clear)
⑧ Cover

Note: The Model B31000 burr hole device package contains the complete burr hole device assembly without the insertion tools or screwdriver.

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Device Specifications

Table 1. Model B32000 and B31000 burr hole devices[a]

Description	Value
Expected lifetime	5 years
Surface area	24.6 cm²
Materials and substances to which the patient can be exposed[b,c,d]	Polyamide, titanium, titanium dioxide

[a] All measurements are approximate.

[b] Discuss any allergies or other intolerances related to the materials and substances with the patient before the procedure.

[c] Tested for category 1A or 1B carcinogenic, mutagenic, or toxic for reproduction (CMR) substances, and endocrine disrupting chemicals (EDC). No known CMRs or EDCs found for the materials and substances listed in this table.

[d] The materials listed are ordered from the most to least amount.

Instructions for use

Inspect the burr hole device for any damage before the procedure.

Preparing the burr hole

1. Identify the burr hole location and prepare the target site.

 ⚠ **Warning:** Do not attach the burr hole device to a cranial suture, to cranial bone that is unstable, damaged, diseased, or less than 4 mm thick. If the burr hole device is attached to this type of bone, the following can result:
 - The lead will not be properly anchored, causing potential lead migration, and intermittent or loss of stimulation.
 - The burr hole closure will be unstable, making it difficult to implant the burr hole device. Unstable burr hole closure increases the risk of infection and may cause potential damage to the dura, possible cerebral spinal fluid (CSF) leakage, and damage to the neural (brain) tissue.

2. Make a scalp incision.
3. Create a 14-mm diameter burr hole.
4. Ensure that the burr hole is clear of debris, for example, no bone chips and edges are smooth.

Securing the base of the burr hole device

1. Position the centering tool (preattached to the base) directly over the burr hole (Figure 2).


***Figure 2.** Placing the centering tool and base over the burr hole.*

2. Hold the base in position against the skull, placing a finger on both sides of the centering tool.
3. Attach the base with the two preinserted bone screws using the screwdriver provided (Figure 3).

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
Figure 3. *Tightening the screws.*

a. Tighten the screws until they are past the legs of the centering tool.

 Caution: Do not overtighten the screws. Overtightening the screws may damage the base of the burr hole device, requiring replacement.

b. Confirm that the screws are fully seated, and that the base is secured before removing the centering tool. The screws should be flush with the base to fully seat them.

 Warning: Do not tighten the screws over an exposed burr hole. Keep the burr hole covered until the screws are completely tightened. Tightening the screws over an exposed burr hole could cause additional trauma to the patient.

4. Lift off the centering tool to remove it and set aside for disposal (Figure 4).


Figure 4. *Removing the centering tool.*

Note: If the lead implant procedure could not be completed, or the burr hole device will not be used to anchor a lead, proceed to "Snapping the cover into the base" on page 13.

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Inserting the support clip into the base


Figure 5. *Support clip features (top view).*

1. Hole
2. Tab 2
3. Lead holder post
4. Static side of the securing mechanism
5. Support clip orientation feature
6. Tab 1
7. Closing notch of the securing mechanism
8. Triangle orientation feature

1. Ensure that the burr hole is clear of all obstructions, including bone. Obstructions will prevent proper insertion of the support clip.
2. Fully insert the metal tip of the insertion tool into the hole located near tab 2 on the support clip (Figure 5).
3. Use the insertion tool to angle the support clip for alignment with the base, and position the static side of the securing mechanism against the lead (Figure 5 and Figure 6).

 Caution: Minimize contact between the lead body and the support clip to prevent changing the position of the implanted lead electrodes. Lead movement may change the position of the lead electrodes, requiring surgical replacement or revision.

4. Insert tab 1 of the support clip at an angle, while aligning the support clip orientation features between the base orientation features (Figure 6).

 Note: If the support clip is placed parallel to the base instead of at an angle, the securing mechanism may not engage properly.
5. Press tab 2 of the support clip into the base, and confirm the support clip is fully connected to the base (Figure 6).


Figure 6. *Inserting the support clip into the base.*

1. Tab 1
2. Support clip orientation feature
3. Lead body
4. Base orientation features
5. Tab 2

Note: A click may be heard or felt to confirm that the support clip is connected to the base. If the click is not heard or felt, visually confirm that all support clip orientation features are connected to the base. Confirm both tab 1 and tab 2 are fully connected.

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6. Use a back-and-forth motion, and gently remove the insertion tool from the support clip.

Closing the support clip and anchoring the lead

1. Keep the lead stylet in place.
2. Insert the plastic tip of the insertion tool into the closing notch on the support clip (Figure 7).


Figure 7. *Closing the securing mechanism and routing the lead.*

① Closing notch
② Static side of the securing mechanism
③ Lead holder posts
④ Locking tab
⑤ Routed lead body

3. Slide the closing notch on the support clip until the securing mechanism closes into the locked position against the lead.

 Note: The securing mechanism must be in the locked position to secure the lead. A click may be heard or felt to confirm that the securing mechanism is locked. If a click is not heard, visually confirm that the locking tab is visible.
4. Gently remove the lead stylet.
5. Select an exit slot for the lead, and route the lead between the lead holder posts near that exit slot (Figure 7).

 △ **Caution:** Minimize pushing and pulling on the lead when routing the lead between the lead holder posts. Lead movement may change the position of the lead electrodes, requiring surgical revision or replacement.

Opening the support clip

1. Locate the opening notch on the closed support clip (Figure 8).


Figure 8. *Opening the support clip.*

① Opening notch on the locking tab
② Routed lead body

Note: The opening notch is only visible when the securing mechanism on the support clip is locked.

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2. Use the plastic tip of the insertion tool, and apply downward pressure on the opening notch. The securing mechanism will disengage and unlock from the locked position.
3. Open the support clip.
4. Reposition the lead as appropriate and repeat the steps for "Closing the support clip and anchoring the lead" on page 12.

Snapping the cover into the base

1. Confirm that the securing mechanism is locked and the lead holder posts are unobstructed.

 Note: If the securing mechanism of the support clip is not locked, the cover may not snap into the base.
2. Align the triangle orientation feature on the cover over the identical triangle orientation feature on the support clip (Figure 9).
3. Align the holes on the cover over the lead holder posts (Figure 9).

 ⚠ **Warning:** Ensure that the lead is routed between the lead holder posts, and the cover is aligned with the base before snapping the cover into the base. An unsecured lead, or misaligned cover may damage the lead and require surgical replacement.


Figure 9. *Aligning the cover to snap into the base.*

① Triangle orientation features align
② Cover legs insert into the base
③ Cover holes align over the lead holder posts

 Note: If the orientation features on the cover are aligned over the support clip features, the cover will connect to the base.
4. Lower the cover and press down to snap all four legs of the cover onto the base (Figure 9).

 Note: Confirm that all four cover legs are inserted into the base.
5. Close and suture the skin flap over the burr hole device, taking care to avoid the lead body, if present.

Removing the cover from the base

⚠ **Warning:** Damage to the cover may occur during removal. If the cover is removed, use a new cover. Using the same cover after it is removed, may result in lead migration and require surgical revision or replacement.

1. Locate the exit slot on the cover that is opposite of the slot containing the routed lead.
2. Insert the plastic tip of the insertion tool into the open exit slot of the cover (Figure 10).

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
Figure 10. *Unsnapping the cover from the base.*

① Open exit slot
② Lead body
③ Lead exit slot

3. Apply upward pressure with the tip of the insertion tool to unsnap the cover from the base.
4. Remove the cover from the base and set it aside for disposal.

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Authorized Representative in the European Community EC REP
Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel. +31-45-566-8000

Europe/Africa/Middle East Headquarters
Medtronic International Trading Sàrl
Route du Molliau 31
Case Postale 84
CH - 1131 Tolochenaz
Switzerland
www.medtronic.eu
Tel. +41-21-802-7000

Asia-Pacific
Medtronic International Ltd.
50 Pasir Panjang Road
#04-51 Mapletree Business City
Singapore 117384
Singapore
Tel. +65-6870-5510

© Medtronic 2020
All Rights Reserved.
M000202C011 Rev A

M000202C011

2021-03-01

English

Trademarks

Medtronic and Medtronic logo are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

Description

The Medtronic Model B31061 SenSight™ connector plug is a component of a neurostimulation system.

Intended purpose

The connector plug fits into an unused socket of a dual channel neurostimulator, or the distal end of an unused DBS extension.

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Neurostimulators: The connector plug is compatible with the following Medtronic neurostimulators:

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Package contents

This package contains 2 connector plugs.

Specifications

Figure 1. Connector plug.

① Shaft
② Ring (for setscrew contact)
③ Handle

All measurements are approximate.

Expected lifetime = 11 years

Shaft diameter = 1.3 mm

Shaft length = 12.8 mm

Overall length = 17.1 mm

Surface area = 1.4 cm²

The material and substance to which the patients can be exposed is polyether urethane.

Discuss allergies or other intolerances related to the material with the patient before the procedure.

Note: This material was tested for category 1A or 1B carcinogenic, mutagenic, or toxic for reproduction (CMR) substances, and endocrine disrupting chemicals (EDC). No known CMRs or EDCs were found in the material used for this component.

Instructions for use

Before opening the package, verify the model number and use-by-date.

Using the connector plug with a neurostimulator

⚠ **Caution:** Always block an unused neurostimulator socket with the connector plug. An open socket may cause pocket stimulation and allow tissue ingress. A socket with ingressed tissue may be unsuitable for future use.

1. If using only 1 extension, insert a connector plug into the unused neurostimulator socket. Hold the connector plug handle and insert it fully into the socket, until the handle meets the neurostimulator header block.
2. Tighten the setscrew for the connector plug. Insert the torque wrench (packaged with the neurostimulator) fully into the self-sealing grommet. Rotate the torque wrench clockwise, until it clicks at least once.

For complete information about implanting the neurostimulator, refer to the appropriate neurostimulator implant manual.

Using the connector plug with an extension

⚠ **Caution:** Always protect an unused extension with the connector plug. If a connector plug was not used to protect the extension, intermittent stimulation may occur, or the extension may be unusable, which would require surgical replacement or revision.

1. Insert the connector plug into the unused extension. Hold the connector plug handle, and insert it fully until the handle meets the distal connector end of the extension.
2. Use the torque wrench (packaged with the extension) to tighten the setscrew for the connector plug.

⚠ **Caution:** Do not tighten the setscrew without stabilizing the setscrew connector block. Failure to stabilize the setscrew connector block can break or damage the extension, resulting in intermittent or loss of stimulation, requiring surgical replacement or revision.

⚠ **Caution:** Ensure that the torque wrench is perpendicular to the setscrew axis, and tighten the setscrew until the torque wrench clicks at least once. Failure to hold the wrench perpendicular while tightening can result in an insufficiently tightened setscrew, or too loose of a setscrew, and require surgical replacement or revision.

3. To remove the connector plug from the extension, expose the connector plug.
4. Retract the setscrew so that the connector plug can be removed.

⚠ **Caution:** Do not remove the connector plug from the extension until the setscrew is retracted. Removing the connector plug before the setscrew is retracted may cause damage to the extension, and require surgical replacement or revision.

5. Remove the connector plug and set it aside for disposal. Refer to the information for prescribers manual for information on component disposal.

For complete information about implanting the extension, refer to the appropriate extension implant manual.

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