ANGIPLAST PVT. LTD

PRODUCT INFORMATION

Product Name & Code	:	Angeltouch IV Cannula with Luer Lock, Wings and Injection Port [IT-31]
	:	Angeltouch IV Cannula with Luer Lock and without Wings and Injection Port [IT-32]
	:	Angeltouch IV Cannula with Luer Lock and Wings and without Injection Port [IT-33]
Product Description	:	The sterile IV cannula is used infusion or transfusion on continuous / intermittent basis and short-term peripheral catheterization of vein to give access for drug delivery. The sterile IV cannula is used for fluid introduced or withdrawn from human circulation system, comprising the channel with two female luer lock connections for infusion of two fluids at the same time.
Product Specification	:	 Bevel type needle for easier less traumatic insert. Thin Wall PTFE catheter with X- Ray opaque line. Flexible and kink resistant catheter is double tapered with rounded tip for smooth introduction and minimum discomfort to the patient. Injection port with one way valve and color coded cap for intermittent and safe administration of drugs during infusion. Flexible wings provide easy and secure fixation and also prevent skin contact with luer lock connector. High transparent flash back chamber for an easy view of returning blood. Blood stopper Luer lock cap for closing the proximal end of the catheter after removing the needle.

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ANGIPLAST PVT. LTD

PRODUCT INFORMATION

Material List & Make	:	**Sr. No.**	**Name of Components**	**No.**	**Material Used**
		1	Cannula	1	Stainless Steel (SS)
		2	Catheter	1	PTFE / FEP
		3	Wings	1	PTFE / FEP
		4	Protective Cap for Catheter	1	Poly-propylene (PP)
		5	Rotating Cap for Luer Lock	1	High Density Polyethylene (HDPE)
		6	Flash Back Chamber	1	PTFE / FEP
		7	Unit Pack Pouch	1	Blister Pack
Instruction For Use	:	**INSTRUCTIONS FOR USE :** - Select proper type and size of I.V. Cannula, always use the smallest size cannula possible to adequately deliver the desired I.V. Fluids. - Carefully select and clean the venipuncture site. - Inspect cannula package for integrity and expiry and then remove cannula from package. - Remove needle cover and grip the cannula from needle hub projection. - Insert cannula into the vein at a low angle and check for blood in flash back chamber to ensure proper venipuncture. - Advance the catheter into vein and simultaneously withdraw the needle gently. - Do not try to reinsert a partially or completely withdrawn needle. - Avoid blood spillage by pressing a finger on the vein above the catheter tip. - Replace needle with either the I.V. Infusion set line or the luer lock plug. - Discard the needle in an appropriate waste container. - Tape down the wings of the cannula to the patient skin. - Cover the puncture site with a sterile dressing. - For I.V. Cannula with injection port drugs can be injected intermittently with the help of a syringe through the injection port after removing the port cap. - For I.V. Cannula with injection port the injection port cap should remain closed when not in use. - Inspect venipuncture site regularly for any reactions and check all connections.			
Contraindications	:	**CONTRAINDICATIONS :** - Administration of highly viscous fluids. - Large blood transfusion. - Not to be used in patients with known hypersensitivity to any of the materials used. - Administration of irritant fluid or drugs.			

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ANGIPLAST PVT. LTD

PRODUCT INFORMATION

Indications	:	**INDICATIONS :** - Infusion of fluids such as solutions, Parenteral nutritions and administration of other drugs. - Maintaining hydration and / or correct dehydration in patients who are unable to take sufficient volume of oral fluids. - Transfusion of blood or blood components.
Warnings	:	**WARNINGS :** • Product is single use only. Re-use of the product can cause Infection & other blood related transfusion diseases like HIV, AIDS & Hepatitis etc. • Read instructions before use. • Check expiry date prior to use. • The product should be used according to the instructions for use. • Ensure before use that product / packing is not tampered or damaged. • Discard after use. Do not recycle, clean or resterile. • Ensure proper disposal of product and packing after use. • Use the product immediately after opening the individual blister packing. • This product should be used by a Doctor / Registered Practitioner or Paramedic. • Store in a cool and dry place, protect from excessive heat or direct sunlight.
Product Use	:	For Single Use Only
Is product Sterilized	:	Yes
Sterilization Method	:	Ethylene Oxide Gas Sterilization
Storage Conditions	:	Do not store in direct sunlight, at extreme temperature or in high humidity. Keep away from getting wet. Store in cool & dry place.

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ANGIPLAST PVT• LTD