Angiplast Pvt. Ltd.

PRODUCT INFORMATION

Product Name	:	Angeltouch Three Way Stop Cock [HIGH PRESSURE / PRESSURE RESISTANT]
Product Code	:	IT-30
Product Description	:	The sterile Three way stop cock is used for infusion medicinal fluid. It comprises the channel with two female luer lock connection for the infusion of two fluid at same time. Three way stop cock manufactured from transparent polycarbonate / PP / HDPE. Three way stop cock is E.T.O. sterile and Pyrogen free. It is also available in non-sterile. It can withstand pressure upto 5-6 bars. It can withstand high cracking resistance. Proper ditting of female luer lock with rotating male luer lock found. It proves the presence of built-in fitting. The filling volume must be not more than 0,30 ml (0,217 ml).
Product Specification	:	 Inside is completely smooth & available with molded lure lock. It is manufactured from non-toxic medical grade material.

Material List& Make :

Sr. No.	Name of Components	No.	Material Used
1	Three Way – Body	1	Poly-carbonate
2	Three Way – Handle	1	High Density Poly-ethylene
3	Protective Cap	1	PP
4	Protective Cap	2	PP
5	Unit Pack Pouch	1	Blister Pack

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Angiplast Pvt. Ltd.

PRODUCT INFORMATION

Instruction For Use	:	**Contraindications:** Not to be used for high viscous fluids and lipids administration. Product should not be used in patients with known hypersensitivity to any of the materials used. Not to be used for blood transfusion. **Instructions for use:** Check visually whether the package is intact. Open packaging and remove stopcock from sterile package. Connect one or two I.V. infusion line as desired to the ports. Remove protector from male portion & hold it upright. Open clamps of I.V. infusion line and displace air entirely from stopcock by rotating HS tap to open position. Rotate tap to close position. Connect it to venipuncture device & secure it with rotating lock. Rotate the tap to open position and adjust the flow of infusion set.
Product Use	:	For Single Use Only
Is product Sterilized	:	Yes
Sterilization Method	:	Ethylene Oxide Gas Sterilization
Storage Conditions	:	Do not store in direct sunlight, at extreme temperature or in high humidity. Keep away from getting wet. Store in cool & dry place.
Packing Detail	:	Inner Box Qty: 50 nos.
		Shipper Box Qty.: 1000 Nos.
		Shipper Box Size: 690 mm x 415 mm x 230 mm

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