Bausch

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TÜVRheinland

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization

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has established and applies a quality management system for medical devices for the following scope:

Design and development, production and distribution of medical devices for dental application for the purpose of articulation and occlusion testing
(see attachment for additional sites included)

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date:	2018-05-05
Certificate Registration No.:	SX 60127886 0001
An audit was performed. Report No.:	21190937 011
This Certificate is valid until:	2021-05-04

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ENGLISH

General description and purpose of the medical device:

Arti-Brux® is a red indicator paint, which is applied with a brush (single use product) to the occlusal surfaces of restorative and prosthetic appliances as well as hard splints. After drying, Arti-Brux® Quick forms a thin and even lacquer coating. The lacquer coating is intended for short-term use (up to max. 24 hours). The lacquer coating is saliva proof and is being rubbed off by occlusal contacts. Arti-Brux® Quick shows occlusal contacts and movement traces of the teeth by means of the removed lacquer layer.

For professional use only.

Indication

Arti-Brux® Quick BK 88 is a red liquid contact color:

- To control active dynamic occlusion on restorative and prosthetic restorations as well as on bite splints.
- For short-term control of occlusion (max. 24h) – Occlusogram.
- For visualization of enamel fissures as well as margin fissures of fillings and crowns.

Contraindication

- Arti-Brux® Quick consists mainly of concentrated alcohol (ethanol 96%). Arti-Brux® should not be applied directly invasively to patients who are addicted to alcohol, pregnant women or children. Also contact with injured mucosa should be prevented.
- Arti-Brux® Quick must not be used in case of iodine allergy.
- Arti-Brux® Quick must not be applied to open dentin or prepared teeth.

Side Effects

- With insufficiently polymerized acrylics or vacuum splints, made of different thermoplastics, discoloration of the plastic may occur.
- If the surfaces to be tested are not properly polished or sand blasted, the lacquer can penetrate into micro cracks of the plastic or ceramic surfaces and cause discoloration.
- When Arti-Brux® Quick is applied, the color can penetrate into enamel cracks or margins of fillings and cause temporary or even permanent discoloration.
- After longer use of Arti-Brux® Quick (up to max. 24 hours) discoloration of the saliva may also occur; the lacquer consists of a natural resin and red food coloring. After use, the mouth should be rinsed with water.

Notes

- Arti-Brux® Quick is a single use product and may only be used within one treatment session to prevent the transmission of germs.
- Disposal: the single use applicator and excessive paint must be disposed of with the usual contaminated waste.

Package content

- Single use applicator, content 0,15ml
- Instructions for use

Application

- The surfaces to be tested have to be cleaned and dried thoroughly.
- Apply the paint to the occlusal surfaces to be tested.
- To control bruxism, especially sleep bruxism, the patient should insert the splint shortly before going to bed.

Effect/Diagnosis

Check the surfaces to be tested after the end of the observation period. The coating is rubbed off by tooth contacts (occlusal surfaces). These contact patterns are used to check the occlusal contacts in static and dynamic. A uniform wear, e.g. on bite splints, is an indicator that the individual splint is evenly loaded. Strong one-sided wear indicates parafunctions (e.g. bruxism) or an incorrectly adjusted splint. The same applies to prosthetic restorations.

Cleaning (of the surfaces to be tested)

After successful inspection, remove the coating from the surface. Appropriate tools are: tooth brush and appropriate tooth paste, alcohol bath or wadding soaked with alcohol (Ethanol 96%). Prosthetic restorations can be cleaned with a toothbrush or standard practice equipment for the professional tooth cleaning.

Disposal Instructions

Arti-Brux® Quick consists of a red food color and a natural resin, which is also approved for food and pharmaceutical applications. All components of the kit can be disposed of with the usual practice waste.

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„Gebrauchsanweisung beachten"
„Follow manufacturers instructions"

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REF Bestellnummer
Reference number

LOT Chargen-Nummer
Batch code (Lot number)

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H225 Highly flammable liquid and vapor
P210 Keep away from heat, hot surfaces, sparks, open flames and other sources of ignition. No smoking.

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H319 Causes severe eye irritation.
P305 IF IN EYES:
P351 Rinse cautiously with water for several minutes.
P338 Remove contact lenses, if present and easy to do. Continue rinsing.

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Date of issue: 2019-12-06

1 General description and purpose of the medical device
Sterilizable, flexible blue silicone sleeves of 45 mm length with an outer diameter of 4.5 mm.
For professional use only.

Indication
Arti-Grip™ sleeves are used to improve the clamping force of Arti-Fol® forceps („Miller"-forceps). The silicone sleeves, pre-cut in 45 mm length, are pulled over one or both clamping surfaces of "Miller" forceps. The slightly adhesive surface of the silicone sleeves results in a stronger clamping effect. Thin occlusion foils and Shimstock foils are thus held securely. Dentists can use the Arti-Grip™ -sleeves, in the patient s mouth as well as on models. Arti-Grip™-sleeves can also be used in the laboratory.

Contraindication
No known contraindications.

Side effects
Possible side effects may be specific allergic reactions. However, given the product's history, there has not been a verified documented report of allergic reactions. Side effects or interactions may occur if the products are used with new, unknown products or materials.

2 Notes
- Before each use the Arti-Grip™-sleeves must be removed from the clamping surfaces of the forceps and subsequently must be cleaned, disinfected and sterilized.
- The Arti-Grip™-sleeves can be steam-sterilized (humid heat, 134°C).
- Arti-Grip™-sleeves can be used together with Arti-Fol®-forceps BK 132, BK 133, BK 142 and BK 144 and Arti-Fol® metallic, Shimstock-film.
- Disposal: before disposing Arti-Grip™-Silicone sleeves, they either have to be sterilized or disposed of with the normal contaminated practice waste.

3 Package content
- 20 pieces Silicone sleeves Arti-Grip™ (non-sterile)
- Instructions for Use

4 Preparation
Arti-Grip™-Silicone sleeves must be cleaned, disinfected and sterilized before each use. This also applies in particular to the initial application after delivery, since the instrument is delivered non-sterile. Please refer to the detailed requirements for reprocessing, described in section 6 Reprocessing.

5 Application
- The Arti-Grip™-sleeves are removed from the packaging in accordance with standard practice hygiene (use of disposable gloves).
- The sleeves should first be moistened with medical alcohol (ethanol or isopropyl alcohol, min. 70%), in order to facilitate to pull them more easily onto the clamping surfaces of the forceps.
- Subsequently, a strip of occlusion testing material (paper or foil) is clamped between the clamping surfaces of the forceps.
- The occlusion test material is inserted (buccal), with the aid of the forceps, between the respective teeth of the upper and lower jaw.
- The static or dynamic occlusion test is being performed.
- The occlusion test material is removed from the mouth and disposed of with the contaminated, normal practice waste.
- The silicone sleeves are removed from the clamping surfaces of the forceps.
- The forceps and the silicone sleeves are prepared for reprocessing.

6 Reprocessing
In principle, the following should be followed: "Hygiene requirements for the processing of medical devices: Recommendation of the Commission for Hospital Hygiene and Infection Prevention (KRINKO) at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM)," Federal Health Bulletin 2012 • 55:1244-1310.

The Arti-Grip™-Silicone sleeves must be reprocessed immediately after each treaᵐent (within a maximum of 2 hours):
- Thoroughly clean the Silicone sleeves by hand, using an appropriate small brush with firm bristles. The Silicone sleeves „Arti-Grip" must also be cleaned from the inside, using a suitable tool.
- Note: The cleaning should be carried out in a water bath, without any additional cleaning agent, below the water surface, to achieve a sufficient cleaning of Arti-Grip™-Silicone sleeves and to avoid protein fixation as well as to protect the environment against contamination with splashed water.
- Rinse the instruments with water (at least drinking water quality).
- Place the parts in a practice-customary cleaning and disinfecting bath. Examples:
 - o Becht Premium Concentrate
 - o Dürr Dental ID 213 instrument disinfection
 - o Pluradent Pluline instrument bath
 - o Schülke & Mayr gigasept® instru AF
 - o Note: "List of disinfectants and disinfecting procedures tested and approved by the Robert Koch Institute" or the VAH list of disinfectants.
 - o Note: To prevent protein fixation, formaldehyde-containing cleaning agents and disinfectants may be used only after adequate cleaning.
 - o Note: Strictly follow the instructions provided by the manufacturer of the cleaning agent or disinfectant. In particular, always follow the required concentration and residence time!
 - o Note: In case of automatic cleaning, the manufacturer's instructions have to be followed strictly!
- Final lavage of Arti-Grip™-Silicone sleeves with water (at least drinking water quality, recommendable: demineralized water with microbiological quality equal to drinking water)
- Drying.
- Visual check for corrosion, damaged surfaces, chipping, damage to shape and contamination. Damaged sleeves shall be discarded (limited number of reconditioning cycles see Section "Reusability"). In case of residual contamination, the entire cleaning procedure with all stages shall be repeated (cleaning, intermittent rinsing, disinfection, final rinsing and drying).
- The Silicone sleeves must be free of any residue and dry before further preparation.
- Silicone sleeves do not require servicing.
- Silicone sleeves should be packed and sealed in adequately sized, single-use sterilization bags that comply with EN 868-2ff ISO 11607 (suitable for steam sterilization). Follow the instructions of the manufacturer of the sterilization bags and sealing machines and the current standard requirements.
- Sterilization must be completed in a validated procedure using moist heat in an autoclave in accordance with DIN EN 13060 Type B and DIN EN 285 and ANSI AAMI ST79. Follow the instructions for use of the autoclave manufacturer.
- Sterilize the Silicone sleeves using moist heat (saturated steam) and a pre-vacuum procedure for **5 minutes at 134°C**.
- After completion of sterilization, the instrument parts must be stored dry and dust-free in the sealed sterilization packaging.
- The recommended storage period for sterile medical products is described in Standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.

7 Reusability
Frequent reconditioning does not have any effect or limit on the Arti-Grip™-Silicone sleeves, as the end of the product lifetime is determined by wear and damage due to use.

The use of damaged or soiled instruments is the responsibility of the user.
In case of non-compliance, all liability is excluded.

8 Disposal Instructions
Before disposing Arti-Grip™-Silicone sleeves, they either have to be sterilized or disposed of with the normal contaminated practice waste.

9 Notification of Incidents
Serious adverse events occurring in connection with the product must be reported to the manufacturer and the competent authority of the Member State.

10 Symbols

Manufacturer

Production Date

Steam sterilizable with saturated steam (autoclaving) at 134°C

„Follow manufacturers instructions"

REF Reference number

UDI Unique Device Identifier

CE Conformity with the relevant EU regulations

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Date of issue: 2020-03-27

General description and intended purpose of the medical product

Arti-Spot® is a contact color for testing the accurate fit of crowns, inlays, onlays, telescope crowns, clasps as well as friction surfaces of attachments and occlusal surfaces (**do not apply in oral cavity**), like gold or ceramics. Arti-Spot® can be applied to the surface to be tested with the enclosed brush (only for use in the laboratory) or with a disposable application brush, leaving a thin film of paint within a few seconds, when the solvent has evaporated.

The thin color skin is destroyed exactly at contact points and thus makes occlusion contacts or points of interferences visible when fitting crowns or telescopic crowns, for example.

For professional use only.

Indication

Arti-Spot® is a liquid contact color, available in the following colors: White (BK 85), Red (BK 86) and Blue (BK 87) for:

- Testing the accurate fit of crowns, inlays, onlays, telescope crowns, clasps as well as friction surfaces of attachments and occlusal surfaces (**do not apply in oral cavity**).

Contra indication

- Arti-Spot® contains solvents and therefore must not be applied directly invasively. The contact color must be applied outside the oral cavity on the surfaces of crowns, inlays, onlays, telescope crowns, clasps and frictions to be tested and these may only be applied to the oral cavity (invasively applied) after the color has dried (without solvent).

Side-effects

- Possible side-effects can be allergic reactions.
- Arti-Spot® can cause discoloration of prosthetic restorations made of plastic. It should only be used on metal or ceramic surfaces.

Notes

- The product must be stored in the sealed bottle at temperatures between 0°C and 40°C.
- The product contains a solvent that can evaporate quickly. After use, the bottle must be closed again immediately.
- The solvent contained in the product is highly flammable. Do not use the product near sources of ignition such as Bunsen burners (see also point Symbols).
- When used on patients, always apply the paint to the parts to be tested with a disposable brush outside the mouth.
- Arti-Spot® can be diluted with concentrated medical alcohol (ethanol 96%) to adjust the color viscosity.
- Arti-Spot® may cause discoloration of prosthetic restorations made of plastic. It should only be used on metal or ceramic surfaces.
- Arti-Spot® can be applied to Bausch Fleximeter-Strips® for combined usage.
- Stains on plastic surfaces or other materials can be removed with a dry paper towel. Color residues can also be removed with hot water or a steam jet. Tenacious paint residues can, if necessary, be dissolved and removed with concentrated alcohol (ethanol 96%).

Package content

- 1 bottle Arti-Spot®, contents 15ml (screw cap with holder for the brush)
- 1 brush
- Quick guide

Preparation

- Shake the bottle briefly before use (mixing balls are audible) so that the paint is mixed evenly; if the paint is too viscous, it can be diluted with medical alcohol (ethanol 96%, undenatured) at any time.

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Instructions for use
Arti-Spot® High-Spot Indicator
BK 85, BK 86, BK 87

Application

- Arti-Spot® is used in accordance with normal practice hygiene.

Color application (and start of application)

- Apply Arti-Spot® sparingly with a brush to the surface to be tested.
- Let the paint dry for approx. 10 seconds.
- Subsequently insert the part to be tested (crown, bridge, abutment of implants, etc.). The color skin is destroyed at the point of contact and the contact point is thus visible.
- Do not reuse the application brush, to avoid transmission of germs!
- After adjustment remove excessive color.

Effect/Diagnosis

After insertion of the surface to be tested into the oral cavity, interferences can be detected on the occlusal surfaces. The color skin is being destroyed on these occlusal surfaces. The base material then shines clearly through and highspots can easily be detected.

Reusability

- It is not permitted to reuse Single Use items (single use application brush). Safe use cannot be guaranteed if these products are used again, since there is an infection risk and / or the safety of the products is not liable.
- The use of damaged or soiled instruments is the responsibility of the user.
- In case of non-compliance, all liability is excluded.

Disposal instructions

Arti-Spot® can be disposed of with the usual practice waste.

Notification of incidents

Serious adverse events occurring in connection with the product must be reported to the manufacturer and the competent authority of the Member State.

Symbols

Manufacturer

Expiry date

„Follow manufacturer’s instructions“

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Permissible storage temperature range

Conformity with the relevant EU regulations

H224 Extremely flammable liquid and vapor
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P403 Store in well-ventilated place.
P235 Keep cool.

H319 Causes serious eye irritation.
P280 Wear protective gloves/clothing and eye/face protection.
P305 IF IN EYES:
P351 Rinse cautiously with water for several minutes.
P338 Remove contact lenses, if present and easy to do. Continue rinsing.
P337 If eye irritation persists:
P313 Get medical advice/attention.
H336 May cause drowsiness or dizziness.
P261 Avoid breathing vapour or spray.
P233 Keep container tightly closed.

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ENGLISH

General description and purpose of the medical device:
Green color paste for marking sore spots on dentures and partial dentures. **For professional use only.**

Indication
Bio-Ink® Flow is indicated for marking sore spots on the gum and subsequent marking on the denture.
Sore Spots also known as decubitus can occur frequently with new dentures or partial dentures. Also a degeneration of the alveolar ridge can cause sore spots as well as unbalanced occlusal forces on the denture. The adaption of the oral mucosa to the new denture often causes inflammatory lesions. Bio-Ink® Flow is a highly viscous paste. By punctiform, laminar or linear application to the mucosa the Palatal Vibrating Line and the mucobuccal fold can be marked. If necessary, the denture can be adjusted.

Contraindication
A direct application onto highly inflamed areas or purulent lesions of the oral mucosa areas should be avoided.

Side Effects
Possible side effects could be allergic reactions.

Notes
• The color paste can be applied directly with the single use cannulas and spread with the single use applicator (e.g. Magic-Brush™) to avoid cross-contamination.
• Alternatively, the paste can be used with a single use applicator (e.g. Magic-Brush™) and then applied to the surface to be tested.
• The mucous membranes, palate and mandible must be dried with clean pulp or cotton roll, previously.
• The set contains components (Magic-Brush™ applicator, single use cannulas and a Multi-Hold™ holder for brush). Except the Multi-Hold™ holder for brush, all components are for single use only and must not be reprocessed.
• The product must be stored in the sealed syringe at temperatures between 0°C und 40°C.

Package content
• 1 x Bio-Ink® Flow Paste, Syringe 1ml
• 10 single use Applicator Needles (cannulas)
• 1 holder for brush (Multi-Hold™ for Magic-Brush™)
• 50 brushes Magic-Brush™ (single use product)
• Instructions for use

Preparation
• Before use, the mucosa membrane, the palate and the mandible must be dried carefully with a pulp or a cotton roll so that the color paste does not run.
• If necessary, clean the holder for the applicator.
• If necessary, fit the applicator into the holder.

Application
Before use make sure that the mucosa, palate and alveolar ridge are dried carefully with a cellulose pulp or a cotton roll so that the color paste does not run. There should be no saliva on the surface to be tested. The denture must also be cleaned and dried carefully beforehand. Insert the applicator (Magic-Brush™) into the provided holder (Multi-Hold™).
The color paste can be applied directly with the disposable cannulas and spread with the disposable applicator (e.g. Magic-Brush™).
Alternatively, the color can be used with the disposable applicator (z.B. Magic-Brush™) and then applied to the surface to be tested.
To facilitate handling, the applicator can be angled at the tip.

Effect/Diagnosis
The color marking of the sore spot will be transferred to the denture when the denture is seated into position. When the denture is removed the spot is visible on the surface of the denture. If necessary the denture can be adjusted by the dentist.

Cleaning (of the surface to be tested)
After control, the color marking can be easily removed with water and a suitable surfactant (dishwashing detergent).

Reprocessing
In principle, the following should be observed: „Requirements for hygiene in the processing of medical devices Recommendation of the Commission for hospital hygiene and infection prevention (KRINKO) at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM)", Bundesgesundheitsblatt 2012 -55: 1244-1310.

• The instrument (holder) must be reprocessed immediately after each treatment (within a maximum of 2 h).
• Thoroughly clean the holder by hand with a small brush with firm bristles.
• Note: The cleansing should be done in a water bowl without further additives below the water surface in order to achieve adequate cleansing of the instrument parts whilst avoiding protein fixation as well as to protect the environment from contamination due to splash water.
• Intermediately rinse the holder with water (minimum drinking water quality).
• Place the instrument in a standard cleaning and disinfection bath. Examples: Becht Premium Konzentrat, Dürr Dental ID 213 Instrument disinfection, Pluradent instrument bath PLULINE, Schülke & Mayr gigasept® instru AF (see also "List of disinfectants and disinfection methods tested and approved by the Robert Koch Institute" or disinfectant list VAH).
• Note: Formaldehyde containing cleaning and disinfecting agents may only be used after sufficient cleaning in order to avoid protein fixation.
• Note: The instructions for use of the manufacturer of the cleaning and disinfecting agent must be strictly observed. In particular, the concentrations to be used and the reaction times must be observed!
• Note: For machine cleaning, the manufacturer's instructions for use must be followed carefully!
• Final rinse of the instrument parts with water (minimum drinking water quality, recommended: demineralized water with microbiological quality according to drinking water).
• Dry.
• Visual inspection for damaged surfaces, chipping and impurities. Damaged instruments must be discarded (For a limited restriction of reuse, see the chapter „Reusability"). If residual contamination still exists, the entire cleaning process must be repeated with all steps (Cleaning, intermediate rinsing, disinfection, final rinsing and drying).
• The instrument must be free of any residue and dry before further processing.
• Maintenance of the instruments is not required.
• The instrument is packaged and sealed in a sufficiently large single sterilization bag according to EN 868-2 ff ISO 11607 (suitable for steam sterilization). Observe the instructions of the manufacturers of the sterilization bags and the sealing machines and the current normative requirements.
• Sterilization must be carried out in a validated process with moist heat in an autoclave according to DIN EN 13060 Type B or DIN EN 285 resp. ANSI AAMI ST79. Observe the manufacturers operating instructions of the autoclave.
• Sterilize the instrument with humid heat (saturated water vapor) using a pre-vacuum procedure for 5 minutes at 134°C.
• The recommended storage life for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.

Reusability
• Frequent reprocessing restricts the life of the plastic instrument. The instrument must be disposed of in case of cracks in the plastic.
• It is not permitted to reuse single use items (applicator (Magic-Brush™) as well as single use cannulas). Safe use cannot be guaranteed if these products are used again, since there is an infection risk and / or the safety of the products is not liable.
• The use of damaged or soiled instruments is the responsibility of the user.
• In case of non-compliance, all liability is excluded.

Disposal Instructions
Bio-Ink®-Flow consists of a green food color and natural waxes, which are approved for food, cosmetics and pharmaceutical applications. All components of the kit can be disposed with the usual practice waste.

Notification of Incidents
Serious adverse events occurring in connection with the product must be reported to the manufacturer and the competent authority of the Member State.

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„Gebrauchsanweisung beachten"
„Follow manufacturers instructions"

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REF
Bestellnummer
Reference number

LOT
Chargen-Nummer
Batch code (Lot number)

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ENGLISH

General description and purpose of the medical device:
Green color paste for marking sore spots on dentures and partial dentures.

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Notes
• The color paste can be applied directly with the disposable cannulas (single use) and spread with the disposable applicator (e.g. Magic-Brush™) to avoid cross-contamination.
• Alternatively, the paste can be used with a single use applicator (e.g. Magic-Brush™) and then applied to the surface to be tested.
• The product must be stored in the sealed syringe at temperatures between 0°C und 40°C.

Package content
• 1 x Bio-Ink® Flow Paste, Syringe 1ml
• 5 single use cannulas
• Instructions for use

Application
Before use make sure that the mucosa, palate and alveolar ridge are dried carefully with a cellulose pulp or a cotton roll so that the color paste does not run. There should be no saliva on the surface to be tested. The denture must also be cleaned and dried carefully beforehand.

Effect/Diagnosis
The color marking of the sore spot will be transferred to the denture when the denture is seated into position. When the denture is removed the spot is visible on the surface of the denture. If necessary, the denture can be adjusted by the dentist.

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Instructions for Use

Y-Holder for Articulating Silk and Occlusion test films BK 140

Manufacturer
Dr. Jean Bausch GmbH & Co. KG
Oskar-Schindler-Str. 4
50769 Cologne - Germany
Tel: +49 221 70936-0
Fax: +49 221 70936-66
E-mail: info@bauschdental.de
Internet: www.bauschdental.de

Date of issue: 2019-12-06

1 General Description and purpose of the product

The Y-holder is a laboratory instrument, made of stainless steel, for easier handling of wide articulating silks or occlusion test foils in the articulator.
The Y-holder is not intended for use on humans and is not a medical device!
The Y-holder is sterilizable.

Indication
Laboratory-Instrument.

Contraindication
Not intended for use on humans!

Side effects
N/A.

2 Notes

- Before use, the Y-holder (instrument) has to be cleaned, disinfected and sterilized, if necessary.
- The instrument can be sterilized with steam (moist heat, 134°C).
- The instrument can be used in combination with articulating silks or occlusion test foils, produced by Dr. Jean Bausch GmbH & Co. KG.
- Disposal: The instrument must be disposed of with the usual laboratory waste.

3 Package content

- 1 Y-Holder (BK 140)
- Instructions for Use

4 Preparation

The instrument should be cleaned before use. This applies in particular to initial use after delivery, as the instrument is delivered non-sterile.

5 Application

- Insert the occlusion test material (e.g. silk or foil) into the holder so that it protrudes over the clamping jaws on all sides.
- Close the clamping jaws by pushing the locking slider on the handle forward.
- Perform the occlusion test on the model in the articulator.
- After the occlusion test has been completed, open the Y-holder by pushing the locking slide back and dispose of the used occlusion test material with the usual laboratory waste.
- If necessary, prepare the Y-holder for reprocessing.

6 Disposal Instructions

The instrument has to be disposed of with the normal laboratory waste.

7 Symbols

Manufacturer

Production Date

„Follow manufacturers instructions“

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Date of issue: 2019-12-05

1 General description and purpose of the medical device
Fleximeter®-Strips are flat, precisely cut silicone tabs in various colors with a disc- and rod-shaped, tapered end and defined thicknesses of 1.0mm, 1.5mm and 2.0mm.
Thickness precisions: 1.0mm ± 0.2mm, 1.5mm ± 0.2mm, 2.0mm ± 0.2mm

For professional use only.

Indication
Fleximeter®-Strips are used for specification of the preparation distance or the preparation height when preparing teeth for restorations (e.g. crowns, bridges and double crowns), and to determine the material thickness needed to adjust bite guard splints. Fleximeter®-Strips, along with the liquid Arti-Spot® or Arti-Spray® occlusion spray, can be used in combination to mark contact points. Dentists can use the Fleximeter®-Strips separately or in pairs, in the patient´s mouth as well as on models. Fleximeter®-Strips can also be used in the laboratory.

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Side Effects
Possible side effects may be specific allergic reactions. However, given the product's history of more than ten years, there has not been a verified documented report of allergic reactions. Side effects or interactions may occur if the products are used with new, unknown products or materials.

2 Notes
- Before each use, Fleximeter®-Strips must be cleaned, disinfected and sterilized.
- Fleximeter®-Strips (instruments) can be steam-sterilized (humid heat, 134°C).
- Fleximeter®-Strips can be used in combination with the products Arti-Spot® or Arti-Spray®.
- Before disposing of Fleximeter®-Strips, they either have to be sterilized or disposed of with the normal contaminated practice waste.

3 Package content
- Fleximeter®-Strips (non-sterile – Type: according to product label)
- Instructions for use

4 Preparation
The instrument must be cleaned, disinfected and sterilized before each use. This also applies in particular to the initial application after delivery, since the instrument is delivered non-sterile. Please refer to the detailed requirements for reprocessing, described in section 6 Reprocessing.

5 Application
- Fleximeter®-Strips are used taking into account standard dental hygiene procedures.
- To accurately determine occlusal measurements, Fleximeter®-Strips are inserted between the affected teeth of the upper and lower jaw.
- Fleximeter®-Strips can be used to change the occlusal spacing.
- To determine the raising of the bite, Fleximeter®-Strips are inserted between the affected teeth of the upper and lower jaw.
- Depending on the respective indication, the entire surface or one end of the Fleximeter®-Strips can be used.
- Fleximeter®-Strips can be coated with a color layer of Arti-Spot® or Arti-Spray®.

Caution:
- If a patient accidentally swallows Fleximeter®-Strips, do not induce vomiting!
- If a patient aspirates Fleximeter®-Strips, immediately take action for removal and call an ambulance if necessary.

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The instruments must be reprocessed immediately after each treatment (within a maximum of 2 hours):
- Thoroughly clean the instruments by hand, using an appropriate small brush with firm bristles.
- Note: The cleaning should be carried out in a water bath, without any additional cleaning agent, below the water surface, to achieve a sufficient cleaning of the Fleximeter-Strips and to avoid protein fixation as well as to protect the environment against contamination with splashed water.
- Rinse the instruments with water (at least drinking water quality).
- Place the parts in a practice-customary cleaning and disinfecting bath. Examples:
 - o Becht Premium Concentrate
 - o Dürr Dental ID 213 instrument disinfection
 - o Pluradent Pluline instrument bath
 - o Schülke & Mayr gigasept® instru AF
 - o Note: "List of disinfectants and disinfecting procedures tested and approved by the Robert Koch Institute" or the VAH list of disinfectants.
 - o Note: To prevent protein fixation, formaldehyde-containing cleaning agents and disinfectants may be used only after adequate cleaning.
 - o Note: Strictly follow the instructions provided by the manufacturer of the cleaning agent or disinfectant. In particular, always follow the required concentration and residence time!
 - o Note: In case of automatic cleaning, the manufacturer's instructions have to be followed strictly!
- Final lavage of Fleximeter-Strips with water (at least drinking water quality, recommendable: demineralized water with microbiological quality equal to drinking water)
- Drying.
- Visual check for corrosion, damaged surfaces, chipping, damage to shape and contamination. Damaged instruments shall be discarded (limited number of reconditioning cycles see Section "Reusability"). In case of residual contamination, the entire cleaning procedure with all stages shall be repeated (cleaning, intermittent rinsing, disinfection, final rinsing and drying).
- The instrument and all parts must be free of any residue and dry before further preparation.
- The instrument does not require servicing.
- All parts of the instrument should be packed and sealed in adequately sized, single-use sterilization bags that comply with EN 868-2ff ISO 11607 (suitable for steam sterilization). Follow the instructions of the manufacturer of the sterilization bags and sealing machines and the current standard requirements.
- Sterilization must be completed in a validated procedure using moist heat in an autoclave in accordance with DIN EN 13060 Type B and DIN EN 285 and ANSI AAMI ST79. Follow the instructions for use of the autoclave manufacturer.
- Sterilize the instrument using moist heat (saturated steam) and a pre-vacuum procedure for **5 minutes at 134°C**.
- After completion of sterilization, the instrument parts must be stored dry and dust-free in the sealed sterilization packaging.
- The recommended storage period for sterile medical products is described in Standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.

7 Reusability
Frequent reconditioning does not have any effect or limit on the use of the instrument, as the end of the product lifetime is determined by wear and damage due to use.

The use of damaged or soiled instruments is the responsibility of the user.

In case of non-compliance, all liability is excluded.

8 Disposal Instructions
Before disposing of Fleximeter®-Strips, they either have to be sterilized or disposed of with the normal contaminated practice waste.

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© 2019 Dr. Jean Bausch GmbH & Co. KG • 50769 Köln • Germany - © Bausch Articulating Papers, Inc. • Nashua, NH • U.S.A.