Box Size: 4.72 x 5.69 Overall Size: 27 x 25 Folds to: 3.375 x 5

COVIDIEN™


Bralon™, Polysorb™, Sofsilk™, Surgidac™

Single Stitch Reload

Surgidac™

Triple Stitch Reload

PT00133901 Page 1 of 2

en

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference to surgical techniques.

This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reuse, reprocessing and/or re-sterilization of this device may create the risk of contaminations, patient infection, permanent impairment or life-threatening injury. Do not reuse, reprocess or resterilize this device.

DESCRIPTION

The Bralon™, Polysorb™, Sofsilk™, and Surgidac™ reloads are available in a single stitch unit with one sutured needle. The Surgidac™ reloads are available in a triple stitch unit with three sutured needles.

Bralon™ braided nylon sutures are inert, nonabsorbable, sterile, surgical sutures composed of the long-chain aliphatic polymers Nylon 6 and Nylon 6.6. Bralon™ sutures are coated uniformly with polydiethylene adipate to enhance handling characteristics, ease of passage through tissue and reduction of capillarity. Bralon™ braided sutures are nonabsorbable nylon surgical sutures which are available either dyed black, with Logwood extract, or undyed (white).

Bralon™ braided nylon sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopoeia (EP).

Polysorb™ sutures are absorbable, sterile, surgical sutures composed of Lactomer™ glycolide/lactide copolymer which is a synthetic polyester composed of glycolide and lactide (derived from glycolic and lactic acids). Polysorb™ sutures are prepared by coating the suture with a mixture of a copolymer of glycolide and lactide and calcium stearoyl lactylate.

Polysorb™ sutures are colored violet to increase visibility and are also available undyed. Polysorb™ sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopoeia (EP) except for minor variations in suture diameters. Such variations are:

Maximum suture oversize in diameter (mm) from USP size	USP size designation (mm)	Maximum average (mm)
4 – 0	0.15 – 0.199	0.050
3 – 0	0.20 – 0.249	0.050
2 – 0	0.30 – 0.339	0.050
0	0.35 – 0.399	0.050

Sofsilk™ silk sutures are nonabsorbable, sterile, non-mutagenic surgical sutures composed of silkworm cocoon silk fibers called fibroin. This protein is derived from the domesticated silkworm species Bombyx mori of the family bombycidae. The silk fibers are treated to remove the naturally occurring wax, gums, bleached and braided to produce Sofsilk™ surgical silk sutures. The braided sutures are coated uniformly with silicone to reduce capillarity, and to increase surface lubricity to enhance handling characteristics, ease of passage through tissue and knot run-down properties.

Sofsilk™ sutures are available colored black with Logwood extract.

Sofsilk™ sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopoeia (EP) for nonabsorbable surgical sutures.

Surgidac™ polyester sutures are nonabsorbable, sterile, surgical sutures composed of Polyethylene terephthalate). They are prepared from fibers of high molecular weight, long-chain, linear polyesters having recurrent aromatic rings as an integral component. The monofilament sutures are uncoated and the braided sutures are available uncoated or coated uniformly with polybutylene adipate to increase surface lubricity, thereby enhancing handling characteristics, ease of passage through tissue, and knot run-down properties. Surgidac™ polyester sutures are inert and elicit only minimal local tissue reaction. They are available dyed with D&C Green No. 6 (colorant ant are also available undyed (white).

Surgidac™ sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopoeia (EP) for nonabsorbable surgical sutures.

NOTE: Clinical evidence has not been established for use in pediatric population.

RELOAD INDICATIONS

The Bralon™, Polysorb™, Sofsil™, and Surgidac™ reloads are intended for use with the Endo Stitch™ and SILS™ Stitch suturing devices. The Endo Stitch™ and SILS™ Stitch suturing devices have application in endoscopic surgery for the placement of interrupted or running sutures in soft tissue.

SUTURE INDICATIONS

Bralon™ braided nylon sutures, Sofsilk™ sutures and Surgidac™ polyester sutures are indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological surgery.

Polysorb™ sutures are indicated for use in soft tissue approximation or ligation and ophthalmic surgery, but not in cardiovascular or neural tissue.

ACTIONS

Bralon™ braided nylon sutures, Sofsilk™ sutures and Surgidac™ polyester sutures elicit a minimal acute inflammatory reaction in tissue, which is followed by gradual encapsulation of the suture by fibrous connective tissue. Bralon™, as with any nylon suture, though not absorbed, progressive hydrolysis of the suture, in vivo, may result in gradual loss of its tensile strength over time.

Polysorb™ - Progressive loss of tensile strength and eventual absorption of Polysorb™ sutures occurs by means of hydrolysis, where the Lactomer™ glycolide/lactide copolymer is broken down by glycolic and lactic acids which are subsequently absorbed and metabolized by the body. Absorption begins as a loss of tensile strength without appreciable loss of mass. Studies indicate tensile strength averages for Polysorb™ sutures are approximately 140% of USP and EP minimum knot strength initially, and approximately 80% at two weeks and 10% at three weeks post implant. Absorption of Polysorb™ sutures is essentially complete between the 56th and 70th day.

CONTRAINDICATIONS

This device is intended for use only as indicated. It is not intended for use when endoscopic techniques generally are contraindicated.

SUTURE CONTRAINDICATIONS

Bralon™ - Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, nylon sutures should not be used where permanent retention of tensile strength is required.

Polysorb™ - Polysorb™ sutures, being absorbable, should not be used where extended approximation of tissue is required.

Sofsilk™ - The use of this suture is contraindicated in patients with known sensitivities or allergies to silk.

Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, Sofsilk™ sutures should not be used where permanent retention of tensile strength is required.

Surgidac™ None known.

WARNINGS AND PRECAUTIONS

1. Ensure that the toggle levers are fully retracted prior to opening the jaws of the device.
2. The reloading button should never be pressed when the instrument is in the body cavity, as this will release the needle.
3. Inspect the application site to ensure hemostasis. Place additional sutures or use electrocautery if necessary to complete hemostasis.
4. Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. A thorough understanding of the operating principles, risks versus benefits, and the hazards involved in utilizing any endoscopic approach is necessary to avoid possible injury to the user and/or patient.
5. Verify mechanical and electrical compatibility of devices from different manufacturers prior to using them together in a procedure.
6. As with any foreign material in the presence of bacterial contamination that may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of contaminated or infected wounds.
7. This device is provided STERILE and is intended for use in a SINGLE procedure only. DISCARD AFTER USE. DO NOT RESTERILIZE.

SUTURE PRECAUTIONS

1. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.
2. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.
3. Skin sutures that must remain in place longer than 7 days may cause irritation and should be snipped off or removed as indicated. Under some circumstances, immobilization by external support may be employed at the discretion of the surgeon.

SUTURE WARNINGS

1. Bralon™ braided nylon sutures, Sofsilk™ sutures, Surgidac™ polyester sutures: As with any foreign body, prolonged contact of any suture with salt solutions may result in calculus formation. Do not resterilize. Discard opened, unused sutures. Users should be familiar with surgical procedures and techniques involving nonabsorbable sutures before employing Bralon™ sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.
2. Polysorb™ - Do not resterilize. Discard opened, unused sutures. Store at room temperature and avoid prolonged exposure to elevated temperatures. In surgery of the urinary tract, Polysorb™ sutures may be employed if such a tract is, or any other, suture with salt solutions, as calculus formation may result. Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing Polysorb™ sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. As any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice must be followed with respect to drainage and closure of contaminated or infected wounds. The use of this suture may be inappropriate in patients with any conditions that, in the opinion of the surgeon, may cause or contribute to delayed wound healing. As this is an absorbable suture material, the use of supplemental support should be considered by the surgeon in closure of the abdomen, chest, joints or other sites subject to expansion or requiring additional support.

ADVERSE REACTIONS

Bralon™ braided nylon sutures, Sofsilk™ sutures, Surgidac™ polyester sutures: Adverse effects associated with the use of this product are: wound dehiscence, gradual loss of tensile strength over time, calculi formation when prolonged contact with salt solutions occurs, enhanced bacterial infectivity, minimal acute inflammatory reaction, and transitory local irritation.

Polysorb™ - Adverse effects associated with this product are: wound dehiscence, failure to provide adequate wound support in sites where expansion, stretching or distention occur, failure to provide adequate wound support in patients with conditions which may delay wound healing, localized irritation when the suture are left in place for greater than 7 days, calculi formation when prolonged contact with salt solutions occurs, enhanced bacterial infectivity, minimal acute inflammatory reaction, and pain, edema and erythema at the wound site.

1 SCHEMA/KEY

A) TRIPLE STITCH RELOAD
B) SUTURE SPOOLS
C) DISPENSER CAP
D) DISPENSER CAP ARROWS
E) NEEDLES
F) RELOAD NOTCH
G) BODY
H) Bralon™, Sofsilk™, and Surgidac™ SINGLE STITCH RELOAD
I) TAB
J) Polysorb™ SINGLE STITCH RELOAD
K) PAPER INSERT

2 LOADING

1. For Polysorb™ reloads, remove the paper insert prior to loading.
 J) Polysorb™ RELOAD
 K) PAPER INSERT
2. To load the Endo Stitch™ suturing device, insert the tip of the open jaws, with the logo facing up, into the single use loading unit and press firmly until fully seated.
3. Prior to loading the SILS™ Stitch instrument, ensure that the distal tip of the instrument is in the straight position and is aligned so that the logo (on the distal tip) and printing (on the handle) are both facing up. Please refer to the SILS™ Stitch instrument instructions for use for device-specific loading and unloading instructions. Squeeze the handles together and pull both toggle levers back until you hear an audible click.

NOTE: The Black Loading Buttons will cover the red indicator boxes.

L) TOGGLE LEVERS
M) BLACK LOADING BUTTONS
N) RED INDICATOR BOXES

4. Prior to removal of the device from the reload cartridge, ensure that the jaws are in a closed position. Press on the dispenser cap (triple stitch reload) or the tab (single stitch reload) located on the top of the cartridge. Remove the suturing device and entire length of suture.

NOTE: Ensure that the needle is fully seated within the jaws of the suturing device by opening the jaws and inspecting the needle.

5. If using the Surgidac™ triple stitch reload, expose a new sutured needle by rotating the dispenser cap until the arrow on the dispenser cap aligns with the reload notch on the body.
6. (I) DISPENSER CAP ARROWS
 (E) RELOAD NOTCH

3 UNLOADING

1. To unload the Endo Stitch™ suturing device and SILS™ Stitch suturing device, remove from the body cavity.

NOTE: Prior to unloading the SILS™ Stitch instrument, ensure the distal tip of the instrument is in the straight position. To facilitate unloading, slide the toggle levers forward.

2. Align toggle levers with the arrows on the body housing, simultaneously squeeze the handles and press the black reloading buttons on the front and back of the device and slide them forward to expose the red indicator boxes.

L) TOGGLE LEVERS
M) BLACK LOADING BUTTONS (front and back)
N) RED INDICATOR BOXES
O) ARROWS
P) 5X7" Stitch
Q) Endo Stitch™

3. Release the handles and remove the needle. The suturing device is now ready to be reloaded.

CAUTION: Do not depress reloading buttons while the suturing device is in the body cavity as this will release the needles.

RELOADING

To reload the Endo Stitch™ and SILS™ Stitch suturing devices, follow the instructions for loading.

AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.

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