
ChM®

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ROTARY REUSABLE SURGICAL INSTRUMENTS DRILL BITS, REAMERS, KIRSCHNER WIRES, CUTTERS, TAPS, TREPHINES, COUNTERSINKS

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1 INDICATIONS

1. Surgical drill bits, reamers, Kirschner wires, cutters, taps, trephines and countersinks are specially designed surgical accessories to help the surgeon cut, drill and thread the bony structures. They belong to a group of rotary surgical instruments used with a drilling machine or other active drive.
2. Surgical and orthopaedic instruments are intended for use only by skilled and trained medical professionals who are familiar with their use and application.

2 DESCRIPTION

1. The unit package contains one piece of the product in non-sterile condition. Clear plastic bags are a typical packaging material. The products may also be supplied as a complete set (*arranged on palettes and placed into specially designed sterilization containers*).
2. The package is equipped with the product label. The label (*as a primary label*) contains, among others:
 1) Logo **ChM** and the address of the manufacturer.
 2) Catalogue number (*REF*), e.g.: 40.XXXX.XXX, and device name and size.
 3) Production batch number (*LOT*), e.g. XXXXXXX.
 4) NON-STERILE sign - indicates non-sterile product.
 5) Information symbols (*described in the footer of this Instructions For Use*).
 6) CE conformity mark and the Notified Body number (*0197*).
3. Depending on the size or type of the product, the following information may be marked on its surface: manufacturer's logo, production batch no. (*LOT*), catalogue no. (*REF*), type of material and device size.

3 MATERIALS

1. The instruments are manufactured from:
 1) Corrosion-resistant steel. The protective layer (*passive layer*) against corrosion is formed on the surface of the device due to high content of chromium.
 2) Implantable steel.
 3) Implantable cobalt alloy.
2. If the material of the device cannot be specified, please contact **ChM** representative.

4 WARNINGS AND PRECAUTIONS

1. Improper, careless and inconsistent with the recommendations provided below handling of the instruments can lead to their chemical, electrochemical or mechanical damage which can adversely affect corrosion resistance and shorten the service life of the devices.
2. Instruments are intended only for specific procedures and must be used strictly according to their intended purpose. Use of instruments not in accordance with their intended purpose may lead to malfunction, accelerated wear and, in consequences, damage to the instrument.
3. The surgeon should be familiar with all components of the device before use and should personally verify if all components and instruments are present before the surgery begins.
4. Before the procedure begins, all instruments should be carefully inspected for their condition and proper functioning. They should be undamaged and without any signs of corrosion. Blades and cutting edges should be sharp and undamaged. Damaged or corroded instruments should be immediately replaced. The use of bent, damaged or corroded instruments is not allowed.
5. The surgeon should verify if the rotary instrument has been properly inserted and tightened into the drive before activation to avoid migration and possible injury.
6. Tissue structures close to the operative site must be protected.
7. Operating room personnel should wear proper eye protection during the surgical procedure requiring the use of rotary instruments and medical drive.
8. Proper drill guides that are available in instrument sets should be used.
9. Firm control of the drive with the instrument attached must be maintained to avoid injury to the patient or operating room personnel.
10. Manual cooling with the use of physiological saline should be applied.
11. Do not apply excessive force when using the instrument – it may lead to its permanent damage and, in consequences, to malfunction of the device.
12. Do not use counterclockwise rotation with an intramedullary reamer. The maximum torque for intramedullary reamers is:

1) 3 Nm for working diameters of ≤ 7.5 mm.
2) 10 Nm for working diameters above 7.5 mm.

13. Avoid excessive speed as it can cause temporary rise in temperature of adjacent bone and surrounding tissues above normal physiological level. It may result in tissue necrosis.
14. Instrument (*e.g. drill bit*) breakage occurs mainly when working around metallic implants (*e.g. fixation plates, intramedullary nails*). To prevent these incidents, ensure that the tip of the drill bit does not have contact with any implant.
15. Collision of the instrument with metal operating equipment, retractor or other device may cause damage that necessitates intraoperative replacement of that instrument.
16. Swarfs and possible metal particles should be carefully filtered by suction.
17. Instruments are subject to constant wear processes. While rare, intraoperative fracture or breakage of the instrument can occur. Instruments which have been subjected to prolonged use or excessive forces are more susceptible to fractures, depending on care taken during surgery and the number of procedures performed. Should breakage occur, the instrument parts must be removed and disposed of immediately in accordance with valid medical facility procedures.
18. In order to confirm the removal of all undesired metal fragments from the surgical field, intraoperative X-Ray examination is recommended.
19. In the case of suspected or documented allergy or intolerance to metallic materials, surgeon shall find out if the patient develops allergic reaction to the instrument material by ordering appropriate tests.
20. Instrument which had contact with tissues or body fluids of another patient cannot be re-used prior to its reprocessing due to a potential risk of cross-infection caused by viruses, bacteria and prions.

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1. Prior to use of a non-sterile device, the following rules apply:
 1) The device must undergo cleaning, disinfection and sterilization procedures.
 2) Effective cleaning is a complicated procedure depending on the following factors: the quality of water, the type and the quantity of used detergent, the technique of cleaning (*manual, automated*), the proper rinsing and drying, the proper preparation of the device, the time, the temperature and carefulness of the person conducting this process, etc.
 3) The hospital facility remains responsible for the effectiveness of the conducted cleaning, packaging and sterilization processes with the use of existing equipment, materials and properly trained personnel.
2. Preparation at the place of use.
 1) Immediately after use, remove from instrument blood and other contaminants with disposable cloth or paper towels. Additionally, it is recommended to rinse the instrument under running water or to place it in the aqueous disinfectant solution. Do not let blood, tissues, body fluids or other biological impurities dry out on the surface of the device.
 2) In order to prevent blood and debris from drying out on the instrument surface, transport the product to the processing area in a closed container or covered with a damp cloth.
 3) In order to avoid contamination during transportation, the dirty instruments should be separated from the clean ones.
3. Preparation for washing and disinfection (*for all methods*).
 1) The used instruments should be reprocessed as soon as possible.
 2) If the instrument can be disassembled, it must be done before cleaning processes.
 3) Rinse under running water and remove surface debris using a disposable cloth, paper towel or plastic brushes (*nylon brushes are recommended*). Particular attention should be paid to openings and places difficult to be cleaned. Very dirty devices should be soaked in an aqueous solution of a detergent or a washing-disinfecting agent, e.g. neodisher® MediClean forte, at temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 4) CAUTION: It is forbidden to use brushes made of metal, bristles or materials which could damage the product.
4. Cleaning and disinfection process.
 1) This Instructions for Use describes two **ChM**-approved cleaning and disinfection methods: manual with ultrasound cleaning and automated method. It is recommended to use automated cleaning and disinfection procedures (*in a washer-disinfector*).
 2) The chosen washing and disinfecting agents must be suitable and approved for use with medical devices. It is important to follow the instructions and restrictions specified by the producer of those cleaning agents. It is recommended to use aqueous solutions of washing-disinfecting agents with a pH value between 10.4 and 10.8. **ChM** used the following materials during the validation process of the described recommendations for cleaning and disinfection. It is allowed to use other materials than those listed below which may also give a comparable effect:
 a) detergent - Dr.Weigert (*producer*) neodisher® MediClean forte (*name of the detergent*);
 b) disinfectant - Dr.Weigert (*producer*) neodisher® Septo Active (*name of disinfectant*).

3) To prevent product damage (*pitting, rust, discoloration*), do not use aggressive cleaning agents (*NaOH, NaOCl*), saline solutions and unsuitable cleaning agents.
4) Where possible, it is recommended to use demineralized water to avoid the formation of spots and stains caused by chlorides and other compounds present in ordinary water.
5) Manual with ultrasound cleaning.
 a) Equipment and materials: a device for ultrasound cleaning, soft, lint-free cloths, plastic brushes, syringes, aqueous solutions of cleaning agent.
 b) Manual cleaning: Initial manual cleaning must be performed prior to ultrasound cleaning.
 c) Rinse under running water until the product is visually clean. Use plastic brushes to remove heavy or large debris.
 d) Soak the product for at least 10 minutes in an aqueous solution of a detergent at temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 e) Rinse the product under cold water for at least 2 minutes, paying particular attention to the holes and places difficult to be cleaned.
 f) Prepare fresh washing solution. Clean the surfaces and gaps of the product, carefully. Use suitable brushes to clean the holes. Clean the product immersed in the solution.
 g) Rinse the product thoroughly under warm running water for at least 2 minutes, paying special attention to the gaps, blind holes, hinges and joints. When cleaning, use brushes and perform multiple reciprocating movements on the surface of the product.
 h) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-h until the product is visually clean.
 i) Ultrasound cleaning: prepare an aqueous cleaning solution at a temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the cleaning agent, in respect of temperature, concentration, exposure time and water quality*). Immerse fully the product in the aqueous cleaning solution and have it washed in ultrasounds for 15 minutes.
 j) Rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned.
 k) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-k until the product is visually clean.
 l) Use demineralized water for final rinsing of the device.
 m) Dry the device thoroughly using disposable, soft, lint-free cloth or compressed air.
 n) Prepare an aqueous solution of disinfecting agent at a temperature of 20+/-2°C using 20g of the agent per 1 liter of water. Immerse the product in the solution, exposure time – 15 min (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 o) After the exposure time, rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned.
 p) The cannulated instruments should be treated using a compressed air or air supplied from the syringe.
 q) Dry the device thoroughly. It is recommended to dry the product in a dryer at a temperature ranging from 90°C to 110°C, for 60 min.
 r) Visually inspect the entire surface of the device.
 s) CAUTION: If the obstruction in the cannula cannot be removed as indicated in the Instructions for Use, the device should be considered at the end of its useful life and should be discarded in accordance with facility procedures and guidelines.
6) The automated method using a washer - disinfector.
 a) Equipment and materials: a washer - disinfector, aqueous solutions of cleaning agent.
 b) Cleaning in the washer-disinfector must be preceded by a manual and ultrasound cleaning, following the procedure described in subsections c-h of paragraph 5.
 c) CAUTION: The equipment used for washing/disinfection should meet the requirements of ISO 15883. Procedure of washing in the washer-disinfector shall be performed according to internal hospital procedures, recommendations of the washer-disinfector manufacturer, and instructions for use prepared by the washing-disinfecting agent manufacturer.
 d) The device should undergo the process of machine washing in the washer-disinfector using the following cycle parameters: (*1*) - pre-washing in cold tap water, duration – 2 min; (*2*) - washing in an aqueous solution of cleaning agent at 55+/-2°C and pH of 10.4 - 10.8, duration – 10 min; (*3*) - rinsing under demineralized water, duration – 2 min; (*4*) - thermal disinfection in demineralised water at 90°C, minimal duration – 5 min; (*5*) - drying at the temperature ranging from 90°C to 110°C, duration - 40 min.

5. Inspection
 1) Each time before re-use and re-sterilization, all medical devices should be inspected.
 2) Generally unmagnifed visual inspection under good light conditions is sufficient.
 3) All parts of the product should be checked for visible dirt and corrosion. Particular attention should be paid to:
 a) Places where dirt can be found, such as joints, latches, etc.
 b) Holes, grooves and gaps the debris could have been pressed into during use, e.g. drill flutes adjacent to the cutting tip.
 4) Each time before re-use and re-sterilization, the functional check of the product should be performed, consisting of:
 a) Verifying the connections in the mating instruments, such as tips, shafts and quick coupling devices.
 b) Verifying the connections in the mating instruments, such as drill bits with drives.
 c) Verifying all rotating devices for straightness (*this can be simply achieved by rolling the device on a flat surface*).
 d) Verifying „*flexible*" *devices (e.g. intramedullary reamers*) for damage to the spiral elements.
 e) Verifying cutting edges for sharpness.
 f) Verifying instruments for damage to material structure (*cracks, dents, peels, etc.*).
 5) Damaged or defective product cannot be approved for further use.
 6) Prior to storage, the instrument must be checked for dryness.
 7) CAUTION:
 a) The **ChM** does not define the maximum number of uses appropriate for re-usable medical instruments. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful and proper use reduces the risk of damage to the product and extends its serviceable life.
 b) The manufacturer does not recommend using any preservatives on medical devices.
6. Packaging
 1) Washed and dried devices shall be stored (*if possible*) in suitable stands placed in special sterilization containers. Separate items should be packed in a packaging intended for the recommended steam sterilization. Sterilization containers, item packaging and packaging process itself have to meet the requirements of ISO 11607 standards. The packaging procedure must be performed in controlled purity conditions. The device must be packed so that during its removal from the packaging, when used, there is no risk for its re-contamination
7. Sterilization
 1) The instrument is washed, disinfected and dried to undergo a sterilization process. The recommended sterilization procedure for non-sterile medical devices validated by **ChM** is autoclave vacuum steam sterilization (*with water vapor under overpressure*) with the following parameters:
 a) temperature: 134°C,
 b) minimum exposure time: 5 min,
 c) minimum drying time: 20 min.
 2) In order to minimize the potential risk of transmission of Creutzfeldt-Jakob disease (*especially for surgical instruments that could come into contact with the central nervous system*), national health authorities in some countries recommend sterilization according to the following parameters:
 a) temperature: 134°C,
 b) minimum exposure time: 18 min,
 c) minimum drying time: 20 min.
 3) CAUTION:
 a) The sterilization process used in the medical facility must be validated and routinely controlled in accordance with the requirements of EN ISO 17665-1.
 b) Sterilization must be effective and in accordance with requirements of the EN 556-1 standard to ensure the required level of guaranteed sterility SAL 10^{-6} (*where SAL stands for Sterility Assurance Level*).
 c) Device must not be sterilized in the packaging in which it was delivered, except specially designed sterilization containers.
 d) The method of sterilization using ethylene oxide, gas plasma and dry heat should not be used, unless the Instructions for Use for the product contains sterilization recommendations using these methods.
 e) The sterilization temperature for plastic products (*PPSU, PEEK, PTFE, silicone*) cannot be higher than 140°C.

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7 COMPATIBILITY

1. **ChM** specialist instrument sets are designed for insertion of **ChM** implants. A specific, illustrated operating technique that describes the proper use of instruments included in the instrument set that is designed for particular implant system, is provided together with such instrument set. It is not allowed to combine **ChM** instruments with products from other manufacturers. The physician bears all responsibility for the use of the **ChM** instruments together with implants and instruments from other manufacturers.

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EN

INSTRUCTIONS FOR USE

Important product information for

REUSABLE ORTHOPAEDIC AND SURGICAL INSTRUMENTS TRIALS

1 INDICATIONS *4*

2 DESCRIPTION *4*

3 MATERIALS *4*

4 WARNINGS AND PRECAUTIONS *5*

5 CLEANING, DISINFECTION, STERILIZATION *5*

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1 INDICATIONS

1. Trials function differently depending on the operated segment of the human skeleton. **ChM** produces the following types of trials:
 1) Trials for spine intervertebral cages
 a) Cages trials are used for the initial assessment of the size of the intervertebral space in order to facilitate the selection of the proper implant size (*the intervertebral cage*) in spine intervertebral stabilization procedures.
 2) Trials for sterno-costal plate
 a) Plate trials are used for initial rough assessment of the size and shape of the patient's chest in order to facilitate the selection of the proper implant size (*sterno-costal plate*), and for the initial determination of the shape of the implant to match the shape of the chest, in the treatment procedures for the chest deformity.
 3) Trials for plates
 a) Plate trials are used for determination of length and/or shape of the implant (*plate*) to match the implantation site and the type of fracture/correction.
 4) Trials for rod
 a) Rod trials are used for initial rough assessment of the size and shape of the rod and to facilitate the selection of the proper size of the spinal rod, in the spinal stabilization procedures using transpedicular screws.
 5) Trials for nails
 a) Nails trials are used for the proper implant selection. They are used to reflect the initial curvature of the bone or initial assessment of the nail length.
 6) Trials for head of radial head prosthesis
 a) Head trials are used as "*duplicates*" of head prosthesis in the initial phase of the surgery. Test fitting of the diameter and height of the head allows for intraoperative selection of proper implant size.
2. Surgical and orthopaedic instruments are intended for use only by skilled and trained medical professionals who are familiar with their use and application.

2 DESCRIPTION

1. The unit package contains one piece of the product in non-sterile condition. Clear plastic bags are a typical packaging material. The products may also be supplied as a complete set (*arranged on palettes and placed into specially designed sterilization containers*).
2. The package is equipped with the product label. The label (*as a primary label*) contains, among others:
 1) Logo **ChM** and the address of the manufacturer.
 2) Catalogue number (*REF*), e.g.: 40.XXXX.XXX, and device name and size.
 3) Production batch number (*LOT*), e.g. XXXXXXX.
 4) NON-STERILE sign - indicates non-sterile product.
 5) CE conformity mark and the Notified Body number (*0197*).
 7) Information symbols (*described in the footer of this Instructions For Use*).
3. Depending on the size or type of the product, the following information may be marked on its surface: manufacturer's logo, production batch no. (*LOT*), catalogue no. (*REF*), trial material, trial size and word „*Trial*", *additionally, the information* „*Do not implant*", which is permanently placed on the surface of the product.

3 MATERIALS

1. Trials are made of: titanium, titanium alloys, stainless steel, aluminum alloys and plastics - mainly PPSU (*Polyphenylsulfone*), PEEK (*Polyetheretherketone*) and silicone. These materials are approved for production of surgical instruments in accordance with applicable procedures.
2. Instruments manufactured from titanium and its alloys are characterized by: high strength, small weight and high resistance to corrosion. Their surfaces are anodized.
3. Devices made of aluminum form protective oxide layer on their surface as a result of electrochemical treatment. The surface may be dyed or stays in natural colour (*dark grey*).
4. Devices made of aluminium with processed layer have good corrosion resistance. However, the contact with strong alkaline cleaning and disinfecting agents, solutions containing iodine or some metal salts, due to chemical interference with the processed aluminium surface, shall be avoided.
5. If the material of the device cannot be specified, please contact **ChM** representative.

4 WARNINGS AND PRECAUTIONS

1. Instruments are intended for use only by skilled and trained medical professionals who are familiar with their use and application.
2. Improper, careless and inconsistent with the recommendations provided below handling of the instruments can lead to their chemical, electrochemical or mechanical damage which can adversely affect corrosion resistance and shorten the service life of the devices.
3. The surgeon should be familiar with all components of the device before use and should personally verify if all components and instruments are present before the surgery begins.
4. Tissue structures close to the operative site must be protected.
5. Collision of the instrument with metal operating equipment, retractor or other device may cause damage that necessitates intraoperative replacement of that instrument.
6. Do not apply excessive force when using the instrument – it may lead to its permanent damage and, in consequences, to malfunction of the device.
7. Deformable trials (*e.g. plate, rod trials*) have limited lifetime which depends on, e.g. the degree of deformation applied, the number of bendings performed, etc. The breakage of the instrument should be considered the end of the product lifetime.
8. Instruments are subject to constant wear processes. While rare, intraoperative fracture or breakage of the instrument can occur. Instruments which have been subjected to prolonged use or excessive forces are more susceptible to fractures, depending on care taken during surgery and the number of procedures performed. Should breakage occur, the instrument parts must be removed and disposed of immediately in accordance with valid medical facility procedures.
9. In order to confirm the removal of all undesired metal fragments from the surgical field, intraoperative X-Ray examination is recommended.
10. Instruments are intended only for specific procedures and must be used strictly according to their intended purpose. Use of instruments not in accordance with their intended purpose may lead to malfunction, accelerated wear and, in consequences, damage to the instrument.
11. Use of trials manufactured by **ChM** together with implants and instruments manufactured by other companies may lead to incorrect readings, damage or destruction of trials or implants/instruments and improper course of surgery and healing process.
12. Instrument which had contact with tissues or body fluids of another patient cannot be re-used prior to its reprocessing due to a potential risk of cross-infection caused by viruses, bacteria and prions.
13. In the case of suspected or documented allergy or intolerance to metallic materials, surgeon shall find out if the patient develops allergic reaction to the instrument material by ordering appropriate tests.

5 CLEANING, DISINFECTION, STERILIZATION

1. Prior to use of a non-sterile device, the following rules apply:
 1) The device must undergo cleaning, disinfection and sterilization procedures.
 2) Effective cleaning is a complicated procedure depending on the following factors: the quality of water, the type and the quantity of used detergent, the technique of cleaning (*manual, automated*), the proper rinsing and drying, the proper preparation of the device, the time, the temperature and carefulness of the person conducting this process, etc.
 3) The hospital facility remains responsible for the effectiveness of the conducted cleaning, packaging and sterilization processes with the use of existing equipment, materials and properly trained personnel.
2. Preparation at the place of use.
 1) Immediately after use, remove from instrument blood and other contaminants with disposable cloth or paper towels. Additionally, it is recommended to rinse the instrument under running water or to place it in the aqueous disinfectant solution. Do not let blood, tissues, body fluids or other biological impurities dry out on the surface of the device.
 2) In order to prevent blood and debris from drying out on the instrument surface, transport the product to the processing area in a closed container or covered with a damp cloth.
 3) In order to avoid contamination during transportation, the dirty instruments should be separated from the clean ones.
3. Preparation for washing and disinfection (*for all methods*).
 1) The used instruments should be reprocessed as soon as possible.
 2) If the instrument can be disassembled, it must be done before cleaning processes.
 3) Rinse under running water and remove surface debris using a disposable cloth, paper towel or plastic brushes (*nylon brushes are recommended*). Particular attention should be paid to openings and places difficult to be cleaned. Very dirty devices should be soaked in an aqueous solution of a detergent or a washing-disinfecting agent, e.g. neodisher® MediClean forte, at temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).

4) CAUTION: It is forbidden to use brushes made of metal, bristles or materials which could damage the product.

4. Cleaning and disinfection process.

1) This Instructions for Use describes two **ChM**-approved cleaning and disinfection methods: manual with ultrasound cleaning and automated method. It is recommended to use automated cleaning and disinfection procedures (*in a washer-disinfector*).
2) The chosen washing and disinfecting agents must be suitable and approved for use with medical devices. It is important to follow the instructions and restrictions specified by the producer of those cleaning agents. It is recommended to use aqueous solutions of washing-disinfecting agents with a pH value between 10.4 and 10.8. **ChM** used the following materials during the validation process of the described recommendations for cleaning and disinfection. It is allowed to use other materials than those listed below which may also give a comparable effect:
 a) detergent - Dr.Weigert (*producer*) neodisher® MediClean forte (*name of the detergent*);
 b) disinfectant - Dr.Weigert (*producer*) neodisher® Septo Active (*name of disinfectant*).
3) To prevent product damage (*pitting, rust, discoloration*), do not use aggressive cleaning agents (*NaOH, NaOCl*), saline solutions and unsuitable cleaning agents.
4) Where possible, it is recommended to use demineralized water to avoid the formation of spots and stains caused by chlorides and other compounds present in ordinary water.
5) Manual with ultrasound cleaning.
 a) Equipment and materials: a device for ultrasound cleaning, soft, lint-free cloths, plastic brushes, syringes, aqueous solutions of cleaning agent.
 b) Manual cleaning: Initial manual cleaning must be performed prior to ultrasound cleaning.
 c) Rinse under running water until the product is visually clean. Use plastic brushes to remove heavy or large debris.
 d) Soak the product for at least 10 minutes in an aqueous solution of a detergent at temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 e) Rinse the product under cold water for at least 2 minutes, paying particular attention to the holes and places difficult to be cleaned.
 f) Prepare fresh washing solution. Clean the surfaces and gaps of the product, carefully. Use suitable brushes to clean the holes. Clean the product immersed in the solution.
 g) Rinse the product thoroughly under warm running water for at least 2 minutes, paying special attention to the gaps, blind holes, hinges and joints. When cleaning, use brushes and perform multiple reciprocating movements on the surface of the product.
 h) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-h until the product is visually clean.
 i) Ultrasound cleaning: prepare an aqueous cleaning solution at a temperature of 40+/-2°C and pH of 10.4 - 10.8 (*follow the information contained in the instructions prepared by the manufacturer of the cleaning agent, in respect of temperature, concentration, exposure time and water quality*). Immerse fully the product in the aqueous cleaning solution and have it washed in ultrasounds for 15 minutes.
 j) Rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned.
 k) Visually inspect the entire surface of the product for debris and impurity. Repeat the steps described in subsections c-k until the product is visually clean.
 l) Use demineralized water for final rinsing of the device.
 m) Dry the device thoroughly using disposable, soft, lint-free cloth or compressed air.
 n) Prepare an aqueous solution of disinfecting agent at a temperature of 20+/-2°C using 20g of the agent per 1 liter of water. Immerse the product in the solution, exposure time – 15 min (*follow the information contained in the instructions prepared by the manufacturer of the agent, in respect of temperature, concentration, exposure time and water quality*).
 o) After the exposure time, rinse the product thoroughly under demineralized water, paying particular attention to the holes and places difficult to be cleaned.
 p) The cannulated instruments should be treated using a compressed air or air supplied from the syringe.
 q) Dry the device thoroughly. It is recommended to dry the product in a dryer at a temperature ranging from 90°C to 110°C, for 60 min.
 r) Visually inspect the entire surface of the device.
 s) CAUTION: If the obstruction in the cannula cannot be removed as indicated in the Instructions for Use, the device should be considered at the end of its useful life and should be discarded in accordance with facility procedures and guidelines.
6) The automated method using a washer - disinfector.

a) Equipment and materials: a washer - disinfector, aqueous solutions of cleaning agent.
b) Cleaning in the washer-disinfector must be preceded by a manual and ultrasound cleaning, following the procedure described in subsections c-h of paragraph 5.
c) CAUTION: The equipment used for washing/disinfection should meet the requirements of ISO 15883. Procedure of washing in the washer-disinfector shall be performed according to internal hospital procedures, recommendations of the washer-disinfector manufacturer, and instructions for use prepared by the washing-disinfecting agent manufacturer.
d) The device should undergo the process of machine washing in the washer-disinfector using the following cycle parameters: (*1*) - pre-washing in cold tap water, duration – 2 min; (*2*) - washing in an aqueous solution of cleaning agent at 55+/-2°C and pH of 10.4 - 10.8, duration – 10 min; (*3*) - rinsing under demineralized water, duration – 2 min; (*4*) - thermal disinfection in demineralised water at 90°C, minimal duration – 5 min; (*5*) - drying at the temperature ranging from 90°C to 110°C, duration - 40 min.

5. Inspection
1) Each time before re-use and re-sterilization, all medical devices should be inspected.
2) All parts of the product should be checked for visible dirt and corrosion. Particular attention should be paid to:
a) Places where dirt can be found, such as joints, latches, etc.
b) Holes, grooves and gaps the debris could have been pressed into during use.
3) Generally unmagnifed visual inspection under good light conditions is sufficient.
4) Each time before re-use and re-sterilization, the functional check of the product should be performed, consisting of:
a) Verifying the connections in the mating instruments, such as tips, shafts and quick coupling devices.
b) Verifying the correct functioning of mechanisms, e.g. screw, ratchet, snap mechanism, etc.
c) Verifying instruments for damage to material structure (*cracks, dents, peels, etc.*).
5) CAUTION:
a) The **ChM** does not define the maximum number of uses appropriate for re-usable medical instruments. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful and proper use reduces the risk of damage to the product and extends its serviceable life.
b) The manufacturer does not recommend using any preservatives on medical devices.
6) Damaged or defective product cannot be approved for further use.
7) When verifying the condition of the rod trials, care must be taken to ensure that the continuity of the silicone coating is maintained. In the case of damage (*e.g. cuts or defects expose the metal core*), the device should be considered at the end of its useful life and should be discarded in accordance with facility procedures and guidelines.
8) Prior to storage, the instrument must be checked for dryness.

6. Packaging
1) Washed and dried devices shall be stored (*if possible*) in suitable stands placed in special sterilization containers. Separate items should be packed in a packaging intended for the recommended steam sterilization. Sterilization containers, item packaging and packaging process itself have to meet the requirements of ISO 11607 standards. The packaging procedure must be performed in controlled purity conditions. The device must be packed so that during its removal from the packaging, when used, there is no risk for its re-contamination

7. Sterilization
1) The instrument is washed, disinfected and dried to undergo a sterilization process. The recommended sterilization procedure for non-sterile medical devices validated by **ChM** is autoclave vacuum steam sterilization (*with water vapor under overpressure*) with the following parameters:
a) temperature: 134°C,
b) minimum exposure time: 5 min,
c) minimum drying time: 20 min.
2) In order to minimize the potential risk of transmission of Creutzfeldt-Jakob disease (*especially for surgical instruments that could come into contact with the central nervous system*), national health authorities in some countries recommend sterilization according to the following parameters:
a) temperature: 134°C,
b) minimum exposure time: 18 min,
c) minimum drying time: 20 min.
3) CAUTION:
a) The sterilization process used in the medical facility must be validated and routinely controlled in accordance with the requirements of EN ISO 17665-1.
b) Sterilization must be effective and in accordance with requirements of the EN 556-1 standard to ensure the required level

of guaranteed sterility SAL 10^{-6} (*where SAL stands for Sterility Assurance Level*).

c) Device must not be sterilized in the packaging in which it was delivered, except specially designed sterilization containers.
d) The method of sterilization using ethylene oxide, gas plasma and dry heat should not be used, unless the Instructions for Use for the product contains sterilization recommendations using these methods.
e) The sterilization temperature for plastic products (*PPSU, PEEK, PTFE, silicone*) cannot be higher than 140°C.

6 STORAGE

1. The devices should be properly stored. When storing surgical instruments, it is recommended that they never be stacked together. It may lead to damage of cutting edges (*nick or dull*) and/or initiation of corrosion centers. Instruments should be stored in a clean and dry room, at room temperature and off the direct sunlight. If possible, instruments should be stored in suitable palettes placed into specially designed sterilization containers.

7 COMPATIBILITY

1. Trials are compatible with adequate implant systems of **ChM** production. It is not allowed to use products manufactured by **ChM** with products from other manufacturers. This is the physician who takes responsibility for the use of instruments from other manufacturers.

If these instructions appear unclear, please contact the manufacturer, who shall provide all required explanations.
Updated INSTRUCTIONS FOR USE are available at the following website: **ifu.chm.eu**

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3. ACCESSORIES 24

3.1. Tourniquets 24

3.2. Mobile stand for tourniquet control UNITS **[30.0040.000]** 24

4. CONSTRUCTION 25

4.1. Control panel 25

4.2. The display main window 27

5. USAGE 28

5.1. EZO-01 control unit activation 28

5.2. Use of the device when battery powered 28

5.3. Occlusion pressure 28

5.4. Occlusion time 28

5.5. Tourniquet inflation 29

5.6. Tourniquet deflating 29

5.7. Menu 29

5.7.1. TESTING 31

6. USE OF EZO-01 CONTROL UNIT DURING THE PROCEDURE 31

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1. INDICATIONS

EZO-01 electronic single tourniquet control unit is used to supply compressed air to tourniquets used to control the blood flow in the extremities in order to produce bloodless surgical field.

The built-in compressor and DC battery allow the functioning of the device without its connection to the compressed air and electrical AC 230V network. Thanks to that, EZO-01 can be used during patient's transportation, power supply AC 230V failure and when there is lack of compressed air network in the facility.

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2. SPECIFICATIONS AND EQUIPMENT

Power supply	100÷240V, 50/60Hz
Battery life	90min
Power consumption	65W
Work factor	Compressed air ***(compressed by built-in compressor)***
Output pressure ***(in tourniquet)***	100÷550mmHg with 5mmHg increments
Setting range of time	1÷99min. with 1min. increments
Dimensions ***(H/W/D)***	21cm x 21cm x 21.5cm
Weight	4kg
Ingress Protection Rating	IP20, EN 60529

The manufacturer reserves the right to structural changes.

Equipment supplied with the device.

	name	Pcs
	Power cord AC 230V, 1.8m	1
	Tourniquet control unit air hose - blue cat. No. 30.0035.010 ***(spiral version, max. length ≈1.7m)*** or 30.0035.1xx * ***(straight version)***	1

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3. ACCESSORIES

EZO-01 unit should only be used with original accessories provided by ChM sp. z o.o.

3.1. TOURNIQUETS

No	Product name	The range of use *(limb circumference in cm)*	Dimensions	Cat. No
1	Single arm tourniquet	25÷40	64x13	30.0009
2	Single tourniquet for children	14÷20	50x6	30.0010
3	Single femoral tourniquet	38÷58	85x14	30.0012
4	Single femoral tourniquet	38÷58	120x13.5	30.0013
5	Single femoral tourniquet	38÷58	140x13.5	30.0008
6	Conical single femoral tourniquet	40÷60	110x11	30.0014
7	Single arm tourniquet long	38÷58	82x8	30.0015
8	Single arm tourniquet	25÷40	62x7	30.0016
9	Single tourniquet for babies	10÷17	30x3	30.0017

3.2. MOBILE STAND FOR TOURNIQUET CONTROL UNITS [30.0040.000]

The stand is used as a mobile platform for the tourniquet control unit, as well as the accessories needed to be used with the device, such as tourniquets, skin protection padding, power cords and other accessories.


Technical parameters	
Dimensions *(H/W/D)*	900 x 570 x 570 mm
The height with the control unit	1100 mm
Weight	9 kg

Tourniquets and stand are additional accessories of the control unit and are offered as a separate accessory. Prior to use, consult Instructions for Use supplied with these products.

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
4. CONSTRUCTION

1. Control panel
2. Handle
3. USB port *(for use only by service personnel)*
4. ON/OFF AC power switch
5. AC power connection socket
6. Fuse drawer
7. Holes for mounting the control unit on the mobile stand 30.0040.000

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
	Control panel elements	
1		ENTER It approves changes introduced to: time, pressure, and in the menu options.
2		PRESSURE Entry into the pressure settings.
3		TIMER Entry into the time settings.
4		CONTROL KNOB Turn of the control knob: -changes the time and pressure ***(after entering the pressure or time settings)***. -navigates the menu options ***(after entering the MENU window)***. Pressing the control knob equals pressing the ENTER button.
5		STAND BY It turns on/off the EZO-01.
6		MENU Entry into the EZO-01 settings window.
7		INFLATE It activates the process of filling up the tourniquet and the work of EZO-01 control unit.
8		DEFLATE It activates the process of emptying the tourniquet and finishes the work of EZO-01 control unit.
9		TOURNIQUET CONNECTION PORT Place where the **ChM** tourniquet is connected to the EZO-01 control unit.

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4.2. THE DISPLAY MAIN WINDOW

The main window is displayed when the control unit is turned on and when the control unit cooperates with the tourniquet.

1. The value of the pressure in the tourniquet
2. The value of the time of occlusion
3. Counting down ↓ / Counting up ↑
4. The control unit in edit mode of tourniquet pressure
5. The control unit in edit mode of occlusion time
6. Mains power / Battery power *(with information on the battery charge level)*

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
5. USAGE

5.1. EZO-01 CONTROL UNIT ACTIVATION

Use AC power cord ***(supplied with the device)*** to connect the EZO-01 to 230V/50Hz network with a protective grounding. Turn the switch located on the back panel of the device ON. Press the STAND BY button located on the control panel. The internal auto-test will be activated. Afterwards, the main window will be displayed.

If the internal auto-test result is negative, an error message will appear on the screen.

Further work with the device will not be possible. Contact the device manufacturer.

Connect the tourniquet ***(using a spiral hose supplied with the** control unit)* to the slot located on the control panel.

5.2. USE OF THE DEVICE WHEN BATTERY POWERED

In the case of: power failure, when the control unit is not plugged in or the switch on the back panel is OFF, the EZO-01 automatically switches into battery power mode. The battery power icon ,will be displayed as well as information on the battery charge level ***(tab. 1)***.

Operating time of the EZO-01 on battery power is dependent on many factors. When the battery charge level goes below 30%, initiate battery charging by connecting the control unit to AC 230V power network ***(the display shows the mains power icon)***.

Icon	Battery charge level
	100÷70%
	70÷50%
	50÷30%
	30÷10%
	10÷0%

Tab.1. Battery charge level

5.3. OCCLUSION PRESSURE

Push the PRESSURE button. The display will show the value of the required pressure and the icon indicating that the unit is in pressure edit mode. By turning control knob, set the desired pressure ***(in the range from 100 to 550 mmHg)***. Confirm the selected value pressing ENTER button or control knob. If the pressure selection is performed before inflating the tourniquet, once approved, the displayed value returns to the value of 000 ***(the actual pressure in the tourniquet)***. If the pressure is edited after inflating the tourniquet, pressure approval automatically changes the pressure in the tourniquet to the selected one.

5.4. OCCLUSION TIME

Push the TIMER button. The display will show the value of the required time and the icon indicating that the unit is in time edit mode. By turning control knob, set the desired time [in the range from 1 to 99 min]. Confirm the selected value pressing ENTER button or control knob. Time edit mode is completed. The time can be changed both, before and after inflating the tourniquet.

The clock is activated automatically when the tourniquet is being inflated. When the selected time has been reached, a beep sound will be emitted. If the tourniquet is not deflated or new time is not set, the timer will start to count up from the 00.

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The tourniquet does not get deflated when the selected occlusion time is reached.

The occlusion time ***(after inflating the tourniquet)*** may be displayed as:

- the increasing value from zero to the set value,the arrow facing up will be displayed ↑
- the decreasing value from the set value to zero, the arrow facing down will be displayed ↓

One of the abovementioned options is selected acc. to p.5.7 of these instructions.

5.5. TOURNIQUET INFLATION

Press and hold pressed for minimum 2s INFLATE button. The control unit will start to fill in the tourniquet, what will be shown on the display as an increasing pressure level up to the selected one. Automatically, the occlusion time measurement will be activated.

Inflating loose, not applied on a limb or rolled into a tight roll tourniquet can cause its damage.

5.6. TOURNIQUET DEFLATING

Press and hold pressed for minimum 2s DEFLATE button. The control unit will start to empty the tourniquet, what will be shown on the display as an decreasing pressure level down to 000 value. Automatically, the occlusion time measurement will be stopped.

5.7. MENU

Changes in other parameters of the EZO-01 unit are introduced from the settings window when pressing the MENU button. For navigation in the settings window and changes introduction, use the control knob. Press MENU button to leave the ***(Options)*** settings window.

Initiation of settings window is not possible during the operation of the control unit with the tourniquet *(after its inflating)*.







Fig. 2. The menu window

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ICON	DESCRIPTION	WINDOW VIEW
	Timer settings. Two options: • Count down *(from the specified value down to zero)*; • Count up *(from zero up to specified value)*	
	Settings of the system clock. Two options: • Show *(show the date and time settings)* • Set clock *(set the date and time)*	
	Language options. Three possible options: • Polish • English • Russian	
	Display settings. Six levels of brightness.	
	Speaker. Key Tone. Two options: • ON • OFF No option to mute the alarm.	
	Information concerning general EZO-01 parameters: • Last pressure value • Last time value • Speaker *(on/off)* • Battery level	
	Autotest. Defining the efficiency of the unit control - tourniquet system. See p. 5.7.1 of these instructions.	
	Data Transmission Mode. Repairs can only be performed by qualified service personnel.	

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5.7.1. TESTING

For testing, use tourniquet that will be used during the procedure. When replacing the tourniquet, the testing should be repeated.

Connect the **ChM** tourniquet using the spiral air hose supplied with the EZO-01 control unit to the tourniquet connection port located on the control panel. Roll up the tourniquet tightly and secure against unrolling. Then enter the Menu and select Autotest. During the testing, the tourniquet will be inflated and the device will check the tightness of the system: control unit – tourniquet. If the system is tight, the testing will end with a PASS sign - the control unit and tourniquet can be used during the procedure. Displayed FAILD sign indicates that the system is not tight. Check the condition of the tourniquet, air hoses, pneumatic connectors and their connection. Eliminate leaks if possible ***(e.g. replace defective tourniquet or air hose)***. Run the test again.

Should it be impossible to remove or detect the leaks, it is strictly forbidden to use the set: control unit-tourniquet for surgery. Only the set with positive test results *(PASS displayed)* can be used for surgery

6. USE OF EZO-01 CONTROL UNIT DURING THE PROCEDURE

Prior to use, read p. 5 USAGE of these instructions.

Before using the EZO-01 control unit in surgical procedure for the first time, the device should be cleaned in accordance with section 9 of this manual.

The EZO-01 control unit can only be used by qualified medical personnel familiar with the operation of the device, such as anesthetists, who control the device outside the operating field.

Prior to each use:

- perform a visual inspection of the product,
- verify the basic functioning of the unit control ***(inflating, deflating, etc.)***,
- verify the battery charge ***(use control unit only with fully charged battery)***,
- perform the testing in accordance with section 5.7.1 of these instructions,
- prepare a replacement unit.

Should any defect appear, do not use the device until the defect is removed or until the defect is deemed to have no negative effect on the proper functioning of the device.

Connect the device to the AC 230V network during the procedure. The device functioning on the battery power during the procedure is allowed only in the case of failure of AC power ***(for more information about working on battery power - see p.5.2 of these instructions)***. Set the occlusion pressure and time. Apply the tourniquet on the operated limb. Make sure that the pneumatic hose has not been bent. Inflate the tourniquet by pressing INFLATE button. If it is necessary to regulate the occlusion pressure or time, these changes can be introduced

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after inflating the tourniquet.
When the occlusion time is reached, a beep sound will be emitted. If the tourniquet is not deflated or new time is not set, the timer will start to count from the 00.

The tourniquet does not get deflated when the selected occlusion time is reached.

Press and hold for minimum 2 sec. DEFLATE button to empty the tourniquet. After emptying the tourniquet ***(the display shows the 000 pressure value)***, remove the tourniquet from the limb, then disconnect it from the unit control. The EZO-01 may be switched off by pressing STAND BY button.

The surgeon decides about the: setting values of the control unit *(pressure, occlusion time, etc.)*, the place the tourniquet is applied, the moment the tourniquet is inflated, deflated and other parameters.

7. GRAPHIC AND SOUND MESSAGES

EZO-01 single tourniquet control unit is equipped with a range of graphic and sound messages that inform the user about relevant parameters or irregularities during the operation. The table ***(Tab. 3.)*** lists the causes of the messages, types of messages and the actions to be taken.

Cause	Message	Action
Low battery level	- graphic information - sound signal	- begin the process of charging the battery, connect the control unit to AC 230V 50Hz network, switch on the rear panel must be in ON position
Leakage in control unit - tourniquet system.	- graphic information - sound signal	- verify the tightness of cords, tourniquets and all connections between the tourniquet and the control unit. Remove the leakage. If the leakage cannot be removed, depending on its place, replace the tourniquet or the control unit.
Occlusion time limit reached	- sound signal	- action depends on a surgeon's decision ***(deflating the tourniquet, continuation of occlusion)***

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8. ELECTROMAGNETIC COMPATIBILITY

8.1. ELECTROMAGNETIC EMISSIONS

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions, CISPR 11	Group 1, Class B	The device uses RF energy only for its internal function. Therefore, its RF emission is very low and it is not likely to cause any interference in nearby electronic equipment.
Harmonic emissions IEC 61000-3-2	Class A	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations / flicker emissions, IEC 61000-3-3	Comply	

8.2. ELECTROMAGNETIC IMMUNITY

Immunity test standard	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge ***(ESD)*** IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be made from wood, concrete or ceramic tile. If floors are covered with a synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2kV for power supply lines ±1kV for in/ out lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV line to line ±2kV line to earth	±1kV line to line ±2kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% Ut ***(0.5 cycle)*** 40% Ut ***(5 cycles)*** 70% Ut ***(25 cycles)*** <5% Ut ***(5 seconds)***	<5% Ut ***(0.5 cycle)*** 40% Ut ***(5 cycles)*** 70% Ut ***(25 cycles)*** <5% Ut ***(5 seconds)***	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires its continuous use even during power outages, it is recommended that the device be connected to an emergency power supply. Ut - mains supply voltage
Power frequency ***(50/60 Hz)*** magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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Immunity test standard	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 800 MHz	V1 = 3V/m	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.17 \sqrt{P}$ $d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts ***(W)*** according to the transmitter manufacturer and d is the recommended separation distance in metres ***(m)***. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a), should be less than the compliance level in each frequency range b). Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 kHz to 2500 MHz	E1 = 3V/m	

a) Field strengths from fixed transmitters, such as base stations for radio telephones ***(cellular/ cordless)*** and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

8.3. SEPARATION DISTANCES

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres ***(m)*** can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts ***(W)*** according to the transmitter manufacturer.

Rated maximum output power of transmitter *(W)*	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P}$[m]	**80 MHz to 800 MHz** $d = 1.17 \sqrt{P}$[m]	**800 MHz to 2.5 GHz** $d = 2.33 \sqrt{P}$[m]
0.01	0.12	0.04	0.1
0.1	0.37	0.11	0.32
1	1.2	0.35	1.0
10	3.7	1.1	3.2
100	12	3.5	10

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The outer surfaces of the device should be cleaned with a soft, moistened ***(not dripping)*** cloth with a mild detergent solution then wiped dry. In the case of contaminating the housing of the device with infectious agents ***(blood, body fluids, etc.)***, disinfection using disinfectants with a neutral pH is recommended. In such cases, soak a soft cloth in disinfectant, wipe the housing and allow it to dry.

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The EZO-01 control unit should NOT be sterilized

10. STORAGE

The device should be stored in a clean and dry room, off the direct sun, at the temperature of 15 ÷ 30°C and humidity of 20 to 70%.

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The non-warranty services are to be handled also by the manufacturer. The manufacturer shall not be liable for repairs ***(and resultant problems)*** carried out by entities other than **ChM** or appointed service points.

To maintain trouble-free and safe operation throughout the life of the device, it is recommended to perform annual technical inspections at the manufacturer's site.

The lifetime of the EZO-01 control unit is 10 years. After this period, the manufacturer *(ChM)* may refuse to perform service repairs or technical inspections.

For servicing, warranty and post-warranty repairing, please contact the manufacturer's representative or **ChM** directly

ChM sp. z o.o.
Lewickie 3b
16-061 Juchnowiec Kościelny, POLAND
tel: +48 85 86 86 100
fax: +48 85 86 86 101
e-mail: chm@chm.eu
www.chm.eu

12. THE LIFETIME AND DISPOSAL

The **ChM** sp. z o.o. does not define the maximum number of uses for the EZO-01 control unit. The useful life of the device depends on many factors including the method and duration of each use, applied cleaning and disinfection and the handling between uses. Nontheless, the lifetime of the EZO-01 control unit is 10 years. After this period, the manufacturer ***(ChM)*** may refuse to perform service repairs or technical inspections.

The device should be disposed of in accordance with the applicable electrical and electronic equipment law. It is recommended that the battery be replaced at least once a year.

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Warnings


Fig. 4. Warnings marked on the casing.

SYMBOL	CONCEPT, MEANING	EXPLANATION
	It is forbidden to dispose of the product into unsorted municipal waste.	The waste devices and other electronic and electrical products should be collected separately and disposed of in accordance with the applicable law on waste electrical and electronic equipment.
	Caution	Read the provided Instructions for Use before operating the device.
	BF Type of Application	The degree of protection against electric shock.
CE 0197	European conformity mark with the number of the notified body	Certifies compliance of the unit structure with the guidelines of the European Community
	Instructions For Use	Before the first use, refer to the Instructions for Use.

Tab. 4. Explanation of the symbols on the date plate

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14. TROUBLESHOOTING

The following table lists some potential problems that may occur while using the device. For each problem, the most likely causes are provided.

If any other problems than the ones listed below have occurred or if the suggested actions have not solved the problem, please contact the manufacturer.
DO NOT perform repairs on your own unless otherwise stated in this manual.

Problem	Possible causes	Solution
Accidental or unintended operations of the device	The device driver wet due to improper clearing	Leave to dry
		Send to **ChM** service
	System error	Send to **ChM** service
	System failure	Send to **ChM** service
	Incomplete information provided in the safety instructions	Report to the manufacturer
	Overheating of the device	Send to **ChM** service
	Strong electromagnetic interference	Turn off all electrical equipment not in use; Change the location of electrical equipment; Plug in electrical equipment to other available outlets
	Exceeding the declared lifetime of the device *(10 years)*	Send to **ChM** service for checking the further usefulness of the product
No response of the device to the set parameters	System failure	Send to **ChM** service
	System hand up	Turn the device off and then back on
	Overheating of the device	Send to **ChM** service
Tourniquet does not inflate	The air hose may be bent or incorrectly connected	Check the hose and the connection
	Damaged tourniquet	Replace the tourniquet for a new one
	Damaged pneumatic air hose	Send to **ChM** service
		Replace the air hose
	System failure	Send to **ChM** service
Tourniquet does not deflate	The air hose may be bent or incorrectly connected	Check the hose and the connection
	Damaged tourniquet	Replace the tourniquet for a new one
	Damaged pneumatic circuit	Send to **ChM** service
		Replace the air hose
	System failure	Send to **ChM** service
Dropping of the device	System failure	Send to **ChM** service
	Display damage	Send to **ChM** service
	Battery damage	Send to **ChM** service
	Appearance of dangerous voltage on the device housing	Send to **ChM** service
	Damage to the power cord	Replace the power cord for the new one
	Damage to the air hose	Replace the air hose for the new one

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Problem	Possible causes	Solution
Unit does not turn ON.	Unit not plugged in and battery fully depleted	Plug in the unit
	The unit is plugged in but the switch is OFF and battery is fully depleted	Turn the switch ON
	Damaged power cord	Replace the power cord for the new one
	Battery failure	Send to **ChM** service
	System failure	Send to **ChM** service
The device does not maintain the set pressure	The air hose may be bent or incorrectly connected	Check the hose and the connection
	Damaged tourniquet	Replace the tourniquet for a new one
	Damaged air hose	Replace the air hose for the new one
		Send to **ChM** service
	System error	Send to **ChM** service
	Exceeding the declared lifetime of the device ***(10 years)***	Send to **ChM** service for checking the further usefulness of the product
Tourniquet cannot be connected to the unit	Damaged air hose	Replace the air hose for the new one
		Send to **ChM** service
	Using tourniquets other than **ChM**'s	Use **ChM** tourniquets
Discoloration and scratches on the surface of the device	Use of detergents other than those recommended by the manufacturer	Consult the Instructions for Use of the unit
	Overheating of the device	Send to **ChM** service
	Exceeding the declared lifetime of the device ***(10 years)***	Send to **ChM** service for checking the further usefulness of the product
Unit switches off	Battery fully depleted	Charge the battery
		Replace the battery for the new one in the manufacture's service
	System error	Send to **ChM** service
	System failure	Send to **ChM** service
	Strong electromagnetic interference	Turn off all electrical equipment not in use; Change the location of electrical equipment; Plug in electrical equipment to other available outlets
Too short battery life	End of battery life	Send to **ChM** service
		Consult the Instructions for Use of the unit
	Failure to follow manufacturer's instructions	Consult the Instructions for Use of the unit
	Exceeding the declared lifetime of the device ***(10 years)***	Send to **ChM** service for checking the further usefulness of the product
Information cannot be read from the device display	Display failure	Send to **ChM** service
	System failure	Send to **ChM** service
	System error	Send to **ChM** service
Unit cannot be set to STANDBY	Pressure sensed in tourniquet	Deflate the tourniquet
	System failure	Send to **ChM** service
Alarm indicating a leaking system	The air hose may be bent or incorrectly connected	Check the hose and the connection
	Damaged tourniquet	Replace the tourniquet for a new one
	Damaged pneumatic air hose	Send to **ChM** service
		Replace the air hose for a new one
	System failure	Send to **ChM** service

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Problem	Possible causes	Solution
Alarm indicating discharged battery	Battery fully depleted	Charge the battery
	Battery failure	Send to **ChM** service
	System failure	Send to **ChM** service
No external power supply icon when the device is plugged in	The device is not plugged in	Plug in the device
	Damaged power cord	Replace the power cord for the new one
	The switch is OFF	Switch ON
	No mains voltage	Check the voltage In the mains
	System failure	Send to **ChM** service
The battery does not charge	The device is not plugged in	Plug in the device
	The unit is plugged in but the switch is OFF and battery is fully depleted	Turn the switch ON
	Battery failure	Send to **ChM** service
	System failure	Send to **ChM** service
Possible electrical shock	Damage to the power cord	Send to **ChM** service
Electrical interference	Occurrence of electrical noise	Turn off all electrical equipment not in use
		Change the location of electrical equipment
		Plug in electrical equipment to other available outlets
The unit has been sterilized	Staff negligence	Send to **ChM** service
	Negligence of Instructions for Use supplied with the product	

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EN SINGLE ADJUSTABLE TOURNIQUET CONTROL UNIT

1. TECHNICAL PARAMETERS 11
2. INDICATIONS FOR USE 12
3. CONNECTION TO THE AIR RESERVOIR OR COMPRESSED AIR MAIN OR NITROGEN 12
4. FUNCTIONING 12
5. INDICATIONS AND OPERATION OF THE TIMER 13
6. ACCESSORIES 14
7. CLEANING AND DISINFECTION 16
8. SERVICE 16

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
Fig.1. Single adjustable tourniquet control unit

1 – indicator of feed pressure presence
2 – manometer of pressure measurements in tourniquet *(mmHg)*
3 – inflate/ deflate switch
4 – knob of pressure regulation in tourniquet
5 – spiral hose to connect the tourniquet
6 – timer
7 – compressed air or nitrogen terminal
8 – battery LR6 1.5V slots

1. TECHNICAL PARAMETERS

Work factor	Compressed air or nitrogen from the air reservoir
Feed pressure	0.2÷1 MPa
Indicator of feed pressure presence	Installed
Initial pressure	0÷500 mmHg
Time measurement	Timer with setting at every 1 second
Sound signaling of time passing	Installed
Electric feeding	2 batteries LR6-1.5V
Overall dimensions (width/depth/height)	150/210/180
Weight	2.5 kg

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2. INDICATIONS FOR USE

Single adjustable tourniquet control unit is used to provide compressed air or nitrogen of precisely controlled pressure to a single tourniquet placed on the patient's limb during the operation in a bloodless field. It enables continued control and regulation of the pressure in the single tourniquet. Due to this, the possibility of the clamp degree regulation is obtained.
Installed timer has sound signal activated when the set time of clamp or release passes.
Single adjustable tourniquet control unit is also applicable for exact control of clamp when collecting blood.

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- set inflate/ deflate switch *(3)* to the deflate position,
- using armored wire provided with the device, connect the external network or a compressed air or nitrogen reservoir with the terminal *(7)* of the single adjustable tourniquet control unit *(Fig. 1)*.

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4. FUNCTIONING

- set the inflate/ deflate switch *(3)* to the deflate position,
- connect single tourniquet to the hose *(5)*,
- roll up the tourniquet in a tight roll or wrap it around stiff roller *(it protects the **tourniquet** from damage during the inflating process)*,
- set the inflate/ deflate switch *(3)* to the inflate position,
- observing the manometer *(2)*, set desired pressure using the knob *(4)*,
- set the inflate/ deflate switch *(3)* to the deflate position,
- place the tourniquet on the patient's limb,
- set desired time of the pressure using the timer buttons *(6)*,
- using the switch *(3)*, inflate the single tourniquet and press START/STOP button to activate the timer *(6)*,
- use switch *(3)* to deflate and inflate the tourniquet when the set time of the pressure passes or when the timer *(6)* stops.

WARNING:
Sound signal or timer stop *(6)* does not signify automatic inflating or deflating of the tourniquet. These procedures shall be activated manually using the switch *(3)*.

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1. As a stopwatch with countdown of a given time value in the scope from 99min. 59sec. down to 0, with acoustic signalization of established time passing.
2. As a stopwatch with counting the time up in the scope between 0 to 99min. 59sec.

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

2. Trolley for tourniquet control units and cylinders.

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
7. CLEANING AND DISINFECTION

The outer surfaces of the device should be cleaned with a soft, moistened *(not dripping)* cloth with a mild detergent solution then wiped dry. In the case of contaminating the housing of the device with infectious agents *(blood, body fluids, etc.)*, disinfection using disinfectants with a neutral pH is recommended. In such cases, soak a soft cloth in disinfectant, wipe the housing and allow it to dry.
Accessories used with the device such as mobile stand for tourniquet control units **[30.0027.000]** and tourniquet control unit holder **[30.0023.000]** should be cleaned and disinfected in accordance with the above-mentioned recommendations.

NOTE:
Do not allow the detergents or disinfectants to enter inside the device housing. Should the detergents or disinfectants enter inside the device housing, do not use the single adjustable tourniquet control unit and send it to the manufacturer for cleaning processes.

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Changes introduced by the user and not accepted by **ChM** may be the basis for warranty annulment.
The warranty does not include product damage due to: negligence, inappropriate usage, usage not in accordance with the intended purpose, mechanical damage.

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It is recommended to perform technical inspection of the device at the manufacturer site at least once every 12 months.

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1. TECHNICAL PARAMETERS 12

2. INDICATIONS FOR USE 13

3. CONNECTION TO THE AIR RESERVOIR OR COMPRESSED AIR MAIN OR NITROGEN 13

4. FUNCTIONING 13

5. INDICATIONS AND OPERATION OF THE TIMER 14

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
Fig.1. Dual adjustable tourniquet control unit

1 - indicator of feed pressure presence
2 - manometers of pressure measurements in tourniquets *(mmHg)*
3 - inflate/ deflate switches
4 - knobs of pressure regulation in tourniquets
5 - spiral hoses to connect the tourniquets
6 - timer
7 - compressed air or nitrogen terminal
8 - Battery LR6 1.5V slots

1. TECHNICAL PARAMETERS

Work factor	Compressed air or nitrogen from the air reservoir
Feed pressure	0.2÷1 MPa
Indicator of feed pressure presence	Installed
Initial pressure	0÷500 mmHg from two separate terminals each smoothly regulated
Time measurement	Timer with setting at every 1 second
Sound signaling of time passing	Installed
Electric feeding	2 batteries LR6-1.5V
Overall dimensions *(width/depth/height)*	200/210/180
Weight	4 kg

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2. INDICATIONS FOR USE

Dual adjustable tourniquet control unit is used to provide compressed air or nitrogen of precisely controlled pressure to a dual tourniquet or two single tourniquets placed on the patient's limb during the operation in a bloodless field. It enables continued and independent control and regulation of the pressure in the dual tourniquet or two single tourniquets. Due to this, the possibility of the clamp degree regulation is obtained. Installed timer has sound signal activated when the set time of clamp or release passes.
Dual adjustable tourniquet control unit is also applicable for exact control of clamp when collecting blood.

3. CONNECTION TO THE AIR RESERVOIR OR COMPRESSED AIR MAIN OR NITROGEN

- set inflate/ deflate switch *(3)* to the deflate position,
- using armored wire provided with the device, connect the external network or a compressed air or nitrogen reservoir with the terminal *(7)* of the dual adjustable tourniquet control unit *(Fig. 1)*.

WARNING:
The output pressure of the external network or a compressed air or nitrogen reservoir must not exceed 1 MPa and must not be lower than 0.2MPa.

4. FUNCTIONING

- set the inflate/ deflate switches *(3)* to the deflate position,
- connect dual tourniquet or two single tourniquets to the hoses *(5)*. If dual tourniquet has marked one section, this mark has to be attributed to one of the hoses *(the red or blue one)*. The colors on the frontal panel are equal to colors of the hoses, which enables the identification of working tourniquet section,
- roll up the tourniquets in a tight roll or wrap them around stiff roller *(it protects the tourniquets from damage during the inflating process)*,
- set the inflate/ deflate switches *(3)* to the inflate position,
- observing the manometers *(2)*, set desired pressure using the knobs *(4)* in particular section of dual tourniquet or in individual single tourniquets,
- set the inflate/ deflate switches *(3)* to the deflate position,
- place the tourniquet*(s)* on the patient's limb,
- set desired time of the pressure using the timer buttons *(6)*,
- using the switches *(3)*, inflate the suitable section of the dual tourniquet or desired single tourniquet and press START/STOP button to activate the timer *(6)*,
- use switches *(3)* to deflate and inflate particular sections of dual tourniquet or individual single tourniquets when the set time of the pressure passes or when the timer *(6)* stops. During such procedures it has to be remembered which of the dual tourniquet part, or which single tourniquet matches red titles and which blue titles on the front panel of the device.

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WARNING:
Sound signal or timer stop *(6)* does not signify automatic inflating or deflating of the tourniquet. These procedures shall be activated manually using the switches *(3)*.

5. INDICATIONS AND OPERATION OF THE TIMER

Installed timer enables work in two modes:

1. As a stopwatch with countdown of a given time value in the scope from 99min. 59sec. down to 0, with acoustic signalization of established time passing.

2. As a stopwatch with counting the time up in the scope between 0 to 99min. 59sec.

Operation in work mode 1

- Pressing the MIN *(1)* and SEC *(2)* button, set on the display the time that shall be measured.
- Pressing the START/STOP *(3)* button, activate counting down.

Fig.2. Timer

WARNING:
Stop the counting by pressing the START/STOP *(3)* button. Press again the button *(3)* to continue the counting down.

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- When the set time is reached, the acoustic signal is activated for 15sec. The earlier activation of the signal is possible also by pressing the START/STOP *(3)* button which causes also the former set time to present on the display.
- Timer zeroing is done by simultaneous pressing the MIN *(1)* and SEC *(2)* buttons *(in any moment)*.

Operation in work mode 2

- Press START/STOP *(3)* button to activate time counting; passing time is presented on the display.
- Press START/STOP *(3)* button to stop the timer.
- Timer zeroing is done by simultaneous pressing the MIN *(1)* and SEC *(2)* buttons *(in any moment)*.


6. ACCESSORIES

1. Mobile stand for tourniquet control units.

- the stand **[30.0027.000]** serves as a mobile platform for the tourniquet control unit. Additionally, it is equipped with a basket for tourniquets and pneumatic hoses.

NOTE:
Mobile stand for tourniquet control units is additional equipment of the control unit and is offered as a separate device.

Installation of the control unit on the stand:

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
2. Trolley for tourniquet control units and cylinders:

- the trolley **[30.0024.000]** serves as a mobile platform for the tourniquet control unit and compressed air cylinder.

NOTE:
Trolley for tourniquet control units and cylinders is additional equipment of the control unit and is offered as a separate device.


Installation of the control unit on the trolley:


3. Tourniquet control unit holder:


- Holder **[30.0023.000]** is used to mount the control unit on the instrumental rails of Eurostandard type.

Installation of the control unit on the rail:

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Accessories used with the device such as mobile stand for tourniquet control units **[30.0027.000]** and tourniquet control unit holder **[30.0023.000]** should be cleaned and disinfected in accordance with the above-mentioned recommendations.

NOTE:
Do not allow the detergents or disinfectants to enter inside the device housing. Should the detergents or disinfectants enter inside the device housing, do not use the dual adjustable tourniquet control unit and send it to the manufacturer for cleaning processes.

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8. SERVICE

When servicing and repairing is concerned, please contact the manufacturer's representative or contact **ChM** directly:

ChM sp. z o.o.
Lewickie 3b
16-061 Juchnowiec Kościelny, POLAND
tel: +48 85 86 86 100
fax: +48 85 86 86 101
e-mail: chm@chm.eu

Changes introduced by the user and not accepted by **ChM** may be the basis for warranty annulment. The warranty does not include product damage due to: negligence, inappropriate usage, usage not in accordance with the intended purpose, mechanical damage.

Accessories and spare parts: filter, initial pressure regulator, electric-pneumatic transmitter, precise pressure regulator, safety valve, release valve, manometer, timer, spiral wire, wire to connect compressed air.

Warranty and non-warranty services and inspections are to be handled solely by experienced personnel with appropriate qualifications. The user is prohibited to replace any of the constituent parts of the device *(except for R6 or LR6 - Fig. 1)*.

The abovementioned services and inspections are performed solely by the device manufacturer: **ChM** sp. z o.o.

It is recommended to perform technical inspection of the device at the manufacturer site at least once every 12 months .

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- 4.1. CONTROL PANEL 6

5. OPERATION 7

- 5.1. CONNECTION TO THE NETWORK OF MEDICAL GASES 7
- 5.2. TOURNIQUET CONNECTION 8
- 5.3. TOURNIQUET INFLATION AND DEFLATION. PRESSURE REGULATION IN TOURNIQUET 8
- 5.4. SETTING TOURNIQUET OCCLUSION TIME 8

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7. INDICATIONS FOR USE

Single adjustable tourniquet control unit is intended to provide compressed air to tourniquets used to control blood flow in extremities in order to produce a bloodless surgical field.

8. TECHNICAL PARAMETERS AND EQUIPMENT

Work factor	Medical compressed air
Feed pressure	400kPa
Output pressure (in tourniquet)	0÷500mmHg
Time measurement	1s ÷ 99min 59s with 1s increments
Dimensions (W/D/H)	20.5cm x 21cm x 16cm
Weight	2.2kg

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	Name	Pcs.
	Tourniquet control unit air hose - blue cat. no. 30.0035.010 *(spiral version, max. length ≈1.7m)* or 30.0035.1xx * *(straight version)*	1
	Compressed air spiral hose cat. no. 30.0037.xxx	1

The hose length and type of quick coupling for connection to the medical gases outlet point acc. to the customer's specification *(table 1)*

Table 1.

Quick coupling connector type *(acc. to ISO 9170-1)*	**Hose cat. no.**
SS 875 24 30 *(AGA Swedish Standard)*	30.0037.1xx*
DIN 13260-2 *(German Standard)*	30.0037.2xx*
ENV 737-6 *(European Standard)*	30.0037.3xx*
UNI 9507 *(Italian Standard)*	30.0037.4xx*
BS 5682-1998 *(British Standard)*	30.0037.5xx*
AS 2896 *(Australian Standard)*	30.0037.6xx*
AFNOR NF S 90-116 *(French Standard)*	30.0037.7xx*

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Tourniquets are additional equipment of the control unit and are offered as separate devices.

9.2. MOBILE STAND FOR TOURNIQUET CONTROL UNITS [30.0040.000]

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
Mobile stand for tourniquet control units is additional equipment of the control unit and is offered as a separate device.

Technical parameters	
Dimensions *(H/W/D)*	900 x 570 x 570 mm
The height with the control unit	1100 mm
Weight	9 kg

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

10.1. CONTROL PANEL

4 - Knob of pressure regulation **5** - Spiral hose terminal to connect the tourniquet **6** - Manometer of pressure measurements in tourniquet **7** - Timer

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
11. OPERATION

11.1. CONNECTION TO THE NETWORK OF MEDICAL GASES

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Place the device on a flat surface that will provide a good grip for the rubber feet, so the uncontrolled movements of the unit are eliminated. In case of working with the unit attached to the mobile stand **[30.0040.000]**, use the brakes located at the wheels to immobilize the stand. Before connecting the tourniquet control unit, turn knob of pressure regulation *(4)* fully counter-clockwise. Use quick coupling connector *(**A**)* to connect the compressed air spiral hose **[30.0037]** with the unit using a medical compressed air terminal located on the rear panel of the device. Connect quick coupling connector *(**B**)* with a network of medical gases – compressed air of 400kPa - the operating rooms are equipped with.

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11.2. TOURNIQUET CONNECTION

Use tourniquet control unit air hose – blue **[30.0035.010]** to connect the tourniquet control unit with tourniquet. Use quick coupling connectors **C** to connect air hose with spiral hose terminal *(5)* located on the unit front panel. Connect threaded quick coupling connector D of air hoses with **ChM** tourniquet.

11.3. TOURNIQUET INFLATION AND DEFLATION. PRESSURE REGULATION IN TOURNIQUET

Use knob of pressure regulation *(4)* to inflate, deflate and regulate the pressure in the tourniquet.

NOTE:
Knob is equipped with protection against accidental changing in the pressure. In order to change the pressure in the tourniquet, pull the knob and then by turning it clockwise or counter-clockwise, change the pressure..

In order to inflate and increase the pressure in the tourniquet, turn the knob *(4)* clockwise and observe the monometer *(6)* where the increase in pressure will be indicated. End the regulation as soon as the desired pressure is achieved.

NOTE:
The maximum pressure is 500mmHg

In order to decrease the pressure in tourniquet, turn the knob *(4)* counter-clockwise. The decrease in pressure will be indicated on the monometer *(6)*. End the regulation as soon as the desired pressure is achieved. In order to empty the tourniquet, turn the knob counter-clockwise fully. The monometer *(6)* should indicate 0mmHg.

11.4. SETTING TOURNIQUET OCCLUSION TIME

The tourniquet control unit has been equipped with a timer *(7)* which allows for controlling the tourniquet occlusion time. In order to set the time of occlusion, press MIN button to increase the number of minutes. Press SEC button to increase the number of seconds. Press START/STOP button to start the countdown. When the set time of occlusion has passed, sound signal is activated.

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
14. THE LIFETIME AND DISPOSAL

The **ChM** sp. z o.o. does not define the maximum number of uses for the EZO-01 control unit. The useful life of the device depends on many factors including the method and duration of each use, applied cleaning and disinfection and the handling between uses. Nontheless, the lifetime of the EZO-01 control unit is 10 years. After this period, the manufacturer (**ChM**) may refuse to perform service repairs or technical inspections.

15. LABELS AND WARNINGS


Nameplate

Explanation of symbols on the nameplate

SYMBOL	MEANING	EXPLANATION
	It is forbidden to throw this product away with household refuse	Worn out devices and other electronic and electrical products must be collected separately and disposed of in compliance with the Act on waste electrical and electronic equipment.
CE 0197	European conformity mark with the number of the notified body	CE mark confirms that the design of the device is compliant with the guidelines of the European Community
	Instructions for Use	Prior to the first use of the device, consult Instructions for Use.

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- 4.1. CONTROL PANEL 14

5. OPERATION 15

- 5.1. CONNECTION TO THE NETWORK OF MEDICAL GASES 15
- 5.2. TOURNIQUET*(S)* CONNECTION 16
- 5.3. TOURNIQUET INFLATION AND DEFLATION. PRESSURE REGULATION IN TOURNIQUETS 16
- 5.4. SETTING TOURNIQUET OCCLUSION TIME 16

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Prior to first use, consult information on how to use the device

1. INDICATIONS FOR USE

Dual adjustable tourniquet control unit is intended to provide compressed air to tourniquets used to control blood flow in extremities in order to produce a bloodless surgical field. The unit is equipped with two independent terminals *(red and blue)*, so the device can work with one or two single tourniquets or one dual tourniquet.

2. TECHNICAL PARAMETERS AND EQUIPMENT

Work factor	Medical compressed air
Feed pressure	400kPa
Number of terminals to provide compressed air to tourniquets	2 independent terminals *(color-coded - red and blue)*
Output pressure (in tourniquets)	0÷500mmHg *(independent regulation for each terminal)*
Time measurement	1s ÷ 99min 59s with 1s increments
Dimensions (W/D/H)	20.5cm x 21cm x 16cm
Weight	3.2kg

The manufacturer reserves the right to introduce design changes.

Equipment supplied with tourniquet control unit

	Name	Pcs.
	Tourniquet control unit air hose - blue cat. no. 30.0035.010 *(spiral version, max. length ≈1.7m)* or 30.0035.1xx * *(straight version)*	1
	Tourniquet control unit air hose - red cat. no. 30.0035.020 *(spiral version, max. length ≈1.7m)* or 30.0035.1xx * *(straight version)*	1
	Compressed air spiral hose cat. no. 30.0037.xxx	1

The hose length and type of quick coupling for connection to the medical gases outlet point acc. to the customer's specification (table 1)

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


4. CONSTRUCTION

1 - Control panel **2** - Handle **3** - Medical compressed air 400kPa terminal **4** - Holes for mounting the unit on the mobile stand **[30.0040.000]**

4.1. CONTROL PANEL

4 - Knob of pressure regulation - red terminal **5** - Knob of pressure regulation - blue terminal **6** - Spiral hose terminal to connect the tourniquet - red **7** - Spiral hose terminal to connect the tourniquet - blue **8** - Manometer of pressure measurements in tourniquet - red terminal **9** - Manometer of pressure measurements in tourniquet - blue terminal **10** - Timer

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
5. OPERATION

5.1. CONNECTION TO THE NETWORK OF MEDICAL GASES

To use the device, it must be supplied with the medical compressed air of a pressure of 400kPa from the hospital medical gases network. Use compressed air spiral hose **[30.0037]** to connect the unit with a network of medical gases.

NOTE:
There are several standards of medical gases network connections *(Table 1)*. Check the compliance of the quick coupling connector standard of the hose B **[30.0037]** with the standard of medical gases outlet point used in the hospital..

Place the device on a flat surface that will provide a good grip for the rubber feet, so the uncontrolled movements of the unit are eliminated. If the case of working with the unit attached to the mobile stand **[30.0040.000]**, use the brakes located at the wheels to immobilize the stand. Before connecting the tourniquet control unit, turn knobs of pressure regulation *(4)* and *(5)* fully counter-clockwise. Use quick coupling connector *(**A**)* to connect the compressed air spiral hose **[30.0037]** with the unit using a medical compressed air terminal located on the rear panel of the device. Connect quick coupling connector *(**B**)* with a network of medical gases – compressed air of 400kPa - the operating rooms are equipped with.

NOTE:
Keep the order of connections performed.

NOTE:
Quick coupling connector *(**B**)* may only be connected to the network of medical gases which is compressed air of 400kPa. Any attempt to connect the tourniquet control unit to the gas outlet other than medical compressed air of 400kPa may result in damage to the quick coupling connector and the device.


When not in use, disconnect the unit from the compressed air network following the procedure opposite to the connection procedure.

NOTE:
The unit when connected to the compressed air network and is not supplying the tourniquet with air, still consumes a small amount of air.

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5.2. TOURNIQUET*(S)* CONNECTION

Use tourniquet control unit air hose – blue **[30.0035.010]** and/or tourniquet control unit air hose – red **[30.0035.020]** to connect the tourniquet control unit with tourniquet*(s)*. For easy identification, air hoses are color-coded in compliance with the colors of the unit terminals. Use quick coupling connectors (***C***) to connect air hoses with spiral hose terminals located on the unit front panel so that the colors of air hoses match the colors of the unit terminals. Connect threaded quick coupling connectors (***D***) of air hoses with two single or one dual tourniquet. In the case of using one single tourniquet, connect it with the red or blue terminal.

NOTE:
In the case of using only one of the terminals, set the pressure in the other one to 0mmHg (turn the knob of the pressure regulation of the unused terminal fully counter-clockwise)

5.3. TOURNIQUET INFLATION AND DEFLATION. PRESSURE REGULATION IN TOURNIQUETS

Use knobs of pressure regulation of red *(4)* and blue *(5)* terminals to inflate, deflate and regulate the pressure in the tourniquets.

NOTE:
Knobs are equipped with protection against accidental changing in the pressure. In order to change the pressure in the tourniquet, pull the knob and then by turning it clockwise or counter-clockwise, change the pressure.

In order to inflate and increase the pressure in tourniquets, turn the knob *(4)* or *(5)*, depending on the regulated terminal, clockwise and observe the corresponding monometer *(8)* or *(9)* where the increase in pressure will be indicated. End the regulation as soon as the desired pressure is achieved.

NOTE:
The maximum pressure for both terminals is 500mmHg.

In order to decrease the pressure in tourniquets, turn the knob *(4)* or *(5)* counter-clockwise. The decrease in pressure will be indicated on the corresponding monometer *(8)* or *(9)*. End the regulation as soon as the desired pressure is achieved. In order to empty the tourniquets, turn the knob *(4)* or *(5)* counter-clockwise fully. The corresponding monometer *(8)* or *(9)* should indicate 0mmHg.

5.4. SETTING TOURNIQUET OCCLUSION TIME

The tourniquet control unit has been equipped with a timer *(10)* which allows for controlling the tourniquet occlusion time. In order to set the time of occlusion, press MIN button to increase the number of minutes. Press SEC button to increase the number of seconds. Press START/STOP button to start the countdown. When the set time of occlusion has passed, sound signal is activated.

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NOTE:
The countdown of the set time of occlusion does not start automatically upon inflating the tourniquet.

To reset the timer settings, press simultaneously both MIN and SEC buttons. Pressing START/STOP button at zeroed settings will cause the timer to start count up the time from 00:00 to the maximum of 99min 59sec. In this mode, the occlusion time cannot be set and the sound signal will not be activated when time has passed.

6. CLEANING AND DISINFECTION

The outer surfaces of the device should be cleaned with a soft, moistened *(not dripping)* cloth with a mild detergent solution then wiped dry. In the case of contaminating the housing of the device with infectious agents *(blood, body fluids, etc.)*, disinfection using disinfectants with a neutral pH is recommended. In such cases, soak a soft cloth in disinfectant, wipe the housing and allow it to dry.

NOTE:
Do not allow the detergents or disinfectants to enter inside the device housing. Should the detergents or disinfectants enter inside the device housing, send the product to the manufacturer for cleaning processes.

7. WARRANTY AND SERVICE

The warranty on the device is granted by the manufacturer for a period of 12 months. The warranty does not include product damage due to: negligence, inappropriate usage, usage not in accordance with the intended purpose, mechanical damage. Warranty services are to be handled solely by the manufacturer - **ChM** or an appointed service point. Unauthorised repairs and modifications will render the warranty void.
The non-warranty services are to be handled also by the manufacturer. The manufacturer shall not be liable for repairs *(and resultant problems)* carried out by entities other than **ChM** or appointed service points.
To maintain trouble-free and safe operation throughout the life of the device, it is recommended to perform annual technical inspections at the manufacturer's site.

The lifetime of the tourniquet control unit is 10 years. After this period, the manufacturer (ChM) may refuse to perform service repairs or technical inspections.

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3. PRECAUTIONS IN USE 31

4. TECHNICAL PARAMETERS AND ACCESSORIES 33

- 4.1. ACCESSORIES SUPPLIED WITH EZO-02 CONTROL UNIT 34
- 4.2. ADDITIONAL ACCESSORIES OF EZO-02 CONTROL UNIT 34

5. CONSTRUCTION 36

- 5.1. Housing components 36
- 5.2. CONTROL PANEL 37
- 5.3. display 38

6. OPERATING INSTRUCTIONS 40

- 6.1. CHANGE OF OCCLUSION PRESSURE 40
- 6.2. CHANGE OF OCCLUSION TIME 40
- 6.3. Tourniquet inflation and deflation 41
- 6.4. Tourniquet tightness testing 42
- 6.5. Changing the basic parameters of EZO-02 control unit 43
- 6.6. ALARM AND INFORMATION SYSTEM 44
- 6.7. Work with the device 47

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1. INTENDED USE

EZO-02 electronic dual tourniquet control unit is a device that generates the pressure in cuffs to produce a bloodless operating field in extremities in surgical procedures. The device maintains in the tourniquets the occlusion pressure defined by the operator, measures the pressure time and monitors the tightness of the system: tourniquet and control unit. Since provided with two independent tourniquet ports, EZO-02 is able to supply pressure to one/ two single cuffs or one dual cuff.

The built-in compressor and DC battery allow for functioning of the device without its connection to the compressed air and electric AC 230V networks.

2. ADVERSE EFFECTS

Prolonged ischemia caused by the use of a tourniquet can result in temporary or permanent damage to tissues, blood vessels and nerves, and in changes in blood clotting.

Excessively inflated tourniquet may cause limb paralysis.

Insufficiently inflated tourniquet may cause passive hyperemia with possible irreversible loss of function.

Intraoperative bleeding may be caused by:

- blood remaining in the limb because of insufficient exsanguination.
- inadequate *(insufficient)* tourniquet pressure applied
- blood entering through the nutrient vessels of the long bones.

3. PRECAUTIONS IN USE

Read the Instructions for Use prior to the first use.

EZO-02 control unit can only be used by qualified medical personnel familiar with the operation of the device and as indicated in this document.

Before each use, visually check the technical condition of the control unit and accessories that will be used with the device *(air hoses, power cord, tourniquets, control unit stand and holder, etc.)*. Should any non-conformances be detected, do not use the device until they are removed or the defective accessory is replaced with a fully functioning ones.

To ensure proper functioning of the device, use only **ChM** original accessories. Using accessories of another manufacturers may result in unplanned functioning of the device resulting in damage to the device and in extreme cases loss of health and even life of a patient.

EZO-02 control unit should be used in the electromagnetic environment described in section 10 of this manual - ELECTROMAGNETIC COMPATIBILITY. Compliance with the guidelines provided therein will allow for proper operation of the device and devices located in its vicinity.

The power cord of the unit should be connected to the power supply socket equipped with protective grounding. Do not connect the device to a network other than AC 230V equipped with protective grounding.

Verify the battery charge status before each use. A fully charged battery will allow the device to work for a long time in the event of 230V AC power outage.

If the power supply gets dirty, before cleaning it, disconnect the device from the 230V AC mains. Clean the device in accordance with the CLEANING AND DISINFECTION instructions.

It is recommended that EZO-02 control unit be cleaned prior to the first use and after each use in accordance with the instructions provided in the chapter *"CLEANING AND DISINFECTION"*. Prior to cleaning, ALWAYS unplug the device from the power supply network to minimize the risk of electric shock.

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DO NOT use the unit in the presence of oxygen, nitrous oxide or a mixture consisting of flammable anesthetics and air.

It is recommended that an alternative *(spare)* device and accessories be available during a surgery since unexpected technical problems can occur.

DO NOT modify any system components yourself. Failure to follow the above-mentioned can result in injury to the medical staff or operator.

DO NOT use the device with a damaged cord or plug.

To ensure proper functioning of the EZO-02 control unit, **ChM** recommends its annual maintenance. Only trained and experienced personnel of the manufacturer can perform such inspections. The manufacturer shall not be liable for any damage resulting from improper operation or maintenance of the product performed by unauthorized persons.

The main power source of the EZO-02 control unit is an external power supply AC 230V. The internal battery power supply should only be used in case of external power supply failure and for patient transportation. Due to the fact that the battery is a consumable element, it is recommended that the battery be replaced at least once a year, or ad hoc when its reduced functional properties are detected.

The EZO-02 electronic dual tourniquet control unit and its accessories cannot be sterilized. This process can destroy the product.

Inflation of a loose, not applied on a limb or tightly rolled up tourniquet may result in its damage.

This is the surgeon that decides on the parameters to be set *(pressure, pressure duration, etc.)*, the location of the cuff to be applied, the inflation and deflation time, and other parameters.

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4. TECHNICAL PARAMETERS AND ACCESSORIES

Power supply	100÷240V, 50÷60Hz
Battery supply	Li-Ion 14.8V
Working time on battery	~ 10 hours
Output pressure *(in tourniquet)*	100-600mmHg with 5mmHg increments
Work factor	Compressed air *(compressed by built-in compressor)*
Mode of operation	Continuous operation
Application part type	BF, EN 60601-1
Ingress Protection Rating	IPX0, EN 60529
Directive 93/42/EEC classification	IIa
Alarm time setting range	1-300min with 1min increments
Type of protection against electric shock	Class I or Internally Powered Equipment (When the unit is operating on backup battery, the type of protection against electric shock changes to internally powered equipment)
EMC	EN 60601-1-2
Applicable standards	EN 60601-1
Dimensions H/W/D	21cm x 21cm x 21.5cm
Height	4.1 kg
Volume of sound signals	45-85 dB

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4.1. ACCESSORIES SUPPLIED WITH EZO-02 CONTROL UNIT

	Name	pcs.
	Power cord AC 230B 1.8 m	1
	Tourniquet control unit air hose – blue Cat. No. 30.0035.010 *(spiral version, max. length ≈1.7m)* or 30.0035.1xx* *(straight version)*	1
	Tourniquet control unit air hose – red Cat. No. 30.0035.020 *(spiral version, max. length ≈1.7m)* or 30.0035.2xx* *(straight version)*	1

* XX – hose length max. 10 m *(e.g. 25 - 2.5m; 05 - 0.5m for 10m - 00)*

4.2. ADDITIONAL ACCESSORIES OF EZO-02 CONTROL UNIT

The elements provided below are additional accessories of the EZO-02 control unit and are offered as a separate accessory. Prior to use, consult Instructions for Use supplied with these products.

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4.2.2. MOBILE STAND FOR TOURNIQUET CONTROL UNITS CAT. NO. 30.0040.000


The stand is used as a mobile platform for the **ChM** family of tourniquet control units, as well as the accessories needed to be used with the device, such as tourniquets, skin protection padding, power cords and other accessories.

Technical parameters	
Dimensions *(H/W/D)*	900 x 570 x 570 mm
The height with the control unit	1100 mm
Weight	9 kg

4.2.3. TOURNIQUET CONTROL UNIT HOLDER – CAT. NO. 30.4351.000


The holder is used to install the unit on supporting rail systems compatible with PN-EN ISO19054.

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
5. CONSTRUCTION

5.1. HOUSING COMPONENTS

1. Control panel *(detailed description in point 5.2 of this manual - Control panel)*

2. USB port *(for use only by service personnel)*

3. Handle - to move the device

4. ON/OFF AC power switch - when the device is plugged in to the 230V AC network and the switch is in OFF position, the unit is disconnected from the mains.

CAUTION: The power switch does not disconnect internal battery power supply.


5. AC 230V power connection socket - place for power cord AC 230V connection

6. Fuse drawer - place for two slow-blow fuses 3.15A 250VAC 5x20mm for short-circuit protection of the device.

7. Mounting holes - used for mounting the control unit on the mobile stand 30.0040.000 or unit holder 30.4351.000 *(For more information, see section 4.2 of this manual)*.

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5.2. CONTROL PANEL

8. LCD display - graphically presents all information related to the EZO-02 operation *(detailed description in section 5.3 of this manual)*

9. Control knob – sets and confirms the parameters and settings of the device *(ENTER)*

10. Unit activation/ deactivation button.

11. MENU activation/ deactivation button.

12. Button for changing the occlusion pressure for the red port

13. Inflation button for the red port

14. Deflation button for the red port

15. Button for changing the occlusion time for the red port

16. Tourniquet red connection port

17. Button for changing the occlusion pressure for the blue channel

18. Inflation button for the blue port

19. Deflation button for the blue port

20. Button for changing the occlusion time for the blue port

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
5.3. DISPLAY

5.3.1. MAIN WINDOW

The main window is displayed after activating the device and while maintaining the occlusion pressure in the tourniquet*(s)*. It presents information about the actual and introduced occlusion pressure and time for the cuff*(s)* connected to the device. The main window also displays icons, warnings and information related to the operation of the control unit.

22. Actual pressure in the red port tourniquet

23. Introduced pressure in the red port tourniquet

24. Actual occlusion time for the red port tourniquet

25. Introduced *(total)* occlusion time for the red port tourniquet

26. Actual pressure in the blue port tourniquet

27. Introduced pressure in the blue port tourniquet

28. Actual occlusion time for the blue port tourniquet

29. Introduced *(total)* occlusion time for the blue port tourniquet

30. Battery charge status icon

	≤ 100%
	≤ 75%
	≤ 50%
	≤ 25%
	≤ 5%

31. 230V AC power icon *(light gray when inactive)*


32. Icon of temporary *(1 minute)* mute of the alarm sound *(light gray when inactive)*

33. Alarm icon *(always with a text message and an audible signal, light gray when inactive)*

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5.3.2. MENU

The menu window is activated after pressing the MENU button on the control panel. Entering menu settings, the parameters of EZO-02 unit control can be modified to suit the operator's preferences. Tightness test of the tourniquet may also be performed.

34. Tightness test of the tourniquet.
35. Screen brightness settings
36. Volume button settings
37. Menu language selection
38. Date and time settings

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6. OPERATING INSTRUCTIONS

6.1. CHANGE OF OCCLUSION PRESSURE

When activated for the first time, the occlusion pressure for both tourniquet ports is set to 200mmHg. The information about the introduced value is displayed in the window [23] for the red port and [27] for the blue port.

Press the PRESSURE button of the port for which the pressure needs to be changed - red [12] or blue [17]. Depending on the selected port, a pulsating mark  will appear in the place of introduced pressure [23] or [27] *(indicating that the device is in pressure editing mode)*, and in place of the actual occlusion pressure [22] or [26], the pulsating pressure value, which will be edited, occurs.

Turn the CONTROL KNOB [9] left or right to change the value of the introduced pressure as required. Afterwards, confirm the selected value by pressing the ENTER button *(pressing the knob)*. The selected value will appear in the introduced pressure window [23] or [27].

For leaving without change confirmation, press the PRESSURE button again *(for the port for which the change was introduced)*.

Change of the occlusion pressure is also possible while the device maintains the pressure in the tourniquet. The procedure for pressure editing is identical to that described above. Having confirmed the new pressure by pressing ENTER *(pressing the knob)*, the device will automatically change the actual occlusion pressure in the tourniquet to the introduced one. To exit without changing the value, press the PESSURE button *(for the port for which the change was introduced)*.

Each time the device is turned on, the occlusion pressure is set to the last value introduced.

6.2. CHANGE OF OCCLUSION TIME

When activated for the first time, the occlusion time for both tourniquet ports is set to 40min. The information about the introduced value is displayed in the window [25] for the red port and [29] for the blue port.

Press the TIMER button of the port for which the time needs to be changed - red [15] or blue [20]. Depending on the selected port, a pulsating mark  will appear in the place of introduced time [25] or [29] *(indicating that the device is in time editing mode)*, and in place of the actual occlusion time [24] or [28], the pulsating time value, which will be edited, occurs.

Turn the CONTROL KNOB [9] left or right to change the value of the introduced time as required. Afterwards, confirm the selected value by pressing the ENTER button *(pressing the knob)*. The selected value will appear in the introduced pressure window [25] or [29].

For leaving without change confirmation, press the TIMER button again *(for the port for which the change was introduced)*.

Change of the occlusion time is also possible while the device maintains the pressure in the tourniquet. The procedure for time editing is identical to that described above. Having confirmed the new time by pressing ENTER *(pressing the knob)*, the device will automatically change the actual occlusion time to the introduced one. To exit without changing the value, press the TIMER button *(for the port for which the change was introduced)*.

CAUTION: **Change introduced to the occlusion time while the device maintains the pressure in the tourniquet does not restart the time counting from the introduced value, but edits the time that was previously set and counted down from.**

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Each time the device is turned on, the occlusion time is set to the last value introduced.

CAUTION: **When the introduced occlusion time is up, completion is signaled graphically and acoustically. The cuff should be deflated** *(p. 6.3 of the manual)* **or the introduced occlusion time must be extended** *(p. 6.2 of the manual)***. Elapsed occlusion time does not deflate the tourniquet.**

6.3. TOURNIQUET INFLATION AND DEFLATION

Press INFLATE [13], [18] button of a given tourniquet port to start inflating the cuff. Upon pressing the INFLATE button, the pressure in the chosen cuff rapidly increases to the set occlusion pressure. The occlusion time countdown for the tourniquet port will start automatically.

The device has been equipped with a protection against uncontrolled deflation of the cuff - for deflating, press and hold pressed the DEFLATE [14], [19] button for 3s. The screen will display the information *„deflating in..."* with 3 seconds counting down. Afterwards, the pressure in the chosen cuff will decrease rapidly to the zero value. Deflation of the cuff automatically stops the occlusion time countdown and resets it to zero.

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
6.4. TOURNIQUET TIGHTNESS TESTING

Turn on the device and go to the menu window by pressing the MENU button [11]. Roll up the tourniquets to be tested tightly and secure against unrolling. Connect it to the red [16] or blue [21] port.

CAUTION: **Launching the menu window is not possible when the control unit maintains the occlusion pressure in the tourniquet.**

Use the CONTROL KNOB [9] and select the CUFFS TESTING window - confirm by pressing ENTER *(the knob)*. Use the KNOB [9] to select the port to which the tested cuff has been connected *(red or blue)*. Confirm by pressing the ENTER button again. The CUFFS TESTING window will appear on the display.

The test starts after pressing the ENTER button. To stop the test, press ENTER again.

After the test, information about the test result is displayed.

CAUTION: **Should the result be negative, a leak was detected in the cuff tested. This cuff must not be used for the procedure! Replace the tourniquet for a new one and restart the test.**

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6.5. CHANGING THE BASIC PARAMETERS OF EZO-02 CONTROL UNIT

- Turn on the device and go to the menu window by pressing the MENU button [11]

CAUTION: **Launching the menu window is not possible when the control unit maintains the occlusion pressure in the tourniquet.**





- Turn the CONTROL KNOB [9] left or right to select the parameter to be changed. Confirm by pressing ENTER. By turning the KNOB [9] again, change the value of the selected parameter. Confirm the change by pressing ENTER. To exit without confirming the change, press the MENU button.
- To enter the main window, press the MENU button.

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



6.6. ALARM AND INFORMATION SYSTEM

The EZO-02 dual tourniquet control unit has been equipped with an alarm and information system informing the user of important parameters or abnormalities occurring during operation of the device.





Tables below list the reasons for a message, the type of a message, and the action to be taken.

Alarm	Messages	Cause	Solution
LOW BATTERY LEVEL	Graphic 	Battery charge status ≤25%	• The device is unplugged from the 230V network. Plug in the device to start the battery charging process. The power supply icon should be activated.
	Sound 		• Power switch [4] on the back panel *(see p 5.1 of this manual)* in the OFF position. Switch it ON to start the battery charging process. The power supply icon should be activated.
BATTERY DISCHARGED	Graphic 	Battery charge status ≤5%	• The device is unplugged from the 230V network. Plug in the device to start the battery charging process. The power supply icon should be activated.
	Sound 		• Power switch [4] on the back panel *(see p 5.1 of this manual)* in the OFF position. Switch it ON to start the battery charging process. The power supply icon should be activated.

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Alarm	Messages	Cause	Solution
SYSTEM LEAKS *(red port)*	Graphic  Sound 	Leakage detected for the red tourniquet port	For the cuff connected to the red port, verify the tightness of pneumatic connections, air hoses and the cuff itself. Remove the leakage. If the leakage cannot be removed, use a different tourniquet or control unit, depending on the defect location.
SYSTEM LEAKS *(blue port)*	Graphic  Sound 	Leakage detected for the blue tourniquet port	For the cuff connected to the blue port, verify the tightness of pneumatic connections, air hoses and the cuff itself. Remove the leakage. If the leakage cannot be removed, use a different tourniquet or control unit, depending on the defect location.

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Alarm	Messages	Cause	Solution
End of Time *(red port)*	Graphic  Sound 	The occlusion time introduced by the operator for the red port tourniquet has elapsed.	Deflate the tourniquet *(p 6.3 of this manual)* or extend the set occlusion time *(p 6.2 of this manual)* for the red port.
End of Time *(blue port)*	Graphic  Sound 	The occlusion time introduced by the operator for the blue port tourniquet has elapsed.	Deflate the tourniquet *(p 6.3 of this manual)* or extend the set occlusion time *(p 6.2 of this manual)* for the blue port.

CAUTION: Audible alarms can be muted for 1 minute by pressing the MENU button. When muted, is activated.

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6.7. WORK WITH THE DEVICE

CAUTION: **Before each use for the procedure, check the technical condition of the device and accessories** *(which will be used with the control unit)* **and the battery charge status. Only technically sound device and accessories may be used for the procedure.**

Turn on the device by pressing the STAND BY button. Use suitable tourniquet ports to install a cuff or cuffs. Perform the cuff tightness test according to section 6.4 of the manual. Positive result of the test means that the device is ready for use. When activated, the last introduced time and occlusion pressure are loaded from the device memory. Apply, as intedned, the tourniquet on a limb, keeping in mind that the pneumatic air hose must be free of bends. Introduce the required pressure and occlusion time. The minimum effective pressure should be determined on the basis of such factors as: the location of the tourniquet applied *(upper or lower limb)*, the patient's systolic blood pressure before surgery and the maximum predicted increase in systolic blood pressure during surgery. Press INFLATE button in the control panel section to which the tourniquet is attached *(red or blue)* to start inflate the cuff. The device will inflate the selected tourniquet to the introduced pressure and will activate the countdown process. If change in pressure or occlusion time is required, these modifications can be introduced when the tourniquet has been inflated. The information on the inflation of the selected cuff will be displayed on the LCD screen.

When the introduced occlusion time has elapsed, text information and an audible signal will be displayed. The tourniquet should be deflated *(p. 6.3 of the manual)* or the introduced occlusion time should be extended *(p. 6.2 of the manual)*.

Immediately, when completely deflated, remove the tourniquet along with the underlying bandages and protective sleeve from the patient.

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7. CLEANING AND DISINFECTION

Before cleaning or disinfection, disconnect the device from 230V network!

The outer surfaces of the device should be cleaned with a soft, moistened *(not dripping)* with a mild detergent solution cloth then wiped dry. In case of contaminating the housing of the device with infectious agents *(blood, body fluids, etc.)*, disinfection using disinfectants with a neutral pH is recommended. In such cases, soak a soft cloth in disinfectant, wipe the housing and allow it to dry.

Do not allow the detergents or disinfectants to enter inside the device housing. Should the detergents or disinfectants enter inside the device housing, switch the device off, disconnect from 230V network and send it to the manufacturer for cleaning processes.

In is forbidden to sterilize the EZO-02 unit control.

8. STORAGE

Store the product in a clean and dry room, off the direct sunlight, at a temperature of 15-30°C and air humidity of 20-70%.

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9. WARRANTY AND SERVICE

The warranty on the device is granted by the manufacturer for a period of 12 months. The warranty does not include product damage due to: negligence, inappropriate usage, usage not in accordance with the intended purpose, mechanical damage. Warranty services are to be handled solely by the manufacturer - **ChM** or an appointed service point. Unauthorised repairs and modifications will render the warranty void.

The non-warranty services are to be handled also by the manufacturer. The manufacturer shall not be liable for repairs *(and resultant problems)* carried out by entities other than **ChM** or appointed service points.

To maintain trouble-free and safe operation throughout the life of the device, it is recommended to perform annual technical inspections at the manufacturer's site

The lifetime of the EZO-02 control unit is 10 years. After this period, the manufacturer (ChM) may refuse to perform service repairs or technical inspections.

For servicing, warranty and post-warranty repairing, please contact the manufacturer's representative or **ChM** directly.

ChM sp. z o.o.
Lewickie 3b
16-061 Juchnowiec Kościelny, POLAND
tel: +48 85 86 86 100
fax: +48 85 86 86 101
e-mail: chm@chm.eu
www.chm.eu

10. THE LIFETIME AND DISPOSAL

The **ChM** sp. z o.o. does not define the maximum number of uses for the EZO-01 control unit. The useful life of the device depends on many factors including the method and duration of each use, applied cleaning and disinfection and the handling between uses. Nontheless, the lifetime of the EZO-02 control unit is 10 years. After this period, the manufacturer (**ChM**) may refuse to perform service repairs or technical inspections.

The device should be disposed of in accordance with the applicable electrical and electronic equipment law.

It is recommended that the battery be replaced at least once a year.

Worn-out lithium-ion battery must be disposed of in accordance with the current law on waste electrical and electronic equipment.

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Problem	Possible causes	Solution
Unit does not turn ON.	Unit not plugged in and battery fully depleted	Plug in the unit
	The unit is plugged in but the switch is OFF and battery is fully depleted	Turn the switch ON
	Damaged power cord	Replace the power cord for the new one
	Battery failure	Send to **ChM** service
	System failure	Send to **ChM** service
The device does not maintain the set pressure	The air hose may be bent or incorrectly connected	Check the hose and the connection
	Damaged tourniquet	Replace the tourniquet for a new one
	Damaged air hose	Replace the air hose for the new one
		Send to **ChM** service
	System error	Send to **ChM** service
	Exceeding the declared lifetime of the device *(10 years)*	Send to **ChM** service for checking the further usefulness of the product
Tourniquet cannot be connected to the unit	Damaged air hose	Replace the air hose for the new one
		Send to **ChM** service
	Using tourniquets other than **ChM**'s	Use **ChM** tourniquets
Discoloration and scratches on the surface of the device	Use of detergents other than those recommended by the manufacturer	Consult the Instructions for Use of the unit
	Overheating of the device	Send to **ChM** service
	Exceeding the declared lifetime of the device *(10 years)*	Send to **ChM** service for checking the further usefulness of the product
Unit switches off	Battery fully depleted	Charge the battery
		Replace the battery for the new one in the manufacture's service
	System error	Send to **ChM** service
	System failure	Send to **ChM** service
	Strong electromagnetic interference	Turn off all electrical equipment not in use; Change the location of electrical equipment; Plug in electrical equipment to other available outlets
Too short battery life	End of battery life	Send to **ChM** service
		Consult the Instructions for Use of the unit
	Failure to follow manufacturer's instructions	Consult the Instructions for Use of the unit
	Exceeding the declared lifetime of the device *(10 years)*	Send to **ChM** service for checking the further usefulness of the product
Information cannot be read from the device display	Display failure	Send to **ChM** service
	System failure	Send to **ChM** service
	System error	Send to **ChM** service
Unit cannot be set to STANDBY	Pressure sensed in tourniquet	Deflate the tourniquet
	System failure	Send to **ChM** service
Alarm indicating a leaking system	The air hose may be bent or incorrectly connected	Check the hose and the connection
	Damaged tourniquet	Replace the tourniquet for a new one
	Damaged pneumatic air hose	Send to **ChM** service
		Replace the air hose for a new one
	System failure	Send to **ChM** service

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




Symbol	Meaning	Description
	Do not discard the product with household waste	Worn devices and other electrical products must be collected separately and disposed of in accordance with the applicable electrical and electronic equipment law.
	Classification Type BF	Degree of protection against electric shock.
CE 0197	European conformity mark with the notified body number	Conform to the European Medical Device Directive 93/42/EEC for class IIa medical devices *(notified body number)*
	Instructions for Use	Prior to first use, consult Instructions for Use
IPxx	Protection degree	Ingress protection degree
	Caution	Read the provided Instructions for Use before operating the device.
	General alarm	Alarms emitted while operating the device
	Manufacturer/ Date of manufacture	The address of the device manufacturer and the date of manufacture of the device.

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13. ELECTROMAGNETIC COMPATIBILITY

The EZO-02 unit is intended for use in an electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such the environment.

13.1. ELECTROMAGNETIC EMISSIONS

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions, CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emission is very low and it is not likely to cause any interference in nearby electronic equipment.
RF emissions, CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions, IEC 61000-3-3	Comply	

13.2. ELECTROMAGNETIC IMMUNITY

Immunity test standard	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge *(ESD)* IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be made from wood, concrete or ceramic tile. If floors are covered with a synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for in/ out lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±1 kV line to line ±2 kV line to earth	
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% Ut *(0.5 cycle)* 40% Ut *(5 cycles)* 70% Ut *(25 cycles)* <5% Ut *(5 seconds)*	<5% Ut *(0.5 cycle)* 40% Ut *(5 cycles)* 70% Ut *(25 cycles)* <5% Ut *(5 seconds)*	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires its continuous use even during power outages, it is recommended that the device be connected to an emergency power supply. Ut - mains supply voltage
Power frequency *(50/60 Hz)* magnetic field IEC 61000-4-8	3A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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Immunity test standard	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	V1 = 3V/m 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d = 1.20 √P d = 1.20 √P 80 MHz to 800 MHz d = 2.30 √P 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts *(W)* according to the transmitter's data provided by the manufacturer and d is the recommended separation distance in meters *(m)*. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a), should be less than the compliance level in each frequency range b). Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 800 MHz 3V/m 800 MHz to 2.5 GHz	E1 = 3V/m 80 MHz to 800 MHz E2 = 3V/m 800MHz to 2.7GHz	

a) Field strengths from fixed transmitters, such as base stations for radio telephones *(cellular/ cordless)* and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b) For the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

13.3. SEPARATION DISTANCES

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters *(m)* can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts *(W)* according to the transmitter's data provided by the manufacturer.

Rated maximum output power of transmitter *(W)*	Separation distance according to frequency of transmitter *(m)*		
	150 kHz to 80 MHz **d = 1.20 √P[m]**	**80 MHz to 800 MHz** **d = 1.20 √P[m]**	**800 MHz to 2.5 GHz** **d = 2.30 √P[m]**
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.30

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