Gebrauchsanweisung – Istruzioni d'uso
Instruction manual – Mode d'emploi
Instrucciones de uso

# Zirkonzahn

# FRESCO, FRESCO DENTINE +, FRESCO DYNAMIC DENTINE, FRESCO LIQUID, FRESCO GEL

CE 1936

GS0118V01-210811

Zirkonzahn srl
Via An der Ahr 7
39030 Gais – ITALY
Tel +39 0474 066 660 – Fax +39 0474 066 661
E-Mail: info@zirkonzahn.com

## Deutsch

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## English

Date of update: 11/08/2021

**Intended use:**
FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE are paste-like ceramic materials for veneering and characterising framework materials made of post-sintered zirconia. They can be used for the manufacture of single crowns, bridges, inlays, onlays, veneers as well as for full contour restorations.
FRESCO LIQUID and FRESCO GEL are used to change the viscosity of the FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE products.
FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE were developed for the use with all Zirkonzahn zirconia products and ICE Stains.

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**Contraindications:**
Should intolerance to the products be suspected, use them only after an allergy screening test has been performed with issue of medical statement of non-allergy conditions. Do not use FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE with severe bruxism.

**Work steps:**
1) Clean and sandblast the post-sintered zirconia structure.
2) Remove ceramic paste from jar. Refresh with FRESCO GEL, if necessary.
3) Then, apply a covering and bubble-free layer of the respective FRESCO product. Ceramic pastes can be applied using a metal tool or a brush for ceramics.
4) Carry out the Fresco firing n. 1 (see firing table FRESCO below).
5) Apply 3D Base Glaze on monolithic structures and carry out glaze firing (see firing table FRESCO below). In the case of slightly reduced teeth, characterisation with stains can be performed directly after the Fresco firing n. 1. After that, the stain firing can be carried out.
6) Carry out glaze and/or stain firing (see firing table FRESCO below), to finalise the monolithic structure in the tooth area (optional).
7) Apply onto the zirconia structure with a metal instrument or brush for ceramics and complete the slightly reduced tooth form. The paste can be diluted with Fresco Liquid, if required.
8) If the pastes are too liquid during application, the applied masses can be dabbed with a handkerchief. After pastes application and modelling, it is recommended to dry the masses slightly in order to work in the final surface structure with a slightly coarser bristle brush.
9) Perform Fresco firing n. 2 (see the firing table FRESCO below). The mass has a slight shrinkage after firing. If the result is not yet satisfactory, the work step can be repeated with the same firing process.
10) Finalise with suitable polishing paste and rubber polishers (optional).
11) Optional: in order to polish the zirconia structure, an additional glaze firing may be carried out. For that, the same parameters of the second firing (glaze and/or stain firing, see firing table FRESCO below) shall be followed, except of the the holding time, which should be 1,5 minutes.

**Caution:**
In order to guarantee reproducibility and colour stability when veneering with FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE, the layer thickness must be selected depending on the colour of the zirconia sub structure. The darker the chosen shade, the more the zirconia structure must be reduced, and the veneering layer increased accordingly. For light dentine shades, however, thin layers ( ≤ 0.4 mm) are already sufficient and therefore a smaller reduction of the zirconia structure is sufficient. This should also be taken into account with any previous Colour Liquid application.

**Firing table FRESCO*:**

| | 1st firing – Fresco firing n. 1 | 2nd Firing – Glaze and/or stain firing | 3rd firing – Fresco firing n. 2 |
|---|---|---|---|
| Temperature T | 820 °C | 800 °C | 790 °C |
| Heating rate | 20-45 °C | 20-45 °C | 20-45 °C |
| Stand-by temperature B | 500 °C | 350 °C | 500 °C |
| Holding time H | 3 min. | 3 min. | 5 min. |
| Closing time S | 8-12 min. | 6-8 min. | 8-12 min. |
| Vacuum start | 450 °C | No vacuum | 450 °C |
| Vacuum end | 820 °C | No vacuum | 789 °C |

*This firing table is applicable for the following products: FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE.

When veneering massive zirconia structures (e.g. full-arch bridges) with FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE, we recommend a closing time of 12 minutes and not to exceed a heating rate of 20-25 °C/min, in order to guarantee a well-sintered ceramic layer. Furthermore, holding times for larger restorations should be extended by 1 minute each.

**Technical data:**

| Product | CTE (*10⁻⁶ K⁻¹) | Type and class according to ISO 6872 |
|---|---|---|
| Fresco | 9.7 | Type 1, Class 1 |
| Fresco Dentine + / Fresco Dynamic Dentine | 9.5 | |

**Storage and stability:**
Products must be stored in tightly sealed original packaging. FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE ceramic products may be stored at max. 35 °C / 95 °F. If stored too cold, the liquids may freeze and must be defrosted before use.

Zirkonzahn srl
Via An der Ahr 7
39030 Gais – ITALY
Tel +39 0474 066 660 – Fax +39 0474 066 661
E-Mail: info@zirkonzahn.com

**Hazard and safety instructions for working with FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE products:**
All FRESCO, FRESCO DENTINE + and FRESCO DYNAMIC DENTINE products are classified as non-hazardous. Therefore, there are no specific safety instructions to be followed.
Observe the safety data sheet and follow general safety and hygiene measures necessary when working with chemicals.

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**N.B.: All medical devices and all products with open packaging are non-returnable**

Gebrauchsanweisung – Istruzioni d'uso
Instruction manual – Mode d'emploi
Instrucciones de uso

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## English

Date of publication: 21/12/2021

**Intended use:**
TEMP BASIC is used for manufacturing provisional crowns and bridges in the anterior and posterior tooth areas.
PRIME and MULTISTRATUM® PRIME are used for manufacturing provisional crowns and bridges in the anterior and posterior tooth area and thanks to their mechanical properties, they are especially suited for large-span bridges.
TEMP PREMIUM FLEXIBLE and MULTISTRATUM® FLEXIBLE are used for manufacturing provisional crowns and bridges in the anterior and posterior tooth areas. Due to the resins' flexibility, slightly diverging insertion directions are being compensated.
DENTURE GINGIVA PRIME and DENTURE GINGIVA FLEXIBLE are used for manufacturing dentures.
THERAPON TRANSPA, TEMP BASIC TRANSPA, TEMP PREMIUM FLEXIBLE TRANSPA and PRIME TRANSPA are designed for the creation of bite splints, orthodontic splints and occlusal splints for bruxism.
SCREW BLANK is used for manufacturing sealing screws for sealing zirconia structures with screw threads.
Once a year, the attending dentist should examine the dental structures and decide whether they can still remain in the mouth of the patient. The attending dentist must inform the patient about this before the beginning of the treatment. After a while, optical and mechanical defects/damages can occur. Should this be the case, Zirkonzahn recommends to replace the restoration with a new one.

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**Contraindications:**
Do not use the products in combination with cloves, nutmeg, pimento and teas containing essential oils.
Should intolerance to the product be suspected, use it only after an allergy screening test has been performed with issue of medical statement of non-allergy conditions.
Avoid contact of the product with liquids that are not specified in the "work steps" section hereunder or with liquids that are not specifically developed for use with this type of resin in the dental field, e.g. lubricants.

**Work steps:**
1) Fix the resin disc in the specific clamping device. Use resin-specific working parameters. Mill dry.
2) Check milled framework for obvious chips or cracks. Only frameworks in perfect condition must be used for the patient.
3) Manual post-processing can be performed. The elaboration is carried out with tungsten-carbide burs. First with rough, subsequently with fine tungsten-carbide bur.
4) Pre-polish with rubber polishers, then polish with pumice stone. The high gloss polishing is carried out with goat-hair brush and Composite Polishing Paste.
5) Sandblast provisional crowns and bridges on the inside before application.
6) Free the final work from dust (steam clean) and clean before insertion.
7) Use eugenol-free cements for the insertion of provisional crowns and bridges. Follow manufacturer's instructions!
8) Use suitable dental bonding agents for applying teeth to denture bases. Follow the manufacturer's instructions!
9) Disinfect the final work with suitable disinfection methods before placement into the patient's mouth.

ATTENTION: Do not overheat the resin when using burs, brushes, steam jets; use only instruments that can ensure appropriate temperature control during each working step.

**Recommended minimum dimensioning for provisional crowns and bridges:**
Wall thickness (cervical): min. 0.7 mm
Wall thickness (occlusal): min. 1.2 mm
Cross section anterior tooth area: min. 12 mm2
Cross section posterior tooth area: min. 16 mm2
Pontics: max. 2

**Safety instructions:**
Grinding dust can irritate eyes, skin or mucous membranes. Please observe the product-specific safety data sheets.

**Technical Data:**

| | E-modulus | Flexural strength | Water sorption | Water solubility |
|---|---|---|---|---|
| **Temp Basic** | 1800 MPa | > 60 MPa | < 25 µg/mm³ | < 2 µg/mm³ |
| **Prime/ Multistratum® Prime/Denture Gingiva Prime** | 3300 MPa | ≥ 105 MPa | < 30 µg/mm³ | < 1 µg/mm³ |
| **Temp Premium Flexible/Screw Blank/Denture Gingiva Flexible** | 2400 MPa | ≥ 105 MPa | < 5 µg/mm³ | < 2 µg/mm³ |
| **Multistratum® Flexible** | 2400 MPa | ≥ 120 MPa | < 5 µg/mm³ | < 2 µg/mm³ |
| **Therapon Transpa** | 2850 MPa | ≥ 110 MPa | < 22 µg/mm³ | < 0.5 µg/mm³ |

**Handling:**
Do not use the product if the packaging is open or damaged upon delivery. Always check the product for impurities and defects. Only use products that are in perfect condition.
All materials can be processed with any Zirkonzahn milling unit using the appropriate tools and milling parameters. Resin discs are already polymerized when supplied, but avoid dropping them nevertheless. Product quality and product durability can be compromised if handled improperly. MULTISTRATUM® FLEXIBLE and MULTISTRATUM® PRIME are resins that have a colour gradient from dentine to enamel. Discs with a diameter of 98 mm are marked with a sticker on the incisal side.
Cleaning in an ultrasonic bath (approx. 40 °C) with gentle cleaning agents.

**Storage and stability:**
Store in original packaging, dust free, dry and at room temperature. Do not expose to physical impact or vibration, avoid contact with liquids, protect discs from contamination.
Keep out of the reach of children.

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## English

Date of update: 10/03/2021

**Intended use:**
Zirkonzahn COLOUR LIQUID PRETTAU® AQUARELL TRANSLUCENT increases the translucency in the incisal area of the framework and is intended for use with the following Zirkonzahn zirconia products: PRETTAU®, PRETTAU® 2 and PRETTAU® 2 DISPERSIVE®. Zirkonzahn COLOUR LIQUID PRETTAU® AQUARELL OPAQUE reduces the translucency and is preferably applied to the inner side of the zirconia structure (e.g. copings, crowns) to cover dark substructures (e.g. abutments). Zirkonzahn COLOUR LIQUID PRETTAU® AQUARELL OPAQUE is intended for use with the following Zirkonzahn zirconia products: PRETTAU® 2, PRETTAU® 2 DISPERSIVE®, PRETTAU® 3, PRETTAU® 3 DISPERSIVE®, PRETTAU® 4 ANTERIOR® and PRETTAU® 4 ANTERIOR® DISPERSIVE®.
When using these products, the instructions for use of the respective zirconia materials must be observed.

**General information:**
The information provided in this instruction manual must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel.
The user is obliged to use the products in accordance with the present instruction manual and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user will be held fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

**Contraindications:**
Should intolerance to the product be suspected, use it only after an allergy screening test has been performed with issue of medical statement of non-allergic conditions.
These products are not intended for use with other materials than the pre-sintered Zirkonzahn zirconia devices mentioned above in the intended use and may not be used for the production of implantable products.
COLOUR LIQUID PRETTAU® AQUARELL TRANSLUCENT may not be used for framework areas, which are exposed to high mechanical stresses.

**Work steps:**
1) Air blow the zirconia structure.
2) Do not moisten structure before colouring.
3) ATTENTION: Shake container well before use!
   In case of storage at too high temperatures, vapours can form inside the bottle and build up pressure. In this case, the container needs to be opened with care under the extractor hood. The liquid can then be used.
4) Colour the areas of the structure, where the translucency should be modified, with a metal-free brush (observe "Intended use" above!).
5) Further colouring of the structure is carried out with COLOUR LIQUID PRETTAU® AQUARELL (observe instructions for use of COLOUR LIQUID PRETTAU®).
6) Let the structure dry under an infrared lamp (drying times: see instructions for use COLOUR LIQUID PRETTAU®).
7) After brushing, close the liquid container firmly.

**Storage and stability:**
Do not store above 25 °C / 77 °F; short-term exceeding of the temperature up to 45 °C / 113 °F is allowed (e.g. during transport). Preferably store the liquids at temperatures below 18 °C/64 °F. If stored too cold, the product may freeze and must be defrosted at room temperature before use.
Do not use after expiration date.

**Hazard and safety instructions for working with Colour Liquid products:**
Colour Liquid products are harmful to health and contact can cause severe skin burns, allergic skin reactions, irritation of the respiratory tract and serious eye damage.
Observe the safety data sheet and follow general safety and hygiene measures necessary when working with chemicals. Wear protective gloves and goggles when working. Avoid inhalation of vapours, mist and aerosols. Wash hands, forearms and face thoroughly after handling.
When not in use, keep container tightly closed and in a well-ventilated place.
If the liquid product gets in contact with eyes, rinse them cautiously with water for several minutes.
If the product is on skin, immediately take off all contaminated clothing and rinse skin with water. Wash contaminated clothing before reuse.
If the product is swallowed, rinse the mouth with water and immediately consult a doctor or hospital. Do not induce vomiting!
If, despite the precautions taken, irritation, rashes or general discomfort occur, seek medical attention.

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## English

Date of publication: 27/06/2022

**Intended use:**
TECNO MED MINERAL, TECNO MED MINERAL DENTINE and TECNO MED MINERAL TISSUE are used for manufacturing reduced crowns and bridges (max. two pontics and min. 13 mm² connection cross-section), telescopic primary crowns, copings as well as frameworks for composite veneered bridges. TECNO MED is designed for manufacturing friction copings on telescopic crowns or attachments.

The materials were designed for a period of five years from the first contact in the patient's mouth. After a period of five years, and then every year, the attending dentist should examine the structure and decide whether it can still remain in the mouth of the patient. The attending dentist must inform the patient about this before the beginning of the treatment. After five years, optical and mechanical defects/damages may occur on the material. Therefore, Zirkonzahn recommends removing the dental structure and replacing it with a new one after the five years.

**General information:**
The information provided in this instruction manual must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel. The user is obliged to use the products in accordance with the present instruction manual and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user is fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

**Contraindications:**
TECNO MED, TECNO MED MINERAL, TECNO MED MINERAL TISSUE and TECNO MED MINERAL DENTINE must not be used for the manufacturing of implants or root posts. TECNO MED MINERAL, TECNO MED MINERAL TISSUE and TECNO MED MINERAL DENTINE must not be used for the manufacturing of structures with more than 2 pontics or conical crowns.
Should intolerance to the products be suspected, use the products of the present instruction manual only after an allergy screening test has been performed with issue of medical statement of non-allergy conditions.
Avoid contact of the product with liquids that are not specified in the "work steps" section hereunder or with liquids that are not specifically developed for use with this type of resin in the dental field, e.g. lubricants.

**Work steps (observe only steps 1-7 and 12-13 for TECNO MED):**
1) Mill resin discs dry, use resin-specific processing parameters.
2) Check milled framework for obvious chips or cracks. Only frameworks in perfect condition must be used for dental restorations.
3) Manual post-processing with tungsten carbide burs can be carried out.
4) Sandblast with aluminium oxide, 50 µm with 2 – 2.5 bar.
5) If the framework is steam cleaned, avoid intense localised heating.
6) Clean structure with ethanol and let it dry.
7) Use Composite Bonding Liquid for Tecno Med 7ml.
8) Slightly polymerise (in order to create a dispersion layer). Duration of the intermediate polymerisation varies depending on the manufacturer of the device. The following serve as guidelines: Spektra 2000 = 90 sec; Spektra LED = 30 sec; Spectra SL 400 = 90 sec; PLC Spectra = 60 sec; Q-PLC = 10 sec.
9) If the polymerised dispersion layer of the Composite Bonding Liquid appears white on the surface, it is over-polymerised and must be removed. Repeat the preceding operation with reduced polymerisation time.
10) Apply opaquer. Observe the manufacturer's indications!
11) Layer with photopolymerising resin (manufacturer of the resin can be selected freely). Follow the general instructions for use of the resin manufacturer.
12) Use of adhesive cement. Follow the general instructions for use of the cement manufacturer.
13) Disinfect the final work with suitable disinfection methods before placement into the patient's mouth.

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**Safety instructions:**
Grinding dust can irritate eyes, skin or mucous membranes. Observe the product-specific safety data sheets.

**Technical Data:**

| | Tecno Med Mineral | Tecno Med Mineral Tissue | Tecno Med Mineral Dentine | Tecno Med |
|---|---|---|---|---|
| Flexural strength | 201 MPa | 189 MPa | 198 MPa | 174 MPa |
| E-modulus | 4900 MPa | 4400 MPa | 5100 MPa | 4300 MPa |
| Water sorption | 7.3 µg/mm³ | 7.2 µg/mm³ | 6.5 µg/mm³ | 4.7 µg/mm³ |
| Water solubility | 1.2 µg/mm³ | 1.2 µg/mm³ | 0.6 µg/mm³ | 1.7 µg/mm³ |

**Handling:**
Do not use the product if the packaging is open or damaged upon delivery. Always check the product for impurities and defects. Only use products that are in perfect condition.
All materials can be processed with any Zirkonzahn milling unit using the appropriate tools and milling parameters.
Do not drop the resin discs. Product quality and product durability can be compromised if handled improperly.

TECNO MED, TECNO MED MINERAL, TECNO MED MINERAL DENTINE and TECNO MED MINERAL TISSUE can be sterilised. If the sterilisation process includes thermic treatment, do not exceed temperatures of 135 °C.

**Storage and stability.**
Store in original packaging, dust free, dry and at room temperature. Do not expose to physical impact or vibration, avoid contact with liquids, protect discs from contamination.

**Information:**
**Manufacturer:** Zirkonzahn srl – Via An der Ahr 7 – 39030 Gais – Italy – www.zirkonzahn.com
**Emergency Phone:** +39 0474 066 660 (Mon. – Fri. 8 am – 6 pm)

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Zirkonzahn®

# ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE, ICE BUILD UP LIQUID, FRESCO GINGIVA, FRESCO ENAMEL, FRESCO LIQUID, FRESCO GEL, FRESCO ESSENTIAL CE 1936

## English

Date of update: 21/12/2021

**Intended use:**
ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE, FRESCO ENAMEL and FRESCO GINGIVA are dental ceramics for veneering and characterising framework materials made of post-sintered zirconia.
ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE and FRESCO ENAMEL can be used for the manufacture of single crowns, bridges, inlays, onlays, veneers as well as for full contour restorations.
FRESCO GINGIVA is a paste-like ceramic material which allows to reproduce the natural gingiva structure in the simplest way.
ICE BUILD UP LIQUID is used for mixing of ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE + and ICE CERAMICS DYNAMIC DENTINE.
FRESCO LIQUID and FRESCO GEL are used to change the viscosity of FRESCO ENAMEL and FRESCO GINGIVA products.
The Set FRESCO ESSENTIAL includes different FRESCO GINGIVA and FRESCO ENAMEL pastes as well as the products FRESCO GEL and FRESCO LIQUID.
ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE, FRESCO ENAMEL and FRESCO GINGIVA have been developed for use with all Zirkonzahn zirconia types and ICE Stains.

**General information:**
The information provided in this instruction manual must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel. The user is obliged to use the products in accordance with the present instruction manual and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user will be held fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

**Contraindications:**
Should intolerance to the products be suspected, use them only after an allergy screening test has been performed with issue of medical statement of non-allergy conditions.
In case of severe bruxism, ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE And FRESCO ENAMEL products should not be used.

**Work steps for use of ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE materials:**
1) Apply a dentine wash (= thin ceramic layer) with a brush on the post-sintered, cleaned, sandblasted and reduced zirconia structure. Mix for the dentine wash the ceramic powder with ICE BUILD UP LIQUID.
2) Carry out wash firing (see firing table ICE CERAMICS below).
3) Apply the different ceramics, which have been mixed with the ICE BUILD UP LIQUID, with a thin brush onto the zirconia structure and design the final tooth anatomy.
4) Carry out dentine firing (see firing table ICE CERAMICS below).
5) Finalise with diamond-burs and diamond-grinding wheels.
6) Wash, dry and, if intended, characterise with stains (see instructions for use ICE Stains)
7) Carry out final stain firing (see firing table ICE CERAMICS below).

**Firing table ICE Ceramics*:**

| | 1st Firing – wash firing | 2nd Firing – dentine firing | 3rd Firing – stain firing/glaze firing |
|---|---|---|---|
| Temperature T | 810 °C | 810 °C | 810 °C |
| Heating rate | 20-35 °C | 20-35 °C | 20-35 °C |
| Stand-by temperature B | 300 °C | 300 °C | 300 °C |
| Holding time H | 5-10 min. | 1-3 min. | 1-2 min. |
| Closing time S | 6-8 min. | 6-8 min. | 6-8 min. |
| Vacuum start | 450 °C | 450 °C | No vacuum |
| Vacuum end | 809 °C | 809 °C | No vacuum |

*This firing table is applicable for the following products: ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE + and ICE CERAMICS DYNAMIC DENTINE.

When veneering massive zirconia structures (e.g. full-arch bridges), we recommend not to exceed a heating rate of 20-25 °C/min and to extend the holding time of the end temperature for at least 2 minutes in order to guarantee a well-sintered ceramic layer.

**Work steps for the use of FRESCO ENAMEL and FRESCO GINGIVA:**
1) Clean and sandblast the post-sintered zirconia structure.
2) Remove the FRESCO ENAMEL or FRESCO GINGIVA paste from the jar. Refresh with FRESCO GEL, if necessary.
3) Then, apply a covering and bubble-free layer of FRESCO ENAMEL on the slightly reduced tooth area. Apply FRESCO GINGIVA 5 covering the area to be veneered. Ceramic pastes can be applied using a metal instrument or a brush for ceramics.
4) Carry out the Fresco firing n. 1 (see firing table FRESCO below).
5) Apply 3D Base Glaze on monolithic structures and carry out glaze firing (see firing table FRESCO below). In the case of slightly reduced teeth, characterisation with stains can be performed directly after the Fresco firing n. 1. After that, the stain firing can be carried out.
6) Carry out glaze and/or stain firing (see firing table FRESCO below), to finalise the monolithic structure in the tooth area (optional).
7) Apply FRESCO ENAMEL onto the zirconia structure with a metal instrument or brush for ceramics and complete slightly reduced tooth form. The paste can be diluted with Fresco Liquid if required.
8) FRESCO GINGIVA 4-6 are mainly applied as a primer. FRESCO GINGIVA 1-3 are used for characterisation of the gingiva. Apply paste with a metal instrument or brush for ceramics onto the zirconia structure.
9) If the pastes are too liquid during application, the applied masses can be dabbed with a handkerchief. After applying and structuring the pastes, it is recommended to dry the masses slightly in order to work in the final surface structure with a slightly coarser bristle brush.
10) Perform the Fresco firing n. 2 (see firing table FRESCO below). The mass has a slight shrinkage after firing. If the result is not yet satisfactory, the work step can be repeated with the same firing process.
11) Finalize with suitable polishing paste and rubber polishers (optional).
12) Optional: in order to polish the zirconia structure, an additional glaze firing can be carried out. For that, the same parameters of the second firing (glaze and/or stain firing, see firing table FRESCO below) shall be followed, except of the the holding time, which should be 1,5 minutes.

**Caution:**
In order to guarantee reproducibility and colour stability when veneering with FRESCO ENAMEL or FRESCO GINGIVA, the layer thickness must be selected depending on the colour of the zirconia substructure. The darker the chosen shade, the more the zirconia structure must be reduced, and the veneering layer increased accordingly. For light dentine shades, however, thin layers (≤ 0.4 mm) are already sufficient and therefore a smaller reduction of the zirconia structure is sufficient. This should also be taken into account with any previous Colour Liquid application.

**Firing table FRESCO*:**

| | 1st firing – Fresco firing n. 1 | 2nd Firing – Glaze and/or stain firing | 3rd firing – Fresco firing n. 2 |
|---|---|---|---|
| Temperature T | 820 °C | 800 °C | 790 °C |
| Heating rate | 20-45 °C | 20-45 °C | 20-45 °C |
| Stand-by temperature B | 500 °C | 350 °C | 500 °C |
| Holding time H | 3 min. | 3 min. | 5 min. |
| Closing time S | 8-12 min. | 6-8 min. | 8-12 min. |
| Vacuum start | 450 °C | No vacuum | 450 °C |
| Vacuum end | 820 °C | No vacuum | 789 °C |

*This firing table is applicable for the following products: FRESCO ENAMEL and FRESCO GINGIVA.

When veneering massive zirconia structures (e.g. full-arch bridges) with FRESCO GINGIVA and FRESCO ENAMEL, we recommend a closing time of 12 minutes and not to exceed a heating rate of 20-25 °C/min, in order to guarantee a well-sintered ceramic layer. Furthermore, holding times for larger restorations should be extended by 1 minute each.

**Technical data:**

| Product | CTE (* 10⁻⁶K -1) | Type and class according to ISO 6872 |
|---|---|---|
| ICE Ceramics / Fresco Enamel | 9.7 | |
| ICE Ceramics Dentine +/ ICE Ceramics Dynamic Dentine / ICE Ceramics Tissue/ Fresco Gingiva | 9.5 | Type 1, Class 1 |

**Storage and stability:**
Products must be stored in tightly sealed original packaging.
FRESCO ENAMEL, FRESCO GINGIVA, FRESCO LIQUID, FRESCO GEL and ICE BUILD UP LIQUID must be stored at max. 35 °C / 95 °F. If stored too cold, the liquids may freeze and must be defrosted before use.

**Hazard and safety instructions for working with ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE, FRESCO GINGIVA and FRESCO ENAMEL:**
All ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE, FRESCO GINGIVA and FRESCO ENAMEL products are classified as non-hazardous. Therefore, there are no specific safety instructions to be followed.
Observe the safety data sheet and follow general safety and hygiene measures necessary when working with chemicals.

It may be possible that not all products included in this instruction manual are available in all countries.

**Information:**
**Manufacturer:** Zirkonzahn srl – Via An der Ahr 7 – 39030 Gais – Italy – www.zirkonzahn.com
**Emergency Phone:** +39 0474 066 660 (Mon. – Fri. 8 am – 6 pm)

**N.B.: All medical devices and all products with open packaging are non-returnable.**

GS0003V01-21122

Zirkonzahn srl
Via An der Ahr 7
39030 Gais – ITALY
Tel +39 0474 066 660 – Fax +39 0474 066 661
E-Mail: info@zirkonzahn.com

CH REP Jan Möstel
Robert-Seidel-Hof 70, 8048 Zürich

Gebrauchsanweisung – Istruzioni d'uso
Instruction manual – Mode d'emploi
Instrucciones de uso

# Zirkonzahn

**ICE CERAMICS, ICE CERAMICS TISSUE, ICE CERAMICS DENTINE +, ICE CERAMICS DYNAMIC DENTINE, ICE BUILD UP LIQUID, FRESCO GINGIVA, FRESCO ENAMEL, FRESCO LIQUID, FRESCO GEL, FRESCO ESSENTIAL** CE 1936

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| Tipologia di Colour Liquid | Da utilizzare con |
|---|---|
| Colour Liquid | Colour Liquid Thinner |
| Colour Liquid Prettau® | Colour Liquid Thinner |
| Colour Liquid Waterbased | Colour Liquid Waterbased Thinner |
| Colour Liquid Prettau® Aquarell | Colour Liquid Waterbased Thinner |
| Colour Liquid Prettau® 2 Aquarell | Colour Liquid Waterbased Thinner |
| Colour Liquid Prettau® 4 Anterior® Aquarell | Colour Liquid Waterbased Thinner |

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## English

Date of update: 27/06/2022

**Intended use:**
Zirkonzahn COLOUR LIQUID is meant for colouring pre-sintered zirconia structures. Zirkonzahn COLOUR LIQUID has been developed for use with the ICE and PRETTAU® zirconia groups. Observe the relative instructions of use when using these products.

| Colour Liquid | For colouring | Characteristics | Drying time* at 70 °C |
|---|---|---|---|
| Colour Liquid | ICE | Acid-based | 40– 60 min. |
| Colour Liquid Prettau® | Prettau® | Acid-based | 20 – 30 min. |
| Colour Liquid Waterbased | ICE | Water-based | 40– 60 min. |
| Colour Liquid Prettau® Aquarell | Prettau® | Water-based | 20– 30 min. |
| Colour Liquid Prettau® 2 Aquarell | Prettau® 2 | Water-based | 20– 30 min. |
| Colour Liquid Prettau® 4 Anterior® Aquarell | Prettau® 4 Anterior® | Water-based | 20– 30 min. |

** The drying time depends on the size of the structure. A longer drying time does not affect the shade result.*

Colour the different zirconia types only with the respective COLOUR LIQUID (see table above).

**General information:**
The information provided in this instruction manual must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel. The user is obliged to use the products in accordance with the present instruction manual and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user will be held fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

**Contraindications:**
Should intolerance to the product be suspected, use it only after an allergy screening test has been performed with issue of medical statement of non-allergic conditions. The product is not intended for use with materials other than pre-sintered zirconia oxide and may not be used for the production of implantable products.

**Work steps:**
**Dipping technique:**
1) Air blow the zirconia structure.
2) Do not moisten structure before colouring.
3) Pour the COLOUR LIQUID from its transport container into a clean and dry glass container for dipping.
4) ATTENTION: Shake container well before use!
   In case of storage at too high temperatures, vapours can form inside the bottle and build up pressure. In this case, the container needs to be opened with care under the extractor hood. The liquid can then be used.
5) If you intend to individually characterise the structure with Colour Liquid intensive colours, this must be carried out before dipping the structure.
6) Dip the structure for 5-10 sec. into the colouring liquids using metal free tweezers.
7) Air-blow.
8) Let the structures dry under an infrared lamp (drying times: see table above)
9) After colouring, pour the liquid back into the transport container and close it firmly.

**Brushing technique:**
1) Air blow the zirconia structure.
2) Do not moisten structure before colouring.
3) ATTENTION: Shake container well before use!
   In case of storage at too high temperatures, vapours can form inside the bottle and build up pressure. In this case, the container needs to be opened with care under the extractor hood. The liquid can then be used.
4) Use a metal-free brush for colouring.
5) Let the structure dry under an infrared lamp (drying times: see table above).
6) After colouring, close the liquid container tightly.

**Dilution:**
The colour effect of COLOUR LIQUID can be adjusted individually. By adding the respective Thinner, the shade can be brightened up.

| Colour Liquid product | Use with |
|---|---|
| Colour Liquid | Colour Liquid Thinner |
| Colour Liquid Prettau® | Colour Liquid Thinner |
| Colour Liquid Waterbased | Colour Liquid Waterbased Thinner |
| Colour Liquid Prettau® Aquarell | Colour Liquid Waterbased Thinner |
| Colour Liquid Prettau® 2 Aquarell | Colour Liquid Waterbased Thinner |
| Colour Liquid Prettau® 4 Anterior® Aquarell | Colour Liquid Waterbased Thinner |

1) Add Thinner to COLOUR LIQUID.
2) Shake or stir the liquids.
3) Further colouring as described above.

**Stabilizer:**
After prolonged use, it is possible that the stabiliser liquid contained in the COLOUR LIQUID evaporates and the zirconia becomes blemished. In this case, 5% of the respective Stabilizer product should be added again to the COLOUR LIQUID.

| Colour Liquid product | Use with |
|---|---|
| Colour Liquid | Colour Liquid Stabilizer |
| Colour Liquid Prettau® | Colour Liquid Stabilizer |
| Colour Liquid Waterbased | Colour Liquid Waterbased Stabilizer |
| Colour Liquid Prettau® Aquarell | Colour Liquid Waterbased Stabilizer |
| Colour Liquid Prettau® 2 Aquarell | Colour Liquid Waterbased Stabilizer |
| Colour Liquid Prettau® 4 Anterior® Aquarell | Colour Liquid Waterbased Stabilizer |

Zirkonzahn srl
Via An der Ahr 7
39030 Gais – ITALY
Tel +39 0474 066 660 – Fax +39 0474 066 661
E-Mail: info@zirkonzahn.com

CH REP Jan Möstel
Robert-Seidel-Hof 70, 8048 Zürich

| Tipologia di Colour Liquid | Da utilizzare con |
|---|---|
| Colour Liquid | Colour Liquid Stabilizer |
| Colour Liquid Prettau® | Colour Liquid Stabilizer |
| Colour Liquid Waterbased | Colour Liquid Waterbased Stabilizer |
| Colour Liquid Prettau® Aquarell | Colour Liquid Waterbased Stabilizer |
| Colour Liquid Prettau® 2 Aquarell | Colour Liquid Waterbased Stabilizer |
| Colour Liquid Prettau® 4 Anterior® Aquarell | Colour Liquid Waterbased Stabilizer |

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**Adding colour pigments:**
Over time, the colour pigments added to the product for a better representation of the colour may evaporate (e.g. if stored in too much light). These colour pigments can be added again using COLOUR LIQUID AQUARELL BIO-PIGMENT (1 cm of liquid in the bottle corresponds to 1-2 drops of COLOUR LIQUID AQUARELL BIO-PIGMENT). Shake COLOUR LIQUID AQUARELL BIO-PIGMENT well before use.

| Aquarell-shade | Bio-pigment |
|---|---|
| A1-D4, 0M1, 0M2, 0M3 | A1 – D4 |
| Cervical A-D | Orange |
| Incisal Blue | Blue |
| Incisal Grey | Grey |
| Incisal Violet / LuNe | Violet |
| Orange 1 | Orange |
| Brown 2 | Brown |
| Tissue A-C | Tissue |

**Storage and stability:**
Do not store above 25 °C / 77 °F; short-term exceeding of the temperature up to 45 °C / 113 °F is allowed (e.g. during transport). Preferably store the liquids at temperatures below 18 °C/64 °F. If stored too cold, the liquid may freeze and must be defrosted at room temperature before use. Do not use after expiration date.
Zirkonzahn COLOUR LIQUID must be stored in a dark and cool environment.

**Hazard and safety instructions for working with Colour Liquid products:**
Colour Liquid products are harmful to health and contact can cause severe skin burns, allergic skin reactions, irritation of the respiratory tract and serious eye damage. Acid-based Colour Liquid products can additionally be corrosive to metals.
Observe the safety data sheet and follow general safety and hygiene measures necessary when working with chemicals. Wear protective gloves and goggles when working. Avoid inhalation of vapours, mist and aerosols. Wash hands, forearms and face thoroughly after handling.
When not in use, keep container tightly closed and in a well-ventilated place.
If the liquid product gets in contact with eyes, rinse them cautiously with water for several minutes.
If the product is on skin, immediately take off all contaminated clothing and rinse skin with water. Wash contaminated clothing before reuse.
If the product is swallowed, rinse the mouth with water and immediately consult a doctor or hospital. Do not induce vomiting!
If, despite the precautions taken, irritation, rashes or general discomfort occur, seek medical attention.

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