Gebrauchsanweisung – Istruzioni d'uso
Instructions for use – Mode d'emploi
Instrucciones de uso

Zirkonzahn®

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Zirkonzahn®

Torque Ratchet Wrench, Screwdriver Handpiece, Adapter Ratchet Wrench, Safety Adapter Ratchet Wrench, Screwdriver, Screwdriver RW, Screwdriver Zirkonzahn MUA, RA Tool Extension

CE

GS0125V01-230209

Zirkonzahn srl
Via An der Ahr 7
39030 Gais – ITALY
Tel +39 0474 066 660 – Fax +39 0474 066 661
E-Mail: info@zirkonzahn.com

CH REP Jan Möstel
Robert-Seidel-Hof 70, 8048 Zürich

English

Date of publication: 09/02/2023

Intended use:
TORQUE RATCHET WRENCH is a reusable manual ratchet tool used with different specially designed Zirkonzahn inserts, such as adapters, screwdrivers and extractors.
The ADAPTER RATCHET WRENCH and the SAFETY ADAPTER RATCHET WRENCH are inserted into the ratchet to provide a firm and secure connection to the different tools, such as the SCREWDRIVER.
The SCREWDRIVER HANDPIECE is a reusable manual tool used with different specially designed Zirkonzahn inserts, such as screwdrivers and extractors.
The SCREWDRIVER and SCREWDRIVER RW are used for manual insertion as well as tightening or loosening of abutments or abutment screws. The SCREWDRIVER (SIAK4xxxLuu) can only be used in combination with one of the adapters or with Zirkonzahn SCREWDRIVER HANDPIECE. The screwdrivers SCREWDRIVER RW and the screwdriver SCREWDRIVER ZIRKONZAHN MUA are used without adapter and can be used alone as well as in combination with the TORQUE RATCHET WRENCH. The screwdrivers of the present instructions for use are available in different versions. Each version is compatible with certain implant systems. Further information about compatibilities can be found in our product catalogue or in our website www.zirkonzahn.com.
The RA TOOL EXTENSION is used to extend the shank of the different insert tools, such as the SCREWDRIVER.

The TORQUE RATCHET WRENCH, the SCREWDRIVER HANDPIECE and the related accessories can be used both in laboratory and patient's mouth.
The TORQUE RATCHET WRENCH can also be used with or without torque limitation. With torque limitation, the ratchet can be used in a torque range from 10 to 35 Ncm. Without torque limitation (when locked), do not exceed a torque of 70 Ncm.
When using the TORQUE RATCHET WRENCH, the SCREWDRIVER HANDPIECE and their accessories, observe and follow the implant and/or abutment manufacturer's instructions.

General information:
The information provided in these instructions for use must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel. The user is obliged to use the products in accordance with the present instructions for use and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user is fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

Contraindications:
Should intolerance to the following alloy components be suspected, use these products only after an allergy screening test has been performed with a medical statement issue of non-allergy conditions to: iron (Fe), manganese (Mn), silicon (Si), chromium (Cr), nickel (Ni), molybdenum (Mo), copper (Cu), niobium (Nb), tantalum (Ta), titanium (Ti), vanadium (V), aluminium (Al) and PEEK.
The users must not modify the products of the present instructions for use, except for those procedures which are expressly allowed. Non-allowed changes result in loss of approval as medical device.

Product components:
The TORQUE RATCHET WRENCH consists of the following components:
1) Ratchet head with opening for different inserts; screwed with threaded joint
2) Housing rod
3) Spring
4) Adjustment sleeve for viewing the set torque
5) Plate
6) Ratchet handle for adjusting the torque

The SCREWDRIVER HANDPIECE consists of the following components:
1) Connection for different tools with right angle connection
2) Screwdriver shank
3) Pivot pin at the rear end (this component is fixed with its own screw in Version 2)

Do not remove the connection from the screwdriver shank. The pivot pin, on the other hand, can be unscrewed from the shank for cleaning and sterilisation. With the new version, first remove the screw using the screwdriver (SIAK4xxxLuu).

The SCREWDRIVER ZIRKONZAHN MUA consists of Holding Screw and Screwdriver Cover.

The RA TOOL EXTENSION and all other SCREWDRIVERs by Zirkonzahn consist of only one component and cannot therefore be disassembled.

The ADAPTER RATCHET WRENCH consists of different pressed parts and must not be disassembled.

The SAFETY ADAPTER RATCHET WRENCH has a specific set screw for fixing the respective tool. For cleaning, this screw can be removed from the SAFETY ADAPTER RATCHET WRENCH using the SCREWDRIVER (SIAK4xxxLuu). Use the SAFETY ADAPTER RATCHET WRENCH only with tightened set screw.

Product models and application:
TORQUE RATCHET WRENCH: the ratchet is used with different inserts. The use depends on the application. Regardless of the use, the TORQUE RATCHET WRENCH can be used with or without torque limitation.
ADAPTER RATCHET WRENCH (SIAK0001) and SAFETY ADAPTER RATCHET WRENCH (SIAK0003): both are adapters for the TORQUE RATCHET WRENCH and are used to insert different Zirkonzahn tools with right angle connection to the ratchet. The possibility of attaching for example dental floss to the SAFETY ADAPTER RATCHET WRENCH and fixing the tool in the adapter can prevent the adapter or tool from falling down, ensuring additional safety.
SCREWDRIVER HANDPIECE: tools with right angle connection, such as the SCREWDRIVER, can be directly inserted into this tool. The firm handle ensures a safe use of the specific tools.
SCREWDRIVER RW: screwdriver in three different lengths for fixing and removing zirconia screws. These screwdrivers can be used alone or in combination with the TORQUE RATCHET WRENCH.
SCREWDRIVER ZIRKONZAHN MUA (SIAK2001): screwdriver for fixing and removing Zirkonzahn Multi Unit Abutments NON HEX. This screwdriver can be used alone or in combination with the TORQUE RATCHET WRENCH.
SCREWDRIVER (SIAK4xxxLuu): screwdriver for fixing and removing screws of different implant systems. These screwdriver versions can only be used in combination with the ADAPTER RATCHET WRENCH, SAFETY ADAPTER RATCHET WRENCH or SCREWDRIVER HANDPIECE. When correctly inserted, the screwdriver makes a click sound. Before use, make sure these screwdrivers are inserted correctly.
The RA TOOL EXTENSION is used to extend the shank of different tools, e.g. the SCREWDRIVER. It may only be used in combination with the ADAPTER RATCHET WRENCH and the SCREWDRIVER HANDPIECE.
The INSTRUMENT CASE UNIVERSAL SCREWDRIVER SET is used to safely store the TORQUE RATCHET WRENCH, the SCREWDRIVER HANDPIECE, the ADAPTER RATCHET WRENCH, the SAFETY ADAPTER RATCHET WRENCH and the SCREWDRIVER. The INSTRUMENT CASE UNIVERSAL SCREWDRIVER SET has been specially conceived to keep the products in it even during sterilisation.

IMPORTANT: before using the INSTRUMENT CASE UNIVERSAL SCREWDRIVER SET for the first time, clean it in a suitable cleaning device.

The list of compatible systems is constantly extended and is available at www.zirkonzahn.com. Not all products are available in all countries.
For all TORQUE RATCHET WRENCH accessories which are not part of these instructions for use, observe the corresponding specific instructions for use.

Before the first use and after each subsequent use on the patient, disassemble, clean and dry the TORQUE RATCHET WRENCH, the SCREWDRIVER HANDPIECE and the SCREWDRIVER ZIRKONZAHN MUA. Immediately before the next use, reassemble the components and proceed with the sterilisation. The SCREWDRIVER, the SCREWDRIVER RW, the ADAPTER RATCHET WRENCH, the SAFETY ADAPTER RATCHET WRENCH, the RA TOOL EXTENSION and the INSTRUMENT CASE UNIVERSAL SCREWDRIVER SET should not be disassembled, but must be cleaned, dried, and sterilised immediately before first use and after each subsequent use on the patient.

Lock the ratchet head to use the ratchet function without torque limitation: turn back the ratchet handle until the ratchet head can move freely. Pull back the housing rod and turn it by 90°. Screw the ratchet handle back in until the torque is at about 35 Ncm.
When using the locked ratchet, do not exceed a torque of 70 Ncm.

1) Loosen the ratchet handle

2) Push back the housing rod
3) Turn the housing by 90°

4) Push the housing rod forward to the blocking position on the ratchet head

5) Fix the housing rod with the ratchet handle

The torque function can be used only when the TORQUE RATCHET WRENCH is not locked. Via ratchet handle, it is possible to set a tightening torque between 10 and 35 Ncm. The marking on the adjusting sleeve shows the value of the set torque on the housing rod. When the set torque is reached, the buckling mechanism is activated to bend the ratchet head. Do not stress the ratchet any further, otherwise the torque would be significantly exceeded with the risk of damaging the TORQUE RATCHET WRENCH.


In case of non-use: turn back the ratchet handle until the ratchet head can move freely.

Sterilisation: turn the ratchet head by 45° so that its longer end lies on the shank and a gap opens allowing the steam to flow through.

Warnings and precautions:
- Do not use the TORQUE RATCHET WRENCH and the SCREWDRIVER HANDPIECE if damaged or with inserts which are not produced by Zirkonzahn or not conceived for the use in combination with these tools.
- The RA TOOL EXTENSION has been specially conceived for use in combination with Zirkonzahn TORQUE RATCHET WRENCH and SCREWDRIVER HANDPIECE. Other combinations, e.g. with contra-angle handpiece, are not allowed.
- The use of non-compatible screwdrivers may cause damage to the screw and possible risks to the patient. Before using the screwdrivers, check whether they are compatible or not.
- Do not use the SCREWDRIVER and SCREWDRIVER RW in combination with any other systems than those specifically indicated. The SCREWDRIVER ZIRKONZAHN MUA has been specially designed for Zirkonzahn Multi Unit Abutment NON HEX and must not be used for any other screws or abutments.
- Make sure the contact surface between screw head and suitable screwdriver is clean of any residues and liquids.
- The products durability depends on how they are used and on the quality of the carried out products maintenance. In any case, before each use of the products described in this instruction for use, carry out a visual inspection for defects.
- The INSTRUMENT CASE UNIVERSAL SCREWDRIVER SET is not a sterile-product. For this reason, keep it in a sealed sterilisation pouch to ensure a proper sterilisation also of the tools it contains. Strictly observe the instructions for use and the durability of the used sterilisation pouch.
- When not using the ratchet, turn the ratchet handle back until the ratchet head can move freely.
- Each ratchet component and in particular the different inserts are very small. Therefore, take the necessary countermeasures to exclude any risk of ingestion. In addition, before each use check the ratchet for possible missing parts. Use the product only in perfect condition.
- Repairs can only be carried out by Zirkonzahn.
- It is recommended to calibrate the TORQUE RATCHET WRENCH once a year. To do that, please contact the Zirkonzahn Support Team at the headquarters in Gais (Tel. +39 0474 066 630 or e-mail: info@zirkonzahn.com).
- Clean and sterilise the components before possible dispatch. During shipment, pack sharp or pointed parts (screwdrivers, drills) in such a way that any risk of injury is excluded.

Use of the TORQUE RATCHET WRENCH:
1) Before each use, check all products for damage and contamination. Use only products in perfect condition. Before use in the patient's mouth, sterilise the products.
2) Before use, insert the SCREWDRIVER into the ADAPTER RATCHET WRENCH or SAFETY ADAPTER RATCHET WRENCH. Make sure the screwdrivers snap in correctly (click sound). When using the SAFETY ADAPTER RATCHET WRENCH, well fix the set screw to guarantee a firm hold of the tool. Use the SCREWDRIVER SIAK (SIAK4xxxLuu) to tighten the screw. In addition, the specific safety concept of the SAFETY ADAPTER RATCHET WRENCH requires the use of a tear-resistant thread, e.g. dental floss. This can be attached through one of the two eyelets at the adapter end and fastened with knots. It is recommended to well fix the other thread end to a finger. When using the ADJUSTING ADAPTER RATCHET WRENCH, please observe the corresponding instructions for use.
3) If for a specific patient case is required a longer screwdriver, the SCREWDRIVER can be inserted into the RA TOOL EXTENSION. This may only be inserted into the TORQUE RATCHET WRENCH with the ADAPTER RATCHET WRENCH.

4) First tighten the screw lightly by hand and then tighten it with the TORQUE RATCHET WRENCH.
5) When using the TORQUE RATCHET WRENCH with torque limitation, set the torque.
6) Before tightening or loosening the screw, make sure the corresponding SCREWDRIVER is correctly seated in the screw head.
7) Tighten carefully, making sure the TORQUE RATCHET WRENCH is held as horizontally as possible. As soon as the buckling mechanism is triggered, release the ratchet.
CAUTION! Risk of breakage: after the buckling mechanism is triggered, do not stress the TORQUE RATCHET WRENCH any further as this can damage the screw and the ratchet itself.
8) Carefully remove from the screw the TORQUE RATCHET WRENCH together with the screw attachment.
9) Immediately after use, immerse the TORQUE RATCHET WRENCH with the used inserts in a suitable cleaning solution to remove any impurities. Observe the cleaning solution instructions.
10) Then clean the TORQUE RATCHET WRENCH and the inserts; sterilise them immediately before each use (see "sterilisation and cleaning").

Use of the SCREWDRIVER HANDPIECE:
The SCREWDRIVER HANDPIECE has no mechanism to regulate or monitor the torque. It may therefore only be used if the torque is not relevant.
Its use is limited to lightly tightening denture screws, removing screws and using the TITANIUM BASE EXTRACTOR.

Before each use, check all products for damage and contamination. Use only products in perfect condition. Before use in the patient's mouth, sterilise the products.

Tighten the rear screw of the SCREWDRIVER HANDPIECE to better hold the screwdriver.

Using the side connection, insert the corresponding tool into the SCREWDRIVER HANDPIECE until it snaps in correctly with an audible click. If for a specific patient case is required a longer screwdriver, the SCREWDRIVER can be inserted into the RA TOOL EXTENSION.

Before tightening or loosening the screw, make sure the corresponding SCREWDRIVER is correctly seated in the screw head as well as the TITANIUM BASE EXTRACTOR in the base. Align the SCREWDRIVER HANDPIECE as perpendicular to the base as possible to avoid possible leverage.

When tightening, apply only very low strength. As soon as resistance is perceived, use the TORQUE RATCHET WRENCH and follow the points above.

When unscrewing the screw or base, increase the applied strength slowly, as an excessive strength could cause a base or screw breakage and/or injury to the patient.

Immediately after use, immerse the SCREWDRIVER HANDPIECE with the used inserts in a suitable cleaning solution to remove any impurities. Observe the instructions of the cleaning solution manufacturer.

Then clean the SCREWDRIVER HANDPIECE and the inserts; sterilise them immediately before each use (see "sterilisation and cleaning").

Store in the INSTRUMENT CASE UNIVERSAL SCREWDRIVER SET only cleaned medical devices.

Important: in any case, observe the implant or abutment manufacturer's instructions regarding the method of screwing and the recommended tightening torque. The tightening torque of Zirkonzahn screws is indicated on the packaging or available at http://www.zirkonzahn.com/assets/files/brochueren/INT-Drehmomenttabelle-fuer-Zirkonzahn-schrauben.pdf.

Measurement accuracy:
The measuring accuracy of the TORQUE RATCHET WRENCH is ensured by the spring used in the calibration and has the following tolerances:

Torque	Tolerance
10 – 15 Ncm	+/- 1.50 Ncm
20 – 25 Ncm	+/- 1.85 Ncm
30 – 35 Ncm	+/- 1.95 Ncm

Sterilisation and cleaning:
ATTENTION: possible incorrect torque limitation due to wrong combination! Spring and adjusting sleeve of the TORQUE RATCHET WRENCH are perfectly matched to each other in each ratchet sample. For this reason, make sure to always use the original two components. Pay particular attention when handling these components to ensure they cannot be confused with similar ones, e.g. by using a different container for each TORQUE RATCHET WRENCH during cleaning.
The INSTRUMENT CASE UNIVERSAL SCREWDRIVER SET is usually not exposed

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to any special impurities. For this reason, this product can also be cleaned mechanically without manual pre-cleaning. The box must be opened and placed in the cleaning machine.
1) Before cleaning, disassemble the TORQUE RATCHET WRENCH, the SCREWDRIVER HANDPIECE (version 1) and the SCREWDRIVER ZIRKONZAHN MUA (see "Disassembly").
2) The SCREWDRIVER HANDPIECE (version 2) can also be cleaned as a whole unit. Remove any visible residue with a paper towel and rinse with water. Immerse all products in warm water (approx. 25 °C) for at least 10 minutes.
3) Immerse the products in a warm (approx. 25 °C) enzymatic cleaning solution. Follow the indications of the cleaning solution (e.g. 2% for 1 minute).
4) Remove from the product all visible residues with a soft plastic brush under cold running water (approx. 10 °C).
5) Clean all channels with an irrigation cannula with 20 ml enzymatic cleaning solution and check for visible impurities. If there are still any, repeat this step.
6) Rinse the components under tap water until no residues of the cleaning solution remain on the device. Clean all lumina and channels with running water and an irrigation cannula.
7) Place the components into a cleaning/disinfection device and follow the manufacturer's instructions.
8) Dry with compressed air and a lint-free, single-use wipe.
9) Check all components for damage and corrosion. If such defects are found, do not use the tool anymore. Repairs may only be carried out by Zirkonzahn.
10) Before use into the patient's mouth, sterilise all components with the following procedure: steam sterilisation using vacuum method (at 132 °C for 4 min.). Reassemble the tools before sterilisation (see "Assembly"). When cleaning the TORQUE RATCHET WRENCH, loosen the spring until the ratchet head can easily move. The ratchet should only be sterilised in the correct position (ratchet head turned by 45°, see illustration). Product sterilisation is user's responsibility.

The rear screw of the SCREWDRIVER HANDPIECE must (version 1) or can (version 2) remain unscrewed during the sterilisation process.

Disassembly - TORQUE RATCHET WRENCH:
1) Remove the insert from the head of the TORQUE RATCHET WRENCH.
2) Unscrew the ratchet handle until it comes off the housing rod.
3) Hold the TORQUE RATCHET WRENCH by its head and pull off the plate, the adjusting sleeve, the housing rod and the spring one after the other.
4) For optimal cleaning (strongly recommended in case of heavy soiling), the ratchet head can also be removed by loosening the screws (hexagon 1.3) and then cleaned. This component can be divided into six small parts: threaded joint, pressure spring, ratchet sleeve, ratchet housing and two screws.

Disassembly - other tools:
When using the ADAPTER RATCHET WRENCH or SAFETY ADAPTER RATCHET WRENCH, remove the inserted tool before cleaning.
The SCREWDRIVER ZIRKONZAHN MUA and the SCREWDRIVER HANDPIECE are disassembled into two components by unscrewing the inner screw.
For the SCREWDRIVER HANDPIECE (version 2), the rear side screw must be removed with the screwdriver (SIAK4xxxLuu). However, the design of this medical device has been changed so that the screwdriver no longer needs to be disassembled for cleaning and sterilisation.

Assembly:
1) After cleaning, sterilising and drying each components, the TORQUE RATCHET WRENCH, the SCREWDRIVER ZIRKONZAHN MUA and the SCREWDRIVER HANDPIECE can be reassembled.
2) When assembling the head of the TORQUE RATCHET WRENCH, make sure that the incision on the ratchet sleeve points in the direction of the arrow on the "IN" side.
3) Hold the TORQUE RATCHET WRENCH by its head and reinsert the housing rod, the spring, the adjusting sleeve and the plate one after the other. Important: observe the insertion direction of the adjusting sleeve, otherwise too low torques will be achieved. The gear-like end of the adjusting sleeve must point towards the ratchet head. Turn in the ratchet handle until is firmly fixed.
4) Check whether the buckling mechanism works in both directions when not locked.

The SCREWDRIVER ZIRKONZAHN MUA and the SCREWDRIVER HANDPIECE are put back together by screwing in the inner screw.

Storage and stability:
The products are delivered in a non-sterile state and have to be stored in dry and dust-free conditions. Retain all information affixed to the original packaging throughout the life cycle of the device in such a way it is identifiable as belonging to the specific medical device.

US Federal Law restricts the sale of these products to the following professions: doctors, dentists and licensed practitioner in the dental field.

Information:
Manufacturer: Zirkonzahn srl – Via An der Ahr 7 – 39030 Gais – Italy – www.zirkonzahn.com
Emergency Phone: +39 0474 066 660 (Mon. – Fri. 8 am – 6 pm)

N.B.: All medical devices and all products with open packaging are non-returnable.

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Zirkonzahn®

Universal Extractor RW 1.6, Universal Extractor 1.6, Universal Extractor RW 2.0, Universal Extractor 2.0, Ti-Base Extractor RW, Ti-Base Extractor A15, Universal Extractor

CE

English

Date of publication: 15/02/2023


Indications:
UNIVERSAL EXTRACTOR RW 1.6, UNIVERSAL EXTRACTOR 1.6, UNIVERSAL EXTRACTOR RW 2.0, UNIVERSAL EXTRACTOR 2.0, TI-BASE EXTRACTOR RW and TI-BASE EXTRACTOR A15 are used to loosen and remove fixed abutments or titanium bases from implants or laboratory analogues. The UNIVERSAL EXTRACTOR can also be used to remove directly screwed superstructures e.g. made with metal, zirconia or resin. Zirkonzahn's extractors are suitable for use both in the laboratory and in the patient's mouth.

General information:
The information provided in these instructions for use must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel. The user is obliged to use the products in accordance with the present instructions for use and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user will be held fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

Contraindications and side effects:
Should intolerance to the product be suspected, use it only after an allergy screening test has been performed with issue of medical statement of non-allergic conditions. The TI-BASE EXTRACTOR and the TI-BASE EXTRACTOR A15 are made with different stainless steel alloys. The TI-BASE EXTRACTOR RW contains nickel. The UNIVERSAL EXTRACTOR is made of Ti-6Al-4V ELI.
The user must not modify the products of the present instructions for use, except that these procedures are expressly allowed. Non-allowed modifications result in loss of approval as medical device.

Product components:
The TI-BASE EXTRACTOR RW and the TI-BASE EXTRACTOR A15 are available as one-piece instruments.

The UNIVERSAL EXTRACTOR RW 1.6 and RW 2.0 consist of three parts that can be screwed together:
1) Shank
2) Locknut
3) Mandrel with handle

ATTENTION: The locknut is firmly screwed to the shank and must not be loosened.

The UNIVERSAL EXTRACTOR 1.6 and 2.0 are available as one-piece instruments.

Product models and application:
TI-BASE EXTRACTOR RW (= EXAAxx(x)yy, except EXAA01xx and EXAA10xx) and TI-BASE EXTRACTOR A15 (= EXABxx(x)yy):
One-piece tools for removing abutments/titanium bases from implants or laboratory analogues. Each TI-BASE EXTRACTOR has been specifically designed for one implant system and must be used only with the corresponding system. The list of compatible systems is constantly extended and is available at www.zirkonzahn.com.
Some TI-BASE EXTRACTORS can also be used as screwdrivers for the screw of the corresponding implant system.

The TI-BASE EXTRACTOR can be used alone or in combination with Zirkonzahn's Torque Ratchet Wrench. The TI-BASE EXTRACTOR A15 can only be used in combination with Zirkonzahn's Adapter Ratchet Wrench/Adjusting Ratchet Wrench/Safety Adapter Ratchet Wrench.

UNIVERSAL EXTRACTOR RW 1.6 and RW 2.0 (= EXAA01xx and EXAA10xx):
Tools composed of three pieces, two of which are firmly screwed. The special design of these extractors makes them suitable for removing abutments/titanium bases and different superstructures without internal threads for different implant systems. Any information about compatibility with a specific implant system must be requested from Zirkonzahn. The UNIVERSAL EXTRACTOR RW 1.6 and RW 2.0 can be used alone or in combination with Zirkonzahn's Torque Ratchet Wrench.

UNIVERSAL EXTRACTOR 1.6 and 2.0 (= EXAA01xx and EXAA10xx):
The special design of these extractors makes them suitable for removing abutments/titanium bases and different superstructures without internal threads for different implant systems. Any information about compatibility with a specific implant system must be requested from Zirkonzahn. The UNIVERSAL EXTRACTOR 1.6 und 2.0 can only be used in combination with Zirkonzahn's Adapter Ratchet Wrench/Adjusting Ratchet Wrench/Safety Adapter Ratchet Wrench.

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Warnings, precautions and safety instructions:
- TI-BASE EXTRACTOR RW and TI-BASE EXTRACTOR A15 products must only be used in combination with the intended implant system.
- During insertion, make sure that the contact surfaces of the implant and the corresponding extractor are clean and free of residues and liquids.
- There is no specific indication regarding the period of time for which the UNIVERSAL EXTRACTOR, TI-BASE EXTRACTOR RW and TI-BASE EXTRACTOR A15 work properly. This depends on how the product is utilised and maintained by the user. In any case, before each use of the products described in this instruction for use, carry out a visual inspection for defects.
- Due to the small size of some product components, when using UNIVERSAL EXTRACTOR, TI-BASE EXTRACTOR and TI-BASE EXTRACTOR A15 special precautions must be taken to avoid ingestion or inhalation of these components by the patient.
- If the products are used in combination with the Adjusting Adapter Ratchet Wrench and the Torque Ratchet Wrench, their specific instructions for use must also be taken into account.
- The use of the UNIVERSAL EXTRACTOR requires skill. It is therefore recommended to practise on a laboratory model first.
- Do not use the UNIVERSAL EXTRACTOR, TI-BASE EXTRACTOR RW and TI-BASE EXTRACTOR A15 if they are damaged or modified.
- The products packaging cannot be autoclaved.

Work steps:
1) Select the suitable extractor and clean, dry and sterilise it before the first use on the patient and after each subsequent use according to the procedure below. Sterilise the product immediately before each use.
2) In the first step, remove the clinical screw of the crown/abutment/titanium base.
3) Remove the structure:
 => For all extractors, make sure that the torque used is very low.
 a) TI-BASE EXTRACTOR: Screw the extractor lightly into the internal thread. At this point, if desired, attach the connection end of the extractor into the ratchet. Screw the extractor until the implant is reached. Turn it further to loosen and remove the abutment/titanium base from the implant/analogue.
 b) TI-BASE EXTRACTOR RW: Before use, attach the extractor to the Adapter Ratchet Wrench/Adjusting Ratchet Wrench/Safety Adapter Ratchet Wrench. Screw the extractor lightly into the internal thread. At this point, if desired, attach the connection end of the extractor into the ratchet. Screw the extractor until the implant is reached. Turn it further to loosen and remove the abutment/titanium base from the implant/analogue.
 c) UNIVERSAL EXTRACTOR RW 1.6 and RW 2.0: Screw firmly the locknut to the shank. Hold the extractor on the locknut either by hand or with suitable instruments. To facilitate this, the two sides are provided with holes. Insert the shank into the hole of the part to be removed until the shank appears at the end of the hole. As the lateral pressure is reduced, it is easier to hold the extractor. Do not bring in further the UNIVERSAL EXTRACTOR: this can lead to complications with certain implant systems. Then, screw the mandrel into the shank as deeply as possible. Loosen and remove the abutment/titanium base from the implant/analogue by turning it further.
 d) UNIVERSAL EXTRACTOR 1.6 and 2.0: Before use, attach the extractor to the Adapter Ratchet Wrench/Adjusting Ratchet Wrench/Safety Adapter Ratchet Wrench. Hold the extractor either by hand or with suitable tools at the lock nut. To facilitate this, both sides are provided with holes. Insert the shank into the hole of the part to be removed until the shank appears at the end of the hole. As the lateral pressure is reduced, it is easier to hold the extractor. Do not bring in further the UNIVERSAL EXTRACTOR: this can lead to complications with certain implant systems. Screw the arbor into the shaft as deeply as possible. Loosen and remove the abutment/titanium base from the implant/analogue by turning it further.

Sterilisation and cleaning:
1) Disassemble the UNIVERSAL EXTRACTOR RW 1.6 and RW 2.0 parts before cleaning and sterilisation. To do so, with the handle turn the shank out of the mandrel. ATTENTION: The locknut must not be loosened from the shank.
2) Remove any visible residue with a paper towel and rinse with water.
3) Immerse the product in warm water (approx. 25 °C) for a minimum of 10 minutes.
4) Immerse the product in a warm (approx. 25 °C) enzymatic cleaning solution. Follow the indications of the cleaning solution (e.g. 2 % for 1 minute).
5) Remove all visible residues from the product with a soft plastic brush under cold running water (approx. 10 °C).
6) Clean lumina and channels with an irrigation canula with 20 ml enzymatic cleaning solution and check for visible impurities. If there are still any, repeat this step.
7) Rinse the product under tap water until no residues of the cleaning solution remain on the device. Clean all lumina and channels with running water and an irrigation canula.
8) Load the product into a cleaning/disinfection system and follow the manufacturer's instructions.
9) Dry with compressed air and lint-free, single-use wipes.
10) Sterilise all parts before inserting them into the patient's mouth with the following procedure: steam sterilisation with vacuum; at 132 °C for 4 minutes. The sterility of the product is the user's responsibility.

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It may be possible that not all products included in these instructions for use are available in all countries.
⚠ US Federal Law restricts the sale of these products to the following professions: doctors, dentists and licensed practitioner in the dental sector.

Disposal:
Medical devices that have already been used on patients have to be considered as contaminated devices and must be disposed of in accordance with local regulations. Non-contaminated products can be disposed of as industrial waste. A separate collection of reusable material (e.g. metal, plastics = packaging) is recommended.

CAUTION – Risk of injury: any product breakage may result in sharp edges and/or pointed corners. These products must be disposed of in suitable containers in accordance with the current regulations.

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Zirkonzahn®

Universal Extractor RW 1.6, Universal Extractor 1.6, Universal Extractor RW 2.0, Universal Extractor 2.0, Ti-Base Extractor RW, Ti-Base Extractor A15, Universal Extractor

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English

Date of publication: 15/02/2023

Indications:
IMPRESSION COPING products are prefabricated tools for transfering the intraoral position of an implant or abutment to a model. They are available for both open (IMPRESSION COPING OPEN TRAY HEX / IMPRESSION COPING OPEN TRAY NON HEX) and closed (IMPRESSION COPING CLOSED TRAY HEX / IMPRESSION COPING CLOSED TRAY NON HEX) impressions.

IMPRESSION COPING and IMPRESSION COPING SCREW are available for different implant systems and can be used only in combination with the specified implant system.

General information:
The information provided in these instructions for use must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel. The user is obliged to use the products in accordance with the present instructions for use and with appropriate hygiene measures as well as to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user will be held fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

Contraindications and side effects:
Should intolerance to the following alloy components be suspected, use it only after an allergy screening test has been performed with a medical statement issue of non-allergy conditions to: titanium (Ti), vanadium (V), aluminium (Al), iron (Fe), manganese (Mn), chromium (Cr), nickel (Ni), nitrogen (N), hydrogen (H), silicon (Si) and carbon (C).
The user must not modify the products described in the present instructions for use, except for expressly specified operations. Non-allowed modifications result in loss of approval as medical device.

Components IMPRESSION COPING CLOSED TRAY HEX / IMPRESSION COPING CLOSED TRAY NON HEX:
1) Sealing cap
2) Impression coping with individual implant connection
3) Screw for impression coping

Components IMPRESSION COPING OPEN TRAY HEX / IMPRESSION COPING OPEN TRAY NON HEX:
1) Impression copings with individual implant connection
2) Screw for impression coping
3) Driver

Component IMPRESSION COPING ZIRKONZAHN MUA SCREW:
1) Driver
2) Impression coping with individual implant connection

Product models and application:
- IMPRESSION COPING CLOSED TRAY HEX / IMPRESSION COPING CLOSED TRAY NON HEX: set for impression taking by means of a closed impression tray. It consists of: impression coping for closed impression, screw and cap. Recommended for patients with narrow mouth opening or pronounced gag reflex as well as for implants placed in narrow spaces, e.g. in the molar region.
- IMPRESSION COPING OPEN TRAY HEX / IMPRESSION COPING OPEN TRAY NON HEX: set for impression taking by means of an open (perforated) impression tray. It consists of: impression coping for open impression, screw and driver. For standard procedures and particularly recommended for the simultaneous impression of several implants/abutments, especially in case of strongly divergent implant/abutment angles.
- IMPRESSION COPING ZIRKONZAHN MUA OPEN TRAY HEX / IMPRESSION COPING ZIRKONZAHN MUA OPEN TRAY NON HEX: set for impression taking by means of an open (perforated) impression tray, consisting of impression coping and driver and specially conceived for Zirkonzahn Multi Unit Abutments. This system works on the same principle of IMPRESSION COPING OPEN TRAY HEX / IMPRESSION COPING OPEN TRAY NON HEX, but the products design is different. Zirkonzahn also offers its standard solutions (open and closed impression) for Zirkonzahn Multi Unit Abutments set.
- IMPRESSION COPING OPEN TRAY PIN: extension for the IMPRESSION COPING SCREW. It prevents impression material from penetrating the IMPRESSION COPING OPEN TRAY HEX / IMPRESSION COPING OPEN TRAY NON HEX. It is included in the set for open impressions and can be also purchased separately.
- IMPRESSION COPING CLOSED TRAY CAP: sealing cap for closing the impression coping in the closed impression system IMPRESSION COPING CLOSED TRAY HEX / IMPRESSION COPING CLOSED TRAY NON HEX. It is included in the set for closed impressions and can be also purchased separately.
- IMPRESSION COPING SCREW: screw for both impression systems (open and closed). It is included in both sets and can be also purchased individually.
- Zirkonzahn impression copings are provided with or without anti-rotation protection.

⚠ Not all models are available in all countries.

Warnings, precautions and safety instructions:
- IMPRESSION COPING products may only be used in combination with the original implant or abutment.
- During insertion, make sure that the contact surface between the impression coping and the implant/abutment is clean and free of residues and liquids.
- There is no specific indication regarding the period of time until the IMPRESSION COPING products work properly. It depends on how the product is used and maintained. Anyway, before each use of the products described in these instructions for use, carry out a visual inspection for defects.
- Do not use damaged or altered IMPRESSION COPING products.
- Due to their small size, when using the IMPRESSION COPING products, special precautions must be taken to avoid ingestion or inhalation by the patient.
- Screws may only be tightened observing torque indications.
- Impression copings must not be rotated when fixed onto the implant/abutment.
- The products packaging cannot be autoclaved.
- Screws and impression copings of the different systems are very similar and can be therefore confused easily. For this reason, it is advisable to assemble the set parts correctly and sterilise them together.

GS0121V01-230215

Zirkonzahn srl
Via An der Ahr 7
39030 Gais – ITALY
Tel +39 0474 066 660 – Fax +39 0474 066 661
E-Mail: info@zirkonzahn.com

CH REP Jan Möstel
Robert-Seidel-Hof 70, 8048 Zürich

Work steps:
1) Before first use and after each subsequent use on the patient, the products must be cleaned, dried and sterilised according to the procedure described below. Sterilisation must take place immediately before each use.
2) OPEN IMPRESSION: insert the IMPRESSION COPING SCREW in the IMPRESSION COPING OPEN TRAY HEX / IMPRESSION COPING OPEN TRAY NON HEX and tighten it onto the implant/abutment using the IMPRESSION COPING OPEN TRAY PIN. Use a Zirkonzahn Screwdriver compatible with the "Zirkonzahn" system (SIAK4001Lxx) to screw in and out the driver. Observe the tightening torque on the product label.
CLOSED IMPRESSION: place the IMPRESSION COPING CLOSED TRAY HEX / IMPRESSION COPING CLOSED TRAY NON HEX on the implants/abutments and fix it with the corresponding IMPRESSION COPING SCREW. Use a Zirkonzahn "Screwdriver White Metal Scanmarker and Impression Coping" (SIAK4002Lxx) to screw in and out the screw. Observe the tightening torque on the product label. After that, secure the IMPRESSION COPING CLOSED TRAY CAP on the impression coping. Once the impression has been taken, the cap remains attached to the impression.
3) Check the exact fit of the impression copings in the implant/abutment. It is recommended to take an X-ray.
4) Apply suitable impression material around the impression coping and close the gap between the impression coping and the gingiva. Fill the impression tray with impression material and place it on the impression coping. Observe the instructions of the respective impression material.
5) If the IMPRESSION COPING CLOSED TRAY HEX / IMPRESSION COPING CLOSED TRAY NON HEX is used, secure the IMPRESSION COPING CLOSED TRAY CAP on the screwed impression coping.
6) SYSTEM WITH OPEN IMPRESSION TRAY: the driver protrudes from the impression tray so that it can be used to detach the impression coping together with the impression tray from the implant/abutment once the impression material has hardened. For a safe loosening of the impression coping, remove any residues of impression material from the driver before it hardens.
7) SYSTEM WITH CLOSED IMPRESSION TRAY: once the material has hardened, detach the sealing cap from the impression coping by gently pulling it upwards. The cap remains in the impression material.
8) Check the impression and make sure the impression coping is correctly placed. Use the impression to make the model.
9) The impression coping can be combined with both Zirkonzahn laboratory and scan analogues. Screw the desired analogue to the impression coping. Use the correct screwdriver and observe the tightening torque indicated on the label.
10) When using laboratory analogues, fill the impression with the desired casting compound. If scan analogues are used, scan the scan analogue screwed to the impression.
11) Finalising: once the casting compound has hardened, the impression coping can be detached from the impression by unscrewing the screw. When using the CLOSED TRAY system, pull down the impression first.
12) In case of digital model production, Zirkonzahn.Software can be used to create a CAM file based on the scan. At this point, the model can be milled.

Sterilisation and cleaning:
1) For sterilisation and cleaning, all components of the different sets must be treated individually. This means that all components must be separated from the impression coping.
2) Remove any visible residue with a paper towel and rinse with water.
3) Immerse the product in warm water (approx. 25 °C) for at least 10 minutes.
4) Immerse the product in a warm (approx. 25 °C) enzymatic cleaning solution. Follow the indications of the cleaning solution (e.g. 2 % for 1 minute).
5) Under running cold water (approx. 10 °C), remove all visible impurities from the product with a soft plastic brush.
6) Clean all channels with an irrigation cannula with 20 ml enzymatic cleaning solution and check for visible impurities. If there are still any, repeat this step.
7) Rinse the product under tap water until no more residues of the cleaning solution remain on the device. Clean all channels under running water and with an irrigation cannula.
8) Place the product into a cleaning/disinfection system and follow the manufacturer's instructions.
9) Dry with compressed air and a lint-free, single-use wipe.
10) Sterilise all parts before inserting them into the patient's mouth with the following procedure: steam sterilisation with vacuum at 132 °C for 4 minutes. Product sterilisation is user's responsibility.

Storage and stability:
The product is delivered in a non-sterile state and has to be stored in the original packaging in dry and dust-free conditions.

It may be possible that not all products included in these instructions for use are available in all countries.
⚠ US Federal Law restricts the sale of these products to the following professions: doctors, dentists and licensed practitioner in the dental field.

Disposal:
Medical devices which have been already used on patients have to be considered as contaminated devices and must be disposed of in accordance with local regulations. Non-contaminated products can be disposed of as industrial waste. A separate collection of reusable material (e.g. metal, plastic = packaging) is recommended.

CAUTION – Risk of injury: any product breakage may result in sharp edges and/or pointed corners. These products must be disposed of in suitable containers and in accordance with the current regulations.

For more information:
For more information, visit our website www.zirkonzahn.com.

Information:
Manufacturer: Zirkonzahn srl – Via An der Ahr 7 – 39030 Gais – Italy – www.zirkonzahn.com
Emergency Phone: +39 0474 066 660 (Mon. – Fri. 8 am – 6 pm)

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English

Date of publication: 17/03/2022

Intended use:
ZZ BASE WMSM (WHITE METAL SCANMARKER) and ZZ BASE WSM (WHITE SCANMARKER) are used to locate the exact position and orientation of the implant or abutment in the patient's mouth.
They can be used for this purpose both in the patient's mouth and on the model. In both cases, the acquired data is transferred into the software via a scanner (an intraoral scanner for scanning in the patient's mouth or a laboratory scanner – like the S300 ARTI SCANNER, S600 ARTI SCANNER or S900 ARTI SCANNER – if scanning on the model).

ZZ BASE WMSM (WHITE METAL SCANMARKER) and ZZ BASE WSM (WHITE SCANMARKER) are available for different implant systems and can be used only in combination with the specified implant system. An overview of the implant systems and components currently stored in our software is available at this webpage: www.zirkonzahn.com/implant-systems. Our virtual libraries are constantly expanded.

General information:
The information provided in this instruction manual must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel.
The user is obliged to use the products in accordance with the present instruction manual and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user is fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

Contraindications:
Should intolerance to the products be suspected, use them only after an allergy screening test has been performed with issue of medical statement of non-allergy conditions.
ZZ BASE WMSM (WHITE METAL SCANMARKER) and ZZ BASE WSM (WHITE SCANMARKER) must not be modified by the user, except for operations that are expressly specified. Non-allowed modifications result in loss of approval as medical device.

Product models:
- ZZ BASE WMSM (WHITE METAL SCANMARKER) is made of a titanium alloy with surface treatment (upper part) as well as of medical titanium alloy Ti-6Al-4V ELI (lower part and screw). The screw is firmly integrated in the scanbody. This scanbody is radiopaque and therefore easily visible on X-ray images.
- ZZ BASE WSM (WHITE SCANMARKER) is made of ceramic-reinforced polyetheretherketone (PEEK). The screw provided with this component is made of medical titanium alloy Ti-6Al-4V ELI.

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Warnings and precautions:
- ZZ BASE WMSM (WHITE METAL SCANMARKER) and ZZ BASE WSM (WHITE SCANMARKER) can only be used in combination with the intended original implant, abutment or analogue.
- During insertion, make sure that the contact surface between the scanmarker and the implant/abutment/analog is clean and free of residues and liquids.
- ZZ BASE WSM (WHITE SCANMARKER) may only be used once in the patient's mouth and only up to 15 times on the model, as the accuracy of fit may be reduced due to the material properties.
- There is no specific indication regarding the period of time for which the ZZ BASE WHITE METAL SCANMARKER work properly. This depends on how the product is utilised and maintained by the user. In any case, before each use of the products described in this instruction for use, carry out a visual inspection for defects.
- Due to the small size of the products, when using ZZ BASE WMSM (WHITE METAL SCANMARKER)/ZZ BASE WSM (WHITE SCANMARKER) special precautions must be taken to avoid ingestion or inhalation of the parts by the patient.
- Screws may only be slightly tightened observing torque indications.
- ZZ BASE WMSM (WHITE METAL SCANMARKER)/ZZ BASE WSM (WHITE SCANMARKER) must not be used if damaged or altered.
- Do not disassemble the ZZ BASE WMSM (WHITE METAL SCANMARKER).
- Do not rotate the ZZ BASE WMSM (WHITE METAL SCANMARKER)/ZZ BASE WSM (WHITE SCANMARKER) when fixed onto the implant/abutment.
- ZZ BASE WMSM (WHITE METAL SCANMARKER) and ZZ BASE WSM (WHITE SCANMARKER) must not be used until they have cooled down after autoclaving. Otherwise, inaccurate scan results may occur due to thermal expansion.
- The products packaging cannot be autoclaved.
- The screws and scanbodies of the individual systems are very similar and can therefore be confused easily. For this reason, it is advisable to assemble the set parts correctly and sterilise them together.

Application in the patient's mouth:
1) Before first use and after each subsequent use on the patient, the ZZ BASE WMSM (WHITE METAL SCANMARKER)/ZZ BASE WSM (WHITE SCANMARKER) and the corresponding screw must be cleaned, dried and sterilised. Immediately before the next use, reassemble the components and proceed with the sterilisation (see „Sterilisation and cleaning").
2) Place the ZZ BASE WMSM (WHITE METAL SCANMARKER)/ZZ BASE WSM (WHITE SCANMARKER) onto the implant or abutment in the patient's mouth. For optimal scan results, align the marking on the ZZ BASE WSM (WHITE SCANMARKER) to the mesio-labial side. There is no specific orientation to consider for scanning the ZZ BASE WMSM (WHITE METAL SCANMARKER).
3) Tighten the screw slightly (observe torque values).
4) Scan with an intraoral scanner (observe the manufacturer's instructions).
5) Remove the ZZ BASE WMSM (WHITE METAL SCANMARKER)/ZZ BASE WSM (WHITE SCANMARKER).

Application on the model:
1) Place the ZZ BASE WMSM (WHITE METAL SCANMARKER)/ZZ BASE WSM (WHITE SCANMARKER) on the model and screw it (observe torque values). For optimal scan results, align the marking on the ZZ BASE WSM (WHITE SCANMARKER) to the mesio-labial side. There is no specific orientation to consider for scanning the ZZ BASE WMSM (WHITE METAL SCANMARKER).
2) Carry out the scanning process with the S300 ARTI SCANNER, S600 ARTI SCANNER or S900 ARTI SCANNER.
3) Remove the ZZ BASE WMSM (WHITE METAL SCANMARKER)/ZZ BASE WSM (WHITE SCANMARKER).

Tightening torque screws:
ZZ BASE WMSM (WHITE METAL SCANMARKER) and ZZ BASE WSM (WHITE SCANMARKER) may only be tightened by hand to avoid deformations, which might lead to reduced fit accuracy. Tightening torques are indicated on the packaging or can be found at http://www.zirkonzahn.com/assets/files/brochueren/INT-Drehmomenttabelle-fuer-Zirkonzahn-schrauben.pdf.

Sterilisation and cleaning:
1) Remove any visible residue with a paper towel and rinse with water.
2) Immerse the product in warm water (approx. 25 °C) for a minimum of 10 minutes.
3) Immerse the product in a warm (approx. 25 °C) enzymatic cleaning solution. Follow the indications of the cleaning solution (e.g. 2 % for 1 minute).
4) Remove all visible impurities from the product with a soft plastic brush under cold running water (approx. 10 °C).
5) Clean lumina and channels with an irrigation canula with 20 ml enzymatic cleaning solution and check for visible impurities. In this case, repeat this step.
6) Rinse the product under tap water until no residues of the cleaning solution remain on the device. Clean all lumina and channels with running water and an irrigation canula.
7) Load the product into a cleaning/disinfection system and follow the manufacturer's instructions.
8) Dry with compressed air and lint-free, single-use wipes.
9) Sterilise all parts before inserting them into the patient's mouth with the following procedure: steam sterilisation with vacuum; at 132 °C for 4 minutes. The product sterilisation is the user's responsibility.

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US Federal Law restricts these products to sale by or on the order of a licensed practitioner in the dental sector.

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