Gebrauchsanweisung – Istruzioni d'uso
Instruction manual – Mode d'emploi
Instrucciones de uso

Zirkonzahn®

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Zirkonzahn®

ZIRCONIA CLEAR 2, PRETTAU® COLOURED, PRETTAU® DISPERSIVE®, PRETTAU® 2 COLOURED, ANATOMIC POSTERIOR PAINTABLE, ANATOMIC POSTERIOR COLOURED, ANATOMIC POSTERIOR DISPERSIVE®, ZIRCONIA CLEAR 4, ICE PREMIUM, ICE ABUTMENT, ICE TRANSLUCENT PLUS COLOURED, ICE TRANSLUCENT PLUS DISPERSIVE®, ICE ECO LINE, ICE PLUS, ANATOMIC PAINTABLE, ANATOMIC COLOURED, ANATOMIC DISPERSIVE®, ZIRCONIA CLEAR 1, ZIRCONIA CLEAR 5, Z-WHITE, ZIRCONIA CLEAR 6, PRETTAU® 3 COLOURED, PRETTAU® 4 ANTERIOR® COLOURED, ANATOMIC CLEAR PAINTABLE, ANATOMIC CLEAR COLOURED, ANATOMIC CLEAR DISPERSIVE®, ZIRCONIA CLEAR 3

CE 0476

GS0116V01-220328

Zirkonzahn srl
Via An der Ahr 7
39030 Gais – ITALY
Tel +39 0474 066 660 – Fax +39 0474 066 661
E-Mail: info@zirkonzahn.com

CH REP Jan Möstel
Robert-Seidel-Hof 70, 8048 Zürich

English

Date of update: 28/03/2022

Intended use:
Zirconia Clear 2, Prettau® Coloured, Prettau® Dispersive®, Prettau® 2 Coloured, Anatomic Posterior Paintable, Anatomic Posterior Coloured, Anatomic Posterior Dispersive®, Zirconia Clear 4, ICE Premium, ICE Abutment, ICE Translucent Plus Coloured, ICE Translucent Plus Dispersive®, ICE Eco Line, ICE Plus, Anatomic Paintable, Anatomic Coloured, Anatomic Dispersive®, Zirconia Clear 1, Zirconia Clear 5, Z-White and Zirconia Clear 6 are intended for the manufacturing of metal-free partial and single crowns, full arch occlusally screwed bridges, inlays, onlays, and veneers as full contour restorations as well as for reduced structures in combination with veneering ceramics. The products have to be categorised as class 5 according to ISO 6872.
Prettau® 3 Coloured, Prettau® 4 Anterior® Coloured, Anatomic Clear Paintable, Anatomic Clear Coloured, Anatomic Clear Dispersive® and Zirconia Clear 3 are destined for the manufacturing of metal-free partial and single crowns, max. 3-unit bridges, inlays, onlays and veneers, as full contour restoration as well as for reduced structures in combination with veneering ceramics and implant superstructures for 3-unit bridges in the anterior and posterior tooth region. The products have to be categorised as class 4 according to ISO 6872.
The products have been developed for use with Colour Liquids, ICE Ceramics, ICE Stains as well as ICE Stains 3D by Enrico Steger . Observe the relative instructions of use when using these products.
The blocks are suitable for all milling units, which are able to process presintered zirconia and which have the proper clamping device for the corresponding block.

⚠ Attention: The pre-sintered structure must never come into contact with water or steam.

General information:
The information provided in this instruction manual must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel. The user is obliged to use the products in accordance with the present instruction manual and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user will be held fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

Contraindications:
Should intolerance to the products be suspected, use them only after an allergy screening test has been performed with issue of medical statement of non-allergy conditions.

Work steps:
1) Fix the zirconia block in the specific clamping device.
2) Start the milling process (the shrinkage factor is indicated on the packaging). The blocks are milled dry and by using a suction unit.
3) Check milled framework for obvious chips or cracks. Only frameworks in perfect condition must be used for dental restorations. In case of unsatisfactory results, the milling burs have to be checked and, if necessary, replaced; the milling unit needs to be cleaned and recalibrated.
4) Separate the structure with a diamond bur from the blank (dry and using a suction unit).
5) Finish structure with tungsten carbide burs.
6) Air blow the structure.
7) Colour structure with the corresponding Colour Liquids (see manual for Colour Liquids).
8) Drying time: min. 45 minutes under an infrared lamp.
9) Place the structure in the sintering furnace suitable for sintering zirconia, close the furnace door and select programme (observe chart of sintering parameters, if necessary the sinter programmes must be updated). Bridge restorations must be sufficiently stabilised to avoid warpage during the sintering process in the furnace. The sintering support in the form of a so-called sintering base should be provided starting from a size of 3 elements. Zirkonzahn has secured the necessary licence for its clients until the end of the patent term.
10) Remove sintered structures from the furnace at room temperature.
11) Check structure on working model for seat and margin fit. Trim with diamond turbine on water cooled turbine if required.
12) Sandblast structure using 120 micron aluminium oxide and a pressure of 4 bar.
13) Clean structure with water, steam-clean.
14) Structure is now ready for veneering.
⚠ Attention: Observe CTE value and follow instructions of the respective ceramics producers!
15) A translucent cement should be used for cementation.
Non-translucent cement could shimmer through the primary structure.
⚠ Attention: Follow the instructions of the respective cement producer.

Dimensioning of connectors:
The products of class 4 according to ISO 6872 are intended for bridges with a max. span of up to 3 units. When manufacturing bridges, the dimensions of the connectors need to be taken into account. It should be 3 x 4 mm in the anterior region and 4 x 4 mm in the posterior region. For all other zirconia types, the dimensions of the connectors can be calculated using the calculation tool on our website (www.zirkonzahn.com/en/calculation-tool).

Chart final temperatures and CTE values:

Material	CTE value	Final sintering temperature
Zirconia Clear 2 Prettau® 2 Coloured Anatomic Posterior Paintable Anatomic Posterior Coloured Anatomic Posterior Dispersive®	~ 10,4*10-6 K-1	1600 °C
Zirconia Clear 4 Prettau® Coloured Prettau® Dispersive®	~ 10,6*10-6 K-1	1600 °C
Prettau® 3 Coloured	~ 10,3*10-6 K-1	
Prettau® 4 Anterior® Coloured Anatomic Clear Paintable Anatomic Clear Coloured Anatomic Clear Dispersive®	~ 10,1*10-6 K-1	
Zirconia Clear 3 ICE Eco Line Zirconia Clear 1 ICE Premium ICE Abutment Z-White Zirconia Clear 6	~ 10,6*10-6 K-1	1500 °C
ICE Translucent Plus Coloured ICE Translucent Plus Dispersive® Anatomic Paintable Anatomic Dispersive® Zirconia Clear 5	~ 10,7*10-6 K-1	
ICE Plus	~ 9,8*10-6 K-1	1500 °C

Chart of sintering parameters – heating rates (valid for all materials)

Weight (incl. sintering stabiliser/dental unit)	Heating rate	Holding time
0 – 2 g	max. 14 °C/min.	1.5 h
2 – 4 g	max. 8 °C/min.	2 h
> 4 g	max. 5 °C/min.	2 h

For the veneering of multi-unit structures you must observe the following parameters during the heating and the cooling cycles:

0 – 1 g	per dental unit	55 °C/min.
1 – 2 g	per dental unit	45 °C/min.
2 – 3 g	per dental unit	35 °C/min.
over 3 g	per dental unit	25 °C/min.

Safety instructions:
Please observe product-specific safety data sheet!

Storage and stability:
Store in a dry place at room temperature. Do not expose to heavy shocks or vibrations. Blocks must not come into contact with liquids. Protect the zirconia discs from dirt.

Avoid any exposure to humidity. Avoid the formation and deposition of dust.

⚠ US Federal Law restricts these products to sale by or on the order of a licensed practitioner in the dental sector.

Information:
Manufacturer: Zirkonzahn srl – Via An der Ahr 7 – 39030 Gais – Italy – www.zirkonzahn.com
Emergency Phone: +39 0474 066 660 (Mon. – Fri. 8 am – 6 pm)

N.B.: All medical devices and all products with open packaging are non-returnable.

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Zirkonzahn®

ZIRCONIA CLEAR 2, PRETTAU® COLOURED, PRETTAU® DISPERSIVE® , PRETTAU® 2 COLOURED, ANATOMIC POSTERIOR PAINTABLE, ANATOMIC POSTERIOR COLOURED, ANATOMIC POSTERIOR DISPERSIVE®, ZIRCONIA CLEAR 4, ICE PREMIUM, ICE ABUTMENT, ICE TRANSLUCENT PLUS COLOURED, ICE TRANSLUCENT PLUS DISPERSIVE®, ICE ECO LINE, ICE PLUS, ANATOMIC PAINTABLE, ANATOMIC COLOURED, ANATOMIC DISPERSIVE®, ZIRCONIA CLEAR 1, ZIRCONIA CLEAR 5, Z-WHITE, ZIRCONIA CLEAR 6, PRETTAU® 3 COLOURED, PRETTAU® 4 ANTERIOR® COLOURED, ANATOMIC CLEAR PAINTABLE, ANATOMIC CLEAR COLOURED, ANATOMIC CLEAR DISPERSIVE®, ZIRCONIA CLEAR 3

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Zirkonzahn®

PRETTAU®, ANATOMIC COLOURED, PRETTAU® 2, PRETTAU® 2 DISPERSIVE®, PRETTAU® 3, PRETTAU® 3 DISPERSIVE®, PRETTAU® 4 ANTERIOR®, PRETTAU® 4 ANTERIOR® DISPERSIVE®, ICE TRANSLUCENT, ICE TRANSLUCENT PLUS

CE 0476

GS0001V01-2203

Zirkonzahn srl
Via An der Ahr 7
39030 Gais – ITALY
Tel +39 0474 066 660 – Fax +39 0474 066 661
E-Mail: info@zirkonzahn.com

CH REP Jan Möstel
Robert-Seidel-Hof 70, 8048 Zürich

English

Date of update: 17/03/2022

Intended use:
Prettau®, Prettau® 2, Prettau® 2 Dispersive®, ICE Translucent, ICE Translucent Plus and Anatomic Coloured are intended for the manufacturing of metal-free partial and single crowns, full arch occlusally screwed bridges, inlays, onlays, and veneers as full contour restorations as well as for reduced structures in combination with veneering ceramics. The products have to be categorised as class 5 according to ISO 6872.
Prettau® 3, Prettau® 3 Dispersive®, Prettau® 4 Anterior® and Prettau® 4 Anterior® Dispersive® are destined for the manufacturing of metal-free partial and single crowns, max. 3-unit bridges, inlays, onlays and veneers, as full contour restoration as well as for reduced structures in combination with veneering ceramics and implant superstructures for 3-unit bridge restorations in the anterior and posterior tooth region. The products have to be categorised as class 4 according to ISO 6872.
The products have been developed for use with Colour Liquids, ICE Ceramics, ICE Stains and ICE Stains 3D. Observe the relative instructions of use when using these products.
The blocks are suitable for all milling units, which are able to process presintered zirconia and which have the proper clamping device for the corresponding block.

⚠ Attention: the pre-sintered structure must never come into contact with water or steam.

General information:
The information provided in this instruction manual must be passed on to any person using the products mentioned therein. The products must only be used by qualified personnel.
The user is obliged to use the products in accordance with the present instruction manual and with appropriate hygiene measures and to verify on his/her own responsibility whether the products are suitable for the individual patient situation. The user will be held fully responsible for the appropriate and correct use of the products. Zirkonzahn assumes no liability for incorrect results in form of direct or indirect damages or any other damages that occur from the use and/or the processing of the products. Any claim for damages (including punitive damages), is limited to the commercial value of the Zirkonzahn products. Independently of this, the user is obliged to report all serious incidents that occur in connection with the products to the competent authority and to Zirkonzahn.

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Work steps:
1) Fix the zirconia block in the specific clamping device.
2) Start the milling process (the shrinkage factor is indicated on the packaging). The blocks are milled dry and by using a suction unit.
3) Check milled framework for obvious chips or cracks. Only frameworks in perfect condition must be used for dental restorations. In case of unsatisfactory results, the milling burs have to be checked and, if necessary, replaced; the milling unit needs to be cleaned and recalibrated.
4) Separate the structure with a diamond bur from the blank (dry and using a suction unit).
5) Finish structure with tungsten carbide burs.
6) Air blow the structure.
7) Colour structure with the corresponding Colour Liquids (see manual for Colour Liquids).
8) Drying time: min. 45 minutes under an infrared lamp.
9) Place the structure in the sintering furnace suitable for sintering zirconia, close the furnace door and select programme (observe chart of sintering parameters, if necessary the sinter programmes must be updated). Bridge restorations must be sufficiently stabilised to avoid warpage during the sintering process in the furnace. The sintering support in the form of a so-called sintering base should be provided starting from a size of 3 elements. Zirkonzahn has secured the necessary licence for its clients until the end of the patent term.
10) Remove sintered structures from the furnace at room temperature.
11) Check structure on working model for seat and margin fit. Trim with diamond turbine on water cooled turbine if required.
12) Sandblast structure using 120 micron aluminium oxide and a pressure of 4 bar.
13) Clean structure with water, steam-clean.
14) Structure is now ready for veneering.
⚠ Attention: observe CTE value and follow instructions of the respective ceramics producers!
15) A translucent cement should be used for cementation.
Non-translucent cement could shimmer through the primary structure.
⚠ Attention: follow the instructions of the respective cement producer.

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For the veneering of bridges of various elements you must observe the following parameters during the heating and the cooling cycles:

0 – 1 g	per dental unit	55 °C/min.
1 – 2 g	per dental unit	45 °C/min.
2 – 3 g	per dental unit	35 °C/min.
über 3 g	per dental unit	25 °C/min.

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