TERUMOBCT

## 10763, 10790

Mirasol® Platelet Disposable Kit for Treatment in Platelet Additive Solution
Jednorázová souprava Mirasol® pro ošetření trombocytů v náhradním roztoku trombocytů
Mirasol® engangskit til trombocytter til behandling i trombocytopbevaringsvæske
Mirasol®-Thrombozyten-Einwegkit zur Behandlung von Thrombozyten in Thrombozyten-Lagerlösung
Αναλώσιμο κιτ αιμοπεταλίων Mirasol® για επεξεργασία σε προσθετικό διάλυμα αιμοπεταλίων
Kit desechable de plaquetas de Mirasol® para el tratamiento en solución aditiva para plaquetas
Dispositif à usage unique Mirasol® Plaquettes pour traitement dans une solution additive plaquettaire
Kit monouso Mirasol® per il trattamento delle piastrine con soluzione additiva per piastrine
Disposable Mirasol® bloedplaatjeskit voor behandeling in bloedplaatjesbewaarvloeistof
Mirasol® engangssett for blodplater til behandling i PAS (blodplatetilsetning)
Jednorazowy zestaw Mirasol® do obróbki płytek krwi w roztworach wzbogacających do płytek (PAS)
Kit descartável para plaquetas Mirasol® para tratamento em solução aditiva de plaquetas
Одноразовый комплект Mirasol® для обработки тромбоцитов в растворе для хранения тромбоцитов
Mirasol® komplet za jednokratnu upotrebu za obradu trombocita u aditivnom rastvoru za trombocite
Mirasol®-engångsset för behandling av trombocyter i tillsatslösning för trombocyter
Trombosit Ek Solüsyonu'nda İşleme İçin Mirasol® Trombosit Tek Kullanımlık Seti

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## Description

The Mirasol® Pathogen Reduction Technology (PRT) System consists of the Mirasol Illuminator, the Mirasol Manager software (optional), and a Mirasol disposable kit. The treatment procedure involves mixing a blood product with Riboflavin Solution in the disposable kit and illuminating the mixture with a calculated light energy dose in the Illuminator.

The Mirasol Platelet Disposable Kits for Treatment in Platelet Additive Solution (PAS) are for use with apheresis Platelets and whole blood-derived Platelets. These sterile disposable kits contain an Illumination/Storage Bag and Riboflavin Solution. The Mirasol DOUBLE Platelet Disposable Kit for treatment in PAS also contains a pre-attached second Platelet Storage Bag designed for splitting products without connecting a separate storage bag. After optional splitting and sampling, Mirasol-treated Platelet products are ready for storage or transfusion.

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## Disposable Kit Overview

Both Mirasol Platelet Disposable Kits for Treatment in PAS consist of:


- Clear outer packaging with a white peel strip for opening.
- One opaque foil pouch containing one bag of Riboflavin Solution (500 µM riboflavin in normal saline) (see Figure 1 and Figure 2).
- One Illumination/Storage Bag with attached tubing, frangible connectors, spike port, and sample bulb (see Figure 1 and Figure 2).

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
The Mirasol DOUBLE Platelet Disposable Kit also contains:

- A Sample Bag
- A second Platelet Storage Bag with attached spike port and sample bulb (see Figure 2).

1 Inlet luer with blue end cap
2 Inlet lines
3 Spike port
4 "3S" set type label
5 Opaque foil pouch
6 Riboflavin Solution Bag
7 Sample bulb
8 Frangible connectors
9 Pre-attached label
10 Illumination/Storage Bag

**Figure 1: Mirasol Platelet Disposable Kit for Treatment in PAS ("single-bag kit")**

1 Inlet luer with blue end cap
2 Inlet lines
3 "3D" set type label
4 Opaque foil pouch
5 Riboflavin Solution Bag
6 Frangible connectors
7 Sample Bag
8 Sample Bag luer with blue end cap
9 Spike ports
10 Illumination/Storage Bag
11 Second Platelet Storage Bag
12 Sample bulb
13 Pre-attached labels

**Figure 2: Mirasol DOUBLE Platelet Disposable Kit for Treatment in PAS**

## Warnings

1. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

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Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
  - Fluid leaks
  - Parts that are warped or deformed
  - Plastics that are brittle and discolored
  - Filters that have reduced filtration capabilities
- Exposure to excessive ethylene oxide (EO) residuals
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

## Cautions

### General

1. Each operator should be thoroughly familiar with these Instructions for Use and the Mirasol Illuminator Operator's Manual.
2. Use aseptic technique throughout all procedures and follow standard blood establishment practices.

### Disposable Kit

1. Use the Mirasol Platelet Disposable Kits for Treatment in PAS described here only for treatment of Platelets in PAS. Different Mirasol Disposable Kits are available for treatment of Platelets in Plasma.
2. Use only Disposable Kits manufactured for the Mirasol PRT system.
3. Do not use the Disposable Kit if it has passed the printed expiration date (YYYY-MM-DD).
4. Carefully inspect for leaks when using the Disposable Kit. If a leak is observed, do not use the Platelet product, or manage the product according to your institution's standard operating procedure (SOP).
5. The Illumination/Storage Set is sterilized with ethylene oxide (EO), which in rare cases has been reported to cause recipient reactions.
6. The light energy dose delivered to the Platelet product is calculated based on the Illumination/Storage Set tare weight configured on the Illuminator, among other factors. Therefore, the configured tare weight must be accurate. Do not remove components from the Illumination/Storage Set (except the Illumination Bag inlet line) prior to illumination unless the correct tare weight has been configured on the Illuminator. See the Mirasol Illuminator operator's manual for instructions on setting tare weight. If the correct tare weight cannot be configured, contact Terumo BCT, Inc.
7. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies

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have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.

8. The Mirasol Illumination Bag has not been qualified for centrifugation.

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#### Blood Components

1. Do not agitate Platelet products during the minimum rest period listed in Table 2.
2. All Platelet products must be illuminated according to the treatment time constraints listed in Table 2. This is to ensure adequate Platelet cell quality after treatment and storage, and to minimize the risk for potential endotoxin build-up to levels that may be harmful to the recipient.
3. After adding Riboflavin Solution and until the time of transfusion, protect the Platelet product from prolonged exposure to direct sunlight or strong ambient light sources. See Table 4 for storage specifications.
4. If the physician transfusing the Platelet product prescribes gamma irradiation, gamma irradiate the Platelet product only after treating with the Mirasol system.
5. If the product contains any platelet aggregates, do not continue with Mirasol system treatment until the platelet aggregates have dissipated.

## Notes

Only Platelet products collected in ACD-A, CPD, or CPDA-1 anticoagulants have been evaluated for use in the Mirasol System.

## Instructions for Use

#### Treating the Platelet Product with the Mirasol System

**Caution:** Use the Mirasol Platelet Disposable Kits described here only for treatment of Platelets in PAS. Different Mirasol Disposable Kits are available for treatment of Platelets in Plasma.

**Note:** The following steps and figures refer to **both** Platelet Disposable Kits for Treatment in platelet additive solution. Unique considerations for the DOUBLE Platelet Disposable Kit are noted.

1. Review the treatment time constraints listed in Table 2 and the storage volume, concentration, and other conditions listed in Table 4.
2. If treating an apheresis Platelet product, allow it to rest according to the rest periods listed in Table 2.

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3. Just prior to use, open the white peel strip on the outer packaging of the Disposable Kit. Remove the Illumination/Storage Set from the white spacer and inspect the set.

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
a. The blue end cap on either the inlet line or the sample bag is unattached.
b. The integrity of the Illumination/Storage Set is compromised for any reason.
c. The frangible connectors are broken.
d. The tubing is severely kinked.
e. There is debris on the outer surface of the Illumination/Storage Bag. Do not clean with disinfection chemicals.
f. The barcodes are illegible.
g. Any pre-attached label is missing or not securely attached.

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4. If you are using the DOUBLE set, permanently seal and disconnect the luer line on the sample bag (see Figure 2, item 8).
5. Sterile connect the Platelet Product Bag to the part of the Illumination/Storage Bag inlet line with the luer and blue end cap, which results in removal of the line with the luer and blue end cap. See Figure 3.

1 Platelet Product Bag
2 Inlet luer with blue end cap
3 Sterile connection on inlet line
4 Illumination/Storage Bag
5 Pre-attached label

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6. Appropriately label the Illumination/Storage Bag according to your facility's standard operating procedure (SOP), using space on the pre-attached label. For the DOUBLE set, label the second Platelet Storage Bag, if required by your facility's SOP.

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7. Just before you are ready to transfer the Riboflavin Solution to the Illumination/Storage Bag, open the foil pouch, remove the Riboflavin Solution Bag, and inspect the bag and solution.

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
8. You can either transfer the Riboflavin Solution and the Platelet product to the Illumination/Storage Bag at the same time, or you can transfer them one at a time.
   - To transfer the Platelet product and the Riboflavin Solution at the same time:
     a. Sterile connect the Riboflavin Solution Bag to the second inlet line on the Illumination/Storage Bag. See Figure 4.
     b. Open both sterile connections and drain the Platelet product and the Riboflavin Solution into the Illumination/Storage Bag.
     c. Permanently seal the inlet line between the Platelet Product Bag and the Y connector, leaving the Riboflavin Solution Bag attached.
   - To transfer the Platelet product and the Riboflavin Solution one at a time:
     a. Open the sterile connection and drain the Platelet product into the Illumination/Storage Bag.
     b. Allowing enough tubing length for connecting to the Riboflavin Solution Bag, permanently seal the inlet line between the Platelet Product Bag and the Y connector.
     c. Sterile connect the Riboflavin Solution Bag to the same inlet line you used to transfer the Platelet product.
     d. Open the sterile connection and drain the Riboflavin Solution into the Illumination/Storage Bag.

1 Platelet Product Bag
2 Riboflavin Solution Bag
3 Sterile connections on inlet lines
4 Y connector
5 Pinch clamp
6 Illumination/Storage Bag

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9. To express residual air from the Illumination/Storage Bag:
   a. Hold the Illumination/Storage Bag so that the ports are facing up, as shown in Figure 4.
   b. Gently apply pressure to the Illumination/Storage Bag to move the Platelets up through the inlet line tubing until they just reach the Riboflavin Solution Bag.
   c. Drain the Platelets from the Riboflavin Solution Bag and inlet line into the Illumination/Storage Bag, closing the clamp on the inlet line before the air interface reaches the Illumination/Storage Bag port.

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
**10.** Permanently seal the inlet line as near to the Illumination/Storage Bag as possible. The sealed inlet line must be less than 2.5 cm (1 inch) long, as shown in Figure 6, for it to properly fit in the Illuminator. Disconnect and discard the Riboflavin Solution Bag.

**1** Sample bulbs
**2** Spike ports
**3** Sealed inlet lines
**4** Frangible connectors
**5** Illumination/Storage Bags
**6** Pinch clamps
**7** Sample Bag
**8** Y connector
**9** Second Platelet Storage Bag

**Figure 6: Single-bag (a) and DOUBLE (b) Illumination/Storage Sets ready for illumination**

**11.** If the Mirasol Illuminator has been configured with a tare weight for a DOUBLE set with some components (the Sample Bag, second Platelet Storage Bag, or both) removed, permanently seal the lines to these components and remove them now.

**12.** Perform the illumination procedure according to the Mirasol Illuminator operator's manual.

**13.** After performing the illumination procedure, break the frangible connector on the spike port. You can break the second frangible connector if access is needed to split the product or to take a sample (see “Sampling the Mirasol-Treated Platelet Product” on page 12).

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**14.** Determine whether you will split the Mirasol-treated Platelet product.

- If you will split the product, see “Splitting the Mirasol-Treated Platelet Product” on page 12.
- If you will not split the product, but are using the DOUBLE set, permanently seal the line between the Y connector and the second Platelet Storage Bag. Remove the second Platelet Storage Bag, leaving the Sample Bag.

**15.** If necessary, take a sample of the Mirasol-treated Platelet product. See “Sampling the Mirasol-Treated Platelet Product” on page 12.

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16. Establish appropriate storage conditions according to Table 4.
    - Until the time of transfusion, protect the treated Platelet product from prolonged exposure to direct sunlight or strong ambient light.
    - The laboratory/blood center should follow applicable standards for product quality acceptance criteria.

## Splitting the Mirasol-Treated Platelet Product

Table 4 shows recommended storage conditions. If required, split the treated product prior to storage.

To split the product using the single-bag set, connect an accessory Platelet Storage Bag to the Illumination/Storage Bag sample bulb outlet line. An accessory Platelet Storage Bag (Catalog #70030) is available from the manufacturer.

To split the product using the DOUBLE set:

1. Close the pinch clamps on the lines to the Sample Bag and the sample bulb (see Figure 6).
2. Break the frangible connector on the outlet port of the Illumination/Storage Bag.
3. Split the Platelet product between the Illumination/Storage Bag and the second Platelet Storage Bag.
4. If desired, express residual air from the Illumination/Storage Set into the Sample Bag.
5. Close the pinch clamp between the Y connector and the second Platelet Storage Bag.
6. If desired, collect a product sample in the Sample Bag. See "Sampling the Mirasol-Treated Platelet Product".
7. Permanently seal the line between the Y connector and the Sample Bag and remove the Sample Bag.

## Sampling the Mirasol-Treated Platelet Product

You can collect a product sample using the Sample Bag on the DOUBLE set, the sample bulb on the single-bag set or the DOUBLE set, or a separate sampling device. See the Mirasol Product Sampling Guidance document for detailed sampling guidance.

**Note:** If you expressed air into the Sample Bag while splitting the product, you may have allowed platelets into the Sample Bag line. In this case, you should use the Sample Bag to take a sample only **before** storage, because the viability of the platelets in the Sample Bag line will degrade during storage.

To collect a product sample in the Sample Bag, complete the following steps:

1. Make sure the pinch clamp between the Y connector and the second Platelet Storage Bag is closed.
2. Break the frangible connector (if it is not already broken) on the outlet line of the Illumination/ Storage Bag.
3. Open the Sample Bag pinch clamp (if it is not already open) and gently apply pressure to the Illumination/Storage Bag to move Platelets up through the line until the desired amount of the mixture is in the Sample Bag.
4. Close the Sample Bag pinch clamp.
5. Permanently seal the line between the Y connector and the Sample Bag and remove the Sample Bag.

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

## Symbols and Certification

### Table 1: Symbols and Certification

| Symbol | Description |
|---|---|
|  | Indicates the package contents. |
|   | Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste. |
|  | Indicates a non-pyrogenic fluid pathway. |
|  | Indicates that the product should not be used if the package is damaged. |
|  | Indicates that the product should not be used if the package is damaged. |
|  | Indicates that the product should not be vented. Do not remove the end caps from the set. |
|   | Indicates the presence of a sterile fluid path. The method of sterilization is ethylene oxide (EO). |
|   | Indicates the presence of a sterile fluid path. The method of sterilization is steam. |
|  | Indicates the medical device manufacturer. |
|   | Indicates the Authorized representative in the European Community. |
|  | Indicates the date after which the medical device is not to be used. |
|   | Indicates the manufacturer's catalog number so that the medical device can be identified. |
|   | Indicates the manufacturer's batch code so that the batch or lot can be identified. |
|  | Indicates a medical device that is not to be resterilized. |
|  | Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. |

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| Symbol | Description |
|---|---|
|  | Indicates that the user should consult the instructions for use. |
|   | Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended. |
|   | Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP). |
|  | Indicates a medical device that needs protection from light sources. |
|  | Indicates a medical device that needs to be protected from moisture. |
|  | Indicates the range of temperatures to which the device can be safely exposed. |
|  | Indicates the product quantity when the quantity is placed in the square. |
|  | Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile). |

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| Characteristic | Apheresis Platelets | Whole-Blood-Derived Platelets |
|---|---|---|
| Anticoagulant | ACD-A | CPD and CPDA-1 |
| Use of platelet additive solution (PAS) | Use only approved platelet additive solutions. Solutions containing phosphate buffer are recommended for storage of Mirasol-treated platelet products. | |
| Platelet volume including PAS (prior to riboflavin solution addition) | 250 mL to 450 mL | |
| Platelet yield | ≥ $2.4 \times 10^{11}$ platelets | |
| Platelet concentration (prior to riboflavin solution addition) | 0.8 to $1.5 \times 10^{6}$ platelets/µL | |
| % Plasma (including anticoagulant) | 30% to 45% | |
| Treatment time constraints | • Rest (no agitation) for a minimum of 2 hours after collection and before treatment.<br>• Treat within 2 to 22 hours of collection. | • Treat within 8 hours of preparing the final (pooled) whole-blood-derived platelet concentrate, but not more than 32 hours after whole blood collection. This time limit assumes that leukoreduction occurs during or after the pooling step. |
| Residual white blood cells | ≤ $250 \times 10^{6}$ per unit of $3.5 \times 10^{11}$ yield | |
| Gamma irradiation | If desired, gamma irradiation may be performed after Mirasol PRT system treatment at any point during the storage period.<br>**Note:** Mirasol PRT system treatment can be used instead of gamma irradiation to prevent transfusion-associated graft-versus-host disease (TA-GVHD). | |
| *The Mirasol PRT system for treating Platelets in PAS has been verified for these incoming product conditions. Certain conditions outside of these ranges may be acceptable but should be evaluated by the user as required. | | |

**Table 3: Disposable Kits***

| Characteristic | Specification: Single-Bag Illumination/Storage Set | Specification: DOUBLE Illumination/Storage Set |
|---|---|---|
| Mirasol Platelet Disposable Kit | Each Disposable Kit consists of one Illumination/Storage Set and one 35 mL bag of 500 µM riboflavin in 0.9% sodium chloride solution. | |
| Illumination/Storage Bag | Citrate vinyl bag | |

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**Table 3: Disposable Kits* (continued)**

| Characteristic | Specification: Single-Bag Illumination/Storage Set | Specification: DOUBLE Illumination/Storage Set |
|---|---|---|
| Second Platelet Storage Bag | N/A | Citrate vinyl bag |
| Illumination/Storage Bag sterility | The Illumination/Storage Set has a sterile and non-pyrogenic fluid pathway, sterilized by ethylene oxide (EO). | |
| Illumination/Storage Bag label | The pre-attached label covers 110 $cm^2$ (17 $in^2$) surface area.<br>This size label must be present during illumination. Additional labels may be applied over the pre-attached label, but they may not cover any bag surface area outside the edges of the pre-attached label. | |
| Tubing dimensions for sterile connections and seals | Outer diameter = 4.06 mm ± 0.08 mm<br>Wall thickness = 0.56 mm ± 0.05 mm | |
| Illumination/Storage Bag and second Platelet Storage Bag unrecoverable volume | 3.4 mL | |
| Sample Bag volume | N/A | 50 mL |
| Sample bulb volume | 1.5 to 4 mL | |
| Riboflavin Solution content | 500 µM riboflavin in 0.9% sodium chloride solution<br>pH 3.5 to 5.1<br>Riboflavin Material Safety Data Sheets are available upon request. | |
| Riboflavin Solution sterility | The Riboflavin Solution is sterile and non-pyrogenic, sterilized by steam. | |
| Riboflavin Solution volume | 35 mL ± 5 mL | |
| Riboflavin Solution raw materials | Riboflavin: USP, EP<br>Sodium chloride: USP, EP<br>Hydrochloric acid: USP, EP<br>Water for injection: USP, EP | |
| Riboflavin Solution time constraints | Treat the Platelet product within 1 hour of removing the Riboflavin Solution Bag from the foil pouch. | |
| Disposable kit storage conditions | Long-term storage range: 2 °C to 30 °C (36 °F to 86 °F)<br>Storage humidity: Up to 75% relative humidity (RH), non-condensing<br>Permitted excursions: Up to 60 °C (140 °F) for up to 2 weeks | |
| Case quantity | 25 kits per case | |
| Case size | 24.9 cm × 35.8 cm × 35.8 cm | 24.9 cm × 43.5 cm × 35.8 cm |

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**Table 3: Disposable Kits* (continued)**

| Characteristic | Specification: Single-Bag Illumination/Storage Set | Specification: DOUBLE Illumination/Storage Set |
|---|---|---|
| Case weight | 5.0 kg | 7.0 kg |

**Table 4: Platelet Product Storage Specifications***

| Characteristic | Storage Conditions Per Bag |
|---|---|
| Storage volume | 150 mL to 435 mL** |
| Storage concentration | 0.70 to 1.5 × $10^6$ platelets/µL |
| Yield | < 5.1 × $10^{11}$ platelets** |
| Storage period and environment | • After adding Riboflavin Solution and until the time of transfusion, protect the Platelet product from prolonged exposure to direct sunlight or strong ambient light sources.<br>• To ensure sufficient protection from light, store the Platelet product in an enclosed Platelet storage cabinet or on any intermediate shelf of an open tabletop Platelet agitator.<br>• **Do not** store the Platelet product on the top shelf of an open tabletop Platelet agitator where it will be exposed to direct ambient light or sunlight.<br>• Store in PAS for a maximum of 7 days after collection at 22 °C ± 2 °C with gentle agitation.<br>**Note:**<br>Any extension of platelet storage time beyond your institution's current limits should be evaluated per Directive [2004/33/EC] and validated according to local blood bank procedures. |

***The Mirasol PRT system for treating Platelets in PAS has been verified for these storage conditions. Certain conditions outside of these ranges may be acceptable but should be evaluated by the user as required.**

****Storage of products at yields between 4.5 and 5.1 × $10^{11}$ platelets or volumes less than 170 mL for a period of greater than three days may show evidence of a decline in some in vitro cell quality parameters. The extent of this decline may vary depending on collection and processing methods used by the site.**

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## Description

The Mirasol® Pathogen Reduction Technology (PRT) System consists of the Mirasol Illuminator, the Mirasol Manager Software (optional), and a Mirasol Disposable Kit. The sterile Disposable Kits contain an Illumination Bag, Storage Bag and Riboflavin Solution.

The Mirasol PRT system is for use with apheresis and whole blood-derived Plasma. The treatment procedure involves mixing the plasma product with Riboflavin Solution in the Illumination Bag and illuminating the mixture with a calculated light energy dose in the Illuminator. The Mirasol-treated Plasma product is then transferred from the Illumination Bag into the Storage Bag. The Mirasol-treated Plasma in the Storage Bag is frozen into Mirasol-treated fresh frozen Plasma (FFP) and is then ready for storage.

## Intended Use

The Mirasol PRT system is intended to reduce the pathogen load and inactivate residual white blood cells in donor Plasma for transfusion.

**Note:** Laboratory studies of photochemical treatment processing with the Mirasol PRT system for Plasma have shown a reduction in the levels of certain viruses, bacteria, and parasites as well as effective inactivation of residual white blood cells. However, no pathogen reduction process has been shown to eliminate all pathogens or exclude the possibility of adverse events due to contaminating white blood cells.

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
## Disposable Kit Overview

The Mirasol Plasma Disposable Kit consists of:

- Clear outer packaging with a white peel strip for opening.
- One opaque foil pouch containing one bag of Riboflavin Solution (500 µM riboflavin in normal saline) (see Figure 1).
- One Illumination/Storage Set with an Illumination Bag, a Storage Bag, attached tubing, frangible connectors, and a sample bulb (see Figure 1).

1 Inlet luer with blue end cap
2 Inlet lines
3 Luer with blue end cap on storage bag
4 Sample bulb
5 Storage bag
6 Illumination bag
7 Frangible connectors
8 "MA" set type label
9 Pre-attached label
10 Opaque foil pouch
11 Riboflavin solution bag

**Figure 1: Mirasol Plasma Disposable Kit**

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   Reuse of a single-use product could result in:

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- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

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### Disposable Kit

1. Use the Mirasol Plasma Disposable Kit described here only for treatment of Plasma. Different Mirasol Disposable Kits are available for treatment of Platelets.
2. Use only Disposable Kits manufactured for the Mirasol PRT system.
3. Do not use the Disposable Kit if it has passed the printed expiration date (YYYY-MM-DD).
4. Carefully inspect for leaks when using the disposable set. If a leak is observed, discard the blood product, or manage the product according to your institution's standard operating procedure (SOP).
5. The Illumination/Storage Set is sterilized with ethylene oxide (EO), which in rare cases has been reported to cause recipient reactions.
6. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.
7. The Mirasol Illumination Bag has not been qualified for centrifugation.

### Riboflavin Solution

1. Riboflavin Solution is sensitive to light. Treat the Plasma within 1 hour of removing the Riboflavin Solution Bag from the foil pouch.
2. Riboflavin Solution is not intended for direct intravenous infusion.
3. Riboflavin Solution can interfere with Plasma protein assays run on certain instruments. The type of instrument and sample dilution should be considered prior to running samples with Riboflavin Solution. Contact Terumo BCT, Inc. for validation instructions specific to the analytical instrument.

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### Blood Components

1. All Plasma products must be illuminated according to the treatment time constraints listed in Table 2, Incoming Product Specifications. This is to ensure adequate Plasma protein quality after treatment and storage, and to minimize the risk for potential endotoxin build-up to levels that may be harmful to the recipient.
2. After adding Riboflavin Solution and until the time of transfusion, protect the Plasma product from prolonged exposure to direct sunlight or strong ambient light sources. See Table 3 for storage specifications.

## Notes

1. Plasma products collected in ACD-A, CPD, CDPA-1 and 4% citrate anticoagulants have been evaluated for use in the Mirasol PRT system.
2. Use of Mirasol-treated FFP for multiple Plasma exchange procedures has not yet been evaluated.

## Instructions for Use

### Treating the Plasma Product with the Mirasol System

**Caution:** Use the Mirasol Plasma Disposable Kit described here only for treatment of Plasma. Different Mirasol Disposable Kits are available for treatment of Platelets.

1. Review the treatment time constraints listed in Table 2, Incoming Product Specifications, and the storage volume and other specifications listed in Table 3, Mirasol-Treated FFP Storage Conditions and Protein Quality.
2. Just prior to use, open the white peel strip on the outer packaging of the Disposable Kit, shown in Figure 1. Remove and inspect the Illumination/Storage Set.

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- **a.** The blue end cap on either the inlet line or the storage bag is unattached.
- **b.** The integrity of the Illumination/Storage Set is compromised for any reason.
- **c.** The frangible connectors are broken.
- **d.** The tubing is severely kinked.
- **e.** There is debris on the outer surface of the Illumination/Storage Bag. Do not clean with disinfection chemicals.
- **f.** The barcodes are illegible.
- **g.** Any pre-attached label is missing or not securely attached.

**Note:** In general, you should use the Riboflavin Solution Bag with the Illumination/Storage Set from the same Disposable Kit. However, if any of the circumstances listed in items c through g in the **Caution** above occurs after Riboflavin Solution has been added to the Plasma product, transfer the Plasma product to a new Illumination/Storage Set from the same lot number. Discard the original Illumination/Storage Set and the unused Riboflavin Solution.

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3. Permanently seal and disconnect the luer line on the Storage Bag using a tubing sealer device. Seal the line as near to the Storage Bag as possible, as shown in Figure 2, Item 4.
4. Sterile connect the Plasma Product Bag to the inlet line with the luer and blue end cap, which results in removal of the line with the luer and blue end cap. See Figure 2.

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
**5.** Appropriately label the Storage Bag according to your facility's Standard Operating Procedure, using space on the pre-attached label.

**Note:** Additional labels may be applied over the pre-attached label on the Storage Bag. Labels may not be applied to the Illumination Bag prior to illumination. The amount of UV energy delivered to the Plasma product is calculated using the surface area of the non-labeled Illumination Bag, among other factors.

**Note:** The set type label should be present.

See Figure 2, Item 7.

1. Plasma product bag
2. Inlet luer with blue end cap
3. Sterile connection on inlet line
4. Sealed luer line
5. Plasma storage bag
6. Pre-attached label
7. Set type label
8. Illumination bag

**Figure 2: Connecting the Plasma Product Bag to the Illumination/Storage Set**

**6.** Just before you are ready to transfer the Riboflavin Solution to the Plasma product, open the foil pouch, remove the Riboflavin Solution Bag, and inspect the bag and solution.

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
7. You can either transfer the Riboflavin Solution and the Plasma product to the Illumination Bag at the same time, or you can transfer them one at a time.
   - To transfer the Plasma product and the Riboflavin Solution at the same time:
     **a.** Sterile connect the Riboflavin Solution Bag to the second inlet line on the Illumination Bag. See Figure 3.
     **b.** Open both sterile connections and drain the Plasma product and the Riboflavin Solution into the Illumination Bag.
     **c.** Permanently seal the inlet line between the Plasma Product Bag and the Y connector, leaving the Riboflavin Solution Bag attached.
   - To transfer the Plasma product and the Riboflavin Solution one at a time:
     **a.** Open the sterile connection and drain the Plasma product into the Illumination Bag.
     **b.** Allowing enough tubing length for connecting to the Riboflavin Solution Bag, permanently seal the inlet line between the Plasma Product Bag and the Y connector.
     **c.** Sterile connect the Riboflavin Solution Bag to the same inlet line you used to transfer the Plasma product.
     **d.** Open the sterile connection and drain the Riboflavin Solution into the Illumination Bag.

1 Plasma product bag
2 Riboflavin solution bag
3 Sterile connections on inlet lines
4 Y connector
5 Pinch clamp on inlet line
6 Storage bag
7 Illumination bag

**Figure 3: Connecting the Riboflavin Solution Bag to the Illumination/Storage Set**


8. Express residual air from the Illumination Bag by completing the following steps:
   - Hold the Illumination Bag so that the ports are facing up, as shown in Figure 3.
   - Gently apply pressure to the Illumination Bag to move the Plasma up through the inlet line tubing until it just reaches the Riboflavin Solution Bag.
   - Drain the Plasma from the Riboflavin Solution Bag and inlet line into the Illumination Bag, closing the clamp on the inlet line before the air interface reaches the Illumination Bag port.

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**Figure 4: Residual air in the illumination bag**


**1** Acceptable amount of residual air

**2** Too much residual air

**9.** Permanently seal the inlet line as near to the Illumination Bag as possible. The sealed inlet line must be less than 2.5 cm (one inch) long, as shown in Figure 5, for proper fit in the Illuminator. Disconnect and discard the Riboflavin Solution Bag.

**10.** Move the clamp between the Illumination Bag and the Storage Bag as near to the Storage Bag as possible (see Figure 5, Item 3.)

**1** Sealed inlet line
**2** Frangible connectors
**3** Pinch clamp
**4** Sample bulb
**5** Illumination bag
**6** Storage bag

**Figure 5: Illumination/Storage Set with sealed inlet line and clamp near Storage Bag**

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11. Perform the illumination procedure according to the Mirasol Illuminator Operator's Manual.
12. If a post-illumination sample is required, break the frangible connector on the outlet port to the sample bulb. See Figure 5. Fill the sample bulb and permanently seal and disconnect the sample bulb from the Illumination Bag.

    If a post-illumination sample is not required, proceed to step 13.
13. Break the frangible connector on the outlet port of the line that connects the Illumination Bag to the Storage Bag. See Figure 5.

**Note:** The broken part of the frangible connector will remain in the Illumination Bag before and after transfer of the Mirasol-treated Plasma to the Storage Bag.

14. Drain the Mirasol-treated Plasma product from the Illumination Bag into the Storage Bag.
15. Express residual air from the Storage Bag by completing the following steps:
    - Hold the Storage Bag so that the ports are facing up.
    - Gently apply pressure to the Storage Bag to move the Plasma up through the Storage Bag inlet line tubing until it just reaches the Illumination Bag.
    - Drain the Plasma from the Illumination Bag and the Storage Bag inlet line into the Storage Bag, closing the clamp on the inlet line before the air interface reaches the Storage Bag port.
16. Permanently seal the Storage Bag inlet line as near to the Storage Bag as possible. Disconnect and discard the Illumination Bag.
17. Store the treated Plasma product.
    - Until the time of transfusion, protect the treated plasma product from prolonged exposure to direct sunlight or strong ambient light sources.
    - Follow the Plasma product storage specifications listed in Table 3, Mirasol-Treated FFP Storage Conditions and Protein Quality.

## Symbols and Certification



**Table 1: Symbols and Certification**

| Symbol | Description |
|---|---|
|  | Indicates the package contents. |
|   | Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste. |
|  | Indicates a non-pyrogenic fluid pathway. |
|  | Indicates that the product should not be used if the package is damaged. |
|  | Indicates that the product should not be used if the package is damaged. |
|  | Indicates that the product should not be vented. Do not remove the end caps from the set. |

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## Table 1: Symbols and Certification (continued)

| Symbol | Description |
|---|---|
|  | Indicates the presence of a sterile fluid path. The method of sterilization is ethylene oxide (EO). |
|  | Indicates the presence of a sterile fluid path. The method of sterilization is steam. |
|  | Indicates the medical device manufacturer. |
|  | Indicates the Authorized representative in the European Community. |
|  | Indicates the date after which the medical device is not to be used. |
|  | Indicates the manufacturer's catalog number so that the medical device can be identified. |
|  | Indicates the manufacturer's batch code so that the batch or lot can be identified. |
|  | Indicates a medical device that is not to be resterilized. |
|  | Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. |
|  | Indicates that the user should consult the instructions for use. |
|  | Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended. |
|  | Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP). |
|  | Indicates a medical device that needs protection from light sources. |
|  | Indicates a medical device that needs to be protected from moisture. |
|  | Indicates the range of temperatures to which the device can be safely exposed. |
|  | Indicates the product quantity when the quantity is placed in the square. |

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| Symbol | Description |
|---|---|
|  | Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile). |

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| Characteristic | Specification |
|---|---|
| Type of Plasma | Apheresis or whole-blood-derived Plasma |
| Anticoagulant | ACD-A or 4% citrate for apheresis<br>CPD or CPDA-1 for whole blood |
| Starting Plasma volume | 170 mL to 360 mL |
| Treatment time constraints | Apheresis Plasma should be treated and frozen to ≤ −30 °C within 6 hours of collection, but is qualified out to 8 hours after collection.<br>Whole blood-derived Plasma must be treated and frozen to ≤ −30 °C within 6 hours of separation. Plasma may be stored in whole blood form at 4 °C ± 2 °C or 22 °C ± 2 °C for up to 18 hours prior to separation.<br>Illumination of previously frozen Plasma is acceptable if the product was initially frozen within 8 hours of collection. Once thawed, the Plasma must be treated and re-frozen within 2 hours. |
| Residual white blood cells | ≤ $250 \times 10^6$ per unit of $3.5 \times 10^{11}$ yield |
| Platelet content | Residual levels qualified to concentrations of ≤ $2.1 \times 10^{12}$ platelets per liter |
| *The Mirasol PRT system for treating Plasma has been verified for these incoming product conditions. Certain conditions outside of these ranges may be acceptable but should be evaluated by the user, as required. | |

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**Table 3: Mirasol-Treated FFP Storage Conditions and Protein Quality**

| Characteristic | Specification |
|---|---|
| Mirasol-treated FFP storage period | Maximum storage period: Qualified out to 2 years from date of collection when stored at ≤ −30 °C. |
| Illumination/storage volumes | 200 mL to 400 mL after addition of Riboflavin Solution and treatment |
| Thawed Mirasol-treated FFP storage period | Mirasol-treated FFP should be used immediately when thawed and is qualified out to 6 hours after thaw. Mirasol-treated FFP is qualified for thawing in a water bath or microwave. |
| Mirasol-treated FFP protein quality | Mirasol-treated FFP meets the Council of Europe guidelines for protein content of FFP. |

**Table 4: Disposable Kits***

| Characteristic | Specification |
|---|---|
| Mirasol Plasma Disposable Kit | Each Disposable Kit consists of one Illumination/Storage Set and one 35 mL bag of 500 µM riboflavin in 0.9% sodium chloride solution. |
| Illumination/Storage Set | Citrate vinyl bag connected to a PVC Universal Blood Storage Bag |
| Illumination/Storage Set sterility | The Illumination/Storage Set has a sterile and non-pyrogenic fluid pathway, sterilized by ethylene oxide (EO). |
| Illumination/Storage Set label | The pre-attached label covers 110 $cm^2$ (17 $in^2$) surface area of the Storage Bag.<br>Additional labels may be applied over the pre-attached label.<br>No labels should be applied to the Illumination Bag prior to illumination. |
| Tubing dimensions for sterile connections and seals | Outer diameter = 4.06 mm ± 0.08 mm<br>Wall thickness = 0.56 mm ± 0.05 mm |
| Illumination Bag unrecoverable volume | 3.4 mL |
| Storage Bag unrecoverable volume | 3.7 mL |
| Sample bulb volume | 1.5 mL to 4.0 mL |
| Riboflavin Solution content | 500 µM riboflavin in 0.9% sodium chloride solution<br>pH 3.5 to 5.1<br>Riboflavin Material Safety Data Sheets are available upon request. |
| Riboflavin Solution sterility | The Riboflavin Solution is sterile and non-pyrogenic, sterilized by steam. |

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**Table 4: Disposable Kits* (continued)**

| Characteristic | Specification |
|---|---|
| Riboflavin Solution volume | 35 ± 5 mL |
| Riboflavin Solution raw materials | Riboflavin: USP, EP<br>Sodium chloride: USP, EP<br>Hydrochloric acid: USP, EP<br>Water for injection: USP, EP |
| Riboflavin Solution time constraints | Treat the Plasma product within 1 hour of removing Riboflavin Solution Bag from foil pouch. |
| Disposable kit storage conditions | Long-term storage range: 2 °C to 30 °C (36 °F to 86 °F)<br>Storage humidity: Up to 75% relative humidity (RH), non-condensing<br>Permitted excursions: Up to 60 °C (140 °F) for up to 2 weeks |
| Case quantity | 25 kits per case |
| Case size | 24.9 cm × 35.8 cm × 35.8 cm |
| Case weight | 6.0 kg |

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This page intentionally left blank.

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## Description

The Mirasol® Pathogen Reduction Technology (PRT) System consists of the Mirasol Illuminator, the Mirasol Manager Software (optional), and a Mirasol Disposable Kit. The sterile disposable kit contains an illumination bag, a storage bag, and riboflavin solution.

The Mirasol PRT system is for use with Whole Blood. The treatment procedure involves mixing the Whole Blood product with riboflavin solution in the illumination bag and illuminating the mixture with a calculated light energy dose in the Illuminator. The Mirasol-treated Whole Blood product is then transferred from the illumination bag into the storage bag.

## Intended Use

The Mirasol Pathogen Reduction Technology (PRT) System for Whole Blood is an extracorporeal blood treatment process that uses ultraviolet light energy and riboflavin (vitamin $B_2$) to reduce the pathogen load and to inactivate white blood cells in Whole Blood for transfusion.

## Indications for Use

- Mirasol PRT system treatment can be used instead of gamma irradiation to prevent transfusion-associated graft-versus-host disease (TA-GVHD).
- Mirasol PRT system treatment has demonstrated the ability to reduce the occurrence of transfusion-related transmission of malaria.
- For additional information regarding the Mirasol-treated Whole Blood product, including in vitro and in vivo data, refer to the Physician's Information Booklet (Terumo BCT part number 777714-013).

**Note:** Laboratory studies of photochemical treatment processing with the Mirasol PRT system for Whole Blood have shown a reduction in the levels of certain viruses, bacteria, and parasites, as well as effective inactivation of residual white blood cells. However, no pathogen reduction process has been shown to eliminate all pathogens or excludes the possibility of adverse events due to contaminating white blood cells.

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   - Product performance issues due to a loss of product integrity, including but not limited to the following:
     - Fluid leaks
     - Parts that are warped or deformed
     - Plastics that are brittle and discolored
     - Filters that have reduced filtration capabilities
   - Exposure to excessive ethylene oxide (EO) residuals

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- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross‑contamination

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2. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.
3. Only Whole Blood products collected in either CPD anticoagulant or CPDA-1 anticoagulant have been evaluated for use with the Mirasol PRT system. Do not use Whole Blood products collected in any other anticoagulant.

## Cautions

For additional information regarding the Mirasol-treated Whole Blood product, including in vitro and in vivo data, refer to the Physician's Information Booklet (Terumo BCT part number 777714-013).

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### Disposable Kit

1. Use the Mirasol Whole Blood Disposable Kit described here only for treatment of Whole Blood. Contact your Terumo BCT representative for information about other Mirasol products.
2. Use only Disposable Kits manufactured for the Mirasol PRT system to ensure compatibility and adherence with operational specifications.
3. Do not use the disposable kit if the kit label, the bag label, or the riboflavin solution bag label has passed the printed expiration date (YYYY-MM-DD).
4. Carefully inspect for leaks when using the disposable kit. If a leak is observed, manage the product according to your facility's standard operating procedure (SOP).
5. The Illumination/Storage Set is sterilized with ethylene oxide (EO), which in rare cases has been reported to cause recipient reactions.
6. The Mirasol Illumination Bag has not been qualified for centrifugation.

### Riboflavin Solution

1. Riboflavin solution is sensitive to light. Add the riboflavin solution to the Whole Blood product within 1 hour of removing the riboflavin solution bag from the foil pouch.
2. Riboflavin solution is not intended for direct intravenous infusion.
3. Do not add more than one bag of riboflavin solution to any one product.

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### Blood Components

1. All Whole Blood products must be illuminated according to the treatment time constraints listed in Table 1 in order to ensure adequate product quality after treatment and storage.
2. After adding riboflavin solution and until the time of transfusion, protect the Whole Blood product from exposure to direct sunlight or strong ambient light sources. See Table 3 for storage specifications.
3. Mirasol PRT system treatment and storage of Whole Blood result in levels of extracellular potassium that are higher than those in non-treated Whole Blood units. Considerations should be made for transfusion recipients who may be susceptible to hyperkalemia.

## Disposable Kit Overview

The Mirasol Whole Blood Disposable Kit consists of:

- Clear outer packaging with a white peel strip for opening
- One opaque foil pouch containing one bag of riboflavin solution (500 µM riboflavin in normal saline) (see Figure 1)
- One illumination/storage set with an illumination bag, a storage bag, and attached tubing and frangible connectors (see Figure 1)


**Figure 1: Mirasol Whole Blood Disposable Kit**

**1** Inlet line to the illumination bag with a luer and a blue end cap
**2** Sealed inlet line
**3** Drain line from the illumination bag to the storage bag
**4** Vent line with a blue end cap on the storage bag
**5** Storage bag
**6** Pre-attached label
**7** Set type label

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**8** Spike port

**9** Frangible connector on the spike port

**10** Frangible connector on the drain line from the illumination bag to the storage bag

**11** Illumination bag

**12** Opaque foil pouch

**13** Riboflavin solution bag

**14** Drain line for the riboflavin solution

**Note:** Callouts Ⓐ through Ⓔ are reference points for the steps in "Treating the Whole Blood Product with the Mirasol PRT System" beginning on page 7.

## Ancillary Equipment

The following ancillary equipment is needed to complete a Mirasol treatment procedure:

- Tubing sealer device
- Sterile connection device
- Hematocrit measurement device

Use this equipment according to your facility's SOP.

## Incoming Product Specifications

**Table 1: Incoming product specifications**

| Characteristic | Specification |
|---|---|
| Anticoagulant | CPD or CPDA-1 |
| Starting Whole Blood volume | **463 mL to 541 mL**, including anticoagulant and donor blood<br>This acceptable range is calculated as follows:<br>(450 mL +5%/−10% donor blood) + (63 mL ± 5 mL AC) |
| Treatment time constraints | Whole Blood must be treated within 24 hours of collection. |
| Illumination hematocrit (after riboflavin solution has been added) | 30% to 45% |
| Storage temperature before Mirasol PRT system treatment | Follow your facility's SOP for Whole Blood products. |

Whole Blood donation volume for Mirasol PRT system treatment is specified as **450 mL +5%/ −10%**. The addition of anticoagulant increases the volume by **63 mL ± 5 mL**. Prior to Mirasol PRT system treatment, **35 mL ± 5 mL** riboflavin solution is added to the Whole Blood unit. The following table shows the low limit and high limit of possible volumes and the weight of each volume.

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**Table 2: Whole Blood volume and weight specifications**

| | Low Limit | Target | High Limit |
|---|---|---|---|
| Volume/weight* of Whole Blood collection (450 mL +5%/−10%) | 405 mL/426 g | 450 mL/474 g | 473 mL/502 g |
| Volume/weight* of Whole Blood unit including anticoagulant (63 mL ± 5 mL) | 463 mL/485 g | 513 mL/539 g | 541 mL/572 g |
| Volume/weight* of Whole Blood unit including anticoagulant and riboflavin solution (35 mL ± 5 mL) | 493 mL/516 g | 548 mL/575 g | 581 mL/613 g |
| *Weight does not include disposable set tare weight. | | | |

## Instructions for Use

### Treating the Whole Blood Product with the Mirasol PRT System

**Caution:** Use the Mirasol Whole Blood Disposable Kit described here only for treatment of Whole Blood. Different Mirasol disposable kits are available for treatment of platelets and plasma.

1. Review the treatment time constraints listed in Table 1 and the storage conditions listed in Table 3.
2. Just prior to use, open the white peel strip on the outer packaging of the disposable kit. Remove and inspect the illumination/storage set.

**Caution:** Do not use the disposable kit under any of the following circumstances:

- **a.** The blue end cap on either the inlet line (Ⓓ) or the storage bag line (Ⓐ) is not attached.
- **b.** The integrity of the illumination/storage set is compromised for any reason.
- **c.** The frangible connectors are broken.
- **d.** The tubing is severely kinked.
- **e.** There is debris on the outer surface of the illumination bag. Do not clean with disinfection chemicals.
- **f.** The barcodes are illegible.
- **g.** The pre-attached label is missing or not securely attached.


**Note:** In general, you should use the riboflavin solution bag with the illumination/storage set from the same disposable kit. However, if any of the circumstances listed in items c through g in the **Caution** above occurs after the riboflavin solution has been added to the Whole Blood product, transfer the Whole Blood product to a new illumination/storage set from the same lot number. Discard the original illumination/storage set and the unused riboflavin solution.

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**3.** Permanently seal and disconnect the vent line on the storage bag (Line Ⓐ) using a tubing sealer device. Seal the line as near to the storage bag as possible, as shown in Figure 2.

**1** Disconnected vent line with a blue end cap
**2** Sealed vent line on the storage bag
**3** Storage bag
**4** Illumination bag


**Figure 2: Removing the vent line**

**4.** Before connecting the whole blood product bag to the illumination bag, verify that the Whole Blood product contains only CPD anticoagulant or CPDA-1 anticoagulant. **If the Whole Blood product does not contain either CPD anticoagulant or CPDA-1 anticoagulant, do not continue with the Mirasol process.** Only Whole Blood products collected in either CPD anticoagulant or CPDA-1 anticoagulant have been evaluated for use with the Mirasol PRT system.

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**5.** Using a sterile connection device, connect the whole blood product bag to the illumination bag inlet line (Line Ⓓ), and remove the part of the line that has a luer and a blue end cap, as shown in Figure 3. **Do not** open the sterile connection at this time. (You will transfer the Whole Blood and the riboflavin solution to the illumination bag at the same time in step 8.)

1 Whole blood product bag (donor blood bag)
2 Sterile connection on the inlet line
3 Disconnected inlet line with a luer and a blue end cap
4 Sealed vent line on the storage bag
5 Storage bag
6 Pre-attached label
7 Set type label
8 Illumination bag

**Figure 3: Connecting the whole blood product bag to the illumination bag inlet line**

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6. Transfer an identification label from the donor whole blood product bag to the storage bag according to your facility's standard operating procedure (SOP), using the space on the pre-attached label. Additional labels may be applied over the pre-attached label on the storage bag.

**Note:** Labels may not be applied to the illumination bag before illumination. The amount of UV energy delivered to the Whole Blood product is calculated using the surface area of the non-labeled illumination bag, among other factors.

**Note:** The following set type label should be present on the illumination bag in the position shown in Figure 3, Item 7:

7. Just before you are ready to transfer the riboflavin solution to the illumination bag, open the foil pouch, remove the riboflavin solution bag, and inspect the bag and solution.

**Caution:** Do not use the riboflavin solution if:

- The integrity of the riboflavin solution bag is compromised for any reason.
- Any crystals or particulates are visible in the riboflavin solution.

**Note:** In general, you should use the riboflavin solution bag with the illumination/storage set from the same disposable kit. However, if either of the two circumstances listed in the **Caution** above occurs, use a new riboflavin solution bag of the same lot number from another disposable kit and discard the unused illumination/storage set.

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
**8.** Transfer the Whole Blood product and the riboflavin solution into the illumination bag at the same time by performing these steps:

**a.** Using a sterile connection device, connect the tubing on the riboflavin solution bag (Line Ⓔ) to the sealed inlet line on the illumination bag (Line ⓒ). See Figure 4.

**b.** Open both sterile connections and drain the Whole Blood product and the riboflavin solution into the illumination bag.


**1** Whole blood product bag (donor blood bag)
**2** Riboflavin solution bag
**3** Sterile connections on the inlet lines
**4** Y connector
**5** Sealed vent line on the storage bag
**6** Storage bag
**7** Illumination bag

**Figure 4: Transferring the Whole Blood product and the riboflavin solution into the illumination bag**

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**c.** Permanently seal and disconnect the inlet line between the empty whole blood product bag and the Y connector, leaving the riboflavin solution bag attached. See Figure 5.

**1** Whole blood product bag (donor blood bag)
**2** Riboflavin solution bag
**3** Sealed Whole Blood product inlet line
**4** Y connector
**5** Storage bag
**6** Illumination bag

**Figure 5: Disconnecting the whole blood product bag**

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9. Express residual air from the illumination bag into the empty riboflavin solution bag by performing these steps:
   - **a.** Hold the illumination bag so that the ports are facing up, as shown in Figure 6.
   - **b.** Gently apply pressure to the illumination bag to move the residual air and Whole Blood up through the inlet line, allowing the Whole Blood to fill the riboflavin solution bag approximately half full.
   - **c.** Close the pinch clamp, and then gently mix the Whole Blood in the riboflavin solution bag to ensure that the solution is homogeneous. Gently tap the ports on the riboflavin solution bag to remove any solution that may have been trapped during mixing. Do not shake or vigorously mix the Whole Blood in the riboflavin solution bag, as this can cause the formation of air bubbles and possible hemolysis.
   - **d.** Open the pinch clamp to drain the Whole Blood from the riboflavin solution bag and the inlet line into the illumination bag, closing the pinch clamp on the inlet line before the air interface reaches the illumination bag port. Gently mix the Whole Blood and the riboflavin solution to obtain a homogeneous solution.


**Note:** It is very important that the product be gently but thoroughly mixed to achieve an accurate hematocrit measurement.

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**e.** Open the pinch clamp on the inlet line and hold the illumination bag so that the ports are facing up, as shown in Figure 6. Gently apply pressure to the illumination bag to move the residual air and Whole Blood up through the inlet line until it reaches halfway between the Y connector and the riboflavin solution bag. Close the pinch clamp on the inlet line.

1 Riboflavin solution bag
2 Inlet line
3 Y connector
4 Pinch clamp
5 Storage bag
6 Illumination bag

**Figure 6: Expressing residual air from the illumination/storage set**

Residual air should now be removed from the illumination bag. The following figure shows acceptable and unacceptable amounts of residual air.

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
**Figure 7: Residual air in the illumination bag**

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
10. Permanently seal the inlet line just below the riboflavin solution bag using a tubing sealer device (see Figure 8).
11. Disconnect and discard the riboflavin solution bag.
12. Permanently seal the inlet line as near as possible to the illumination bag, using a tubing sealer device. The sealed inlet line must be less than 2.5 cm (1 in) long for proper fit into the Illuminator. Disconnect the tubing segment between the illumination bag and the riboflavin solution bag. (This segment is approximately 18 cm to 30 cm [8 in to 12 in] long. See Figure 8.) This segment will be used to collect a hematocrit sample.

**1** Disconnected riboflavin solution bag
**2** Y connector
**3** Tubing segment for the hematocrit sample
**4** Pinch clamp, moved as near as possible to the storage bag
**5** Storage bag
**6** Pre-attached label
**7** Set type label
**8** Remaining segment of the sealed inlet line (less than 2.5 cm/1 in)
**9** Illumination bag

**Figure 8: Sealing and disconnecting the riboflavin solution bag and the hematocrit sample segment**

13. Extract a hematocrit sample from the segment of inlet line tubing that was disconnected.
    - Extract the sample using, for example, a plain tube or a microcentrifuge tube; **do not** use a tube that contains EDTA or AC, which will dilute the sample and change its hematocrit.
    - Follow your facility's SOP first to mix the sample so that it is homogeneous, and then to measure the hematocrit.

**Note:** Always use the hematocrit of the pre-illumination Whole Blood (including riboflavin solution). If the hematocrit of the blood donor or of the collected Whole Blood is used, the blood product will receive more light energy than necessary, and the quality of the Whole Blood may be decreased.

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
14. Move the pinch clamp on the line between the storage bag and the illumination bag (Line Ⓑ) so it is as near as possible to the storage bag (see Figure 8).
15. Perform the illumination procedure according to the instructions in the Mirasol PRT System for Whole Blood Illuminator operator's manual.
16. After the illumination procedure, remove the illumination/storage set from the Illuminator, hang the illumination bag, and break the frangible connector on the outlet port of the line that connects the illumination bag to the storage bag (Line Ⓑ). See Figure 9.

**Note:** The broken part of the frangible connector will remain in the illumination bag before and after the Mirasol-treated Whole Blood is transferred to the storage bag.

17. Drain the treated Whole Blood product into the storage bag.

1 Illumination bag
2 Broken frangible connector
3 Drain line from the illumination bag (also known as "storage bag inlet line" and Line Ⓑ)
4 Pinch clamp
5 Storage bag

**Figure 9: Draining the Whole Blood product**

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
**18.** Express residual air from the storage bag by performing these steps:

**a.** Hold the storage bag so that the ports are facing up, as shown in Figure 10.

**b.** Gently apply pressure to the storage bag to move the residual air and Whole Blood up through the storage bag inlet line (Line Ⓑ) until it just reaches the illumination bag.

**c.** Close the pinch clamp on the inlet line to hold the blood in the tubing.


**1** Illumination bag
**2** Drain line from the illumination bag to the storage bag
**3** Pinch clamp
**4** Storage bag spike port
**5** Storage bag

**Figure 10: Expressing residual air from the storage bag**

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**19.** Permanently seal the storage bag inlet line (Line Ⓑ) near the illumination bag, and then disconnect and discard the illumination bag (see Figure 11).

**1** Disconnected illumination bag
**2** Drain line from the illumination bag to the storage bag
**3** Sealed storage bag inlet line
**4** Pinch clamp
**5** Storage bag

**Figure 11: Sealing and disconnecting the illumination bag**

**20.** Remove the pinch clamp from the storage bag inlet line (Line Ⓑ), and make tubing segments according to your facility's SOP.

**21.** Label the product and release it to inventory according to your facility's SOP.

- Until the time of transfusion, protect the treated Whole Blood product from exposure to direct sunlight or strong ambient light sources.
- Follow the Whole Blood product storage specifications listed in Table 3.

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## Storage Conditions

**Table 3: Mirasol-treated Whole Blood storage conditions and quality***

| Characteristic | Specification |
|---|---|
| Time Between Treatment and Storage | Place the treated Whole Blood product into storage within 45 minutes.<br>Until the time of transfusion, protect the treated Whole Blood product from exposure to direct sunlight or strong ambient light sources. |
| Maximum Storage Period | • Whole Blood units that are held at controlled room temperature or below prior to treatment on the Mirasol PRT system within 24 hours of collection have been verified for storage for up to 14 days at 1 °C to 6 °C.<br>• Whole Blood units that are held at controlled room temperature or below prior to treatment on the Mirasol PRT system within 8 hours of collection have been verified for storage for up to 21 days at 1 °C to 6 °C. |
| Storage Temperature | 1 °C to 6 °C |
| Transport Temperature | 1 °C to 10 °C |
| Leukoreduction | Optional. Follow your facility's SOP. |
| *The Mirasol PRT system for Whole Blood has been verified for these product storage conditions. Certain conditions outside of these ranges may be acceptable but should be evaluated by the user as required. | |

## Mirasol Whole Blood Disposable Kit Specifications

**Table 4: Mirasol Whole Blood disposable kit**

| Characteristic | Specification |
|---|---|
| Mirasol Whole Blood Disposable Kit | Each disposable kit consists of one illumination/storage set and one 35 mL bag of 500 µM riboflavin in 0.9% sodium chloride solution. |
| Illumination/storage set | Citrate vinyl bag connected to a PVC universal blood storage bag. |
| Illumination/storage set sterility | The illumination/storage set has a sterile and non-pyrogenic fluid pathway, sterilized by ethylene oxide (EO). |
| Illumination/storage set label | The pre-attached label covers a 110 $cm^2$ (17 $in^2$) surface area of the storage bag.<br>Additional labels may be applied over the pre-attached label.<br>**Note: Do not apply any labels to the illumination bag prior to illumination.** |
| Tubing dimensions for sterile connections and seals | Outer diameter = 4.06 mm ± 0.08 mm<br>Wall thickness = 0.56 mm ± 0.05 mm |

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**Table 4: Mirasol Whole Blood disposable kit (continued)**

| Characteristic | Specification |
|---|---|
| Riboflavin solution content | 500 µM riboflavin in 0.9% sodium chloride solution<br>pH 3.5 to 5.1<br>Riboflavin Material Safety Data Sheets are available upon request. |
| Riboflavin solution sterility | The riboflavin solution is sterile and non-pyrogenic, sterilized by steam. |
| Riboflavin solution volume | 35 mL ± 5 mL |
| Riboflavin solution raw materials | Riboflavin: USP, EP<br>Sodium chloride: USP, EP<br>Hydrochloric acid: USP, EP<br>Water for injection: USP, EP |
| Disposable kit storage conditions | Long-term storage range: 2 °C to 30 °C (36 °F to 86 °F)<br>Storage humidity: Up to 75% relative humidity (RH), non-condensing<br>Permitted excursions: Up to 60 °C (140 °F) for up to 2 weeks |
| Case quantity | 25 kits per case |
| Case size | 24.9 cm × 35.8 cm × 35.8 cm |
| Case weight | 6.0 kg |

## Symbols and Certification


**Table 5: Symbols and certification**

| Symbol | Description |
|---|---|
| CE 2797 | Indicates the product was manufactured in accordance with Annex II of the European Council Directive 93/42/EEC, as amended. |
| EC REP | Indicates the Authorized representative in the European Community. |
|  | Indicates the medical device manufacturer. |
|  | Indicates the package contents. |
|  | Indicates that the user should consult the instructions for use. |
|  | Indicates the product quantity when the quantity is placed in the square. |

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**Table 5: Symbols and certification (continued)**

| Symbol | Description |
|---|---|
|   | Indicates the manufacturer's catalog number so that the medical device can be identified. |
|   | Indicates the manufacturer's batch code so that the batch or lot can be identified. |
|  | Indicates that the product should not be used if the package is damaged. |
|  | Indicates the date after which the medical device is not to be used. |
|  | Indicates a non-pyrogenic fluid pathway. |
|   | Indicates the presence of a sterile fluid path. The method of sterilization is ethylene oxide (EO). |
|   | Indicates the presence of a sterile fluid path. The method of sterilization is steam. |
|  | Indicates the range of temperatures to which the device can be safely exposed. |
|  | Indicates that the product should not be vented. Do not remove the end caps from the set. |
|   | Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP). |
|  | Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. |
|   | Indicates a medical device that is not to be resterilized. |
|   | Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste. |
|  | Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile). |

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## Description

The Mirasol® Pathogen Reduction Technology (PRT) System consists of the Mirasol Illuminator, the Mirasol Manager software (optional) and a Mirasol disposable kit. The treatment procedure involves mixing a blood product with Riboflavin Solution in the Illumination/Storage Bag and illuminating the mixture with a calculated light energy dose in the Illuminator.

The Mirasol Platelet Disposable Kits for Treatment in Plasma are for use with apheresis Platelets and whole-blood-derived Platelets. These kits may also be used to treat Plasma-reduced platelet concentrate (PPC) in plasma with the addition of Platelet Additive Solution (PAS) following treatment. These sterile disposable kits both contain an Illumination/Storage Bag and Riboflavin Solution. The Mirasol DOUBLE+ Platelet Disposable Kit for Treatment in Plasma also contains a Platelet Additive Solution (PAS) Assembly that allows PAS to be added to the treated PPC products through a sterile barrier filter and a pre-attached second Platelet Storage Bag designed for splitting products without connecting a separate storage bag. After optional splitting and sampling, Mirasol-treated Platelet products are ready for storage or transfusion.

**Note:** The Mirasol Platelet Disposable Kits for Treatment in Plasma described here are for treatment of Platelets in plasma only, with an optional addition of PAS following treatment. For more information about additional Mirasol Platelet Disposable Kits, contact your Terumo BCT representative.

## Intended Use

The Mirasol PRT system is intended to reduce the pathogen load and inactivate residual white blood cells in donor Platelet concentrates for transfusion.

**Note:** Laboratory studies of photochemical treatment processing with the Mirasol PRT System for Platelets have shown a reduction in the levels of certain viruses, bacteria, and parasites as well as effective inactivation of residual white blood cells. However, no pathogen reduction process has been shown to eliminate all pathogens or exclude the possibility of adverse events due to contaminating white blood cells.

Mirasol PRT system treatment can be used instead of gamma irradiation to prevent transfusion-associated graft-versus-host disease (TA-GVHD).

Mirasol PRT system treatment in combination with leukoreduction can be used instead of cytomegalovirus (CMV) testing to prevent transfusion-transmitted CMV infection.

## Contraindications

There are no known contraindications for use of this product.

## Disposable Kit Overview

Both Mirasol Platelet Disposable Kits for Treatment in Plasma consist of:


- Clear outer packaging with a white peel strip for opening.
- One opaque foil pouch containing one bag of Riboflavin Solution (500 µM riboflavin in normal saline) (see Figure 1 and Figure 2).
- One Illumination/Storage Bag with attached tubing, frangible connectors, spike port, and sample bulbs (see Figure 1 and Figure 2).

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
The Mirasol DOUBLE+ Platelet Disposable Kit for Treatment in Plasma also contains:

- A PAS assembly, which consists of a sample bag and sterile barrier filter.
- A second platelet storage bag with attached spike port and sample bulb (see Figure 2).

**1** Inlet luer with blue end cap
**2** Inlet lines
**3** Spike port
**4** "LT" set type label
**5** Opaque foil pouch
**6** Riboflavin solution bag
**7** Sample bulb
**8** Frangible connectors
**9** Pre-attached label
**10** Illumination/storage bag

**Figure 1: Mirasol Platelet Disposable Kit for Treatment in Plasma ("single-bag kit")**

**1** Inlet luer with blue end cap
**2** Inlet lines
**3** Spike ports
**4** "AS" set type label
**5** Opaque foil pouch
**6** Riboflavin solution bag
**7** Sample bulb
**8** Frangible connectors
**9** Pre-attached labels
**10** Illumination/storage bag
**11** Second platelet storage bag
**12** PAS assembly

**Figure 2: Mirasol DOUBLE+ Platelet Disposable Kit for Treatment in Plasma**

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## Warnings

1. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re‑sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re‑sterilized.

   Reuse of a single‑use product could result in:

   - Product performance issues due to a loss of product integrity, including but not limited to the following:
     - Fluid leaks
     - Parts that are warped or deformed
     - Plastics that are brittle and discolored
     - Filters that have reduced filtration capabilities
   - Exposure to excessive ethylene oxide (EO) residuals
   - Viral infections such as hepatitis or human immunodeficiency virus (HIV)
   - Bacterial infections
   - Cross‑contamination

   Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

## Cautions

### General

1. Each operator should be thoroughly familiar with these Instructions for Use and the Mirasol Illuminator operator's manual.
2. Use aseptic technique throughout all procedures and follow standard blood establishment practices.

### Disposable Kit

1. Use the Mirasol Platelet Disposable Kits for Treatment in Plasma described here only for treatment of Platelets in plasma.
2. Use only Disposable Kits manufactured for the Mirasol PRT system.
3. Do not use the Disposable Kit if it has passed the printed expiration date (YYYY-MM-DD).
4. Carefully inspect for leaks when using the Disposable Kit. If a leak is observed, do not use the Platelet product, or manage the product according to your institution's standard operating procedure (SOP).
5. The Illumination/Storage Set is sterilized with ethylene oxide (EO), which in rare cases has been reported to cause recipient reactions.
6. The light energy dose delivered to the Platelet product is calculated based on the Illumination/ Storage Set tare weight configured on the Illuminator, among other factors. Therefore, the configured tare weight must be accurate. Do not remove components from the Illumination/Storage

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Set (except the Illumination Bag inlet line) prior to illumination unless the correct tare weight has been configured on the Illuminator. See the Mirasol Illuminator operator's manual for instructions on setting tare weight. If the correct tare weight cannot be configured, contact Terumo BCT, Inc.

7. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.
8. The Mirasol Illumination Bag has not been qualified for centrifugation.

#### Riboflavin Solution

1. Riboflavin Solution is sensitive to light. Treat the Platelet product within 1 hour of removing the Riboflavin Solution Bag from the foil pouch.
2. Riboflavin Solution is not intended for direct intravenous infusion.
3. Samples containing Riboflavin Solution cannot be measured for $pO_2$ levels when measured on the Radiometer ABL® 700 Series Blood Gas Analyzer. The riboflavin interferes with this instrument's $pO_2$ measurement.

#### Blood Components

1. Do not agitate Platelet products during the minimum rest period listed in Table 2.
2. All Platelet products must be illuminated according to the treatment time constraints listed in Table 2. This is to ensure adequate Platelet cell quality after treatment and storage, and to minimize the risk for potential endotoxin build-up to levels that may be harmful to the recipient.
3. After adding Riboflavin Solution and until the time of transfusion, protect the Platelet product from prolonged exposure to direct sunlight or strong ambient light sources. See Table 4 for storage specifications.
4. If the physician transfusing the Platelet product prescribes gamma irradiation, gamma irradiate the Platelet product only after treating with the Mirasol PRT system.
5. If the product contains any platelet aggregates, do not continue with Mirasol PRT system treatment until platelet aggregates have dissipated.

## Notes

Only Platelet products collected in ACD-A, CPD, or CPDA-1 anticoagulants have been evaluated for use in the Mirasol PRT system.

## Instructions for Use

#### Treating the Platelet Product with the Mirasol PRT System

**Caution:** Use the Mirasol Platelet Disposable Kits for Treatment in Plasma described here only for treatment of Platelets in plasma. Different Mirasol Disposable Kits are available for treatment of Platelets in Platelet Additive Solution (PAS).

**Note:** The following steps and figures refer to both Mirasol Platelet Disposable Kits for Treatment in Plasma. Unique considerations for the DOUBLE+ Platelet Disposable Kit for Treatment in Plasma are noted.

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1. Review the treatment time constraints listed in Table 2 and the storage volume, concentration, and other specifications listed in Table 4.
2. If treating an apheresis Platelet product, allow it to rest according to the rest periods listed in Table 2.
3. Just prior to use, open the white peel strip on the outer packaging of the Disposable Kit. Remove the Illumination/Storage Set from the white spacer and inspect the set.


**Caution:** Do not use the Disposable Kit under any of the following circumstances:

- **a.** The blue end cap on the inlet line is unattached.
- **b.** The integrity of the Illumination/Storage Set is compromised for any reason.
- **c.** The frangible connectors are broken.
- **d.** The tubing is severely kinked.
- **e.** There is debris on the outer surface of the Illumination/Storage Bag. Do not clean with disinfection chemicals.
- **f.** The barcodes are illegible.
- **g.** Any pre-attached label is missing or not securely attached.

**Note:** In general, you should use the Riboflavin Solution Bag with the Illumination/Storage Set from the same Disposable Kit. However, if any of the circumstances listed in items c through g in the Caution above occurs after Riboflavin Solution has been added to the Platelet product, transfer the Platelet product to a new Illumination/Storage Set from the same lot number. Discard the original Illumination/Storage Set and the unused Riboflavin Solution.

**Note:** Use sterile connection devices and tubing sealer devices according to standard blood establishment practice.

4. Sterile connect the Platelet Product Bag to the part of the Illumination/Storage Bag inlet line with the luer and blue end cap, which results in removal of the line with the luer and the blue end cap. See Figure 3.

1. Platelet product bag
2. Inlet luer with blue end cap
3. Sterile connection on inlet line
4. Illumination/storage bag
5. Pre-attached label

**Figure 3: Connecting the Platelet Product Bag to the Illumination/Storage Bag**

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5. Appropriately label the Illumination/Storage Bag according to your facility's standard operating procedure (SOP), using space on the pre-attached label. For the DOUBLE+ set, also label the second Platelet Storage Bag, if required by your facility's SOP.

**Note:** Additional labels may be applied over the pre-attached label on the Illumination/Storage Bag prior to illumination, but they must not cover any bag surface area outside the edges of the pre-attached label. The amount of UV energy delivered to the Platelet product is calculated using the surface area of the pre-attached label, among other factors.

6. Just before you are ready to transfer the Riboflavin Solution to the Illumination/Storage Bag, open the foil pouch, remove the Riboflavin Solution Bag, and inspect the bag and solution.


**Caution:** Do not use the Disposable Kit if:

- The integrity of the Riboflavin Solution Bag is compromised for any reason.
- Any riboflavin crystals or particulates are visible in the Riboflavin Solution.

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7. You can either transfer the Riboflavin Solution and the Platelet product to the Illumination/Storage Bag at the same time, or you can transfer them one at a time.
   - To transfer the Platelet product and the Riboflavin Solution at the same time:
     **a.** Sterile connect the Riboflavin Solution Bag to the second inlet line on the Illumination/ Storage Bag. See Figure 4.
     **b.** Open both sterile connections and drain the Platelet product and the Riboflavin Solution into the Illumination/Storage Bag.
     **c.** Permanently seal the inlet line between the Platelet Product Bag and the Y connector, leaving the Riboflavin Solution Bag attached.
   - To transfer the Platelet product and the Riboflavin Solution one at a time:
     **a.** Open the sterile connection and drain the Platelet product into the Illumination/Storage Bag.
     **b.** Allowing enough tubing length for connecting to the Riboflavin Solution Bag, permanently seal the inlet line between the Platelet Product Bag and the Y connector.
     **c.** Sterile connect the Riboflavin Solution Bag to the same inlet line you used to transfer the Platelet product.
     **d.** Open the sterile connection and drain the Riboflavin Solution into the Illumination/ Storage Bag.

1 Platelet product bag
2 Riboflavin solution bag
3 Sterile connections on inlet lines
4 Illumination/storage bag
5 Y connector
6 Pinch clamp

**Figure 4: Connecting the Riboflavin Solution Bag to the Illumination/Storage Bag**

8. To express residual air from the Illumination/Storage Bag:
   **a.** Hold the Illumination/Storage Bag so the ports are facing up, as shown in Figure 4.
   **b.** Gently apply pressure to the Illumination/Storage Bag to move the Platelets up through the inlet line tubing until they just reach the Riboflavin Solution Bag.
   **c.** Drain the Platelets from the Riboflavin Solution Bag and inlet line into the Illumination/ Storage Bag, closing the clamp on the inlet line before the air interface reaches the Illumination/Storage Bag port.

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The following figure shows acceptable and unacceptable amounts of residual air.



**Figure 5: Residual air in the illumination bag**

**1** Acceptable amount of residual air

**2** Too much residual air

**9.** Permanently seal the inlet line as near to the Illumination/Storage Bag as possible. The sealed inlet line must be less than 2.5 cm (1 inch) long, as shown in Figure 6, for it to properly fit in the Illuminator. Disconnect and discard the Riboflavin Solution Bag.

**1** Sample bulbs
**2** Spike ports
**3** Sealed inlet lines
**4** Frangible connectors
**5** Illumination/storage bags
**6** Pinch clamps
**7** PAS assembly
**8** Y connector
**9** Second platelet storage bag

**Figure 6: Single-bag (a) and DOUBLE+ (b) Illumination/Storage Sets ready for illumination**

**10.** If the Illuminator has been configured with a tare weight for a DOUBLE+ set with some components (the PAS Assembly, second Platelet Storage Bag, or both) removed, permanently seal the lines to these components and remove them now.

**11.** Perform the illumination procedure according to the Mirasol Illuminator operator's manual.

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12. After performing the illumination procedure, break the frangible connector on the spike port. You will break the second frangible connector when access is needed to add PAS, split the product, or take a sample.

**Note:** The broken parts of the frangible connectors will remain in the Illumination/Storage Bag before and after transfusion of the Platelet product.

13. Establish appropriate storage conditions according to specifications listed in Table 4. If necessary, add PAS to the Mirasol-treated Platelet product. See "Adding PAS to a Mirasol-Treated PPC Product".
    - Until the time of transfusion, protect the treated platelet product from prolonged exposure to direct sunlight or strong ambient light sources.
    - The laboratory/blood center should follow appropriate standards for product quality acceptance criteria.
14. If necessary, take a sample of the Mirasol-treated Platelet Product. See "Sampling the Mirasol-Treated Platelet Product Following PAS Addition" on page 12.

#### Adding PAS to a Mirasol-Treated PPC Product

Use only approved Platelet Additive Solutions. Solutions containing phosphate buffer are recommended for storage of Mirasol-treated platelet products.

**Note:** Platelet Additive Solution must be added within 2 hours of completion of the illumination process.

Add PAS to the Mirasol-treated PPC product according to the storage concentration specifications listed in Table 4 and according to your facility's SOP. Split product as required according to the specifications in Table 2.

To add PAS to the single bag set, complete the following steps:

1. Sterile connect the PAS bag to the sample bulb outlet line.
2. Break the frangible connector on the sample bulb outlet line.
3. Transfer the appropriate amount of PAS into the Illumination/Storage Bag. A minimum of 135 mL of PAS is recommended.
4. Permanently seal the sample bulb outlet line as near as possible to the Illumination/Storage Bag, but allowing enough tubing length for future sampling needs.

If required, split the treated product prior to storage. To split the product using the single-bag set, connect an accessory Platelet Storage Bag to the Illumination/Storage Bag sample bulb outlet line. An accessory Platelet Storage Bag (Catalog #70030) is available from the manufacturer.

To add PAS to the DOUBLE+ Set, complete the following steps:

1. Close the Sample Bag and sample bulb pinch clamps.
2. Connect the PAS bag to the Sterile Barrier Filter inlet line.
3. Transfer the appropriate amount of PAS through the Sterile Barrier Filter.
4. If desired, express residual air from the Illumination/Storage Set into the Sample Bag.
5. If desired, collect a product sample in the Sample Bag. See Sampling the Mirasol-Treated Platelet Product Following PAS Addition.

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
**Figure 7: PAS Assembly**

1. Sterile Barrier Filter inlet Line
2. Sterile Barrier Filter
3. Sample Bag
4. Pinch clamps
5. Upper Y connector
6. Lower Y connector

6. Permanently seal the PAS Addition Line between the upper and lower Y connectors, and remove and discard the PAS Assembly.
7. Gently mix the PAS and the treated Platelet product in the Illumination/Storage Bag.
8. If necessary, split the Platelet product between the Illumination/Storage Bag and the second Platelet Storage Bag.

## Sampling the Mirasol-Treated Platelet Product Following PAS Addition

You can collect a product sample using the sample bag or sample bulb on the DOUBLE + set, the sample bulb line on the single bag set, or a separate sampling device. See the Mirasol Product Sampling Guidance document for detailed sampling guidance.

- To collect a product sample in the Sample Bag after adding PAS, but before storage, complete the following steps:
  1. Gently mix the PAS and the Mirasol-treated Platelet product.
  2. Make sure the Filter pinch clamp is closed.
  3. Open the Sample Bag pinch clamp and gently apply pressure to the Illumination/Storage Bag to move Platelets up through the PAS Addition Line until the desired amount of the mixture is in the Sample Bag. If product resides in the PAS Addition Line, flush the line several times before collecting the sample.
  4. Close the Sample Bag pinch clamp.
  5. Permanently seal the PAS Addition Line between the upper and lower Y connectors and remove the PAS Assembly.

To collect a product sample in the sample bulb before or after storage, follow standard blood establishment practices. If product resides in the sample bulb line, flush the line several times before collecting the sample.



To collect a product sample in a separate sampling device, connect it to the appropriate location using a sterile connection device.

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**Table 1: Symbols and Certification (continued)**

| Symbol | Description |
|---|---|
|  | Indicates that the user should consult the instructions for use. |
| CE 2797 | Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended. |
|   | Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP). |
|  | Indicates a medical device that needs protection from light sources. |
|  | Indicates a medical device that needs to be protected from moisture. |
|  | Indicates the range of temperatures to which the device can be safely exposed. |
|  | Indicates the product quantity when the quantity is placed in the square. |
|  | Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile). |

## Return of Used Product

If for any reason the product must be returned to Terumo BCT, Inc., a returned goods authorization (RGA number) is required prior to shipping.

**It is the responsibility of the healthcare institution to adequately prepare and identify the product for return shipment.**

Shipping instructions, including an RGA number and instructions for cleaning, packaging, and labeling, may be obtained from the Terumo BCT Quality Assurance Department. Please contact your local representative for information regarding returned goods and product complaints.

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## Additional Information

**Table 2: Incoming Product Specifications***

<table>
<tr><th>Characteristic</th><th colspan="2">Specification</th></tr>
<tr><th></th><th>Standard Platelets in Plasma</th><th>Plasma-Reduced Platelets for Final Storage in PAS</th></tr>
<tr><td>Anticoagulant</td><td>Apheresis: ACD-A<br>Whole-Blood-Derived: CPD and CDPA-1</td><td>Apheresis: ACD-A<br>Whole-Blood-Derived: CPD and CDPA-1</td></tr>
<tr><td>Platelet volume (prior to Riboflavin Solution addition)</td><td>170 mL to 360 mL</td><td>90 mL to 360 mL</td></tr>
<tr><td>Platelet concentration (prior to Riboflavin Solution addition)</td><td>0.80 to 2.10 × $10^6$ Platelets/µL</td><td>1.75 to 3.40 × $10^6$ Platelets/µL</td></tr>
<tr><td rowspan="2">Treatment time constraints</td><td>Apheresis:<br>• Rest (no agitation) for a minimum of 2 hours after collection and before treatment.<br>• Treat within 2 hours to 22 hours of collection.</td><td>Apheresis:<br>• Rest (no agitation) for a minimum of 2 hours after collection and before treatment.<br>• Treat within 2 hours to 18 hours of collection.</td></tr>
<tr><td colspan="2">Whole-Blood-Derived:<br>• Treat within 8 hours of preparing the final (pooled) whole-blood-derived platelet concentrate, but not more than 32 hours after whole blood collection. This time limit assumes that leukoreduction occurs during or after the pooling step.</td></tr>
<tr><td>Residual white blood cells</td><td colspan="2">≤ 250 × $10^6$ per unit of 3.5 × $10^{11}$ yield</td></tr>
<tr><td>Gamma irradiation</td><td colspan="2">If desired, gamma irradiation may be performed after Mirasol PRT system treatment at any point during the storage period.<br>Note:<br>Mirasol PRT system treatment can be used instead of gamma irradiation to prevent transfusion-associated graft-versus-host disease (TA-GVHD).</td></tr>
<tr><td colspan="3">*The Mirasol PRT system for treating Platelets in Plasma has been verified for these incoming product conditions. Certain conditions outside of these ranges may be acceptable but should be evaluated by the user, as required.</td></tr>
</table>

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**Table 3: Disposable Kits***

| Characteristic | Specification | |
|---|---|---|
| | **Single-Bag Illumination/Storage Set** | **DOUBLE+ Illumination/Storage Set** |
| Mirasol Platelet Disposable Kit for Treatment in Plasma | Each Disposable Kit consists of one Illumination/Storage Set and one 35 mL bag of 500 µM riboflavin in 0.9% sodium chloride solution. | |
| Illumination/Storage Bag | Citrate vinyl bag | |
| Second Platelet Storage Bag | N/A | Citrate vinyl bag |
| Illumination/Storage Set sterility | The Illumination/Storage Set has a sterile and non-pyrogenic fluid pathway, sterilized by ethylene oxide (EO). | |
| Illumination/Storage Bag label | The pre-attached label covers 110 $cm^2$ (17 $in^2$) surface area. This size label must be present during illumination. Additional labels may be applied over the pre-attached label, but they may not cover any bag surface area outside the edges of the pre-attached label. | |
| Tubing dimensions for sterile connections and seals | Outer diameter = 4.06 mm ± 0.08 mm<br>Wall thickness = 0.56 mm ± 0.05 mm | |
| Illumination/Storage Bag and Second Platelet Storage Bag unrecoverable volume | 3.4 mL | |
| Sterile barrier filter membrane pore size | N/A | 0.2 µm |
| Sample Bag volume | N/A | 50 mL |
| Sample Bulb volume | 1.5 mL to 4 mL | |
| Riboflavin Solution content | 500 µM riboflavin in 0.9% sodium chloride solution<br>pH 3.5 to 5.1<br>Riboflavin Material Safety Data Sheets are available upon request. | |
| Riboflavin Solution sterility | The Riboflavin Solution is sterile and non-pyrogenic, sterilized by steam. | |
| Riboflavin Solution volume | 35 mL ± 5 mL | |
| Riboflavin Solution raw materials | Riboflavin: USP, EP<br>Sodium chloride: USP, EP<br>Hydrochloric acid: USP, EP<br>Water for injection: USP, EP | |
| Riboflavin Solution time constraints | Treat the Platelet product within 1 hour of removing the Riboflavin Solution Bag from the foil pouch. | |

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**Table 3: Disposable Kits* (continued)**

| Characteristic | Specification | |
|---|---|---|
| | **Single-Bag Illumination/Storage Set** | **DOUBLE+ Illumination/Storage Set** |
| Disposable Kit storage conditions | Long-term storage range: 2 °C to 30 °C (36 °F to 86 °F)<br>Storage humidity: Up to 75% relative humidity (RH), non-condensing<br>Permitted excursions: Up to 60 °C (140 °F) for up to 2 weeks | |
| Case quantity | 25 kits per case | |
| Case size | 24.9 cm × 35.8 cm × 35.8 cm | 24.9 cm × 43.5 cm × 35.8 cm |
| Case weight | 5.0 kg | 7.0 kg |

**Table 4: Platelet Product Storage Conditions***

| Characteristic | Specification | |
|---|---|---|
| | **Standard Platelets in Plasma** | **Plasma-Reduced Platelets for Final Storage in PAS** |
| Use of Platelet Additive Solution (PAS) | Not applicable | PAS must be added within 2 hours of completion of the illumination process. |
| Total storage volume and yield | If total volume > 400 mL and/or total Platelets per bag > $5.1 \times 10^{11}$**, treated product must be divided within 24 hours of treatment and stored in 2 separate storage bags. | |
| Storage volume per bag | 170 mL to 400 mL | 250 mL to 400 mL |
| Storage concentration | 0.70 to $2.1 \times 10^6$ Platelets/µL | 0.70 to $1.5 \times 10^6$ Platelets/µL after addition of PAS |

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**Table 4: Platelet Product Storage Conditions* (continued)**

<table>
<tr><th>Characteristic</th><th colspan="2">Specification</th></tr>
<tr><th></th><th>Standard Platelets in Plasma</th><th>Plasma-Reduced Platelets for Final Storage in PAS</th></tr>
<tr><td rowspan="3">Storage period and environment</td><td colspan="2">• After adding Riboflavin Solution and until the time of transfusion, protect the Platelet product from prolonged exposure to direct sunlight or strong ambient light sources.<br>• To ensure sufficient protection from light, store the Platelet product in an enclosed Platelet storage cabinet or on any intermediate shelf of an open tabletop Platelet agitator.<br>• Do not store the Platelet product on the top shelf of an open tabletop Platelet agitator where it will be exposed to direct ambient light or sunlight.</td></tr>
<tr><td>Store for a maximum of 5 days after collection at 22 °C ± 2 °C with gentle agitation.</td><td>Store for a maximum of 7 days after collection at 22 °C ± 2 °C with gentle agitation.</td></tr>
<tr><td colspan="2">Note:<br>Any extension of platelet storage time beyond your institution’s current limits should be evaluated per Directive [2004/33/EC] and validated according to local blood bank procedures.</td></tr>
<tr><td colspan="3">*The Mirasol PRT System for treating Platelets in Plasma has been verified for these storage conditions. Certain conditions outside of these ranges may be acceptable but should be evaluated by the user, as required.<br>**Storage of products at yields between 4.5 and 5.1 × $10^{11}$ yield for a period greater than 3 days may show evidence of a decline in some in vitro cell quality parameters. The extent of this decline may vary depending on collection and processing methods used by the site.</td></tr>
</table>

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©2016, 2021 Terumo BCT, Inc.
Terumo BCT Biotechnologies, LLC
©2016, 2019 Terumo BCT, Inc.