Medtronic

Mo.Ma Ultra™

Cerebral Protection Device
Устройство за церебрална защита
Zařízení pro cerebrální ochranu
Anordning til beskyttelse af hjernen
Zerebralprotektionssystem
Συσκευή προστασίας εγκεφάλου
Dispositivo de protección cerebral
Tserebraalkaitseseade
Aivojen suojauslaite
Dispositif de protection cérébrale
Uređaj za zaštitu mozga
Cerebrális védőeszköz
Perangkat Perlindungan Serebral
Dispositivo di protezione cerebrale
Церебральды қорғау құрылғысы
뇌 보호 장치
Galvos smegenų apsaugos įtaisas
Cerebrālā aizsargierīce
Уред за церебрална заштита
Instrument voor cerebrale bescherming
Enhet for beskyttelse av hjernen
Urządzenie zabezpieczające naczynia mózgowe
Dispositivo de proteção cerebral
Dispozitiv de protecție cerebrală
Система церебральной защиты
Zariadenie na cerebrálnu ochranu
Pripomoček za cerebralno zaščito
Uređaj za cerebralnu zaštitu
Cerebralt skydd
Serebral Koruma Cihazı
Система церебрального захисту
Thiết bị bảo vệ não

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni per l'uso • Қолдану жөніндегі нұсқаулар • 사용 지침 • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання • Hướng dẫn sử dụng

CE
2797

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
EN
1 Y-piece with hemostatic valve and flexible extension tubing (17 cm length)
2 Exit port of the working channel
3 Hollow mandrel, 0.035 in (0.89 mm) guidewire compatible, to insert inside the working channel during Mo.Ma Ultra cerebral protection device introduction
4 Shaft with 9 Fr [8 Fr] outer diameter
5 3-way stopcock
6 Balloon inflation/deflation ports with attached 1-way stopcocks (11).
7 Distal balloon (ECA), with a central radiopaque marker, occlusion range up to Ø 6 mm
8 Proximal balloon (CCA), with a central asymmetric marker, occlusion range up to Ø 13 mm
9 Distal tip valve (“fish-mouth-valve”)
10 3x filter baskets
11 Two 1-way stopcocks for connection to balloon inflation/deflation port (6).
12 30 mL self-locking syringe
13 T-Safety Connector

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1 Description

There are two different Mo.Ma Ultra cerebral protection device versions:

- 9 Fr Mo.Ma Ultra double balloon with a 6 Fr working channel compatible with 9 Fr introducer sheath
- 8 Fr Mo.Ma Ultra double balloon with a 5 Fr working channel compatible with an 8 Fr introducer sheath

All information related to the 8 Fr version will be indicated in brackets […]

The Mo.Ma Ultra cerebral protection device consists of a three-lumen shaft (one working channel and two inflation-deflation lumen) integrating two compliant balloons and a handle, which is located at the proximal end. The compliant balloons can be independently inflated: the proximal balloon in the common carotid artery (CCA) and the distal balloon in the external carotid artery (ECA). The balloons are mounted on the distal portion of the shaft with the distal balloon close to the tip and the proximal balloon at a distance of 60 mm from the distal balloon (balloon markers relative distance), just before the exit port of the working channel. The two compliant balloons can be inflated up to 13 mm (proximal CCA balloon) and 6 mm (distal ECA balloon). A hollow mandrel is included in the packaging and is intended to be inserted inside the Mo.Ma Ultra cerebral protection device's working channel to improve trackability during device insertion over a 0.035 in (0.89 mm) guidewire. The system is compatible with 9 Fr [8 Fr] introducer sheath, its usable length (distance measured from the handle to the tip of the catheter) is 950 mm and its total length (distance measured from the hemostatic valve to the tip of the catheter) is 1,110 mm (please refer to the device drawing in *Figure 1*).

The working channel is a 0.083 in (2.10 mm) [0.069 in (1.76 mm)] internal diameter fully usable lumen, with a length of 1,045 mm (distance measured from the hemostatic valve to the exit port), which serves as guiding catheter and as blood aspiration lumen. The exit-port of the working channel is located between the two balloons allowing easy access to the internal carotid artery (ICA).

Two radiopaque markers, one for each balloon, allow correct and prompt device positioning. In addition the asymmetry of the proximal marker allows correct and easy orientation of the exit-port of the working channel towards the ICA in order to facilitate guidewire introduction.

A one way tip-valve located at the distal end of the device prevents eventual backflow from the ipsilateral ECA to the brain through the ICA.

2 Contents

One Mo.Ma Ultra cerebral protection device, one hollow mandrel, one hemostatic valve with 3-way stopcock and extension line, three 40 µm cell filters, one 30 mL syringe with fix male Luer for the purging procedure, one T-safety connector, and two 1-way stopcocks.

3 Intended purpose

The intended purpose of the Mo.Ma Ultra cerebral protection device is to capture and remove debris that becomes dislodged during interventional procedures in the carotid arteries.

3.1 Intended patient population

The Mo.Ma Ultra cerebral protection device is intended for patients in need of embolic protection during carotid angioplasty and/or carotid stenting procedures.

This device is not intended for pediatric use.

3.2 Indications for use

The Mo.Ma Ultra cerebral protection device is indicated to be used in patients eligible for carotid angioplasty and/or stenting with stenosis involving the ICA and/or the carotid bifurcation and reference diameter of ECA from 3 to 6 mm and reference diameter of CCA from 5 to 13 mm.

3.3 Contraindications

The Mo.Ma Ultra cerebral protection device is contraindicated for use in patients showing one of the following criteria

- Severe disease of the ipsilateral CCA
- Intracranial tumors, aneurysms or severe intracranial stenosis distal to the target lesion
- Relevant acute neurological events occurring in the last 5 days before scheduled intervention and diagnosed by neurological assessment, such as computed tomography (CT) or cranial magnetic resonance imaging (MRI)
- Inability to respond to external questions and stimuli and to exert a pressure with the contralateral hand
- Patient suffering from dementia

Procedure related

- Severe chronic renal failure (Creatinine values >2.5 mg/dL)
- Contraindication for antiplatelets and/or anticoagulation therapy

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- Allergy to contrast media
- Severe peripheral vascular disease preventing femoral access, hemorrhagic or hyper-coagulable status and/or inability to obtain hemostasis at the site of the femoral puncture
- Inability to accept a temporary pacemaker

3.4 Clinical benefits

The Mo.Ma Ultra cerebral protection device is indicated to be used in patients eligible for carotid angioplasty and/or stenting with stenosis involving the ICA and/or the carotid bifurcation and a reference diameter of the ECA from 3 mm to 6 mm and a reference diameter of the CCA from 5 mm to 13 mm.

The Mo.Ma Ultra cerebral protection device is a proximal cerebral protection device with a double-occlusion balloon system designed to contain and remove all sizes of debris that can dislodge during interventional procedures in the carotid arteries. The Mo.Ma Ultra cerebral protection device allows for proximal embolic protection to be established prior to crossing a carotid lesion. The positive impact of the Mo.Ma Ultra cerebral protection device is to provide protection from stroke.

Note: After the EUDAMED website is launched, the Summary of Safety and Clinical Performance (SSCP) can be found at https://ec.europa.eu/tools/eudamed using the Basic UDI-DI: 0763000B000050179.

Table 1. The Mo.Ma Ultra cerebral protection device clinical benefits

Clinical benefit	**Study rate (n/N) [95% CI][a]**
Reduction in incidence of major and minor stroke and TIA[b] within 30 days of procedure	
Published study showing the benefit	
Cremonesi 2015	1% (4/393) [0.3, 2.6]
Dyuzhikov 2013	0% (0/170) [0.0, 2.1]
Montorsi 2016	0% (0/61) [0.0, 5.9]
Montorsi 2020	1.9% (1/52) [0.0, 10.3]
Palombo 2016	0% (0/35) [0.0, 10.0]
Lee 2015	0% (0/32) [0.0, 10.9]
Bijuklic 2012	0% (0/31) [0.0, 11.2]
Technical success (ability to use the cerebral protection device as well as conjunctive devices when described)	
Published study showing the benefit	
Cremonesi 2015	100% (393/393) [87.2, 100]
Montorsi 2020	100% (52/52) [93.2, 100.0]
Palombo 2016	97.1% (34/35) [85.1, 99.9]
Lee 2015	93.8% (30/32) [79.2, 99.2]
Atchaneeyasakul 2016	100% (27/27) [87.2, 100.0]

[a] N = total number of patients; n = number of patients with stroke or technical success; CI = confidence interval

[b] TIA = transient ischemic attack

3.5 Intended users

The Mo.Ma Ultra cerebral protection device should only be used by physicians who have adequate training and experience performing stenting and/or angioplasty of the carotid arteries.

This device should only be used in a clinical setting in a sterile environment.

3.6 Device performance characteristics

The Mo.Ma Ultra cerebral protection device allows for proximal embolic protection. The large diameter working channel with the exit port serves as a guide catheter enabling use of interventional devices for aspiration of debris. The compliant balloons allow for inflation in the ECA and CCA . When inflated, the balloons establish blood-flow suspension, including antegrade flow from the CCA and retrograde flow from the ECA. The spacing of the balloons as marked by radiopaque markers, and the asymmetry of the proximal radiopaque marker, allow correct and easy orientation of the exit port of the working channel toward the ICA in order to facilitate access to the lesion. A one-way tip valve located at the distal end of the device prevents eventual backflow from the ipsilateral ECA to the brain through the ICA.

4 Warnings

- The following substances are known carcinogenic, mutagenic, or toxic for reproduction (CMR) or endocrine disrupting materials and substances to which patients can be exposed and have been identified to be present above the > 0.1% weight / weight threshold:

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 - Di(2-ethylhexyl)phthalate (DEHP, CAS 117-81-7) has been found to be in the device's clear extension tubing at a proportion of > 12% weight / weight.
 - 1-Methyl-2-Pyrrolidone (CAS 872-50-4) has been found to be in the polyamide balloon inflation lumen at a proportion of > 0.1% weight / weight.
- The Mo.Ma Ultra cerebral protection device contains Nitinol. Nitinol is a nickel-titanium alloy which can cause an allergic reaction in some people. However, because the substance does not contact the blood in procedures involving the Mo.Ma Ultra cerebral protection device, the risk of an allergic reaction is minimal.
- Perform balloon purging as described in this Instruction for Use, before inserting the Mo.Ma Ultra cerebral protection device inside the patient.
- Avoid the positioning of the Mo.Ma Ultra cerebral protection device without the hollow mandrel.
- Exercise care during handling and avoid acute bends of the device before and during balloon purging procedure.
- When inflating the occlusive balloon(s), inflation control must be made by angiographic visual estimation of balloon cylindrical shape deformation (not by pressure!).
- After inflating the balloon(s) immediately perform angiographic check of blood flow blockage as described in the present Instruction for Use, and check patient's clinical tolerance to occlusion right after.
- Should patient intolerance to occlusion occur during the procedure, immediately remove all debris performing syringe blood aspirations and deflate the proximal (CCA) balloon right after.
- Before deflating the occlusion balloon(s) always verify that no more debris are retrieved in the aspirated blood.
- Appropriate antiplatelet and/or anticoagulant therapy should be administered to the patient, as determined by the physician in accordance with standard protocols for carotid stenting. ACT should be maintained > 250 s through the intervention.
- ICA lesion must not be crossed by guidewires or any other interventional catheters before the balloon(s) have been inflated and before having checked that blood flow has been effectively blocked.
- When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation.
- Do not place the occlusion balloon(s) into highly calcified vessel segments of the carotid vessels.
- Do not manipulate the Mo.Ma Ultra cerebral protection device in an inflated state.
- If resistance occurs during manipulation do not force or continue. The reason of resistance must first be ascertained by fluoroscopy, road mapping or DSA before the device is moved backwards or forwards.
- Use only a mixture of contrast medium and saline solution to fill the balloon(s) (50/50 - 30/70). Never use air or any gaseous medium or pure contrast dye to inflate the balloon(s).
- Do not use with Lipiodol or Ethiodol contrast media, or other such contrast media, which incorporate the components of these agents.
- Do not expose the Mo.Ma Ultra cerebral protection device to organic solvents, e.g. alcohol, acetone.
- Use the catheter by date specified on the package label.

5 Precautions

- Carefully inspect the sterile package and the device before use. Do not use the device if the packaging or the device is damaged.
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Allergic reactions to contrast medium should be identified before percutaneous transluminal angioplasty.
- The general technical requirements for catheter insertion must be observed at all times. This includes purging the balloons, flushing the components with sterile, isotonic saline solution prior to use and the usual prophylactic, systemic heparinization.
- Catheter applications vary and the technique must be selected on the basis of the patient's condition and the experience of the interventionalist.
- Confirm the compatibility of other devices (guidewires, balloon dilatation catheters, stent delivery systems, etc.) with the Mo.Ma Ultra cerebral protection device before use.

6 Potential complications / adverse effects

The complications that may result from a carotid balloon dilatation and stenting procedure, aided by a proximal flow blockage cerebral protection device, include but are not limited to:

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- Arterio-venous fistula
- Pseudoaneurysm
- Local infections

Procedure related:

- Bradycardia
- Hypotension
- Carotid artery spasm
- Dissection of the carotid arteries
- Air emboli
- Cerebrovascular accident (Stroke (ischemic, hemorrhagic), TIA)
- Acute myocardial infarction (unstable angina)
- Intravascular stent migration

Angiography related:

- Hyper-/ Hypotension
- Pain and tenderness
- Arrhythmias
- Sepsis/infection
- Systemic embolization
- Endocarditis
- Short-term hemodynamic deterioration
- Death
- Drug reactions
- Allergic reaction to contrast medium
- Pyrogenic reaction

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

7 Preparation technique

1. Prior to use, carefully examine the entire Mo.Ma Ultra cerebral protection device for defects or damage during shipment.
2. Connect the hemostatic working channel valve to the proximal port of the connector (1).
3. Flush the working channel with saline solution through the hemostatic valve and extension tubing.
4. Flush the hollow mandrel with saline solution.
5. Insert the tip of the mandrel through the hemostatic valve and into the working channel. Advance the mandrel completely.
6. Close the hemostatic valve (by rotating the upper cap) in order to block the mandrel in the correct position.

7.1 Purging

Prepare the supplied 30 mL syringe filled with 10 mL contrast/saline mix (50/50 - 30/70).

1. Prepare in sequence both the ECA distal balloon inflation line and the CCA proximal balloon inflation line as follows: screw the syringe to the inflation line.
2. Keeping the syringe downward, draw vacuum and wait for about 15 seconds, until no air bubbles raise from the device.
3. Keeping the tip downward get back to ambient pressure (release the syringe plunger).
4. Detach the syringe from the device. Repeat the previous steps for the remaining inflation line.

Make sure that a drop of solution raise from the luer port as proof of purging correctly done.

The device is now purged

7.2 Device insertion, positioning, orientation

1. After a standard retrograde, femoral approach, dilate the puncture site for insertion of a 9 Fr [8 Fr] (about 25 cm) introducer sheath.
2. Under fluoroscopic control, using standard diagnostic technique, engage the ECA with a 0.035 in (0.89 mm) super-stiff guidewire.

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3. Carefully introduce the proximal section of the 0.035 in (0.89 mm) guidewire into the white distal tip valve (9.“fish-mouth-valve”), which is mounted on the distal tip of the Mo.Ma Ultra cerebral protection device. After engaging the distal tip of the mandrel, continue advancing the guidewire all through the mandrel length.
4. Make sure that the guidewire enters the tip of the hollow mandrel, by which the working channel is occluded.
5. Advance the Mo.Ma Ultra cerebral protection device (including mandrel) under fluoroscopy over the wire through the 9 Fr [8 Fr] introducer until the distal marker (on the ECA balloon) is positioned in the ECA close to the carotid bifurcation and max. 1.5 cm distally from the bifurcation ostium.

 Warning: Place the ECA balloon sufficiently deep into the ECA, so that the distal balloon (ECA balloon) cannot slip back into the CCA. However, do not place the ECA balloon deeper than 1.5 cm (measured from the ostium into the ECA) to avoid CCA proximal balloon interference with the stent placement.
6. The proper orientation of the working channel exit port towards the ICA can be checked using the central asymmetric marker on the proximal balloon. The marker should point towards the ICA ostium.
7. Once the Mo.Ma Ultra cerebral protection device is in place, open the hemostatic valve, remove the mandrel and, as soon as back bleeding occurs, close the hemostatic valve, leaving the 0.035 in (0.89 mm) guidewire in place.

8 Balloon inflation, flow blockage

Caution: Use only the supplied 30 mL syringe with fix male Luer and T-safety connector for inflation of the balloons. The proper inflation is NOT pressure depending, but to be angiographically controlled.

Do not use the Mo.Ma Ultra cerebral protection device in ECA reference diameters > 6 mm and in CCA reference diameters > 13 mm.

1. Prepare the supplied 30 mL syringe filled with contrast/saline mix (50/50 – 30/70) and connect it to T-safety connector.
2. Join together the T-safety connector with the one-way stopcock.
3. Connect the assembly to the syringe and, keeping the tip upward flush it.
4. Ensure that no air bubbles remain in the syringe and in the T-safety connector.
5. Join the assembly with the Mo.Ma Ultra cerebral protection device inflation line previously prepared.
6. Screw the assembly and make it secure connected with the device. Screwing, handle the assembly by the one way stopcock.
7. With the one way stopcock open gently inflate the balloon.
8. In case of overpressure from the T-safety connector few drops of water will escape.
9. As soon as balloon inflation is completed, turn the one way stopcock from the present open-position 90 degrees to the closed-position.
10. With the one way stopcock closed, detach the syringe and the T-safety connector.
11. Inject 5 mL contrast dye through the working channel to ensure the proper occlusion of the ECA.
12. After proper positioning of the distal balloon in the ECA, retrieve the 0.035 in (0.89 mm) guidewire from the ECA.
13. Repeat the same inflation procedure from number 1 to 11 for the proximal (CCA) balloon.

9 Check of flow blockage

1. After the balloon(s) have been inflated, connect a syringe filled with contrast/saline mix (50/50 – 30/70) to the side-port of the working channel.
2. First try to aspirate slowly 10 mL of blood into the syringe.

 Note: This is a prophylactic safety test, checking for the ability to aspirate blood after the postdilatation procedure from the ICA. In case it will not be possible to aspirate blood through the working channel, the reason in most of the cases will be the suction of the vessel wall onto the working channel exit port.
3. Under fluoroscopic control slowly inject 5 mL of contrast/saline mix and verify that contrast stagnation occurs at the carotid bifurcation.

 Note: If the total flow blockage is achieved in the CCA and ICA, the contrast medium will stagnate in the occluded segment and can be used as a kind of “road map” for the adequate positioning of any stent or dilation balloon catheter.

 Caution: Ensure that the injection of contrast medium with a maximum of 5 mL is performed in a slow and controlled way. Do not inject ANY more liquids during the following procedure.
4. Now the treatment of the carotid artery can be performed under full cerebral protection.

10 Check of patient tolerance

- Perform a continuous, very close monitoring of the patient, when initiating the CCA occlusion. Observe the patient and talk to him or her in order to verify his mental and physical tolerance of the unilateral flow blockage. A whistling toy should be available to verify patient’s ability to squeeze it with his or her contralateral hand.
- Control and maintain patient systemic blood pressure in order to avoid hypotension and hypertension.

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- Should the patient manifest immediate or delayed intolerance after inflation of the balloons, the operator must first aspirate blood until no more debris appear on the microporus filter to prevent cerebral embolism immediately followed by deflation of the proximal balloon (CCA) in order to restore cerebral blood flow. Please refer to the section “DEBRIS REMOVAL”.
- At this point operator should wait a couple minutes until the patient has recovered from symptoms. Then re-inflate the proximal balloon in order to check persistence of intolerance. In this case it is the operator’s responsibility to decide either to interrupt the procedure or to use an alternative protection device (i.e. Filter).

Caution: Should the patient not tolerate the occlusion of the CCA, the operator must first aspirate blood to prevent cerebral embolism immediately followed by deflation of the proximal balloon (CCA) in order to restore cerebral blood flow.

11 Carotid treatment procedure

- To treat the diseased ICA, the operator can use all material of choice, which is compatible with 6 Fr [5 Fr] introducer sheaths. The working channel of the Mo.Ma Ultra cerebral protection device has an inner fully usable lumen diameter of 0.083 in (2.10 mm) [0.069 in (1.76 mm)] and is used as guiding catheter or sheath for the advancement and retrieval of guidewires, pre- and postdilatation balloon catheters and carotid stent delivery systems.
- Carotid angioplasty and stenting are to be performed according to standard technique and are left to operator’s experience and skills.
- It is mandatory to leave the ICA guidewire across the treated lesion, while the Mo.Ma Ultra cerebral protection device is in place.

12 Debris removal

Caution: At the end of the procedure, normally after postdilating the stented lesion, perform debris removal by blood aspiration, while the Mo.Ma Ultra cerebral protection device with fully inflated balloons is still in place.

Before performing blood aspiration verify that the hemostatic valve in the Mo.Ma Ultra cerebral protection device connector is completely closed.

1. Connect a 20 mL empty syringe to the 3-way stop cock of the side port (5.)
2. Properly close the hemostatic valve (1.)
3. Open the stop-cock and start slowly aspirating blood.

 Caution: During aspiration, carefully monitor the patient due to any occurring symptoms of intolerance.

 Note: This procedure can be done at any stage during the intervention but it is a must at least after postdilatation of the stented lesion.
4. Slowly aspirate at least 60 mL of blood (three 20 mL syringes). Check the last syringe for debris, by flushing the blood content through the supplied microporous filter. If still debris are visible in the filter basket, aspirate additional 20 mL syringes of blood, until no more debris can be detected.

Warning: If no blood can be aspirated through the working channel, carefully advance a suitable guiding catheter through the working channel, slightly distal from the working channel exit port. Aspirate now 60 mL through this guiding catheter.

13 Balloon deflation

1. Prior to deflation remove the T-safety connector from the 30 mL inflation syringe. Alternatively, any other syringe can be used for deflation.
2. When the absence of any residual debris has been ensured, deflate the distal balloon (ECA) under fluoroscopic control by connecting the syringe to the distal inflation port and apply negative pressure.
3. Deflate the proximal balloon (CCA) under fluoroscopy by connecting the syringe to the proximal inflation port and by applying negative pressure.
4. Perform control angiography of the carotid artery.

14 System withdrawal

1. Retract every guidewire out of the working channel.
2. Gently retrieve the Mo.Ma Ultra cerebral protection device under fluoroscopic control.

Caution: IMPORTANT NOTE ON BALLOON INFLATION

CCA and ECA balloons are highly compliant type and are designed to achieve atraumatic vessel occlusion by changing their shape from circular to cylindrical thus featuring a wide balloon-vessel contact area plus giving the operator an immediate angiographic feedback that vessel occlusion has been achieved. For this reason a pressure guided operation of the balloons must not be performed. Volumes needed to occlude carotid vessels depend upon vessel diameters. In vitro testing has shown that for the distal balloon, in case of the maximum considered diameter of 6 mm, the volume needed to achieve cylindrical balloon shape and thereby effective occlusion, should not exceed 0.3 mL. In the same way, for the proximal balloon, in case of the maximum considered diameter of 13 mm, the volume needed to achieve cylindrical balloon shape and thereby effective occlusion, should not exceed 1.5 mL.

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Common carotid vessel size require less volumes to achieve occlusion than the above mentioned quantities. Therefore during inflation the operator must carefully watch the fluoroscopic image in order to detect the cylindrical shape deformation of the balloons.

15 How supplied

The Mo.Ma Ultra cerebral protection device is supplied sterile (via ethylene oxide) and intended for single use only. The Mo.Ma Ultra cerebral protection device is sterilized by ethylene oxide gas. It will remain sterile as long as the packaging remains unopened and undamaged. It is supplied with a single pouch. Use product by the date printed on the label.

Caution: Do not use if the inner package is open or damaged.

16 Storage

Store the device in a dry area away from sunlight.

Do not store the device where it is directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used before the Use-by date on the package label.

17 Disposal

Dispose of the device in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

18 Disclaimer of warranty

Note: This disclaimer of warranty does not apply in Australia.

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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© 2023 Medtronic
Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States or in other countries.
Medtronic, Inc.
Medtronic B.V.