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English

Product name Intra-oral wound dressing

Intended Purpose
It is used to protect oral wound by forming a physical barrier on the wound.

Features & benefits (Performance Characteristic)
This product is non-eugenol protective dressing material applied to intraoral wounds.
1. Advanced concept in intraoral dressing
2. Protect oral wounds from food and cigarette smoke.
3. Aid in hemostasis
4. Protect suture thread from tongue irritation
5. Strong adhesion using hydrophilic polymer
6. Safe and easy to use. Easy to cut into different shapes/sizes
7. Self adhesive with saliva
8. Reduce irritation sore in mouth
9. Mint flavor refreshes the mouth after oral surgery

Indication
1. Protect wounds after oral surgery and treatment
2. Promote healing after oral surgery and treatment
3. Canker sore
4. Protect suture thread from tongue irritation
5. Protect sore from orthodontia
6. Protect fluorine varnish or sensitive teeth

How to Use
• Preparation
1. Check product contamination, packing damage, and expiration date.
2. Please carefully read the instructions before use.
(Some bubble in the product are from the manufacturing process and safe)

• Instructions
1. Irrigate wound with sterile or saline solution..
2. Cut product into a proper shape and size.
3. Slightly remove moisture with gauze on the wound. An excessive exudation may reduce its adhesive strength and its use time.
4. Remove transparent release paper and apply on wound.
5. Gently press product for 5 to 10 seconds while product adheres to wound.
(According to circumstances, please press repeatedly)

Contra-indications
1. Do not use if the product appears damaged or contaminated before use.
2. Do not use on deep wounds with excessive bleeding.

CAUTIONS
⚠ This medical device is for single use. If reused, the product may cause infection, cross infection, and inflammation. This product shall be disposed in accordance with disposal regulation of each CA for fear of infection.
1. Please carefully read the instructions before use.
2. If you experience any swelling, allergic reaction, upset stomach or any signs of infection, consult with your doctor.
3. If the product is removed forcibly, the wound may come off. The attachment side automatically falls off when it dissolved completely in the mouth.
4. Remove the dropped protection layer from the wound.
5. Pay particular attention if it is being used on young children under 5 years old.
6. The time of attachment may be reduced when stimulated by tongue or food.
• A notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established;

Storage
1. Store indoors (1℃~35℃) in a clean place.
2. Avoid fire, direct sunlight, and high temperature environment.
3. Keep out of reach of children.
• This product may be hardened according to the storage and management method. Although its adhesion lasts at least 4 hours, the attachment method and the environment may result in longer or shorter adhesion times.
• Some bubble on the surface of the strip are from the manufacturing process and are safe.

Applicable patient's profile
Patients who need a dressing for protecting oral wounds.
Pay particular attention if product is being used on young children under 5 years old.
Device user : Person 6 years old or Older
For devices intended for use by lay persons, the circumstances in which the user should consult a healthcare professional;

TBM Corporation B-301, 313, Cheomdangwagi-ro, Buk-gu, Gwangju, 61008, Korea Tel:+82-62-971-2845 Fax:+82-62-971-2815
website:www.tbmkorea.com (User manual can be downloaded from the homepage.)
EC REP Emergo Europe B.V. Prinsessegracht 20, 2514 AP The Hague, The Netherlands

Model Number OB23, OB53 **Doc. No** TF-ORA-IFU, Rev 0 **Date of Issue of the IFU** 2020-02
* This product shall be disposed according to disposal regulation of each CA.

Symbol Description

REF	Model Number	LOT	Lot number		Manufacturer
	Date of manufacture		Use by		Do not reuse
	Caution		Storage and Transport Temperature limit		Keep away from sunlight
	Consult instructions for use		Do not use if package is damaged	EC REP	Authorized representative in the European community
CE	CE mark	GMP	Certificate of GMP	UDI Barcode	UDI Barcode (Issuing Agency: GS1)

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Emergo Europe B.V. Prinsessegracht ٢٥١٤ ,٢٠ AP The Hague
The Hague ,Emergo Europe B.V. Prinsessegracht ٢٥١٤ ,٢٠ AP