Concept Medical

MagicTouch SCB

SIROLIMUS COATED BALLOON CATHETER

EN	Sirolimus Coated Balloon Catheter Instructions for Use	Page **04**
FR	Cathéter à ballon enduit de sirolimus Instructions d'utilisation	Page **10**
DE	Sirolimus beschichteter Ballonkatheter Gebrauchsanweisung	Page **18**
IT	Catetere a palloncino rivestito con Sirolimus Istruzioni per l'uso	Page **25**
ES	Catéter con Balón Recubierto de Sirolimus Instrucciones de uso	Page **33**
PT	Cateter de balão revestido com Sirolimus Instruções de uso	Page **40**
RU	Баллонный катетер, покрытый сиролимусом Инструкции по использованию	Page **47**
TR	Sirolimus Kaplı Balon Kateter Kullanım Talimatları	Page **55**

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ENGLISH

1.0 PRODUCT DESCRIPTION

MagicTouch - SCB Sirolimus Coated Balloon Catheter (MagicTouch - SCB) is a device/drug combination product consisting of two components: 1). a standard PTCA balloon catheter 2). coated with encapsulated sirolimus formulation. Generic name of MagicTouch - SCB is "Sirolimus Coated Balloon Catheter" and GMDN Code is 62218.

The coating of MagicTouch - SCB utilizes Nanolute technology to deliver polymer-free sirolimus encapsulated in a Generally Recognized as Safe (GRAS) phospholipid excipient. The amphiphilic properties of the phospholipid allow the Sirolimus to remain encapsulated and protected from degradation for a time sufficient to ensure treatment of the target lesion. Coating is on circumferential surface of the balloon length ensures homogeneous drug delivery.

The use of Nanolute technology provides low in-transit drug loss resulting in reduced shedding of particulate matter that can contribute to downstream embolization. This technology also allows for better drug diffusion and drug retention throughout the artery, as well as enhanced bioavailability due to deeper penetration of the drug into the artery tissue.

1.1 DEVICE COMPONENT DESCRIPTION

MagicTouch - SCB is indicated for gaining lumen diameter in Coronary artery disease by performing Percutaneous Transluminal Coronary Angioplasty, after appropriate lesion preparation in native coronary arteries with reference vessel diameter 1.50 mm to 4.00 mm in diameter and up to 36 mm lesion length.

The balloon catheter is a monorail (Rapid Exchange) PTCA catheter with a semi-compliant inflatable balloon at the distal part and an atraumatic, tapered tip to aid in crossing the tight stenosis. The catheter is double lumen in the distal part, one lumen is used for inflation and deflation of the balloon, and the second lumen allows access to the distal part of the catheter for guide wire insertion. The maximum recommended guide wire diameter is 0.014". A luer lock fitting (Hub) allows the connection with an inflation device. The internal tubing of the balloon has two radio opaque markers (Platinum - Iridium)to provide visual reference points for balloon positioning relative to the stenosis within the vessel. The catheter working length is 140cm and is available in diameters between 1.50 mm and 4.00mm, and lengths of 10, 15, 20, 25, 30, 35 and 40 mm as shown in Table 1.


Figure 1: Representative diagram of MagicTouch - SCB

Table 1: Device Component Description

Brand Name	MagicTouch - SCB
Available Lengths (mm)	10, 15, 20, 25, 30, 35, 40
Available Diameters (mm)	1.50, 2.00, 2.25, 2.50, 2.75, 3.00, 3.50, 4.00
Drug Component	Sirolimus Drug (also known as Rapamycin)
Excipient	Lipoid E 80 (Egg phospholipid)
Drug Dose	1.27µg/mm^2
Delivery System Workable Length	140cm (1400mm)
Balloon Inflation Pressure	Nominal Inflation Pressure: 8 bar Rated Burst Pressure: 14 bar*
Guide wire compatibility (max.)	(max.) 0.014"
Sheath Compatibility	6F (2.00mm)
Guiding Catheter Inner Diameter	6F (1.98mm) (inner lumen ≥ 0.078")
Catheter Shaft Outer Diameter	Proximal 1.95 F (0.65 +0.025−0.005 mm)/ Distal 2.67F(0.89± 0.02 mm)/ Intermediate 2.55F (0.85 ± 0.02 mm)
Catheter Configuration	Rapid Exchange (Monorail)
Radiopaque Marker Bands	90% Pt / 10% Ir alloy
Sterilization	Ethylene Oxide Gas

**Note: Do not exceed Rated Burst Pressure (RBP)*

1.2 DRUG COMPONENT DESCRIPTION

Active Pharmaceutical Ingredient (API): Sirolimus is a macrocyclic lactone produced by Streptomyces hygroscopius.

Sirolimus appearance: White to off-white powder.

Sirolimus solubility: Freely soluble in chloroform, acetone and acetonitrile and insoluble in water and melts at approximately 183°C -185°C.

Sirolimus Molecular formula: $C_{51}H_{79}NO_{13}$

Molecular weight: 914.2 g/mol

CAS Registry no.: 53123- 88-9

Chemical name: (3S, 6R, 7E, 9R, 10R, 12R, 14S, 15E, 17E, 19E, 21S, 23S, 26R, 27R, 34aS) - 9, 10, 12, 13, 14, 21, 22, 23, 24, 25, 26, 27, 32, 33, 34, 34a - hexadecahydro - 9, 27 - dihydroxy - 3 - [(1R) - 2 - [(1S, 3R, 4R) - 4 - hydroxy - 3 - methoxycyclohexyl] - 1 - methylethyl] - 10, 21 - dimethoxy - 6, 8, 12, 14, 20, 26 - hexamethyl - 23, 27 - epoxy - 3H - pyridol [2, 1 - c] [1, 4] oxaazacyclohentriacontine - 1, 5, 11, 28, 29 (4H, 6H, 31H) - pentone.


Figure 2: Structural Formula of Sirolimus

Table 2: Product Sizes, Reference and Drug Content

Product Reference	Balloon Diameter (mm)	Balloon Length (mm)	Drug Content (µg)
CMN15010	1.50	10	59.82
CMN20010	2.00	10	79.76
CMN22510	2.25	10	89.73
CMN25010	2.50	10	99.70
CMN27510	2.75	10	109.66
CMN30010	3.00	10	119.63
CMN35010	3.50	10	139.57
CMN40010	4.00	10	159.51
CMN15015	1.50	15	89.73
CMN20015	2.00	15	119.63
CMN22515	2.25	15	134.59
CMN25015	2.50	15	149.54
CMN27515	2.75	15	164.50
CMN30015	3.00	15	179.45
CMN35015	3.50	15	209.36
CMN40015	4.00	15	239.27
CMN15020	1.50	20	119.63
CMN20020	2.00	20	159.51
CMN22520	2.25	20	179.45
CMN25020	2.50	20	199.39
CMN27520	2.75	20	219.33
CMN30020	3.00	20	239.27
CMN35020	3.50	20	279.15
CMN40020	4.00	20	319.02
CMN15025	1.50	25	149.54
CMN20025	2.00	25	199.39
CMN22525	2.25	25	224.31
CMN25025	2.50	25	249.24
CMN27525	2.75	25	274.16
CMN30025	3.00	25	299.09
CMN35025	3.50	25	348.93
CMN40025	4.00	25	398.78
CMN15030	1.50	30	179.45
CMN20030	2.00	30	239.27
CMN22530	2.25	30	269.18

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2.0 INDICATION FOR USE

MagicTouch - SCB is indicated for treatment of Coronary Artery Disease by performing Percutaneous Transluminal Coronary Angioplasty, after appropriate lesion preparation. MagicTouch - SCB is therefore indicated to dilate the diseased segment(s) in coronary artery by delivering drug in coronary artery or graft with lesion lengths of 6.00 mm to 36.00 mm and reference vessel diameter from 1.50 mm to 4.00 mm.

The device is intended for use in following conditions:

- Patients with de-novo lesion;
- Patients with In-stent restenosis of previous stent implantation;
- Patients with disease in small vessels where stent implantation is either not feasible or involving higher risk of restenosis;
- Patients with bifurcation lesions.

3.0 CONTRAINDICATIONS

MagicTouch - SCB is contraindicated in the following types of patients:

- Patients with a hypersensitivity to Sirolimus drug or its structural activity related compounds.
- Patients with a known hypersensitivity to excipients with phospholipid or related origins.
- Lesions not suitable for interventional treatment.
- Patients whose diseased segment cannot be pre-dilated or prepared before drug coated balloon treatment.
- Patients who cannot receive recommended antiplatelet or anticoagulation therapy.

4.0 WARNINGS

- Contents supplied are sterilized using Ethylene Oxide (EO) process.
- Ensure that the product package is intact, and the inner package housing device has not been opened or damaged as this may indicate the sterile barrier has been broken.
- Do not use the product if sterile barrier is opened or opened prior to intended use.
- Inspect the product before use for any damage during removal from the protective catheter holder and protective cap. Return to company using proper handling precautions if found to be damaged for further inspection and analysis.
- The use of product carries the risks associated with coronary artery interventions, including subacute thrombosis, vascular complications or bleeding events.
- Patients with known hypersensitivity to Sirolimus drug and excipient (Phospholipids) may suffer an allergic response to this device.
- Patients who are unlikely to comply with the recommended Antiplatelet therapy should not receive this product.
- The product shall be discarded in a proper way after use as it may have residual drug (Sirolimus) on its surface.
- MagicTouch - SCB is intended for single use only.
- **Do not** re-sterilize and re-use. Reuse or re-sterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Reuse or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness and death. Concept Medical is not be responsible for any direct, incidental or consequential damages resulting from re-sterilization or reuse.
- In case of procedures where the patient requires a stent implantation, the stenting must be done prior to device treatment at the target lesion.
- **Do not** exceed rated burst pressure (RBP). Use of correct pressure monitoring device is recommended to prevent over pressurization. Inflation in excess of the RBP may cause the balloon to rupture.
- **Do not** use MagicTouch - SCB after the "Use-by Date".
- When MagicTouch - SCB is in the body, it should be manipulated only under sufficient and/or high-quality fluoroscopy. Prior to withdrawing the catheter from the lesion, the balloon must be completely deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- **Do not** use air or any gaseous medium to inflate MagicTouch - SCB. Use only the recommended inflation medium. Handle the product with caution to avoid any damage to the balloon coating or folded balloon. Avoid touching the balloon while loading it on guidewire. Always use dry Sterile gloves while loading MagicTouch - SCB on Guidewire.
- Take appropriate precautions to prevent or reduce clotting when catheter is used. Flush or rinse ancillary product, if needed (except MagicTouch - SCB) entering the vascular system with sterile isotonic saline or a similar solution via the guide wire access port prior to use. Consider the use of systemic heparinization.

5.0 PRECAUTIONS FOR USE

5.1 GENERAL PRECAUTIONS:

- Only physicians who have received adequate training in Interventional cardiology shall perform this treatment.
- The treatment shall be performed at hospitals where emergency coronary bypass graft surgery can be readily performed.
- Prior to Percutaneous Coronary Intervention (PCI), if a dental or surgical procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventional cardiologist and patient should carefully consider whether a drug coated balloon and its associated recommended antiplatelet therapy is the appropriate PCI treatment choice.
- Do not expose system to organic solvents such as alcohol, saline and detergents as it may impair with product performance in treatment.
- Preparation of the lesion i.e. **Pre-dilatation of the lesion before the treatment with MagicTouch - SCB is mandatory.**
- To minimize the possible introduction of air into the system, it's imperative that, before proceeding, careful attention is paid to the maintenance of tight catheter connections and through aspiration and flushing of the system.
- Never advance the Drug delivery catheter without the guide wire extending from the tip.
- Never attempt to move the guide wire when the balloon is inflated. Do not advance the catheter against significant resistance. The cause of resistance should be determined via fluoroscopy and appropriate remedial action taken.
- MagicTouch - SCB is not intended for stent expansion.
- MagicTouch - SCB is not intended for pre-dilation of the lesion.
- Use of MagicTouch - SCB in patients and lesions outside intended use (off label) is not recommended.
- In case of procedures where the patient requires a stent implantation, the Stenting must be done prior to MagicTouch - SCB treatment at the target lesion.

 Note: The use of drug coated balloons for treatment of coronary artery disease conditions is a relatively new concept and studies have been performed to evaluate results in In-stent restenosis, small vessels and de novo lesions with data comparable with use of stents. Drug release using MagicTouch - SCB provides longer drug retention and bio-available and thereby prevent re-occurrence. However, where the product is used outside the specific indications, the outcomes may differ from the previously documented studies with drug coated balloons.

 Using drug coated balloon treatment in patients and lesions outside labelled indications may have an increased risk of adverse events, including thrombosis, embolization, dissection, myocardial infarction or death.

5.2 HANDLING PRECAUTIONS

- Special care must be taken not to handle or in any way disrupt the coating on balloon. This is most important while removing

the catheter from the packaging, placing it over the guide wire and advancing it through large bore rotating haemostatic valve and guiding catheter hub.

- Balloon manipulation, trying to refold may damage the coating, contamination or fragmentation of coating on the balloon.

5.3 DRUG INTERACTION

Several drugs are known to affect metabolism of Sirolimus, and other drug interactions may be inferred from known metabolic effects. Sirolimus is known to be a substrate for both cytochrome P450 IIIA4 (CYP3A4) and P-glycoprotein as per published papers. Consideration should be given to the potential drug interaction when deciding for MagicTouch - SCB treatment in a patient who is taking a drug that could interact with Sirolimus. The effect of drug interactions on safety and effectiveness of MagicTouch - SCB has not been determined.

5.4 IMMUNE SUPPRESION POTENTIAL

Sirolimus, the active pharmaceutical ingredient of MagicTouch - SCB, is an immunosuppressive agent that is also available in oral formulations. When multiple MagicTouch - SCB treatment is performed in patients, systemic concentration of Sirolimus may reach immunosuppressive levels momentarily, especially in patients who also have hepatic insufficiency or who are taking drugs that inhibit CYP3A4 or P-glycoprotein. This possibility should be considered for such patients, particularly if they are also taking oral Sirolimus, other immunosuppressive drugs or are otherwise at risk of immunosuppression.

5.5 MULTIPLE INFLATION OF DEVICE

- Multiple inflation of MagicTouch - SCB has not been evaluated in humans.
- It is recommended that treatment using **single inflation of device with 60 seconds** at nominal pressure **or in case of multiple inflations total inflation time should be summing up to 60 seconds** duration and is suggested at nominal pressure without removing the balloon catheter during the procedure.

5.6 BRACHYTHERAPY

The safety and effectiveness of MagicTouch - SCB in patients with previous brachy therapy has not been established. Vascular brachy therapy is being studied for treatment of in-stent restenosis in stents; however, there is no specific evidence of safety and efficacy of treatment. MagicTouch - SCB and vascular brachy therapy alters the arterial biology and combined vascular response of these two therapies have not been established.

5.7 MAGICTOUCH - SCB INFLATION PRECAUTIONS

- MagicTouch - SCB is a drug delivery device to deliver drug upon appropriate inflation to the vessel.
- The vessel shall be pre-dilated with an appropriately sized balloon.
- Do not prepare or pre-inflate the balloon prior to treatment other than as directed. Use the balloon purging technique.
- Guiding catheters should have lumen sizes that are suitable to accommodate the balloon across the lesion. Wrong sizing may lead to damage of balloon coating.
- Inflating the balloon above its Rated burst pressure (RBP) may cause rupture in the target lesion and may lead to procedure requiring surgical treatment. Always maintain inflation pressure below the RBP.
- Inflating a balloon in a tight lesion may lead to dissection of the vessel distal/proximal to the lesion and may cause acute closure of the vessel leading to additional intervention (CABG, further dilatation, placement of additional stents or other intervention).
- Do not expand the balloon if it is not properly positioned at the target lesion.
- Expansion of balloon has a potential to compromise side branch flow, hence caution needs to be taken when inflation is done.
- When treating multiple lesions, the distal lesion should be stented initially if required, followed by MagicTouch - SCB treatment. This can be followed by proximal lesion treatment. Procedure in this order reduces the chances of disrupting the coating surface of balloon catheter.
- Inflation pressure of balloon should always be monitored. Do not exceed Rated burst pressure as listed in the compliance chart on label. Use of pressures higher than those specified on the product label may result in ruptured balloon with possible intima damage and dissection.
- Do not attempt to pull a partially inflated or inflated balloon back through the guiding catheter, as serious injury to intima layer may occur. Remove the device after complete deflation through the catheter. Use proper retrieval methods for the device as it may cause trauma at the insertion site and may lead to bleeding, haematoma or pseudo aneurysm.
- **Ensure treatment device length more than or equal to lesion length with addition of 2 mm on proximal and distal end for complete coverage**.

5.8 MAGICTOUCH - SCB REMOVAL PRECAUTIONS

- Should any unusual resistance be felt at any time during either lesion access or removal of the balloon, the entire system should be removed.
- While removing the system, advance the guidewire into the coronary anatomy as far as distally possible. Tighten the rotating haemostasis valve to secure the MagicTouch - SCB with the guiding catheter; then removing the Balloon catheter system and guiding catheter as a single unit.
- Failure to follow this procedure or to apply excessive force during withdrawal can potentially result in damage to the anatomy and cause any other vascular complication.
- If it is necessary, retain the guidewire in position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.
- Note: After use, this product may be a biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable government regulations.

5.9 POST PROCEDURE MEDICATIONS AND PRECAUTIONS

- Care must be taken when the MagicTouch - SCB is being withdrawn after treatment. Post procedure treatment of Dual antiplatelet therapy (DAPT) is recommended for a period of minimum 3 months. Patients who require early discontinuation of antiplatelet therapy (e.g. secondary to active bleeding) should be monitored carefully.

6.0 DRUG INFORMATION

6.1 MECHANISM OF ACTION

Sirolimus is a potent inhibitor of T-Lymphocyte activation, Smooth Muscle Cell and Endothelial cell proliferation in response to stimulation by cytokines, growth factors and antigens. In cells, Sirolimus binds to the immunophilin, FK Binding protein-12 (FKBP-12). Sirolimus- FKBP-12 complex binds to and inhibits the activation of the Mammalian Target of Rapamycin (mTOR), a key regulatory kinase. This inhibition suppresses cytokine-driven cell proliferation, inhibiting the progress of the cell cycle from the G_1 to S phase.

6.2 PHARMACOKINETIC OF MAGICTOUCH - SCB

In a recent GLP animal study performed by Dr. Renu Virmani's laboratory (CV Path, Gaithersburg, MD) MagicTouch - SCB sirolimus pharmacokinetic was quantified over a period of time in coronary artery segment (treated, untreated proximal and untreated distal), non-target tissues, myocardium segment, nontarget organs and whole blood.

6.2.1 TREATED CORONARY ARTERY SEGMENTS

In the **treated coronary artery segments**, the C_{max} of Sirolimus was measured at 1 hour after balloon deployment (mean ± standard deviation: 1959 ng/g ± 641 ng/g). Two weeks after balloon deployment, the Sirolimus concentration had dropped to 168.1± 76.9 ng/g (mean ± standard deviation). Thirty (30) days after balloon deployment, Sirolimus concentrations were still above the lower limit of quantification in all (6 out of 6) coronary artery samples with a mean concentration of 47.7± 38.3 ng/g. One hundred and twenty (120) days after balloon deployment, Sirolimus was still quantifiable in 1 out of 6 treated coronary artery samples (4.0 ng/g).


Figure 3: Pharmacokinetic analysis of mean Sirolimus concentration-time data

6.2.2 WHOLE BLOOD

Sirolimus was quantifiable in EDTA **whole blood** as early as 3 minutes after deployment of the Sirolimus-coated balloons. The maximum concentration (C_{max}) in EDTA whole blood of a mean 29.8 ng/mL was reached after 15 minutes. Twenty-four hours after balloon deployment, the mean EDTA whole blood concentration had dropped to 1.70 ± 0.93 ng/ml. No quantifiable Sirolimus concentrations were found in any of the

animals after 30 days and above following balloon deployment. The Sirolimus apparent half-life was a mean 5.7 hours and the mean residence time was 3.7 hours. The $AUC_{0\text{-}T}$ over the observation period was a mean 197.8 ng/ml·h.

Table 3: Pharmacokinetic analysis of Sirolimus concentration-time data.

$T_{1/2}$	T_{max}	C_{max}	T_{last}	C_{last}	$AUC_{0\text{-}T}$	MRT
[h]	[h]	[ng/mL]	[h]	[ng/mL]	[ng/mL·h]	[h]
5.7	0.25	29.8	24	1.70	197.8	3.7

Abbreviations: $T_{1/2}$: apparent half-life, T_{max}: time to maximum concentration, C_{max}: maximum concentration, T_{last}: time to last quantifiable concentration, C_{last}: last quantifiable concentration, $AUC_{0\text{-}T}$ area-under-the-concentration time curve over the observation period, MRT: mean residence time


Figure 4: Sirolimus concentration in EDTA whole blood concentration-time data.

6.2.3 UNTREATED CORONARY ARTERY SEGMENTS

The Sirolimus concentrations in the **untreated proximal coronary artery segments** were markedly lower than those in the treated artery tissues. The C_{max} of a mean (± standard deviation) 192.7 ng/g ± 115.4 ng/g was reached after 1 hour. Sirolimus was still quantifiable in the proximal untreated coronary arteries in 4 out of 5 proximal coronary artery samples at 30 days after balloon deployment (mean ± standard deviation: 30.9 ng/g ± 25.5 ng/g). Sirolimus was still quantifiable in 2 out of 5 of the untreated proximal coronary artery segments after 120 days (7.5 ng/g and 9.3 ng/g).

Sirolimus in the **untreated distal coronary artery segments** was reached after 1 hour and the mean concentration was 1107 ng/g (± standard deviation: 650.8 ng/g). After 30 days, Sirolimus was still quantifiable in 1 out of 5 untreated distal coronary artery segments (19.4 ng/g). There was no quantifiable Sirolimus concentration in any of these samples on study days 60 and later.

6.2.4 MYOCARDIUM SEGMENTS

In all **myocardium segments**, the highest C_{max} of Sirolimus was measured at 1 hour after balloon deployment. The maximum concentrations in the myocardium samples ranged from a mean 2541 ng/g ± 745.2 ng/g in the samples collected adjacent to the left anterior descending coronary artery (LAD) to 520.8 ng/g ± 866.6 ng/g in the samples collected in the apical posterior region. In the samples collected adjacent to the LAD, after one (1) day, the mean Sirolimus tissue concentration had dropped to 249.2 ng/g ± 121.3 ng/g and after three (3) days to a mean 17.6 ng/g ± 1.4 ng/g. The Sirolimus concentrations in the other myocardium regions followed a similar time-concentration pattern. After 30 and 60 days, no Sirolimus could be quantified in any of the myocardium samples.

6.2.5 NON-TARGET ORGANS

In the **non-target tissues**, the highest Sirolimus concentration was found in the lung with 459.4 ng/g after 1 hour after balloon deployment. The organ with the second highest concentrations was the spleen with 125.4 ng/g after 1 hour. The highest Sirolimus concentration in the liver was 81.5 ng/g, also reached after 1 hour. Among the non-target tissues, the maximum Sirolimus concentrate was lowest in the bone (rib) with the peak concentration at 5.1 ng/g. In all of the non-target tissues, Sirolimus was not quantifiable any longer on day 30 and later.


Figure 5: Pharmacokinetic analysis of mean Sirolimus concentration-time data in non-target organs

Results from a pharmacokinetic study evaluating the MagicTouch - SCB showed a rapid decrease in blood concentration of sirolimus over the first 24 hours after treatment, while tissue concentrations of sirolimus remained detectable for much longer up to 120 days. Drug availability of long duration allows prevention of restenosis in treated segments.

6.3 DRUG INTERACTION FOLLOWING TREATMENT OF MAGICTOUCH - SCB

Drug interaction studies have not been performed with MagicTouch - SCB because of limited systemic exposure to Sirolimus. However, consideration should be given to the potential for both local and systemic drug interactions in the vessel wall when deciding to implant MagicTouch - SCB in a subject taking a drug with the known interaction with Sirolimus.

Sirolimus is metabolized by cytochrome P450 3A4 (CYP3A4) in the gut wall and liver and undergoes efflux from enterocytes of the small intestine by p-glycoprotein (P-gp). Absorption and subsequent elimination of systemically absorbed Sirolimus may be influenced by drugs that affect CYP3A4 and P-gp.

Medications that are strong inhibitors of CYP3A4 and P-gp may increase Sirolimus levels while inducers of CYP3A4 and P-gp may reduce Sirolimus metabolism in-vivo.

Sirolimus when prescribed as oral medication may interact with the drugs/foods listed below:

- P-gp Inhibitors such as but not limited to Cyclosporine, Digoxin
- Antibiotics such as but not limited to Ciprofloxacin, Ofloxacin
- Glucocorticoids
- CYP3A4 Inhibitors such as but not limited to Ketoconazole, Fluconazole
- CYP3A4 Inducers such as but not limited to Rifampin, Phenytoin, Dexamethasone
- Herbal preparations (St. John's wort- Hypericum perforatum)
- Amphotericin B such as but not limited to Abelcet, Amphocin
- Cimetidine such as but not limited to Tagamet
- vitamins
- Vaccination- sirolimus is potent immunosuppressant hence it may affect response to vaccination making the vaccination less effective. For some period after receiving MagicTouch - SCB, use of live vaccines such as, but not limited to measles, mumps, BCG, yellow fever, varicella, Ty21a Typhoid should be avoided.
- Grapefruit/grapefruit juice: Grapefruit/grapefruit juice reduces CYP3A4 mediated metabolism of sirolimus.
- Medicine to lower your cholesterol or triglycerides
- Medicine for high blood pressure or heart problems
- Anti-seizure medicine
- Medicines used to treat stomach acid, ulcers, or other gastrointestinal problems
- Drugs that may increase sirolimus blood concentration:
 - Antifungal agents: such as but not limited to Itraconazole, Clotrimazole
 - Gastrointestinal prokinetic agents such as but not limited to metoclopramide
 - Macrolide antibiotics: such as but not limited to Troleandomycin
 - Calcium channel blockers: such as but not limited to Diltiazem, verapamil
 - Other drugs: Danazol, HIV Protease Inhibitor (Indinavir)
- Drugs that may decrease sirolimus level:
 - Anticonvulsants such as but not limited to Phenobarbitol
 - Antibiotics such as but not limited to Rifabutin

6.4 MUTAGENESIS, CARCINOGENICITY AND REPRODUCTIVE TOXICOLOGY

- In vivo micronucleus test, Chromosomal aberration test and Bacterial reverse mutation (Ames) test conducted on demonstrated that it is non-mutagenic.
- Carcinogenicity and Reproductive Toxicology: No long-term studies have been performed to evaluate the carcinogenic or genotoxic potential of the drug sirolimus coated balloons and there are no adequate and well-controlled studies published in women and men intending to have children or pregnant women. However, considering MagicTouch - SCB as non-implantable device, and limited exposure of sirolimus to vessel wall, properties of sirolimus alone was found to be non-carcinogenic. Therefore, we do not know at the present time whether sirolimus (Rapamune, Rapamycin) is safe to be used in any kind of treatment during pregnancy or in patients intending to have children.

7.0 POTENTIAL ADVERSE EVENTS

Possible complications linked to the use of the balloon catheter during the procedure:

- Dissection or perforation of the coronary artery

- Injury or rupture of the coronary artery
- Total occlusion
- Thrombosis
- Arterial spasm
- Ventricular fibrillation
- Disturbance of cardiac conductibility
- Embolism
- Infection and or pain at the access site
- Haematoma
- Drug reaction to antiplatelet agent/anticoagulation agent/contrast medium
- Emboli, distal
- Nausea & vomiting
- Vascular complication which may requires vessel repair.
- Alopecia
- Anaemia
- Blood product transfusion
- Gastrointestinal Symptoms
- Hepatic Enzyme Changes
- Histological changes in vessel wall, including inflation, cellular damage or necrosis.
- Myalgia/Arthralgia
- Peripheral Neuropathy
- These complications can directly result in the patient's death.

Possible complications that could occur following an angioplasty procedure with balloon catheter, on short and medium term:

- Restenosis of the dilated artery
- Unstable angina
- Acute myocardial infarction
- Disturbance of cardiac conductibility
- Bleeding complications or haematoma

These complications can directly result in the patient's death.

8.0 INDIVIDUALIZATION OF TREATMENT

- Selection of patients for MagicTouch - SCB treatment should be carefully done looking at the risks and benefits to the patient. Stent is generally avoided in patients with high risk of bleeding and MagicTouch - SCB may provide an alternative treatment strategy in such cases.
- Premorbid conditions that increase the risk of a poor initial result and the risks of emergency referral for bypass surgery (diabetes mellitus, renal failure, and severe obesity) should be reviewed. Multivariable patient conditions and outcomes modelling suggest that treatment assignment/strategy remained an independent predictor of clinical and angiographic outcomes even after adjusting for other baseline and procedural confounding variables.

9.0 PATIENT COUNSELLING INFORMATION

Physicians should consider the following in counselling patient about this product:

- Discuss the risks associated with MagicTouch - SCB balloon dilation treatment
- Discuss the risks associated with a Sirolimus-coated treatment
- Discuss the risks/benefits issues for this particular patient
- Discuss alteration to current lifestyle immediately following the procedure and over the long term.

10.0 CLINICAL DATA OF MAGICTOUCH - SCB

Extensive clinical evidence is available in more than 2,200 patients across the globe. Patients enrolled in all clinical studies belong to different countries demonstrating promising performance of the device in patients with different ethnicities across the globe.

Brief descriptions of all Clinical studies are described in the table below:

Sr. No.	Study Name	Indication (s)	Study Location (s)	N Patients, Lesions Enrolled	Results
1	**NANO-LUTE Study**	In-stent restenosis Small vessel disease Bifurcation lesions	India	450 patients, 479 lesions ISR: 211 De-novo SV: 197	**At 2-years:** **Overall population, N=408**[1] MACE: 4.2% TLR: 3.2% TV-MI: 0.2% Cardiac Death: 0.7% **Unpublished updated data*:** At 2-years, N=450 MACE: 5.6% TLR: 4.2% TV-MI: 0.4% Cardiac Death: 0.9% At 3-years, N=422 MACE: 7.1% TLR: 5.0% TV-MI: 0.5% Cardiac Death: 1.7% **De-novo small vessels** **Unpublished updated data*:** At 2-years, N=197 MACE: 3.6% TLR: 2.5% TV-MI: 0% Cardiac Death: 1.0% At 3-years, N=190 MACE: 4.7% TLR: 3.2% TV-MI: 0% Cardiac Death: 1.6% **In-stent restenosis** **Unpublished updated data*:** At 2-years, N=211 MACE: 8.5% TLR: 6.6% TV-MI: 0.9% Cardiac Death: 0.9% At 3-years, N=191 MACE: 11% TLR: 7.9% TV-MI: 1.0% Cardiac Death: 2.1%
2	**FASICO Registry**	In-stent restenosis Small vessel disease De-novo	Italy	32 patients, 34 lesions (12 patients with DES-ISR)	**Outcomes at 6.9 ± 1.7 months**[2] **No. of patients =32,** MACE: 9.4% TLR: 9.4% Cardiac Death: 0% MI: 0%
3	**FASICO NATIVES**	De-novo lesions	Italy	20 patients, 20 lesions	**Outcomes at 6.6 ± 2.5 months**[3] **No. of patients = 20** Lumen gain: 0.82 ± 0.38Late lumen loss: 0.09 ± 0.34
4	**EAST-BOURNE**	De-novo In-stent restenosis	Italy, Malaysia, India, Spain, Greece & Ireland	1471 Patients De-novo = 55% DES-ISR = 36% BMS-ISR = 9 %	**Rate of target vessel revascularization at 12 months**[4] **Overall population, N=605** MACE: 5.8% TLR: 2.5% MI:2.3% Death: 1% BARC: 0.3% **ISR, N=278** MACE: 8.6% TLR: 5.4% MI:1.8% Death: 1.4% BARC: 0%
5	**Real-World UK Retro-spective study**	De-novo In-stent restenosis	UK	288 patients 373 lesions ISR lesions: 140 De-novo lesions: 233	**Clinical outcomes at 6-months**[5] **Overall population, N=288** Death:3.4% Cardiac death:1.7% Target vessel-MI: 3.4% TLR per lesion: 12% MACE : 10% **De-novo population, N=186** Death:3% Cardiac death: 1.7% Target vessel-MI: 2% TLR per lesion:9% MACE:6% **ISR population, N=102** Death: 4.3% Cardiac death: 2.5% Target vessel-MI:6% TLR per lesion:17% MACE: 17.6%

6	**First-In-Man Assessment of the Novel MagicTouch - SCB Sirolimus Drug Coated Balloon Catheter for the Treatment of Coronary Lesions**	In-stent restenosis	Brazil	30 patients 26 Patients with QCA 23 Patients with IVUS	**Clinical outcomes up to 12 months*** **30 Days:** Death: 3.3% MI: 0% TLR: 0% MACE: 3.3% **6- Months** Death: 6.9% MI: 0% TLR: 3.4% MACE: 6.9% **12 Months** Death: 10.7% MI: 0% TLR: 7.1% MACE: 17.9% **QCA at 6-months, N =26:** Late Lumen Loss: 0.29 ± 0.29 **IVUS Results, N=23** NIH Volume (mm^3): 050 ± 045 % of NIH obstruction: 019 ± 011 Lumen Volume (mm^3): 218 ± 172

**Data on Concept Medical file.*

1. El-Mokdad R, di Palma G, Cortese B. Long-term follow-up after sirolimus-coated balloon use for coronary artery disease. Final results of the Nanolutè study. Catheterization and Cardiovascular Interventions. 2020 Mar 19.
2. Immediate and short-term performance of a novel sirolimus-coated balloon during complex percutaneous coronary interventions. The FAtebenefratelliSIrolimusCOated-balloon (FASICO) registry." Cardiovascular Revascularization Medicine 18.7 (2017): 487-491.
3. Cortese, Bernardo, et al. "Angiographic performance of a novel sirolimus-coated balloon in native coronary lesions: the FAtebenefratelliSIrolimusCOated NATIVES prospective registry." Journal of Cardiovascular Medicine 20.7 (2019): 471-476.
4. EASTBOURNE interim analysis-the largest trial on DCB ever, Bernardo Cortese, presented at Indialive 2020.

 Mid-term clinical outcomes from use of Sirolimus coated balloon in coronary intervention; data from real world population." Catheterization and Cardiovascular Interventions (2020).

11.0 HOW SUPPLIED

- STERILE: This device is sterilized with ethylene oxide gas and is non-pyrogenic. Do not use if the package is opened or damaged. Device is for one use only (single use). Do not re-sterilize.
- CONTENTS: One (1) MagicTouch - SCB (Rapid exchange delivery system).
- STORAGE: Temperature limitation: 8 to 25°C, store at dry place. Keep away from direct sunlight.

12.0 OPERATOR'S MANUAL FOR USE

MagicTouch - SCB is supplied and packaged in sterile condition. The product is packed in a box with identification labels on each packing. There are two stages of packing for maintaining the integrity and sterility of device - first pack is a Tyvek pouch which is sealed and stored in a secondary laminated aluminium foil pouch for protection. It is important that integrity of packing is verified before usage of device.

12.1 ACCESS TO STERILE PACK

Open the product box pack at tear line provided at the top. This provides access to the secondary aluminium foil pouch in which the sterile package is kept. Tear open the laminated aluminium foil pouch at the cut provided on the side. The inner most Tyvek pouch holds sterile MagicTouch - SCB. MagicTouch - SCB is stored in a coiled hoop. Open the Tyvek pouch and pass or drop the sterile balloon catheter into the sterile field using aseptic technique.

12.2 INSPECTION PRIOR TO USE

Before opening, carefully inspect the package, and check for damage to sterile barrier. Prior to using the device, carefully remove the system from package and inspect it for any bends, kinks or other damage. Do not use the device if any damage to the sterile package or balloon system is noted.

12.3 MATERIAL REQUIRED

Quantity	Material
N/A	Appropriate guiding catheter(s)
2-3	10-20 cc syringes
1,000 u /500 cc	Sterile Heparinized Normal Saline (HepNS)
1	0.014" (0.36 mm) diameter guidewire
1	Rotating haemostatic valve with an appropriate internal diameter (min. I.D. of 0.096" [2.4 mm])
N/A	Contrast diluted 1:1 with normal saline
1	Inflation device(With Luer Fitting)
1	Stopcock (3-way minimum)
1	Torque device
1	Guidewire Introducer
N/A	Appropriate anticoagulation and anti-platelet drugs

12.4 PREPARATION OF MAGICTOUCH - SCB

PRECAUTION

AVOID any manipulation of the drug coated balloon catheter during preparation as this may disrupt the integrity of coating on the device.

DO NOT wipe the distal part having balloon using any media as it may damage the coating.

ANY advancement of the drug coated balloon catheter should be in deflated condition and with use of guidewire in the coronary artery.

12.5 MAGICTOUCH - SCB SYSTEM PREPARATION

12.5.1 INITIAL PREPARATION OF MAGICTOUCH - SCB WITH INFLATION DEVICE

- Prepare the inflation device with diluted contrast medium.(60% Contrast Diluted 1:1 with HepNS in the Inflation device).
- Inflation Device Should be a 20 ml or 30 ml disposable device capable of producing a maximum pressure of 30-35 ATM/bars, fitted with a threaded plunger assembly with lock / release bar, a flexible high-pressure extension tube, and a three-way medium pressure stopcock.
- The inflation device should have a pressure gauge either analog or digital to measure the inflation pressure which would aid in visibility of the inflation pressure.
- For treatment of occluded vessels, contrast visualization of the distal vessel to confirm position of guidewire within the lumen is recommended.
- Attach the inflation device to the 3-way stopcock.
- Attach stopcock to the balloon inflation hub.
- Open the stopcock to balloon catheter system and leave inflation device system on neutral position.

12.5.2 ATTACHING THE INFLATION DEVICE TO MAGICTOUCH - SCB

- Create a fluid-fluid connection between prepared MagicTouch - SCB and the inflation device extension tube by connecting the luer connectors securely.
- Squeeze the trigger and pull back on the handle to apply a vacuum to the balloon.

12.5.3 MAGICTOUCH - SCB INFLATION AND DEFLATION USING INFLATION DEVICE:

- To inflate the balloon, squeeze or push (as the case may be) the trigger of the Inflation Device allowing the plunger to return to the resting position. Release grip on the trigger, locking the plunger into position. To increase pressure, rotate the handle clockwise until the desired inflation pressure is reached. The lock mechanism maintains the pressure.

 NOTE: Loss of pressure may indicate a leak in the system

12.6 MAGICTOUCH - SCB DELIVERY PROCEDURE

1. Prepare the vascular access site according to standard practice. Required access to the lesion must be done using an appropriate guide wire (recommended 0.014").
2. **Pre-dilate the lesion with a PTCA balloon catheter**. Limit the longitudinal length of pre-dilatation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the MagicTouch - SCB.
3. Maintain neutral pressure on the inflation device. Open the rotating haemostatic valve as widely as possible.
4. Backload the delivery system onto the proximal portion of the guidewire while maintaining the guidewire position across the target lesion.
5. Advance the drug coated balloon over the guidewire to the target lesion. Use the radiopaque balloon markers to position the balloon across the lesion; perform check angiography shoot to confirm the position of the balloon. Close the rotating

haemostatic valve at the right position of balloon catheter in lesion to lock it in place.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the drug coated balloon before inflation, the entire system should be removed as a single unit.

See Precautions - 5.9 Removal precautions for removing balloon catheter system.

12.7 MAGICTOUCH - SCB DEPLOYMENT PROCEDURE

1. Before deployment, reconfirm the correct position of the drug coated balloon relative to the target lesion via the radiopaque balloon markers.
2. Attach the inflation device (only partially filled with contrast media) to a three-way stopcock and apply negative pressure to purge the balloon of air.
3. Turn the stopcock on the catheter to the off position and purge the inflation device of air. Close the side port of the stopcock.
4. Under fluoroscopic visualization, inflate the balloon to at least the nominal pressure to deliver drug from the balloon surface, but do not exceed the labelled rated burst pressure. Optimal treatment and drug delivery require the balloon to be in full contact with the artery wall, with the internal diameter matching the size of the reference vessel diameter. Higher pressures may be required for certain lesions but Rated burst pressure should always not be exceeded. Refer to balloon compliance chart for further details on diameters and sizes.
5. Fully cover the entire lesion and dilated balloon area (including dissections) with the MagicTouch - SCB, allowing for adequate coverage into healthy tissue proximal and distal to the lesion.
6. If more than one MagicTouch - SCB is needed to cover the lesion area, adequately overlap balloon treatment, considering overlapping of the first treated segment. Ensure no gap between Drugs coated balloon treatments by visual estimation and including overlap to account for any geographical miss.
7. Keep MagicTouch - SCB inflated for a **minimum 60 seconds** (single inflation) at nominal pressure or **in case of multiple inflations the total inflation time should be summing up to 60 seconds** at nominal pressure without removing the balloon catheter.
8. Deflate the balloon by pulling a vacuum with the inflation device. Make certain that the balloon is fully deflated before attempting to move the catheter.
9. Confirm that vessel dilation is adequately done by angiographic injection through the guiding catheter.
10. Further dilation of treatment area should only be done if a scaffolding support of stent is planned at the site. A larger balloon may be used if the required diameter of the vessel is not achieved.

 NOTE: It is recommended that the guidewire and/or the balloon catheter remain across the lesion until the procedure is complete. Contrast media have different viscosities and may affect the inflation/deflation time.

12.8 MAGICTOUCH - SCB REMOVAL PROCEDURE

1. Ensure that the balloon is fully deflated.
2. While maintaining the guidewire position and negative pressure on the inflation device, withdraw the drug coated balloon delivery system.

 Note: Should unusual resistance be felt at any time during either lesion access or removal of the drug coated balloon catheter before inflation, the entire system should be removed as a single unit. See Precautions - 5.9 Removal precautions for removing balloon catheter system.
3. Repeat angiography to assess the treatment area.

12.9 MAGICTOUCH - SCB COMPLIANCE CHART

Refer to labelling.

13.0 STERILISATION AND STORAGE CONDITIONS

MagicTouch - SCB has been sterilized using Ethylene Oxide gas (EtO) and is non-pyrogenic. The device is for single use only. Storage temperature limitation: 8 to 25° C. Store at dry place. Keep away from direct sunlight.

14.0 PATIENT INFORMATION

In addition to these Instructions for Use booklet, the following patient specific information regarding the MagicTouch - SCB is available in product box:

A Patient Treatment Card includes specific information about MagicTouch - SCB (Product name, Size, Reference Number, Serial Number, LOT Number, & Expiry Date) and contact details for patient and his treating physician. It is always recommended that patients treated with MagicTouch - SCB keep this card in their possession and present it to their treating physician before going for any other treatment / surgical intervention.

15.0 DISCLAIMER OF WARRANTY

Concept Medical B.V. does not provide any expressed or implied warranty, including without limitation, any implied warranty of merchantability, fitness for a particular purpose or non-infringement, with respect to the product(s) described in this publication and Concept Medical B.V. hereby expressly disclaims any such warranties.

Descriptions or specifications in Concept Medical B.V. printed matter, including this publication, are meant solely to generally describe the relevant product at the time of manufacture and do not constituteany express warranties.

No person has the authority to bind Concept Medical B.V. to any representation or warranty.

Under no circumstances shall Concept Medical B.V. be liable for any medical or other expenses or direct, indirect, incidental, consequential, punitive or exemplary damages with respect to the products described herein, including without limitation, in respect of any use, defect, failure or malfunction, whether based upon warranty, contract, tort or other legal theory, even if Concept Medical B.V. has been advised or should have been aware of the possibility of such damages, other than as expressly provided by specific law without possibility of exclusion.

The exclusions, disclaimers and limitations contained in this section shall be interpreted in such manner as to be effective and valid under applicable law, but in case any one or more of the provisions contained herein shall, for any reason, be held to be invalid, illegal or unenforceable in any respect, such provision shall be ineffective to the extent, but only to the extent, of such invalidity, illegality or unenforceability without invalidating the remainder of this section or any other provisions hereof.

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1. El-Mokdad R, di Palma G, Cortese B. Long-term follow-up after sirolimus-coated balloon use for coronary artery disease. Final results of the Nanolute study. Catheterization and Cardiovascular Interventions. 2020 Mar 19.
2. Immediate and short-term performance of a novel sirolimus-coated balloon during complex percutaneous coronary interventions. The FAtebenefratelliSIrolimusCOated-balloon (FASICO) registry." Cardiovascular Revascularization Medicine 18.7 (2017): 487-491.
3. Cortese, Bernardo, et al. "Angiographic performance of a novel sirolimus-coated balloon in native coronary lesions: the FAtebenefratelliSIrolimusCOated NATIVES prospective registry." Journal of Cardiovascular Medicine 20.7 (2019): 471-476.
4. EASTBOURNE interim analysis-the largest trial on DCB ever, Bernardo Cortese, presented at Indialive 2020.Mid-term clinical outcomes from use of Sirolimus coated balloon in coronary intervention; data from real world population." Catheterization and Cardiovascular Interventions (2020).

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Hiç kimsenin Concept Medical B.V.'i herhangi bir temsil veya garantiye yasal olarak bağlayıcı olma yetkisi yoktur.