ANDOCOR

ARTERIAL CANNULA
CANULE ARTÉRIELLE
ARTERIELLE KANÜLE
ARTERIËLE CANULE
CÁNULA ARTERIAL
CANNULA PER ARTERIE
АРТЕРИАЛНА КАНЮЛА
ARTERIEKANYLE
ΑΡΤΗΡΙΑΚΗ ΚΑΝΟΥΛΑ
ARTÉRIÁS KANÜL
ARTERIJOS KANIULĖ
KANIULA TĘTNICZA
CÂNULA ARTERIAL
ARTÄRKANYL
ARTERIJSKA KANILA

INSTRUCTIONS FOR USE
MODE D'EMPLOI
GEBRAUCHSANWEISUNG
GEBRUIKSAANWIJZING
INSTRUCCIONES DE USO
ISTRUZIONI PER L'USO
ИНСТРУКЦИИ ЗА УПОТРЕБА
BRUGSANVISNING
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
HASZNÁLATI ÚTMUTATÓ
NAUDOJIMO INSTRUKCIJA
INSTRUKCJA UŻYCIA
INSTRUÇÕES DE UTILIZAÇÃO
BRUKSANVISNING
UPUTSTVO ZA UPOTREBU

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ARTERIAL CANNULA

Description
Andocor arterial cannulae are soft cannulae with encapsulated steel wire reinforcement or without wire reinforcement. The cannula consists of a curved, angled or beveled tip that is permanently attached to a PVC tapered body or wirewound PVC body.

Indications for Use
These cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass. This device is intended to be used by a cardiovascular surgeon and is intended to be used on patients undergoing cardiopulmonary bypass surgery.

Contraindications
This device is not designed for use except as indicated.

Warnings
Read all warnings, precautions and instructions carefully prior to use. Store in a dry place, below a temperature of 40 °C. Avoid bending or kinking the cannula. This cannula is for short-term, temporary use only. Check that connectors, luer caps or vented caps are secured before use.
Do not clamp on reinforced tubing areas.

Precautions
Proper surgical procedures and techniques are the responsibility of the medical professional. Exercise care during the connection of the cannula and the perfusion line to insure that no air is trapped in the perfusion line. Placement location and orientation of the sutures in the aorta is critical to proper positioning of the cannula. Once the cannula is firmly secured with the sutures, it is difficult to change the position of the cannula.
Do not exceed recommended flow.

Instructions for use
Place a purse string suture in the anterior aspect of the distal ascending aorta at the desired location. Place a tourniquet over this purse string. The tourniquet on the inner purse string is tightened and secured to the cannula by tying a securing ligature around the cannula and tourniquet at a convenient site. Position the cannula and arterial line so that the cannula has a gentle curve and enters the aorta perpendicularly. Arterial cannulae allow oxygenated blood to be returned to the patient's arterial system. Certain arterial cannulae are manufactured with suture rings. These rings are provided as a fixation point for suturing and act to retain the tip of the cannula within the vessel.
Contains di(2-ethylhexyl) phthalate(DEHP).

Disposable
Single use only.
Do not reuse or resterilize this device. Reuse or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device. Sterile and non-pyrogenic.

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AORTIC ROOT CANNULA

Description
The Aortic Root Cannula and Cardioplegia consists of a soft, flexible, thin wall tip permanently attached to a flexible tube. A female luer fitting is permanently attached to the distal end of the tube. The introducer is packaged within the cannula body. The Aortic Root Cannula with Vent Line consists of a soft, flexible, thin wall tip permanently attached to a flexible Y-type fitting.

Indications for use
This cannula is intended for use during cardiopulmonary bypass for delivering cardioplegia solutions and venting of the heart. The cannula may also be used to aspirate air from the aorta at the conclusion of the cardiac procedure.
This device is intended to be used by a cardiovascular surgeon and is intended to be used on patients undergoing cardiopulmonary bypass surgery.

Contraindications
This device is not intended for use except as indicated above.

Warnings
Read all warnings, precautions and instructions carefully prior to use.
Store in a dry place, below a temperature of 40°C.
If the introducer needle is removed from the cannula BEFORE placement into the aorta, use caution when reinserting the introducer to prevent cutting or gouging the inner surface of the cannula or tip.
DO NOT re-insert the introducer needle into the cannula once removed AFTER placement into aorta.

Instructions for use
Place two purse string sutures using pledgets at the desired cannulation site.
Draw the sutures through tourniquets.
Close the clamp on the vent line.
Insert the cannula into the aorta in the center of the area enclosed by the purse string sutures.
Caution: Cannula must be fully inserted with side holes inside the aorta to avoid risk of false lumen within the aortic wall. During infusion of the cardioplegia solution, care should be taken to avoid manipulations of the heart and aorta that could displace the cannula.
Contains di(2-ethylhexyl) phthalate (DEHP).

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MULTI-STAGE VENOUS CANNULA

Indications for Use
This cannula is intended for use in venous drainage via the right atrium and inferior vena cava simultaneously during cardiopulmonary bypass surgery up to six hours or less.
This device is intended to be used by a cardiovascular surgeon and is intended to be used on patients undergoing cardiopulmonary bypass surgery.

Contraindications
This device is not intended for use except as indicated above. This cannula is contraindicated for use in procedures where the chambers of the heart will be opened or in cases with septal defects.

Warnings
Read all warnings, precautions and instructions carefully prior to use.
Store in a dry place, below a temperature of 40°C.
The introduction of air into the right heart may result in the immediate total loss of venous drainage to the extracorporeal circulation. Position the cannula and venous line tubing carefully to avoid excessive tension on the right atrium. Do not bend the cannula too acutely or excessive tension on the right atrium may result. It is suggested that the cannula and venous line tubing exit from the chest toward the head when the cannula tip is placed in the inferior vena cava.

Instructions for use
Select a cannula size which will provide adequate venous return rates for a given size patient.
Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the cannula onto the sterile field utilizing an aseptic technique.
Verify that the occluder is firmly engaged in the distal tip of the cannula.
Contains di(2-ethylhexyl)phthalate (DEHP).

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SINGLE STAGE VENOUS RETURN CANNULA

Description
Venous return cannulae from Andocor consist of soft PVC tubing with or without encapsulated steel wire reinforcement, terminating in a tapered tip on the distal end. They are offered in straight and right angled models. Multiple tip variations are available: bullet tip, lighthouse tip and 90° curved metal tip.

Indications for use
Venous return cannulae are intended for cannula drainage from the superior and inferior vena cava during cardiopulmonary bypass surgery up to six hours or less.
This device is intended to be used by a cardiovascular surgeon and is intended to be used on patients undergoing cardiopulmonary bypass surgery.

Contraindications for use
This device is not designed for use except as indicated above.
Single tube venous drainage from the right atrium should not be used in procedures where air from the procedure site may enter the cannula.

Warnings and precautions
Read all warnings, precautions and instructions carefully prior to use.
Store in a dry place, below a temperature of 40°C.
This device is designed, intended, and distributed for single use only.
Use the largest cannula possible, consistent with good surgical technique.
Wire reinforced cannulae should be clamped in the non-reinforced section located at the connector end since clamping of the reinforced section may produce permanent cannula deformation.
Proper surgical procedures and techniques are the responsibility of the medical profession.

Instructions for use
Remove the device from the packaging, using the aseptic technique. If the package is damaged, do not use the device as sterility may have been compromised.
Using the appropriate surgical technique, insert the cannula into the cannulation site and secure in place.
Back flush the device to remove air from the cannula tubing.
Carefully connect the proximal end of the device to the venous drainage line of the extracorporeal circuit, hereby preventing entry of air into the circuit, and initiate cardiopulmonary bypass.
Position the cannula so that it has a gentle curve, in order to prevent kinking of the cannula tubing.
After completion of the bypass procedure, remove the cannula and close the cannulation site using the appropriate surgical technique.

Disposable
Single use only.
Do not reuse or resterilize this device. Reuse or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device.
Sterile and non-pyrogenic.

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LEFT HEART VENT CATHETER

Product Description
The Left Heart Vent Catheter consists of a flexible tube with a perforated distal segment.
A small tapered tip is attached to the perforated distal end to aid in the insertion of the catheter across the mitral valve. A malleable stainless steel introducer is provided with the catheter.
Depth markings on the tube indicate insertion depth.

Indications for use
This catheter is intended for use in venting the left heart during cardiopulmonary bypass surgery up to six hours or less.
This device is intended to be used by a cardiovascular surgeon and is intended to be used on patients undergoing cardiopulmonary bypass surgery.

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Warnings
Read all warnings, precautions and instructions carefully prior to use.
Store in a dry place, below a temperature of 40 °C.

Instructions for use
Inspect the package and product for damage and expiration date. If undamaged and unexpired open the package and transfer the catheter onto the sterile field utilizing an aseptic technique.
Insert the catheter at least past the last set of holes on the distal tip of the catheter or at least to the first line on the catheter. During insertion, palpate the left ventricle to insure proper positioning of the catheter across the mitral valve.
After the procedure remove the vent catheter from the heart and tie the purse string suture.

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PERICARDIAL SUMP

Description
The Pericardial Sump consists of a soft PVC tubing and features a weighted metal tip. Certain models feature a smooth metal tip, while others are constructed with a wirewound tip which ends in a rounded metal tip. The devices terminates in a 1/4" connector on the proximal end.

Indications for use
The Pericardial Sump is intended for the removal of excess fluid from the surgical field.
This device is intended to be used by a cardiovascular surgeon and is intended to be used on patients undergoing cardiopulmonary bypass surgery.

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Warnings
Read all warnings, precautions and instructions carefully prior to use.
Store in a dry place, below a temperature of 40°C.
Proper surgical procedures and techniques are the responsibility of the medical professional. Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience and the type of procedure employed.
The model with a wirewound tip should not be placed through the leaflets of a valve as damage to the valve may occur.
This device contains stainless steel, do not use in an MRI environment.

Instructions for use
Remove the device from the packaging, using the aseptic technique. If the package is damaged, do not use the device as sterility may have been compromised. Flush the device. Transfer the device to the surgical field, using the appropriate surgical technique. Connect the proximal end of the device to the suction line.
Contains di(2-ethylhexyl) phthalate (DEHP).

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RETROGRADE CARDIOPLEGIA CANNULA

Description
The Retrograde Cardioplegia Cannula comprises a cannula body with a bullet tip, terminating in a female luer lock on the proximal end. A pressure monitoring line is a part of the Retrograde Cardioplegia Cannula body, which starts at the bullet tip and ends at a 3-way stopcock. A solid guide wire stylet with handle is equipped with a preformed curve for insertion.

Indications for use
The Retrograde Cardioplegia Cannula is intended to be connected to the cardioplegia line for delivery of cardioplegia solutions and blood into the patient's coronary sinus during cardiopulmonary bypass procedures.
This device is intended to be used by a cardiovascular surgeon and is intended to be used on patients undergoing cardiopulmonary bypass surgery.

Contraindications
This device is not designed for use except as indicated above.

Warnings
Read all warnings, precautions and instructions carefully prior to use.
Store in a dry place, below a temperature of 40°C.
This equipment is only intended for short term use (< 6 hours). Appropriate use of this device must be considered by the user, based on medical training, experience, the type of procedure employed and the benefits and risks related to the use of the device for each patient. Insert the catheter before initiating full bypass, to avoid air entering into right atrium or coronary sinus. Use caution during introduction of the catheter and manipulation of the heart while the catheter is in place, in order to minimize trauma to the coronary sinus.

Instructions for use
Remove the device from the packaging, using the aseptic technique. If the package is damaged, do not use the device as sterility may have been compromised. Do not bend the stylet into an acute angle. Place a 1 cm purse string in the mid to the upper right atrium and make a small atriotomy through the purse string. Reposition the cannula if resistance occurs. Do not force the cannula into the coronary sinus in any case. To insert, control the cannula with the left hand while the right index finger is placed at the junction of the IVC and AV groove inferiorly. The tip of the cannula should be placed against the finger and then guided and positioned in the coronary sinus.
The stopcock of the pressure monitoring line must be attached to a pressure monitoring system, make sure that all air has aspirated from the pressure monitoring line. Unclamp the cannula body and connect it to the primed cardioplegia delivery line, make sure that no air is present within the tubing or cannula. Now open the pressure monitoring stopcock.
After crossclamping the aorta, gradually start the retrograde cardioplegia flow. Observe pressure in the coronary sinus. Red coloration of the coronary veins on the anterior surface of the heart should be observed and, likewise, blue coloration of the coronary arteries if blood cardioplegia is used.
Contains di(2-ethylhexyl) phthalate (DEHP).

Disposable
Single use only. Do not reuse or resterilize this device. Reuse or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device.
Sterile and non-pyrogenic.

For further information and/or in case of a complaint, malfunction or any serious incident that occurs in relation to the use of this device, please contact ANDOCOR using the contact details provided in these instructions for use.

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GAS DIFFUSER

Description
The Gas Diffuser consists of a PVC tubing which can be connected to a CO2 or other gas source, with a bacterial filter and a well-defined porous plug which will diffuse the gas as widely and precise as possible.

Indications for use
The gas diffuser is intended for intraoperative wound ventilation with CO2 in order to avoid air embolisms and bacteriological contamination during open heart sugery.
This device is intended to be used by a cardiovascular surgeon and is intended to be used on patients undergoing cardiopulmonary bypass surgery.

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Warnings
Read all warnings, precautions and instructions carefully prior to use.
Store in a dry place, below a temperature of 40°C.
Contains di(2-ethylhexyl) phthalate (DEHP).

Instructions for use
Preparations
1. Ensure that a CO2 source is available with the capacity of 10L/minute via a flow meter for the whole of the estimated period of the operation.
2. The flow meter must be connected to or embody a pressure regulator. The working pressure for the CO2-gas must be between 3,5 and 4.5 kPa.
3. Place the body part into the lower part of the operation wound. Bend the narrow hose so that the body part is located about 5cm below the edge of the wound.
4. Secure the hose.
5. Connect the open part of the wide hose to the flow meter.

Starting and supplying the flow
1. Start a flow of CO2 at 10L/minute for at least 1 minute before the heart and the major blood vessels are to be opened.
2. Maintain a supply of 10L until the heart and the major blood vessels have been closed. Then switch off the flow.
3. If the body becomes completely soaked in fluid, briefly squeezing the body will restore its function.
4. Avoid using a suction rate that exceeds the gas flow rate of CO2. Otherwise there is a risk that room air may enter the operation cavity.

After use
1. Dispose of the gas diffuser in the same way as blood contaminated materials.
2. The gas diffuser is intended for single use only and must not be cleaned and reused.

Disposable
Single use only. Do not reuse or resterilize this device. Reuse or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device. Sterile and non-pyrogenic.

Warranty and limitation of liability
1. The manufacturer warrants that the gas diffuser has been manufactured in accordance with its stated specifications and in compliance with good manufacturing practice, along with regula tory requirements.
2. The manufacturer will replace at no cost a gas diffuser that has a manufacturing defect, on condition that it is detected and reported within 6 months after delivery and that the product lot number can be stated.
3. No liability can be accepted for damage to personnel or property caused by the wrongful application or storage of the product. Liability claims will not be accepted if the gas diffuser is applied in spite of damaged packaging or the sterility has been expired.

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