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EN Operating Instructions
Lunos® Prophy Powder Gentle Clean

The following variants of „Prophy Powder Gentle Clean“ are available:

Spearmint (REF: CP2630A2250), contains limonene
Orange (REF: CP2630A2250), contains limonene
Berry (REF: CP2650A2250), contains limonene
Neutral (REF: CP2610A2250), contains no flavouring agents

Product description
“Prophy Powder Gentle Clean“ is a medical device. This product is only intended for professional dental application and must only be used by trained specialists.
“Prophy Powder Gentle Clean“ is designed for supragingival application using commercial powder jet devices. Please note the operating instructions of your powder jet device before use of the prophy powder.
The formulation is free of animal ingredients, gluten, lactose, fructose and microplastics. This product does not have a cariogenic effect.

Indication
This product is recommended for:
- Supragingival cleansing of the tooth surface
 - Removal of discolouration and plaque
- Preparation and cleaning of the tooth surface before:
 - Adhesive attachment of orthodontic brackets
 - Application of fluoride and bleaching preparations
 - Insertion of composite restorations, fixed dentures, e.g. inlays, onlays, crowns, veneers etc.
 - Application of sealing materials

This prophy powder cannot be used to remove tartar.
This prophy powder cannot be used to prepare enamel.
“Prophy Powder Gentle Clean“ is suitable for children and adults. Do not swallow.

Contraindications
Patients suffering from chronic airway diseases (e.g. bronchitis or asthma) must not be treated with a powder-water jet device. Inhalation of the aerosol powder jet can lead to difficult breathing. Likewise, the product should not be used on non-intact enamel, hypermineralised areas or on patients with dry mouth.
In certain cases, flavouring agents present in the powder may cause skin irritation or allergic reactions. For this reason, patients with known allergies to flavouring agents should be treated with “Prophy Powder Gentle Clean Neutral“.

Warning
Do not swallow. Avoid contact with the eyes. In case of contact with the eyes, rinse with plenty of water. All serious events that occur in connection with the product must be reported immediately to the manufacturer and the competent authority of your member state. The product is intended for use on humans.

Precautions
Use appropriate protective equipment when working with “Prophy Powder Gentle Clean“.
The patient and operator need to wear suitable eye protection, e.g. safety glasses.
The operator needs to wear suitable protective equipment in the form of mouth and nose protection and suitable protective gloves.
Use a suitable spray mist aspiration cannula, e.g. the Dürr Dental Prophy Cannula.
“Prophy Powder Gentle Clean“ may only be used by appropriately trained dental specialists.
Make sure that the powder jet device is not positioned directly at fixed prostheses or fillings, as the restorations may be damaged.
Do not direct the spray nozzle at the gingiva during the treatment. Direct exposure of the soft tissue may damage the gingiva. The powder jet nozzle must not be aimed directly into the gingival pocket, as there is a risk of emphysema arising in the tissue and/or of embolism in the case of air entering the blood circulation.
To protect sensitive skin, please apply Vaseline to the patient's lips before commencing the treatment.
After completion of treatment, clean the patient's face with a moist disposable cloth, e.g. the Lunos® Prophy Wellness Towel.
Please advise your patient that the ingestion of food or beverages (tea, coffee or other food items) within the first two to three hours after completion of the treatment can cause renewed discolouration of the teeth.
Clean the surfaces of the dental unit with a moist cloth following the use of “Prophy Powder Gentle Clean“.
Please make sure to use “Prophy Powder Gentle Clean“ in powder jet devices for supragingival application only.

Preparation
Please remove all residual powder from the powder chamber of the powder jet device before refilling it with “Prophy Powder Gentle Clean“.
Shake the closed bottle before removing any “Prophy Powder Gentle Clean“. Remove the sealing foil and ensure that it is undamaged. For this purpose, completely remove the lid from the bottle, pull off the seal and then fully screw the lid back onto the bottle.
It is essential that you make sure that the powder chamber is completely dry before filling it. There must be no residual moisture in the chamber, e.g. from reprocessing measures.
When filling the powder chamber with powder, you need to make sure not to exceed the maximum filling level of the powder chamber.
To minimise the amount of dust generated when you fill the chamber, it is recommended that you use the lever to compress the bottle slightly.

Storage and shelf life
Please make sure that the bottles are carefully closed again immediately after removing the powder so that the prophy powder does not absorb any moisture.
Store in enclosed, well-ventilated rooms away from moisture and bright light between 0 °C and 25 °C at 10 % to 75 % relative humidity. Keep out of direct sunlight. When the storage conditions are observed, the shelf life is 36 months in the unbroken container and 12 months after opening.

Disposal
Dispose of in accordance with the local official regulations.

Composition
Prophy Powder Gentle Clean Neutral: Trehalose, excipients
Prophy Powder Gentle Clean Orange: Trehalose, excipients, orange flavouring (contains limonene)
Prophy Powder Gentle Clean Spearmint: Trehalose, excipients, spearmint flavouring (contains limonene)
Prophy Powder Gentle Clean Wild Berry: Trehalose, excipients, berry flavour (contains limonene)

UK-REP: DÜRR DENTAL (Products) UK Ltd.
14 Linnell Way – Telford Way Industrial Estate
Kettering, Northants NN16 8PS, United Kingdom ☎ +44 1536526740

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Prophy Powder Perio Combi (REF: CPZ640A1950)

Product description

"Prophylaxis Powder Perio Combi" is a medical device. This product is only intended for professional dental applications and must be used exclusively by qualified personnel! "Prophylaxis Powder Perio Combi" is designed for sub- and supragingival application in commercial prophy devices. Please note the operating instructions of the prophy powder jet device before use of the prophylaxis powder. The formulation is free of animal ingredients, gluten, lactose, fructose and microplastics.

Indication

This product is recommended for:

- Removal of periodontal biofilm in the scope of periodontal maintenance therapy (recall therapy), with a suitable subgingival nozzle, e.g. the Perio Tip from Dürr Dental
- Cleaning implant surfaces to support treatment of mucositis or peri-implantitis after completed initial therapy
- Gentle removal of plaque at the sulcus and particularly sensitive dental surfaces
- Removal of soft deposits/biofilm in supragingival regions
- Removal of plaque on orthodontic appliances (brackets)

Prophylaxis Powder Perio Combi is suitable for children and adults. This product does not have a cariogenic effect. Do not swallow.

Contraindication

Patients suffering from chronic airway diseases (e.g. bronchitis or asthma) must not be treated with a powder-water jet device since the powder present in the aerosol can lead to difficult breathing. Likewise, the product should not be used on non-intact enamel, dentin defects, cracked teeth or dental restorations. The treatment of deep periodontal pockets can lead to bacteraemia. Therefore, suitable measures and restrictions must be undertaken for the treatment of risk patients (endocarditis, pregnant and lactating women, diabetes, immunodeficiency (neutropenia, agranulocytosis, diabetes, hemophilia), patients on treatment (radiation therapy, chemotherapy, antibiotic therapy)).

Warning

Do not swallow. Avoid contact with the eyes. In case of contact with the eyes, rinse with plenty of water. All serious events that occur in association with the product must be immediately reported to the manufacturer and the competent authority of your member state. The product is intended for use on human subjects.

Precautions

Use appropriate protective equipment when working with "Prophylaxis Powder Perio Combi". The patient and operator need to wear suitable eye protection, e.g. protective goggles. The operator must be equipped with respiratory protection in the form of mouth and nose protection and suitable protective gloves. Use a suitable spray mist aspiration cannula, e.g. the Dürr Dental Prophylaxis Cannula. "Prophylaxis Powder Perio Combi" may only be used by appropriately trained dental specialists. Please make sure not to direct the jet of the powder jet device directly at fixed dental restorations and fillings as this may damage the restorations. Do not direct the spray nozzle directly at the gingiva during the treatment. Direct exposure of the soft tissue may damage the gingiva. Also there is a risk of emphysema and the risk of embolism from air being introduced into blood circulation. To protect sensitive skin, please apply Vaseline to the patient's lips before commencing the treatment. After completion of treatment, clean the patient's face with a moist disposable cloth, e.g. the Lunos® Prophylaxis Wellness Cloth. Please advise your patient that the ingestion of food or beverages (tea, coffee or other stains) within the first two hours after completion of the treatment can cause renewed discoloration of the teeth. Clean the surfaces of the dental unit with a moist cloth following the use of the prophylaxis powder.

Preparation

Please remove all residual powder from the powder chamber of the powder jet device before refilling it with "Prophylaxis Powder Perio Combi". Shake the closed bottle before removing any "Prophylaxis Powder Perio Combi". Open the bumper-proof seal before removing the powder. For this purpose, take the lid off the bottle, pull off the seal and replace the lid. Fill the powder chamber of the powder jet device with the bottle. It is essential to ensure that the powder chamber is completely dry before filling it. There must be no residual moisture in the chamber, e.g. from reprocessing. Immediately after filling the powder chamber with powder, you need to make sure not to exceed the maximum filling level of the powder chamber.

To minimize the amount of dust generated when you fill the chamber, it is recommended that you use the rocker and compress the bottle slightly.

Storage and shelf life

Please make sure that the bottles are carefully closed again immediately after removing the powder so that the prophylaxis powder does not absorb any moisture. Store unopened, well-ventilated rooms with a maximum ambient temperature of 30°C and 25°C at 10% to 75% relative humidity. Keep out of direct sunlight. When the storage conditions are observed, the shelf life is 36 months in the sealed container and 12 months after opening.

Disposal

Dispose of in accordance with the local official regulations.

Composition

Prophy Powder Perio Combi: Trehalose, excipients

UK-REP: DÜRR DENTAL (Products) UK Ltd.
Unit 2 Network Centre, Industrial Estate
Kettering, Northants NN16 8PS, United Kingdom ☎ +44 1536526740

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