Amplatzer™ Amulet™ Left Atrial Appendage Occluder

Instructions for Use

Abbott

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Amplatzer™ Amulet™ Left Atrial Appendage Occluder

EN: Instructions for Use

CAUTION: Read all instructions carefully. Failure to follow these instructions, warnings, and precautions may lead to patient injury or device damage.

Device Description

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is a transcatheter, self-expanding device intended for use in preventing thrombus embolization from the left atrial appendage (LAA). The device is constructed from a nitinol mesh and consists of a lobe and a disc connected by a central waist. Polyester patches are sewn into both the lobe and disc to facilitate occlusion. The lobe has stabilizing wires to improve device placement and retention. The device has threaded screw attachments at each end for connection to the delivery and loading cables. The device has radiopaque markers at each end and at the stabilizing wires that permit visibility during fluoroscopy.

Refer to the figures and tables in Appendix A: Supplemental Information for more information about the device. Table 1 identifies device and delivery sheath dimensions; Table 2 lists sizing information. The table below identifies what Figure 1 and Figure 2 illustrate.

Figure 1. Occluder Components		**Figure 2. Delivery Components**	
1	Proximal end screw	9	Loader with device
2	Waist of device	10	Loader hub
3	Lobe	11	14 Fr flush adaptor (for sizes 28 mm–34 mm)
4	Distal end screw	12	Hemostasis valve
5	Marker bands	13	Delivery cable
6	Stabilizing wires	14	Delivery cable vise
7	Platinum thread	15	14 Fr sheath adaptor (for sizes 16 mm–25 mm)
8	Disc	16	Loading cable vise
		17	Loading cable

Intended Purpose

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is a percutaneous, transcatheter device intended to occlude the left atrial appendage.

Indication for Use

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is indicated to prevent thrombus embolization from the left atrial appendage (LAA) in patients who have nonvalvular atrial fibrillation.

Performance Characteristics

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder consists of the following:

- LAA occluder
- Delivery cable
- Loader
- Loading cable
- Loading cable vise
- Hemostasis valve
- Delivery cable vise
- 14 French sheath adaptor (16 mm–25 mm sizes only)
- 14 French flush adaptor (28 mm–34 mm sizes only)

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is intended to be delivered with the Amplatzer™ Amulet™ Delivery Sheath or the Amplatzer™ Steerable Delivery Sheath.

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is a self-expanding device made primarily from a nitinol wire mesh. The occluder construction includes a lobe and disc connected by a central waist. Eight unique size configurations are available to facilitate the occlusion of a variety of left atrial appendage (LAA) anatomies. The occluder is pre-loaded in a loader and pre-attached to a delivery cable (both loader and delivery cable are included in the package). A loading cable and loading cable vise are included, should the implant require reloading.

The loader reduces the size of the LAA occluder, allowing introduction of the occluder into the delivery sheath. The 14 French flush adaptor is attached to the 14 French loader on devices with lobe sizes 28 mm, 31 mm, and 34 mm to allow connection of the loader to the hemostasis valve. The hemostasis valve, pre-attached to a three-way stopcock, allows flushing of the occluder within the loader and controls blood backflow during the procedure. The 14 French sheath adaptor is included in devices with lobe sizes 16 mm, 18 mm, 20 mm, 22 mm, and 25 mm to facilitate the optional connection of the 13 French loader to a 14 French Amplatzer™ Amulet™ Delivery Sheath or Amplatzer™ Steerable Delivery Sheath.

Limitations

The implantation of the device can be limited by the presence of a previously implanted closure device within the interatrial septum, such as an atrial septal defect occluder or a patent foramen ovale occluder.

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Device Lifetime

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is expected to remain implanted permanently. The device lifetime for the Amplatzer™ Amulet™ Left Atrial Appendage Occluder is at least 10 years, as demonstrated by bench testing of the occluder.

Intended Clinical Benefit

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is a percutaneous, transcatheter device intended to close the LAA in patients with nonvalvular atrial fibrillation. The following clinical benefits are expected to be achieved:

- A reduced risk of an ischemic stroke due to thrombus embolization (an ischemic stroke is one caused by a blocked artery) compared to patients not taking oral anticoagulants (OAC)[1]
- A reduced risk of major bleeding events compared to patients taking OAC[2, 3]
- A minimally invasive transcatheter procedure that results in a high implant-success rate and a low rate of procedure-related major adverse events[1, 4]
- A high rate of immediate, complete closure of the LAA, which makes it possible for patients to be discharged without OAC[1, 4]

Contraindications

The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is contraindicated for patients:

- with the presence of intracardiac thrombus
- with active endocarditis or other infections producing bacteremia
- where placement of the device would interfere with any intracardiac or intravascular structures

Warnings

- Physicians must be prepared to deal with urgent situations such as pericardial effusion or device embolization that can require removal of the device.
- Remove embolized devices. Do not remove an embolized device unless the device is fully recaptured inside a delivery sheath.
- This device should be used only by physicians who are trained in standard transcatheter techniques. The physician should determine which patients are candidates for procedures that use this device.
- The use of post-procedure anticoagulation therapy may be associated with an increased potential for a late pericardial effusion. Physicians should monitor for signs and symptoms of pericardial effusion and obtain appropriate imaging when indicated. Physicians should also consider routine echocardiography to screen for pericardial effusion.
- The Amplatzer™ Amulet™ LAA Occluder consists of a nickel-titanium alloy, which is generally considered safe. However, in vitro testing has demonstrated that nickel is released from this device for a minimum of 60 days. Patients who are allergic to nickel may have an allergic reaction to this device, especially those with a history of metal allergies. Certain allergic reactions can be serious; patients should be instructed to notify their physicians immediately if they suspect they are experiencing an allergic reaction such as difficulty breathing or inflammation of the face or throat. Some patients may also develop an allergy to nickel if this device is implanted.
- Do not use this device if the sterile package is open or damaged.
- The device was sterilized with ethylene oxide and is for single use only. Do not re-use or resterilize this device. Attempts to re-use or resterilize this device can cause a malfunction, insufficient sterilization, or harm to the patient.
- DO NOT use the Amplatzer™ Amulet™ Left Atrial Appendage Occluder after the Use-by date stated on the package label.
- If the device is retracted while it is in the sheath, the device and the sheath must both be removed and replaced. Failure to replace both the device and the sheath may result in sheath and/or device malfunction.
- If the device is retracted farther than the radiopaque markers (fully recaptured), the device and the delivery sheath must both be removed and replaced. Failure to replace both the device and the sheath may result in sheath or device malfunction.

Precautions

- The physician should exercise clinical judgment in situations that involve the use of antithrombotic drugs before, during, and/or after the use of this device.
- The physician should exercise caution if implanting a device in a patient who has an implantable cardioverter defibrillator (ICD) or pacemaker and leads.

1. Hildick-Smith D, Landmesser U, Camm AJ, et al. “Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study.” European Heart Journal 41, no. 30 (2020): 2894-2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, et al. “Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation.” JACC Cardiovascular Interventions 14, no. 1 (2021): 69-78.
3. Chen S, Julian Chun K.R., Ling Z, et al. “Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update.” Journal of Cardiovascular Development and Disease 8, no. 6 (2021): doi: 10.3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Epub ahead of print. PMID:34459659.

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- The physician should have the guidewire in the left upper pulmonary vein when making exchanges in the left atrium.
- Ensure that the vasculature is adequate for the sheath size being selected.
- The physician should exercise caution if performing ablation at or near the implant site after the device is implanted.
- Use standard interventional cardiovascular catheterization techniques when using Amplatzer™ products.

Precautions for Special Populations

- Pregnancy: The safety and effectiveness of this occluder has not been established during pregnancy. Fluoroscopic X-ray guidance is used during placement of the device. The risk of increased X-ray exposure for patients who are pregnant must be weighed against the potential benefits of this technique.
- Nursing mother: The safety and effectiveness of this occluder has not been established in lactating mothers. There has been no quantitative assessment for the presence of leachables in breast milk.
- Pediatric Population: The safety and effectiveness of this occluder has not been established in a pediatric population.

Storage and Handling Information

There are no special storage requirements beyond what is typical for medical devices.

This Instructions for Use is recyclable.

Summary of Safety and Clinical Performance

A summary of the safety and clinical performance for this device is available at https://ec.europa.eu/tools/eudamed: This is the SSCP location after the launch of the European Database on Medical Devices/Eudamed. Search for the device using the basic UDI-DI: 5415067AMP1000DQ.

Patient Identification Card

A patient identification card is included with product documentation. Hospital personnel should complete the front side of the Patient Identification Card by printing the patient's name, implant date, along with the physician's name and healthcare facility information and phone number in the designated areas. A peel-away device label is provided in the sterile packaging, which includes device information (model number, lot number, UDI number, and device name). Apply the peel-away device label to the back of the card in the designated area. If a peel-away label is not available, write the information in the designated area on the back of the card. In procedures where multiple implantable devices are implanted, include information for all implants in the designated area.

NOTE: Write the device lot number in the SN/LOT field.

After the card is completed, separate the top portion from the Patient Identification Card by tearing along the perforated line. Provide the completed Patient Identification Card to the patient. Instruct the patient to carry the Patient Identification Card at all times and to present the implant card to hospital staff in emergency situations. To obtain a replacement card if a patient loses or damages their card, contact Abbott Technical Support. Once the replacement card is received, handwrite the required information in the respective areas and provide to the patient.

Product Materials

The LAA occluder is composed of the following materials, which are intended to come into contact with tissue:[5]

Table 1. Product Composition

Component	Material Description	Substance	Substance Concentration within each Component (% w/w)
Distal and proximal patch	Polyester sheet	Polyethylene terephthalate	90–100%
Yarn	Polyester Filament Yarn and Staple Fiber	Polyester Blend	100%
Device Body and Stabilizing Wires	Nitinol Wire	Nickel	<57%
		Titanium	<43%
		Cobalt	<0.1%
		Iron	<0.1%
		Carbon	<0.1%
		Oxygen	<0.1%
		Niobium	<0.1%
		Copper	<0.1%
		Chromium	<0.1%
		Hydrogen	<0.1%
		Nitrogen	<0.1%
Marker Bands	80% Platinum / 20% Iridium alloy	Platinum	80%
		Iridium	20%
Radiopaque Cable	90% Platinum / 10% Iridium alloy	Platinum	90%
		Iridium	10%

5. All information included in this table was either self-reported in vendor declarations or reported in vendor certification of conformance documents.

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Table 1. Product Composition (Continued)

Component	Material Description	Substance	Substance Concentration within each Component (% w/w)
End Screw	316L Stainless Steel	Iron	<69%
		Chromium	<17%
		Nickel	<11%
		Molybdenum	<3%
		Manganese	<2%
		Silicon	<1%
		Nitrogen	<0.1%
		Phosphorus	<0.1%
		Carbon	<0.1%
		Sulfur	<0.1%
		Cobalt	<0.4%

CAUTION: One or more components of this device may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt: Chemical Abstracts Service (CAS) No. 7440-48-4; EC No. 231-158-0.

Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

MR Conditional

Non-clinical testing has demonstrated that the Amplatzer™ Amulet™ Left Atrial Appendage Occluder is MR Conditional. A patient with the Amplatzer™ Amulet™ Left Atrial Appendage Occluder can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0 Tesla (3.0T)
- Maximum spatial gradient field of 19 T/m (1900 G/cm)
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than or equal to 4ºC after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 20 mm from the Amplatzer™ Amulet™ Left Atrial Appendage Occluder when imaged with a gradient echo pulse sequence in a 3.0T MR system.

Potential Adverse Events

Potential adverse events that can occur during or after this procedure include, but are not limited to:

Air embolism
Airway trauma
Allergic reaction/toxic effect due to contrast media, medication, or metal
Anemia
Anesthesia reaction (nausea, vasovagal reaction, confusion/altered mental status/other)
Arrhythmia
Atrial septal defect
Bleeding
Cardiac arrest
Cardiac tamponade
Chest pain or discomfort
Congestive heart failure
Death
Device embolism
Device erosion
Device malfunction
Device malposition
Device migration
Device-related thrombus
Fever
Foreign body embolization
Hematuria
Hypertension/hypotension
Infection
Multi-organ failure
Myocardial infarction
Perforation
Pericardial effusion
Pleural effusion
Pericarditis
Renal failure/dysfunction
Respiratory failure
Seizure
Significant residual flow
Stroke/transient ischemic attack
TEE complications
Thrombocytopenia
Thromboembolism (including peripheral and pulmonary thromboembolism)
Thrombus formation
Valvular regurgitation/insufficiency
Vascular access site injury (hematoma, pseudoaneurysm, arteriovenous fistula, groin pain, or other)
Vessel trauma/injury

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Directions for Use

Materials recommended for use with the Amplatzer™ Amulet™ Left Atrial Appendage Occluder

- 0.035-inch guidewire
- 35-cc syringe
- Delivery sheath with inner diameter appropriate for device size (see Table 1 in Appendix A: Supplemental Information)

NOTE: The Amplatzer™ Amulet™ Left Atrial Appendage Occluder is intended for use with an Amplatzer™ Amulet™ Delivery Sheath or an Amplatzer™ Steerable Delivery Sheath.

NOTE: Do not oversize or undersize the delivery sheath. Follow the recommended sizes listed in Table 1 in Appendix A: Supplemental Information.

Procedure

1. Evaluate the patient using transesophageal echocardiography (TEE) or intracardiac echocardiography (ICE) to rule out the presence of intracardiac thrombus (including left atrial appendage thrombus) and presence of pericardial effusion.
2. Prepare the patient for a standard transcatheter procedure, potentially including administration of oral antiplatelet medication. Effective anticoagulation therapy should be maintained during and after the procedure as determined by the physician. After a transseptal puncture is performed, maintain a recommended activated clotting time (ACT) of 250-350 seconds until the procedure is complete.
3. Perform a transseptal puncture.
4. Gain access to the left atrium and place the guidewire into the left upper pulmonary vein.

 CAUTION: Use caution when advancing the guidewire.
5. Use angiography, TEE (preferably 3D), and/or pre-procedural cardiac CT to measure the left atrial appendage, including the depth of the left atrial appendage (shown as W in Table 2 of Appendix A: Supplemental Information) and the maximum width of the orifice (shown as U in the graphic in Table 2 of Appendix A). Image the left atrial appendage until it is clearly visible.
 - Identify and measure the left atrial appendage at the landing zone (defined as about 10 mm-12 mm from the orifice) for the device lobe (shown as V in Table 2 in Appendix A: Supplemental Information) to determine the appropriate device size to occlude the left atrial appendage. Confirm the orifice measurement (U) is less than the disc size of the selected device.
 - See Table 2 in Appendix A to determine the appropriate device size to occlude the left atrial appendage. Procedures guided by ICE should be preceded by a TEE (preferably 3D) or cardiac CT to measure the left atrial appendage.

 WARNING: Do not implant the device if the measurements of the left atrial appendage do not fall within the sizing chart in Table 2 in Appendix A.

 NOTE: Ensure the measurement system is appropriately calibrated.
6. Prepare the delivery sheath for use according to the manufacturer's instructions for use.
7. Prepare the Amplatzer™ Amulet™ Left Atrial Appendage Occluder and delivery assembly for use:
 - Inspect the sterile pouch and verify that it is unopened and undamaged. Do not use the device or delivery assembly if the sterile barrier has been compromised.
 - Gently open the sterile packaging and inspect the components for damage. Do not use damaged or kinked components.
 - Remove the delivery cable and device from packaging.
8. If desired, slide the delivery cable vise over the delivery cable. Tighten the delivery cable vise onto the delivery cable.
9. Connect the hemostasis valve to the loader hub and turn the rotating luer to tighten.
10. Advance the tip of the device to the tip of the loader.
11. Ensure the device is attached to the delivery cable by rotating the delivery cable/delivery cable vise clockwise until resistance is felt.
12. Fill a 35-cc syringe with sterile saline. Connect the syringe to the hemostasis valve. Submerse device in the sterile saline. Flush sterile saline through the loader and simultaneously tighten the hemostasis valve onto the delivery cable. Fill the syringe and flush sterile saline through the loader two more times.

 NOTE: During additional flushing ensure air does not enter the system.

 NOTE: The loader tubing may be massaged to help remove air.

 NOTE: During steps 12–16 make sure the tip of the loader is submerged in sterile saline at all times so air does not enter the loader.
13. Fully retract the syringe plunger while holding the syringe body. Allow time for sterile saline to completely fill the syringe.
14. Advance the syringe plunger to force sterile saline through the loader.

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NOTE: During this process, a small amount of air may accumulate in the syringe. Ensure the syringe is held upright and no air is introduced back into the loader.

15. Flush the device again using steps 13 and 14.

 NOTE: If air is visible in the loader, additional flush steps may be performed as described in steps 12-15 or the device may be removed from the loader and disconnected from the delivery cable. Reload the device by performing steps 34–43.

 WARNING: Do not retract the device.

16. Retract the plunger to fill the syringe with sterile saline to accommodate fluid-to-fluid connection with the delivery sheath.
17. Advance the dilator and delivery sheath over the guidewire to the entrance of the left atrial appendage. Retract the dilator until it is inside the sheath tip. Position the sheath tip slightly beyond the landing zone of the left atrial appendage.

 WARNING: Use caution when advancing the dilator and sheath to prevent trauma.

18. Slowly remove the dilator from the delivery sheath. Remove the guidewire from the patient.
19. If a 14 Fr delivery sheath was selected for use with a 16 mm–25 mm device, tightly connect the 14 Fr sheath adaptor to the delivery sheath hub.

 WARNING: Proximal adaptor end is only compatible with supplied loader fitting. Do not flush or aspirate with any other accessories.

20. Allow blood backflow to purge all air from the delivery sheath assembly by holding the hub below the plane of the patient's body.
21. Immediately connect the loader to the delivery sheath in a fluid-to-fluid fashion and tighten the rotating luer.
22. Slightly loosen the hemostasis valve. Introduce the device into the delivery sheath and advance the device to the tip of the delivery sheath. Look for air bubbles when advancing. If air is present, allow backflow to purge air from the delivery sheath assembly.

 WARNING: Aspirating is not recommended because air can be introduced into the system.

 WARNING: Do not rotate the delivery cable when advancing the device through the delivery sheath.

23. Hold the delivery cable in place while retracting the delivery sheath to partially expose the lobe (see Figure 3). If necessary, adjust the position of the sheath and lobe for alignment.
24. Use fluoroscopy and/or echocardiography to guide the device into the left atrial appendage.

 WARNING: If the device needs to be retracted after this point, the device and sheath must both be removed and replaced.

 NOTE: Delivery cable connection can be checked by rotating the delivery cable/delivery cable vise clockwise until resistance is felt.

25. Continue to deploy the lobe by pushing the cable forward and/or pulling the sheath back. Stop when the lobe is fully deployed within the left atrial appendage at the intended landing zone.

 CAUTION: Do not advance the delivery cable or the sheath after the lobe is fully deployed.

26. The device lobe should be perpendicular with the axis of the left atrial appendage at the intended landing zone (see Figure 4).
27. Maintain a slight tension on the delivery cable while retracting the delivery sheath to expose the disc (see Figure 4). The device disc should cover the orifice. Do not leave a gap between the orifice and the disc (see Figure 5: Arrows in this figure identify possible gaps).
28. Confirm proper device placement using echocardiography and fluoroscopy. Indications of proper device placement include:
 - At least 2/3 of the device lobe should be distal to the left circumflex artery on echocardiography. (see Figure 6).
 - The device lobe should be slightly compressed and have good apposition to the left atrial appendage wall.
 - The orientation of the device lobe must be in line with the axis of the intended landing zone in the left atrial appendage (see Figure 4, Figure 5, and Figure 6).
 - The disc must be separated from the lobe (see Figure 4).
 - The disc will have a concave shape.
29. If device repositioning is required:
 - Partially recapture the device by holding the delivery cable in place while readvancing the delivery sheath until the platinum thread marker is aligned with the distal edge of the delivery sheath's marker band (see Figure 3). Do not retract the device any farther to prevent damaging the delivery sheath.

 WARNING: If the device is retracted farther than the radiopaque markers, do not readvance the device or perform injections. The device and the delivery sheath must both be removed and replaced.

 - Reposition and redeploy the device (see steps 23–28).

 NOTE: The device can be partially recaptured and redeployed a maximum of three times. If the device position is still unsatisfactory, then remove and replace both the device and the sheath.

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30. Perform echocardiography and fluoroscopy to confirm that the device is in place.
31. When proper device placement is confirmed the device may be released:
 - Detach the device by turning the delivery cable/delivery cable vise counter-clockwise.

 CAUTION: Do not exceed eight rotations. If difficulty occurs while attempting to release the device, stop and release the built up tension. Advance the cable until it is perpendicular with the device and attempt to release.
32. Remove the delivery cable and delivery sheath from the patient.
33. Perform echocardiography and fluoroscopy to confirm that the device remains in place.

Device Reloading Procedure

34. Insert loading cable through the distal loader tip.
35. Attach distal end screw to loading cable.
36. Rotate the loading cable vise clockwise until resistance is felt, ensuring that the loading cable is secured to the device.
37. Immerse the device and hub end of loader in sterile saline to remove air from both components. Device may be massaged to help remove air.
38. Pull the loading cable vise until the lobe is fully retracted within the loader, but stop before the disc is loaded.
39. Connect the delivery cable to the exposed proximal end screw of the device. Grasp the loader and device and rotate the delivery cable clockwise to thread the device onto the delivery cable. Stop rotating the delivery cable when resistance is felt.
40. Immerse the device in sterile saline.
41. Pull the loading cable vise until the device is loaded into the proximal end of the loader.

 NOTE: Do not pull the device to the distal end of the loader.
42. Remove the loading cable by rotating the loading cable vise counterclockwise.
43. Return to Step 9.

Post-procedure Instructions

- Monitor the patient overnight. Perform a transthoracic echocardiogram (TTE) to make sure the device is in the correct position and no pericardial effusion is present before the patient is discharged.
- Aspirin (or an alternate antiplatelet agent) is recommended for patients for six months post-implant. The decision to continue antiplatelet therapy (such as aspirin alone) after six months is at the discretion of the physician.
- Clopidogrel (or an alternate antiplatelet) is recommended for patients and prescription should follow routine standard of care.
- Oral anticoagulation is not recommended unless residual flow around the device is >5 mm.
- Appropriate endocarditis prophylaxis is recommended.
- Instruct the patient when to seek medical attention.
- Provide routine standard of care follow-up, including echocardiography, for evaluation of residual shunt and adverse events (including thrombus formation and pericardial effusion).

Reporting Device Incidents

If, in the course of using this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Disposal

Amplatzer™ Amulet™ Left Atrial Appendage Occluder components should be appropriately classified as biohazards and disposed of in compliance with applicable facility procedures and local and country laws and regulations.

Warranty

Abbott Medical warrants to buyer that, for a period equal to the validated shelf life of the product, this product shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Abbott Medical's obligation under this warranty is limited to replacing or repairing at its option, at its factory, this product if returned within the warranty period to Abbott Medical and after confirmed to be defective by the manufacturer.

EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY AS TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

See the Terms and Conditions of Sale for further information.

Australian Warranty

This section is only provided in English.

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This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or Australian Sponsor.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical at the address provided with the product and after Abbott Medical confirms that the goods are defective.

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical’s liability for the breach of the term or guarantee is limited to, at Abbott Medical’s option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED.

See the Terms and Conditions of Sale for further information.

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1. Hildick-Smith D, Landmesser U, Camm AJ, et al. „Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study.“ European Heart Journal 41, no. 30 (2020): 2894–2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, et al. „Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation.“ JACC Cardiovascular Interventions 14, no. 1 (2021): 69–78.
3. Chen S, Julian Chun K.R., Ling Z, et al. „Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update.“ Journal of Cardiovascular Development and Disease 8, no. 6 (2021): doi: 10.3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Epub ahead of print. PMID:34459659.

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1. Hildick-Smith D, Landmesser U, Camm AJ, et al. “Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study.” European Heart Journal 41, no. 30 (2020): 2894-2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, et al. “Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation.” JACC Cardiovascular Interventions 14, no. 1 (2021): 69-78.
3. Chen S, Julian Chun K.R., Ling Z, et al. “Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update.” Journal of Cardiovascular Development and Disease 8, no. 6 (2021): doi: 10.3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Elektronische Veröffentlichung vor dem Druck. PMID:34459659.

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1. Hildick-Smith D, Landmesser U, Camm AJ, et al. “Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study.” European Heart Journal 41, no. 30 (2020): 2894-2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, et al. “Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation.” JACC Cardiovascular Interventions 14, n.º 1 (2021): 69-78.
3. Chen S, Julian Chun K.R., Ling Z, et al. “Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update.” Journal of Cardiovascular Development and Disease 8, n.º 6 (2021): doi: 10.3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Epub ahead of print. PMID:34459659.

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1. Hildick-Smith D, Landmesser U, Camm AJ, et al. "Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study." European Heart Journal 41, no. 30 (2020): 2894-2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, et al. "Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation." JACC Cardiovascular Interventions 14, no. 1 (2021): 69-78.
3. Chen S, Julian Chun K.R., Ling Z, et al. "Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update." Journal of Cardiovascular Development and Disease 8, no. 6 (2021): doi: 10.3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Epub ahead of print. PMID:34459659.

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1. Hildick-Smith D, Landmesser U, Camm AJ, et al. “Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study.” European Heart Journal 41, no. 30 (2020): 2894-2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, et al. “Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation.” JACC Cardiovascular Interventions 14, no. 1 (2021): 69-78
3. Chen S, Julian Chun K.R., Ling Z, et al. “Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update.” Journal of Cardiovascular Development and Disease 8, no. 6 (2021): doi: 10.3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Epub ahead of print. PMID:34459659.

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1. Hildick-Smith D, Landmesser U, Camm AJ, et al. “Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study.” European Heart Journal 41, no. 30 (2020): 2894–2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, et al. “Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation.” JACC Cardiovascular Interventions 14, no. 1 (2021): 69–78.
3. Chen S, Julian Chun K.R., Ling Z, et al. “Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update.” Journal of Cardiovascular Development and Disease 8, no. 6 (2021): doi: 10.3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Epub ahead of print. PMID:34459659.

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1. Hildick-Smith D, Landmesser U, Camm AJ. i in. „Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study.” European Heart Journal 41, nr 30 (2020): 2894–2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, i in. „Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation” JACC Cardiovascular Interventions 14, nr 1 (2021): 69–78.
3. Chen S, Julian Chun K.R., Ling Z, et al. „Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update.“ Journal of Cardiovascular Development and Disease 8, no. 6 (2021): doi: 10.3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Publikacja elektroniczna poprzedzająca wydanie wersji drukowanej. PMID:34459659.

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1. Hildick-Smith D, Landmesser U, Camm AJ, et al. “Left Atrial Appendage Occlusion with the Amplatzer Amulet Device: Full Results of the Prospective Global Observational Study.” European Heart Journal 41, no. 30 (2020): 2894-2901.
2. Nielsen-Kudsk JE, Kasper K, Damgaard D, et al. “Clinical Outcomes Associated With Left Atrial Appendage Occlusion Versus Direct Oral Anticoagulation in Atrial Fibrillation.” JACC Cardiovascular Interventions 14, no. 1 (2021): 69-78.
3. Chen S, Julian Chun K.R., Ling Z, et al. “Comparison of Left Atrial Appendage Occlusion versus Non-Vitamin-K Antagonist Oral Anticoagulation in High-Risk Atrial Fibrillation: An Update.” Journal of Cardiovascular Development and Disease 8, no. 6 (2021): doi: 10,3390/jcdd8060069
4. Lakkireddy D, Thaler D, Ellis CR, Swarup V, Sondergaard L, Carroll J, Gold MR, Hermiller J, Diener HC, Schmidt B, MacDonald L, Mansour M, Maini B, O'Brien L, Windecker S; Amulet IDE Investigators. AMPLATZER™ AMULET™ LEFT ATRIAL APPENDAGE OCCLUDER VERSUS WATCHMAN™ DEVICE FOR STROKE PROPHYLAXIS (AMULET IDE): A RANDOMIZED CONTROLLED TRIAL. Circulation. 2021 Aug 30. doi: 10.1161/CIRCULATIONAHA.121.057063. Epub ahead of print. PMID:34459659.

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™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third party trademark, which is property of its respective owner.
© 2022 Abbott. All Rights Reserved.
Abbott Medical UK Limited