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1 Device description


1.1 Specifications

The Avalus bioprosthesis, Model 400, consists of a polyester covered base frame and trileaflet support frame structure shaped to the geometry of a trileaflet valve. The base frame and trileaflet support frame are injection-molded using Polyetheretherketone (PEEK) material. The base frame PEEK material is impregnated with barium sulphate to allow for radiographic visualization. The leaflets are laser cut from bovine tissue that have been cross-linked in buffered glutaraldehyde. The leaflets are inserted between the trileaflet support frame and the base frame, and then all components are securely sutured together. A sewing ring, fabricated from polyester cloth, is integrated into the inflow base frame cover to allow the user to suture and seat the Avalus bioprosthesis in the supra-annular position. Sewing markers are located on the sewing ring in the mid sinus area of each cusp to provide guidance for even spacing of the implant sutures (*Figure 1*). A disposable holder is attached to the outflow of the valve to facilitate implantation. The aortic bioprosthesis (valve) is treated with an alpha amino oleic acid (AOA) antimineralization process that has been shown to mitigate porcine leaflet calcification in animal studies.

The disposable holder is designed to fit the reusable Medtronic valve handle, Model 7420. The holder features a single cut point to remove the holder from the valve. Double-ended sizers are used to select the appropriate valve size. The barrel end represents the valve orifice. The replica end imitates the prosthesis geometry.

Table 1. Available sizes and sewing ring diameters (nominal values in millimeters)

Size	19 mm	21 mm	23 mm	25 mm	27 mm	29 mm
(1) Stent diameter (TAD)	19	21	23	25	27	29
(2) Internal orifice diameter	17.5	19.5	21.5	23.5	25.5	27.5
(2a) Internal orifice diameter[a]	18	20	22	24	26	28
(3) External sewing ring diameter	27.0	29.0	31.0	33.0	36.0	38.0
(4) Valve profile height	13.0	14.0	15.0	16.0	17.0	18.0
(5) Aortic protrusion	11.0	12.0	13.0	14.0	15.0	16.0

Note: TAD—tissue annulus diameter

[a] Measurement shows stent frame including tissue (2) and stent frame excluding tissue (2a)

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Table 2. Heart valve materials

Component	Materials in contact with the body	Surface area[a]
Tissue	Processed bovine pericardium	35 cm^2
Frame/Polymer stent	Polyetheretherketone (PEEK), impregnated with barium sulfate	18 cm^2
Cloth	Polyester	60 cm^2
Suture	Polyester and ultra-high-molecular-weight polyethylene (UHMWP)	7 cm^2

[a] Surface area is based on the material in the largest size valve.

The aortic bioprosthesis (valve) is treated with an alpha amino oleic acid (AOA) antimineralization process that has been shown to mitigate porcine leaflet calcification in animal studies. The Avalus heart valve is stored in a 0.2% glutaraldehyde solution. As part of the procedure, it is rinsed prior to implantation, reducing glutaraldehyde residuals to safe levels.

1.2 Avalus system and accessories

The Medtronic valve handle, Model 7420, is used to aid in the implantation of the bioprosthesis. Use only Avalus sizers, Model 7400S, to determine the appropriate Avalus bioprosthesis size.

Caution: Do not use the sizers or handles until they have been thoroughly cleaned and sterilized. Refer to the appropriate instructions for use for further instructions.

Caution: Do not use other manufacturers' valve sizers, or sizers from another Medtronic prosthesis, to size the Avalus bioprosthesis.

1.3 Expected device lifetime

The device lifetime of an Avalus bioprosthesis varies from patient to patient. At some point, the valve may need to be replaced. The expected lifetime of an Avalus bioprosthesis is projected to be a minimum of ten years for many patients based on clinical evidence and testing according to industry standards.

2 Intended purpose

The intended purpose of the Avalus bioprosthesis is to restore unidirectional blood flow through the aortic annulus by replacing a diseased or damaged heart valve.

2.1 Indications for use

The Avalus bioprosthesis is indicated to replace damaged or diseased aortic valves with moderate or greater stenosis or regurgitation as assessed by echocardiography, consistent with EACTS guidelines.

2.2 Contraindications

Patients who are not candidates for surgical AVR per contemporary EACTS guidelines.

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2.3 Patient target group

The Avalus bioprosthesis is intended for adult patients with an indication for bioprosthetic aortic valve replacement per contemporary EACTS guidelines.

2.4 Intended users

Device implantation should only be performed by medical professionals trained in aortic valve replacement techniques. It is the responsibility of the hospital to ensure that the individuals assigned to operate this device are thoroughly trained in its use and alert to potential problems.

Clinicians implanting surgical valves should be experienced in aortic valve replacement and review the provided instructions for use before surgery. Medtronic does not require that specific training be provided to surgeons, but does offer aortic valve instruction and training materials to surgeons and medical professionals. Medtronic trains field personnel to ensure they understand aortic valve replacement before providing clinicians information on aortic valve replacement.

Prior to implant, patients may receive educational training materials on valve replacement therapy.

2.5 Expected clinical benefits

The clinical benefits to patients from Avalus bioprosthesis implantation include the following:

1. Improved quality of life through the relief of symptoms related to severe aortic stenosis or regurgitation from a native or failed surgical valve.
2. Positive impact on hemodynamic function as evaluated by echocardiogram.
3. An increased chance of survival.

The Summary of Safety and Clinical Performance (SSCP) for this device can be found at https://ec.europa.eu/tools/eudamed. Search for the SSCP using the device name and the following elements as applicable: device model, reference number, and catalog number.

3 Warnings and precautions

3.1 Warnings

- **Read all warnings, precautions, and instructions for use carefully before use. Failure to read and follow all instructions or failure to observe all stated warnings could cause serious injury to the patient.**
- This device was designed for single use only. Do not reuse, reprocess, or resterilize this device. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination, which could result in patient injury, illness, or death.
- Use aseptic technique in all procedures. Failure to follow aseptic techniques could create a risk of contamination of the device, which could result in patient injury, illness, or death.

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- As with any implanted medical device, there is potential for patient immunological response, including an allergic response. Care should be exercised in patients with hypersensitivities to the device materials outlined in Section 1.1, Specifications.
- Check the shipping temperature indicator inside the carton. If the shipping temperature indicator window is black, the valve is not suitable for clinical use.
- Do not use the valve in any of the following circumstances:
 - The valve has been dropped, damaged, or mishandled in any way
 - Particulate matter is present on the valve
 - The Use-by date has elapsed
 - All tamper strips are damaged
 - The serial number tag does not match the number on the container label
 - The glutaraldehyde storage solution does not completely cover the valve or fluid is leaking from the packaging
 - The valve shows signs of deterioration or discoloration
 - The sewing ring marker is not visible
- Do not handle the tissue portion of the valve with instruments. Extreme care must be taken to prevent damage to the delicate valve tissue. Even a minor perforation could enlarge over time to significantly impair valve function.
- Do not place the nonsterile exterior of the valve container in the sterile field.
- Do not add antibiotics or any other substance to the storage solution or the rinse solution. This could result in an inflammatory response or regurgitation.
- Do not oversize. Implanting too large a valve may cause procedural complications such as subannular obstruction, deformation of the valve, coronary occlusion, or damage to surrounding tissues. These complications could lead to prolonged procedure time, compromised hemodynamics, reduced valve life, or prosthetic valve explant.
- Do not undersize. Implanting too small a valve may cause stenosis or paravalvular leak.
- Do not use the Avalus sizers, Model 7400S, or the Medtronic valve handle, Model 7420, until they have been thoroughly cleaned and sterilized. Refer to the appropriate instructions for use for further instructions.
- Do not use other manufacturers' valve sizers, or sizers for another Medtronic prosthesis, to size the valve.
- Do not allow the valve to dry. To prevent the valve tissue from drying, periodically irrigate the valve with sterile, normal saline during implantation.

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- Carefully consider these potential hazards when selecting an appropriate valve substitute:
 - Calcific degeneration could cause accelerated deterioration of the valve in patients with altered calcium metabolism (for example, chronic renal failure, hyperparathyroidism).
 - Calcification may occur earlier in children, adolescents, or young adults.
 - Premature calcification may also occur in older adults who accept a biologic prosthesis.
 - Patients with a bioprosthesis that require chronic anticoagulation are at additional risk of bleeding.
 - Stenosis and regurgitation of the bioprosthesis may occur in patients with coagulation disorders such as AT3 deficiency.
 - Paravalvular leak is more likely to occur in patients with aneurysmal aortic or degenerative conditions, cystic medial necrosis, or Marfan's syndrome.
- Do not implant a valve without completing proper rinsing procedures. Patient exposure to excessive glutaraldehyde may cause infection or inflammatory response.
- Avoid contact of the valve or the rinse solution with towels, linens, or other sources that may transfer particulate matter to the leaflet tissue.
- Use caution when placing sutures through the sewing ring to avoid laceration of the leaflet tissue. If the leaflet tissue is damaged, the valve must be explanted and replaced.
- Do not excessively bend the bioprosthesis commissures when tying knots.
- Do not handle the leaflet tissue.

3.2 Precautions

- Avoid prolonged contact with the valve storage solution. Glutaraldehyde could cause irritation of the eyes, nose, skin, and throat if continued exposure occurs. Avoid prolonged exposure to, or breathing of, the chemical vapor. Use only with adequate ventilation. If skin contact occurs, immediately flush the affected area with copious amounts of water for 10 to 15 minutes. If eye contact occurs, flush the eye with water for 15 minutes and seek immediate medical attention.
- Dispose of waste substances in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

4 Potential adverse events

Adverse events potentially associated with the implantation and use of a bioprosthetic heart valve, or with the use of related accessories during the implantation procedure, include the following:

- Angina
- Aortic tissue damage
- Cardiac dysrhythmias
- Embolism
- Endocarditis

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- Heart failure
- Hemolysis
- Hemolytic anemia
- Hemorrhage, anticoagulant/antiplatelet-related
- Immunological response, including allergic response
- Inflammatory reaction
- Infection other than endocarditis
- Leak, transvalvular or paravalvular
- Myocardial infarction
- Nonstructural dysfunction (leaflet entrapment/impingement, obstructive pannus ingrowth, suture dehiscence, inappropriate sizing)
- Pericardial effusion
- Pleural effusion
- Prosthesis regurgitation
- Prosthesis stenosis
- Stroke
- Structural deterioration (calcification, leaflet tear)
- Tamponade
- Valve thrombosis

These complications could lead to:

- Reoperation
- Explant of the bioprosthesis
- Permanent disability or organ damage
- Death

If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

5 Instructions for use

5.1 Handling and preparation instructions

Proper size selection is critical to heart valve replacement. Use Avalus sizers, Model 7400S, to select the appropriately sized Avalus bioprosthesis. For further information, refer to the instructions for use.

The exterior of the container and lid are nonsterile. The valve and all packaging components inside the container are sterile and must be handled within the sterile operative field.

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1. Prepare 2 rinse basins each containing 500 mL of sterile, normal saline solution.
2. Examine the tamper strips to verify that the container has not been damaged or previously opened. Do not use the valve if all of the tamper strips are damaged.
3. Turn the lid counterclockwise and open the container (*Figure 2*).
4. Hold the retainer with the thumb and index finger. Slowly lift the valve out of the container, allowing for drainage of the glutaraldehyde storage solution (*Figure 3*).
5. Verify that the serial number on the retainer matches the serial number on the container lid, shelf carton, and Patient Registration Form. Record the serial number in the patient's record using the stickers provided on the Patient Registration Form.

 Caution: Failure to record the serial number accurately could result in a loss of traceability.
6. Hold the retainer upright. Remove the retainer cap by turning it counterclockwise using the thumb and index finger (*Figure 4*). The holder will be visible.
7. Insert a sterile handle into the holder. Rotate the handle clockwise into the threaded opening of the holder until resistance is felt (*Figure 5*).
8. Pull upward on the handle to remove the valve from the retainer (*Figure 6*).

5.2 Rinse procedure

1. Place the entire valve and holder in 1 of the rinse basins (*Figure 6*).
2. Use the handle to continually agitate the valve for a minimum of 30 seconds.
3. Gently squeeze the sewing ring to remove any residual glutaraldehyde. Do not touch the tissue portion of the valve.

 Warning: Take care not to damage the valve when squeezing. Distorting or damaging the valve could lead to regurgitation.
4. Place the valve in the second basin. Keep the valve in the basin until required by the surgeon.

5.3 Aortic valve implantation

The Avalus bioprosthesis is designed for implantation in the supra-annular position. Use Avalus sizers, Model 7400S, to visualize placement of the sewing cuff above the annulus and to confirm placement and fit of the valve in the supra-annular space.

The choice of surgical technique is left to the discretion of the individual surgeon when implanting in the supra-annular position. Any modified technique should be in accordance with *Chapter 3, Warnings and precautions, page 13*. In general, follow these steps:

1. Orient the valve so that the stent posts are pointing up and the coronary ostia are not obstructed by the stent and sewing ring.

 Warning: Incorrect positioning of the valve could result in serious injury or death.
2. Irrigate the valve periodically with sterile, normal saline to prevent drying of the delicate tissue.
3. Position the valve in the annulus and place the sutures in the sewing ring.

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Warning: Take care when lowering the valve into place with the handle. Pushing too hard on the handle could cause damage to the annulus.

4. Suture the valve in place using an appropriate suture technique. Refer to *Chapter 3, Warnings and precautions, page 13* to avoid potential problems.
5. Cut the sutures at the single cut point indicated by the white arrow to release the holder (*Figure 7*). Gently pull the handle away from the valve to remove the holder.

 Warning: Failure to remove the valve holder could result in serious injury, including stenosis or congestive heart failure.
6. Remove the holder from the handle and discard the holder.
7. Tie all knots and trim suture tails.

 Warning: Do not tie the sutures into the valve holder sutures.

 Warning: Take care to cut the sutures close to the knots. Exposed suture tails can damage leaflet tissue.

 Warning: Verify that no suture remnants are present, removing them before completing the valve implantation. Allowing suture remnants to remain on the valve could lead to patient embolism.

6 Individualization of treatment

Except where contraindicated, during the initial stages after implantation, patients receiving a bioprosthetic heart valve should be maintained on anticoagulant therapy.

Consider long-term anticoagulant or antiplatelet therapy in patients with a dilated left atrium, a history of thromboembolic events, or a cardiac rhythm of atrial fibrillation or atrial flutter.

7 Patient counseling information

In accordance with local regulations, healthcare providers should review the instructions for use for applicable information to be shared with the patient. A patient implant card, which contains identifying information about the implanted device, is included in the device package. After device implant, complete the patient implant card and provide it to the patient before they are discharged. Healthcare providers should communicate the following instructions to their patients:

- Always carry their implant card with them.
- Access additional information about their device on the website that is listed on their patient implant card.

 Note: If the patient is unable to access the website, the physician must provide the information from the website to the patient.
- Always inform any healthcare personnel that they have an implanted device before any procedure has begun.
- Contact their healthcare provider if they notice any new or changing symptoms.

Continuous medical follow up (at least annually) is recommended so that complications related to the valve can be diagnosed and properly managed.

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Patients with a bioprosthetic valve are at risk for bacteremia (for example, undergoing dental procedures). Advise the patient about prophylactic antibiotic therapy.

8 How supplied

8.1 Packaging

The Avalus bioprosthesis is chemically sterilized and is supplied sterile in a buffered 0.2% glutaraldehyde solution. Sterility is compromised if the glass jar-and-lid container is opened or damaged. The outside of the container is not sterile. Do not place it in the sterile field.

8.2 Storage

Store the Avalus bioprosthesis between 5°C and 25°C (41°F and 77°F). Refrigeration is not required, and freezing could damage the valve. Room temperature storage up to 25°C (77°F) is satisfactory. Do not expose the valve to sunlight or other ultraviolet light sources, or placed where significant temperature fluctuations could occur.

Maintain appropriate inventory control so that bioprostheses with earlier Use-by dates are implanted first to avoid expiration dates.

8.3 Return of explanted bioprostheses

Contact a Medtronic sales representative to request a product return kit to return explanted bioprostheses for analysis and disposal. If a kit is not available, place the explanted bioprosthesis in a container filled with glutaraldehyde or 10% buffered formalin immediately after explantation. For further return instructions, contact a Medtronic sales representative.

9 Patient information

9.1 Registration information

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the European Union.

A Patient Registration Form is included in each device package. After implantation, please complete all requested information. Return the original form to the Medtronic address indicated on the form. Provide the patient with the temporary identification card before discharge.

An Implanted Device Identification Card is provided to the patient. This card contains the name and telephone number of the patient's physician, as well as information that medical personnel would require in the event of an emergency.

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10 Postoperative information

10.1 Magnetic resonance imaging (MRI) compatibility

The Avalus bioprosthesis Model 400, is magnetic resonance (MR) safe. The valve contains no metal and, therefore, poses no known hazards in all MR environments.

MR scanning at 3.0 T and 1.5 T can be performed immediately after implantation. The valve will not cause any harm to the patient when exposed to MR scanning.

10.2 Image artifact

MR image quality could be compromised if scanning on or near the position of the valve. If necessary, optimize the MR imaging parameters for the presence of the valve.

11 Disclaimer of warranty

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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© 2022 Medtronic
Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (“TM*”) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
Medtronic, Inc.
Medtronic B.V.