Medtronic

Avalus™ 7400S

Sizers
Сепаратори
Měrky
Sizere
Obturatoren
Μετρητές μεγέθους
Medidores
Kalibraatorid
Sovittimet
Calibreurs
Alati za određivanje veličine
Méretezők
Pengukur
Misuratori
Өлшегіштер
Kalibratoriai
Izmēra noteicēji
Калибратори
Sizers
Målere
Przymiary
Calibradores
Calibradores
Calibroare
Измеритель диаметра клапана Avalus
Meradlá
Naprave za določanje velikosti
Regulatori veličine
Storleksmätare
Boyut belirleyiciler
Вимірювачі
Dụng cụ định cỡ

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Petunjuk Penggunaan • Istruzioni per l'uso • Қолдану жөніндегі нұсқаулар • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по применению • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання • Hướng dẫn sử dụng

! USA **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

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Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States or in other countries.

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1 Device description

1.1 Specifications

The Avalus™ sizers, Model 7400S, are manufactured using transparent polysulfone and nitinol to enable direct observation of the bioprosthetic valve profile within the patient's native valve annulus. Each sizer has a barrel and replica head connected by a malleable nitinol wire (*Figure 1*).

Sizers are provided for each size of the Avalus™ bioprosthesis. The size is indicated on the handle and barrel head.

The size is equal to the center of the valve's trileaflet support frame and the outer diameter of the barrel head (*Figure 2*).

The replica head represents the valve profile with the sewing cuff folded upward in its implanted state when the sutures are placed close to the stent (*Figure 2*). The valve orifice area is represented by the inner diameter of the replica head. The stent posts of the replica head represent the geometry of the valve stent post outer diameter, height, and width.

1.2 Avalus system and accessories

Table 1. Parts and accessories available separately

Model number	Description
7420	Medtronic valve handle
400	Avalus bioprosthesis
T7400	Avalus accessory tray

2 Intended purpose

The intended purpose of the Avalus sizers 7400S is to facilitate the implantation of the Avalus bioprosthesis.

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2.1 Indications for use

The Avalus sizers are intended to aid the surgeon in the selection of the optimal size of the Avalus bioprosthesis.

2.2 Contraindications

The Avalus sizers are contraindicated for use with heart valve products other than the Avalus bioprosthesis.

2.3 Patient target group

This device is intended for patients who are candidates for a tissue prosthetic valve based on a physician evaluation.

2.4 Intended users

Device implantation should only be performed by medical professionals trained in aortic valve replacement techniques. It is the responsibility of the hospital to ensure that the individuals assigned to operate this device are thoroughly trained in its use and alert to potential problems.

Clinicians implanting surgical valves should be experienced in aortic valve replacement and review the provided instructions for use before surgery. Medtronic does not require that specific training be provided to surgeons, but does offer aortic valve instruction and training materials to surgeons and medical professionals. Medtronic trains field personnel to ensure they understand aortic valve replacement before providing clinicians information on aortic valve replacement.

Prior to implant, patients may receive educational training materials on valve replacement therapy.

2.5 Expected clinical benefits

The clinical benefits to patients from Avalus bioprosthesis implantation include the following:

1. Improved quality of life through the relief of symptoms related to severe aortic stenosis or regurgitation from a native or failed surgical valve.
2. Positive impact on hemodynamic function as evaluated by echocardiogram.
3. An increased chance of survival.

The Summary of Safety and Clinical Performance (SSCP) for this device can be found at https://ec.europa.eu/tools/eudamed. Search for the SSCP using the device name and the following elements as applicable: device model, reference number, and catalog number.

3 Warnings and precautions

3.1 Warnings

- **Read all warnings, precautions, and instructions for use carefully before use. Failure to read and follow all instructions or failure to observe all stated warnings could cause serious injury to the patient.**
- Carefully inspect the device before each use for cracks or flaws that could result from sterilization, handling, or general use. Ensure the size is legible. Failure to discard a device that shows signs of cracking, corrosion, or degradation could result in patient injury.
- The reusable sizers are provided nonsterile and must be thoroughly cleaned and sterilized before each use. Failure to follow these Instructions for Use or hospital-approved cleaning and sterilization instructions could result in patient infection.
- If the sizers fall outside of the sterile field, clean and resterilize them before use. Follow these Instructions for Use or hospital-approved cleaning and sterilization instructions.
- Do not oversize. Implanting too large a valve may cause procedural complications such as subannular obstruction, deformation of the valve, coronary occlusion, or damage to surrounding tissues. These complications could lead to prolonged procedure time, compromised hemodynamics, reduced valve life, or prosthetic valve explant.
- Do not undersize. Implanting too small a valve may cause stenosis or paravalvular leak.
- Do not use other manufacturers' sizers or sizers from other Medtronic heart valve products to size the Avalus bioprosthesis. Other sizers may not indicate the appropriate Avalus bioprosthesis size.
- Do not use solvent-based cleaners (for example, acetone or toluene) when cleaning valve sizers or accessory trays. Use water-soluble surfactant-type soaps or mild detergents as cleaning agents. The use of solvent-based cleaners could result in patient infection or damage to the sizers.
- **If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.**

3.2 Precautions

- Dispose of waste substances in accordance with applicable laws, regulations, and hospital procedures, including those regarding biohazards, microbial hazards, and infectious substances.

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4 Instructions for use

Use both ends of the sizer to determine the bioprosthesis size. Use the barrel head to determine the patient's annulus size and optimal valve size. Use the replica head to confirm the fit of the Avalus bioprosthesis and adjust the valve size selection as necessary.

1. Determine the patient's annulus size by finding the barrel head that fits comfortably in the patient's annulus without deforming the anatomy. Use the barrel end to determine the optimal valve size for the patient's annulus size.
2. Use the corresponding replica end to check for interference with the coronary ostia above the valve and whether the patient's aortic root allows this valve size to be implanted. If the potential exists for valve prosthesis interference with coronary ostia or surrounding structures, tissue obstruction to valve prosthesis inflow, or distortion of the valve prosthesis stent, consider reorientation of the valve, selecting a smaller valve size, or addressing the obstruction or interference surgically.
3. Confirm the appropriate valve size by finding the replica sizer head that fits comfortably in the patient root in the intended supra-annular position without obstruction or deforming the anatomy.

5 How supplied

5.1 Packaging

The sizers are supplied nonsterile. Thoroughly clean and sterilize the sizers before each use.

5.2 Available sizes

The sizers are provided in a complete set with the following sizes: 17 mm[1], 19 mm, 21 mm, 23 mm, 25 mm, 27 mm, and 29 mm.

5.3 Storage

The sizers should be kept in a clean, cool, and dry environment.

6 Processing and reprocessing

The following instructions have been validated by Medtronic for preparing this medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials, and personnel in the processing facility, achieves the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the following instructions should be properly evaluated for effectiveness and potential adverse consequences.

Warning: If within a European member state, consider specific national regulations related to reprocessing.

Warning: Failure to follow reprocessing instructions can decrease device durability and device safety.

Warning: Thoroughly clean the device to eliminate soil before sterilization. Failure to follow these Instructions for Use or hospital-approved cleaning and sterilization instructions could result in patient infection or damage to accessories.

- It is recommended that devices are reprocessed as soon as reasonably practical following use.
- Do not allow contaminated devices to dry before implementing any cleaning procedure.
- After use, in order to prevent organics from drying, place the device in a container with an enzymatic solution or enzymatic spray per the manufacturer's instructions until reprocessing can be performed.
- To ensure containment for safe transportation, hospital staff should reference their own protocols, including those regarding protection of the medical device, environment, health care personnel, and mode of transportation. Ensure adequate personal protective equipment (PPE) is used during the handling of medical devices.
- Examine the components for cracking or degradation before use. Do not process if components show signs of cracking or degradation.

6.1 Preparation before cleaning

Disassembly of simple assemblies is necessary to allow for more complete cleaning and sterilization per these Instructions for Use. Disassembly should not require any mechanical tooling.

6.2 Point of use (within a maximum of 1 hour)

Wipe down the device to remove blood and surgical debris. Reprocess as soon as possible following use. If reprocessing cannot be performed within an hour of use, submerge the device in a pH-neutral enzymatic solution (for example, ENZOL) or spray with a pH-neutral enzymatic spray per the manufacturer's instructions until reprocessing can be performed.

6.3 Cleaning and disinfection for countries within the EU

Warning: Do not use solvent-based cleaners.

[1] The 17 mm size is not available within the EU market.

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6.3.1 Automated cleaning and disinfection

Automated cleaning pretreatment

1. Immerse the devices in a covered container with an enzymatic detergent (for example, ENZOL) using warm tap water 22 °C to 43 °C (72 °F to 110 °F).
2. Clean each component thoroughly using a soft bristle brush (for example, nylon toothbrush).

 Note: Take care to clean areas where soil may be impacted or shielded from the cleaning process.
3. Rinse the devices in warm running tap water 22 °C to 43 °C (72 °F to 110 °F) for a minimum of 2 minutes or until all traces of the cleaning solution are removed.
4. Visually inspect the devices for any remaining soil.
5. Repeat the above steps, if needed.
6. Allow the devices to dry.
7. Place the device into a perforated instrument tray suitable for use in the washer-disinfector. Ensure the device is placed in such a manner that will protect it from unwanted motion or potential mechanical damage during the automated cycles.

 Warning: Do not clean in the accessory tray.

 Place the perforated instrument tray with the device into the washer-disinfector.

Automatic cleaning and disinfection

Using a washer-disinfector complying with the requirements of the ISO 15883 series, perform the process steps (*Table 2*).

Caution: Do not use instrument lubricant on the devices.

Table 2. Automatic washer-disinfector cleaning parameters

Process step	Time (min)	Time (min) for inactivation of prions	Temperature	Water quality	Cleaning agent
1 - Pre-wash	1	3	<35 °C (95 °F)	Cold tap water	N/A
2 - Alkaline wash	5	10	≥55 °C (≥131 °F)	Hot tap water or purified water. For prion inactivation, deionized or reverse osmosis (DI/RO) water	Alkaline detergent[a] diluted per the manufacturer's specifications
3 - Neutralizing wash	1	2	≥43 °C (≥109 °F)	Hot tap water or hot purified tap water	Neutralizer[b] diluted as required to meet neutralization needs
4 - Intermediate rinse	2	2	≥60 °C (≥140 °F)	Hot tap water or hot purified tap water	N/A
5 - Thermal rinse (disinfection)	5	5	≥90 °C (≥194 °F)	Hot purified water[c]	N/A
6 - Drying	15	15	High setting or ≥82 °C (≥180 °F)	N/A	N/A

[a] Tested with neodisher SeptoClean™, 10 mL/L (1%) and 1 mL/L (0.1%). pH of ≥11.5 (1% neodisher SeptoClean™) required in order to inactivate prions. Equivalent cleaner may be used. Ensure hospital procedures and regulations are followed.

[b] Tested with neodisher Z™, 1 mL/L (0.1%).

[c] Purified water with endotoxin levels less than 10 EU/mL and water microbial levels less than 10 CFU/mL.

Caution: Do not omit the neutralizing wash step (Step 3) of the automatic washing and disinfection process. Excluding this step could decrease the usable lifetime of the device.

6.3.2 Manual cleaning and disinfection

Manual cleaning

Medtronic has established the following manual cleaning procedure.

1. **Rinse:** Remove gross soil using running tap water for a minimum of 30 seconds (water temperature 25 °C - 40 °C (77 °F - 104 °F)).

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2. **Soak:** Submerge the device in an alkaline[2] or pH-neutral enzymatic[3] bath diluted as specified by the manufacturer's Instructions for Use and soak for a minimum of 1 minute (water temperature 25 °C - 40 °C (77 °F - 104 °F)).
3. **Scrub:** While submerged, scrub the device with a soft nylon bristle general instrument brush for a minimum of 30 seconds (water temperature 25 °C - 40 °C (77 °F - 104 °F)). For hard-to-reach crevices on the sizers, use a 7 mm soft nylon bristle brush. For hard-to-reach crevices on the sterilization tray, use a 12 mm soft nylon bristle brush.
4. **Rinse:** Rinse the devices under running tap water for a minimum of 30 seconds (water temperature 25 °C - 40 °C (77 °F - 104 °F)).
5. **Inspection:** Visually inspect the device for any remaining soil. If needed, repeat the cleaning process from Step 2.
6. **Final rinse:** Rinse the devices under purified water (DI or RO) for a minimum of 30 seconds[4].
7. **Dry:** Dry the device with medical grade compressed air and a soft, lint-free cloth until completely dry.

Manual disinfection

Warning: Chemical disinfection may reduce the device's useful lifetime.

1. Chemical disinfectant: Submerge the device in a chemical[5] disinfection bath diluted as specified by the manufacturer's Instructions for Use.
2. Completely dissolve in lukewarm water by stirring. Solution is ready to use after 15 minutes.
3. Fully submerge the devices and soak in the working solution for 15 minutes. All surfaces must be completely wetted with the disinfectant solution. Air bubbles must be removed.
4. Using DI water, rinse off completely the chemical disinfection working solution for a minimum of 3 minutes[6].
5. Allow the devices to dry completely, using a soft, lint-free cloth if necessary.

6.4 Cleaning for countries outside the EU

Warning: Do not use solvent-based cleaners.

Note: Disinfection is not required because of terminal sterilization.

6.4.1 Automated cleaning

Automated cleaning pretreatment

1. Immerse the devices in a covered container with an enzymatic detergent (for example, ENZOL) using warm tap water 22 °C to 43 °C (72 °F to 110 °F).
2. Clean each component thoroughly using a soft bristle brush (for example, nylon toothbrush).
 Note: Take care to clean areas where soil may be impacted or shielded from the cleaning process.
3. Rinse the devices in warm running tap water 22 °C to 43 °C (72 °F to 110 °F) for a minimum of 2 minutes or until all traces of the cleaning solution are removed.
4. Visually inspect the devices for any remaining soil.
5. Repeat the above steps, if needed.
6. Allow the devices to dry.
7. Place the device into a perforated instrument tray suitable for use in the washer-disinfector. Ensure the device is placed in such a manner that will protect it from unwanted motion or potential mechanical damage during the automated cycles.
 Warning: Do not clean in the accessory tray.
 Place the perforated instrument tray with the device into the washer-disinfector.

Automatic washer

Medtronic used the Prolystica™ family of cleaners according to the manufacturer's recommendations to validate the automated cleaning process. It is the responsibility of the processor to ensure that the processing is conducted in accordance with a validated method. Any deviation by the processor from these recommendations should be evaluated. The cleaning method may be manual or automated per these Instructions for Use or an equivalent validated method.

1. Place the devices in an automated washer (for example, HAMO™ LS-1000 washer/disinfector or washer/disinfector meeting ISO 15883-1). Avoid contact between devices while loading the washer.

[2] neodisher SeptoClean™ or equivalent may be used in accordance with the detergent label instructions.

[3] ENZOL™ enzymatic or equivalent may be used in accordance with the detergent label instructions.

[4] Endotoxin levels less than 10 EU/mL and water microbial levels less than 10 CFU/mL.

[5] neodisher SeptoActive or equivalent may be used in accordance with the disinfectant label instructions. For neodisher SeptoActive, use 20 g/L (2.0 %). pH range is approximately 7.6 – 7.9.

[6] Endotoxin levels less than 10 EU/mL and water microbial levels less than 10 CFU/mL.

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Caution: Do not use instrument lubricant on the devices.

Treatment	Time (min)	Temperature	Water quality	Cleaning agent
Enzymatic wash	04:00	55 °C-60 °C (131 °F-140 °F)	tap water	Prolystica Ultra Concentrate Enzymatic Cleaner
Wash	02:00	45 °C-50 °C (113 °F-122 °F)	tap water	Prolystica Ultra Concentrate Neutral
Rinse	02:00	65 °C-70 °C (149 °F-158 °F)	tap water	not applicable
Dry	15:00	70 °C-80 °C (158 °F-176 °F)	not applicable	not applicable

2. Visually inspect each component for any remaining soil and moisture. If soil remains, repeat the process.
3. Dry the devices with filtered, compressed air or a lint-free wipe, if needed.

6.4.2 Manual cleaning

Medtronic has established the following manual cleaning procedure:

1. Submerge the devices in an enzymatic detergent (for example, ENZOL™) diluted with warm tap water 22 °C to 43 °C (72 °F to 110 °F) according to the manufacturer's recommendations for a minimum of 5 minutes.
2. Remove visible soil with a soft-bristled brush (for example, nylon toothbrush).
3. Rinse the devices with warm running tap water at 22 °C to 43 °C (72 °F to 110 °F) for a minimum of 1 minute.
4. Dry with a clean, lint-free wipe.
5. Visually inspect each component for any remaining soil or moisture. If soil remains, repeat the cleaning process. If wet, use a lint-free wipe to dry.

6.5 Sterilization and resterilization for all geographies

- A standard packaging material (for example, a sterilization wrap) meeting requirements in ISO 11607-1 may be used. Use standard packaging material to double wrap the instrument. Ensure that the pack is large enough to contain the instrument without stressing the seals. The validation performed by Medtronic utilized the appropriate accessory tray for each device.
- When sterilizing multiple devices in 1 autoclave cycle, ensure that the sterilizer's maximum load is not exceeded. Refer to the sterilizer equipment's instructions for use. Parameters below have been validated to provide a 10^-6 sterility assurance level (SAL) in accordance with requirements in ISO 17665-1. In accordance with this validation, users performing sterilization processes should also follow the recommended process standard ISO 17665-1.
- Examine the devices for cracking or degradation before use. Do not use if components show signs of cracking or degradation.
- Use a water source which meets the equipment manufacturer's specifications for steam quality, according to EN 285.

6.5.1 Steam sterilization

Table 3. Sterilization cycle parameters for the United States and its territories

Cycle type	Gravity displacement	Prevacuum (dynamic-air-removal)	
Temperature	132 °C (270 °F)	132 °C (270 °F)	135 °C (275 °F)
Exposure time	15 minutes	4 minutes	3 minutes
Dry time	15 minutes	30 minutes	16 minutes

For medical facilities located outside the United States and its territories: Some non-U.S. health care authorities recommend sterilization with parameters that will minimize the potential risk of transmitting Creutzfeldt-Jakob Disease (CJD). This recommendation is especially important for surgical instruments that could come into contact with the central nervous system.

Table 4. Sterilization cycle parameters for medical facilities outside the United States and its territories

Cycle type	Gravity displacement[c]	Prevacuum (dynamic-air-removal)	Prevacuum (dynamic-air-removal) for CJD[a]
Temperature	134 °C (273 °F)	134 °C (273 °F)	134 °C (273 °F)
Exposure time	20 minutes	4 minutes	18 minutes
Dry time[b]	30 minutes	30 minutes	30 minutes

[a] Medtronic recommends incinerating devices that have directly contacted patients suspected or confirmed with Transmissible Spongiform Encephalopathy (TSE)/CJD diagnosis. *NHS Estates HTM 2010 Parts 4 & 6: Appendix 2, Items contaminated with TSE Agents* and *WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies* refer to a TSE decontamination cycle using a steam autoclave at a

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temperature of 134 °C to 137 °C (273 °F to 279 °F) for a single cycle of 18 minutes or multiple cycles totaling 18 minutes (for example, six 3-minute cycles).

[b] The minimum dry times were validated using sterilizers with vacuum drying capabilities. Drying cycles using ambient atmospheric pressure may require longer dry times. Refer to the sterilizer manufacturer's recommendations.

[c] Do not use gravity displacement within the CE Market.

6.6 Reuse limitations

- The device is intended to be stored indefinitely. The material and construction of these devices is suited to handle multiple sterilization cycles without impacting safe use or functionality.
- End of life is determined by excessive wear or damage from normal use.
- See the EXAMINATION section of this document below to determine if the device is at the end of its useful life.

EXAMINATION – Devices must always be examined before and after cleaning and before use. Examination should be thorough and must include a visual inspection of working surfaces. Never use or process instruments or accessories that are damaged or nonfunctional.

Visual inspection – Ensure that the following are true:

- No cracks, deformation, corrosion, or distortion are present on any part of the device.
- All device markings are legible, including the device name and size.
- All parts are present and free of damage and deterioration.

If the device passes the visual inspection as described in these Instructions for Use, the instrument can be expected to perform as intended.

7 Disclaimer of warranty

The warnings contained in the product labeling provide more detailed information and are considered an integral part of this disclaimer of warranty. Although the product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental, or consequential damages caused by any use, defect, failure, or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort, or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this disclaimer of warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this disclaimer of warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this disclaimer of warranty did not contain the particular part or term held to be invalid.

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© 2023 Medtronic
Medtronic, Inc.
Medtronic B.V.