Medtronic

SimuPlus™ 7700FR, 7700FB

Flexible Annuloplasty Ring and Band • Гъвкави анулопластични пръстен и лента • Flexibilní anuloplastický prstenec a proužek • Fleksibel annuloplastikring og -bånd • Flexibler Anuloplastie-Ring / Flexibles Anuloplastie-Band • Εύκαμπτος δακτύλιος και εύκαμπτος ανοικτός δακτύλιος δακτυλιοπλαστικής • Anillo y banda flexibles para anuloplastia • Painduv annuloplastika rõngas ja võru • Taipuisa anuloplastiarengas ja -nauha • Anneau et bande pour annuloplastie flexibles • Fleksibilni prsten i fleksibilna traka za anuloplastiku • Rugalmas anuloplasztikai gyűrű és pánt • Anello ed emianello flessibili per anuloplastica • Lankstus anuloplastikos žiedas ir lankas • Elastīgais anuloplastijas gredzens un aploce • Флексибилен прстен и обрач за анулопластика • Flexibele annuloplastiekring en -band • Fleksibelt annuloplastikkring og -bånd • Elastyczny pierścień zamknięty i pierścień otwarty do annuloplastyki • Anel e banda de anuloplastia flexível • Anel e banda de anuloplastia flexíveis • Inel și bandă de anuloplastie flexibile • Гибкое кольцо и полукольцо для аннулопластики • Flexibilný anuloplastický prstenec a pásik • Upogljiv obroč in trak za anuloplastiko • Savitljivi prsten i traka za anuloplastiku • Flexibel anuloplastikring och flexibelt anuloplastikband • Esnek Anüloplasti Halkası ve Bandı • Гнучке кільце та півкільце для анулопластики

Instructions for Use • Инструкции за употреба • Návod k použití • Brugsanvisning • Gebrauchsanweisung • Οδηγίες χρήσης • Instrucciones de uso • Kasutusjuhend • Käyttöohjeet • Mode d'emploi • Upute za upotrebu • Használati útmutató • Istruzioni per l'uso • Naudojimo instrukcijos • Lietošanas pamācība • Упатство за употреба • Gebruiksaanwijzing • Bruksanvisning • Instrukcja użytkowania • Instruções de utilização • Instruções de utilização • Instrucțiuni de utilizare • Инструкция по эксплуатации • Pokyny na používanie • Navodila za uporabo • Uputstva za upotrebu • Bruksanvisning • Kullanım Talimatları • Інструкція з використання

! USA **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

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Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (“TM*”) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.

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English

1 Device description

The SimuPlus flexible annuloplasty ring Model 7700FR (*Figure 1*) and the SimuPlus flexible annuloplasty band Model 7700FB (*Figure 2*) consist of a flexible braided polyester fabric with a 3 mm cross-section. The ring or band may be implanted in the mitral annulus to stabilize the annulus following valve repair. The band may be implanted in the tricuspid annulus to stabilize the annulus following valve repair. The ring and band are marked at three points by colored sutures. Two markers correspond to the trigones (mitral valve) or septal commissures (tricuspid valve), and a central marker identifies the centerline of the ring or band. The individual ring or band size (26 mm to 42 mm in 2 mm increments) refers to the circumference between the green trigone markers on the ring or band. Silicone markers, impregnated with barium sulfate and tungsten, enable radiographic visualization. The ring and band are mounted to a holder to ease positioning and implantation.

The Medtronic SimuPlus annuloplasty ring and band holder incorporates a chordal reference line to assist the surgeon during chordal repair procedures. A laterally stretched suture is present across the holder opening above the typical area of leaflet coaptation. This suture helps approximate the level of the annular plane. If the surgeon decides to include artificial chords in the valve repair, the chordal reference line may be useful when determining the level at which to tie the artificial chord knots. The chordal reference line does not require any additional steps to remove. It is removed with the holder after the retention sutures are cut.

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2 Indications

The SimuPlus flexible annuloplasty ring and band are indicated for use in patients undergoing surgery for diseased or damaged mitral valves. The band is indicated for the same use for tricuspid valves. The SimuPlus flexible annuloplasty ring and band provide support for the mitral annulus or tricuspid annulus (band only) and restrict expansion of the annulus.

3 Contraindications

- severe, generalized, or localized bacterial endocarditis
- heavily calcified valves
- greatly dilated annulus (not reducible by standard techniques)
- severe valvular dysfunction (not correctable by standard techniques)
- valvular retraction with severely reduced mobility
- congenital malformations with lack of valvular tissue

4 Warnings

- For single use only. Do not resterilize the ring or band.
- Only surgeons who have received adequate training to determine whether incompetent, stenotic, or diseased heart valves are capable of being repaired or replaced should use this device.
- Only surgeons who have received appropriate training in valve repair, including ring and band implant and sizing techniques, should use this device.
- Correct annuloplasty ring and band sizing is an important element of a successful valve repair. Undersizing the ring or band can result in valve stenosis. Oversizing the ring or band can result in valve regurgitation.
- During implantation, place horizontal, interrupted or continuous sutures into the annulus following the line of leaflet insertion.
- In the tricuspid position, do not place sutures in the coronary sinus, right coronary artery, AV node or other conduction tissue, as this may result in impairment of the cardiac conduction system.
- It is necessary to secure the ends of the ring and band into fibrous trigone tissue to maintain permanent attachment of the ring or band.
- Only surgeons with appropriate training and experience should include artificial chords in a heart valve repair procedure. The level at which to tie the chords should be determined by the surgeon's judgment, based on the surgeon's training, level of experience, and the characteristics of the patient and lesion.
- Correct placement and anchoring of artificial chords are important elements of successful valve repair. Improper placement and

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anchoring of artificial chords can result in leaflet tethering or regurgitation.

- Determining the correct length of artificial chords is an important element of successful valve repair. Artificial chords that are too long can result in mitral valve prolapse. Artificial chords that are too short can result in restricted leaflet motion.
- Suture knots must be securely tied. Loose knots and long suture tails may be a source for hemolysis, thrombosis, or thromboembolism.
- Do not cut the ring or band, as resultant loose threads can be a source of hemolysis, thrombosis, or thromboembolism.
- Do not cut the chordal reference line before removing the holder from the ring or band. Loose threads may fall into the patient anatomy and cause embolism if not retrieved, or tissue damage may occur during retrieval.
- Intraoperative or postoperative echocardiography should be used to evaluate the effectiveness of the valve repair. Minimizing regurgitation and preventing systolic anterior motion (SAM) are important elements of an effective repair.
- Surgeons who use annuloplasty rings or bands should be current on all anticoagulation regimens.
- When postoperative anticoagulant therapy is used, the patient's anticoagulation status should be carefully monitored.
- Patients with intra-atrial thrombi or a giant left atrium may benefit from long-term anticoagulation therapy.
- The surgeon may desire that patients in atrial fibrillation remain on anticoagulation therapy until sinus rhythm is established.

5 Precautions

- Do not use cutting edge needles, as they may damage the annuloplasty device potentially leading to ring or band dehiscence, ring or band fracture, and possible mitral or tricuspid regurgitation.
- Take care not to damage the annuloplasty ring or band during handling.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.

6 Potential adverse events

While infrequent, certain complications have been reported when using annuloplasty rings or bands. These include the following:

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- uncorrected or recurrent regurgitation
- stenosis
- ring or band dehiscence
- hemolysis (even with mild regurgitation)
- low cardiac output
- heart block
- systolic anterior motion (SAM) and left ventricular outflow tract obstruction (LVOTO)
- damage to coronary arteries
- endocarditis
- thrombosis
- thromboembolism
- anticoagulant-related hemorrhage
- ring or band fracture
- leaflet perforation
- bleeding diathesis

The potential for these complications should be considered when selecting the most beneficial surgical procedure for each patient.

To avoid or minimize occurrence of these adverse events, the annuloplasty repair, including sizing and implantation, should be conducted in accordance with the methods described in these instructions for use by surgeons with appropriate training and experience in valve repair.

7 Individualization of treatment

To allow for healing and incorporation of the annuloplasty ring or band by host tissue, regardless of cardiac rhythm, postoperative anticoagulation therapy should be considered for at least six weeks following surgery.

8 Patient counseling information

8.1 Patient dental care

Patients with annuloplasty rings or bands who undergo dental or other potentially bacteremic procedures must be considered for prophylactic antibiotic therapy.

8.2 MRI conditions

The SimuPlus flexible annuloplasty ring or band is MR Safe. This device poses no known hazards in any magnetic resonance (MR) environment.

The device will not cause any harm to the patient when exposed to MR scanning immediately after implantation. MRI at 1.5 Tesla and 3.0 Tesla and

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other field strengths may be performed immediately following the implantation of the device.

The presence of other implants or medical circumstances of the patient may require lower limits on some or all of the above parameters.

9 How supplied

9.1 Packaging

The SimuPlus flexible annuloplasty ring and band Model 7700FR and Model 7700FB are available in the following sizes: 26 mm, 28 mm, 30 mm, 32 mm, 34 mm, 36 mm, 38 mm, 40 mm, and 42 mm. The products are nonpyrogenic. They are sterilized using irradiation and supplied in double-aseptic transfer pouches contained in an outer carton. The packaging system is designed to ease placement of the device into the sterile field. The ring and band assemblies are sterile if the pouches are undamaged and unopened. The outer surfaces of the outer pouch are nonsterile and must not be placed in the sterile field.

Warning: Carefully inspect all product packaging for damage or defects before use. Do not use the product if the unit container has been opened or damaged, if any device damage is observed, or if the integrity of the sterile barrier has been compromised. Contact your Medtronic representative about how to return the device.

9.2 Storage

Store the product in its original packaging, including the outer shelf carton, in a clean, cool, and dry area to protect the product and minimize the potential for contamination.

The sterility and nonpyrogenicity of the SimuPlus flexible annuloplasty ring and band are validated to remain unaffected until the use-by date identified on the shelf carton, provided the pouches are not opened or damaged.

10 Instructions for use

10.1 Sizing

The SimuPlus sizers are used to select the proper SimuPlus flexible annuloplasty ring or band size for repair of the mitral or tricuspid valve. Sizers are reusable but must be cleaned and sterilized by autoclave (steam) before each use.

Proper ring and band size selection is an important part of valvular annuloplasty to help restore proper function. For size selection, use a SimuPlus polysulfone sizer Model 7700PS or silicone sizer Model 7700FS. The polysulfone sizers interface with the Medtronic annuloplasty handle Model 7686, provided separately (*Figure 3*).

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Handles and sizers are provided non-sterile and must be cleaned and sterilized before use.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the ring or band. Other ring and band sizers may not indicate the appropriate SimuPlus flexible annuloplasty ring and band size.

10.2 Mitral annulus sizing

Determine the proper SimuPlus flexible annuloplasty ring or band size to implant.

1. Lower the sizer onto the valve annulus and align the trigone markers on the sizer with the trigones (*Figure 4*).
2. Determine the appropriate mitral ring or band size by locating the sizer with a trigone marker spacing that best corresponds to the intertrigonal distance.

Note: The outer periphery of the sizer corresponds to the centerline of the ring or band and indicates where the implantation sutures will be placed on the valve annulus.

Other considerations for sizing include the patient disease state and physician surgical technique.

10.3 Tricuspid annulus sizing

Determine the proper SimuPlus flexible annuloplasty band size to implant.

1. Lower the sizer onto the valve annulus and align the commissure markers on the sizer with the anteroseptal and posteroseptal commissures of the septal leaflet (*Figure 5*).
2. Determine the appropriate tricuspid band size by locating the sizer with a commissure marker spacing that best corresponds to the circumferential distance between the commissures.

Note: The outer periphery of the sizer corresponds to the centerline of the band and indicates where the implantation sutures will be placed on the valve annulus.

Other considerations for sizing include the patient disease state and physician surgical technique.

10.4 Handling and preparation instructions

1. Open the box and remove the product literature, Patient Registration Form, and envelopes.
2. Remove the double-aseptic transfer pouch containing the device.

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3. Inspect the pouches ensuring they have not been opened or damaged. The ring and band assemblies are sterile as long as the inner pouch has not been compromised. If the inner pouch is damaged, do not implant the device.
4. If the outer pouch is damaged, the exterior surface of the inner pouch may not be sterile.
5. Open the outer transfer pouch and, while still holding the bottom of the outer pouch, pass the inner pouch into the sterile field.
6. The inner pouch should be opened only in the sterile field.

10.5 Device implantation

1. Remove the appropriate size SimuPlus flexible annuloplasty ring or band from the sterile package using aseptic technique.
2. Remove the serial number identification tag by cutting the retaining suture and record the serial number in the patient's record. Verify that the serial number matches the serial number on the Patient Registration Form.

 Warning: The serial tag must be removed from the ring or band for proper function. Do not cut or tear the ring or band fabric during removal of the serial tag.
3. Snap the handle into the handle flange on the holder.

10.6 Mitral annulus suture placement

Ring: Place sutures in each trigone, approximately 4 mm apart. Place additional interrupted sutures, approximately 4 mm apart, in the anterior and posterior portion of the mitral annulus (*Figure 6*, number 1).

Band: Place sutures in each trigone, approximately 4 mm apart. Place additional interrupted sutures, approximately 4 mm apart, in the posterior portion of the mitral annulus (*Figure 6*, number 2).

Warning: Avoid placing sutures in the circumflex coronary artery.

Note: Pledgets may be used for trigone sutures and mid-posterior annulus suture to reduce the possibility of ring or band dehiscence.

Note: If performing chordal repair, the artificial chords should be placed using conventional surgical repair techniques before placing annular sutures through the annuloplasty ring.

10.7 Tricuspid suture placement

Place sutures in the septal commissures, approximately 4 mm apart. Place additional interrupted sutures, approximately 4 mm apart, in the septal, anterior, and posterior portion of the tricuspid annulus (*Figure 6*, number 3).

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Warning: Do not place sutures in the coronary sinus, right coronary artery, AV node, bundle of HIS, or other conduction tissue.

10.8 Ring and band suture placement

1. Orient the ring or band holder assembly on the annulus.
2. Pass sutures through the ring or band, approximately 2 mm to 4 mm apart, entering at the bottom of the ring or band and exiting the periphery of the ring or band (*Figure 7*).
3. Place sutures completely around the ring or band.

10.9 Holder and handle removal

10.9.1 Holder and handle removal without chordal repair

1. Use the handle to push the ring or band down onto the valve annulus while pulling back on the sutures.
2. Use a scalpel to cut the 1 retention suture on the flange (*Figure 8*).

 Warning: Do not cut the ring or band fabric while cutting the flange retention suture.
3. Remove the flange and handle and dispose of the flange.
4. Securely tie all knots around the ring or band and trim all excess suture tails (*Figure 9*).
5. Use a scalpel to cut the 2 retention sutures on the holder in the areas indicated (*Figure 10* and *Figure 11*).

 Warning: Do not cut the ring or band fabric while cutting the holder retention suture.

 Warning: Do not cut the chordal reference line before removing the holder from the ring or band. Loose threads may fall into the patient anatomy and cause embolism if not retrieved, or tissue damage may occur during retrieval.
6. Remove the holder from the annulus.

 Warning: For proper valve function, the holder must be removed from the ring or band at the end of the procedure. Do not leave the holder attached to the ring or band under any circumstances.

 Note: The holder posts can be grasped with forceps or other instrument for ease of removal (*Figure 12*).
7. Dispose of the holder.
8. Test the completed valve repair for competency (*Figure 13*).

10.9.2 Holder and handle removal with chordal repair

1. Pull the artificial chords through the ring or band holder.
2. Use the handle to push the ring or band down onto the mitral valve annulus while pulling back on the sutures.

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3. Use a scalpel to cut the 1 retention suture on the handle flange (*Figure 8*).

 Warning: Do not cut the ring or band fabric while cutting the flange retention suture.

4. Remove the flange and handle and dispose of the flange.
5. Securely tie all knots around the ring or band and trim all excess suture tails (*Figure 9*).
6. Tie the artificial chordal knots.

 Note: The surgeon should determine the appropriate artificial chord length based on factors such as patient anatomy, disease state, leaflets being repaired, and the surgeon's experience and training in surgical valve repair. The chordal reference line is a reference feature that can be used to compare the artificial chord length to the level of the annular plane.

 Warning: Ensure correct placement and anchoring of artificial chords. Improper placement and anchoring of artificial chords can result in leaflet tethering or regurgitation.

 Warning: Determine the correct length of the artificial chords. Artificial chords that are too long can result in mitral valve prolapse or regurgitation. Artificial chords that are too short may result in restricted leaflet motion or regurgitation.

 Warning: When tying the chordal knots, avoid displacing the papillary muscles, which may result in artificial chords that are too short.

7. Use a scalpel to cut the 2 retention sutures on the holder in the areas indicated (*Figure 10* and *Figure 11*).

 Warning: Do not cut the ring or band fabric while cutting the holder retention suture.

 Warning: Do not cut the chordal reference line before removing the holder from the ring or band. Loose threads may fall into the patient anatomy and cause embolism if not retrieved, or tissue damage may occur during retrieval.

8. Remove the holder from the annulus.

 Warning: For proper valve function, the holder must be removed from the ring or band at the end of the procedure. Do not leave the holder attached to the annuloplasty ring or band under any circumstances.

 Note: The holder posts can be grasped with a forceps or other instrument for ease of removal (*Figure 12*).

9. Dispose of the holder.
10. Test the completed valve repair for competency (*Figure 13*).

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10.10 Accessories

Use one of the SimuPlus polysulfone or silicone sizers to determine the appropriate device size.

Use only handle Model 7686 to interface with the holder and polysulfone sizers.

For detailed information on accessory use and sterilization, refer to the instructions for use for SimuPlus accessories (sizers and handle).

Warning: Do not use annuloplasty sizers or handles from other manufacturers or other Medtronic annuloplasty products to size the SimuPlus flexible annuloplasty ring or band.

10.11 Sterilization

The ring and band are provided sterile (irradiation) on the holder and must not be resterilized. Rings and bands that have been damaged or contaminated from patient contact should not be used.

11 Registration Information

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information.

A patient registration form is included in each device package.

After implantation, please complete all requested information.

The serial number may be found on the package and on the identification tag attached to the device. Return the original form to the Medtronic address indicated on the form and provide the temporary identification card to the patient prior to discharge.

12 Disclaimer of Warranty

The following disclaimer of warranty applies to United States customers only:

ALTHOUGH THE SIMUPLUS FLEXIBLE ANNULOPLASTY RING AND BAND, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY

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ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

13 Disclaimer of Warranty

The following disclaimer of warranty applies to customers outside the United States:

ALTHOUGH THE SIMUPLUS FLEXIBLE ANNULOPLASTY RING AND BAND, HEREAFTER REFERRED TO AS "PRODUCT," HAVE BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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English

1 Device description

The SimuForm semi-rigid annuloplasty ring Model 7800RR (*Figure 1*) consists of a MP35N wire stiffener in the posterior segment, running from trigone to trigone. The ring stiffener is enclosed within a close-coiled MP35N spring that passes around the circumference of the annuloplasty ring. The spring is covered by a thin silicone sheath. Braided polyester fabric is used to cover and form the body of the ring. The ring has two green markers to indicate the anterior and posterior trigones. A green demarcation suture runs around the upper face of the ring. The individual ring size (24 mm to 40 mm in 2 mm increments) refers to the inner circumference between the green trigone markers on the ring. The internal

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spring and stiffener provide radiographic visualization around the circumference of the ring.

The Medtronic SimuForm annuloplasty ring holder incorporates a chordal reference line to assist the surgeon during chordal repair procedures. A laterally stretched suture is present across the ring holder opening above the approximate area of leaflet coaptation. This suture helps approximate the level of the annular plane. If the surgeon decides to include artificial chords in the valve repair, the chordal reference line may be useful when determining the level at which to tie the artificial chord knots. The chordal reference line does not require any additional steps to remove. It is removed with the holder after the retention sutures are cut.

2 Indications

The SimuForm semi-rigid annuloplasty rings are for use in patients undergoing surgery for diseased or damaged mitral valves. The SimuForm semi-rigid annuloplasty ring provides support for the mitral annulus and restricts expansion of the annulus.

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4 Warnings

- For single use only. Do not resterilize the ring.
- Only surgeons who have received adequate training to determine whether incompetent, stenotic, or diseased heart valves are capable of being repaired or replaced should use this device.
- Only surgeons who have received appropriate training in valve repair, including ring implant and sizing techniques, should use this device.
- Correct annuloplasty ring sizing is an important element of a successful valve repair. Undersizing the ring can result in valve stenosis. Oversizing the ring can result in valve regurgitation.
- Care should be taken that sutures are not placed in atrial tissue, as this may result in impairment of the cardiac conduction system.
- It is necessary to secure the ends of the ring into fibrous trigone tissue to maintain permanent attachment of the ring.

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- Only surgeons with appropriate training and experience should include artificial chords in a heart valve repair procedure. Determining the suitable level at which to tie the chords is the surgeon's judgement, based on the surgeon's training, level of experience, and the characteristics of the patient and lesion.
- Correct placement and anchoring of artificial chords are important elements of successful valve repair. Improper placement and anchoring of artificial chords can result in leaflet tethering or regurgitation.
- Determining the correct length of artificial chords is an important element of successful valve repair. Artificial chords that are too long can result in mitral valve prolapse. Artificial chords that are too short can result in restricted leaflet motion.
- Suture knots must be securely tied. Loose knots and long suture tails may be a source for hemolysis, thrombosis, or thromboembolism.
- Do not cut the ring, as resultant loose threads can be a source of hemolysis, thrombosis, or thromboembolism.
- Do not cut the chordal reference line before removing the holder from the ring. Loose threads may fall into the patient anatomy and cause embolism if not retrieved, or tissue damage may occur during retrieval.
- Do not alter or deform the ring to conform to annular anatomy as this could lead to ring fracture and possible mitral regurgitation.
- Intraoperative or postoperative echocardiography should be used to evaluate the effectiveness of the valve repair. Minimizing regurgitation and preventing systolic anterior motion (SAM) are important elements of an effective repair.
- Surgeons who use annuloplasty rings should be current on all anticoagulation regimens.
- When postoperative anticoagulant therapy is used, the patient's anticoagulation status should be carefully monitored.
- Patients with intra-atrial thrombi or a giant left atrium may benefit from long-term anticoagulation therapy.
- The surgeon may desire that patients in atrial fibrillation remain on anticoagulation therapy until sinus rhythm is established.
- Do not squeeze the ring with sharp instruments as this may damage the surface of the stiffener, which may result in ring fracture and possible mitral regurgitation.
- The ring is indicated for use in the mitral position only. Use in the tricuspid position may result in ring fracture, conduction system damage, tricuspid regurgitation, or ring dehiscence.

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5 Precautions

- Do not use cutting edge needles, as they may damage the annuloplasty device potentially leading to ring dehiscence, ring fracture, and possible mitral regurgitation.
- Take care not to damage the annuloplasty ring during handling.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.

6 Potential adverse events

While infrequent, certain complications have been reported when using annuloplasty rings. These include the following:

- uncorrected or recurrent regurgitation
- stenosis
- ring dehiscence
- hemolysis (even with mild regurgitation)
- low cardiac output
- heart block
- systolic anterior motion (SAM) and left ventricular outflow tract obstruction (LVOTO)
- damage to coronary arteries
- endocarditis
- thrombosis
- thromboembolism
- anticoagulant-related hemorrhage
- ring fracture
- leaflet perforation
- bleeding diathesis

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7 Individualization of treatment

To allow for healing and incorporation of the annuloplasty ring by host tissue, regardless of cardiac rhythm, postoperative anticoagulation therapy should be considered for at least six weeks following surgery.

8 Patient counseling information

8.1 Patient dental care

Patients with annuloplasty rings who undergo dental or other potentially bacteremic procedures must be considered for prophylactic antibiotic therapy.

8.2 MRI conditions

Nonclinical testing has demonstrated the SimuForm semi-rigid ring is MR Conditional. A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T or 3 T
- Maximum spatial field gradient of 2500 gauss/cm (25 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Operating Mode).

Under the scan conditions defined above, the SimuForm semi-rigid ring is expected to produce a maximum temperature rise of less than 2.2°C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the device extends approximately 20 mm from the SimuForm semi-rigid ring when imaged with a gradient echo pulse sequence and 3 T MRI system.

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9 How supplied

9.1 Packaging

The SimuForm semi-rigid annuloplasty ring Model 7800RR is available in the following sizes: 24 mm, 26 mm, 28 mm, 30 mm, 32 mm, 34 mm, 36 mm, 38 mm, and 40 mm. The products are nonpyrogenic. They are sterilized using irradiation and supplied in double-aseptic transfer pouches contained in an outer carton. The packaging system is designed to ease placement of the device into the sterile field. The ring assemblies are sterile if the pouches are undamaged and unopened. The outer surfaces of the outer pouch are nonsterile and must not be placed in the sterile field.

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The sterility and nonpyrogenicity of the SimuForm semi-rigid annuloplasty ring are validated to remain unaffected until the use-by date identified on the shelf carton, provided the pouches are not opened or damaged.

10 Instructions for use

10.1 Sizing

The SimuForm polysulfone sizers are used to select the proper SimuForm semi-rigid annuloplasty ring size for repair of the mitral valve. Sizers are reusable but must be cleaned and sterilized by autoclave (steam) before each use.

Proper ring size selection is an important part of valvular annuloplasty to help restore proper function. Use a SimuForm polysulfone sizer set (Model 7800PS) and annuloplasty handle (Model 7686) for size selection (*Figure 2*).

Handles and sizers are provided nonsterile and must be cleaned and sterilized before use.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the ring. Other ring sizers may not indicate the appropriate SimuForm semi-rigid annuloplasty ring size.

10.2 Mitral annulus sizing

Determine the proper SimuForm semi-rigid annuloplasty ring size to implant.

1. Lower the sizer onto the valve annulus and align the trigone markers on the sizer with the trigones (*Figure 3*).
2. Determine the appropriate mitral ring size by locating the sizer with a trigone marker spacing that best corresponds to the intertrigonal distance.

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Note: The outer periphery of the sizer corresponds to the outer diameter of the ring stiffener and indicates where the implantation sutures will be placed on the valve annulus.

Other considerations for sizing include the patient disease state and the physician surgical technique.

10.3 Handling and preparation instructions

1. Open the box and remove the product literature, Patient Registration Form, and envelopes.
2. Remove the double-aseptic transfer pouch containing the device.
3. Inspect the pouches ensuring they have not been opened or damaged. The ring assembly is sterile as long as the inner pouch has not been compromised. If the inner pouch is damaged, do not implant the device.
4. If the outer pouch is damaged, the exterior surface of the inner pouch may not be sterile.
5. Open the outer transfer pouch and, while still holding the bottom of the outer pouch, pass the inner pouch into the sterile field.
6. The inner pouch should be opened only in the sterile field.

10.4 Device implantation

The SimuForm semi-rigid annuloplasty ring assembly may be used with or without the handle (annuloplasty handle, Model 7686, provided separately).

Caution: The SimuForm semi-rigid annuloplasty ring is designed for implantation only with interrupted-suture techniques.

Warning: The handle must be cleaned and sterilized before use.

1. Remove the appropriate sized SimuForm semi-rigid annuloplasty ring from the sterile package using aseptic technique.
2. Remove the serial number identification tag by cutting the retaining suture and record the serial number in the patient's record. Verify that the serial number matches the serial number on the Patient Registration Form.

 Warning: The serial tag must be removed from the ring for proper function. Do not cut or tear the ring fabric during removal of the serial tag.
3. Align and engage a sterile handle into the handle flange of the holder. The thin section of the handle is malleable to permit alignment of the ring holder with the annulus.

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10.5 Mitral annulus suture placement

The implantation sutures should be placed outside the green demarcation marker to avoid the internal stiffener.

Place sutures in each trigone, approximately 4 mm apart. Place additional interrupted sutures, approximately 4 mm apart, in the anterior and posterior portion of the mitral annulus (*Figure 4*).

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Note: Pledgets may be used for trigone sutures and mid-posterior annulus suture to reduce the possibility of ring dehiscence.

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10.6 Ring suture placement

1. Orient the ring holder assembly on the annulus.
2. Pass the sutures through the ring, approximately 2 mm to 4 mm apart, entering at the bottom of the ring and exiting at the periphery of the ring (*Figure 5*).
3. Place the sutures completely around the ring.

10.7 Holder and handle removal

10.7.1 Holder and handle removal without chordal repair

1. Use the handle to push the ring down onto the valve annulus while pulling back on the sutures.
2. Use a scalpel to cut the 1 retention suture on the flange (*Figure 6*).

 Warning: Do not cut the ring fabric while cutting the flange retention suture.
3. Remove the flange and handle and dispose of the flange.
4. Securely tie all knots around the ring and trim all excess suture tails (*Figure 7*).
5. Use a scalpel to cut the 2 retention sutures on the holder in the areas indicated (*Figure 8*).

 Warning: Do not cut the ring fabric while cutting the holder retention sutures.

 Warning: Do not cut the chordal reference line before removing the holder from the ring. Loose threads may fall into the patient anatomy and cause embolism if not retrieved, or tissue damage may occur during retrieval.
6. Remove the holder from the annulus.

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Warning: For proper valve function, the holder must be removed from the ring at the end of the procedure. Do not leave the holder attached to the annuloplasty ring under any circumstances.

Note: The holder posts can be grasped with a forceps or other instrument for ease of removal (*Figure 9*).

7. Dispose of the holder.
8. Test the completed valve repair for competency (*Figure 10*).

10.7.2 Holder and handle removal with chordal repair

1. Pull the artificial chords through the ring holder.
2. Use the handle to push the ring down onto the mitral valve annulus while pulling back on the sutures.
3. Use a scalpel to cut the 1 retention suture on the flange (*Figure 6*).

 Warning: Do not cut the ring fabric while cutting the flange retention suture.
4. Remove the flange and handle and dispose of the flange.
5. Securely tie all knots around the ring and trim all excess suture tails (*Figure 7*).
6. Tie the artificial chordal knots.

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 Warning: When tying the chordal knots, avoid displacing the papillary muscles, which may result in artificial chords that are too short.
7. Use a scalpel to cut the 2 retention sutures on the holder in the areas indicated (*Figure 8*).

 Warning: Do not cut the ring fabric while cutting the holder retention sutures.

 Warning: Do not cut the chordal reference line before removing the holder from the ring. Loose threads may fall into the patient anatomy and cause embolism if not retrieved, or tissue damage may occur during retrieval.

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8. Remove the holder from the annulus.

 Warning: For proper valve function, the holder must be removed from the ring at the end of the procedure. Do not leave the holder attached to the annuloplasty ring under any circumstances.

 Note: The holder posts can be grasped with a forceps or other instrument for ease of removal (*Figure 9*).
9. Dispose of the holder.
10. Test the completed valve repair for competency (*Figure 10*).

10.8 Accessories

Use one of the SimuForm polysulfone sizers to determine the appropriate device size.

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For detailed information on accessory use and sterilization, refer to the instructions for use for SimuForm accessories (sizers and handle).

Warning: Do not use annuloplasty sizers or handles from other manufacturers or other Medtronic annuloplasty products to size the SimuForm semi-rigid annuloplasty ring.

10.9 Sterilization

The ring is provided sterile (irradiation) on the holder and must not be resterilized. Rings that have been damaged or contaminated from patient contact should not be used.

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ALTHOUGH THE SIMUFORM SEMI-RIGID ANNULOPLASTY RING, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

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ALTHOUGH THE SIMUFORM SEMI-RIGID ANNULOPLASTY RING, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable

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law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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Tricuspid Band

1. Device description

The Tri-Ad™ 2.0 Adams tricuspid band, Model 900SFC, consists of an MP35N polished and formed wire stiffener fitted with 2 end caps. The band is covered with a braided polyester fabric. The stiffener runs from the posteroseptal commissure to the anterior segment distal to the anteroposterior commissure. The stiffener helps remodel and stabilize the enlarged portion of the tricuspid annulus (Figure 1). The band has 3 radial green markers to locate the following 3 areas of the band:

- The end of the stiffener section at the posteroseptal commissure.
- The anteroposterior commissure.
- The end of the stiffener section in the anterior segment of the band.

A green demarcation suture runs around the middle of the band in the stiffener region. The band size refers to the dimension of the band where the implantation sutures are placed, except in the stiffener region. The internal stiffener, and the radiopaque core in both fully flexible ends of the band, provide radiographic visualization along the entire circumference of the band. The band is available in the following sizes: 26 mm, 28 mm, 30 mm, 32 mm, 34 mm, and 36 mm.

2. Indications for use

The Tri-Ad 2.0 Adams tricuspid band is indicated for the reconstruction or remodeling of pathological tricuspid valves. The band provides support and restricts expansion of the tricuspid annulus.

3. Contraindications

- Severe, generalized, or localized bacterial endocarditis
- Heavily calcified valves
- Greatly dilated annulus (not reducible by standard techniques)
- Severe valvular dysfunction (not correctable by standard techniques)
- Valvular retraction with severely reduced mobility
- Congenital malformations with lack of valvular tissue

4. Warnings and precautions

4.1. Warnings

- Do not use this device without adequate training in band implant and sizing techniques.
- Determine whether an incompetent, stenotic, or diseased heart valve can be repaired or replaced.
- Correct annuloplasty band sizing is an important element of a successful valve repair. Undersizing the band can result in valve stenosis, band dehiscence, or band fracture. Oversizing the band can result in valve regurgitation or band fracture.
- Do not place sutures in the coronary sinus, right coronary artery, AV node, or other conduction tissue. Sutures in these locations can result in impairment of the cardiac conduction system.
- Suture knots must be securely tied. Loose knots and long suture tails can be a source for hemolysis, thrombosis, or thromboembolism.
- Do not cut the band. Cutting the band can result in loose threads that can be a source of hemolysis, thrombosis, or thromboembolism. In addition, resultant jagged edges of the stiffener can lead to tissue trauma.
- Do not alter or deform the band to conform to annular anatomy. Altering or deforming the band can result in band fracture, valve regurgitation, or stenosis.
- To avoid damage to the surface of the stiffener, do not squeeze the band with sharp instruments. Squeezing the band with sharp instruments can result in band fracture or tricuspid regurgitation.
- Intraoperative or postoperative echocardiography can be used to evaluate the effectiveness of the valve repair. Minimizing regurgitation is an important element of an effective repair.
- Ensure that current anticoagulation regimens are used.
- Carefully monitor the patient's anticoagulation status when postoperative anticoagulant therapy is used.
- Patients with intra-atrial thrombi or a giant left atrium could benefit from long-term anticoagulation therapy.
- Patients in atrial fibrillation can remain on anticoagulation therapy until sinus rhythm is established.
- The band is indicated for use in the tricuspid position only. Use in the mitral position could result in band fracture, mitral regurgitation, or band dehiscence.

4.2. Precautions

- Do not use cutting edge needles. Damage to the annuloplasty band can result and lead to band dehiscence, band fracture, or tricuspid regurgitation.
- Take care not to damage the annuloplasty band during handling.

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- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device, which could result in patient injury, illness, or death.

5. Potential adverse events

While infrequent, certain complications have been reported when using annuloplasty bands. Reported complications include the following adverse events:

- Uncorrected or recurrent regurgitation
- Stenosis
- Band dehiscence
- Hemolysis (even with mild regurgitation)
- Low cardiac output
- Heart block
- Damage to coronary arteries
- Endocarditis
- Thrombosis
- Thromboembolism
- Anticoagulant-related hemorrhage
- Band fracture
- Leaflet perforation
- Bleeding diathesis

Consider the potential for these complications when selecting the most beneficial surgical procedure for each patient.

Use the sizing and implantation methods described in these Instructions for Use to avoid or minimize occurrence of these adverse events. Ensure that the physician has appropriate training and experience in valve repair.

6. Instructions for use

6.1. Sizing

Select the appropriate band size to help restore proper function of the valve. Use the Tri-Ad™ 2.0 Adams sizer set, Model 7900PS, and the Medtronic annuloplasty handle, Model 7686, to select the proper band size for repair of the tricuspid valve. The reusable sizers and handle are provided nonsterile and must be thoroughly cleaned and sterilized before use. Refer to the appropriate instructions for use for further information.

Align and insert the handle into the snap cavity of the sizer. The thin section of the handle is malleable to align the sizers with the valve annulus.

Warning: Do not use other manufacturers' annuloplasty sizers or sizers from other Medtronic annuloplasty products to size the Tri-Ad 2.0 Adams tricuspid band. Other sizers may not indicate the appropriate band size.

Place an initial suture near the posteroseptal commissure approximately at the level of the coronary sinus (Figure 2).

To determine the appropriate size band, match the outer portion or periphery of the sizer to the location of the implantation suture in the annulus.

1. Lower the sizer onto the valve annulus and align the sizer notches with the commissures of the tricuspid valve.
2. Gently extend the anterior leaflet and cover its surface with the selected sizer (Figure 3).

The sizer that conforms most closely to the remodeled annulus, indicates the appropriate tricuspid band size to be considered for implant. The notches on the sizer represent the 3 commissures of the tricuspid valve. The sizers have a small hole that can be used to help with sizer removal by tying a length of 2-0 (or larger) suture to the sizer through this hole.

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6.2. Handling and preparation instructions

1. Open the box. Remove the product literature and Patient Registration Form.
2. Remove the double-aseptic transfer pouch containing the band mounted to a holder.
3. Inspect the outer pouch. If the outer pouch is damaged, the exterior surface of the inner pouch may not be sterile.
4. Open the outer transfer pouch, and while still holding the bottom of the outer pouch, pass the inner pouch into the sterile field.
5. Inspect the inner pouch. The band assembly is sterile as long as the inner pouch has not been compromised. If the inner pouch is damaged, do not implant the band.
6. Open the inner pouch only in the sterile field.

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6.3. Device implantation

The band assembly may be used with or without the Medtronic annuloplasty handle, Model 7686, provided separately.

1. Thoroughly clean and sterilize the handle prior to each use.
2. With the printed side of the holder facing up, align and engage the sterile handle into the snap fit cavity of the holder. The handle must be inserted into the side of the holder with the snap fit cavity. The thin section of the handle is malleable to align the device with the annulus.
3. Remove the serial number identification tag. Record the serial number in the patient's record. Verify that the serial number matches the serial number on the Patient Registration Form.

Warning: The serial number tag must be removed from the band for proper function. Do not cut or tear the fabric of the band during removal of the serial number tag.

6.4. Tricuspid annulus suture placement

Place interrupted sutures, approximately 4 mm apart, along the annulus ending at the midpoint of the anterior leaflet (Figure 4). Place approximately 8 to 10 sutures in the valve annulus.

Warning: Do not place sutures in the following locations:

- Coronary sinus
- Right coronary artery
- Atrioventricular (AV) node
- Bundle of His, or other conduction tissue.

Note: Pledgets can be used with annular sutures to reduce the possibility of band dehiscence.

6.5. Band suture placement

1. Orient the band assembly on the tricuspid annulus.
2. Entering at the bottom of the band, pass the sutures through the band approximately 2 mm to 4 mm apart. Exit on the top of the band. To avoid the internal stiffener, place the implantation sutures outside the green demarcation marker.

6.6. Holder and handle removal

1. Push the band down onto the valve annulus while pulling back on the sutures.
2. Use a scalpel to cut the 2 retention sutures on the holder in the areas indicated (Figure 5). Dispose of the holder.

 Warning: Do not cut the band fabric while cutting the holder retention suture.

 Warning: For proper valve function, the holder must be removed from the band at the end of the procedure. Do not leave the holder attached to the annuloplasty band under any circumstances.
3. Securely tie all knots around the band and trim all excess sutures (Figure 6).
4. Test valvular competency.

7. Individualization of treatment

Consider postoperative anticoagulation therapy for at least 6 weeks to allow for healing and incorporation of the annuloplasty band by host tissue, regardless of cardiac rhythm.

8. Patient counseling information

Patients with annuloplasty bands are at increased risk for endocarditis due to bacteremia (for example, when undergoing dental procedures). Advise patients about prophylactic antibiotic therapy.

Encourage patients to carry the Implanted Device Identification Card, provided by Medtronic, with them at all times.

9. How supplied

9.1. Packaging

A single annuloplasty band assembly (device and holder) is packaged in sterile, double-aseptic transfer pouches. The packaging system is designed to ease placement of the assembly into the sterile field. The assembly is sterile if the pouches are undamaged and unopened. The outer surfaces of the outer pouch are nonsterile and must not be placed in the sterile field.

9.2. Available sizes

The Tri-Ad 2.0 Adams tricuspid band, Model 900SFC, is available in the following sizes for the tricuspid position: 26 mm, 28 mm, 30 mm, 32 mm, 34 mm, and 36 mm.

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The sterility and nonpyrogenicity of the Tri-Ad 2.0 Adams tricuspid band is validated to remain unaffected until the Use-by date identified on the shelf carton, provided the pouches are not opened or damaged.

10. Patient information

10.1. Registration information

Note: Patient registration does not apply in countries where patient privacy laws conflict with providing patient information, including countries from the European Union.

A patient registration form is included in each device package. After implantation, please complete all requested information.

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An Implanted Device Identification Card is provided to the patient. The card contains the name and telephone number of the patient's physician, as well as information that medical personnel would require in the event of an emergency.

11. Postoperative information

11.1. MRI safety information 

Nonclinical testing has demonstrated the Tri-Ad 2.0 Adams tricuspid band is MR Conditional. A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:

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Under the scan conditions defined above, the Tri-Ad 2.0 Adams tricuspid band is expected to produce a maximum temperature rise of less than 2.2°C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the device extends approximately 20 mm from the Tri-Ad 2.0 Adams tricuspid band when imaged with a gradient echo pulse sequence and 3T MRI system.

12. Accessories

Use the Tri-Ad 2.0 Adams sizer set, Model 7900PS, to determine the appropriate band size. If using a handle, only use the Medtronic annuloplasty handle, Model 7686, to interface with the holder and sizers.

Caution: Do not use other manufacturers' annuloplasty sizers, or sizers for another Medtronic annuloplasty product, to size the Tri-Ad 2.0 Adams tricuspid band.

Caution: Do not use the accessories until they have been thoroughly cleaned and sterilized. Refer to the appropriate accessory instructions for use for further instructions.

13. Disclaimer of warranty

THE FOLLOWING DISCLAIMER OF WARRANTY APPLIES TO UNITED STATES CUSTOMERS ONLY:

ALTHOUGH THE TRI-AD 2.0 ADAMS TRICUSPID BAND, MODEL 900SFC. HEREAFTER REFERRED TO AS "PRODUCT," HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held to be illegal, unenforceable, or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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ALTHOUGH THE TRI-AD 2.0 ADAMS TRICUSPID BAND, MODEL 900SFC, HEREAFTER REFERRED TO AS "PRODUCT," HAS BEEN CAREFULLY DESIGNED, MANUFACTURED, AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

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© 2020 Medtronic
Trademarks may be registered and are the property of their respective owners.
Medtronic B.V.
Medtronic Australasia Pty Ltd
Medtronic of Canada Ltd
M973138A001 Rev. AB