P257202311 ASTE RU

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1) ESTELITE ASTERIA is a light-cured, radiopaque composite resin for use in anterior and posterior restorations. It is indicated for all carious lesion classes, including minimal or no-preparation type restorations. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049. ESTELITE ASTERIA contains 82% by weight (71% by volume) silica-zirconia filler and composite filler. A high filler load offers decreased polymerization shrinkage. Every inorganic filler contained in ESTELITE ASTERIA is a spherical filler (mean particle size: 200 nm, particle size range: 100 to 300 nm) that enables excellent gloss retention and wear resistance. ESTELITE ASTERIA contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), 1,6-bis(methacryl-ethyloxycarbonylamino)trimethyl hexane (UDMA), Triethylene glycol dimethacrylate (TEGDMA), Mequinol, Dibutyl hydroxyl toluene, and UV absorber.
2) ESTELITE ASTERIA incorporates Radical-Amplified Photopolymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time (90 sec.). Please see the table depicting the relationship between curing time and increment depth (■INDICATIONS FOR FILLING AND CURING).
3) ESTELITE ASTERIA is provided in SYRINGE.

■SHADE

ESTELITE ASTERIA contains 7 Body shades and 5 Enamel shades. The ESTELITE ASTERIA shading system is designed to mimic a wide shade range of the natural dentition by layering the appropriate shades selected from 2 categories: Body and Enamel.

Component	Shade	Tip
Body	A1B, A2B, A3B, A3.5B, A4B, B3B, BL	Body shades are designed for reconstructing the dentin layer. The Body shades should cover all enamel margins except the incisal area.
Enamel	Natural Enamel (NE)	NE is recommended to restore translucency in the incisal area in most cases.
	White Enamel (WE)	WE is recommended for the proximal wall. WE is suggested as an alternative to NE in whiter cases.
	Yellow Enamel (YE)	YE is designed to mimic discolored enamel.
	Trans Enamel (TE)	TE is the most translucent in ESTELITE ASTERIA. This shade is suggested as an alternative to NE in highly translucent cases.
	Occlusal Enamel (OcE)	OcE is recommended for the occlusal surface. OcE has exceptional sculptability to shape occlusal cusps and fissures.

■INDICATIONS

- Direct anterior and posterior restorations including occlusal surfaces
- Direct bonded composite veneer
- Diastema closure
- Repair of porcelain/composite

■CONTRAINDICATIONS

ESTELITE ASTERIA contains methacrylic monomers. DO NOT use ESTELITE ASTERIA for patients allergic to or hypersensitive to methacrylic and related monomers or to any of the other ingredients.

■PRECAUTIONS

1) DO NOT use ESTELITE ASTERIA for any purposes other than those listed in these instructions. Use ESTELITE ASTERIA only as directed herein.
2) ESTELITE ASTERIA is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
3) DO NOT use ESTELITE ASTERIA if the safety seals are broken or appear to have been tampered with.
4) If ESTELITE ASTERIA causes an allergic reaction or hypersensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl or latex) at all times when handling ESTELITE ASTERIA to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE ASTERIA comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of ESTELITE ASTERIA with eyes, mucosal membrane, skin and clothing.
 - If ESTELITE ASTERIA comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE ASTERIA comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If ESTELITE ASTERIA comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) ESTELITE ASTERIA should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of ESTELITE ASTERIA, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments and brushes with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
11) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

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ESTELITE ASTERIA is designed to be cured by either a halogen, LED or Plasma Arc light with a wavelength of 400-500 nm. Be sure to light-cure ESTELITE ASTERIA extra-orally and check the time needed for complete hardening of ESTELITE ASTERIA with your light-curing unit before performing the bonding procedure. The following table summarizes the relationship between curing time and increment depth.

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1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 "depth of cure" of ISO4049.
2) Curing time must be at least 20 seconds.

■STORAGE

1) Store ESTELITE ASTERIA at temperatures between 0 - 25°C (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use ESTELITE ASTERIA after the indicated date of expiration on the syringe package.

■DISPOSAL

To safely dispose of excess ESTELITE ASTERIA, extrude unused portion from SYRINGE and light-cure before disposal.

■CUSTOM SHADE GUIDE

See the instructions for ESTELITE ASTERIA Custom Shade Guide prior use.

■CLINICAL PROCEDURE

1. Cleaning
Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste, and then rinse with water.

2. Shade Selection
Select the appropriate shade of ESTELITE ASTERIA using the custom shade guide prepared in advance.
- Complete the shade selection within 5 minutes; teeth become whiter when dehydrated.
- Lightness (color value) is the primary importance on shade selection.
- In the case of whitened teeth, select the shade a few weeks later after the completion of whitening. Whitened teeth tend to become slightly darker with the elapse of time.
- Please see the table depicting the shade and tip (■SHADE).

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4. Cavity Preparation
Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), bevels assist in erasing demarcations between the cavity margins and the restoration, thereby enhancing both esthetics and retention.
- Scalloped (undulant) bevel may be desirable depending on the esthetics and retention.
- In the case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.
- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

5. Pulp Protection
Protect the pulp appropriately, if the cavity is in close proximity to the pulp. In the case of pulp exposure, apply calcium hydroxide for pulp capping. DO NOT USE EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing of ESTELITE ASTERIA.

6. Bonding System
Apply a bonding system according to its manufacturer's instructions.
- If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

7. Dispensing
Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately.
- DO NOT apply unnecessary force to the syringe immediately after removing from the refrigerator.

8. Filling and contouring
Fill the cavity incrementally. Increments should not exceed the indicated curing depth (please refer to the aforementioned table).
- In case that characterization is required, tints (such as ESTELITE COLOR) can be used. Apply the tint according to manufacturer's instructions.
- DO NOT mix ESTELITE ASTERIA with other brands of composite resin to avoid incomplete cure or entrapment of air bubbles. In order to avoid air bubbles entrapment, DO NOT mix ESTELITE ASTERIA with another shade of the paste.

9. Curing
Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing light tip within a distance of 2 mm from the increment.
- If other brands of composite resins are layered over the cured composite, follow its manufacturer's instructions.

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Light-cured Characterizing Material

ESTELITE COLOR

ENGLISH

CAUTION: Read all instructions, precautions and notes before using.

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1) ESTELITE COLOR is a low viscosity, light-cured tint for individualized shade characterizations of direct and indirect resin restorations. ESTELITE COLOR is applied below or in-between the layers of the composite restoration.
2) ESTELITE COLOR contains 58% by weight of silica-zirconia filler and silica-titania filler. The monomer matrix contains Bisphenol A polyethoxy methacrylate (Bis-MPEPP), Triethylene glycol dimethacrylate (TEGDMA), and 1,6-bis(methacryl-ethyloxycarbonylamino)trimethyl hexane (UDMA).

■SHADE

ESTELITE COLOR contains 9 Tints and 4 Opaquers.
Tints: Blue, Gray, Lavender, Ochre, Dark Brown, White, Yellow, Red, Clear
Opaquers: Pink Opaque, Low Chroma Opaque, Medium Chroma Opaque, High Chroma Opaque
Low Chroma Opaque (LCO): for shade range A1, B1, C1 and Bleached teeth
Medium Chroma Opaque (MCO): for shade range A2, B2, D2, and C2
High Chroma Opaque (HCO): for shade range A3, B3, D3, C3 and darker

■INDICATIONS

- Individualized shade characterizations of direct and indirect resin restorations.
- Masking of discoloration.

■CONTRAINDICATIONS

ESTELITE COLOR contains methacrylic monomers. DO NOT use ESTELITE COLOR for patients allergic to or hypersensitive to methacrylic and related monomers.

■PRECAUTIONS

1) DO NOT use ESTELITE COLOR for any purposes other than those listed in these instructions. Use ESTELITE COLOR only as directed herein.
2) ESTELITE COLOR is designed for sale and use by licensed dental care professionals only. It is not designed for sale or use by non-dental care professionals.
3) DO NOT use ESTELITE COLOR if the safety seals are broken or appear to have been tampered with.
4) If ESTELITE COLOR causes an allergic reaction or hypersensitivity, discontinue use immediately.
5) Use examination gloves (plastic, vinyl or latex) at all times when handling ESTELITE COLOR to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE COLOR comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of ESTELITE COLOR with eyes, mucosal membrane, skin and clothing.
 - If ESTELITE COLOR comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE COLOR comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If ESTELITE COLOR comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) ESTELITE COLOR should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of ESTELITE COLOR, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.

■STORAGE


1) Store ESTELITE COLOR at a temperature between 0 - 25℃ (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use ESTELITE COLOR after the indicated date of expiration on the syringe/package.

■DISPOSAL

To safely dispose of excess ESTELITE COLOR, extrude unused portion from the syringe and light-cure before disposal.

■ATTACHING THE SYRING TIP

1) Hold the nozzle and turn the cap counterclockwise for removal. [1]
2) Hold the nozzle and place the Syringe Tip onto it. Turn the Syringe Tip clockwise until it is securely locked. [2]
 - Accessories and Syringe Tip are designed for ESTELITE COLOR and PALFIQUE ESTELITE LV. Avoid using Syringe Tips from other brands with ESTELITE COLOR.
 - Syringe Tips are also sold separately.
 - Use gloves when attaching the Syringe Tip.

■CLINICAL PROCEDURE

Individualized shade characterizations of direct resin restorations

1. **Cavity preparation & Bonding system**
 Prepare the cavity and rinse with water. Apply a bonding system according to manufacturer's instructions.
 - If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.
2. **Dentin shade resin composite**
 Based on the multilayer approach, a dentin shade resin composite is applied over the bonded surface. If required, the mamelons or fissures are sculpted with a fine tipped instrument at this stage. Then the dentin is light cured according to manufacturer's instructions.
3. **Dispensing the ESTELITE COLOR**
 Extrude the small amount of ESTELITE COLOR on the Disposable Dish. If required, mix ESTELITE COLOR with other shades of ESTELITE COLOR. After extruding the material, wipe the nozzle with gauze and replace the cap immediately.
 - Tips for shade selection
 - Incisal edge: Blue (youthful look), Gray (mature look), Lavender (elderly look)
 - Increasing saturation: Ochre (fossae areas), Dark Brown (pits and fissures: stained craze lines)
 - White spots: White
 - Color adjustment: Yellow, Ochre, Red (mixed in varying proportions)
 - Chroma adjustment: Clear (mixed with another shade)
 - Neutralizing grayness: Pink Opaque (mixed with LCO, MCO, or HCO)
 Tints can be mixed in varying proportions to achieve the correct hue, value and chroma, as needed. Opaquers can be mixed in varying proportions to achieve the correct value and chroma, as needed.
 - Protect the dispensed material from ambient light before the application using a Light Shield Cover.
 - DO NOT mix ESTELITE COLOR with other brands of products to avoid incomplete cure.
4. **Applying the ESTELITE COLOR**
 Using a fine tipped instrument or a fine tipped brush, apply ESTELITE COLOR sparingly. The maximum thickness is 0.5 mm (Dark Brown: 0.05 mm).
5. **Curing**
 Use a light curing unit with a wavelength range of 400 to 500 nm (peak: 470 nm, intensity ≥ 500 mW/cm²). Light-cure for at least 30 seconds, ensure the curing light tip does not exceed 2 mm from the surface.
6. **Enamel shade resin composite**
 An enamel shade resin composite is applied over the colored surface. Then the enamel is light cured; finish and polish according to manufacturer's instructions. ESTELITE COLOR must be placed below or in-between the layers of the composite restoration.

Masking of discoloration

1. **Cavity preparation & Bonding system**
 Prepare the cavity and rinse with water. Apply a bonding system according to manufacturer's instructions.
 - If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.
2. **Dispensing the ESTELITE COLOR**
 Extrude the small amount of ESTELITE COLOR on the Disposable Dish. If required, mix ESTELITE COLOR with other shades of ESTELITE COLOR. After extruding the material, wipe the nozzle with gauze and replace the cap immediately.
 - Tips for shade selection
 - Masking dark areas: LCO, MCO, HCO
 - Protect the dispensed material from ambient light before the application using a Light Shield Cover.
 - DO NOT mix ESTELITE COLOR with other brands of products to avoid incomplete cure.
3. **Applying ESTELITE COLOR**
 Using a fine tipped brush or other suitable instruments, apply the appropriate Opaquer of ESTELITE COLOR evenly over the discolored areas to elevate the color value. The maximum thickness is 0.5 mm for each increment.
 - DO NOT apply ESTELITE COLOR directly onto the tooth without application of bonding agent.
 - Tips for thickness of Opaquer
 The required average thickness of Opaquer is 0.35 mm to mask the discolored areas adequately.
4. **Curing**
 Use a light curing unit with a wavelength range of 400 to 500 nm (peak: 470 nm, intensity ≥ 500 mW/cm²). Light-cure for at least 30 seconds, ensure the curing light tip does not exceed 2 mm from the surface.
5. **Resin composite**
 A resin composite is applied over the material. Then the composite is light cured; finish and polish according to manufacturer's instructions. ESTELITE COLOR must be placed below or in between the layers of the composite restoration.

Individualized shade characterizations of indirect resin restorations

Use indirect composite resin system according to manufacture's instructions. During the fabrication, ESTELITE COLOR is placed between the layers of resin-based restorative material according to the procedure described for direct restorations (direct staining technique).

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■ PRODUCT DESCRIPTION AND GENERAL INFORMATION

1) OMNICHROMA is a light-cured, radiopaque composite resin for use in anterior and posterior restorations and is indicated for all carious classes. OMNICHROMA is a single shade material. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049.
2) OMNICHROMA contains 79% by weight (68% by volume) of spherical silica-zirconia filler (mean particle size: 0.3 μm, particle size range: 0.2 to 0.4 μm) and composite filler. OMNICHROMA contains 1,6-bis(methacryl-ethyloxycarbonylamino) trimethyl hexane (UDMA), Triethylene glycol dimethacrylate (TEGDMA), Mequinol, Dibutyl hydroxyl toluene and UV absorber.
3) OMNICHROMA is available in either SYRINGE or Pre Loaded Tip (PLT).

■ INDICATIONS

For use as a tooth shade resin material in dental procedures, such as:
- Direct anterior and posterior restorations
- Direct bonded composite veneer
- Diastema closure
- Repair of porcelain/composite

■ CONTRAINDICATIONS

OMNICHROMA contains methacrylic monomers and UV absorbers. DO NOT use OMNICHROMA for patients allergic to or hypersensitive to methacrylic and related monomers, UV absorbers or any of the other ingredients.

■ PRECAUTIONS

1) DO NOT use OMNICHROMA for any purpose other than those listed in these instructions. Use OMNICHROMA only as directed herein.
2) OMNICHROMA is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
3) DO NOT use OMNICHROMA if the safety seals are broken or appear to have been tampered with.
4) If OMNICHROMA causes an allergic reaction or hypersensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl or latex) at all times when handling OMNICHROMA to avoid the possibility of allergic reactions from methacrylic monomers.
 Notes: Certain substance/materials may penetrate through examination gloves. If OMNICHROMA comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of OMNICHROMA with eyes, mucosal membrane, skin and clothing.
 - If OMNICHROMA comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If OMNICHROMA comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If OMNICHROMA comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) OMNICHROMA should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of OMNICHROMA, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments and brushes with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
11) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.
12) Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.
13) DO NOT mix OMNICHROMA with other brands of resin composite to avoid losing its color adjustment.

■ INDICATIONS FOR FILLING AND CURING

OMNICHROMA is designed to be used with either a halogen or LED curing-light with a wavelength of 400-500 nm. Be sure to light-cure OMNICHROMA extra-orally and check the time needed for complete hardening of OMNICHROMA with your light-curing unit before performing the bonding procedure. The below table illustrates the relationship between curing time and increment depth.

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1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 "depth of cure" of ISO 4049.

■ SPECIAL NOTES FOR THE USE OF PLT

1) PLTs are designed for single patient use only. Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.
2) Dispensers are not provided for OMNICHROMA. Use a dispenser that fits the PLT of OMNICHROMA. For proper use and proper disinfection, see manufacturer's instructions.
3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

■ STORAGE

1) Store OMNICHROMA at temperature between 0 - 25°C (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use OMNICHROMA after the indicated date of expiration on the syringe or PLT package.

■ DISPOSAL

To safely dispose of excess OMNICHROMA, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

■ CLINICAL PROCEDURE

Determine if this product is suitable for application to the case before use. The desirable esthetic outcome may not be obtained depending on the case (e.g. a case in which no tooth structure remains on the lingual surface, a case in which the underlying tooth is extremely discolored, a case in which the color of the underlying tooth is to be changed to another color, etc.).

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2. Isolation
A rubber dam is the preferred method of isolation.

3. Cavity Preparation
Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II), bevels and chamfers assist in ensuring demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.
- Scalloped bevel could be desirable depending on the esthetics and retention.
- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent, followed by a metal primer if metal is exposed on the fractured surfaces according to manufacturer's instructions.

4. Pulp Protection
Glass ionomer or a calcium hydroxide lining should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing OMNICHROMA. OMNICHROMA is not indicated for direct pulp capping.

5. Bonding System
Apply a bonding system according to manufacturer's instructions.
- If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

6. Dispensing
6-1. PLT
Please read ■ **SPECIAL NOTES FOR THE USE OF PLT** before use.
- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.
- Extrude the paste directly into the cavity or indirectly on the mixing pad.

6-2. SYRINGE
Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately.
- DO NOT apply unnecessary force to the syringe immediately after removing from the refrigerator.

7. Filling and Contouring
Fill the cavity incrementally. Increments should not exceed the indicated curing depth (please refer to the aforementioned table).
- Use OMNICHROMA BLOCKER under OMNICHROMA for the direct restoration of extensive Class III or IV cavity. OMNICHROMA can appear semi-translucent when OMNICHROMA is not adequately surrounded by a cavity wall. If applied in direct restorations of extensive Class III or IV cavity, the cured OMNICHROMA may appear dark.
- In the case of repair of porcelain / composite, the application of an opaquer to the fractured surfaces as an initial layer is needed to hide a metal portion if it exists. Then, place OMNICHROMA into the fractured surfaces.
- In the case of masking tooth discoloration, apply OMNICHROMA BLOCKER or a dental color masking material such as an opaquer on the discolored surface to mask the discoloration, then place OMNICHROMA.

8. Curing
Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing light tip within a distance of 2 mm from the increment.

9. Finishing
Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

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1) OMNICHROMA BLOCKER is a light-cured, radiopaque, composite resin. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049. OMNICHROMA BLOCKER contains 82% by weight (71% by volume) of silica-zirconia filler and composite filler. A high filler load offers decreased polymerization shrinkage. All inorganic filler contained in OMNICHROMA BLOCKER is a spherical filler (mean particle size: 0.2 μm, particle size range: 0.1 to 0.3 μm). The monomer matrix contains Bis-GMA and Triethylene glycol dimethacrylate.
2) OMNICHROMA BLOCKER is provided either in SYRINGE or Pre Loaded Tip (PLT).
3) OMNICHROMA BLOCKER is designed for use as a thin layer at the lingual cavity wall of extensive class III and IV restorations. OMNICHROMA BLOCKER can mask slight discoloration or reconstruct a highly opaque tooth.

■ INDICATIONS

For use as a tooth shade resin in dental procedures.

■ CONTRAINDICATIONS

OMNICHROMA BLOCKER contains methacrylic monomers. DO NOT use OMNICHROMA BLOCKER for patients allergic to or hypersensitive to methacrylic and related monomers.

■ PRECAUTIONS

1) DO NOT use OMNICHROMA BLOCKER for any purpose other than those listed in these instructions. Use OMNICHROMA BLOCKER only as directed herein.
2) OMNICHROMA BLOCKER is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
3) DO NOT use OMNICHROMA BLOCKER if the safety seals are broken or appear to have been tampered with.
4) If OMNICHROMA BLOCKER causes an allergic reaction or oversensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl, or latex) at all times when handling OMNICHROMA BLOCKER to avoid the possibility of allergic reactions from methacrylic monomers.
 Note: Certain substances/materials may penetrate through examination gloves. If OMNICHROMA BLOCKER comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of OMNICHROMA BLOCKER with eyes, mucosal membrane, skin and clothing.
 - If OMNICHROMA BLOCKER comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If OMNICHROMA BLOCKER comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If OMNICHROMA BLOCKER comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) OMNICHROMA BLOCKER should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of OMNICHROMA BLOCKER, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.
11) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

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OMNICHROMA BLOCKER is designed to be used with either a halogen or LED curing-light with a wavelength of 400-500 nm. Be sure to light-cure OMNICHROMA BLOCKER extra-orally and check the time needed for complete hardening of OMNICHROMA BLOCKER with your light-curing unit before performing the bonding procedure. The below table illustrates the relationship between curing time and increment depth.

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■ SPECIAL NOTES FOR THE USE OF PLT

1) PLTs are designed for single patient use only. Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.
2) Dispensers are not provided for OMNICHROMA BLOCKER. Use a dispenser that fits the PLT of OMNICHROMA BLOCKER. For proper use and proper disinfection, see manufacturer's instructions.
3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

■ STORAGE

1) Store OMNICHROMA BLOCKER at a temperature between 0 - 25°C (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use OMNICHROMA BLOCKER after the indicated date of expiration on the syringe or PLT package.

■ DISPOSAL

To safely dispose of excess OMNICHROMA BLOCKER, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

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2. Shade Selection
Select OMNICHROMA BLOCKER as the initial layer under OMNICHROMA, with direct restorations of extensive Class III or IV cavity using OMNICHROMA.

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4. Cavity Preparation
Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II), bevels and chamfers assist in erasing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.
- Scalloped bevel could be desirable depending on the esthetics and retention.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent, followed by a metal primer if metal is exposed on the fractured surfaces according to manufacturer's instructions.

5. Pulp Protection
Glass ionomer or a calcium hydroxide lining should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of bonding agent and/or OMNICHROMA BLOCKER. OMNICHROMA BLOCKER is not indicated for direct pulp capping.

6. Bonding System
Apply a bonding system according to manufacturer's instructions.
- If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.

7. Dispensing
7-1. PLT
Please read ■ SPECIAL NOTES FOR THE USE OF PLT before use.
- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.
- Extrude the paste directly into cavity or indirectly on the mixing pad.

7-2. SYRINGE
Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. After dispensing, turn the handle counterclockwise by a half to full turn to release residual pressure inside the syringe and re-cap the syringe immediately after dispensing.
- DO NOT apply unnecessary force to the syringe immediately after removing from the refrigerator.

8. Filling
Fill the cavity with OMNICHROMA BLOCKER as an initial layer. Increments should not exceed the indicated depth (please refer to the aforementioned table).
- DO NOT mix OMNICHROMA BLOCKER with other brands of resin composite to avoid incomplete cure or voids.

9. Curing
Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing-light tip within a distance of 2 mm from the increment.
- Follow the filling and curing instructions of OMNICHROMA to be layered over the cured OMNICHROMA BLOCKER.

10. Finishing
Shape and polish the restoration according to instructions for OMNICHROMA used as the surface layer of the restoration.

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1) ESTELITE Σ QUICK is a light-cured, radiopaque, submicron filled composite resin for use in anterior and posterior restorations and is indicated for all carious classes. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049. ESTELITE Σ QUICK contains 82% by weight (71% by volume) of silica-zirconia filler and composite filler. A high filler load offers decreased polymerization shrinkage. All inorganic filler contained in ESTELITE Σ QUICK is a spherical submicron filler (mean particle size: 0.2 μm, particle size range: 0.1 to 0.3μm) that is designed to deliver excellent gloss retention and wear resistance. The monomer matrix contains Bis-GMA and Triethylene glycol dimethacrylate.
2) ESTELITE Σ QUICK implements Radical-Amplified Photopolymerization initiator technology (RAP technology). Curing time for ESTELITE Σ QUICK can be shortened by one-third using a light curing unit with a wavelength range of 400 to 500nm (peak: 470nm) compared to our conventional resin composite. Please see the table depicting the relationship between curing time and increment depth (■INDICATIONS FOR FILLING AND CURING).
3) RAP technology facilitates a longer working time for ESTELITE Σ QUICK.
4) ESTELITE ΣQUICK is provided either in SYRINGE or Pre Loaded Tip(PLT).

■SHADE

A1, A2, A3, A3.5, A4, A5, B1, B2, B3, B4, C1, C2, C3, OA1, OA2, OA3, OPA2, BW (Bleach White), WE (White Enamel) and CE (Clear Enamel/Incisal).
- Opalescent shades (OA1, OA2, OA3) have adequate opacity designed for blocking out the dark shine that comes through the oral cavity (Class III and Class IV restoration). The opalescent shades also serve as a dentin layer for the multi-shade layering technique. Shades OA1, OA2 and BW closely match those of primary teeth. However, they are not designed for masking metal or a darkly stained area such as tetracycline stained.
- OPA2 is designed for masking slight stain or reconstructing a highly opaque tooth fuse as a thin layer at the lingual cavity wall of class III and IV).
- Please check with your local distributor for available shades.

■INDICATIONS

- Direct anterior and posterior restorations including the occlusal surface and all carious classes
- Direct bonded composite veneer
- Diastema closure
- Repair of porcelain/composite

■CONTRAINDICATIONS

ESTELITE Σ QUICK contains methacrylic monomers. DO not use ESTELITE Σ QUICK for patients allergic to or hypersensitive to methacrylic and related monomers.

■PRECAUTIONS

1) DO NOT use ESTELITE Σ QUICK for any purpose other than those listed in these instructions. Use ESTELITE Σ QUICK only as directed herein.
2) ESTELITE Σ QUICK is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
3) DO NOT use ESTELITE Σ QUICK if the safety seals are broken or appear to have been tampered with.
4) If ESTELITE Σ QUICK causes an allergic reaction or oversensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl, or latex) at all times when handling ESTELITE Σ QUICK to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE Σ QUICK comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of ESTELITE Σ QUICK with eyes, mucosal membrane, skin and clothing.
- If ESTELITE Σ QUICK comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
- If ESTELITE Σ QUICK comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration.
- If ESTELITE Σ QUICK comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
- Instruct the patient to rinse his mouth immediately after treatment.
7) ESTELITE Σ QUICK should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of ESTELITE Σ QUICK, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.

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Be sure to light-cure ESTELITE Σ QUICK extraorally and check the time needed for complete hardening of ESTELITE Σ QUICK with your light curing unit before performing the bonding procedure. The following table illustrates the relationship between curing time and increment depth.

Relationship between curing time and increment depth

Light type	Intensity (mW/cm²)	Curing time (seconds)	Increment depth (mm)[1] A1, A2, A3, A3.5, B1, B2, B3, C1, C2, OA1, OA2, CE(Inc.), BW, WE	Increment depth (mm)[1] A4, A5, B4, C3, OA3, OPA2
Halogen	400	10	1.9	1.6
		20	2.2	1.9
	600	10	2.0	1.7
		20	2.5	2.0
	800	5	1.8	1.6
		10	2.2	1.9
LED	260	10	1.7	1.4
		20	2.0	1.7
	900	5	1.7	1.5
		10	2.0	1.8
Plasma Arc	950	3	1.9	1.5

1) Increment depth was determined on the basis of test results performed in accordance with section 7.10 "depth of cure" of ISO4049.

■SPECIAL NOTES FOR THE USE OF PLT

1) PLTs are designed for single patient use only. Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.
2) Dispensers are not provided for ESTELITE Σ QUICK. Use a dispenser that fits the PLT of ESTELITE Σ QUICK. For proper use and proper disinfection, see manufacturer's instructions.
3) Apply the composite at normal room temperature (18 - 30°C / 62 - 84°F). Material that is too cold might be difficult to extrude.
4) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

■STORAGE

1) Store ESTELITE Σ QUICK at a temperature between 0 - 25°C (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use ESTELITE Σ QUICK after the indicated date of expiration on the syringe or PLT package.

■DISPOSAL

To safely dispose of excess ESTELITE Σ QUICK, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

■CLINICAL PROCEDURE

1. Cleaning
Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.
2. Shade Selection
Select the appropriate shade of ESTELITE Σ QUICK prior to placing the rubber dam using the shade guide provided for ESTELITE Σ QUICK.
- Complete the shade selection within 5 minutes, desiccated teeth are lighter than wet or moist teeth.
- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade.
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

3. Isolation
A rubber dam is the preferred method to isolate the tooth from contamination.
4. Cavity Preparation
Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in erasing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.
- Scalloped (undulant) bevel could be desirable depending on the esthetics and retention.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

5. Pulp Protection
Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE Σ QUICK.
6. Bonding system
Apply the adhesive system according to the manufacturer's instructions.
- If you are using self-curing adhesive systems or dual-curing systems, make sure that they can be used with light-curing composite materials.

7. Dispensing
7-1. PLT
Please read ■SPECIAL NOTES FOR THE USE OF PLT before use.
- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.
- Extrude the paste directly into cavity or indirectly on the mixing pad.

7-2. SYRINGE
Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. Re-cap the syringe immediately after dispensing.
- DO NOT apply unnecessary force to the syringe immediately after removing from the refrigerator.

8. Filling and contouring
Fill the cavity incrementally. Increments should not exceed the indicated depth (please refer to the aforementioned table).
- DO NOT mix ESTELITE Σ QUICK with other brands of resin composite or other shades of ESTELITE Σ QUICK to avoid incomplete cure or voids.

9. Curing
Light-cure each increment for at least the indicated time (please see aforementioned table), keeping the curing light tip within a distance of 2 mm from the increment.
- If other brands of composite resins are layered over the cured composite, follow their instructions.

10. Finishing
Shape and polish the restoration. For finishing, use fine finishing diamond points. To erase the demarcation between the composite and the natural tooth rotate a 12 fluted carbide bur slowly along the margins without water. Use finishing metal strips or vinyl strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool, such as PoGo™ (DENTSPLY/Caulk), Sof-Lex™ (3M ESPE), Identoflex® HiLuster Dia Polishers (Kerr), or D-FINE™ Hybrid Diamond (COSMEDENT, Inc.).
* PoGo, Sof-Lex, Identoflex and D-FINE are not registered trademarks of Tokuyama Dental Corporation.

CAUTION : Federal law restricts this device to sale by or on the order of a dental professional.

IMPORTANT NOTE: The manufacturer of ESTELITE Σ QUICK is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.
Product specifications of ESTELITE Σ QUICK are subjected to change without notification. When the product specification change, the instruction and precautions may change also.

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1) ESTELITE POSTERIOR is a light-cured, radiopaque composite resin for use in posterior restorations. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049. ESTELITE POSTERIOR contains 84% by weight (70% by volume) of silica-zirconia filler (mean particle size: 2 µm, particle size range: 0.1 to 10µm). A high filler load offers decreased polymerization shrinkage. The monomer matrix contains Bis-GMA, Triethylene glycol dimethacrylate (TEGDMA) and Bisphenol A polyethoxy methacrylate (Bis-MPEPP).
2) ESTELITE POSTERIOR implements Radical-Amplified Photopolymerization initiator technology (RAP technology). Curing time for ESTELITE POSTERIOR can be shortened by one-third using a light curing unit with a wavelength range of 400 to 500nm (peak: 470nm) compared to our conventional resin composite. Please see the table depicting the relationship between curing time and increment depth (■ INDICATIONS FOR FILLING AND CURING).
3) RAP technology facilitates a longer working time for ESTELITE POSTERIOR.
4) ESTELITE POSTERIOR is provided in SYRINGE.

■SHADE

PCE, PA1, PA2, PA3.
- Every shade provided for ESTELITE POSTERIOR has adequate translucency designed for posterior restoration.
- PCE serves as a top layer for the multi-shade layering technique.
- Tokuyama Dental Corporation does not provide ESTELITE POSTERIOR Shade Guide. Please make a mockup shade guide by using ESTELITE POSTERIOR.

■INDICATIONS

- Direct posterior restorations including the occlusal surface
- Repair of porcelain/composite

■CONTRAINDICATIONS

ESTELITE POSTERIOR contains methacrylic monomers. DO not use ESTELITE POSTERIOR for patients allergic to or hypersensitive to methacrylic and related monomers.

■PRECAUTIONS

1) DO NOT use ESTELITE POSTERIOR for any purpose other than those listed in these instructions. Use ESTELITE POSTERIOR only as directed herein.
2) ESTELITE POSTERIOR is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
3) DO NOT use ESTELITE POSTERIOR if the safety seals are broken or appear to have been tampered with.
4) If ESTELITE POSTERIOR causes an allergic reaction or oversensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl, or latex) at all times when handling ESTELITE POSTERIOR to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE POSTERIOR comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of ESTELITE POSTERIOR with eyes, mucosal membrane, skin and clothing.
 - If ESTELITE POSTERIOR comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE POSTERIOR comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration.
 - If ESTELITE POSTERIOR comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) ESTELITE POSTERIOR should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of ESTELITE POSTERIOR, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.

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Be sure to light-cure ESTELITE POSTERIOR extraorally and check the time needed for complete hardening of ESTELITE POSTERIOR with your light curing unit before performing the bonding procedure. The following table illustrates the relationship between curing time and increment depth.

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■STORAGE

1) Store ESTELITE POSTERIOR at a temperature between 0 - 25°C (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use ESTELITE POSTERIOR after the indicated date of expiration on the syringe or package.

■DISPOSAL

To safely dispose of excess ESTELITE POSTERIOR, extrude unused portion from SYRINGE and light-cure before disposal.

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2. Shade Selection
Select the appropriate shade of ESTELITE POSTERIOR prior to placing the rubber dam.
- Complete the shade selection within 5 minutes; desiccated teeth are lighter than wet or moist teeth.
- Lightness (color value) is primary importance on shade selection.
- When wavering in shade selection, select the lighter (high value) shade.
- In case of whitened teeth select the shade a few weeks after the completion of whitening. Whitened teeth may become slight darker.

3. Isolation
A rubber dam is the preferred method to isolate the tooth from contamination.

4. Cavity Preparation
Prepare the cavity and rinse with water. Add chamfers to the margins of posterior preparations (class I, II) because chamfers assist in diminishing demarcations between the cavity margins and restoration, thereby enhancing both esthetics and retention.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning, rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent according to manufacturer's instructions.

5. Pulp Protection
Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of bonding agent and/or ESTELITE POSTERIOR.

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7. Dispensing
Remove the syringe cap. Extrude the paste onto a mixing pad by turning the handle clockwise. Re-cap the syringe immediately after dispensing.
- DO NOT apply unnecessary force to the syringe immediately after removing from the refrigerator.

8. Filling and contouring
Fill the cavity incrementally using a suitable placing instrument. Increments should not exceed the indicated depth (please refer to the aforementioned table).
- DO NOT mix ESTELITE POSTERIOR with other brands of resin composite or other shades of ESTELITE POSTERIOR to avoid incomplete cure or entrapment of air bubbles.

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10. Finishing
Shape and polish the restoration. For finishing, use fine finishing diamond points. To diminish the demarcation between the composite and the natural tooth, rotate a 12 fluted carbide bur slowly along the margins without water. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tool, such as PoGo™ (DENTSPLY/Caulk), Sof-Lex™ (3M ESPE), Identoflex® HiLuster Dia Polishers (Kerr), or D-FINE™ Hybrid Diamond (COSMEDENT, Inc.).
* PoGo, Sof-Lex, Identoflex and D-FINE are not registered trademarks of Tokuyama Dental Corporation.

IMPORTANT NOTE: The manufacturer of ESTELITE POSTERIOR is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.
Product specifications of ESTELITE POSTERIOR are subjected to change without notification. When the product specification change, the instruction and precautions may change also.

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PALFIQUE ESTELITE PASTE SYRINGE

ENGLISH

Read all information, precautions, notes and materials before using PALFIQUE ESTELITE PASTE.

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1) PALFIQUE ESTELITE PASTE is a light cured submicron filled composite resin. PALFIQUE ESTELITE PASTE contains 82 % by weight (71 % by volume) of submicron (0.2μm) Silica-zirconia spherical filler and composite filler. The monomer matrix contains Bis-GMA and Triethylene glycol dimethacrylate.
2) PALFIQUE ESTELITE PASTE is radiopaque. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049.
3) PALFIQUE ESTELITE PASTE is provided in SYRINGE.

■SHADE

A1, A2, A3, A3.5, B2, C2, Incisal (Inc.), OA2 and OA3
- Tokuyama Dental Corporation does not provide PALFIQUE ESTELITE PASTE Shade Guide. Please make a mockup shade guide by using PALFIQUE ESTELITE PASTE.

■INDICATIONS

Direct anterior and posterior restorations including the occlusal surface and all carious classes.

■CONTRAINDICATIONS

PALFIQUE ESTELITE PASTE contains methacrylic monomers. DO not use PALFIQUE ESTELITE PASTE for patients allergic to or hypersensitive to methacrylic and related monomers.

■CLINICAL APPLICATION

1) Clean and dry the cavity preparation. If pulpal protection is required, use a lining material. DO NOT use Eugenol based product. Eugenol based products will inhibit curing of the bonding agent and/or composite resin.
2) Apply the adhesive system according to the manufacturer's instructions.
 - If you are using self-curing adhesive systems or dual-curing systems, make sure that they can be used with light-curing composite materials.
3) Place an adequate amount of PALFIQUE ESTELITE PASTE on the mixing pad and fill the cavity. Light irradiation should be for a minimum of 30 seconds for each layer (≤ 2 mm).
4) Shape and polish the restoration.
 - Instruct the patient adequately on the maintenance of the restoration after treatment.

■STORAGE

1) Store PALFIQUE ESTELITE PASTE at a temperature between 0-23°C (32-73°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use PALFIQUE ESTELITE PASTE after the indicated date of expiration.

■DISPOSAL

To safely dispose of excess PALFIQUE ESTELITE PASTE, extrude unused product from syringe with filling instruments and light cure before disposal.

■PRECAUTIONS

1) Use PALFIQUE ESTELITE PASTE only as directed herein.
2) PALFIQUE ESTELITE PASTE is designed for sale and use by licensed dental care professions only. It is not designed for sale and nor is it suitable for use by non-dental care professionals.
3) DO NOT use PALFIQUE ESTELITE PASTE for any purpose other than dental instructions as noted herein.
4) DO NOT use this product if the safety seals are broken or appear to have been tampered with.
5) If PALFIQUE ESTELITE PASTE causes an allergic reaction or oversensitivity, discontinue its use immediately.
6) The use of examination gloves (plastic, vinyl, latex) is recommended at all times to avoid the possibility of allergic reactions from methacrylic monomers. If PALFIQUE ESTELITE PASTE comes in contact with the examination gloves, immediately dispose of the examination gloves and wash hands thoroughly with water.
7) AVOID CONTACT WITH THE EYES. In case of contact with eyes, immediately flush eyes with water and seek medical attention.
8) When using a light-curing unit, protective eye shields, glasses or goggles should be worn.
9) DO NOT mix PALFIQUE ESTELITE PASTE with other brands of composite resin to avoid incomplete cure or air bubble contamination. In order to avoid air bubble contamination, DO NOT mix PALFIQUE ESTELITE PASTE with another shade of the paste.
10) DO NOT use Eugenol materials with PALFIQUE ESTELITE PASTE.
11) In the event that PALFIQUE ESTELITE PASTE contacts the skin, clean the skin immediately with an alcohol soaked cotton swab.
12) Protect pulp with a lining material, if the pulp is exposed.
13) PALFIQUE ESTELITE PASTE should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
14) DO NOT apply unnecessary force to the syringes, especially immediately after taking out of a refrigerator.
15) DO NOT leave PALFIQUE ESTELITE PASTE unnecessarily within the reach of patients and children to avoid being swallowed by mistake.
16) Clean the filling instruments with alcohol after use.

The user and/or patient should report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure the product is suitable for an appropriate application before use.
Specifications are subjected to change without notification.

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1) ESTELITE BULK FILL Flow is a low viscosity, light-cured, radiopaque composite resin. This low stress flowable material can be placed 4mm increments. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type 1 and Class 2 (Group 1) material by ISO 4049.
2) ESTELITE BULK FILL Flow incorporates Radical-Amplified Photopolymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time.
3) ESTELITE BULK FILL Flow contains 70% by weight (56% by volume) of spherical silica-zirconia filler and composite filler. The composite filler offers reduced polymerization shrinkage stress due to its round shape. The spherical filler contained in ESTELITE BULK FILL Flow (mean particle size: 200nm, particle size range: 100 to 300nm) facilitates excellent gloss retention and wear resistance. ESTELITE BULK FILL Flow contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate (TEGDMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), Mequinol, Dibutyl hydroxyl toluene and UV absorber.
4) ESTELITE BULK FILL Flow is provided either in SYRINGE or Pre Loaded Tip (PLT).

■SHADE

U (Universal), A1, A2, A3, B1.
- Please check with your local distributor for available shades.

■INDICATIONS

- Direct anterior and posterior restorations
- Cavity lining
- Blocking out cavity undercuts before fabricating indirect restorations
- Repair of porcelain/composite

■CONTRAINDICATIONS

ESTELITE BULK FILL Flow contains methacrylic monomers. DO NOT use ESTELITE BULK FILL Flow for patients allergic to or hypersensitive to methacrylic and related monomers or to any of the other ingredients.

■PRECAUTIONS

1) DO NOT use ESTELITE BULK FILL Flow for any purposes other than those listed in these instructions. Use ESTELITE BULK FILL Flow only as directed herein.
2) ESTELITE BULK FILL Flow is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
3) DO NOT use ESTELITE BULK FILL Flow if the safety seals are broken or appear to have been tampered with.
4) If ESTELITE BULK FILL Flow causes an allergic reaction or hypersensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl, or latex) at all times when handling ESTELITE BULK FILL Flow to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE BULK FILL Flow comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of ESTELITE BULK FILL Flow with eyes, mucosal membrane, skin and clothing.
 - If ESTELITE BULK FILL Flow comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE BULK FILL Flow comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If ESTELITE BULK FILL Flow comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) ESTELITE BULK FILL Flow should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of ESTELITE BULK FILL Flow, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments and brushes with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
11) ESTELITE BULK FILL Flow is available in semi-translucent shades. If applied in cavity Class III or IV direct restorations, the cured ESTELITE BULK FILL Flow may appear dark. Avoid using ESTELITE BULK FILL Flow for the direct restoration of Class III or IV when opacity is paramount.
12) Light-cure ESTELITE BULK FILL Flow extra orally to confirm the time needed for complete hardening of ESTELITE BULK FILL Flow with your light-curing unit prior to doing the restoration.
13) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

■CURING TIME AND DEPTH

ESTELITE BULK FILL Flow enables a 4mm depth of cure (U shade) by light-curing for 10 seconds with Halogen or LED light-curing unit (Light intensity: 800mW/cm^2 or more, wavelength: 400-500nm). The depth of cure was determined based on test results performed in accordance with section 7.10 "depth of cure" of ISO 4049.

■SPECIAL NOTES FOR THE USE OF PLT

1) PLTs are designed for single patient use only. Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.
2) Dispensers are not provided for ESTELITE BULK FILL Flow. Use a dispenser that fits the PLT of ESTELITE BULK FILL Flow. For proper use and proper disinfection, see manufacturer's instructions.
3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

■STORAGE

1) Store ESTELITE BULK FILL Flow at a temperature between 0 - 25°C (32 - 77°F). Shelf life at this temperature is 3 years.
2) AVOID direct exposure to light and heat.
3) DO NOT use ESTELITE BULK FILL Flow after the indicated date of expiration on the syringe or PLT package.

■DISPOSAL

To safely dispose of excess ESTELITE BULK FILL Flow, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

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2. Shade Selection
Prior to isolation of tooth, select the appropriate shade of ESTELITE BULK FILL Flow using shade card provided for ESTELITE BULK FILL Flow.

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
4. Cavity Preparation
Prepare the cavity and rinse with water.
- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent, followed by a metal primer if metal is exposed on the fractured surfaces according to manufacturer's instructions.

5. Pulp Protection
Calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing ESTELITE BULK FILL Flow. ESTELITE BULK FILL Flow is not indicated for direct pulp capping.

6. Bonding system
Apply the adhesive system according to the manufacturer's instructions.
- If you are using self-curing adhesive systems or dual-curing systems, make sure that they can be used with light-curing composite materials.

7. Dispensing
7-1. PLT
Please read ■SPECIAL NOTES FOR THE USE OF PLT before use.
- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.

7-2. SYRINGE
1) Hold the nozzle and turn the cap counterclockwise for removal. [1]
2) Hold the nozzle and place the Dispensing Tip onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

- Accessories and Dispensing Tip are exclusively designed for ESTELITE BULK FILL Flow. Avoid using Dispensing Tips from other brands with ESTELITE BULK FILL Flow.
- Dispensing Tips are also sold separately.
- Use uncontaminated gloves when handling the syringe and attaching the Dispensing Tip.

8. Filling and contouring
Cavity can be directly filled from the Dispensing Tip/PLT. Increments should not exceed 4mm at a time.
- DO NOT apply excessive force to extrude the material from the syringe or PLT immediately after removal from the refrigerator.
- DO NOT mix ESTELITE BULK FILL Flow with other brands of resin composite or other shades of ESTELITE BULK FILL Flow to avoid incomplete cure or entrapment of air bubbles.
- In case of syringe, after extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.
- This product is provided non-sterile and is not intended to be sterilized by the user.

9. Curing
Light-cure each increment for at least the indicated time (please see the curing recommendations below). Light-cure each layer (up to 4mm increment) using Halogen light-curing or LED light-curing unit with a wavelength of 400-500 nm. Be sure the light-curing unit with an intensity of 800mW/cm^2 or more. Hold the light guide tip as close to the restorative as possible during light exposure.

CURING RECOMMENDATIONS

Shade	Increment Depth	Light Output	Curing Time
All shades	4mm	≥800mW/cm^2	10 sec.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

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Resin-based Dental Restorative Material

ESTELITE UNIVERSAL FLOW

ENGLISH

Carefully read all information, precautions, notes and materials before using ESTELITE UNIVERSAL FLOW.

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1) ESTELITE UNIVERSAL FLOW is a low viscosity, light-cured, radiopaque composite resin for use in anterior and posterior restorations. Direct placement of ESTELITE UNIVERSAL FLOW into a prepared cavity allows for easy handling. It has radiopacity equal to or greater than 1mm aluminum, and is classified as a Type I and Class 2 (Group 1) material by ISO 4049.
2) ESTELITE UNIVERSAL FLOW incorporates Radical-Amplified Photopolymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time. Please see the table depicting the relationship between curing time and increment depth (**■INDICATIONS FOR FILLING AND CURING**).
3) ESTELITE UNIVERSAL FLOW comes in three viscosities, High, Medium and Super Low. Each type contains spherical silica-zirconia filler and composite filler as indicated in below table. The composite filler offers reduced polymerization shrinkage stress due to its round shape. The spherical filler contained in ESTELITE UNIVERSAL FLOW (mean particle size: 200 nm, particle size range: 100 to 300 nm) facilitates excellent gloss retention and wear resistance.

Type	Filler Loading
High	69% by weight (55% by volume)
Medium	71% by weight (57% by volume)
Super Low	70% by weight (56% by volume)

4) ESTELITE UNIVERSAL FLOW contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), Triethylene glycol dimethacrylate (TEGDMA), 1,6-bis(methacryl-ethyloxycarbonylamino)trimethyl hexane (UDMA), Mequinol, Dibutyl hydroxyl toluene and UV absorber.
5) ESTELITE UNIVERSAL FLOW is available in either SYRINGE (L-SYRINGE or SYRINGE (Medium only)) or Pre Loaded Tip (PLT (Medium only)).

■SHADE

Available shade in each type is as below.

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- Opalescent shades (OPA2, OPA3, OPA4) have adequate opacity designed for blocking out the dark shine that comes through the oral cavity (Class III and IV restoration). The opalescent shades also are designed for masking slight stain or reconstructing a highly opaque tooth.
- Please check with your local distributor for available shades.

■INDICATIONS

- Direct anterior and posterior restorations
- Cavity base or liner
- Blocking out cavity undercuts before fabricating indirect restorations
- Repair of porcelain/composite

■CONTRAINDICATIONS

ESTELITE UNIVERSAL FLOW contains methacrylic monomers and UV absorber. DO NOT use ESTELITE UNIVERSAL FLOW for patients allergic to or hypersensitive to methacrylic and related monomers, to UV absorber or to any of the other ingredients.

■PRECAUTIONS

1) DO NOT use ESTELITE UNIVERSAL FLOW for any purpose other than those listed in these instructions. Use ESTELITE UNIVERSAL FLOW only as directed herein.
2) ESTELITE UNIVERSAL FLOW is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
3) DO NOT use ESTELITE UNIVERSAL FLOW if the safety seals are broken or appear to have been tampered with.
4) If ESTELITE UNIVERSAL FLOW causes an allergic reaction or hypersensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl or latex) at all times when handling ESTELITE UNIVERSAL FLOW to avoid the possibility of allergic reactions from methacrylic monomers or UV absorber. Note: Certain substance/materials may penetrate through examination gloves. If ESTELITE UNIVERSAL FLOW comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of ESTELITE UNIVERSAL FLOW with eyes, mucosal membrane, skin and clothing.
 - If ESTELITE UNIVERSAL FLOW comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If ESTELITE UNIVERSAL FLOW comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
 - If ESTELITE UNIVERSAL FLOW comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) ESTELITE UNIVERSAL FLOW should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of ESTELITE UNIVERSAL FLOW, do not leave it unsupervised within the reach of patients and children.
9) Clean the placement instruments and brushes with alcohol after use.
10) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
11) This product is designed to be used at room temperature (18 - 30°C / 62 - 84°F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

■INDICATIONS FOR FILLING AND CURING

ESTELITE UNIVERSAL FLOW is designed to be cured by either a halogen, LED, or Plasma Arc light with a wavelength of 400-500 nm. Be sure to light-cure ESTELITE UNIVERSAL FLOW extra-orally and check the time needed for complete hardening of ESTELITE UNIVERSAL FLOW with your light-curing unit before performing the bonding procedure. The following table summarizes the relationship between curing time and increment depth.

Relationship between curing time and increment depth:

Light type	Intensity (mW/cm²)	Curing time (seconds)	Increment depth (mm) 1)		
			CE, BW, A1, A2, A3, A3.5, A4, B3	A5	OPA2 OPA3 OPA4
Halogen	600	10	1.8	1.6	1.3
		20	2.3	2.0	1.6
	800	10	2.0	1.7	1.5
		20	2.4	2.1	1.7
	1200	5	1.8	1.6	1.3
		10	2.2	2.0	1.6
LED	600	10	1.7	1.5	1.2
		20	2.1	1.8	1.4
	800	10	1.9	1.7	1.3
		20	2.2	2.0	1.6
	1200	10	2.0	1.8	1.5
	2300	3	1.7	1.5	1.2
Plasma Arc	1200	3	1.7	1.5	1.1

1) Increment depth was determined on the basis of test results performed in accordance with "depth of cure" of ISO 4049.

■SPECIAL NOTES FOR THE USE OF PLT

1) PLTs are designed for single patient use only. Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.
2) Dispensers are not provided for ESTELITE UNIVERSAL FLOW. Use a dispenser that fits the PLT of ESTELITE UNIVERSAL FLOW. For proper use and proper disinfection, see manufacturer's instructions.
3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

■STORAGE

1) Store ESTELITE UNIVERSAL FLOW at temperatures between 0 - 25°C (32 - 77°F).
2) AVOID direct exposure to light and heat.
3) DO NOT use ESTELITE UNIVERSAL FLOW after the indicated date of expiration on the syringe or PLT package.

■DISPOSAL

To safely dispose of excess ESTELITE UNIVERSAL FLOW, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

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2. Shade Selection
Prior to isolation of tooth, select the appropriate shade of ESTELITE UNIVERSAL FLOW using shade guide provided for ESTELITE UNIVERSAL FLOW.

3. Isolation
A rubber dam is the preferred method of isolation.

4. Cavity Preparation
Prepare the cavity. In case where no cavity preparation has been made (caries-free cervical defects), clean the tooth surface with a rubber cup and a fluoride-free cleaning paste. Rinse thoroughly with water.
- In the case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid etch for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling agent, followed by a metal primer if metal is exposed on the fractured surfaces according to manufacturer's instructions.


5. Pulp Protection
Calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing ESTELITE UNIVERSAL FLOW. ESTELITE UNIVERSAL FLOW is not indicated for direct pulp capping.

6. Bonding System
Apply the adhesive system according to the manufacturer's instructions.
- If you are using self-curing adhesive systems or dual-curing systems, make sure that they can be used with light-curing composite materials.

7. Dispensing


7-1. PLT (Medium)
Please read **■SPECIAL NOTES FOR THE USE OF PLT** before use.
- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.

7-2. L-SYRINGE (High, Medium and Super Low)
- Hold the nozzle and turn the cap counterclockwise for removal. [1]
- Hold the nozzle and place the Dispensing Tip (metal tip) onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

- Accessories and Dispensing Tip are exclusively designed for the L-SYRINGE of ESTELITE UNIVERSAL FLOW. Use the attached Dispensing Tips (metal tips).
- Dispensing Tips (TOKUYAMA Dispensing Tip) are also sold separately.
- Use uncontaminated gloves when handling the syringe and attaching the Dispensing Tip.

7-3 SYRINGE (Medium)
- Hold the nozzle and turn the cap counterclockwise for removal. [1]
- Hold the nozzle and place the Dispensing Tip (plastic tip) onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

- Accessories and Dispensing Tip are exclusively designed for the SYRINGE of ESTELITE UNIVERSAL FLOW. Use the attached Dispensing Tips (plastic tips).
- Dispensing Tips (TOKUYAMA Plastic Tip) are also sold separately.
- Use uncontaminated gloves when handling the syringe and attaching the Dispensing Tip.

8. Filling and contouring
Cavities can be directly filled from the Dispensing Tip/PLT, or indirectly filled with an instrument after extruding the paste on the mixing pad. Increments should not exceed the indicated curing depth at a time (please refer to the aforementioned table).
- DO NOT apply excessive force to extrude the material from the syringe or PLT immediately after removal from the refrigerator.
- DO NOT mix ESTELITE UNIVERSAL FLOW with other brands of resin composite or other shades of ESTELITE UNIVERSAL FLOW to avoid incomplete cure or entrapment of air bubbles.
- In case of syringe, after extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.
- This product is provided non-sterile and is not intended to be sterilized by the user.

9. Curing
Light-cure each increment for at least the indicated time (please see aforementioned table).
- If other brands of composite resins are layered over the cured composite, follow their instructions.

10. Finishing
Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use.
Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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ESTECEM II *PLUS*

ENGLISH

Read all information, precautions and notes before using.

■PRODUCT DESCRIPTION AND GENERAL INFORMATION

1. ESTECEM II PLUS is a dual-cure (light and/or self-cure), radiopaque, adhesive resin cement system with excellent handling, esthetic and adhesive properties to teeth and all prosthetic materials. The main components of the ESTECEM II PLUS KIT are the PASTE (A/B) and Tokuyama Universal Bond II (A/B). The PASTE is dispensed in an automix syringe. Tokuyama Universal Bond II promotes the adhesion of PASTE to tooth structure. Tokuyama Universal Bond II also enhances the adhesion of PASTE to all prosthetic materials.
2. The PASTE contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate (TEGDMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), Peroxide, Camphorquinone and Silica-Zirconia Filler (Filler Load 74 wt%).
3. Tokuyama Universal Bond II contains Phosphoric acid monomer, Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate (TEGDMA), 2-Hydroxyethyl methacrylate (HEMA), MTU-6 (thiouracil monomer), Silane coupling agent, Peroxide, Borate catalyst, Acetone, Ethanol and Purified water.

■INDICATIONS

- Cementation of crowns, bridges, inlays, and onlays made of glass/oxide ceramics (porcelain, zirconia and alumina), metals/alloys (precious and non-precious) and resin materials including inorganic filler (composite materials)
- Repair of fractured porcelain fused to metal crowns and all ceramic restorations
- Cementation of veneers
- Cementation of adhesion bridges
- Cementation of metal or resin cores, metal or glass-fiber posts

■CONTRAINDICATIONS

ESTECEM II PLUS contains methacrylic monomers, peroxide, organic solvents and acids. DO NOT use ESTECEM II PLUS for patients allergic to or hypersensitive to methacrylic monomers, related monomers, peroxide, organic solvents and acids.

■PRECAUTIONS

1) Use Tokuyama Universal Bond II as primer for the ESTECEM II PLUS PASTE. DO NOT use other primers or adhesives as primer for the adhesive surfaces.
2) DO NOT use ESTECEM II PLUS for any purpose other than those listed in these instructions. Use ESTECEM II PLUS only as directed herein.
3) ESTECEM II PLUS is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
4) DO NOT use the ESTECE II PLUS KIT, PASTE, or Tokuyama Universal Bond II if the safety seals are broken or appear to have been tampered with.
5) If ESTECEM II PLUS causes an allergic reaction or oversensitivity, discontinue its use immediately.
6) Use examination gloves (plastic, vinyl or latex) at all times when handling ESTECEM II PLUS to avoid the possibility of allergic reactions from methacrylic monomers.
 Note: Certain substances/materials may penetrate through examination gloves. If the PASTE, or Tokuyama Universal Bond II comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
7) Avoid contact of the PASTE, or Tokuyama Universal Bond II with eyes, mucosal membrane, skin and clothing.
 - If the PASTE, or Tokuyama Universal Bond II comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If the PASTE, or Tokuyama Universal Bond II comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is seated. Tokuyama Universal Bond II affected areas may whiten from protein coagulation, but such whitening should disappear within 24 hours. If such whitening does not disappear within 24 hours, immediately contact a physician, and the patient should be so advised.
 - If the PASTE or Tokuyama Universal Bond II comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
8) The PASTE or Tokuyama Universal Bond II should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
9) To avoid the unintentional ingestion of the PASTE or Tokuyama Universal Bond II, do not leave it unsupervised within the reach of patients and children.
10) DO NOT expose Tokuyama Universal Bond II or the vapor to open flame because Tokuyama Universal Bond II is flammable.
11) DO NOT use mixing wells supplied by other manufacturers.
12) To avoid cross infection or lowering bond strength, DO NOT reuse the disposable applicator and the disposable mixing well. Clean the mixing well thoroughly with alcohol after each use.
13) DO NOT use ESTECEM II PLUS directly to the cavity being in close proximity to the pulp. Pulp protection with glass ionomer lining or calcium hydroxide is recommended. DO NOT use EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit the curing of Tokuyama Universal Bond II.
14) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
15) Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.
16) DO NOT mix the PASTE or Tokuyama Universal Bond II with other brands of paste, primers or adhesives.
17) The PASTE or Tokuyama Universal Bond II should be brought to room temperature before using. After removing from the refrigerator, allow the product to sit for 20 minutes or until it reaches room temperature. Return the PASTE to the refrigerator after using.

■PRECAUTIONS FOR PASTE

1) Use the Cement Mixing Tip or Mini Mixing Tip S included in this KIT when extruding the mixed PASTE.
2) Use a clean instrument or explorer for the application or removal of the PASTE.
3) To avoid unexpected hardening of the PASTE, after each use remove the mixing-tip from the automix syringe. Using a clean cloth or tissue wipe the excess paste from the tip of syringe prior to recapping the syringe using the cap provided with the syringe at the time of purchase.

■PRECAUTIONS FOR TOKUYAMA UNIVERSAL BOND II

1) To avoid splashing the liquid, the bottles of Tokuyama Universal Bond II should not be opened immediately after taking out of the refrigerator.
2) Use the mixing well or disposable mixing well provided when dispensing and mixing of Tokuyama Universal Bond II. If Tokuyama Universal Bond II is not dispensed into the rubber mixing it may be difficult to remove the cured Tokuyama Universal Bond II from mixing well. The cured Tokuyama Universal Bond II can be removed by inverting the rubber.
3) Some materials and medicaments (hemostatic material) inhibit setting/adhesion of Tokuyama Universal Bond II. DO NOT use products that contain:
 - eugenol,
 - aluminum chloride,
 - ferric sulfate,
 - calcium sulfate,
 - aluminum sulfate,
 - diammine silver fluoride [molecular formula: $Ag(NH_3)_2F$]
4) Tokuyama Universal Bond II does not adhere to tooth structure immediately after applications of materials that contain:
 - hydrogen peroxide (oxydol),
 - sodium hypochlorite

 We advise avoiding the use of these materials on the day of cementation, however these materials can be used for root canal treatment as the inhibition caused by these materials usually evaporates within five days.
5) Be sure to hold the bottles of Tokuyama Universal Bond II A and B vertically when dispensing. Do not hold the bottles diagonally or horizontally while dispensing; this could cause the Tokuyama Universal Bond II flow back into the nozzle and become contaminated.

■STORAGE

1) Store the ESTECEM II PLUS KIT, PASTE under refrigeration at temperatures between 0 to 10℃(32 to 50℉).
2) Store the Tokuyama Universal Bond II at temperatures between 0 to 25 ℃ (32 to 77℉).
3) Keep away from heat, direct sunlight, sparks and open flames.
4) DO NOT use after the indicated expiration date on the syringe/bottle/package.

■DISPOSAL

- To dispose of the PASTE, extrude unused paste from the automix syringe, mix and allow paste to harden before disposal.
- Unused Tokuyama Universal Bond II is to be absorbed into an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

■CLINICAL PROCEDURES

1. Tooth Preparation

Remove contamination from the tooth surface such as plaque, temporary cement residue, oil from materials to test crown fit, oil mist from a hand-piece, saliva, blood and exudates fluids by brushing with flour pumice, ultrasonic scaling or cleaning the surface of the tooth with alcohol. Thoroughly rinse and dry the surface.
- If the adherent surface involves uncut enamel, apply etching agent to the uncut enamel. See section "Selective-enamel-etch technique".
- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.

2. Etching of the tooth structure (optional)

(a.) Selective-enamel-etch technique
Apply phosphoric acid gel only to the enamel surrounding the margin of prepared cavity and leave phosphoric acid gel in place for 10-15 seconds. Rinse the etched surface thoroughly (at least 15 seconds) with water, and then dry with mild air.

(b.) Total etch technique
Apply phosphoric acid gel onto the prepared enamel first, and then the dentin. The etchant should be left to react on the enamel and dentin for 10-15 seconds. Rinse thoroughly (at least 15 seconds) with water and dry with mild air or with cotton pellets; do not overdry.

3. Restoration Preparation

Follow the procedure shown below, depending on the material to be restored.

1) Ceramics other than Porcelain, and Composite materials
 Roughen the interior of the restoration by sandblasting (0.1 to 0.2MPa), air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.
2) Porcelain
 DO NOT roughen the surface; Prepare the interior of the restoration, providing a clean and fresh surface for adhesion. Thoroughly rinse and dry the surface.
3) Metal restorations
 Roughen the area by sandblasting (0.3 to 0.5MPa), air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.
 - When using the sandblaster or air abrasion, use at least 30 to 50 micron grid of alumina powder.
 - If unable to roughen or prepare the adhesive surface, etch the area for adhesion with phosphoric acid then thoroughly rinse and dry.
 - The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol soaking. Clean the restoration after adjustment is complete.
 1) Gypsum and investment
 2) Silicone oil from fit check material
 3) Oil mist from handpiece
 4) Plaque, tartar, saliva, blood and exudate fluids
 5) Temporary restoration and temporary cement
 - If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

4. Dispensing of Tokuyama Universal Bond II

Allow Tokuyama Universal Bond II to reach room temperature prior to use. Open the bottle by removing the cap and dispense one drop each of Tokuyama Universal Bond II A and B into mixing well or disposable mixing well. Mix thoroughly with a disposable applicator until the mixed bonding agent turns green.
- Dispense Tokuyama Universal Bond II B IMMEDIATELY after dispensing A. Tokuyama Universal Bond II A may turn pink due to solvent volatilization if allowed to sit too long. If Bond A turns pink, dispose of it and re-dispense fresh Tokuyama Universal Bond II A and B.
- Complete the application within 1 minute of dispensing when using the mixing well since Tokuyama Universal Bond II contains volatile solvents. After mixing, the color of Tokuyama Universal Bond II changes gradually based on the reaction of the catalyst of Tokuyama Universal Bond II, this is normal and does not indicate a problem with Tokuyama Universal Bond II.
- Complete the application within 3 minutes of dispensing when using the disposable mixing well.
- DO NOT touch the nozzle to the mixing well while dispensing. DO NOT interchange the bottle caps of Tokuyama Universal Bond II A and B. If Tokuyama Universal Bond II A and B come into contact while dispensing, Tokuyama Universal Bond II may cure in the nozzle.
- Wipe off excess Tokuyama Universal Bond II from the tip of the nozzle and the inside of the cap prior to closing.
- Immediately after dispensing, tightly replace the bottle caps of Tokuyama Universal Bond II A and B.
- Do not over tighten the bottle caps.

5. Restoration Pretreatment with the Tokuyama Universal Bond II

1) Application
 Using the disposable applicator, apply the mixed Tokuyama Universal Bond II on the surface to be bonded.
 - Ensure Tokuyama Universal Bond II covers all surfaces where the PASTE will be applied.
 - DO NOT over apply. One application of Tokuyama Universal Bond II is sufficient. Too thick of a layer may compromise the fit of the restoration.
 - DO NOT allow saliva, blood, oil, water etc. to come in contact with the surface to be bonded. If the treated surface is contaminated, thoroughly rinse with water, dry and re-apply fresh Tokuyama Universal Bond II.
 - DO NOT wait too long after the application of Tokuyama Universal Bond II. Air dry within 30 seconds after application, to ensure the appropriate film thickness.
2) Air dry
 Using an oil-free air/water syringe, apply mild air to the surface. Use a vacuum aspirator to prevent spatter of the Tokuyama Universal Bond II.
 - Be sure to dry the surface after Tokuyama Universal Bond II is applied. Incomplete drying may result in poor adhesion.

6. Tooth (Adherent surface) Pretreatment with the Tokuyama Universal Bond II

1) Application
 Using the disposable applicator, apply the mixed Tokuyama Universal Bond II on the surface to be bonded.
 - Ensure Tokuyama Universal Bond II covers all surfaces where the PASTE will be applied. Even if the adherent surface includes materials other than tooth structure such as ceramics, porcelain, composite, and/or metal, no extra primers are needed.
 - Dispensed Tokuyama Universal Bond II can be used for both the restoration and tooth pretreatment if used within 1 minute using the mixing well or 3 minutes if using the disposable mixing well.
 - If the time after dispensing is over 1 minute (mixing well), or over 3 minutes (disposable mixing well), DO NOT use dispensed Tokuyama Universal Bond II due to the possibility of the poor coating properties or insufficient adhesion. If the time after dispensing exceeds the limit mentioned above (1 minute or 3 minutes), dispense again one drop each of Tokuyama Universal Bond II A and B into the mixing well or disposable mixing well and mix thoroughly with a disposable applicator.
 - If the Tokuyama Universal Bond II pools on the cavity floor or a cavosurface angle and is too thick to air thin, remove the excess with a paper point or a new disposable applicator.
 - Be careful not to allow excess Tokuyama Universal Bond II to flow into the subgingival. If Tokuyama Universal Bond II flows into the subgingival, thoroughly rinse and dry the area.
 - DO NOT wait too long after the application of Tokuyama Universal Bond II; this will lead to the possibility of the restoration not seating properly with cementation. Air dry within 30 seconds after application, to ensure the appropriate film thickness.
2) Air dry
 Using an oil-free air/water syringe, apply weak air continuously to the Tokuyama Universal Bond II surface until the runny Tokuyama Universal Bond II stays in the same position without any movement. Finish by using mild air to the surface. Use a vacuum aspirator to prevent spatter of the Tokuyama Universal Bond II.
 - If saliva, blood, or other fluids contaminate the applied Tokuyama Universal Bond II, thoroughly rinse with water, dry and re-apply fresh Tokuyama Universal Bond II.
 - DO NOT rinse the applied Tokuyama Universal Bond II with water except in unintentional contamination.
 - If accidental spattering occurs, it may cause the tissue to whiten or possible allergic reaction.

7. Cementation and Final Bonding

1) Preparing the syringe and mixing tips
 Before each use, remove the cap from automix syringe, and squeeze out a small quantity of paste A and B, making sure equal amounts are being dispensed through the two outlets of the syringe.
 - How to remove the cap: turn the cap 1/4 of a turn counterclockwise to align the projections on the cap with the grooves in the syringe. Holding the base of the cap, remove the cap by twisting and pressing downward.
 - DO NOT allow PASTE A to come into contact with PASTE B, there is a possibility of hardening of the PASTE at the tip of syringe. Wipe the excess PASTE from the tip of syringe with a clean cloth or tissue before using a new mixing tip.
 - If equal amounts of PASTE A and B are not dispensed, there is a possibility of poor polymerization. Ensure equal amounts are being dispensed through the two outlets of the syringe prior to attaching a mixing tip.
 - How to attach the mixing tip: align the projections on the adapter with the grooves in the syringe as closely as possible, and turn the mixing tip 1/4 of a turn clockwise to lock it in place.
 - After use, remove the mixing tip from automix syringe and replace the cap of the automix syringe.
2a) Cementation of Restoration
 Attach a Cement Mixing Tip to the syringe, squeeze PASTE (A and B) from the syringe and apply the mixed PASTE to the inner surface of the restoration. Place the restoration on the tooth with firm pressure. Any excess PASTE at the margin can be removed with a brush tip.
 - DO NOT apply the PASTE directly to the prepared tooth or the cavity treated with Tokuyama Universal Bond II to avoid incorrect seating of the restoration.
 - To keep the appropriate film thickness, seat the restoration as soon as possible after the application of the PASTE (within 1 minute after dispensing), if not there is a possibility that the restoration will not seat properly with cementation.
2b) Placing the core or post
 Attach a Mini Mixing Tip S and a Nozzle for Mini Mixing Tip to the syringe for application in the root canal. Apply the mixed paste over the entire adherent surface of the core or post, or the entire tooth surface within the cavity. Place the core or post quickly into the cavity, slightly vibrating it to prevent air bubbles from entering the root canals.
 - If the paste is applied directly into the cavity, place the core or post into the cavity within 30 seconds after application of the PASTE.
 - DO NOT use Lentulo spirals to insert the cement in the root canal as this can excessively accelerate setting.

8. Removal of Excess PASTE.

There are two methods for removal of excess PASTE.
a) When using the light-cure method
 Once the restoration is seated, light-cure the excess PASTE with a dental curing light for 2 to 4 seconds. The gel state of the excess PASTE can be easily removed.
 - Excess light curing will lead to difficult removal of the excess cement.
 - The Tokuyama Universal Bond II layer may appear fluorescent under the dental curing light irradiation.
b) When using the self-cure method
 The gel state of the excess cement can be easily removed within 1 to 3.5 minutes after seating.
 - Waiting more than 3.5 minutes after seating will lead to difficult removal of the excess cement.
 - The setting time of the PASTE may change (approximately 20 seconds at the maximum) depending on the oral temperature in each clinical case.

9. Final Adjustment

1) Final PASTE hardening
 a) In case of translucent restoration materials such as ceramics or composite materials
 Irradiate the light to the applied PASTE on the seated restoration for 20 seconds or more. When the adherent area is too large to be completed by a single light-cure, multiple cures are recommended.
 - Confirm that the light-curing unit has sufficient intensity (>500mW/cm2) before using.

 Note that using a cracked light guide will lower the intensity.
 b) In case of non translucent restoration materials such as metal
 Light-cure along margins for 20 seconds or more, then allow the PASTE to set for 8 minutes with the patient applying normal occlusal pressure to the restoration.
2) Finishing
 Check patient's bite, if necessary. Re-polish the restoration with material of choice.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use.
Specifications are subject to change without notice. When the product specifications change, the instructions and precautions may also change.

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Protective Sealant for Treatment of Hypersensitive Dentin

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1. TOKUYAMA SHIELD FORCE PLUS is a one-component, single application, light-cured protective sealant designed for the treatment of hypersensitive dentin. TOKUYAMA SHIELD FORCE PLUS effectively infiltrates and seals exposed dentinal tubules to inhibit the flow of fluid through the tubules which is the cause of sensitivity. At the same time the durable coat of cured TOKUYAMA SHIELD FORCE PLUS helps prevents loss of dentin due to abrasion and/or erosion of exposed cervical dentin.
2. A light-curing unit with a wavelength range of 400 to 500 nm (peak: 470 nm) may be used for curing TOKUYAMA SHIELD FORCE PLUS.
3. TOKUYAMA SHIELD FORCE PLUS contains phosphoric acid monomer, bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), triethylene glycol dimethacrylate, 2-hydroxyethyl methacrylate (HEMA), camphorquinone, alcohol and purified water. The pH level immediately after dispensing is approximately 2.0.
4. TOKUYAMA SHIELD FORCE PLUS is dispensed in a bottle.

■INDICATIONS

TOKUYAMA SHIELD FORCE PLUS is indicated for:
- Treatment of hypersensitive dentin,
- Reduction of abrasion and erosion of exposed cervical dentin,
- Alleviation and/or prevention of tooth sensitivity after tooth preparation for direct and indirect restorations.

■CONTRAINDICATIONS

1) TOKUYAMA SHIELD FORCE PLUS contains methacrylic monomers, organic solvents and acids. DO NOT use TOKUYAMA SHIELD FORCE PLUS for patients allergic to or hypersensitive to methacrylic monomers, related monomers, organic solvents and acids.
2) TOKUYAMA SHIELD FORCE PLUS is NOT indicated for:
 - Treatment of subgingival area where the gingival reattachment is expected,
 - Direct pulp capping,
 - Root canal sealing.
3) Not to be used if the patient suffers from severe gingival inflammation (especially conditions which can lead to tissue sloughing, e.g. acute necrotizing ulcerative gingivitis (ANUG)).

■PRECAUTIONS

1) DO NOT use TOKUYAMA SHIELD FORCE PLUS for any purposes other than those listed herein. Use TOKUYAMA SHIELD FORCE PLUS only as directed.
2) TOKUYAMA SHIELD FORCE PLUS is designed for safe and use by licensed dental care professionals only. It is not designed for sale nor is suitable for use by non-dental care professionals.
3) DO NOT use TOKUYAMA SHIELD FORCE PLUS if the safety seals are broken or appear to have been tampered with.
4) If TOKUYAMA SHIELD FORCE PLUS causes an allergic reaction, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl or latex) at all times when handling TOKUYAMA SHIELD FORCE PLUS to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If TOKUYAMA SHIELD FORCE PLUS comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of TOKUYAMA SHIELD FORCE PLUS with eyes, mucous membrane, skin and clothing.
 - If TOKUYAMA SHIELD FORCE PLUS comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If TOKUYAMA SHIELD FORCE PLUS comes in contact with the mucous membrane, wipe the affected area immediately, and thoroughly flush the area with water. Affected areas may be whitened from protein coagulation, but such whitening tissue should disappear within 24 hours. If such whitening does not disappear after 24 hours, immediately contact a physician, and the patient should be so advised.
 - If TOKUYAMA SHIELD FORCE PLUS comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
 - Instruct the patient to rinse his/her mouth immediately after treatment.
7) TOKUYAMA SHIELD FORCE PLUS should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of TOKUYAMA SHIELD FORCE PLUS, do not leave it unsupervised within the reach of patients and children.
9) DO NOT expose TOKUYAMA SHIELD FORCE PLUS or its vapor to open flame because TOKUYAMA SHIELD FORCE PLUS is flammable.
10) To avoid cross infection, DO NOT reuse the disposable applicator included in the TOKUYAMA SHIELD FORCE PLUS package.
11) Clean the dispensing well included in the TOKUYAMA SHIELD FORCE PLUS package thoroughly with alcohol after each use.
12) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.

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1) Some materials and medicaments (hemostatic solutions) inhibit curing of TOKUYAMA SHIELD FORCE PLUS for an extended period even after meticulous cleansing with water. DO NOT USE products which contain:
 - Eugenol
 - Hydrogen peroxide
 - Sodium hypochlorite
 - Diammine silver fluoride [molecular formula: $Ag(NH_3)_2F$]
 - Phenols such as parachlorophenol, 2-methoxyphenol, phenol
 - Aluminum chloride
 - Ferric sulfate
 - Aluminum sulfate
 - Epinephrine
2) Application of phosphoric acid followed by application of TOKUYAMA SHIELD FORCE PLUS may reduce treatment efficacy.

■STORAGE

1) Store under refrigerated conditions at temperatures between 0 to 10°C (32 to 50°F).
2) Keep away from heat, direct sunlight, sparks, and open flames.
3) DO NOT use after the indicated expiration date on the bottle/package.

■DISPOSAL

Unused TOKUYAMA SHIELD FORCE PLUS is to be absorbed into an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

■CLINICAL PROCEDURES

1. Treatment of dentinal hypersensitivity

1.1. Cleaning

Clean tooth surface and remove plaque and saliva with a cotton pellet. A rubber cup with a fluoride-free paste under local anesthesia may also be used to clean.

1.2. Isolation

The rubber dam is the preferred method for isolation.

1.3. Drying

Dry the tooth by using a blotting technique or an air syringe technique.
- DO NOT desiccate the hypersensitive dentin. Desiccation can aggravate the hypersensitivity.
- Glass ionomer lining or calcium hydroxide should be applied if the area to be treated is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of TOKUYAMA SHIELD FORCE PLUS.

1.4. Dispensing

Open the bottle cap and dispense one or two drops of sealant into the Dispensing Well. Close the bottle cap tightly immediately after dispensing.
- Wipe off the excess sealant on the tip of nozzle before closing.
- Do not mix the sealant with other materials.

1.5. Application

Apply generous amounts of the sealant to the area to be treated using the disposable applicator. Leave undisturbed for 10 seconds or more. If the excess sealant pools on the area, blot the excess with a new disposable applicator.
- Protect the dispensed sealant and the inserted applicator from ambient light before the application using a light blocking plate.
- Complete the application within 5 minutes after dispensing because TOKUYAMA SHIELD FORCE PLUS contains a volatile alcohol.
- If saliva, blood, or other fluids contaminate the applied area, thoroughly rinse the area with water, dry and re-apply the sealant.
- Do not rinse the applied sealant with water except in unintentional contamination.

1.6. Air Dry

Using an oil-free air/water syringe, apply weak air flow continuously to the sealant surface until the runny sealant stays in the same position without any movement (usually for 5 seconds). Finish by using a strong air flow for 5 seconds or more. Use a vacuum aspirator to prevent spatter of the sealant.
- If spattering accidentally occurs, it may cause the tissue to whiten or a possible allergic reaction.
- To prevent spattering, refer to the following tips:
 1) To avoid inadvertent strong air;
 a) Be sure to begin with a weak air flow outside the mouth.
 b) Direct the air flow at the applied area.
 2) Extending the distance of the air/water syringe from the tooth reduces the air flow. The mirror reflection technique is also useful in reducing the air flow.

1.7. Light-cure

Keep the curing light tip within a distance of 2 mm from the surface. Then light-cure the surface for 10 seconds or more. If the surface is too large, divide the area into segments and light-cure each segment individually.
- Confirm that the light-curing unit has sufficient intensity (>300 mW/cm²) before using. Note that using a cracked light guide will lower the intensity.

1.8. Check for excess material

Check for excess material in the gingival sulcus, and remove the excess with a scaler or periodontal curette. Instruct the patient to rinse his/her mouth immediately after treatment.

1.9. Test for sensitivity

If a single application of TOKUYAMA SHIELD FORCE PLUS does not sufficiently relieve the hypersensitivity or sustained pain persists, the treatment can be repeated before alternative treatments are indicated. If hypersensitivity persists after a second application, seriously consider other sources of pain, such as traumatic occlusion, pulpal pathology.
- If additional restorations are required due to the esthetic reason, restore according to procedures in step "2.9. Adhesive treatment" and "2.10. Composite resin filling".

1.10. Patient recall

Depending on individual oral hygiene, wear of the TOKUYAMA SHIELD FORCE PLUS may require re-application approximately every six months.

2. TOKUYAMA SHIELD FORCE PLUS sealant under direct restorations to alleviate and/or prevent tooth sensitivity after tooth preparation

2.1. Cleaning

Thoroughly clean the tooth surface using a rubber cup with a fluoride-free paste then rinse with water.

2.2. Isolation

The rubber dam is the preferred method for isolation.

2.3. Cavity preparation

Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in masking demarcations between the cavity margins and the restoration, enhancing both esthetics and retention.

2.4. Drying

Dry the cavity by using a blotting technique or an air syringe technique.
- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- Glass ionomer lining or calcium hydroxide should be applied if the area to be treated is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of TOKUYAMA SHIELD FORCE PLUS.

2.5. Dispensing

Open the bottle cap and dispense one or two drops of sealant into the Dispensing Well. Close the bottle cap tightly immediately after dispensing.
- Wipe off the excess sealant on the tip of nozzle before closing.
- Do not mix the sealant with other materials.

2.6. Application

Apply generous amounts of the sealant to the area to be treated using the disposable applicator. Leave undisturbed for 10 seconds or more. If the excess sealant pools on the area, blot the excess with a new disposable applicator.
- Protect the dispensed sealant and the inserted applicator from ambient light before the application using a light blocking plate.
- Complete the application within 5 minutes after dispensing because TOKUYAMA SHIELD FORCE PLUS contains a volatile alcohol.
- If saliva, blood, or other fluids contaminate the applied area, thoroughly rinse the area with water, dry and re-apply the sealant.
- Do not rinse the applied sealant with water except in unintentional contamination.

2.7. Air Dry

Using an oil-free air/water syringe, apply weak air flow continuously to the sealant surface until the runny sealant stays in the same position without any movement (usually for 5 seconds). Finish by using a strong air flow for 5 seconds or more. Use a vacuum aspirator to prevent spatter of the sealant.
- If spattering accidentally occurs, it may cause the tissue to whiten or a possible allergic reaction.
- To prevent spattering, refer to the following tips:
 1) To avoid inadvertent strong air;
 a) Be sure to begin with a weak air flow outside the mouth.
 b) Direct the air flow at the applied area.
 2) Extending the distance of the air/water syringe from the tooth reduces the air flow. The mirror reflection technique is also useful in reducing the air flow.

2.8. Light-cure

Keep the curing light tip within a distance of 2 mm from the surface. Then light-cure the surface for 10 seconds or more. If the surface is too large, divide the area into segments and light-cure each segment individually.
- Confirm that the light-curing unit has sufficient intensity (>300 mW/cm²) before using. Note that using a cracked light guide will lower the intensity.

2.9. Adhesive treatment

Immediately apply adhesive as described in manufacturer's directions for use.
- If using a total etch adhesive system, DO NOT etch the tooth surface prior to the use of TOKUYAMA SHIELD FORCE PLUS.

2.10. Composite resin filling

Place the restorative material over the cured adhesive according to manufacturer's directions for use.

3. TOKUYAMA SHIELD FORCE PLUS sealant under indirect restorations to alleviate and/or prevent tooth sensitivity after tooth preparation

3.1. Temporary Restoration

3.1.1. Tooth Preparation

Following tooth preparation, core build-up and fabrication of the provisional restoration, clean enamel and dentin with water spray and then air dry. Apply a rubber dam or other suitable isolation technique.
- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- Glass ionomer lining or calcium hydroxide should be applied if the area to be treated is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of TOKUYAMA SHIELD FORCE PLUS.

3.1.2. Dispensing

Open the bottle cap and dispense one or two drops of sealant into the Dispensing Well. Close the bottle cap tightly immediately after dispensing.
- Wipe off the excess sealant on the tip of nozzle before closing.
- Do not mix the sealant with other materials.

3.1.3. Application

Apply generous amounts of the sealant to the area to be treated using the disposable applicator. Leave undisturbed for 10 seconds or more. If the excess sealant pools on the area, blot the excess with a new disposable applicator.
- Protect the dispensed sealant and the inserted applicator from ambient light before the application using a light blocking plate.
- Complete the application within 5 minutes after dispensing because TOKUYAMA SHIELD FORCE PLUS contains a volatile alcohol.
- If saliva, blood, or other fluids contaminate the applied area, thoroughly rinse the area with water, dry and re-apply the sealant.
- Do not rinse the applied sealant with water except in unintentional contamination.

3.1.4. Air Dry

Using an oil-free air/water syringe, apply weak air flow continuously to the sealant surface until the runny sealant stays in the same position without any movement (usually for 5 seconds). Finish by using a strong air flow for 5 seconds or more. Use a vacuum aspirator to prevent spatter of the sealant.
- If spattering accidentally occurs, it may cause the tissue to whiten or a possible allergic reaction.
- To prevent spattering, refer to the following tips:
 1) To avoid inadvertent strong air;
 a) Be sure to begin with a weak air flow outside the mouth.
 b) Direct the air flow at the applied area.
 2) Extending the distance of the air/water syringe from the tooth reduces the air flow. The mirror reflection technique is also useful in reducing the air flow.

3.1.5. Light-cure

Keep the curing light tip within a distance of 2 mm from the surface. Then light-cure the surface for 10 seconds or more. If the surface is too large, divide the area into segments and light-cure each segment individually.
- Confirm that the light-curing unit has sufficient intensity (>300 mW/cm²) before using. Note that using a cracked light guide will lower the intensity.

3.1.6. Dismissal

Wipe the oxygen inhibition layer with alcohol-soaked cotton swab or gauze.

3.1.7. Impression Taking

Take impression with impression material according to manufacturer's instructions.

3.1.8. Temporary Cementing

Seat the temporary/provisional restoration on the tooth with temporary cement (water-based cement) according to manufacturer's instructions. Temporary stopping (gutta-percha) may also be used for temporary filling.
- Do not use resin-based material (cement/filling) because it may adhere to the cured TOKUYAMA SHIELD FORCE PLUS.

3.2. Permanent Restoration

3.2.1. Cleaning of Cavity/Abutment

Remove temporary restoration, residual temporary cement and hard deposits from the tooth by ultrasonic scaler or other cleaning instruments. Thoroughly clean all contaminants from the tooth surface to be bonded. Rinse with water and dry in the usual manner.
- Elements listed below, which inhibit the setting of resin cement, should be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol, citric acid, EDTA, or the application of phosphoric acid for 2 to 3 seconds before cementation.
 1) Silicone oil from fit check material
 2) Oil mist from handpiece
 3) Saliva, blood and exudates

3.2.2. Cementing

Cement the restoration with resin cement according to its manufacturer's instructions.
- The use of multistep resin cements are recommended for cementing. The use of self-adhesive resin cements are not recommended due to potential lowered adhesion of the cement.

The manufacturer of the TOKUYAMA SHIELD FORCE PLUS is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for an appropriate application before use. Product specifications of the TOKUYAMA SHIELD FORCE PLUS are subject to change without notification. When the product specifications are changed, the instructions and precautions may also be changed.

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Self-Cured Dental Universal Adhesive

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1. Tokuyama Universal Bond II is a two-component self-cured dental adhesive system for both direct and indirect restorations that can be used with Self-etch, Selective-enamel-etch and Total-etch techniques. As a universal adhesive, Tokuyama Universal Bond II has been designed to be fully compatible with light-cured, self-cured and dual-cured composite materials.
2. Tokuyama Universal Bond II improves the bond strength of polymerizable resin material (adhesive resin cement, acrylic resin and composite resin) to indirect restorative materials such as glass-ceramics (porcelain), oxide-ceramics (zirconia and alumina), metals (precious and non-precious) and resin materials including inorganic filler.
3. Tokuyama Universal Bond II contains Phosphoric acid monomer, Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate (TEGDMA), 2-Hydroxyethyl methacrylate (HEMA), Thiouracil monomer (MTU-6), Silane coupling agent, Peroxide, Borate catalyst, Acetone, Ethanol and Purified water.

■INDICATIONS

- Direct anterior and posterior restorations with light-curing, dual-curing, and self-curing composite materials
- Intraoral repair of composite restorations, porcelain fused to metal, metal, and all-ceramic restorations without an additional primer
- Cementation of indirect restorations and veneers when combined with light-cure, dual-cure, and self-curing resin cements
- Bonding of core build-up materials made of core build-up materials
- Bonding of denture resin to metal base, clasp or attachment
- Repair of denture with metal base, clasp or attachment
- Bonding of opaque resin to a metal base in the fabrication of resin-faced crowns

■CONTRAINDICATIONS

Tokuyama Universal Bond II contains methacrylic monomers, peroxide, organic solvents and acids. DO NOT use Tokuyama Universal Bond II for patients allergic to or hypersensitive to methacrylic monomers, related monomers, peroxide, organic solvents and acids.

■PRECAUTIONS

1) DO NOT use Tokuyama Universal Bond II for any purpose other than those listed in these instructions. Use Tokuyama Universal Bond II only as directed herein.
2) Tokuyama Universal Bond II is designed for safe and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
3) DO NOT use Tokuyama Universal Bond II if the safety seals are broken or appear to have been tampered with.
4) If Tokuyama Universal Bond II causes an allergic reaction or oversensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl or latex) at all times when handling Tokuyama Universal Bond II to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If Tokuyama Universal Bond II comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of Tokuyama Universal Bond II with eyes, mucosal membrane, skin and clothing.
 - If Tokuyama Universal Bond II comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If Tokuyama Universal Bond II comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the completing the restoration. Affected areas may whiten from protein coagulation, however, whitening should disappear within 24 hours. If whitening does not disappear within 24 hours, immediately contact a physician, and the patient should be so advised.
 - If Tokuyama Universal Bond II comes into contact with the skin or clothing, immediately saturate the area with an alcohol-soaked cotton swab or gauze and flush with water.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) Tokuyama Universal Bond II should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of Tokuyama Universal Bond II, do not leave it unsupervised within the reach of patients and children.
9) DO NOT expose Tokuyama Universal Bond II or its vapor to open flame because Tokuyama Universal Bond II is flammable.
10) To avoid cross-infection or lowering bond strengths, DO NOT reuse the disposable applicator and the disposable mixing well. Clean the mixing well thoroughly with alcohol after each use.
11) DO NOT use Tokuyama Universal Bond II directly to the cavity base in close proximity to the pulp. Pulp protection with glass ionomer lining or calcium hydroxide is recommended. DO NOT use EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit curing of Tokuyama Universal Bond II.
12) Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.
13) DO NOT use mixing wells supplied by other manufacturers.
14) Use the mixing well or disposable mixing well provided when dispensing and mixing of Tokuyama Universal Bond II. If Tokuyama Universal Bond II is not dispensed into the rubber, it may be difficult to remove the cured Tokuyama Universal Bond II from mixing well. The cured Tokuyama Universal Bond II can be removed by inverting the rubber.
15) DO NOT mix Tokuyama Universal Bond II with other brands of primers or adhesives.
16) To avoid splashing the liquid, the bottle of Tokuyama Universal Bond II should not be opened immediately after taking out of the refrigerator. Tokuyama Universal Bond II should be brought to room temperature before using. After removing from the refrigerator, allow the product to sit for 20 minutes or until it reaches room temperature.
17) Do not use Tokuyama Universal Bond II with resin materials for which special adhesives or primers are specified.

■PRECAUTIONS FOR MEDICAMENTS AND MATERIALS

1) Some materials and medicaments (hemostatic material) inhibit setting/adhesion of Tokuyama Universal Bond II.
 DO NOT USE products that contain:
 - eugenol,
 - aluminum chloride,
 - ferric sulfate,
 - calcium sulfate,
 - aluminum sulfate,
 - diammine silver fluoride [molecular formula: $Ag(NH_3)_2F$]
2) Tokuyama Universal Bond II does not adhere to tooth structure immediately after applications of materials that contain:
 - hydrogen peroxide (oxydol),
 - sodium hypochlorite

We advise avoiding the use of these materials on the day of cementation, however these materials can be used for root canal treatment as the inhibition caused by these materials usually evaporates within five days.

■STORAGE

1) Store Tokuyama Universal Bond II at temperatures between 0 to 25°C (32 to 77°F).
2) Keep away from heat, direct sunlight, sparks and open flames.
3) DO NOT use after the indicated expiration date on the bottle/package.

■DISPOSAL

Unused Tokuyama Universal Bond II is to be absorbed into an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

■CLINICAL APPLICATION FOR DIRECT RESTORATIONS AND INTRAORAL REPAIR OF RESTORATIONS WITH COMPOSITE RESIN

1. Cleaning
Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

2. Isolation
The rubber dam is the preferred method of isolation.

3. Cavity preparation
Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II) because bevels and chamfers assist in ensuring demarcations between the cavity margins and the restoration, thereby enhancing both esthetics and retention.
- If the adhesive technique involves uncut enamel, apply etching agent to the uncut enamel. See section "Selective-enamel-etch technique".

4. Drying
Dry the cavity by using a blotting technique or an air syringe technique.
- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- The substances listed below, which inhibit curing of Tokuyama Universal Bond II, should be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol, citric acid, or the application of phosphoric acid for 2 to 3 seconds before application:
 1) Oil mist from handpiece,
 2) Saliva, blood and exudates

5. Etching of the tooth structure (optional)
(a.) Selective-enamel-etch technique
Apply phosphoric acid gel only to the enamel surrounding the margin of prepared cavity and leave phosphoric acid gel in place for 10-15 seconds. Rinse the etched surface thoroughly (at least 15 seconds) with water, and then dry with mild air.
(b.) Total etch technique
Apply phosphoric acid gel onto the prepared enamel first, and then the dentin. The etchant should be left to react on the enamel and dentin for 10-15 seconds. Rinse thoroughly (at least 15 seconds) with water and dry with mild air or with cotton pellets; do not overdry.

6. Preparing the adhesive surface during repair
6.1 Ceramics other than Porcelain, and Composite materials
Bevel the margin if necessary. Roughen the area using air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.
6.2 Metal restorations
Roughen the area using air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.
6.3 Porcelain
Bevel the margin if necessary. DO NOT roughen the surface. Polish the area to prepare a clean and fresh surface for adhesion. Thoroughly rinse and dry the surface.
- When using air abrasion, use at least 30 to 50 micron grid of alumina powder.
- The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol soaking.
 1) Oil mist from handpiece
 2) Plaque, tartar, saliva, blood and exudate fluids
 - If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

7. Dispensing of Tokuyama Universal Bond II
Allow Tokuyama Universal Bond II to reach room temperature prior to use. Open the bottle by removing the cap and dispense one drop each of Tokuyama Universal Bond II A and B into the mixing well or disposable mixing well. Mix thoroughly with a disposable applicator until the mixed bonding agent turns green.
- Dispense Tokuyama Universal Bond II A and B IMMEDIATELY after dispensing. A Tokuyama Universal Bond II A may turn pink due to solvent volatilization if allowed to sit too long. If Bond A turns pink, dispose of it and re-dispense fresh Tokuyama Universal Bond II A and B.
- Complete the application within 1 minute of dispensing when using the mixing well since Tokuyama Universal Bond II contains volatile solvents. After mixing, the color of Tokuyama Universal Bond II changes gradually based on the reaction of the catalyst of Tokuyama Universal Bond II, this is normal and does not indicate a problem with Tokuyama Universal Bond II.
- Complete the application within 3 minutes of dispensing when using disposable mixing well.
- Be sure to hold the bottles of Tokuyama Universal Bond II A and B vertically when dispensing. Do not hold the bottles diagonally or horizontally while dispensing; this could cause the Tokuyama Universal Bond II flow back into the nozzle and become contaminated.
- DO NOT touch the nozzle to the mixing well while dispensing. DO NOT interchange the bottle caps of Tokuyama Universal Bond II A and B. If Tokuyama Universal Bond II A and B come into contact while dispensing, Tokuyama Universal Bond II may cure in the nozzle.
- Wipe off excess Tokuyama Universal Bond II from the tip of the nozzle and the inside of the cap prior to closing.
- Immediately after dispensing, tightly replace the bottle caps of Tokuyama Universal Bond II A and B.
- Do not over tighten the bottle caps.

8. Application of Tokuyama Universal Bond II
Using the disposable applicator, apply the mixed Tokuyama Universal Bond II to the entire adhered surface(s). Take care not to omit any surfaces where the Tokuyama Universal Bond II should be applied.
- In the case of intraoral repair of restorations, if the adherent surface includes tooth structure and restorations (ceramics, porcelain, metals or composite materials): apply the mixed Tokuyama Universal Bond II to the entire adherent surface at once. No extra primers are needed.
- If the applied Tokuyama Universal Bond II pools on the cavity floor or a concave surface angle and is too thick to air thin, remove the excess with a paper point or a new disposable applicator.
- Be careful not to allow excess Tokuyama Universal Bond II to flow into the subgingival. If Tokuyama Universal Bond II flows into the subgingival, thoroughly rinse and dry the area.
- If saliva, blood, or other fluids contaminate the applied Tokuyama Universal Bond II, thoroughly rinse the cavity with water, dry and re-apply fresh Tokuyama Universal Bond II.
- DO NOT rinse the applied Tokuyama Universal Bond II with water except in unintentional contamination.
- DO NOT blow the applied Tokuyama Universal Bond II with air because this will result in a thicker film.

9. Air Dry
Using an oil-free air/water syringe, apply weak air continuously to the Tokuyama Universal Bond II surface until the shiny Tokuyama Universal Bond II stays in the same position without any movement. Finish by using mild air to the surface. Use a vacuum aspirator to prevent splatter of the Tokuyama Universal Bond II.
- If saliva, blood, or other fluids contaminate the applied Tokuyama Universal Bond II, thoroughly rinse with water, dry and re-apply fresh Tokuyama Universal Bond II.
- Do not rinse the applied Tokuyama Universal Bond II with water except in unintentional contamination.
- If accidental spattering occurs, it may cause the tissue to whiten or possible allergic reaction.

10. Composite resin
Restore with a light-cured, self-cured or dual-cured composite resin according to its manufacturer's instructions. The overfilled composite resin should be thoroughly finished and polished.
- The Tokuyama Universal Bond II layer may appear fluorescent under the dental curing light irradiation.

■CLINICAL APPLICATION FOR CEMENTATION OF INDIRECT RESTORATIONS

1. Preparing the adhesive surface of the restoration
1.1 Ceramics other than Porcelain, and Composite materials
Roughen the interior of the restoration by sandblasting (0.1 to 0.2MPa), air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.
1.2 Metal restorations
Roughen the area by sandblasting (0.3 to 0.5MPa), air abrasion or grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.
1.3 Porcelain
DO NOT roughen the surface. Polish the area to prepare a clean and fresh surface for adhesion. Thoroughly rinse and dry the surface.
- When using the sandblaster or air abrasion, use at least 30 or 50 micron grid of alumina powder.
- If unable to roughen or polish the adhesive surface, etch the area for adhesion with phosphoric acid then thoroughly rinse and dry.
- The substances listed below inhibit adhesion and must be removed from the adhesive surface by mechanical cleaning or alcohol soaking. Clean the restoration after adjustment is complete:
 1) Gypsum and investment
 2) Silicone oil from fit check material
 3) Oil mist from handpiece
 4) Plaque, tartar, saliva, blood and exudate fluids
 5) Temporary restoration and temporary cement
 If necessary, protect the surface not being bonded with metal separator, paraffin film, Vaseline, etc.

2. Dispensing of Tokuyama Universal Bond II
Allow Tokuyama Universal Bond II to reach room temperature prior to use. Open the bottle by removing the cap and dispense one drop each of Tokuyama Universal Bond II A and B into the mixing well or disposable mixing well. Mix thoroughly with a disposable applicator until the mixed bonding agent turns green.
- Dispense Tokuyama Universal Bond II A and B IMMEDIATELY after dispensing. A Tokuyama Universal Bond II A may turn pink due to solvent volatilization if allowed to sit too long. If Bond A turns pink, dispose of it and re-dispense fresh Tokuyama Universal Bond II A and B.
- Complete the application within 1 minute of dispensing when using the mixing well since Tokuyama Universal Bond II contains volatile solvents. After mixing, the color of Tokuyama Universal Bond II changes gradually based on the reaction of the catalyst of Tokuyama Universal Bond II, this is normal and does not indicate a problem with Tokuyama Universal Bond II.
- Complete the application within 3 minutes of dispensing when using disposable mixing well.
- Be sure to hold the bottles of Tokuyama Universal Bond II A and B vertically when dispensing. Do not hold the bottles diagonally or horizontally while dispensing; this could cause the Tokuyama Universal Bond II flow back into the nozzle and become contaminated.
- DO NOT touch the nozzle to the mixing well while dispensing. DO NOT interchange the bottle caps of Tokuyama Universal Bond II A and B. If Tokuyama Universal Bond II A and B come into contact while dispensing, Tokuyama Universal Bond II may cure in the nozzle.
- Wipe off excess Tokuyama Universal Bond II from the tip of the nozzle and the inside of the cap prior to closing.
- Immediately after dispensing, tightly replace the bottle caps of Tokuyama Universal Bond II A and B.
- Do not over tighten the bottle caps.

3. Application of Tokuyama Universal Bond II
Using the disposable applicator, apply the mixed Tokuyama Universal Bond II on the surface to be bonded.
- Ensure Tokuyama Universal Bond II covers all surfaces where the cement will be applied.
- DO NOT over apply. One application of Tokuyama Universal Bond II is sufficient. Too thick of a layer may compromise the fit of the restoration.
- DO NOT allow saliva, blood, oil, water, etc. to come in contact with the surface to be bonded. If the treated surface is contaminated, thoroughly rinse with water, dry and re-apply fresh Tokuyama Universal Bond II.
- DO NOT wait too long after the application of Tokuyama Universal Bond II. Air dry within 30 seconds after application to ensure the appropriate film thickness.

4. Air dry
Using an oil-free air/water syringe, apply mild air to the surface. Use a vacuum aspirator to prevent spatter of the Tokuyama Universal Bond II.

5. Cementing
Cement with a resin cement according to its manufacturer's instructions.
- DO NOT apply resin cement paste directly to the prepared tooth or the cavity treated with Tokuyama Universal Bond II to avoid improper seating of the restoration.
- Curing may be accelerated depending upon the resin cement used. To keep the appropriate film thickness, seat the restoration as soon as possible after the application of the resin cement (within 1 minute after dispensing), if not there is a possibility that the restoration will not seat properly with concentration.
- The Tokuyama Universal Bond II layer may appear fluorescent under the dental curing light irradiation.

■CLINICAL APPLICATION FOR BONDING OF ACRYLIC RESIN TO METAL BASE, CLASP OR ATTACHMENT

1. Preparing the adhesive surface
Roughen the area by sandblasting (0.1 to 0.2MPa), or by grinding with a diamond bur using a slow speed to prepare the surface for adhesion. Thoroughly rinse and dry the surface.

2. Dispensing, Application and Air dry
Use the same procedures as described in CLINICAL APPLICATION FOR CEMENTATOIN OF INDIRECT RESTORATIONS.

3. Application, curing, trimming and finishing of acrylic resin
Proceed with the procedure according to the manufacturer's instructions for the use of polymerizable acrylic resin.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the user to ensure that this product is suitable for its appropriate application before use.
Specifications are subject to change without notice. When the product specifications change, the instructions and precautions may also change.

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One Component Self-Etching Light-Cured Dental Adhesive

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1. TOKUYAMA BOND FORCE II is a one-component, one-coat application, self-etching, light-cured, dental adhesive system with excellent handling and adhesive properties and exceptional marginal integrity toward cut/uncut enamel and dentin when used in combination with light- or dual-cured composite materials.
2. A light-curing unit having the camphorquinone (CQ) wavelength range (peak: 470nm, spectrum: 400 to 500nm) can be used for curing TOKUYAMA BOND FORCE II.
3. TOKUYAMA BOND FORCE II contains Phosphoric acid monomer, Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate, 2-Hydroxyethyl methacrylate (HEMA), Camphorquinone, alcohol and purified water. The pH level immediately after dispensing is approximately 2.8.
4. TOKUYAMA BOND FORCE II is available in three convenient dispensing systems, Bottle, Pen-Type and Unit Dose. Please check with your local distributor for availability.

■INDICATIONS

Bonding of light- or dual-cured composite material to:
- cut/uncut enamel,
- cut/uncut dentin,
- fractured porcelain/composite repair.

■CONTRAINDICATIONS

TOKUYAMA BOND FORCE II contains methacrylic monomers, organic solvents and acids. DO NOT use TOKUYAMA BOND FORCE II for patients allergic to or hypersensitive to methacrylic monomers, related monomers, organic solvents and acids.

■PRECAUTIONS

1) DO NOT use TOKUYAMA BOND FORCE II for any purpose other than those listed in these instructions. Use TOKUYAMA BOND FORCE II only as directed herein.
2) TOKUYAMA BOND FORCE II is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor is it suitable for use by non-dental care professionals.
3) DO NOT use TOKUYAMA BOND FORCE II if the safety seals are broken or appear to have been tampered with.
4) If TOKUYAMA BOND FORCE II causes an allergic reaction or oversensitivity, discontinue its use immediately.
5) Use examination gloves (plastic, vinyl or latex) at all times when handling TOKUYAMA BOND FORCE II to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If TOKUYAMA BOND FORCE II comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
6) Avoid contact of TOKUYAMA BOND FORCE II with eyes, mucosal membrane, skin and clothing.
 - If TOKUYAMA BOND FORCE II comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
 - If TOKUYAMA BOND FORCE II comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed. Affected areas may whiten from protein coagulation, but such whitening should disappear within 24 hours. If such whitening does not disappear within 24 hours, immediately contact a physician, and the patient should be so advised.
 - If TOKUYAMA BOND FORCE II comes into contact with the skin or clothing, immediately saturate the area with alcohol-soaked cotton swab or gauze.
 - Instruct the patient to rinse his mouth immediately after treatment.
7) TOKUYAMA BOND FORCE II should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
8) To avoid the unintentional ingestion of TOKUYAMA BOND FORCE II, do not leave it unsupervised within the reach of patients and children.
9) DO NOT expose TOKUYAMA BOND FORCE II or the vapor to open flame because TOKUYAMA BOND FORCE II is flammable.
10) To avoid cross infection or lowering bond strength, DO NOT reuse the disposable applicator. Clean the dispensing well thoroughly with alcohol after each use.
11) DO NOT use TOKUYAMA BOND FORCE II directly to the cavity being in close proximity to the pulp. Pulp protection with glass ionomer lining or calcium hydroxide is recommended.
12) When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
13) Be aware that when a restoration is chipped from misaligned occlusion or bruxism (clenching, grinding or tapping), the repaired restoration may chip again.
14) DO NOT mix TOKUYAMA BOND FORCE II with other brands of primers or adhesives.
15) TOKUYAMA BOND FORCE II should be brought to room temperature before using if it is stored under refrigeration.
16) If using the Pen dispensing system,
 - DO NOT push the button when the cap is closed.
 - DO NOT push the button when the nozzle is facing up.
 - DO NOT hold the Pen diagonally or horizontally while dispensing.
 - DO NOT use the Pen directly in the oral cavity.
 - DO NOT dismantle the Pen.

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1) Some materials and medicaments (hemostatic material) inhibit the adhesion of TOKUYAMA BOND FORCE II for an extended period even after careful meticulous cleansing with water. DO NOT USE products which contain:
 - eugenol,
 - diammine silver fluoride [molecular formula: $Ag(NH_3)_2F$],
 - phenols such as parachlorophenol, guaiacol, phenol,
 - aluminum chloride,
 - ferric sulfate,
 - aluminum sulfate,
 - epinephrine.
2) TOKUYAMA BOND FORCE II does not adhere to tooth structure immediately after applications of materials which contain:
 - hydrogen peroxide (oxydol),
 - sodium hypochlorite.

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■STORAGE

1) Store TOKUYAMA BOND FORCE II at a temperature between 0 to 25 °C (32 to 77 °F). Storage under inappropriate condition will accelerate the degradation of TOKUYAMA BOND FORCE II.
2) If not in use for a prolonged period of time, store in a cool dark place.
3) Keep away from heat, direct sunlight, sparks and open flames.
4) DO NOT use after the indicated expiration date on the bottle/package.

■DISPOSAL

Unused TOKUYAMA BOND FORCE II is to be absorbed into an inert absorbent material such as gauze or cotton, and disposed of in accordance with local regulations.

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3. Cavity preparation
Prepare the cavity and rinse with water. Add bevels to the enamel margins of anterior preparations (class III, IV, V), as well as chamfers to the margins of posterior preparations (class I, II). Bevels and chamfers assist in erasing demarcations between the cavity margins and the restoration, enhancing both esthetics and retention.
- If there is an adherent surface including uncut enamel, apply etching agent to the uncut enamel.
- In case of porcelain/composite repairs, roughen the surface with a bur or a diamond point to prepare the area for adhesion; apply phosphoric acid for cleaning; rinse thoroughly with water; air dry thoroughly and treat with a silane coupling reagent according to its manufacturer's instructions.

4. Drying
Dry the cavity by blotting or with an air syringe.
- DO NOT desiccate the vital tooth. Desiccation can lead to post-operative sensitivity.
- The substances listed below inhibit curing and should be removed from the tooth surface by thoroughly cleaning the tooth surface with alcohol, citric acid, or the application of phosphoric acid for 2 to 3 seconds before application:
 1) Oil mist from handpiece,
 2) Saliva, blood and exudates.


5. Pulp Protection
Glass ionomer lining or calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT USE EUGENOL BASED MATERIALS to protect the pulp as these materials will inhibit the curing of TOKUYAMA BOND FORCE II.

6. Dispensing

6-1. Unit Dose
Twist open the Unit Dose container and dip the disposable applicator.
- When opening the cap, make sure to verify the top and bottom of the container.
- To avoid cross infection, do not reuse the material on other patients.

6-2. Bottle
Open the Bottle by removing the cap and dispense one or two drops of adhesive into the dispensing well. Immediately after dispensing, tightly replace the Bottle cap.
- Be sure to hold the Bottle of TOKUYAMA BOND FORCE II vertically when dispensing. Do not hold the Bottle diagonally or horizontally while dispensing; this could cause the adhesive to flow back into the nozzle and become contaminated.
- Wipe off excess adhesive from the tip of the nozzle prior to closing.

6-3. Pen
Open the Pen by removing the cap when the cap is facing up and dispense one or two drops of adhesive into the dispensing well by holding the Pen vertically and pushing on button repeatedly until adhesive starts to flow (3 or 2 pushes, * This depends on how the button is pushed and the position that the Pen is being held). Immediately after dispensing, tightly replace the Pen cap.
- Be sure to hold the Pen vertically when dispensing. Do not hold the Pen diagonally or horizontally while dispensing; this could cause the adhesive to flow back into the nozzle and become contaminated.
- Wipe off excess adhesive from the tip of the nozzle prior to closing.

Indicator:
The indicator allows users to monitor the remaining amount of adhesive. Black Line moves closer to the E as the amount of material decreases. When the line nears "E", this indicates that the Pen is almost empty. This is not an exact reading but only an indication of the remaining material.

7. Application
Using the disposable applicator, apply TOKUYAMA BOND FORCE II on the surface to be bonded and wait 10 seconds or more.
- Ensure TOKUYAMA BOND FORCE II covers all surfaces where the composite resin will be applied.
- Protect the dispensed adhesive and the inserted applicator from ambient light before the application.
 Unit Dose: keep the applicator in the Unit Dose container.
 Bottle, Pen: use a light blocking plate.
- Complete the application within 5 minutes after dispensing because TOKUYAMA BOND FORCE II contains a volatile alcohol.
- Be careful not to allow excess TOKUYAMA BOND FORCE II to flow into the subgingival. If TOKUYAMA BOND FORCE II flows into the subgingival, thoroughly rinse and dry the area.
- If the adherent surface include ceramics, porcelain or composite, prepare the surface according to the appropriate procedures prior to the application of TOKUYAMA BOND FORCE II.

8. Air Dry
Using an oil-free air/water syringe, apply mild air continuously to the TOKUYAMA BOND FORCE II surface until the runny TOKUYAMA BOND FORCE II stays in the same position without any movement (usually for 5 seconds). Use a vacuum aspirator to prevent spatter of the adhesive.
- If accidental spattering occurs, it may cause the tissue to whiten or possible allergic reaction.
- If the adhesive pools on the cavity floor or the cavosurface angle and is too thick to air thin, blot the excess with a new disposable applicator before mild air application.
- If saliva, blood, or other fluids contaminate the applied adhesive, thoroughly rinse the cavity with water, dry and re-apply new adhesive.
- Do not rinse the applied adhesive with water except in unintentional contamination.

9. Light-cure
Light-cure the surface for 10 seconds or more, keeping the curing light tip within a distance of 2 mm from the surface. If the cavity is too large or the light tip is too distant (e.g. MOD), divide the area into segments and light-cure each segment individually.
- Confirm that the light-curing unit has sufficient intensity (>300mW/cm²) before using. Note that using a cracked light guide will lower the intensity.

10. Light-cured Composite
Restore with light-cured composite resin according to its manufacturer's instructions. The composite resin must be finished and polished.
- When dual-cured composite resins are placed into a cavity, the first increment must be light-cured using a layering technique.
- Do not use self-cured composite resins here because the Phosphoric acid monomer contained in TOKUYAMA BOND FORCE II may interfere with the curing of the self-cured resins, leading to their premature detachment.

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Fig. 10

ENGLISH

ESTELITE ASTERIA Custom Shade Guide is a shade guide kit for making your own custom shade guides. ESTELITE ASTERIA Custom Shade Guide consists of the following items.

Item	Material	Quantity
A. Holder	ABS (Acrylonitrile-Butadiene-Styrene) resin	1
B. Shade Guide Handle	Polypropylene	15
C. Mold	ABS (Acrylonitrile-Butadiene-Styrene) resin	1

■HOW TO USE

Fig. 1. Select the Mold.
The Mold has two types of tooth shapes. The left side molds are wide and shallow. These are suitable for Enamel shades. The right side molds are narrow and deep. These are suitable for Body shades.
Fig. 2. Flip up the cover cap of the Mold.
Fig. 3. Before setting the Shade Guide Handle into the appropriate slot, fill the root portion* of tooth shaped mold with ESTELITE ASTERIA resin paste using an instrument to anchor the Handle to the resin.
Fig. 4. Insert the Handle with its flat side up until it seats in completely and is flush with the Mold surface.
Fig. 5. Fill up the remaining portion of the tooth Mold with ESTELITE ASTERIA resin paste.
Fig. 6. Remove excess paste and flatten the surface.
Fig. 7. Close the cover cap and press down firmly to squash the resin paste.
Fig. 8. Light irradiate for 10 sec. (A4B shade: 20 sec.) through the cover cap.
Fig. 9. Re-open the cover cap and remove the Handle from the Mold. Trim side resin flashes, as needed.
Fig. 10. Set the finished Handle (Individual Shade) in the Holder.

■NOTICE FOR USE OF THE ESTELITE ASTERIA CUSTOM SHADE GUIDE

1. Keep the ESTELITE ASTERIA Custom Shade Guide away from direct sunlight to prevent a change in color.
2. Refer to the following table when cleaning the ESTELITE ASTERIA Custom Shade Guide.

	Cleaning method
Correct	- wash with soap and water - ultrasonic cleaning while immersed in water - immerse in a 3% glutaraldehyde solution - wipe with a cotton swab soaked in an alcohol (70%) solution
Incorrect	- autoclave [1] - dry heat sterilization [1] - immerse in an alcohol (70%) solution [2] - immerse in a sodium hypochlorite solution [3] - immerse in a hydrogen peroxide solution [3]

1) Shade Guide Handle, Holder and Mold will melt.
2) Holder and Mold will whiten and crack.
3) Individual Shade may discolor.

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