Bentley

BeGraft
peripheral

Peripheral Stent Graft System

en	Peripheral Stent Graft System
bg	Система с периферна присадна протеза
cs	Systém periferního stentgraftu
de	Peripheres Stentgraftsystem
da	Perifert stentgraftsystem
es	Sistema de endoprótesis periférica recubierta
et	Perifeerne stentsiiriku süsteem
fi	Perifeerinen stenttisiirrejärjestelmä
fr	Système d'endoprothèse périphérique couverte
el	Σύστημα περιφερικού μοσχεύματος στεντ
hr	Periferni stent graft sustav
hu	Perifériás sztentgraft rendszer
it	Sistema di endoprotesi periferica
lt	Periferinio stento įvedimo sistema
lv	Perifēro stentu implantātu sistēma
nl	Perifeer stentgraftsysteem
no	Perifert stentgraftsystem
pl	System stentgraftu obwodowego
pt	Sistema de endoprótese periférica coberta
ro	Sistemul de endoproteză periferică
ru	Стент-графт периферийный с системой доставки
sk	Systém periférneho stentgraftu
sl	Sistem perifernih žilnih opornic
sv	Periferstentgraftsystem
tr	Periferal Stent Greft Sistemi

PRINT DATE 02.07.2020. Printed documents must be compared with the valid revision before use!
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

en Caution **bg** Внимание **cs** Pozor **de** Achtung **da** Forsigtig **es** Precaución **et** Ettevaatust! **fi** Huomio **fr** Mise en garde **el** Προσοχή **hr** Oprez **hu** Figyelem! **it** Attenzione **lt** Dėmesio! **lv** Uzmanību! **nl** Pas op **no** Forsiktig **pl** Przestroga **pt** Atenção **ro** Atenție **ru** Внимание! **sk** Pozor **sl** Pozor **sv** Försiktig **tr** Dikkat

en Consult instructions for use **bg** Направете справка с инструкциите за употреба **cs** Nahlédněte do návodu k použití **de** Gebrauchsanweisung beachten **da** Se brugsanvisningen **es** Consulte las instrucciones de uso **et** Vt kasutusjuhiseid **fi** Lue käyttöohjeet **fr** Consulter le mode d'emploi **el** Συμβουλευτείτε τις οδηγίες χρήσης **hr** Pogledajte upute za uporabu **hu** Lásd a használati útmutatót **it** Consultare le istruzioni per l'uso **lt** Skaitykite naudojimo instrukcijas **lv** Skatīt lietošanas instrukcijas **nl** Raadpleeg de gebruiksaanwijzingen **no** Se bruksanvisningen **pl** Patrz instrukcja użycia **pt** Consultar as instruções de utilização **ro** Consultați instrucțiunile de utilizare **ru** Ознакомьтесь с инструкцией по применению **sk** Pozri návod na použitie **sl** Glejte navodila za uporabo **sv** Läs bruksanvisningen **tr** Kullanım için talimatlara başvurun

MD **en** Medical Device **bg** Медицинско изделие **cs** Zdravotnický prostředek **de** Medizinprodukt **da** Medicinsk udstyr **es** Producto sanitario **et** Meditsiiniseade **fi** Lääkinnällinen laite **fr** Dispositif médical **el** Ιατροτεχνολογικό προϊόν **hr** Medicinski proizvod **hu** Orvostechnikai eszköz **it** Dispositivo medico **lt** Medicinos priemonė **lv** Medicīniska ierīce **nl** Medisch hulpmiddel **no** Medisinsk utstyr **pl** Wyrób medyczny **pt** Dispositivo médico **ro** Dispozitiv medical **ru** Медицинское изделие **sk** Zdravotnícka pomôcka **sl** Medicinski pripomoček **sv** Medicinteknisk produkt **tr** Tıbbi cihaz

en Manufacturer **bg** Производител **cs** Výrobce **de** Hersteller **da** Producent **es** Fabricante **et** Tootja **fi** Valmistaja **fr** Fabricant **el** Κατασκευαστής **hr** Proizvođač **hu** Gyártó **it** Fabbricante **lt** Gamintojas **lv** Ražotājs **nl** Fabrikant **no** Tilvirker **pl** Producent **pt** Fabricante **ro** Producător **ru** Изготовитель **sk** Výrobca **sl** Proizvajalec **sv** Tillverkare **tr** Üretici

en Date of manufacture **bg** Дата на производство **cs** Datum výroby **de** Herstellungsdatum **da** Fremstillingsdato **es** Fecha de fabricación **et** Valmistamiskuupäev **fi** Valmistuspäivä **fr** Date de fabrication **el** Ημερομηνία κατασκευής **hr** Datum proizvodnje **hu** Gyártás dátuma **it** Data di fabbricazione **lt** Pagaminimo data **lv** Izgatavošanas datums **nl** Productiedatum **no** Fremstillingsdato **pl** Data produkcji **pt** Data de fabrico **ro** Data fabricației **ru** Дата изготовления **sk** Dátum výroby **sl** Datum proizvodnje **sv** Tillverkningsdatum **tr** Üretim tarihi

en Number of units **bg** Брой елементи **cs** Počet jednotek **de** Anzahl der Einheiten **da** Antal enheder **es** Cantidad de unidades **et** Ühikute arv **fi** Yksiköiden määrä **fr** Nombre d'unités **el** Αριθμός μονάδων **hr** Broj jedinica **hu** Az egységek darabszáma **it** Numero di unità **lt** Vienetų skaičius **lv** Vienību skaits **nl** Aantal eenheden **no** Antall enheter **pl** Liczba zestawów **pt** Número de unidades **ro** Numărul de unități **ru** Количество **sk** Počet jednotiek **sl** Število enot **sv** Antal enheter **tr** Birim sayısı

en Keep away from sunlight **bg** Да не се излага на пряка слънчева светлина **cs** Chraňte před slunečním zářením **de** Vor Sonnenlicht schützen **da** Beskyttes mod sollys **es** Mantener alejado de la luz directa del sol **et** Hoida otsese päikesevalguse eest **fi** Pidettävä poissa auringonvalosta **fr** Conserver à l'abri des rayons du soleil **el** Προστατέψτε από την άμεση έκθεση σε ηλιακό φως **hr** Ne izlažite sunčevoj svjetlosti **hu** Napfénytől távol tartandó **it** Tenere lontano dalla luce diretta del sole **lt** Saugokite nuo saulės šviesos **lv** Sargāt no saules gaismas **nl** Buiten direct zonlicht bewaren **no** Må ikke utsettes for sollys **pl** Chronić przed światłem słonecznym **pt** Manter afastado da luz solar **ro** A nu se expune direct la razele solare **ru** Защищать от воздействия прямых солнечных лучей **sk** Nevystavujte slnečným lúčom **sl** Hranite na temnem **sv** Skyddas mot solljus **tr** Güneş ışığından uzak tutun

en Keep dry **bg** Да се съхранява на сухо **cs** Uchovávejte v suchu **de** Trocken aufbewahren **da** Opbevares tørt **es** Mantener en un lugar seco **et** Hoida kuivas **fi** Pidettävä kuivana **fr** Conserver au sec **el** Φυλάσσετε σε ξηρό χώρο **hr** Držati na suhom **hu** Szárazon tartandó **it** Tenere asciutto **lt** Laikykite sausai **lv** Uzglabāt sausā vietā **nl** Koel bewaren **no** Hold tørr **pl** Chronić przed wilgocią **pt** Conservar seco **ro** A se păstra într-un loc uscat **ru** Защищать от влаги **sk** Uchovávajte v suchu **sl** Hranite na suhem **sv** Förvaras torrt **tr** Kuru yerde saklayın

en Use-by date **bg** Валидно до **cs** Datum spotřeby **de** Verwendbar bis **da** Anvendes inden **es** Fecha de caducidad **et** „Kõlblik kuni" kuupäev **fi** Viimeinen käyttöpäivä **fr** Date limite d'utilisation **el** Ημερομηνία λήξης **hr** Upotrijebiti do **hu** Felhasználhatóság **it** Data di scadenza **lt** Naudokite iki **lv** Izlietot līdz **nl** Uiterste gebruiksdatum **no** Siste forbruksdato **pl** Termin ważności **pt** Prazo de validade **ro** A se folosi înainte de **ru** Годен до **sk** Dátum spotreby **sl** Rok uporabnosti **sv** Utgångsdatum **tr** Son kullanma tarihi

en Do not use if package is damaged **bg** Не използвайте, ако опаковката е повредена **cs** Nepoužívejte, je-li balení poškozeno **de** Bei beschädigter Verpackung nicht verwenden **da** Må ikke anvendes, hvis pakken er beskadiget **es** No utilizar si observa que el embalaje está dañado **et** Mitte kasutada juhul, kui pakend on vigastatud **fi** Ei saa käyttää, jos pakkaus on vahingoittunut **fr** Ne pas utiliser si l'emballage est endommagé **el** Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία έχει υποστεί ζημιά **hr** Ne upotrebljavate ako je pakiranje oštećeno **hu** Ne használja, ha a csomagolás sérült **it** Non usare se la confezione è danneggiata **lt** Nenaudokite, jei pakuotė yra pažeista **lv** Nelietot, ja iepakojums ir bojāts **nl** Niet gebruiken indien de verpakking is beschadigd **no** Må ikke brukes dersom pakningen er skadet **pl** Nie należy używać zestawu, jeżeli opakowanie jest uszkodzone **pt** Não utilizar se a embalagem estiver danificada **ro** A nu se folosi dacă ambalajul a fost deschis sau distrus **ru** Не использовать, если упаковка повреждена **sk** Nepoužívajte, ak je balenie poškodené **sl** Ne uporabi, če je embalaža poškodovana **sv** Använd inte produkten om förpackningen är skadad **tr** Ambalajı hasarlıysa kullanmayın

en Single use **bg** За еднократна употреба **cs** Pro jednorázové použití **de** Nur zum Einmalgebrauch **da** Må ikke genbruges **es** No reutilizar **et** Ühekordseks kasutamiseks **fi** Kertakäyttöinen **fr** Utilisation unique **el** Προϊόν μίας χρήσης **hr** Jednokratna uporaba **hu** Nem szabad újrafelhasználni **it** Monouso **lt** Naudokite vieną kartą **lv** Nelietot atkārtoti **nl** Voor eenmalig gebruik **no** Engangsbruk **pl** Wyrób jednorazowego użytku **pt** Utilização única **ro** De unică folosință **ru** Не использовать повторно **sk** Jednorazové použitie **sl** Za enkratno uporabo **sv** Engångsbruk **tr** Tek kullanımlıktır

2 STERILIZE **en** Do not resterilize **bg** Забранено повторно стерилизиране **cs** Nesterilizujte opakovaně **de** Nicht erneut sterilisieren **da** Må ikke resteriliseres **es** No reesterilizar **et** Mitte resterili-seerida **fi** Ei saa steriloida uudelleen **fr** Ne pas restériliser **el** Μην επαναποστειρώνετε **hr** Nemojte ponovno sterilizirati **hu** Nem szabad újrasterilizálni **it** Non risterilizzare **lt** Nesterilizuokite pakartotinai **lv** Nesterilizēt atkārtoti **nl** Niet opnieuw steriliseren **no** Må ikke resteriliseres **pl** Nie sterylizować ponownie **pt** Não reesterilizar **ro** Nu resterilizați **ru** Не стерилизовать повторно **sk** Nesterilizujte opakovane **sl** Ne sterilizirajte ponovno **sv** Sterilisera inte om produkten **tr** Yeniden sterilize etmeyin

STERILE EO **en** Sterilized using ethylene oxide **bg** Стерилизирано спомоща на етиленов оксид **cs** Sterilizováno etylenoxidem **de** Sterilisiert mit Ethylenoxid **da** Steriliseret med ethylenoxid **es** Esterilizar mediante óxido de etileno **et** Steriliseeritud etüleenoksiidiga **fi** Steriloitu eteenioksidilla **fr** Stérilisé à l'oxyde d'éthylène **el** Αποστειρωμένο με οξείδιο του αιθυλενίου **hr** Sterilizirano etilen oksidom **hu** Etilénoxiddal sterilizálva **it** Metodo di sterilizzazione: ossido di etilene **lt** Sterilizuota etileno oksidu **lv** Sterilizēts ar etilēna oksīdu **nl** Gesteriliseerd met gebruik van ethyleenoxide **no** Sterilisert med etylenoksid **pl** Sterylizowano tlenkiem etylenu **pt** Esterilizado em óxido de etileno **ro** Sterilizat cu oxid de etilen **ru** Стерилизовано окисью этилена **sk** Sterilizované etylénoxidom **sl** Sterilizirano z etilenoksidom **sv** Steriliserad med etylenoxid **tr** Etilen oksit kullanılarak sterilize edilmiştir

OTW **en** Over the Wire **bg** Техника „Over-the-wire" **cs** Přes drát **de** Over-the-Wire **da** Over-the-wire **es** Sobre alambre **et** Juhtetraat **fi** Vaijerikäyttöinen **fr** Par fil-guide **el** Με συρμάτινο οδηγό **hr** Preko žice **hu** Sztent a vezetődróton típusú rendszer **it** Over the Wire **lt** Per laidą **lv** Over the Wire tehnoloģija **nl** Voor een wire **no** Over vaieren **pl** Typu over-the-wire **pt** Coaxial **ro** Peste fir **ru** Доставка по проводнику **sk** Pomocou drôtu **sl** Prek žice **sv** Över ledaren **tr** Tel Üstü

REF **en** Catalogue number **bg** Каталожен номер **cs** Katalogové číslo **de** Artikelnummer **da** Katalognummer **es** Número de catálogo **et** Katalooginumber **fi** Luettelonumero **fr** N° de référence **el** Αριθμός καταλόγου **hr** Kataloški broj **hu** Katalógusszám **it** Numero di catalogo **lt** Katalogo numeris **lv** Kataloga numurs **nl** Catalogusnummer **no** Katalognummer **pl** Numer katalogowy **pt** Número de catálogo **ro** Numărul de catalog **ru** Номер по каталогу **sk** Katalógové číslo **sl** Kataloška številka **sv** Katalognummer **tr** Katalog numarası

PRINT DATE 02.07.2020: Printed documents must be compared with the valid revision before use!
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

LOT **en** Lot number **bg** Номер на партида **cs** Číslo šarže **de** Fertigungslosnummer **da** Lotnummer **es** Número de lote **et** Partii number **fi** Eränumero **fr** Numéro de lot **el** Αριθμός παρτίδας **hr** Broj serije **hu** Tételszám **it** Numero di lotto **lt** Serijos numeris **lv** Partijas numurs **nl** Lotnummer **no** LOT-nummer **pl** Numer partii **pt** Número de lote **ro** Număr lot **ru** Номер партии **sk** Číslo šarže **sl** Številka serije **sv** Lotnummer **tr** Lot numarası

REC GW **en** Recommended Guide Wire **bg** Препоръчителен водач **cs** Doporučený vodicí drát **de** Empfohlener Führungsdraht **da** Anbefalet guidewire **es** Alambre guía recommendado **et** Soovitatav juhttraat **fi** Suositeltava ohjauslanka **fr** Fil-guide recommandé **el** Συνιστώμενος συρμάτινος οδηγός **hr** Preporučena žica vodilica **hu** Ajánlott vezetődrót **it** Filo guida consigliato **lt** Rekomenduojama pravedėjo viela **lv** Ieteicamā vadītājstīga **nl** Aanbevolen voerdraad **no** Anbefalt ledevaier **pl** Zalecany prowadnik **pt** Fio-guia recomendado **ro** Ghidajul recomandat **ru** Екомендуемый проводник **sk** Doporučený vodiaci drôt **sl** Priporočena vodilna žica **sv** Rekommenderad ledare **tr** Önerilen Kılavuz Tel

REC MIN INTRO **en** Recommended Minimum Sheath Introducer **bg** Препоръчва се водач с минимална обшивка **cs** Doporučený minimální zavaděč pouzdra **de** Empfohlene Mindestgröße der Einführschleuse **da** Anbefalet min. sheatintroducer **es** Introductor de vaina mínimo recomendado **et** Soovitatud väikseim sisestushülss **fi** Pienin suositeltu sisäänvientiputki **fr** Introducteur minimum recommandé **el** Συνιστώμενος ελάχιστος εισαγωγέας θηκαριού **hr** Preporučena minimalna uvodnica **hu** A bevezetőhüvely ajánlott minimális mérete **it** Valore consigliato minimo per l'introduttore a guaina **lt** Rekomenduojamas mažiausias movos obturatorius **lv** Ieteicamais minimālais apvalka ievadītājs **nl** Aanbevolen minimale sheath introducer **no** Anbefalt minste størrelse innføringshylse **pl** Zalecany minimalny introduktor z koszulkę **pt** Bainha Introdutora Mínima Recomendada **ro** Intubator minim recomandat **ru** Рекомендуемый минимальный размер гильзы интродьюсера **sk** Odporúčané minimálne zavádzacie puzdro **sl** Priporočena najmanjša uvajalna cevka **sv** Rekommenderad minsta hylsinförare **tr** Önerilen Minimum Sheath İntroduser

P **en** Pressure **bg** Налягане **cs** Tlak **de** Druck **da** Tryk **es** Presión **et** Rõhk **fi** Paine **fr** Pression **el** Πίεση **hr** Tlak **hu** Nyomás **it** Pressione **lt** Slėgis **lv** Spiediens **nl** Druk **no** Trykk **pl** Ciśnienie **pt** Pressão **ro** Presiune **ru** Давление **sk** Tlak **sl** Tlak **sv** Tryck **tr** Basınç

NP **en** Nominal Pressure **bg** Номинално налягане **cs** Jmenovitý tlak **de** Nominaldruck **da** Nominelt tryk **es** Presión nominal **et** Nimirõhk **fi** Nimellispaine **fr** Pression nominale **el** Ονομαστική πίεση **hr** Nominalni tlak **hu** Névleges nyomás **it** Pressione nominale **lt** Nominalus spaudimas **lv** Nominālais spiediens **nl** Nominale druk **no** Nominelt trykk **pl** Ciśnienie nominalne **pt** Pressão nominal **ro** Resiunea nominală **ru** Номинальное давление **sk** Menovitý tlak **sl** Nominalni tlak **sv** Nominellt tryck **tr** Nominal basınç

RBP **en** Rated Burst Pressure **bg** Необходимо налягане за разширение **cs** Jmenovitý tlak roztržení **de** Nenn-Berstdruck **da** Nominelt sprængningstryk **es** Presión nominal de rotura **et** Hinnanguline lõhkemisrõhk **fi** Nimellinen puhkeamispaine **fr** Pression d'éclatement nominale **el** Ονομαστική πίεση ρήξης **hr** Maksimalni dopušteni tlak **hu** Névleges repedési nyomás **it** Pressione massima di rottura **lt** Nominalus balionėlio sprogimo spaudimas **lv** Aprēķinātais plīšanas spiediens **nl** Nominale barstdruk **no** Nominelt bristetrykk **pl** Nominalne ciśnienie rozerwania **pt** Pressão nominal de rebentamento **ro** Presiunea de spargere estimată **ru** Давление разрыва баллона **sk** Menovitý tlak roztrhnutia **sl** Nominalni razpočni tlak **sv** Nominellt bristningstryck **tr** Nominal patlama basıncı

STENT **en** Expanded Stent Outer Diameter **bg** Външен диаметър на разгънат стент **cs** Vnější průměr prodlouženého stentu **de** Außendurchmesser des entfalteten Stents **da** Ydre diameter på udvidet stent **es** Diámetro exterior de la endoprótesis expandida **et** Laiendatud stendi välisläbimõõt **fi** Laajennetun stentin ulkohalkaisija **fr** Diamètre externe du stent déployé **el** Εξωτερική διάμετρος εκπτυγμένου στεντ **hr** Vanjski promjer proširenog stenta **hu** Kiterjesztett sztent külső átmérője **it** Diametro esterno dello stent espanso **lt** Išplėsto stento išorinis skersmuo **lv** Izplesta stenta ārējais diametrs **nl** Buitendiameter van geëxpandeerde stent **no** Ytre diameter på utvidet stent **pl** Zewnętrzna średnica rozprężonego stentu **pt** Diâmetro Exterior do Stent Expandido **ro** Diametrul exterior extins al stentului **ru** Наружный диаметр расширенного стента **sk** Vonkajší priemer roztiahnutého stentu **sl** Zunanji premer razširjene žilne opornice **sv** Ytterdiameter för utvidgad stent **tr** Genişletilmiş Stent Dış Çapı

STENT **en** Labelled Stent Length **bg** Дължина на стента, означена върху етикета **cs** Délka stentu uvedená na štítku **de** Auf dem Etikett angegebene Stentlänge **da** Anført stentlængde **es** Longitud de la endoprótesis que se indica en la etiqueta **et** Stendi pikkus, mis on märgitud sildile **fi** Etikettiin merkitty stentin pituus **fr** Longueur du stent indiqué sur l'étiquette **el** Μήκος επισημασμένου στεντ **hr** Označena dužina stenta **hu** A címkén feltüntetett sztenthosszúság **it** Lunghezza dello stent riportata in etichetta **lt** Uzlīmē norādītais stenta garums **lv** Etiketėje nurodytas stento ilgis **nl** Gelabelde stentlengte **no** Lengde på merket stent **pl** Długość stentu podana na etykiecie **pt** Comprimento do stent marcado **ro** Etiketėje nurodytas stento ilgis **ru** Длина стента, указанная на этикетке **sk** Dĺžka stentu na označení **sl** Oznaka za dolžino stenta **sv** Angiven stentlängd **tr** İşaretli Stent Uzunluğu

USAB CAT **en** Usable Catheter Length **bg** Използваема дължина на катетър **cs** Použitelná délka katétru **de** Nutzbare Katheterlänge **da** Anvendelig kateterlængde **es** Longitud utilizable del catéter **et** Kasutatav kateetri pikkus **fi** Käytettävä katetrin pituus **fr** Longueur de cathéter utilisable **el** Ωφέλιμο μήκος καθετήρα **hr** Upotrebljiva dužina katetera **hu** A katéter hasznos hosszúsága **it** Lunghezza utilizzabile del catetere **lt** Naudingas kateterio ilgis **lv** Izmantojamais katetra garums **nl** Bruikbare katheterlengte **no** Brukbar kateter lengde **pl** Użytkowa długość cewnika **pt** Comprimento Útil do Cateter **ro** Lungimea utilizabilă a cateterului **ru** Используемая длина катетера **sk** Použiteľná dĺžka katétra **sl** Uporabna dolžina katetra **sv** Användbar kateterlängd **tr** Kullanılabilir Kateter Uzunluğu

en Non-pyrogenic **bg** Негоримо **cs** Nepyrogenní **de** Nicht pyrogen **da** Ikke-pyrogen **es** Apirógeno **et** Mittepürogeenne **fi** Ei-pyrogeeninen **fr** Apyrogène **el** Μη πυρετογόνο **hr** Nepirogeno **hu** Nem pirogén **it** Non pirogenico **lt** Nepirogeniškas **lv** Nepirogēns **nl** Niet pyrogeen **no** Ikke-pyrogen **pl** Produkt niepirogenny **pt** Apirogénico **ro** Apirogen **ru** Апирогенно **sk** Nepyrogénne **sl** Nepirogeno **sv** Icke-pyrogen **tr** Apirojendir

MR **en** MR Conditional **bg** Безопасно за ЯМР при определени условия **cs** Podmíněně použitelné v prostředí MR **de** Bedingt MRT-tauglich **da** MR-betinget **es** Compatibilidad condicional con la RM **et** MR-tingimuslik **fi** Ehdollisesti MR-yhteensopiva **fr** Compatibilité conditionnelle à la résonance magnétique (RM) **el** Υπό συνθήκες συμβατό με MR **hr** Sigurno za uporabu s MR pri točno propisanim uvjetima **hu** MR-kompatibilis **it** Compatibilità RMI condizionata **lt** Santykinai saugus naudoti MR aplinkoje **lv** Saderīgs ar MR **nl** MR-veilig **no** MR-betinget **pl** Produkt bezpieczny w określonych warunkach wykonywania badania MRI **pt** Condicional para RM **ro** Condiționat(ă) de RM **ru** Условно МР-совместимый **sk** Podmienečne bezpečné v prostredí MR **sl** Slikanje MR dovoljeno pri določenih pogojih **sv** MR-villkorat **tr** MR koşullu

PRINT DATE 02.07.2020: Printed documents must be compared with the valid revision before use!
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

en	**English**	Instructions for use	**5**
bg	**Български**	Инструкции за употреба	**7**
cs	**Česky**	Návod na použití	**10**
de	**Deutsch**	Gebrauchsanweisung	**12**
da	**Dansk**	Brugsanvisning	**15**
es	**Español**	Instrucciones de uso	**18**
et	**Eesti**	Kasutusjuhised	**20**
fi	**Suomi**	Käyttöohjeet	**23**
fr	**Français**	Mode d'emploi	**26**
el	**Ελληνικά**	Οδηγιεσ χρησησ	**29**
hr	**Hrvatski**	Upute za uporabu	**31**
hu	**Magyar**	Használati útmutató	**34**
it	**Italiano**	Istruzioni per l'uso	**37**
lt	**Lietuviškai**	Naudojimo instrukcijos	**39**
lv	**Latviešu**	Lietošanas instrukcijas	**42**
nl	**Nederlands**	Gebruiksaanwijzingen	**45**
no	**Norsk**	Bruksanvisning	**47**
pl	**Polski**	Instrukcja użytkowania	**50**
pt	**Português**	Instruções de utilização	**53**
ro	**Română**	Instrucţiuni de utilizare	**56**
ru	**Русский**	Инструкция по применению	**58**
sk	**Slovenčina**	Návod na použitie	**61**
sl	**Slovenščina**	Navodila za uporabo	**64**
sv	**Svenska**	Bruksanvisning	**67**
tr	**Türkçe**	Kullanım talimatları	**69**

PRINT DATE 02.07.2020: Printed documents must be compared with the valid revision before use!
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

English

Note:

These instructions are intended for all personnel who will handle the "BeGraft Peripheral Stent Graft System" and should be read carefully before using this product!

BeGraft Peripheral Stent Graft System

1 Product Description

The BeGraft Peripheral Stent Graft System consists of

- a balloon-expandable cobalt-chromium (L605) stent, covered with an ePTFE graft material, which is
- pre-mounted on the balloon of an over the wire (OTW) stent delivery system.

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Two radiopaque markers located underneath the balloon identify the BeGraft Peripheral Stent Graft position and fluoroscopically mark the working length of the balloon.

2 How Supplied

Contents per Package:

- One "BeGraft Peripheral Stent Graft System"
- One Compliance Card
- One Implant Information Card
- One Instruction for Use

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Storage: Store in a dry, dark and cool place. Do not use the BeGraft System after the "use by" date (indicated on the product label).

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3 Indications

The BeGraft Peripheral Stent Graft System is indicated for intraluminal chronic placement in iliac and renal arteries for:

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- Any use of the BeGraft Peripheral Stent Graft System for procedures other than those indicated in these instructions is not recommended.
- It is explicitly advised that the product is not designed and tested to be used beyond the approved indications (e.g. in the central circulatory system or in overlapping configuration). The manufacturer shall not be liable if the product is used beyond the approved indications. In this case, the treatment is at the sole responsibility of the treating physician.
- The long-term effectiveness of the BeGraft Peripheral Stent Graft has not been fully established.
- Do not use the system after the "use by" date (indicated on the product label).
- Do not expose this device to organic solvents, ionizing radiation or ultraviolet light.
- The BeGraft Peripheral Stent Graft System should only be used by physicians trained in interventional techniques and percutaneous transluminal angioplasty.
- Assure that the labelled stent dimensions are suitable for the intended procedure.
- Do not use if the package is open or damaged.
- Employ aseptic techniques during removal of the product from the package and during use.
- Prior to use, please carefully inspect the product to verify that it has not been damaged during shipment (e.g. bends, kinks or other damage) and visually ensure that the crimped BeGraft Stent is positioned between the balloon markers.
- Do not touch the crimped BeGraft Peripheral Stent as this could decrease the retention force of the stent graft on the balloon.
- Administer appropriate anticoagulation / antiplatelet treatment before insertion of the BeGraft Peripheral Stent Graft System. Continue the anticoagulation / antiplatelet treatment after the implantation according to current versions of national guide lines and international standards.
- For balloon inflation it is mandatory to use a pressure monitoring device; optimal is the use of an inflation device with a calibrated manometer.
- Only use diluted contrast medium. Do not use air or any gaseous medium to inflate the balloon.
- The inflation pressure must not exceed the Rated Burst Pressure (RBP) indicated on the product label and compliance card.
- Do not over tighten the haemostatic valve, as lumen constriction may occur, possibly impairing inflation/deflation of the balloon.

PRINT DATE 02.07.2020: Printed documents must be compared with the valid revision before use!
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

- Do not attempt to straighten a catheter if the shaft is bent or kinked, as it may break. Instead, prepare a new system.
- Do not use the BeGraft Delivery System for post dilatation of the BeGraft Peripheral Stent.
- The BeGraft Peripheral Stent Graft System is intended to perform as a system. The stent graft should not be removed for use in conjunction with other dilatation catheters, nor should the BeGraft Delivery System be used in conjunction with other stents.
- Only insert the BeGraft Peripheral Stent Graft System under fluoroscopic guidance.
- Use the recommended minimum introducer sheath size (refer to product label), as the integrity of the stent graft could be impaired (dislodgement of the stent graft and / or the ePTFE graft material), causing an embolism in the worst case.
- Only manipulate the BeGraft Peripheral Stent Graft System forward and backward when the balloon has been fully deflated under vacuum.
- Prior to the deployment of the BeGraft Peripheral Stent Graft, utilize fluoroscopy to verify that the BeGraft Peripheral Stent Graft has been properly positioned. Do not deploy the BeGraft Peripheral Stent Graft if it is not properly positioned.
- Once partly or fully deployed, the BeGraft Peripheral Stent Graft cannot be repositioned or removed.
- Under-expansion of the BeGraft Peripheral Stent Graft may result in migration.
- If the balloon of the BeGraft Delivery System gets caught in the BeGraft Peripheral Stent Graft during removal of the balloon, it should be removed immediately by surgery.
- A balloon burst during implantation may cause the migration or deformation of the BeGraft Peripheral Stent Graft. In such a case, the Delivery System should be retracted carefully, leaving the guide wire in place. A new balloon catheter with a small profile and semi-compliant balloon material should be carefully inserted over the guide wire and the BeGraft Peripheral Stent Graft can be expanded completely.
- Placing the BeGraft Peripheral Stent Graft across a bifurcation / sidebranch will compromise the blood flow and impact future diagnostic or therapeutic procedures.
- A stent previously implanted proximally to the lesion/ defect can be deformed or dislodged by contact, e.g. with the guide wire or the BeGraft Delivery System.
- Re-crossing a partially or fully deployed BeGraft Peripheral Stent Graft with adjunct devices must be performed with extreme caution to avoid BeGraft Stent damage or deformation.
- The use of laser devices, mechanical atherectomy and hyperthermia are not recommended in the stented region.

6.1 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the **BeGraft Peripheral Stent Graft is MR Conditional**. A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- static magnetic field of 1.5 Tesla and 3 Tesla, with
- maximum spatial field gradient of 5,100 G/cm (51 T/m)
- maximum force product of 91,000,000 G^2/cm (91 T^2/m)
- theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the BeGraft Peripheral Stent Graft is expected to produce a maximum temperature rise of less than

- **2.7°C** (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of ≈ 1.4°C (2 W/kg, 1.5 Tesla)
- **2.1°C** (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of ≈ 1.0°C (2 W/kg, 3 Tesla)

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In non-clinical testing, the image artifact caused by the device extends approximately 12.4 mm from the BeGraft Peripheral Stent Graft when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

7 Procedure

7.1 Inspection Procedure Prior To Use

Prior to using the BeGraft Peripheral Stent Graft System, carefully remove the system from the package and inspect the catheter for bends, kinks and other damage. Do not use the system if any defects are noted.

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7.5 Stent Graft Delivery Procedure

1. Wipe the exposed guide wire with heparinized saline to remove residual blood or contrast medium.
2. Fully open the haemostatic valve. Maintain neutral pressure on the inflation device.
3. Backload the Delivery System onto the proximal portion of the guide wire while maintaining guide wire position across the target lesion / defect.
4. For insertion of the BeGraft Stent Graft into the appropriately sized introducer sheath, it is recommended to carefully hold the stent graft at the proximal end (with your fingers) when pushing it into the haemostatic valve and to closely observe if the position of the stent graft on the balloon is affected. If this is the case, this device shall not be used and the unit should be returned to Bentley InnoMed GmbH.

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5. Advance the Delivery System over the guide wire to the target lesion / defect. Utilize the radiopaque balloon markers to position the BeGraft Peripheral Stent Graft across the lesion / defect; perform angiography to confirm the BeGraft Peripheral Stent Graft position. If applicable tighten the haemostatic valve.

Note: If you notice during the process of moving the Delivery System into position that the BeGraft Peripheral Stent Graft has moved on the balloon, do not deploy the BeGraft Peripheral Stent Graft.

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6. The BeGraft Peripheral Stent Graft is now ready to be deployed.

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1. Maintain negative pressure to allow the balloon to remain fully deflated during removal.

PRINT DATE 02.07.2020: Printed documents must be compared with the valid revision before use!
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

2. Withdraw the Delivery System through the sheath with the guide wire remaining across the lesion / defect.

Caution: Should unusual resistance be felt at any time during removal of the Delivery System post BeGraft Peripheral Stent Graft implantation, the introducer sheath and Delivery System should be removed as a single unit. In this case, position the proximal balloon marker just distal to the tip of the introducer sheath, then advance the guide wire in the anatomy as far distally as safely possible and secure the Delivery System to the introducer sheath; then remove the introducer sheath and Delivery System as a single unit. Failure to follow these steps and / or applying excessive force to the Delivery System can potentially result in loss or damage of the Delivery System components.

3. Confirm optimal wall apposition of the BeGraft Stent Graft using standard angiography techniques. If the primary angiographic result is not optimal, you can further expand the BeGraft Stent Graft with the help of an appropriate balloon catheter.

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English

Note:

These instructions are intended for all personnel who will handle the "BeGraft Peripheral Plus Stent Graft System" and should be read carefully before using this product!

BeGraft Peripheral Plus Stent Graft System

1 Product Description

The BeGraft Peripheral Plus Stent Graft System consists of

- a balloon-expandable stent graft (composed of two cobalt-chromium (L605) stents assembled one above the other, each covered with ePTFE graft material),
- pre-mounted on the balloon of an over the wire (OTW) stent delivery system.

The Delivery System is compatible with 0.035" (0.89 mm) guide wires and is available in usable lengths of 75 cm and 120 cm.

Two radiopaque markers located underneath the balloon identify the BeGraft Peripheral Plus Stent Graft position and fluoroscopically mark the working length of the balloon.

2 How Supplied

Contents per Package:

- One "BeGraft Peripheral Plus Stent Graft System"
- One Compliance Card
- One Implant Information Card
- One Instruction for Use

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Storage: Store in a dry, dark and cool place. Do not use after the "use by" date (indicated on the product label).

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3 Indications

The BeGraft Peripheral Plus Stent Graft System is indicated for intraluminal chronic placement in iliac and renal arteries for:

- Restoring and improving the patency and
- Treating aneurysms, acute perforations, acute ruptures and fistulas.

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1. When PTA is technically not possible (e. g. not feasible to access the lesion / defect with guide wire or balloon catheter)
2. Aneurysms immediately adjacent to the site of stent implantation
3. Stenosis distal to the site of stent implantation
4. Lesions in or adjacent to essential collateral(s)
5. Lesions in locations subject to external compression
6. Heavily calcified lesions resistant to PTA
7. Patients with diffuse distal disease resulting in poor stent outflow
8. Patients with a history of coagulation disorders
9. Patients with aspirin allergy or bleeding complications and patients unable or unwilling to tolerate anticoagulant / antiplatelet therapy and / or non-responders to anticoagulant / antiplatelet therapy
10. Fresh thrombus formation
11. Patients with known hypersensitivity to the stent material (L605) and / or PTFE

5 Potential Adverse Events

Potential adverse effects include, but may not be limited to the following:

- Abscess
- Acute or subacute closure of stent graft
- Allergic reaction to contrast medium
- Allergic reaction to the implanted stent graft material (L605) and/or PTFE
- Angina pectoris / abdominalis
- Arterial spasm
- Arteriovenous fistulae
- Arrhythmia
- Bleeding or hematoma at the puncture site or from anticoagulation/antiplatelet medication (may require transfusion)
- Death
- Distal emboli (air, particle, thromboembolic emboli) with potential limb loss
- Dislodgement of stent graft and/or loss of ePTFE-cover causing distal embolisation with potential limb loss
- Drug reactions with respect to medication
- Emergency vascular surgery
- Exacerbation of pre-existing conditions
- Hematuria
- Hemoptysis
- Hemorrhage or hematoma
- Hemorrhagic stroke
- Hypo / hypertension
- Infection, sepsis
- Inflammation, fever
- Injury, dissection, perforation, rupture, sudden and total occlusion of the vessel or other damages to the indicated arteries
- Intimal flap, tear or trauma
- Ischemia
- Myocardial infarction
- Nausea and vomiting
- Pseudoaneurysm formation
- Renal insufficiency / failure
- Restenosis of treated segment
- Seizure
- Short-term hemodynamic deterioration
- Thrombosis / embolisation / occlusion (stent graft or other)
- Tissue reaction, necrosis

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- Any use of the BeGraft Peripheral Plus Stent Graft System for procedures other than those indicated in these instructions is not recommended.
- It is explicitly advised that the product is not designed and tested to be used beyond the approved indications (e.g. in the central circulatory system or in overlapping configuration). The manufacturer shall not be liable if the product is used beyond the approved indications. In this case, the treatment is at the sole responsibility of the treating physician.
- The long-term effectiveness of the BeGraft Peripheral Plus Stent Graft has not been fully established.
- Do not use the system after the "use by" date (indicated on the product label).
- Do not expose this device to organic solvents, ionizing radiation or ultraviolet light.
- The BeGraft Peripheral Plus Stent Graft System should only be used by physicians trained in interventional techniques and percutaneous transluminal angioplasty.
- Assure that the labelled stent dimensions are suitable for the intended procedure.
- Do not use if the package is open or damaged.
- Employ aseptic techniques during removal of the product from the package and during use.
- Prior to use, please carefully inspect the product to verify that it has not been damaged during shipment (e.g. bends, kinks or other damage) and visually ensure that the crimped BeGraft Stent Graft is positioned between the balloon markers.
- Do not touch the crimped BeGraft Peripheral Plus Stent Graft as this could decrease the retention force of the stent graft on the balloon.
- Administer appropriate anticoagulation / antiplatelet treatment before insertion of the BeGraft Peripheral Plus Stent Graft System. Continue the anticoagulation / antiplatelet treatment after the implantation according to current versions of national guide lines and international standards.
- For balloon inflation it is mandatory to use a pressure monitoring device; optimal is the use of an inflation device with a calibrated manometer.
- Only use diluted contrast medium. Do not use air or any gaseous medium to inflate the balloon.
- The inflation pressure must not exceed the Rated Burst Pressure (RBP) indicated on the product label and compliance card.
- Do not over tighten the haemostatic valve, as lumen constriction may occur, possibly impairing inflation/deflation of the balloon.

PRINT DATE 29.06.2020: Printed documents must be compared with the valid revision before use!
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

- Do not attempt to straighten a catheter if the shaft is bent or kinked, as it may break. Instead, prepare a new system.
- Do not use the BeGraft Delivery System for post dilatation of the BeGraft Peripheral Plus Stent Graft.
- The BeGraft Peripheral Plus Stent Graft System is intended to perform as a system. The stent graft should not be removed for use in conjunction with other dilatation catheters, nor should the BeGraft Delivery System be used in conjunction with other stents.
- Only insert the BeGraft Peripheral Plus Stent Graft System under fluoroscopic guidance.
- Use the recommended minimum introducer sheath size (refer to product label), as the integrity of the stent graft could be impaired (dislodgement of the stent graft and / or the ePTFE graft material), causing an embolism in the worst case.
- Only manipulate the BeGraft Peripheral Plus Stent Graft System forward and backward when the balloon has been fully deflated under vacuum.
- Prior to the deployment of the BeGraft Peripheral Plus Stent Graft, utilize fluoroscopy to verify that the BeGraft Peripheral Plus Stent Graft has been properly positioned. Do not deploy the BeGraft Peripheral Plus Stent Graft if it is not properly positioned.
- Once partly or fully deployed, the BeGraft Peripheral Plus Stent Graft cannot be repositioned or removed.
- Under-expansion of the BeGraft Peripheral Plus Stent Graft may result in migration.
- If the balloon of the BeGraft Delivery System gets caught in the BeGraft Peripheral Plus Stent Graft during removal of the balloon, it should be removed immediately by surgery.
- A balloon burst during implantation may cause the migration or deformation of the BeGraft Peripheral Plus Stent Graft. In such a case, the Delivery System should be retracted carefully, leaving the guide wire in place. A new balloon catheter with a small profile and semi-compliant balloon material should be carefully inserted over the guide wire and the BeGraft Peripheral Plus Stent Graft can be expanded completely.
- Placing the BeGraft Peripheral Plus Stent Graft across a bifurcation / sidebranch will compromise the blood flow and impact future diagnostic or therapeutic procedures.
- A stent previously implanted proximally to the lesion/ defect can be deformed or dislodged by contact, e.g. with the guide wire or the BeGraft Delivery System.
- Re-crossing a partially or fully deployed BeGraft Peripheral Plus Stent Graft with adjunct devices must be performed with extreme caution to avoid BeGraft Stent damage or deformation.
- The use of laser devices, mechanical atherectomy and hyperthermia are not recommended in the stented region.

6.1 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the **BeGraft Peripheral Plus Stent Graft is MR Conditional**. A patient with this device can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- static magnetic field of 1.5 Tesla and 3 Tesla, with
- maximum spatial field gradient of 5,700 G/cm (57 T/m)
- maximum force product of 102,000,000 G^2/cm (102 T^2/m)
- theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the BeGraft Peripheral Plus Stent Graft is expected to produce a maximum temperature rise of less than

- **2.9°C** (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of ≈ 1.5°C (2 W/kg, 1.5 Tesla)
- **2.1°C** (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of ≈ 0.9°C (2 W/kg, 3 Tesla)

after 15 minutes of continuous scanning.

The effect of heating in the MRI environment for overlapping stent grafts or stent grafts with fractured struts is not known.

In non-clinical testing, the image artifact caused by the device extends approximately 14.6 mm from the BeGraft Peripheral Plus Stent Graft when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

7 Procedure

7.1 Inspection Procedure Prior To Use

Prior to using the BeGraft Peripheral Plus Stent Graft System, carefully remove the system from the package and inspect the catheter for bends, kinks and other damage. Do not use the system if any defects are noted.

7.2 Accessories

- Introducer sheath in recommended minimum size (refer to product label)
- Two to three 10 - 20 ml syringes
- 1000 u / 500 ml Heparinized Normal Saline (HepNS)
- One 0.035" (0.89 mm) guide wire of appropriate length
- 60 % contrast (e.g. with 300 mg Iodine / ml) diluted 1:1 with normal saline
- One inflation device
- One three-way stopcock
- One torque device (if applicable)

7.3 Guide Wire Lumen Flush Procedure

1. In order to avoid loosening of the crimped stent graft from the balloon, hold the protector at its orange marker with one hand and the delivery system with the other hand, then pull the protector off.
2. Verify that the stent graft is located between the radiopaque balloon markers.
3. Attach the syringe with HepNS to the guide wire port.
4. Flush until fluid exits the distal tip.

7.4 Delivery System Preparation Procedure

1. Prepare the inflation device / syringe with diluted contrast medium.
2. Attach the inflation device / syringe to the stopcock; attach to the inflation port.
3. With the tip down, orient the Delivery System vertically.
4. Open the stopcock to the Delivery System; pull negative pressure for 30 seconds; release to neutral pressure for contrast fill.
5. Close the stopcock to the Delivery System; purge the inflation device / syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled.

Warning: If it is not possible to create or to maintain a vacuum, check the connections. If you cannot find the cause, do not use this BeGraft Delivery System under any circumstances and return the unit to Bentley InnoMed GmbH.

7. If a syringe was used, attach a prepared inflation device to the stopcock.
8. Open the stopcock to the Delivery System.
9. Leave on neutral pressure.

7.5 Stent Graft Delivery Procedure

1. Wipe the exposed guide wire with heparinized saline to remove residual blood or contrast medium.
2. Fully open the haemostatic valve. Maintain neutral pressure on the inflation device.
3. Backload the Delivery System onto the proximal portion of the guide wire while maintaining guide wire position across the target lesion / defect.
4. For insertion of the BeGraft Peripheral Plus Stent Graft into the appropriately sized introducer sheath, it is recommended to carefully hold the stent graft at the proximal end (with your fingers) when pushing it into the haemostatic valve and to closely observe if the position of the stent graft on the balloon is affected. If this is the case, this device shall not be used and the unit should be returned to Bentley InnoMed GmbH.

Warning: If a too small introducer sheath is used, the integrity of the stent graft could be impaired (dislodgement of the stent graft and / or the ePTFE graft material), causing an embolism in the worst case.

5. Advance the Delivery System over the guide wire to the target lesion / defect. Utilize the radiopaque balloon markers to position the BeGraft Peripheral Plus Stent Graft across the lesion / defect; perform angiography to confirm the BeGraft Peripheral Plus Stent Graft position. If applicable tighten the haemostatic valve.

Note: If you notice during the process of moving the Delivery System into position that the BeGraft Peripheral Plus Stent Graft has moved on the balloon, do not deploy the BeGraft Peripheral Plus Stent Graft.

Caution: Do not attempt to pull an unexpanded Stent Graft back through the introducer sheath. If an unexpanded Stent Graft needs to be removed, the introducer sheath and the Delivery System should be removed as single unit (see “7.7 Delivery System Removal Procedure”).

6. The BeGraft Peripheral Plus Stent Graft is now ready to be deployed.

7.6 Stent Graft Deployment Procedure

Caution: Refer to the product label for expanded stent outer diameter, nominal deployment pressure (NP) and rated burst pressure (RBP).

PRINT DATE 29.06.2020: Printed documents must be compared with the valid revision before use!
Регистрация МИ в Росздравнадзоре
www.nevacert.ru | info@nevacert.ru

1. Inflate the balloon to the specified nominal deployment pressure (NP). Confirm complete expansion of the stent graft / balloon fluoroscopically. Do not exceed the RBP.
2. Deflate the balloon by pulling negative pressure on the inflation device. Ensure that the balloon is fully deflated before removal of the Delivery System.

7.7 Delivery System Removal Procedure

1. Maintain negative pressure to allow the balloon to remain fully deflated during removal.
2. Withdraw the Delivery System through the sheath with the guide wire remaining across the lesion / defect.

Caution: Should unusual resistance be felt at any time during removal of the Delivery System post BeGraft Peripheral Plus Stent Graft implantation, the introducer sheath and Delivery System should be removed as a single unit. In this case, position the proximal balloon marker just distal to the tip of the introducer sheath, then advance the guide wire in the anatomy as far distally as safely possible and secure the Delivery System to the introducer sheath; then remove the introducer sheath and Delivery System as a single unit. Failure to follow these steps and / or applying excessive force to the Delivery System can potentially result in loss or damage of the Delivery System components.

3. Confirm optimal wall apposition of the BeGraft Peripheral Plus Stent Graft using standard angiography techniques. If the primary angiographic result is not optimal, you can further expand the BeGraft Peripheral Plus Stent Graft with the help of an appropriate balloon catheter.

Note: Once used, this product is potentially biohazardous. Handle and dispose of in accordance with current medical practices and applicable local and/or national laws and regulations.

8 Implant Information Card

An Implant Information Card is supplied with this product. This card must be carefully filled out and the peel-off label of the respective implant, indicating the Manufacturer / Product Name / REF# / Lot#, should be applied to the indicated area on the Implant Information Card (up to 4 peel-off labels can be applied). The patient should carry this Implant Information Card with him and show it to medical personnel in case of further treatments.

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10 Absence of Defectiveness

Bentley InnoMed GmbH manufactured, sterilised and packaged this product with the due care described in the respective laws and regulations for Medical Devices. Each manufacturing lot is individually tested and inspected before it is released for distribution.

As a result of biological differences in individuals, no product is 100% efficient under all circumstances. Therefore it is not possible to reliably predict product efficiency during use, nor is it possible to rule out unsatisfactory results. The product shall not be used if damaged during transport, storage or prior to use. Resterilisation or reuse of the product is prohibited. In case any difficulties arise during the use of the products of Bentley InnoMed GmbH, please send an e-mail to complaint@bentley.global. If the safety of users or patients is affected or if they have been harmed, the statutory notification obligations and deadlines shall be observed.

In the event of complaints or product failure, all components (including accessories required for the intervention such as guide wire, sheath, etc.) must be kept and stored. The manufacturer should be informed immediately after this and the relevant components returned to them.

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
Bentley

BeYond
venous

Venous Self-Expanding Stent System

en	Venous Self-Expanding Stent System
bg	Венозна саморазширяваща се стентова система
cs	Venózní samoexpandibilní systém stentu
de	Venöses Selbstexpandierendes Stentsystem
da	Selvekspanderende venestentsystem
es	Sistema de stent autoexpandible venoso
et	Venoosne iselaienev stendisüsteem
fi	Itselaajeneva laskimostenttijärjestelmä
fr	Système de stent veineux auto-expansible
el	Σύστημα φλεβικού αυτοεκπτυσσόμενου στεντ
hr	Sustav samoširećeg venskog stenta
hu	Vénás öntáguló sztentrendszer
it	Sistema di stent venoso auto-espandente
lt	Savaime išsiskleidžiančio veninio stento sistema
lv	Venozā pašizpletošā stenta sistēma
nl	Veneus zelfexpanderend stentsysteem
no	Venøst selvekspanderende stentsystem
pl	System żylnego stentu samorozprężalnego
pt	Sistema de Stent Venoso Auto-expansível
ro	Sistemul de stent venos autoexpandabil
ru	Система самораскрывающегося венозного стента
sk	Systém samoexpandujúceho venózneho stentu
sl	Sistem samoraztegljive venske žilne opornice
sv	Venöst självexpanderande stentsystem
tr	Kendiliğinden Açılabilen Venöz Stent Sistemi

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
English

Note:

These instructions are intended for all personnel who will handle the "BeYond Venous Self-Expanding Stent System" and should be read carefully before using this product!

BeYond Venous Self-Expanding Stent System

1 Product Description

The "BeYond Venous Self-Expanding Stent System" consists of

- a "BeYond Venous Self-Expanding Stent", loaded in an over-the-wire (OTW), 0.035" compatible, stent delivery system with a useable length of 100 cm and a shaft diameter of 3.3 mm (10 Fr).

The stent is a self-expanding, flexible, slotted-tube nitinol stent and achieves its unconstrained diameter upon deployment into the target vessel. The stent is available in 5 diameters (Ø 10, 12, 14, 16 and 18 mm) each in lengths of 60, 80, 100, 120 and 150 mm.

One radiopaque marker (located at the catheter tip) fluoroscopically identifies the distal end of the loaded stent. The proximal end of the loaded stent can be identified by the end of the metallic middle tube.

2 How Supplied

Contents per Package:

- One "BeYond Venous Self-Expanding Stent System" (placed in a blister with lid)
- One Implant Information Card
- One Instructions for Use

Accessories as listed in chapter 7.2 are not included.

Storage: Store in a dry, dark and cool place. Do not use the "BeYond Venous Self-Expanding Stent System" after the "Use by" date (indicated on the product label).

Non Pyrogenic.

Sterile: Sterilized with ethylene oxide (EO). Do not use if the package is open or damaged.

Note: This single use device cannot be reused on another patient, as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and/or chemical characteristics introduced under conditions of repeated use, cleaning, and/or re-sterilization may compromise the integrity of the design and/or materials, leading to contamination due to narrow gaps and/or spaces and diminished safety and/or performance of the device. Absence of original labelling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and/or user.

3 Indications

The "BeYond Venous Self-Expanding Stent System" is indicated for endovascular treatment of patients with symptomatic obstructions of the femoral vein or iliac vein, e.g.

- acute iliofemoral deep vein thrombosis,
- post-thrombotic syndrome from prior iliofemoral deep vein thrombosis,
- chronic venous insufficiency from non-thrombotic iliofemoral vein compression such as the May Thurner syndrome,
- chronic venous insufficiency from non-thrombotic iliofemoral vein compression such as compression from ectatic or tortuous iliac arteries and
- chronic venous insufficiency from cancer-related iliofemoral vein compression.

4 Contraindications

Contraindications include, but may not be limited to:

1. When PTA is technically not possible (e.g. not feasible to access the lesion with guide wire or balloon catheter)
2. Total venous occlusion without possibility of sufficient dilation for subsequent passage of the "BeYond Venous Self-Expanding Stent System"
3. Patients with known hypersensitivity to the stent material (nickel titanium alloy)
4. Patients with uncontrolled coagulation disorders
5. Patients with allergy to anticoagulants or antiplatelet agents or bleeding complications and patients unable or unwilling to tolerate anticoagulant/antiplatelet therapy and/or non-responders to anticoagulant/antiplatelet therapy
6. Patients with active bleeding events or with an increased risk of bleeding complications which cannot be treated with anticoagulants or antiplatelet agents
7. Patients with sepsis or acute bacterial infections to avoid the risk of stent-related infections
8. Patients with ongoing intravenous drug abuse to avoid the risk of stent-related infections or thrombosis

Common contraindications regarding PTA of vascular lesions are also applicable.

5 Potential Adverse Events

Potential adverse events include, but may not be limited to the following:

- Access site complications including: bleeding, pain, pseudoaneurysm, hematoma, nerve or vessel damage, infection or arteriovenous fistula formation
- Acute thrombosis in cases of lack of perfusion (e.g. spasmus or other occlusions) of the passage vessels and insufficient anticoagulation
- Cerebrovascular dysfunction and/or stroke
- Death
- Embolism/pulmonary embolism
- Emergency vessel surgery
- Inflammation, fever, septicemia, bacteremia
- Hemorrhage
- Hypersensitivity/allergic reaction (e.g. to stent or catheter material and/or contrast medium)
- Hypotension or Hypertension
- Infection
- Myocardial infarction (MI) and/or other cardiovascular dysfunctions (e.g. ischemia, angina, arrhythmia)
- Organ failure (e.g. renal failure)
- Pain
- Pericardial tamponade
- Ischemic events (perfusion limitations) in lateral branches of arteries, veins or distal vessels
- Pseudoaneurysm formation
- Respiratory arrest (apnoea)
- Restenosis
- Retinal or retroperitoneal hemorrhage related to anticoagulant therapy
- Reperfusion damage including cardiac or pulmonary dysfunction caused by aprupt change to haemodynamics, including cardiac arrest and pulmonary failure or pulmonary embolism
- Thrombophlebitis/thrombosis of the stented vein/stent occlusion
- Tissue damage (e.g. necrosis)
- Vasoconstriction (e.g. vasospasm)
- Vascular damage including dissection, perforation and rupture

6 Warnings and Precautions

The medical techniques and procedures described below do not entail all the PTA procedures acknowledged in medicine.

- Any use of the "BeYond Venous Self-Expanding Stent System" for procedures other than those indicated in this Instructions for Use is not recommended.
- The "BeYond Venous Self-Expanding Stent System" should only be used by physicians trained in interventional techniques and percutaneous transluminal angioplasty.
- Insert the "BeYond Venous Self-Expanding Stent System" under fluoroscopic guidance.
- If a long lesion needs to be stented consider using the longest available stent rather than overlapping stents. If multiple stents are placed in an overlapping fashion, they should be of similar composition (nitinol) and the overlapping zone should be at least 10 to 15 mm.
- Administer appropriate anticoagulation/antiplatelet treatment before insertion of the "BeYond Venous Self-Expanding Stent System". Continue the anticoagulation/antiplatelet treatment after the implantation according to current versions of national guide lines and international standards.
- Placing the "BeYond Venous Self-Expanding Stent" across a bifurcation/sidebranch will compromise the blood flow and impact future diagnostic or therapeutic procedures.
- A stent previously implanted proximally to the lesion can be deformed or dislodged by contact, e.g. with the guide wire or the "BeYond Venous Self-Expanding Stent System".
- The use of laser devices, mechanical atherectomy and hyperthermia are not recommended in the stented region.
- Do not use the system after the "Use by" date (indicated on the product label).
- Prior to use, carefully inspect the product to verify that it has not been damaged during shipment (e.g. bends, kinks or other damages).
- Do not use if the package is open or damaged.
- Employ aseptic techniques during removal of the product from the package and during use.
- Do not expose this device to organic solvents, ionizing radiation or ultraviolet light.
- Use the recommended minimum introducer sheath size (refer to label!).
- Ensure that the selected stent dimensions are suitable for the intended procedure. Undersizing the stent diameter can lead to stent migration.
- Do not attempt to straighten the catheter if the shaft is kinked. Instead, use a new product.
- Do not use the product if flushing cannot be performed properly.
- Prior to deployment of the "BeYond Venous Self-Expanding Stent", utilize fluoroscopy to verify that the loaded stent within the delivery system has been properly positioned, otherwise do not deploy the stent.
- Ensure to keep the introducer sheath in place during stent deployment.
- For stent deployment, never push the proximal handle ⑨ forward! This could lead to improper deployment of the stent resulting in major complications (e.g. loss of stent integrity).
- Should unusual resistance be felt at any time during lesion access or removal of the catheter, the introducer sheath and catheter should be removed as a single unit. If possible, retain the guide wire position for subsequent vessel access. Applying excessive force to the delivery system might lead to loss or damage of the stent and/or the components of the "BeYond Venous Self-Expanding Stent System".
- If exessively high force is needed to start the deployment of the stent and the stent is still completely inside the delivery system remove the delivery system and replace with a new unit. If deployment of the stent has already started, fully deploy the stent.
- Repositioning of the stent is not possible! Once deployment has started, the stent must be fully deployed, since the outer tubing of the shaft cannot be re-advanced over the partly deployed stent.
- Re-crossing a "BeYond Venous Self-Expanding Stent" with adjunct devices must be performed with extreme caution to avoid damage, deformation or displacement of the stent.
- Refer to balloon catheter's instructions for use for information, warnings and precautions concerning the pre-/post-dilation procedure.
- The safety and effectiveness of this device for use in the arterial system have not been established.
- The long-term outcomes following repeat dilation of endothelialized stents are unknown.
- The long-term effectiveness of the "BeYond Venous Self-Expanding Stent" has not been fully established.

6.1 MRI Safety and Compatibility

Non-clinical testing has demonstrated that the **"BeYond Venous Self-Expanding Stent"**, in single and in overlapped configuration, is **MR Conditional** as defined in ASTM F2503. A patient with this device can be scanned safely immediately after implantation under the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0 Tesla, with
- Maximum spatial field gradient of 128 T/m (12,800 G/cm)
- Maximum force product of 231 T²/m
- Theoretically estimated maximum whole body averaged specific absorption rate (WBA-SAR) of 2 W/kg for 15 minutes of continuous scanning (Normal Operating Mode).

The "BeYond Venous Self-Expanding Stent" shows a magnetically induced **displacement force** (ASTM F2052- 14) of 2% (1° deflection) and **torque** (ASTM F2213-06) of 1% and should not migrate in this MRI environment.

3.0 Tesla RF-induced heating in non-clinical testing on overlapped configuration (performed according to ASTM F2182- 11) produced a maximum in-vitro temperature rise after 15 minutes of **2.7°C** when scaled to a WBA-SAR of 2 W/kg. Tests were conducted in a 3.0 T Siemens Magnetom Trio with "Numaris/4 syngo MR B17" software. These calculations do not take into consideration the cooling effects of the blood flow.

1.5 Tesla RF-induced heating in non-clinical testing on overlapped configuration (performed according to ASTM F2182- 11) produced a maximum in-vitro temperature rise after 15 minutes of **7.1°C** when scaled to a WBA-SAR of 2 W/kg. Tests were conducted in a 1.5 T Philips Intera with "Release 12.6.1.4, 2012-05-22" software. These calculations do not take into consideration the cooling effects of the blood flow.

RF-induced heating general: The effects of MRI on overlapping stents (other than the tested worst case configuration of 130 mm) or stents with fractured struts are unknown. No other RF heating testing than 1.5 and 3 Tesla was performed.


MR-artefacts: In non-clinical testing, the image artefact caused by the device extends approximately 6.7 mm from the stent when imaged with a gradient echo pulse sequence and a 3 Tesla MR system. Therefore, it may be necessary to optimize the MR imaging parameters in the presence of "BeYond Venous Self-Expanding Stents".

7 Procedure

7.1 Inspection Procedure Prior To Use

Prior to using the "BeYond Venous Self-Expanding Stent System", carefully remove the system from the tray and inspect the catheter for bends, kinks and other damage. Do not use the system if any defects are noted. Verify that the safety lock ⑦ is still tightened. Do not use if the stent ③ is already partially deployed. Ensure that there is no gap between tip ① and shaft ④. If there is a gap, open the safety lock ⑦ (approx. two rotations), move the distal handle ⑤ carefully towards the tip until the gap is closed and tighten the safety lock again.

Warning: Do not attempt to straighten the catheter if the shaft is kinked. Instead, use a new system.

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7.2 Accessories

In addition to the "BeYond Venous Self-Expanding Stent System", the following standard materials may also be required to facilitate delivery and deployment of the stent:

- Introducer sheath of minimally 10 Fr (≥ 3.3 mm)
- Guide wire of maximally 0.035" (≤ 0.89 mm) and of appropriate length
- Sterile syringe for flushing (≥ 5 ml)
- 1000 u/500 ml Heparinized Normal Saline (HepNS)
- Intravenous contrast medium diluted with normal, sterile saline
- Appropriate balloon dilation catheter for pre- and post-dilation
- Appropriate anticoagulation drugs

7.3 Pre-Dilation

Consider pre-dilation of the target vessel and pre-dilate if necessary. Remove the balloon catheter, leaving the guide wire in place.

7.4 Stent Size Selection

Determine the length of the target lesion to identify the appropriate length of stent(s) required. Ensure that the stent is long enough to fully cover the lesion thereby extending the lesion at both sides sufficiently.
Identify the reference vessel diameter proximally and distally to the lesion (after potential pre-dilation). For appropriate sizing refer to the "Stent Diameter Selection Table" below.

Stent Diameter Selection Table

Reference Vessel Diameter	Unconstrained Stent Outer Diameter
8 - 9 mm	10 mm
10 - 11 mm	12 mm
12 - 13 mm	14 mm
14 - 15 mm	16 mm
16 - 17 mm	18 mm

Warning: If a long lesion needs to be stented consider using the longest available stent rather than overlapping stents. If multiple stents are placed in an overlapping fashion, they should be of similar composition (nitinol) and the overlapping zone should be at least 10 to 15 mm.
Warning: Ensure that the selected stent dimensions are suitable for the intended procedure. Undersizing the stent diameter can lead to stent migration.

7.5 Guide Wire Lumen and Delivery System Flushing Procedure

Attach the syringe with HepNS to the guide wire port on the proximal handle ⑨. Flush until HepNS exits the distal tip ①.
Attach the syringe with HepNS to the manifold ⑥. Flush until HepNS exits at the gap between tip ① and shaft ④.

7.6 Stent Delivery Procedure

The intervention must be performed under fluoroscopic guidance.
Wipe the exposed guide wire with HepNS to remove residual blood or contrast medium.
Wipe the usable length portion of the shaft with gauze soaked with HepNS.
Warning: Placing the "BeYond Venous Self-Expanding Stent" across a bifurcation/sidebranch will compromise the blood flow and impact future diagnostic or therapeutic procedures.
Backload the delivery system onto the guide wire while maintaining the guide wire's position across the target lesion.
Advance the delivery system over the guide wire as straight as possible and hold the introducer sheath while introducing the delivery system. In order to verify the position of the loaded stent relative to the lesion, perform X-ray and utilize the radiopaque marker ② at the tip ① and the end of the metallic (radiopaque) middle tube ③. In order to fully cover the lesion, the stent should extend the lesion sufficiently on both sides.
Note: This defines the location of the implanted stent.
The stent is now ready to be deployed and the positions of all components should not be changed any more.

7.7 Stent Deployment Procedure

During stent deployment ensure to keep the introducer sheath in place. Before deployment of the stent, it is important to straighten the shaft ④ of the delivery system as much as possible without further moving the device.
Fully unscrew the safety lock ⑦ by turning it counter-clockwise for at least two rotations. To keep the already positioned stent exactly in place, carefully move the proximal handle ⑨ forward under fluoroscopic control until the end of the metallic middle tube ③ reaches the stent. In order to deploy the stent, hold the proximal handle ⑨ stationary and pull back the distal handle ⑤ in a slow, smooth and consistent fashion until the proximal handle is reached. Use fluoroscopic control in order to verify the stent is fully deployed.

Note: The deployment force might be higher than expected.
Warning: If excessively high force is needed to start the deployment of the stent and the stent is still completely inside the delivery system remove the delivery system and replace with a new unit. If deployment of the stent already started, fully deploy the stent.
Warning: For stent deployment, never push the proximal handle ⑨ forward! This could lead to improper deployment of the stent resulting in major complications (e.g. loss of stent integrity).
Warning: Repositioning of the stent is not possible! Once deployment has started, the stent must be fully deployed, since the outer tubing of the shaft cannot be re-advanced over the partly deployed stent.

7.8 Delivery System Removal Procedure

Withdraw the delivery system through the sheath with the guide wire remaining across the lesion. Check for delivery system integrity after removal from the patient.
Note: The delivery system can be removed without the need to recapture the tip with the outer tubing of the shaft.
Note: Should unusual resistance be felt at any time during lesion access or removal of the delivery system post stent implantation, the introducer sheath and catheter should be removed together as a unit.
Confirm optimal wall apposition of the stent using standard angiography techniques. If the stent is incompletely expanded at any point along the lesion, post-dilation with a balloon catheter can be performed at the discretion of the physician. Select an appropriately sized balloon catheter and post-dilate the stent and the lesion with conventional technique.
Warning: The inflation diameter of the balloon used for post-dilation shall not exceed the unconstrained stent diameter.
Note: Refer to balloon catheter's instructions for use for information, warnings and precautions concerning the post-dilation procedure.
Note: Once used, this product is potentially biohazardous. Handle and dispose of in accordance with current medical practices and applicable local and/or national laws and regulations.

8 Implant Information Card

An Implant Information Card is supplied with this product. This card must be carefully filled out and the peel-off label of the respective implant(s), indicating the Manufacturer/Product Name/REF#/Lot#, should be applied to the indicated area on the Implant Information Card (up to 4 peel-off labels can be applied). The patient should carry this Implant Information Card with him and show it to medical personnel in case of further treatments.

9 References

The physician should consult up-to-date literature on current medical practice of PTA.

10 Absence of Defectiveness

Bentley InnoMed GmbH manufactured, sterilised and packaged this product with the due care described in the respective laws and regulations for Medical Devices. Each manufacturing lot is individually tested and inspected before it is released for distribution.
As a result of biological differences in individuals, no product is 100% efficient under all circumstances. Therefore it is not possible to reliably predict product efficiency during use, nor is it possible to rule out unsatisfactory results. The product shall not be used if damaged during transport, storage or due to improper handling by the operator. Resterilisation or reuse of the product is prohibited. Bentley InnoMed GmbH shall not be liable for any damages arising from any off-label use.

Important Information:
If more than 2 years have elapsed between the revision date stated on the back side of this instruction and today's date, please contact the manufacturer in order to find out if any changes have occurred.

For further information please contact:

Bentley InnoMed GmbH
Lotzenäcker 3
72379 Hechingen
Deutschland (Germany)
T + 49 (0) 7471 984 995 10
F + 49 (0) 7471 984 995 9
info@bentley.global
www.bentley.global
© 2019 Bentley InnoMed GmbH
All rights reserved.

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