Tyvek® Pouch with STERRAD™ Chemical Indicator (Heat Seal (125XX) and Self Seal (123XX)) Tyvek® Roll with STERRAD™ Chemical Indicator (124XX)

ADVANCED STERILIZATION PRODUCTS, INC.

33 Technology Drive, Irvine, CA 92618, USA
+1(888)-STERRAD www.asp.com

EC REP ASP, The Netherlands BV
BIC 1 5657 BX Eindhoven, The Netherlands

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Tyvek® Pouch with STERRAD™ Chemical Indicator (Heat Seal (125XX) and Self Seal (123XX)) Tyvek® Roll with STERRAD™ Chemical Indicator (124XX)

en

INSTRUCTIONS FOR USE

INDICATIONS FOR USE:

Tyvek® Pouches and Rolls with STERRAD™ Chemical Indicator are intended to be used to enclose medical devices that are to be terminally sterilized in the STERRAD Sterilization Systems and Cycles, in a single or double pouch configuration, and to indicate, by color change, that the pouch has been exposed to sterilant. After completion of the sterilization process, the pouch/roll is intended to maintain sterility of the enclosed device until used.

The pouches and rolls are printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2014 [Type I]) that changes from red to the color (or lighter) indicated on the comparator bar included on the shelf pack when exposed to hydrogen peroxide vapor during processing in the STERRAD Sterilization Systems.

The pouches and rolls can be used in the STERRAD Sterilization Systems and Cycles shown in Table 1.

Table 1: Indications for Use

STERRAD / Cycle		Single-channel stainless steel lumens	Other lumens	Load requirements
STERRAD 100S		≥1 mm ID x ≤125 mm long ≥2 mm ID x ≤250 mm long ≥3 mm ID x ≤400 mm long	Metal and nonmetal ≥6 mm ID x ≤310 mm long	≤10 lumens per load
STERRAD NX™ with or without ALLClear™ Technology	STANDARD	≥1 mm ID x ≤150 mm long ≥2 mm ID x ≤400 mm long	N/A	≤10 lumens per load ≤10.7 lbs
	ADVANCED	≥1 mm ID x ≤500 mm long	N/A	≤10 lumens per load ≤10.7 lbs
STERRAD 100NX with or without ALLClear™ Technology	STANDARD	≥0.7 mm ID x ≤500 mm long	N/A	≤10 lumens per load ≤21.4 lbs**
	EXPRESS	N/A	N/A	≤10.7 lbs on bottom shelf only
	DUO	N/A	N/A	Accessory devices (e.g., cameras and light cords) only 2 cameras maximum

**The validation testing for this lumen size was conducted using a maximum of 10 lumens in the USA; 40 lumens for markets outside the USA. Your loads should not exceed the validated maximum number of lumens.

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DESCRIPTION:

1. The pouches and rolls are constructed from two films: a transparent, non-permeable film, and an opaque, permeable film of Tyvek material.
2. Real-time testing supports maintenance of post-processing pouch/roll sterility for a period of 12 months.
3. The Product Codes for all Tyvek Pouches and Rolls with STERRAD Chemical Indicator are listed in Table 2.

Table 2: Product Codes for Tyvek Pouches and Rolls with STERRAD Chemical Indicator

Product Type	Product Codes
Tyvek Self-Seal Pouches with STERRAD Chemical Indicator (123XX)	12320, 12326, 12335, 12332, 12342, 12340, 12348, 12356
Tyvek Rolls with STERRAD Chemical Indicator (124XX)	12407, 12410, 12415, 12420, 12425, 12435, 12442, 12450
Tyvek Heat-Seal Pouches with STERRAD Chemical Indicator (125XX)	12521, 12526, 12532, 12541, 12548, 12543*, 12544*, 12557*, 12558*, 12559*, 12560*

**Indicates Product Codes exclusively for markets outside the USA.*

DIRECTIONS FOR USE:

1. Place medical devices(s) through an open end of pouch/roll. If using a self-seal pouch, peel off self-seal tape and seal the open end of the pouch. If using a heat-seal pouch, place in heat sealer following directions listed in the "NOTE" directly following. For roll, cut to desired length and seal open ends following heat sealer instructions. **NOTE:** When heat sealing the Tyvek pouch/roll for the first time, the operator should establish the seal temperature/dwell time range for their specific machine according to the peel strength required when opening the package. Typical conditions for heat sealer vary by machine, but are typically between 121°C-135°C (250°-275°F), and generally require 1 to 3 seconds of seal time. Seal temperature/time used should be recorded for future reference. Consult heat sealer instructions for use for proper operation.
2. Refer to the comparator bar on the shelf pack for the proper post-processing color change. The chemical indicator changes from red to the color (or lighter) indicated on the comparator bar when exposed to hydrogen peroxide vapor during processing in the STERRAD Sterilization Systems.
3. Correct color change of the chemical indicator bar indicates exposure to hydrogen peroxide, but does not indicate sterility.
4. The pouches and rolls are printed with a chemical indicator bar which is a process indicator (ISO 11140-1:2014 [Type I]) that changes from red to the color (or lighter) indicated on the comparator bar included on the shelf pack when exposed to hydrogen peroxide vapor during processing in the STERRAD Sterilization Systems. Real-time testing has demonstrated stability of the post-processing chemical indicator color change for up to 6 months.

STORAGE AND DISPOSAL:

1. Store at ambient temperature, 15°C-30°C (59°F-86°F), under dry conditions. Avoid exposing the chemical indicator side of pouch/roll to light.
2. Before and after processing, store away from STERRAD Sterilizer or sources of hydrogen peroxide, acids, alkalis, or volatile antimicrobial agents such as glutaraldehyde or formaldehyde.
3. You may dispose this product in accordance with your National, Federal, State and Local regulations or per healthcare facility's policies and procedures.

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⚠ **CAUTION:**

1. After processing, the chemical indicator will not be suitable for use as permanent record.
2. If the chemical indicator bar has not completely changed from red to color (or lighter) indicated by comparator bar, introduction of hydrogen peroxide, an essential part of the STERRAD Sterilization cycle, may not have taken place.
3. Chemical indicator bar offer an additional way to verify processing in the sterilization cycle and should be used in addition to, not in place of, the biological indicator. A color change on the Chemical indicator bar does not signify sterilization; they solely confirm that the indicator has been exposed to the hydrogen peroxide.
4. Do not use after expiration date.
5. Single use only. Do not reuse as the pouches and rolls will not function properly if reused. Should a cycle cancellation occur, repackage the devices using new pouches and rolls when restarting the cycle.
6. Do not use if pouch or roll is damaged.
7. Users should perform their own determination for their heat sealer parameters.
8. Tyvek pouches may not be appropriate for some delicate devices requiring protective packaging such as flexible endoscopes. Always refer to the device manufacturer's instructions for use for special packing instructions.

Tyvek is a DuPont registered trademark.

For a glossary of symbols used, go to eifu.asp.com

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© ASP 2021