Amplatzer™ Multifenestrated Septal Occluder - “Cribriform”

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Amplatzer™ Multifenestrated Septal Occluder - “Cribriform”

EN: Instructions for Use

CAUTION: Read all instructions carefully. Failure to follow these instructions, warnings, and precautions may lead to patient injury or occluder damage.

Device Description

The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is a self-expanding occluder made from braided nitinol designed for occluding multi-fenestrated atrial septal defects. The occluder has two discs, the left atrial disc and right atrial disc that are linked together by a short connecting flexible waist. The discs are designed to cover multiple fenestrations in the atrial septum and the waist is intended to be placed through the most central defect in relationship to the multiple fenestrations. The occluder is filled with polyester patches that are securely sewn to the right atrial disc and the left atrial disc by a polyester thread. The occluder consists of radiopaque marker bands at each end of the device, one attached to the left atrial disc and one attached to the end screw on the right atrial disc. The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is sterile and for single use.

Refer to the figure and tables in Appendix A: Supplemental Information for more information about the occluder. Occluder size selection and delivery system compatibility and dimensions are provided in Table 1 and Table 2, respectively. The following components are illustrated in Figure 1:

Figure 1 in Appendix A: Supplemental Information

1 Left Atrial Disc
2 Waist
3 Right Atrial Disc
4 Marker Band
5 End Screw and Marker Band

Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” Performance Characteristics

The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is a self-expanding occluder made of braided nitinol wire. The occluder has a waist and two retention discs that stabilize the occluder across the atrial septum and designed to cover multiple fenestrations in the atrial septum. Polyester patches are securely sewn within each disc to promote occlusion. The braided nitinol wire and polyester patches promote occlusion and tissue ingrowth. Radiopaque marker bands provide visual guidance under fluoroscopy, aiding in precise occluder positioning. The end screw allows the attachment, recapture, and release of the occluder, enabling precise placement within the anatomy.

Intended Purpose

The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is a device intended for the occlusion of a multifenestrated atrial septal defect.

Indications

The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is indicated for patients with an ostium secundum atrial septal defect.

Target Patient Population

Patients with a multifenestrated ostium secundum atrial septal defect with clinical evidence of right ventricular (RV) volume overload (such as 1.5:1 degree of left-to-right shunt or RV enlargement on echocardiography)

Intended User

Physicians trained in standard transcatheter techniques

Intended Clinical Benefit

The following benefits are expected to be achieved:

- Minimally invasive alternative to surgical closure associated with lower complication rates
- Reduced hospital stay compared to surgery
- Immediate atrial septal defect closure with subsequent alleviation of associated morbidities

Summary of Safety and Clinical Performance

A summary of the safety and clinical performance (SSCP) for this occluder is available at https://ec.europa.eu/tools/eudamed. This is the SSCP location after the launch of the European Database on Medical Devices/Eudamed. Search for the occluder using the Basic UDI-DI: 5415067AMP1100DV.

Limitations

No known limitations

Contraindications

The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is contraindicated for:

- patients known to have extensive congenital cardiac anomaly, which can only be adequately repaired by way of cardiac surgery.
- patients with active endocarditis or other systemic infections requiring treatment.
- patients who are unable to tolerate intra-procedural anticoagulation or post-procedural anti-platelet therapy.
- patients whose vascular access site is inadequate to accommodate appropriate sheath size.
- patients known to have demonstrated intracardiac thrombi (especially left atrial or left atrial appendage thrombi).
- patients whose size or condition (for example, too small for transesophageal echocardiography [TEE] probe, catheter size, vasculature size, active infection) would cause the patient to be a poor candidate for cardiac catheterization.
- patients for whom the radius of the occluder is greater than the distance from the central defect to the aortic root or superior vena cava.

Warnings

- The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” consists of a nickel-titanium alloy, which is generally considered safe. However, in vitro testing has demonstrated that nickel is released from this occluder for a minimum of 60 days. Patients who are allergic to nickel may have an allergic reaction to this occluder, especially those with a history of metal allergies. Certain allergic reactions can be serious; patients should be instructed to notify their physicians immediately if they suspect they are experiencing an allergic reaction such as difficulty breathing or inflammation of the face or throat. Some patients may also develop an allergy to nickel if this occluder is implanted. For more information, see Product Materials.
- DO NOT use the Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” after the Use-by date stated on the package label.
- Do not use this occluder if the sterile package is open or damaged.
- This occluder was sterilized with ethylene oxide and is for single use only. Do not reuse or resterilize this occluder. Attempts to resterilize this occluder can cause a malfunction, insufficient sterilization, or harm to the patient.

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- Physicians must be prepared to deal with urgent situations such as pericardial effusion or occluder embolization that can require the removal of the occluder.
- Remove embolized occluders. Do not remove an embolized occluder through intracardiac structures unless the occluder is fully recaptured inside a catheter.
- Do not use a power injection syringe to inject contrast solution through the sheath.
- Implantation of this device may not replace the need for anticoagulation therapy in patients with an atrial septal defect and paradoxical embolism.
- The use of transthoracic, transesophageal, or intracardiac echocardiographic imaging (TTE, TEE, or ICE) is required before, during, and after implantation. Imaging is mandatory as an aid in positioning the sheath through the most suitable defect to achieve complete closure of the defect. If TEE is used, the patient's esophageal anatomy must be adequate for the placement and manipulation of the TEE probe.
- Do not release the device from the delivery cable if the device does not conform to its original configuration, or if the device position is unstable, or if the device interferes with any adjacent cardiac structure, such as superior vena cava (SVC), pulmonary vein (PV), mitral valve (MV), coronary sinus (CS) or aorta (AO). Recapture the device and redeploy. If still unsatisfactory, recapture the device and replace with a new device.
- Placement of the Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” may impact future cardiac interventions, for example, transseptal puncture and mitral valve repair.

Precautions

- Use standard transcatheter techniques when using Amplatzer™ products.
- This device should be used only by physicians who are trained in standard transcatheter techniques. The physician should determine which patients are candidates for procedures that use this device.
- The physician should exercise clinical judgment in situations that involve the use of anticoagulants or antiplatelet drugs before, during, and/or after the use of this device.
- Prolonged procedures may result in increased exposure to anesthesia, contrast media, and/or radiation.
- Patients should be fully heparinized throughout the procedure using adequate dosing so as to keep the active clotting time (ACT) greater than 200 seconds prior to device insertion.
- Do not use contrast power injection with delivery sheath.
- Store in a dry place.

Patient Selection

- Certain patients may be at higher risk for complications such as tissue erosion and device embolization. If higher-risk patients have devices implanted, closer follow-up is warranted (see Post-Procedure Care). Higher-risk patients include the following:
 - Patients with deformation of the device at the aortic root.
 - Patients with high defects (minimal aortic and superior rims).
 - Patients with less than a 9 mm distance from the central defect to the aortic root or superior vena cava orifice.

Special Populations

- Pregnancy: The safety and effectiveness of this occluder has not been established during pregnancy. Fluoroscopic X-ray guidance is used during placement of the occluder. The risk of increased X-ray exposure for patients who are pregnant must be weighed against the potential benefits of this technique.
- Nursing mothers: The safety and effectiveness of this occluder has not been established in lactating mothers. There has been no quantitative assessment for the presence of leachables in breast milk.

MR Conditional

Non-clinical testing has demonstrated that the Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is MR Conditional. A patient with the Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” can be safely scanned in an MR system under the following conditions:

- Static magnetic field of 1.5 Tesla (1.5T) or 3.0 Tesla (3.0T)
- Maximum spatial gradient field of 19 T/m (1900 G/cm)
- Maximum MR system reported, whole-body averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode)

Under the scan conditions defined above, the occluder is expected to produce a maximum temperature rise of less than or equal to 4°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the occluder extends radially up to 1.1 cm and 1.3 cm (respectively) from the occluder when imaged with a gradient echo pulse sequence in a 1.5T MR system and a gradient echo pulse sequence in a 3.0T MR system.

Potential Adverse Events

Potential adverse events that may occur during or after a procedure placing this Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” include, but are not limited to:

- Air embolus
- Allergic reaction / toxic effect due to contrast media, medication, or metal
- Anesthesia reaction (nausea, vasovagal reaction, confusion/altered mental state or other)
- Apnea
- Arrhythmia
- Arteriovenous fistula
- Bleeding
- Brachial plexus injury
- Cardiac perforation
- Cardiac tamponade
- Cardiovascular injury
- Chest pain
- Death
- Device embolization
- Device malformation
- Device migration
- Dissection
- Endocarditis
- Groin pain
- Headache/migraine
- Heart Block
- Hematoma
- Hemolysis
- Hypertension
- Hypotension
- Infection
- Myocardial infarction
- Perforation of vessel or myocardium
- Pericardial effusion
- Peripheral embolism
- Peripheral pulse loss
- Phrenic nerve injury
- Pleural effusion
- Pseudoaneurysm
- Residual shunt
- Stroke/transient ischemic attack
- Thromboembolism
- Thrombus

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- Erosion
- Fever
- Foreign material embolic event
- Valve damage
- Valvular insufficiency
- Vascular access site injury

Compatibility

The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is compatible with 45° models of the Amplatzer™ TorqVue™ Delivery System, the Amplatzer™ Trevisio™ Intravascular Delivery System, and the Amplatzer™ TorqVue™ Exchange System.

Device Lifetime

The Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is expected to remain implanted permanently. The device lifetime for the Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” is at least 10 years, as demonstrated by bench testing of the occluder.

Patient Identification Card

A patient identification card is present within the packaging of the occluder. Hospital personnel should complete the front side of the Patient Identification Card by printing the patient’s name, implant date, physician’s name, healthcare facility information, and phone number in the designated areas. A peel-away device label is provided in the sterile packaging, which includes occluder information (model number, lot number, UDI number, and device name). Apply the peel-away device label to the back of the card in the designated area. If a peel-away label is not available, write the information in the designated area on the back of the card.

NOTE: Write the occluder lot number in the SN/LOT field.

After the card is completed, separate the top portion from the Patient Identification Card by tearing along the perforated line. Provide the completed Patient Identification Card to the patient. Instruct the patient to carry the Patient Identification Card at all times and to present the card to hospital staff in emergency situations. To obtain a replacement card if a patient loses or damages their card, contact Abbott Medical Technical Support. Once the replacement card is received, handwrite the required information in the respective areas and provide the new card to the patient.

Product Materials

The following product materials are intended to come into contact with tissue:

Component	Material description	Substance	Maximum substance concentration in each component (% w/w)	Maximum substance mass based on largest occluder (milligrams)
Braided Nitinol Wire	Nitinol	Nickel	<57	600.7
		Titanium	<43	453.2
		Carbon	<0.1	1.1
		Cobalt	<0.1	1.1
		Copper	<0.1	1.1
		Chromium	<0.1	1.1
		Hydrogen	<0.1	1.1
		Iron	<0.1	1.1
		Niobium	<0.1	1.1
		Nitrogen	<0.1	1.1
		Oxygen	<0.1	1.1
End Screw	316L / 316 LVM Stainless Steel	Iron	<69	22.9
		Chromium	<18	6.0
		Nickel	<14	4.7
		Molybdenum	<3	1.0
		Manganese	<2	0.7
		Silicon	<1	0.4
		Nitrogen	<0.1	0.1
		Carbon	<0.1	0.1
		Phosphorus	<0.1	0.1
		Sulfur	<0.1	0.1
		Cobalt	<0.4	0.2
Marker Band	Platinum / Iridium Alloy	Platinum	80	20.8
		Iridium	20	5.2
Polyester Thread	Polyester	Polyester Blend	100	2.4
Polyester Patch	Polyester	Polyethylene Terephthalate	100	33.8

CAUTION: One or more components of this occluder may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight (w/w): Cobalt: Chemical Abstracts Service (CAS) No. 7440-48-4; EC No. 231-158-0.

Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

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Directions for Use

All figures and tables referenced in this section are located in Appendix A: Supplemental Information.

Materials recommended for use with the Amplatzer™ Multifenestrated Septal Occluder - "Cribriform":

- 30–60 cc luer-lock syringe
- 0.035-inch guidewire

Pre-Procedure Care

It is recommended that at least 24 hours prior to the procedure that the patient receives aspirin (81 mg or 325 mg) or an alternative antiplatelet/anticoagulant if the patient has an aspirin intolerance.

General instructions for the Amplatzer™ Multifenestrated Septal Occluder - "Cribriform" are provided below:

Procedure

1. Prepare the patient for a standard transcatheter procedure. Maintain a recommended activated clotting time (ACT) of greater than 200 seconds prior to device insertion. Effective anticoagulation therapy should be maintained during the procedure as determined by the physician.
2. Puncture the femoral vein and perform a standard right-heart catheterization.

 CAUTION: When placing the occluder using a compatible delivery system, refer to the instructions for use provided with the delivery system.
3. Perform an angiogram to demonstrate the atrial defect. Catheterize the left atrium using a 45° LAO position and cranial angulation of 35°–45°. Inject contrast medium into the right upper lobe pulmonary vein.
4. Introduce the guidewire to gain access to the left atrium and then into the left upper pulmonary artery.
5. Occluder size and placement (the fenestration in which the occluder will be placed) is based on the location of the fenestrations. The use of transesophageal echocardiography or intra cardiac echocardiography is mandatory in positioning the sheath in the most centrally located fenestration. Once the most centrally located defect has been crossed, measurements should be made to determine the distance to the outer rim of the furthest fenestration and the occluder size selected should be 2:1 of this measurement. The occluder selected must completely cover all fenestrations.
6. Transesophageal echocardiography (TEE) or similar imaging equipment (in other words, intracardiac echocardiography) is required to measure the distance from the most centrally located defect to the aortic root, and the superior vena cava orifice.
7. Size the occluder such that the radius of the discs will not exceed the lesser of the two measurements from the centrally located defect to the aortic root and the superior vena cava orifice. See Table 1 in Appendix A: Supplemental Information.

 NOTE: In Table 1, Column A contains the measurements for the shortest distance from defect to aortic root or the distance from defect to superior vena cava orifice. Column B contains the suggested device size.
8. Prepare the delivery system for use according to the delivery system's instructions for use.
9. Insert the dilator into the delivery sheath and tighten the rotating luer.
10. Advance the dilator and delivery sheath over the guidewire, through the defect, and into the left upper pulmonary vein. Verify the correct position of the sheath by a test hand injection of contrast medium or by echocardiography.

 CAUTION: If strong resistance is met during manipulation, discontinue the procedure and determine the cause of resistance before proceeding. If the cause of the resistance cannot be determined, withdraw all components.
11. Attach the hemostasis valve to the proximal end of the loader and flush with sterile saline. Slowly remove the dilator to prevent ingress of air. Allow blood backflow to purge all air from the system.
12. Slowly remove the guidewire and flush the delivery sheath with sterile saline.
13. Prepare the occluder for use:
 - Inspect the outer sterile pouch and verify that it is unopened and undamaged. Do not use the occluder if the sterile barrier is open or damaged.
 - Gently open the outer, then the inner, pouches.
14. Advance the delivery cable through the hemostasis valve and loader.
15. Attach the device to the distal end of the delivery cable by rotating the device clockwise approximately five turns until the device is secure. Then rotate the device counterclockwise 1/8 of one turn to facilitate release of the device.
16. Immerse the device and loader in sterile saline and pull back on the delivery cable to retract the device inside the loader. Flush the occluder via the side arm of the hemostasis valve.
17. Attach the loader to the delivery sheath and tighten the rotating luer to lock the components together.

 CAUTION: Attach the loader firmly to the sheath to ensure that there are no gaps between the inner surfaces of the two components.
18. Advance the delivery cable and occluder through the delivery sheath until the occluder reaches the tip of the delivery sheath. Do not rotate the delivery cable.

 CAUTION: Do not advance the delivery cable and occluder if resistance is felt.
19. Under fluoroscopic and echocardiographic guidance, deploy the left atrial disc and part of the connecting waist and pull the occluder gently against the atrial septum, which can be felt and also observed by echocardiography.
20. Maintain a slight tension on the delivery cable while retracting the delivery sheath approximately 5–10 cm to deploy the right atrial disc.
21. Confirm correct placement. Use fluoroscopy and echocardiography to confirm that the occluder is in place and evaluate for residual shunt or valve insufficiency.

 WARNING: Do not release the occluder from the delivery cable if the occluder does not conform to its original configuration or if the occluder position is unstable or interferes with any adjacent cardiac structure (such as SVC, PV, MV, CS, AO). Under these circumstances, advance the delivery sheath to recapture the occluder and redeploy. If the occluder position is still unsatisfactory, recapture the occluder and replace it with a new occluder or abort the procedure.

 NOTE: The occluder can be deployed a maximum of two times.
22. Release the occluder:
 - Attach the plastic vise to the delivery cable by tightening the screw on the vise.
 - Release the occluder by rotating the vise counterclockwise as indicated by the arrow.
 - In the unlikely event that this should not be possible, advance the sheath against the right atrial disc to secure the occluder, which will facilitate detachment.
23. Remove the delivery cable and delivery sheath from the patient.

Post-Procedure Care

- A transthoracic echocardiogram should be performed before discharge to assess for device position and residual shunt.
- Patients with any observed pericardial effusion following the occluder implantation should be closely monitored with serial echocardiogram performed until the pericardial effusion resolves.

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- Clinical follow-up with a cardiologist and echocardiograms is recommended at implant, one day post-implant, pre-discharge, and again at one week, one month, six months, and 12 months post-implant. Clinical follow-up with a cardiologist annually thereafter is also recommended.
- Patients should be educated to seek immediate medical attention that includes an echocardiogram, if they develop signs or symptoms of hemodynamic instability such as chest pain, arrhythmia, fainting, or shortness of breath.
- Patients should be instructed to avoid strenuous activity for a minimum of one month post-device implant or as directed by their physician. Strenuous activities may lead to the increased risk of adverse events including erosion. Patients should be reminded that if they experience any symptoms of shortness of breath or chest pain at any time, and especially after strenuous activity, they should seek medical care immediately.
- Patients should take appropriate endocarditis prophylaxis for six months following occluder implantation. The decision to continue endocarditis prophylaxis beyond six months is at the discretion of the physician.
- Patients should be treated with anti-platelet/anticoagulation therapy (such as aspirin) for six months post-implant. The decision to continue anti-platelet/ anticoagulation therapy beyond six months is at the discretion of the physician.

Storage and Handling Information

There are no special storage requirements beyond what is typical for medical devices.

Reporting Device Incidents

If, in the course of using this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Disposal

Amplatzer™ Multifenestrated Septal Occluder - “Cribriform” components should be appropriately classified as biohazards and disposed of in compliance with applicable facility procedures and local and country laws and regulations.

Warranty

Abbott Medical warrants to buyer that, for a period equal to the validated shelf life of the product, this product shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer’s instructions for use and shall be free from defects in materials and workmanship. Abbott Medical’s obligation under this warranty is limited to replacing or repairing at its option, at its factory, this product if returned within the warranty period to Abbott Medical and after confirmed to be defective by the manufacturer.

EXCEPT AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY AS TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

See the Terms and Conditions of Sale for further information.

For Australia Only

This section is only provided in English.

This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or Australian Sponsor.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical at the address provided with the product and after Abbott Medical confirms that the goods are defective.

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical’s liability for the breach of the term or guarantee is limited to, at Abbott Medical’s option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED. See the Terms and Conditions of Sale for further information.

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™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third party trademark, which is property of its respective owner.
© 2023 Abbott. All Rights Reserved.
Abbott Medical UK Limited