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pressing the compound from the syringe into the dermis.
8. If significant resistance is encountered when pushing the plunger, the injection needle may be moved slightly to allow easier placement of the material or it may be necessary to change the injection needle. Needle jams are more likely with use of needles smaller than 27gauge.
9. Advance the needle into the subdermis to the starting location. Carefully push the plunger of the FACETEM syringe to start the injection and slowly inject the FACETEM material in linear threads while withdrawing the needle. Continue placing additional lines of material until the desired level of correction is achieved.
10. Apply slow continuous even pressure to the syringe plunger to inject the implant as you withdraw the needle. The implant material should be completely surrounded by soft tissue without leaving globular deposits. The injected area may be massaged as needed to achieve even distribution of the implant.

STORAGE

FACETEM should be stored at a controlled room temperature between 1 °C and 25 °C.
The expiration date, when stored in these temperatures, is three years from the date of manufacture. Do not use if the expiration date has been exceeded.

DISPOSAL

Used and partially used syringes and injection needles could be biohazardous and should be handled and disposed of in accordance with facility medical practices and local, state or federal regulations.

FURTHER INFORMATION

If further information is needed or required, please contact the manufacturer or local representatives.

PACKING UNIT

Volume	Contents
0.3mL	1 syringe / 1 box
0.8mL	2 syringes / 1 box
1.0mL	1 syringe / 1 box
1.5mL	1 syringe / 1 box

MANUFACTURED BY

CGBIO
CG Bio Co., Ltd.
B1, 1F, 2F, C-dong, 29, Jeyakdanji-ro, Hyangnam-eup, Hwaseong-si, Gyeonggi-do, Korea
Tel. +82-2-550-8300 Fax. +82-31-352-0135

EXPLANATION OF SYMBOLS	
LOT	BATCH CODE
REF	CATALOGUE NUMBER
	MANUFACTURER
	DATE OF MANUFACTURE
	USE BY
2	DO NOT REUSE
2	DO NOT RESTERILIZE
	DO NOT USE IF PACKAGE IS DAMAGED
	TEMPERATURE LIMITATION
	CONSULT INSTRUCTIONS FOR USE
STERILE	STERILIZED USING STEAM OR DRY HEAT

MADE IN KOREA
CSF-109-1/AUG 2021

Facetem™

Natural Face, Natural System

Manual / English

CGBIO
CG Bio Co., Ltd.
B1, 1F, 2F, C-dong, 29, Jeyakdanji-ro, Hyangnam-eup, Hwaseong-si, Gyeonggi-do, Korea

Manual

Facetem™ *Natural Face, Natural System* English

PRODUCT DESCRIPTION

FACETEM is a sterile, latex-free, non-pyrogenic, semi-solid, cohesive subdermal implant. Facetem contains microparticles of biocompatible, biodegradable, synthetic calcium hydroxylapatite (CaHA) suspended in a gel carrier of sterile water for injection, glycerin and sodium carboxymethylcellulose. Facetem (0.3mL, 0.8mL, 1.0mL and 1.5mL) has a CaHA particle size range of 25-45 microns and should be injected with a 25-27 gauge needle.

INTENDED USE / INDICATIONS FOR USE

FACETEM is indicated for subdermal implantation for the plastic and reconstructive surgery of deep dermal and sub-dermal soft tissue augmentation of the facial area such as Nasolabial folds, cheeks, pre-jowl sulcus and acne scar.

CONTRAINDICATIONS

FACETEM is contraindicated for patients with
- Severe allergies manifested by a history of anaphylaxis or history or presence of multiple severe allergies; and
- Bleeding disorders; and
- Known hypersensitivity to any of the components of FACETEM.

WARNINGS

FACETEM is contraindicated for patients with
- Use of FACETEM in any person with active skin inflammation or infection in or near the treatment area should be deferred until the inflammatory or infectious process has been controlled.
- Injection procedure reactions to FACETEM have been observed consisting mainly of short-term (i.e.,< 7 days) bruising, redness and swelling.
- Special care should be taken to avoid injection into the blood vessels. An introduction into the vasculature may occlude the vessels and could cause infarction or embolism.
- Do not overcorrect (overfill) a contour deficiency because the depression should gradually improve within several weeks as the treatment effect of FACETEM occurs.
- The safety and effectiveness of FACETEM for use in the lips has not been established.

PRECAUTIONS

The calcium hydroxylapatite (CaHA) particles of FACETEM are radiopaque and are clearly visible on CT Scans and may be visible in standard, plain radiography. Patients should be informed of the radiopaque nature of FACETEM, so that they can inform their primary care health professionals as well as radiologists.
FACETEM should only be used by health care providers with expertise in the correction of volume deficiencies after fully familiarizing themselves with the product, the product educational materials and the entire package insert.
FACETEM is packaged for single patient use. Do not re-sterilize. Do not use if package is opened or damaged. Do not use if the syringe end cap or syringe plunger is not in place.
Long-term safety and effectiveness of FACETEM beyond one year have not been investigated in clinical trials.
The safety of FACETEM in patients with increased susceptibility to keloid formation and hypertrophic scarring has not been studied.
As with all transcutaneous procedures, FACETEM injection carries a risk of infection. Standard precautions associated with injectable materials should be followed.
Safety of FACETEM for use during pregnancy, in breastfeeding females or in patients under 20 years has not been established.
Patients who are using medications that can prolong bleeding, such as aspirin or warfarin, may, as with any injection, experience increased bruising or bleeding at the injection site.
Universal precautions must be observed when there is a potential for contact with patient body fluids.
After use, treatment syringes and needles may be potential biohazards. Handle accordingly and dispose of in accordance with accepted medical practice and applicable local, state and federal requirements.
The patient should be informed that he or she should minimize exposure of the treated area to extensive sun or heat exposure for approximately 24 hours after treatment or until any initial swelling and redness has resolved.
Safety and effectiveness of FACETEM in the periorbital area has not been established.
No studies of interactions of FACETEM with drugs or other substances or implants have been conducted.

ADVERSE EVENTS / POST-MARKETING SURVEILLANCE

- Implant and injection site reactions, mostly non-serious events, have been reported. These include: discoloration, bruising, swelling, mass formation, erythema, pain, scarring and ischemia. Most instances of discoloration including hyperpigmentation, sometimes described as a blue or brown color and ranging from mild to severe, have occurred within the same day as treatment but have also occurred up to 6 months post treatment.
- These events typically resolve within a few days. Implant and/or injection site bruising, swelling, erythema and pain generally occurred on the same day as treatment usually resolving within 1 to 4 weeks.
- Serious adverse events have been extremely rarely reported. The most commonly reported serious adverse events were hypersensitivity, and implant and/or injection site swelling, ischemia and discoloration.

*Adverse Events descripted above are common observations of similar products from multiple manufacturers.
**Adverse reactions should be reported to CG Bio Co., Ltd. at 82-2-550-8300, or local representatives.

INDIVIDUALIZATION OF TREATMENT

Before treatment, the patient's suitability for the treatment and the patient's need for pain relief should be assessed. The outcome of treatment with FACETEM will vary between patients. In some instances, additional treatments may be necessary depending on the size of the defect and the needs of the patient.

DIRECTIONS FOR USE

General

- The following are required for the percutaneous injection procedure:
 FACETEM syringe(s)
- 25-27 gauge needle with Luer lock fittings

1. Prepare patient for percutaneous injection using standard methods. The treatment injection site should be marked and prepared with a suitable antiseptic. Local or topical anesthesia at the injection site should be used at the discretion of the physician.
2. Prepare the syringes of FACETEM and the injection needle(s) before the percutaneous injection. A new injection needle may be used for each syringe, or the same injection needle may be connected to each new syringe.
3. Remove PET tray from the carton. Open the PET tray and remove the syringe.
4. Remove the Luer syringe cap from the distal end of the syringe prior to attaching the needle. The syringe of FACETEM can then be twisted onto the Luer lock fitting of the needle. The needle must be tightened securely to the syringe and primed with FACETEM. If excess FACETEM is on the surface of the Luer lock fittings, it will need to be wiped clean with sterile gauze. Slowly push the syringe plunger until FACETEM extrudes from the end of the needle. If leakage is noted at the Luer fitting, it may be necessary to tighten the needle, or to remove the needle and clean the surfaces of the Luer fitting or, in extreme cases, replace both the syringe and the needle.
5. Locate the initial site for the implant. Scar tissue may be difficult or impossible to treat. Avoid, if possible, passing through these tissue types when advancing the injection needle.
6. The amount injected will vary depending on the site and extent of the restoration or augmentation desired. FACETEM should be injected subdermally.
7. Insert needle with bevel down at approximately a 30± angle to the skin. Needle should be injected subdermally, placing the needle with an angle of 30-45 degrees. The needle is introduced longitudinally into the deep dermis to create a tunnel. CaHA is then injected into the created space by simultaneously slowly withdrawing the syringe and

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