CRE8™

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IC0713.137.102 Issued 202204 Rev.13

CiD
Carbostent & Implantable Devices

Manufactured by: CiD S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of Alvimedica group
www.alvimedica.com

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ENGLISH

CRE8™
Sirolimus eluting ℰCarbofilm™ coated coronary stent on rapid exchange balloon catheter

1. DESCRIPTION

The CRE8™ device consists of a Sirolimus eluting coronary stent firmly held on the distal end of a semicompliant balloon catheter.
The coronary stent is a flexible implantable device that can be expanded using a PTCA catheter.
The stent is made of cobalt chromium alloy (L605) and is coated with ℰCarbofilm™, a thin carbon film with a high-density turbostratic structure substantially identical with that of the pyrolytic carbon used for mechanical cardiac valve discs. Coating the substrate with ℰCarbofilm™ provides it with the bio- and haemocompatible characteristics proper of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.
The outer surface of the stent has dedicated grooves, fully coated with ℰCarbofilm™, for containing the pharmaceutical formulation Amphilimus™, which is composed of the drug Sirolimus and a mixture of long chain fatty acids. The specific drug dosage for the vessel surface treated is approximately 0.9 µg/mm², which corresponds to a minimum dose of 50 µg on the smaller stent (2.25x8mm) and a maximum dose of 395 µg on the larger stents (4.0x38 mm and 3.5x46 mm).
Two radio-opaque platinum markers at either end of the stent enable accurate positioning on the lesion to be treated.
The rapid-exchange type **balloon catheter** provides a safe means for delivering the coronary stent to the lesion to be treated.
The distal portion of the catheter has a hydrophilic coating and consists of two lumens: one for balloon inflation and deflation, the other for guidewire advancement and retraction.
Two radio-opaque markers, located outside the usable length of the balloon, enable accurate placement across the stenosis.
The proximal portion of the catheter, a stainless steel hypotube, contains the lumen for balloon inflation and deflation.
Two depth indicators, placed 90 and 100 cm from the distal end, indicate the end of the balloon catheter from the guiding catheter, respectively in cases of brachial or femoral/radial approach.
The proximal end of the catheter has a female Luer lock for attachment to an inflation device.
The Manufacturer produces the CRE8™ device and carries out all quality controls both during manufacture and on the finished product, in accordance to the appropriate good manufacturing standards.

How supplied:
CRE8™ device is supplied sterile, non-pyrogenic and individually packed in a foil pouch that must **not be placed into a sterile field**.
NOTE: The foil pouch is the only sterile barrier.
The manufacturer sterilization uses an ethylene oxide and CO_2 mixture.
Sterility and chemical stability are guaranteed while the package remains intact and correctly stored, until the expiry date printed on the packaging (EXPIRY DATE).

Contents:
- One Cre8™ device
- One compliance table

Storage:
Store in a cool, dry place, away from sunlight, at a temperature of 25 °C; excursion permitted to 30 °C.

2. INTENDED USE

The coronary stent is indicated for the treatment of stenotic lesions in coronary arteries and coronary artery bypasses to maintain vessel patency.
Randomized clinical trials have shown that drug-eluting stents can significantly reduce angiographic late loss (late lumen loss), binary restenosis and repetition of target lesion revascularization.

3. INDICATIONS

The stent is indicated for improving coronary luminal diameter in patients with symptomatic ischemic heart disease due to de novo and restenotic lesions of native coronary arteries.

4. CONTRAINDICATIONS

Use of the stent is contraindicated in the following situations:
- Lesions considered untreatable with PTCA or other interventional techniques
- Disorders/allergies that limit the use of anti-platelet and/or anticoagulant therapy
- Severe allergy to the contrast medium or drugs used during the procedure
- Lesions on a vessel with reference diameter < 2.25 mm
- Patients with known hypersensitivity or allergies to Sirolimus, fatty acids (such as stearic acid, palmitic acid, behenic acid).

5. WARNINGS

- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient with great care before making procedural choices.
- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician should refer to updated literature for information regarding the risks and advantages of the various procedures before choosing the procedure to adopt.
- **In patients with known hypersensitivity or allergies to metal components of the stent, the Cre8™ stent should be used only if potential benefits of the stent outweigh potential risks.**
- There are no adequate or well-controlled studies in pregnant women for this product. The Cre8™ stent should be used during pregnancy only if potential benefits of the stent outweigh potential risks.
- The device is intended for single use only. Do not reuse, reprocess or re-sterilize. This could compromise device performance and cause contamination risks of the device and infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The CRE8™ device should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce alterations.
- The stent should not be touched. This is most important during catheter removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
- The CRE8™ device is designed to perform as a system. Do not use its components separately.
- The CRE8™ device is indicated for use in combination with PTCA. The delivery catheter is not a coronary dilatation catheter; it should be used for stent deployment only.
- Do not use the CRE8™ device if its packaging has been opened or damaged, incorrectly stored, or if the stated “Expiry date” has passed. In such cases product performances and sterility are not guaranteed.
- Do not use the CRE8™ if the shaft of the catheter has kinks or bends; in such cases do not try to straighten the catheter.
- The stent should not come into contact with liquids prior to preparation and positioning. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline solution, the contact time should be limited (maximum one minute).
- A heart surgery team should be available for a possible intervention.
- Follow the manufacturer’s instructions when using accessories (guiding catheter, guidewire, haemostatic valve).
- The CRE8™ device should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
- The traditional procedure for stent deployment requires pre-dilatation of the lesion.
- Recent scientific literature describes procedures performed without pre-dilatation.
- The stent should be implanted in the target lesion using its delivery catheter.
- If resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the CRE8™ device and the guiding catheter as a single unit. Applying excessive force and/or incorrectly handling the system may lead to stent deployment or damage to the delivery catheter.
- If the patient has multiple lesions in a single vessel, it is advisable to treat first the distal lesion, then the proximal lesions. This order will reduce the need to cross the proximal stent during implantation of the distal stent and therefore reduces the risk of displacing the proximal stent.
- Implantation of a coronary stent may cause dissection of the vessel distal and/or proximal to implantation site and may also cause acute occlusion of the vessel, making an additional intervention necessary (CABG, further dilatation, placement of additional stents, or other procedures).
- Do not pressurise the system until the stent is positioned across the lesion to be treated.
- Never use air or any other gas to inflate the balloon.
- During catheter inflation do not exceed the Rated Burst Pressure.
- Do not retract the catheter at the end of the procedure until its balloon is fully deflated.
- Do not attempt to reposition a partially expanded stent. Attempted repositioning may cause severe vessel damage.
- If the stent is lost in the coronary vessel, recovery procedures can be initiated. The procedures can, however, cause injury to coronary vessels and/or to the vessel access site.
- Do not attempt to clean or re-sterilize devices that have been in contact with blood and organic tissues. Used devices are to be disposed of as dangerous medical waste with risk of infection.

6. PRECAUTIONS

- Before the angioplasty procedure, examine the catheter to check that it operates correctly (no kinks, bends or other damage) and ensure that the dimensions are correct for the specific procedure in which it is to be used.
- Use of a gauged inflation device is strongly recommended.
- The Cre8™ device should be used only by physicians trained specifically to perform percutaneous transluminal coronary angioplasty (PTCA) and coronary stent implantation.
- In cases of acute myocardial infarction, implantation of the stent should be carefully assessed by the operator due to the risk of acute thrombosis.
- Where vessels are excessively tortuous and there is also proximal atherosclerosis, the catheter may be difficult to advance. In such cases incorrect handling could lead to dissection or rupture of the coronary vessel.
- Placement of a stent has the potential to compromise side branch patency.
- An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter.
- To avoid damaging the stent, use extreme care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.

7. MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that CoCr stent product lines are MR Conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:
- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T).
- Maximum spatial field gradient of 2,890 G/cm (28.90 T/m)

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MRI Heating
Under the scan conditions defined below, CoCr stent product lines in a single stent configuration are expected to produce a maximum temperature rise of less than 2.4 °C in a 1.5 T MRI system and 2.9 °C in a 3 T MRI system.

	1.5 T	3 T
MR system reported, whole body averaged SAR	*2.9 W/kg*	*2.9 W/kg*
Calorimetry measured values, whole body averaged SAR	*2.1 W/kg*	*2.7 W/kg*
Highest temperature change	*2.4 °C*	*2.9 °C*

MR Artifact
In non-clinical testing, the image artifact caused by the device extends approximately 1.0 cm from CoCr stents when imaged with a gradient-echo pulse sequence in a 3 T MRI system.

8. DRUG REGIME
Studies in clinical literature indicate the need to administer appropriate anticoagulant, anti-platelet and vasodilation therapy for successful stent implantation.
Antiplatelet treatment during the period after the procedure should be used per the guidelines from current ESC or ACC/AHA/SCAI for PCI.
While physicians should adhere to DAPT duration defined into the current ESC (2017 and subsequent updating) or ACC/AHA/SCAI Guidelines for PCI, the decision to interrupt or discontinue DAPT is the responsibility of the treating physician, taking into consideration the individual patient's condition (e.g. non-deferrable surgery). Guidelines details on indications and modalities for DAPT interruption are available in the ESC or ACC/AHA/SCAI organizations' websites (https://www.escardio.org/ https://www.acc.org/ https://www.heart.org/ http://www.scai.org).
In case unanticipated interruption or discontinuation of DAPT is required any time after one month following Cre8™ products family implantation, one-year clinical data of Cre8™ products family indicate low stent thrombosis rates and no observed increased risk for stent thrombosis.

9. INTERACTIONS WITH DRUGS
Although no specific clinical data are available, certain drugs, such as Tacrolimus, which act through the same binding protein (FKBP), may interfere with the efficacy of Sirolimus. Drug interaction studies have not been performed. Sirolimus is metabolized by CYP3A4. Strong CYP3A4 inhibitors may cause greater exposure to Sirolimus up to levels associated with systemic effects, especially in the case of multiple implant. Systemic exposure to Sirolimus if the patient is taking concomitant systemic immunosuppressive drugs must also be considered.

10.POTENTIAL ADVERSE EVENTS
Possible adverse effects of PTCA include, but are not limited to:
- Access site haematoma or pain
- Acute myocardial infarction
- Allergic reaction (to contrast medium, drugs used during the procedure or the materials composing the stent)
- Arterial spasm
- Arteriovenous Fistula
- Artery aneurysm or pseudoaneurysm
- Cardiac arrhythmia
- Death
- Dissection, perforation, rupture of the artery
- Distal embolism
- Embolization (air, tissue, device or thrombus)
- Haemorrhage
- Hypo/hypertension
- Infections
- Restenosis of the vessel
- Stent occlusion
- Thrombosis (acute, subacute or late)
- Unstable angina
- Ventricular fibrillation
- Vessel occlusion.

11.DIRECTIONS FOR USE
Auxiliary material (not included in the box)
- Guiding catheter(s) with a 5F diameter (inner diameter 1.47 mm) or larger
- Syringes
- Guidewires 0.014 inches in diameter (0.356 mm) or less
- Inflation device
- Introducer

Model
Each Cre8™ device is identified by a model code and a lot number; available product codes are listed in Table 1. The code is composed of the letters IC, followed by two letters identifying the type of device, two or three digits indicating the nominal diameter of expansion of the stent and a further two digits indicating the length of the stent.
The batch number allows traceability of all information regarding device manufacture and system control in the Quality Assurance archives of the Manufacturer.
To facilitate traceability of the device at the user end, the product code is printed on the adhesive labels included on each box; these labels can be attached to the implant patient's medical chart.

Stent + delivery catheter preparation
During preparation of the system, take special care not to let to any liquid come into contact with the stent. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline solution, the contact time should be limited (maximum one minute).
The stent should be implanted in the target lesion using its delivery catheter.
After having examined the package for any damage, remove the CRE8™ device and bring it into a sterile area.
a) Check that the delivery catheter does not present kinks, bends or other damage.
b) Carefully slide the protective cover off the stent by grasping the cover at the distal end do not touch the part of the cover over or proximal to the stent. Check that the stent is intact and well centered on the balloon.
Damage to the CRE8™ device may impair its performance. If the stent is out of place or damaged, do not use it.
c) Flush guidewire lumen with heparin/saline mixture. **Caution: do not place the stent into contact with liquids while flushing the guidewire lumen.**
d) Prepare the inflation device according to the manufacturer's instructions with dilute contrast medium.
e) Purge any air from the balloon the stent is mounted on as follows:
1) Fill the inflation device/syringe with 4 ml of contrast medium diluted (50% contrast medium and 50% sterile saline solution);
2) After connecting the inflation device/syringe to the delivery catheter Luer connector, point the distal tip of the catheter (balloon) downwards vertically;
3) Apply negative pressure and aspirate for at least 30 seconds. **Allow the pressure to rise gradually back to neutral** as the system fills with contrast medium.
4) Disconnect the inflation device/ syringe from the luer connection and remove all air.
5) Without introducing air, repeat steps 2 and 3 and aspirate for 10-15 seconds until air bubbles no longer appear.

Stent insertion
The guiding catheters with a 5F diameter (inner diameter 1.47 mm) or larger are compatible for use with the CRE8™ device.
Coronary guidewires 0.014 inches in diameter (0.356 mm) or less are compatible for use with the CRE8™ device. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.
a) Keep the delivery catheter the stent is crimped on **at ambient pressure during insertion** (neutral pressure).
b) Flush the exposed portion of the guidewire with a heparin/saline mixture to remove traces of blood and contrast medium.
c) Ensure that the haemostasis valve is completely open before introducing the system.
d) Thread the distal tip of the catheter over the proximal end of the guidewire until the guidewire exits the catheter body at the guidewire exit port.
e) Slowly advance the dilatation catheter until the appropriate depth marker aligns with the haemostatic valve hub. This indicates that the balloon is exiting the guiding catheter.
f) Under fluoroscopy, slowly advance the system over the guidewire to allow the blood the retrograde filling of the guiding catheter; keep advancing the system until the stent reaches the treatment site.
g) When advancing the delivery system into the target vessel, make sure that the stent and the delivery balloon move as a single unit. This should be done by observing, under fluoroscopy, the position of the balloon radiopaque markers in relation to the stent radiopaque markers.
If resistance is encountered at any time, do not force the system: withdraw the CRE8™ device and the guiding catheter as a single unit. Applying excessive force and/or incorrectly handling the system may lead to stent deployment or damage to the delivery catheter.

Stent implantation and expansion
a) Position the stent and, with the help of the **balloon's radiopaque markers**, make sure it is correctly positioned in relation to the lesion to be treated and completely covering it.
b) Expand the stent by slowly inflating the balloon to the nominal diameter. Table 2 lists the diameters of the balloon as inflation pressure varies (column I), for the seven nominal diameter groups: 2.25 mm (column II), 2.5 mm (column III), 2.75 mm (column IV), 3.0 mm (column V), 3.5 mm (column VI), 4.0 mm (column VII) and 4.5 mm (column VIII). Values on dark background refer to pressures above the Rated Burst Pressure*. The stent expands at a minimum balloon inflation pressure of about 5 atm.
c) **Stent elastic recoil is between 2 and 7%* depending on stent model and expansion diameter**.
* NOTE: These values are obtained from the results of in-vitro testing
d) **Deflate the balloon and check angiographically that the stent has fully expanded. If necessary, re-inflate the balloon for optimal implantation**.
e) Before withdrawing the delivery catheter, maintain negative pressure for a minimum of 30 seconds, as per routine PTCA procedure. Ensure the balloon is fully deflated before and remove the catheter. A device with larger and longer balloon requires a longer deflation time.
f) If post-dilatation of the stent is required, a PTCA balloon can be used. Use extreme care:
not to expand the 2.25 mm nominal diameter stent over 2.55 mm.
not to expand the 2.50 – 2.75 mm nominal diameter stents over 3.05 mm.
not to expand the 3.00 – 3.50 mm nominal diameter stents over 3.85 mm.
not to expand the 4.00 – 4.50 mm nominal diameter stent over 5.05 mm.
Optimal expansion requires the stent to be in full contact with the arterial wall, so that the stent diameter equals the diameter of the reference vessel. ENSURE THAT THE STENT IS FULLY EXPANDED.

12. LIABILITY AND WARRANTY
The Manufacturer guarantees that this device has been designed, manufactured and packaged with the greatest care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforestated conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not

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eliminate completely, the risks related to the use of the product.
The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its intended use, as mentioned in other sections of this instructions booklet.
Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.
In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects.
Under such circumstances, the customer shall deliver the device to The Manufacturer, which reserves the right to examine the allegedly faulty device and establish whether it actually presents manufacturing defects. The warranty consists exclusively in the replacement of the defective device with another device of the same type or equivalent, from the same Manufacturer.
The warranty applies under the condition that the device is returned correctly packed to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.
On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for return of the defective device.
The Manufacturer declines any responsibility for cases of negligence to observe the methods of use and the precautions stated in this instructions booklet and for cases of use of the device after the use before date printed on the package. Furthermore, the Manufacturer declines any responsibility related to the consequences arising from the treatment choices and the methods of use or the application of the device; the Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device, or the choice of implantation technique used by the operator.
The agents and the representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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CRE8™

SIROLIMUS ELUTING
CARBOFILM™ COATED CORONARY STENT
ON RX BALLOON CATHETER


TABLE 1		Stent length							
		8mm	12mm	16mm	20mm	25mm	31mm	38mm	46mm
Nominal stent diameter	2.25 mm	**ICLI22508**	**ICLI22512**	**ICLI22516**	**ICLI22520**	**ICLI22525**	**ICLI22531**	-	-
	2.5 mm	**ICLI2508**	**ICLI2512**	**ICLI2516**	**ICLI2520**	**ICLI2525**	**ICLI2531**	**ICLI2538**	**ICLI2546**
	2.75 mm	**ICLI27508**	**ICLI27512**	**ICLI27516**	**ICLI27520**	**ICLI27525**	**ICLI27531**	**ICLI27538**	**ICLI27546**
	3.0 mm	**ICLI3008**	**ICLI3012**	**ICLI3016**	**ICLI3020**	**ICLI3025**	**ICLI3031**	**ICLI3038**	**ICLI3046**
	3.5mm	**ICLI3508**	**ICLI3512**	**ICLI3516**	**ICLI3520**	**ICLI3525**	**ICLI3531**	**ICLI3538**	**ICLI3546**
	4.0mm	**ICLI4008**	**ICLI4012**	**ICLI4016**	**ICLI4020**	**ICLI4025**	**ICLI4031**	**ICLI4038**	-
	4.5mm	-	**ICLI4512**	**ICLI4516**	**ICLI4520**	**ICLI4525**	**ICLI4531**	-	-

TABLE 2		Balloon diameter (mm)						
		2.25	**2.50**	**2.75**	**3.00**	**3.50**	**4.00**	**4.50**
Inflation pressure (atm)	6	2.08	2.31	2.54	2.77	3.26	3.74	4.15
	7	2.13	2.37	2.61	2.85	3.35	3.84	4.27
	8	2.19	2.43	2.68	2.93	3.43	3.92	4.38
	NP 9	**2.25**	**2.50**	**2.75**	**3.00**	**3.50**	**4.00**	**4.50**
	10	2.31	2.55	2.83	3.08	3.58	4.09	4.61
	11	2.36	2.61	2.90	3.15	3.66	4.17	4.68
	12	2.40	2.65	2.95	3.20	3.71	4.25	4.74
	13	2.44	2.70	3.01	3.25	3.77	4.32	4.80
	14	2.48	2.74	3.06	3.29	3.82	4.37	4.86
	15	2.51	2.77	3.09	3.32	3.86	4.41	4.92
	16	2.55	2.80	3.12	3.36	3.90	4.46	4.98
	17	2.59	2.83	3.15	3.39	3.94	4.50	5.04
	RPB 18	2.62	2.86	3.19	3.43	3.98	4.54	5.11
	19	2.65	2.89	3.23	3.47	4.03	4.59	5.19
	20	2.68	2.93	3.27	3.51	4.08	4.65	5.27
	21	2.71	2.96	3.30	3.55	4.12	4.69	5.33
	22	2.74	2.99	3.34	3.60	4.18	4.74	5.40

These values result from in vitro testing

NP Nominal Pressure
RBP Rated Burst Pressure Do not exceed
Pressure above Rated Burst Pressure

If post dilatation is required, a PTCA balloon can be used. Use extreme care:
not to expand the 2.25 mm nominal diameter stents over 2.55 mm
not to expand the 2.50 mm and 2.75 mm nominal diameter stents over 3.05 mm
not to expand the 3.00 mm and 3.50 mm nominal diameter stents over 3.85 mm
not to expand the 4.00 mm and 4.50 mm nominal diameter stents over 5.05 mm

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