*safe*CLINITUBES capillary tubes

989-890. 201409L.

en	INSTRUCTIONS
bg	ИНСТРУКЦИИ
cs	INSTRUKCE
da	BRUGSANVISNING
de	GEBRAUCHSANLEITUNG
el	ΟΔΗΓΊΕΣ
es	INSTRUCCIONES
et	KASUTUSJUHISED
fi	OHJEET
fr	MODE D'EMPLOI
hr	UPUTE
hu	UTASÍTÁSOK
it	ISTRUZIONI
ja	インストラクション
lt	INSTRUKCIJOS
lv	INSTRUKCIJAS
nl	GEBRUIKSAANWIJZING
no	INSTRUKSJON
pl	INSTRUKCJE
pt	INSTRUÇÕES
ro	INSTRUCŢIUNI
ru	ИНСТРУКЦИИ
sk	NÁVOD
sr	INSTRUKCIJA
sv	BRUKSANVISNING
tr	KULLANIM TALİMATLARI
zh	指示

en: Sensitive to light. Store in a dark place.
bg: Светлочувствителен. Съхранявайте на тъмно.
cs: Citlivý na světlo. Skladujte na tmavém místě.
da: Følsom over for lys. Opbevares mørkt.
de: Lichtempfindlich, im Dunkeln aufbewahren.
el: Ευαίσθητο στο φως. Αποθηκεύστε σε σκοτεινό μέρος.
es: Sensible a la luz. Almacenar en lugar oscuro.
et: Tundlik valgusele. Säilitada pimedas.
fi: Herkkä valolle. Säilytettävä pimeässä.
fr: Sensible à la lumière. Conserver à l'obscurité.
hr: Osjetljivo na svjetlost. Čuvati na tamnom
hu: Fényre érzékeny. Tárolja sötétben.
it: Fotosensibile. Conservare al buio.
ja: 感光性です。暗所で保存してください。
lt: Jautrus šviesai. Laikyti tamsioje vietoje.
lv: Jūtīgs pret gaismu. Uzglabāt tumšā vietā.
nl: Lichtgevoelig, niet blootstellen aan direct zonlicht. Opslag in een donkere ruimte.
no: Lysfølsom, lagres mørkt.
pl: Wrażliwe na światło.
pt: Sensível à luz. Armazenar num local escuro.
ro: Sensibil la lumină. A se păstra în spații

ru: Светочувствителен. Хранить в темном месте.
sk: Citlivé na svetlo. Skladujte na tmavom mieste.
sr: Osetljivo na svetlost. Čuvati na tamnom mestu.
sv: Ljuskänslig. Förvara mörkt.
tr: Işığa karşı hassastır. Karanlıkta saklayınız.
zh: 对光敏感。请避光保存。

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23-10-2014 10:31:58

en

Intended use: Heparinized capillary tubes intended for blood sampling.

WARNING – *Risk of infection*
Only lot authorized control and work with blood samples. Make sure gloves are worn.

Sampling procedure using mixing wires and end caps

1. Loosely put an end cap on one end of the capillary tube.
2. Put a mixing wire into the capillary tube and let it slide to the same end as the end cap.
3. Select the sampling site and increase the circulation in the region.
4. Apply antiseptic to the sampling site and let it dry completely prior to collection.
5. Puncture the skin so that drops of blood form quickly.
6. Fill the capillary tube without touch the outside of the blood drop so that air does not enter the capillary tube.
7. Tighten the end cap. Put an cap on the other end of the capillary tube and tighten.
8. Gently invert the capillary tube to mix the blood and anticoagulant.
Note: For the best results, mix the blood with the heparin from one end cap to the other a minimum of ten times. Do not shake the capillary tube. Do not remove it during the capillary tube.
9. For best results, analyze the sample immediately. If storage is unavoidable, store the capillary tube horizontally at room temperature and analyze within 10 min.
Note: Sample deterioration with increased storage time (greater variability of gasses and of thb-measurements).
10. If stored, mix the sample again as described in step 8.
11. Slide the mixing wire to the end of the capillary tube opposite to that from which the blood is to be aspirated. Do **not** remove the mixing wire.
12. Gently remove the end caps.
13. Analyze the sample.
Note: For more information, refer to the instructions for use or operator's manual for your analyzer.

Sampling procedure without using mixing wires and end caps
Please refer to points 3, 4, 5, 6 and 13 in the process above.

WARNING – *Risk of infection*
Discard the waste as infectious waste [1].
[1] CLSI document GPS5-A2, Clinical Laboratory Waste Management

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Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS.
© 2014 Radiometer Medical ApS. All Rights Reserved.