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INSTRUCTIONS FOR USE

July 1, 2024

Product name: Sterile single use dental needles

Brand name: HOGEN SPITZE

Instructions for use

• **Device name:** Sterile single use dental needles


STERILE EO CE 1639

• **Intended purpose:** Sterile dental injection needles for single use which are used in dental cartridge syringe for injection of dental local anesthetics.
• **Patient target groups:** Pediatric to elderly patients of all genders requiring injection of a local anesthetic for dental treatment.
• **Intended users:** Dentists

• **Indications:** For dental treatment
• **Contra-indications:** Metal hypersensitivity, History of untoward reaction to injection in the past
• **Undesirable side-effects:** Stainless steel contains trace amounts of nickel and chromium, therefore may cause a reaction to those with hypersensitive skin.

• **Performance characteristic:** Injection of dental local anesthetics.

• **Method of use**

[1] Break the hermetic seal (Twist cap to open – Do not bend or snap open).
[2] Connect the needle on the dental cartridge syringe.
[3] Remove the protective cap to use.

• **Compatible syringe:** Use the appropriate dental cartridge syringe below according to the thread type of product.
- Metric threaded dental cartridge syringe
- Inch threaded dental cartridge syringe

• **Warnings**
- Do not use before checking the instructions for use
- Do not use if seal is broken
- Do not reuse, risk of patient cross-infection
- Do not use if use by date has expired
- Do not bend the needles. These can break the needle and cause injury
- Dispose of used needles immediately in a sharps disposal container according to the your local regulations

• **Residual risks:**
- Use of unintended user may cause patient injury.

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- Use of needle to patients with a metal hypersensitivity and history of untoward reaction to injection in the past may cause patient injury.
- Reuse of needle may cause re-or cross infection.
- Use of expired product may cause infection.
- Use of expired product may not be able to achieve the intended performance, which can lead to patient injury.
- Use of improper syringe may cause leak local anesthetics.
- Improperly dispose of used needle may cause needle stick injury, which can lead to re-or cross infection.
- Use of needle with a broken seal may cause infection.
- Use without checking the instructions for use may cause injury to the patient due to omission of safety-related information.
- Bending needles may cause needle breakage, which can lead to patient injury.

• **Storage condition:** Room temperature

• **Note:** If any serious incident that has occurred in relation to the medical device, please report it to the manufacturer and/or authorized representative and competent authority.

C-K DENTAL IND.CO.,LTD. 88-50, Sudo-ro, Bucheon-si, Gyeonggi-do, Republic of Korea Tel: +82-32-679-1204, Fax: +82-32-679-1207 Email: ck-dental@hanmail.net	**EC REP OBELIS S.A.** Bd. Général Wahis, 53, 1030 Brussels, Belgium Tel: +32-2-732-59-54, Fax: +32-2-732-60-03 Email: mail@obelis.net

Date of issue and version:	**2023-06 / R02-1**

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C-K DENTAL IND. CO. LTD.