Instruction for Use

Product Name: Intraoral Scanner

Model: IOS-11

Version No	Content/modification	File No.	Date
1.0	First Edition	TF-IOS-05	2023-05-08

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Product photos in kind prevail, Runyes reserve all the right for the final explanation.

Version Number: 2023-05-08

Runyes® 3DS

Intraoral Scanner

Operation Manual

Technical Manual

Ningbo Runyes Medical Instrument Co., Ltd

Ningbo Runyes Medical Instrument Co., Ltd

After sale service: Ningbo Runyes Medical Instrument Co., Ltd

Address: 032 Building, No.456, Tonghui Road, Jiangbei Investment & Pioneering Park C, 315033, Ningbo, China

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User Guide

System software and hardware upgrades are initiated through Runyes only. No software or hardware should be added to or deleted from the Runyes systems without prior approval of Runyes. Doing so may result in damage to the system and will void the product warranty.

The user must read and comply with all safety, warning, and instructional labels on the Runyes products.
Ensure your Runyes products are properly maintained through periodic maintenance. If you suspect equipment malfunction or failure, discontinue using the products and contact local distributor or call Runyes after-sales hotline immediately. Do not attempt to make any repairs on Runyes products.

Product Name:	Intraoral Scanner
Product Model:	IOS-11
Serial Number:	
Date of Manufacture:	

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►Chapter 4 Software Installation 12

4.1. Install software 12

4.2. The image acquisition process of Intraoral Scanner 14

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►Chapter 5 Maintenance....24

5.1. Visual Inspection24

5.2. Cleaning, sterilization, and disinfection....24

►Chapter 6 Common problems and solutions....27

►Chapter 7 Product application precautions.28

►Chapter 8 Circuit connection chart29

►Chapter 9 Technical Specification....30

9.1. Scanner technical Specification....30

9.2. Equipment operating environment：30

9.3. Equipment storage, transportation environment:30

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►Chapter 1 Introduction

Thank you for your trust in us. We hope this product will satisfy you completely.
It is recommended that you read this manual before installation and use, which is beneficial to eliminate and reduce the risk of wind hazard to you and your patients due to disoperation.

1.1. Components

Scanner host, scanner head, base, power adapter, USB control box, data cable, Dongle, and U disk with software.

1.2. Intended use

Using optical scanning method to collect the 3D geometric data of the tooth in the patient's mouth and provide the 3D digital model for CAD/CAM denture design and processing.

1.3. Intend User

The scanner should be operated by a healthcare professional educated and competent to perform dental intraoral image acquisition.

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No

1.7. Classification

IEC 60601-1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2 Medical Electrical Equipment PART 1-2: General Requirements for Basic Safety and Essential Performance Collateral Standard: Electromagnetic Disturbances Requirements and Tests
EN ISO 780 Packaging - Distribution packaging - graphical symbols for handling and storage of packages

- Safety classification: type B medical device equipment.
- Applied parts: scanner head.
- Not belong to the category AP or APG device.
- Operation mode: Continuous.
- Degree of protection against harmful ingress of water: IPX0.

1.8. Working Principle

The intraoral scanner uses the principle of structured light for three-dimensional reconstruction. The scanning system includes two subsystems: projection and imaging. The projection system projects the structural light onto the target object. The imaging system captures the corresponding image and analyzes the three-dimensional information from the stripe image. Then, according to the three-dimensional characteristics, a series of three-dimensional surfaces are matched and fused to obtain a complete three-dimensional object model.

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1.9. Specification

According to the area of Intraoral Scanner, it is classified as follows:

Specification	Scope of reconstruction	Pixel
I0S-11	14 x 14 x 15 mm	1024 X 768 pixels


1.10. Components of Scanner

No.	Components	Quantity	Remarks
1.	Scanner host	1	Without Scanner head
2.	Scanner head	3	Size 1: 78.9mm×19.8mm×15.8mm
		1	Size 2: 79.2mm×16.2mm×12.8mm
3.	base	1	
4.	Data cable	1	
5.	USB control box	1	USB 2.0/3.0
6.	Power adapter	1	INPUT:100-240V~ 50/60Hz 500mA OUTPUT: DC 12V , 1.5A
7.	U disk with software	1	
8.	Dongle	1	
9.	Operational Manual	1	

Please check the following table to see if the parts are complete before the installation of the Intraoral Scanner. If the individual content does not match and you cannot install the system, please contact the local distributor or agent for support.

1) Scanner host：

The scanner's 3D information is reconstructed from 2D images taken by the scanner. The scanner is mainly composed of optical and imaging systems, and the head has protective glass protection.



2) Scanner head：

Made of medical plastic and optical prism, it is installed into the head of the scanner for scanning. The scanner head can be sterilized in accordance with the prescribed method.

The scanner head comes in two sizes, as shown in the following pictures:

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
3) Base:

Scanner host base bracket, when the host is placed on the base in the scanning state, it will enter the standby state.

4) USB control box with data cable:

It is used to transmit the output signal from the scanner host to the computer, and connect the scanner host and the computer.

5) Data Cable:

Used to transmit the output signal from the scanner host to the USB control box.

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It is used to transmit output signals from the base to the computer, connecting the base and the computer.

Dongle	
U disk with software	

7) Power Adapter：

Medical certification power adapter.

INPUT: 100-240V~ 50/60Hz 500mA

OUTPUT: 12V⎓, 1.5A

CAUTION The Intraoral Scanner is integrated with the Power Adapter and cannot be replaced at will.

1.11. Icons

Symbol	Description
	Caution: please refer to the instructions in the user manual.
	Express important instructions. Failure or damage system or other property if not observed.
	Represents useful information and how to use our software.
	Represents warnings and safety instructions. Failure to comply may be a serious threat and injury to the patient and operator

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	Refer to instruction manual/booklet for safety instructions.
	Type B applied part
	Stand-by
	Serial Number
	Date of manufacture
	Manufacturer Manufacture's name and address
	Unique Device Identifier
	Fragile Contents of the transport package are fragile therefore it shall be handled with care
	This way up
	Keep away from rain Transport package shall be kept away from rain
	Package symbol, stacking limit by number, it shows the maximum number of identical transport package which may be stacked on the bottom one, where “n” is the limiting number.

1.12. Expected User Profile

Staffs who are engaged in the diagnosis and treatment of dental diseases in dental hospitals and clinics, and have been trained in the operation of this product.

Operators should be familiar with the use of personal computers and related programs to make it easy to use features on a personal computer.

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►Chapter 2 Safety

♦2.1 General warning

Read the instructions to familiarize yourself with the device before putting it into service.

Read the warnings and the safety instructions carefully.

No modification of this equipment is allowed.

The equipment can only be sent back to our company for repair.

♦2.2 General Safety

The user is responsible for the obligations associated with the installation and operation of the equipment.

Only qualified technicians can use and maintain this equipment.

Turn off the machine and unplug the power cord before cleaning and disinfecting.

Use only original accessories and accessories supplied or recommended by the manufacturer.

Do not connect to devices that are not compatible with the system.

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►Chapter 3 Hardware Installation

3.1. What to do before use

Remove the Intraoral Scanner from the original packaging and check that the device received matches your reservation. All accessories are intact and undamaged.

You need to install the driver in the PC which to identify the Intraoral Scanner.

3.2. Installation notes

1. When installing, please handle with care and minimize the distortion and pull of the wires. Do not tread or roll the wires.
2. Place the base on a flat and stable position to prevent falling.
3. Do not drop the handle and base on the ground so as not to cause irreparable damage.
4. In order to avoid disturbing images, don't let the system approach the strong magnetic field and avoid the static emission source.
5. Although the electromagnetic interference of this product is low, it does not guarantee that the operation does not affect the surrounding equipment. If there is interference, please keep the product away from the interference device.
6. Power adapter plug is intended to be used as the isolation means from supply mains, do not position the Intraoral Scanner so that it is difficult to operate the disconnection device.

3.3. PC Specification

We cannot guarantee the work of the Intraoral Scanner and the pirated Microsoft Windows. So please use the legitimate version of Microsoft Windows 10.

*Operating system

Microsoft Windows 10 64 home edition or above

CPU: Intel i7-6700 quad-core/compliance processor or above;

*Recommend: Intel i7-8700 or above (Laptop i7-8750H or above)

Video card and video memory: NVIDIA GeForce 960GT or above graphics card, more than 2G;

*Recommend: NVIDIA GeForce 1060GT or above graphics card, more than card, more 6G

Memory：16G or higher;

*Recommend: 256G solid state drive or 128G solid state drive +1T mechanical hard disk

Display resolution：1920x1080;

USB port： USB2.0/3.0

Additional equipment connected to the Intraoral Scanner must comply with the current IEC or ISO standards (e.g. IEC60950-1 for data processing equipment). Furthermore all configurations shall comply with the requirements for medical electrical systems (see clause 16 of IEC60601-1).

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Anybody connecting additional equipment to the Intraoral Scanner configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements.

When more than one device is connected to use, the accumulation of leakage current may cause the danger of security. If in doubt, consult your local representative or the technical service department.

NOTE Before installing the software and the server, you must close the Windows system's firewall and the anti-virus software with firewall function to make sure the software can be installed and run properly.

NOTE If you need to install additional software on your computer, only install internationally recognized programs.

3.4. Installation steps


Fig. 3.4.1

1. View Figure 3.4.1, connect the data cable 3 to the host 1 of the oral digital impression machine, and tighten the cable lock shell 2.

Pay attention to insert the data cable 3 into the rear port of the host computer 1 first, and then tighten the cable lock shell 2, pay attention to the direction of the buckle in Figure 3.4.2.

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
Fig. 3.4.2

2. The other end of the data cable 3 is connected to the USB control box 4.

3. Connect the USB cable interface 6 on the USB control box 4 to the computer USB port.

4. The adapter is connected to the power supply, and the adapter cable 5 is connected to the USB control box 4. The indicator light of the USB control box 4 is on and the whole machine is powered on.

5. Pick up the main unit 1 of the oral digital impression machine, press the button 7, the button indicator light is on, indicating that it can be scanned.

Pay attention to install the scanning head before use, and insert and unplug it in parallel with the direction of the card slot as shown in Figure 3.4.3, do not tilt it, and the protection head is inserted and unplugged in the same way.


Fig. 3.4.3

Button indicator:

- green light self-test/scanning
- Blu-ray standby
- blue light flashing
- USB connection failed



6. When the scanning is over, put the main unit 1 of the oral digital impression machine into the base 8, and the main unit will enter the standby mode and stop the scanning work.

7. cut off the connection between the adapter and the power supply or disconnect the connection between the adapter line 5 and the USB control box 4, and the whole machine is powered off.

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
3.5. Cable Connection


1. Connect the scanner host to data cable.

2. Connect the data cable to the USB control box

3. Plug USB control box to PC's USB port.

4. Connect the power adapter to USB control box.

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4.1. Install software

Step 1:

Insert U disk to the computer, manually operate installation.

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
Step 3:

Enter the installation wizard, click "Next".

Step 4:


Read carefully the license terms before installation, if you agree with the agreement, please click "I Agree", if not click "Cancel".

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
Step 5:

When "Dental Viewer’s install shield pops up, click “Next"

Step 6:

Installation process, click" Finish" when it is finished.

4.2. The image acquisition process of Intraoral Scanner

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Warning：Surface high temperature

The Intraoral Scanner base has a preheating device for preheating the scanner head to reduce fogging of the scanner head lens during scanning. The temperature of the heating sheet can reach about 60 degrees. Please do not touch it with your hands to prevent burns.

Description: image brightness.

The image is automatically adjusted to the image brightness, always with the best brightness. The ambient lighting should be as dark as possible. Avoid any external light source entering. Turn off the oral lights.

Important: do not use hemostatic cotton in impression meter.

Do not use hemostatic cotton or other objects that may affect the image in the scanning area, otherwise it will result in an error.

- **Ready to Start Scanning:**
 a) Make sure all the connections be ready;
 b) Enter the patient information, enter the scanning interface, and select the dental arch to be scanned;
 c) Take the scanner host from the base, press the button, when the light becomes green means it is ready to start scanning.

- **Scanner head:**

Important: The camera should be cleaned and disinfected after each use. Please follow the instructions of cleaning and disinfection to avoid cross-infection between patients.

There are four consecutive steps of images capturing:

- occlusal surface
- cheek side
- Lingual Side
- Proximal surface

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A. Occlusal surface scanning


Important: The distance between the radiating window of the impression meter and the surface to be measured is required. The distance must be kept between 0 and 15mm (best distance: 5mm). The camera should not be on the teeth or gums. If the distance is too long, data cannot be collected.

B. Buccal side scanning

The scanner head is located near the adjacent tooth of the prepared tooth.

1. Turn the impression meter to the cheek from 45 degrees to 90 degrees.(max degree)
2. Move the imprinting device over the back of the tooth, through the entire buccal side.

C. Lingual side scanning

The scanner head is located near the adjacent tooth of the prepared tooth.

3. Rotate the scanner from 90 degrees on the cheek to the side of the tongue 45 degrees to the maximum 90 degrees.
4. Move the scanner in the near direction over the front teeth, passing through the entire tongue.

D. Mesial side scanning

Move the scanning tip towards the direction of mesial side by tilting the scanning head, inorder to take the good image of adjacent teeth.

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• Please remove the soft tissue.

• Remove the tooth silver from the activity, and make the tooth silver distance between 1 and 2mm.


The next button is completed, and the calculation stage is optimized. After the optimization calculation, the final generated 3d image will be displayed. Please check the correct output. If part of the image is missing, click “step up” and continue scanning.

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4.3. Running Intraoral Scanner


Step 1:

Open the shortcut of "Runyes 3DS" on the desktop to enter the software interface

Step 2


Click the "New Patient" icon in the red box

Step 3


Enter the "patient name" and other relevant information, click Add

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Step 4

Click the "New Case" icon in the red box

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
Step 6

Click "Scan" in the lower right corner

Step 7

The main function of the main interface

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
1. System time and software version information
2. Process steps
3. Patient name
4. Scan object switching
5. 3D scan results show
6. Function key area: including "True Color", "Lock the Scanning View", "Al Scan", "Metal Scan", "Eraser", "Cut", "Lock Scan Area", "Delete", "Preview", " Reset View", "HD Camera"
7. Scanner Status/Video

Step 8

Switch the scan object to "upper jaw" and scan to obtain a 3D model of the patient's upper jaw.

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Step 9:


Switch the scan object to "lower jaw" and scan to obtain a 3D model of the patient'slower jaw.

Step 10:

Switch the scan object to "occlusion", and ask the patient to bite up and down in a normal state. Scan the teeth and move the scanner head up and down to scan part of the upper and lower teeth, and then the software will automatically occlude the upper and lower teeth.

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
Step 11:

Click "Check" to check whether there is any problem in the final 3D model and whether it meets the design requirements, etc.

Step 12:

Click "Export", select the required file type, then click "OK", you can find the file under the corresponding path.

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►Chapter 5 Maintenance

- When used by patients, do not maintain or disinfect the scanner inside the mouth.
- There are no user repairable or replaceable components in the intraoral scanner.
- You must return to the original factory to open and repair the scanner inside the entrance.

5.1. Visual Inspection

5.1.1. Check the surface of the machine (including power adapter).

If there is any stain or smudge, it should be cleaned in time.

If there is damaged/distorted enclosure, do not use.

Visually check the connectivity of all cables. Reconnect the cable that has been loosened or disconnected.

5.1.2. Please clean the plug of power adapter at least once a year. Too much dust on plug may cause the fire.

5.2. Cleaning, sterilization, and disinfection

5.2.1. Scanning head

a) **Cleaning:** Separate the scanning head and scanner body. First, clean the scanning head and lens with clean water (running water) and alcohol (75% concentration medical alcohol). Pay attention to removing dirt, oil stains, spots and other marks on the lens, and dry the scanning head and lens with a soft cloth (non-woven fabric) and dust-free cotton swab, and do a good job of drying.

b) Sterilization and disinfection:

A. Pressure steam sterilization

Place the scanning head into a specialized sterilization bag and seal it. Use a pressure steam sterilizer to sterilize it in the following two ways:

- Select 134 ℃ sterilization time for more than 3.5 minutes;
- Select 121 ℃ sterilization time for more than 15 minutes;

After sterilization, it is necessary to store according to the sensing requirements. Before use, it is necessary to check whether the scanning head and scanning head lenses are intact. If it needs to be used immediately after sterilization, it should be taken out and left to stand for more than 30 minutes to cool to room temperature before installation and operation.

Attention: The scanning head is a consumable item that can be sterilized 100 times in a high-temperature sterilization environment. When the scanning head shell ruptures, the contact spring falls off, the scanning head lens cracks, severe spots, and dirt cannot be cleaned, it needs to be disposed of.

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B. Soaking disinfection

Soak the scanning head completely in CIDEXOPA solution (0.55% Phthalaldehyde) for more than 5 minutes. After completion, dry the scanning head and scanning head lenses with a soft cloth (non-woven fabric) and a dust-free cotton swab, and use immediately to prevent secondary contamination of the scanning head. Before use, check whether the scanning head and scanning head lenses are intact.

CAUTION Attention: When soaking and disinfecting, the scanning head needs to be vertically placed in the solution, and it must be dried when taken out.

WARNING Attention: The scanning head needs to be sterilized and disinfected for the first use, and it also needs to be sterilized and disinfected for different patients in the future.

5.2.2 Scanner host

a) Cleaning and disinfection of scanner host


Wipe and clean the product shell with a small amount of soapy water dipped in a regular cotton cloth. After cleaning, soap and water should be removed, leaving no residue on the surface, and wiped dry with a clean, dry cotton cloth.

Dip a small amount of 75% medical alcohol onto a soft cloth (non-woven fabric) to wipe the surface of the scanner body. After a certain period of time, air dry it naturally or use another clean and dry soft cloth (non-woven fabric) to dry the residual alcohol.

The above suggestion is to clean and disinfect once a day.

CAUTION Attention: Do not use cleaning materials that can damage the surface of the casing for wiping.

Do not allow liquid to enter the interior of the equipment, causing mechanical damage. Pay special attention to the areas shown in the below picture.



b) Cleaning and disinfection of the front protection window of the scanner host

Gently wipe the surface of the protective window with a small amount of 75% medical alcohol dipped in a dust-free cotton swab to remove dirt, oil stains, spots, and other marks on the lens. Then use another clean and dry dust-free cotton swab to dry the surface of the protective window.

CAUTION Attention: The protective window is a precision optical component that must be carefully cleaned and disinfected, and be careful not to allow excess liquid to flow into other locations when wiping.

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Note1: The scanning head has a reflector made of glass material. Pay attention to the scanning angle and force to prevent damage caused by glass breakage;

Note2: The scanner head is not allowed to contact the scanner's protection window to insert the scanner head in the direction of the scanner body. Do not tilt. Carefully insert the scanner head until the slot is in place.

Note3: If the detection is different from normal, please go to the local distributor to further explain the product problem.

Note that the scanner head is placed vertically into the disinfectant solvent, and do not allow the solvent to enter the scanner head.

Necessary maintenance and precautions:

Operation	Cycle	Responsible person
Visual Inspection [5.1 1]	daily	operators
Visual Inspection [5.1 2]	A year	operators
Maintain Scanner	daily	operators
Cleaning and Disinfection [5.2.1]	daily	operators
Cleaning and Disinfection [5.2.2]	after each use	operators

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►Chapter 6 Common problems and solutions

6.1. After pressing the power button, the device cannot be turned on.

1. Check whether the adapter is in good contact and the power indicator is normal.
2. The equipment power indicator light, the computer cannot identify the equipment.
3. Check whether USB and the platform are plugged in properly.
4. Check whether the platform and computer USB plug is good.

6.2. The computer has identified the device, but the application has no image.

1. The device has dormant function, check whether the equipment is in hibernation state, shake and wake up.
2. Restart the application and reconnect with the device.
3. Turn off the device power and restart the device.

6.3. During the use of the device, an interrupt does not work.

1. The device has the function of thermal protection, and after the protection function, please turn off the device. After 5-10 minutes, the ambient temperature should be re-opened. During the use of the equipment, the camera fog.
2. Use after removing the fog, or put in the platform to heat the lens and continue to use.
3. This electronic equipment waterproof grade is IPX0, do not spray directly, immerse in various liquids.

6.4. During the use of the device, the USB image is slowed down or disconnected.

1. Check the USB cable connection see if it is good.
2. Check Whether the equipment is used for a long time leads to higher internal temperature.
3. If the above issue occurs frequently, it may be that the USB internal data cable is broken. Please contact the local after-sales service staff.

6.5. During the use of the equipment, the image noise is large or the image color is distorted.

1. This equipment belongs to the precision instrument, should be handled gently. Heavier shocks can cause unpredictable damage to internal hardware.
2. Internal software calibration of the equipment.
3. If the above issue occurs, please contact the local after-sales service staff.

6.6. During the use of the equipment, the scanner head lens produces fog, it cannot be scanned normally.

Check if the heater chip on the base fails or is unheated.

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►Chapter 7 Product application precautions.

1. The product is precision optical instrument equipment, need careful operation and maintenance; Accidental collision, falling may cause optical component damage and deviation, and ultimately affect the scan result. If any of the above occurs, please contact the customer service staff in time.

2. Use and maintenance of the scanner, it is necessary to pay attention to protect the leakage optical components, such as scanning prism, scanner ontology protection glass, etc., to avoid erosion corrosive liquid or sharp objects scratching the lens surface, if you have the above case, please contact customer service.

3. When scanner head sterilization/sterilization, it is necessary to avoid the infiltration of liquid from the back end of the scanner head. If there is infiltration and contamination of the inner surface of the lens, it is necessary to clean and blow dry with pure (99.9%) alcohol.


4. The scanner is used to pay attention to the joint, and the extra force makes the connection unstable or even disconnects. Also please be careful not to bend the wire to avoid damage.

5. When the scanner is not used, please remove the scanner head and cover the protective sleeve and close the base power. When not in use for a long time, pull out the power adapter.

6.7. 6. The scanner emits visible light when scanning. Please do not illuminate the eyes with scanning light, otherwise it will cause short-term discomfort.

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Camera resolution: 1.3MP CMOS

Pixel size: 4.8 μ m

Oral scanning light source: LED (without laser radiation)

Export development data format: OBJ/STL/PLY

Light radiant intensity: ≤ 100Mw/cm ^{2}

Size of oral digital impression instrument (excluding scanning head): 200 × 58 × 36mm

Data cable of Scanner (connecting USB control box to Scanner host): 1.9m

USB connection cable (host to computer connection): 0.7mF

Connection method: USB2.0/3.0

9.2. Equipment operating environment：

a} environment temperature：10 ℃- 40℃

b) Relative humidity：≤ 85%;

c} atmospheric pressure：86KPa~106KPa;

9.3. Equipment storage, transportation environment:

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►Chapter 10 Waste disposal

In order to reduce the burden on the environment, recyclable parts should be sent to the recycling center after removing the hazardous materials.

All components and components containing hazardous substances shall be disposed of in accordance with law and environmental provisions. When dealing with waste products, they must be protected from harm.

Also, local directive shall be followed.

△Recyclable ▲unrecyclable

Part	Main material	Recyclable material	Disposal center	Separation of harmful substances
Cover Metal	ABS aluminum	△ △		△
Circuit board		▲		
Wire	copper	△		
packing	paper	△		
Other			△	

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►Chapter 11 Electromagnetic compatibility

For this equipment, special precautions shall be taken for electromagnetic compatibility (EMC) and shall be installed and used in accordance with the EMC information specified in this manual. Portable and mobile radio frequency communications equipment may have an impact on the equipment.

Except as internal components of the spare parts for sale cable (transducer), using the specified attachment and cable (transducer) may lead to increase of equipment or system launch or immunity to reduce equipment or system should not be put to use with other devices or close to, if use must be close to or stacked, validation should be observed in the use of a configuration under normal operation.

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The following cables must be used to meet the requirements of electromagnetic emission and anti-interference:

Name of cable	Length
Power Cable	1. 9m
Cable (connection between base and scanner)	1. 9m
Cable (connection between scanner and computer)	1. 2m

The fundamental nature can be used for image acquisition.

Name	Description
Image acquisition	When you turn on the power, start the software and move the handle, the image display box on the software can display the image normally.

Guidance and manufacture's declaration - electromagnetic emissions		
The Intraoral Scanner is intended for use in the electromagnetic environment specified below. The customer or the user of the Intraoral Scanner should assure that it is used in such an electromagnetic environment:		
Emission Test	Compliance	Electromagnetic environment - guidance
RF emission CISPR 11	Group 1	The Intraoral Scanner uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The Intraoral Scanner is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emission IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emission IEC 61000-3-3	applied	

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Guidance and manufacturer's declaration - electromagnetic immunity.			
The Intraoral Scanner is intended for use in the electromagnetic environment specified below. The customer or the user of the Intraoral Scanner should assure that it is used in such an electromagnetic environment:			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV Differential mode ±2 kV common mode	±1 kV Differential mode ±2 kV common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95% dip in UT) for 0.5 cycle 40 % UT (60% dip in UT) for 5 cycles 70 % UT (30% dip in UT) for 25 cycles	<5 % UT (>95% dip in UT) for 0.5 cycle 40 % UT (60% dip in UT) for 5 cycles 70 % UT (30% dip in UT) for 25 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Intraoral Scanner requires continued operation during power man's interruption, it is recommended that the Intraoral Scanner be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: UT is the a.c. mains voltage prior to application of the test level.			

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<table>
<tr><td colspan="4">Guidance and manufacturer's declaration - electromagnetic immunity.</td></tr>
<tr><td colspan="4">The Intraoral Scanner is intended for use in the electromagnetic environment specified below. The customer or the user of the Intraoral Scanner should assure that it is used in such an electromagnetic environment:</td></tr>
<tr><td>Immunity test</td><td>IEC 60601 The test level</td><td>Compliance level</td><td>Electromagnetic environment - guidance</td></tr>
<tr><td>Conducted RF
IEC 61000-4-6</td><td>3 Vrms
150 kHz to 80 MHz</td><td>3 Vrms</td><td rowspan="2">Portable and mobile RF communications equipment should be used no closer to any part of the Intraoral Scanner, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Recommended separation distance
d = 1.2√P
d = 1.2√P 80 MHz to 800 MHz
d = 2.3√P 800 MHz to 2.5 GHz
Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey [a], should be less than the compliance level in each frequency range.[b]
Interference may occur in the vicinity of equipment marked with the following symbol:
</td></tr>
<tr><td>Radiated RF
IEC 61000-4-3</td><td>3 V/m
80 MHz to 2.5 GHz</td><td>3 V/m</td></tr>
<tr><td colspan="4">Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</td></tr>
</table>

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[a] Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitter transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Intraoral Scanner is used exceeds the applicable RF compliance level above, the Intraoral Scanner should be observed to verify normal operation. If abnormal performance observed, additional measures may be necessary, such as reorienting or relocating the Intraoral Scanner.

[b] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Intraoral Scanner

The Intraoral Scanner is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Intraoral Scanner can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Intraoral Scanner as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter /W	Separation distance according to frequency of transmitter /m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0. 12	0.12	0.23
0.1	0. 38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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