Abbott

Emboshield NAV6™

Embolic Protection System

CE 2797 UKCA 0086

en English:.....5
de Deutsch:.....8
fr Français :.....11
es Español:.....15
it Italiano:.....18
pt Português:.....22
sv Svenska:.....25
nl Nederlands:.....28
da Dansk:.....31
el Ελληνικά:.....35
pl Polski:.....38
hu Magyar:.....42
cs Česky:.....45
tr Türkçe:.....48
sk Slovensky:.....51
no Norsk:.....54
bg Български:.....57
ro Română:.....61
ru Русский:.....64
fi Suomi:.....68
et Eesti keel:.....71
lv Latviešu:.....74
lt Lietuvių:.....77
sl Slovenščina:.....81
uk Українська:.....84
ar العربية:.....88
hr Hrvatski:.....91
sr Srpski:.....94
mk Македонски:.....97

PPL2137364 (2024-01-02)

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INSTRUCTIONS FOR USE

Table of Contents

1.0 DEVICE DESCRIPTION
- 1.1 BareWire™ Filter Delivery Wire
- 1.2 Delivery Catheter
- 1.3 Filtration Element
- 1.4 RX Retrieval Catheter

2.0 HOW SUPPLIED
3.0 INTENDED PURPOSE
4.0 INTENDED USER AND USER TRAINING
5.0 INTENDED CLINICAL BENEFIT
6.0 INDICATIONS
7.0 CONTRAINDICATIONS
8.0 WARNINGS
9.0 PRECAUTIONS
10.0 POTENTIAL ADVERSE EVENTS
11.0 REPORTING OF ADVERSE EVENTS
12.0 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE
13.0 EXPECTED LIFETIME OF THE DEVICE
14.0 DEVICE DISPOSAL
15.0 PATIENT SELECTION AND TREATMENT
16.0 PATIENT COUNSELING AND PATIENT INFORMATION
17.0 CLINICIAN USE INFORMATION
- 17.1 Materials Required
- 17.2 General Preparation and Use (Aseptic Technique)
 - 17.2.1 Sizing Information
 - 17.2.2 Filter Loading
 - 17.2.3 Using Contrast Media
 - 17.2.4 Preparation for Deployment
 - 17.2.5 Device Deployment
 - 17.2.6 Device Retrieval

18.0 PRODUCT INFORMATION DISCLOSURE

This Instructions for Use pertains to the Emboshield NAV6™ Embolic Protection System and associated BareWire™ Filter Delivery Wire.

Table 1: Basic UDI-DI

Product Name	Basic UDI-DI
Emboshield NAV6 Embolic Protection System	8717648EPS0001YC
BareWire Filter Delivery Wires	8717648EPS0002YE

This Instructions for Use (IFU) as per European Regulation (EU) 2017/745 can be found on the Abbott eLabeling website at vascular.elFU.abbott.

1.0 DEVICE DESCRIPTION

The Emboshield NAV6™ Embolic Protection System is a temporary percutaneous transluminal filtration system designed to be used as a guide wire and to capture embolic material released during an interventional procedure within the carotid artery or the arterial vasculature of the lower extremity. The system consists of the following components:

1.1 BareWire™ Filter Delivery Wire

The BareWire Filter Delivery Wire (**Figure 1**) is a 0.014" (0.36 mm) polytetrafluoroethylene (PTFE) coated stainless steel wire with a 3 cm platinum / nickel radiopaque distal tip section. The Filter Delivery Wires in the BareWire range have zebra stripes on the proximal end for ease of identification.


Figure 1

The Emboshield NAV6 Embolic Protection System is supplied pre-loaded with the BareWire™ Workhorse Filter Delivery Wire (190 cm). An alternative range of BareWire Filter Delivery Wires are also available for use with the Emboshield NAV6 Embolic Protection System (**Table 2**). Only the BareWire Filter Delivery Wires shown in **Table 2** are intended for use with the Emboshield NAV6 System.

Table 2: BareWire Filter Delivery Wires

Description	Length	Product Code
BareWire Distal Access	190 cm	22444-19
BareWire Workhorse	190 cm	22445-19
BareWire Workhorse	315 cm	22445-31
BareWire Support	190 cm	22446-19

1.2 Delivery Catheter

The RX Delivery Catheter (**Figure 2**) usable length is 135 cm. The crossing profile is either 0.036" (0.91 mm) or 0.041" (1.04 mm), depending on Filtration Element size. A pull handle is used to deploy the loaded Filtration Element from the pod. Two pairs of silver indicator bands are provided along the catheter shaft. A proximal pair of indicator bands indicates the catheter tip position during advancement through the guide catheter / introducer sheath, and a distal pair of indicator bands indicates the proximity of the RX exit port during catheter retraction. A radiopaque marker band is positioned proximal to the pod. An enamel coating is applied to the shaft of the Delivery Catheter.


Figure 2

1.3 Filtration Element

The Filtration Element (**Figure 3** and **Table 3**) consists of a nylon membrane with a hydrophilic coating and an internal nitinol frame. There are two proximal entry ports and multiple distal perfusion pores. There are proximal and distal marker bands. The Filtration Element is available in two sizes that fit a range of vessel sizes: Small, for use in 2.5 mm - 4.8 mm diameter vessels, and Large, for use in 4.0 mm - 7.0 mm diameter vessels. To aid in visualization, the Filtration Element frame is also radiopaque.


Figure 3

Table 3: Filtration Element

Filtration Element	Dim (A) cm	Dim (B) cm	Dim (C) cm	Vessel Size mm	Unconstrained Filtration Element (FE) Diameter mm
Small	1.9	0.2	0.15	2.5 - 4.8	5.0
Large	2.3	0.2	0.15	4.0 - 7.0	7.2

1.4 RX Retrieval Catheter

The RX Retrieval Catheter (**Figure 4**) has a 139 cm usable length, and a distal radiopaque expansile tip. The maximum outer diameter of the distal tip is 0.067" (1.70 mm). A handle grip is situated at the proximal end. Two pairs of silver indicator bands are provided along the catheter shaft. A proximal pair indicates the catheter tip position during advancement through the guide catheter / introducer sheath, and a distal pair indicates the proximity of the RX exit port during catheter retraction.


Figure 4

2.0 HOW SUPPLIED

Sterile: This device is sterilized by ethylene oxide. Non-pyrogenic.

Contents: Each Emboshield NAV6™ Embolic Protection System contains two pouches with the following content:

Contents of Filtration System primary pack (pouch)	RX Delivery Catheter Filtration Element (FE) Loading Funnel* Syringe BareWire™ Filter Delivery Wire Introducer Tool Torque Device
Contents of RX Retrieval Catheter primary pack (pouch)	RX Retrieval Catheter Torque Device Syringe Tip

* Provided as part of Filtration System Pack Tray

3.0 INTENDED PURPOSE

The Emboshield NAV6™ Embolic Protection System is intended to capture and remove embolic material during interventional procedures.

4.0 INTENDED USER AND USER TRAINING

The intended users of the subject devices are physicians and catheterization lab staff who have been appropriately trained in interventional procedures of carotid and / or lower extremity arteries, and who have received appropriate training on the principles, clinical applications, complications, side effects and hazards commonly associated with percutaneous transluminal angioplasty (PTA), atherectomy and stenting procedures in carotid, and / or lower extremity arteries.

5.0 INTENDED CLINICAL BENEFIT

The following are the intended clinical benefits of the Emboshield NAV6™ Embolic Protection System:

- Low rates of stroke for carotid indication, and
- Low rates of distal embolization for lower extremity indication

6.0 INDICATIONS

The system is designed to reduce distal embolization while maintaining blood flow during interventional procedures.

The system is indicated for use during percutaneous transluminal angioplasty, and stenting procedures for patients with atherosclerotic disease in carotid arteries.

The system is indicated for use during percutaneous transluminal angioplasty, atherectomy and stenting procedures for patients with atherosclerotic disease in the lower extremity arteries.

7.0 CONTRAINDICATIONS

The Emboshield NAV6™ Embolic Protection System is contraindicated for use in:

- Patients in whom anticoagulant and / or antiplatelet therapy is contraindicated.
- Patients with a known allergy or hypersensitivity to device materials (nickel-titanium alloy) or contrast medium, who cannot be adequately premedicated.
- Patients with uncorrected bleeding disorders.
- Inability to cross the lesion with BareWire™ Filter Delivery Wire.
- Diffusely diseased vessels where there is no disease-free section in which to deploy the Filtration Element.
- Insufficient straight section of vessel distal to the lesion to permit Filtration Element deployment.

8.0 WARNINGS

- Refer to instructions supplied with all interventional devices to be used with the Emboshield NAV6™ Embolic Protection System for their intended uses, contraindications, and potential complications.
- The Emboshield NAV6 Embolic Protection System is supplied sterile. Do not use if the package has been opened or is damaged. Carefully inspect the system components prior to use to verify that they have not been damaged and that the size, shape and condition are suitable for the procedure for which they are to be used. A device or access device that is kinked or damaged in any way should not be used.
- The Emboshield NAV6 device can only be used with the BareWire™ Filter Delivery Wire. Use of the device with any guide wire other than the BareWire Filter Delivery Wire will lead to loss of the Filtration Element during the procedure or an inability to retrieve the Filtration Element.
- This device is designed and intended for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and / or delivery system and / or lead to device failure, which may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and / or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious diseases(s) from one patient to another. Contamination of the device and / or delivery system may lead to injury, illness or death of the patient.
- Keep dry. Keep away from sunlight. Do not use the product after the "Use by" Date specified on the label.
- Use of the device should be restricted to physicians trained to the specifics of the device and to the Instructions for Use. Physicians must be knowledgeable of the current medical literature and appropriately trained on the principles, clinical applications, complications, side effects and hazards commonly associated with carotid and / or lower extremity interventional procedures.
- Overstretching of the artery may result in rupture and life-threatening bleeding.
- Appropriate antiplatelet, anticoagulant and, if necessary, vasodilator therapy must be used during the procedure. Anticoagulant therapy sufficient to maintain an Activated Clotting Time of at least 250 seconds for the duration of the procedure is recommended.
- Adequate guide catheter or sheath support is required for the introduction of the RX Delivery Catheter, RX Retrieval Catheter and all interventional devices to be used during the procedure.
- Maintain a snug seal between the device and the hemostatic valve during catheter insertion. Failure to observe this may result in air being drawn into the access device through the hemostatic valve. Device insertion should be performed slowly to minimize the risk of air entrainment. It is therefore recommended that flushing of contrast media (or other fluids) is performed before or after insertion of the catheter, but not while the catheter is within the access device.
- Do not advance any component of the Emboshield NAV6 Embolic Protection System against significant resistance. The cause of any resistance should be determined via fluoroscopy and remedial action taken.
- Torquing the BareWire Filter Delivery Wire against resistance may cause BareWire Filter Delivery Wire damage and / or BareWire Filter Delivery Wire tip separation.
- Perform all exchanges slowly to prevent air embolism or trauma to the artery.
- The Delivery Catheter should not be pulled back quickly in the event that resistance is felt during catheter advancement. Rapid withdrawal of the Delivery Catheter in the presence of resistance may result in premature deployment of the Filtration Element.
- In the event of complications surgical intervention may be required.
- A high-pressure contrast injection may cause Filtration Element movement in the vessel.
- Do not attempt to move the Filtration Element after deployment and prior to the start of retrieval.
- Do not torque the RX Delivery and RX Retrieval Catheters.
- Allow and maintain adequate distance between the Emboshield NAV6 Filtration Element and other interventional devices (including carotid stent systems and balloons) to avoid engagement.
- Use of the stent system or the embolic protection system in the following pre-existing conditions may result in serious patient harms:
 - Hemoglobin (Hgb) < 8 gm/dl (unless on dialysis), platelet count < 50,000, INR > 1.5 (irreversible), or heparin-associated thrombocytopenia.
 - Evidence of a stroke within the previous 30 days.

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- History of ipsilateral stroke with fluctuating neurologic symptoms within 1 year.
- Patients with total occlusion of the target vessel.
- Patients with evidence of intraluminal thrombus thought to increase the risk of plaque fragmentation and distal embolization.
- Any condition that precluded proper angiographic assessment or made percutaneous arterial access unsafe, (e.g., morbid obesity, sustained systolic blood pressure > 180 mmHg).
- Known cardiac sources of emboli.
- Patients with highly calcified lesions resistant to PTA.
- Atherosclerotic disease involving adjoining vessels (i.e., the aortic arch or ostial common carotid artery) precluding safe placement of the guiding catheter or sheath.
- Patients with severe vascular tortuosity.
- Patients with in-stent restenosis.
- The safety and effectiveness of concurrent treatment of lesions in patients with bilateral carotid artery disease have not been established.

- The safety and effectiveness of the Emboshield NAV^6 Embolic Protection System has not been established in pregnant patients and patients under the age of 18.
- Caution should be used if pre-dilating the lesion without embolic protection as this may increase the risk of an adverse outcome.
- The BareWire Filter Delivery Wire should be kept moistened throughout the procedure.
- The Filtration Element can only be positioned proximal to the radiopaque distal section of the Filter Delivery Wire.
- To reduce the potential for the liberation of emboli during lesion crossing, the device should be carefully manipulated and not advanced against resistance.
- If the Filtration Element moves into the stented vessel segment prior to retrieval, DO NOT RETRIEVE. Use the Retrieval Catheter to gently maneuver the Filtration Element distally until it is situated in an unstented portion of vessel. Retrieval should then proceed.
- Maintain proper guiding catheter / sheath support throughout the procedure. Ensure that there is enough distance between the proximal tip of the Filtration Element and the most distal tip of any interventional device to be introduced over the Filter Delivery Wire to avoid engagement. The tip of a balloon catheter or a stent delivery system or an atherectomy device should not contact the Filtration Element. Failure to maintain adequate distance could result in inadvertent Filtration Element movement and Stent Delivery System tip / Filtration Element entanglement and / or Filtration Element / Stent entanglement if guide catheter or sheath prolapse occurs.
- In patients requiring the use of antacids and / or H2-antagonists before or immediately after stent placement, oral absorption of antiplatelet agents (e.g., aspirin) may be adversely affected.
- The minimum guide catheter or sheath internal diameter is printed on the package label. Do not use a smaller guide catheter or sheath than indicated.
- For proper positioning of the filter basket, the vessel distal to the lesion should have an absence of excessive tortuosity and be of adequate length (approximately 4 cm distal to the lesion and proximal to the petrous portion of the vessel).
- Persons with known history of allergies to any of the components of the Emboshield NAV^6 Embolic Protection System listed below may suffer an allergic reaction to this device. Prior to use on the patient, the patient should be counselled on the materials contained in the device and a thorough history of allergies must be discussed. This device contains nickel-titanium alloy, platinum-tungsten alloy, platinum-iridium alloy, polyimide, polytetrafluoroethylene, polyethylene, polyethylene naphthalene, polyamide, gold, polyurethane, polyvinylpyrrolidone (PVP) coating, stainless steel, tungsten, fluorinated ethylene propylene, polypropylene, polyether block amide, silicone, and perfluoropolyether.

9.0 PRECAUTIONS

- The device must only be flushed using the 5-ml syringe and flushing tip provided.
- Use with fixed (passive) hemostatic valves is not recommended.
- For best device performance, the guide wire exit notch should remain within the guiding catheter or sheath.
- The delivery system is not designed for use with power injection. Use of power injection may adversely affect device performance.
- Care must be used when removing the filter basket through a newly deployed stent to maintain filter basket integrity and to avoid disrupting the stent geometry.
- To optimize system performance, it is imperative that careful attention is paid to the preparation of the system using the preparation techniques specified in the Instructions For Use.
- The size of the Filtration Element should be selected to match the diameter of the vessel at the intended site of deployment. The clinician should take into consideration the slight increase in vessel diameter that may occur after treatment of the lesion.
- Precautions to prevent or reduce clotting should be taken when any interventional device is used. Flush or rinse all devices entering the vascular system with heparinised normal saline or alternative anticoagulant, prior to use.
- The Emboshield NAV^6™ Embolic Protection System must be used with a guide catheter or introducer sheath to maintain adequate support for the BareWire™ Filter Delivery Wire throughout the procedure.
- The maintenance of blood flow through the device should be observed throughout the procedure by the use of contrast injection.
- If excessive debris is collected in the Filtration Element such that distal perfusion of contrast media is significantly reduced or no contrast media is perfusing past the filter, the Emboshield NAV^6 Embolic Protection System may have reached its maximum capacity to contain emboli. Remove and replace the Emboshield NAV^6 Embolic Protection System. Otherwise, it may be difficult to completely recover all embolic debris and the potential for thrombus formation may increase.
- Venous access should be available during carotid stenting in order to manage bradycardia and / or hypotension by either pharmaceutical intervention or placement of a temporary pacemaker, if needed.
- Removal of the BareWire Filter Delivery Wire with the Emboshield NAV^6 Filtration Element through any interventional devices other than the Emboshield NAV^6 RX Retrieval Catheter has not been tested.
- The minimum expanded stent internal diameter required for retrieval of a large embolic load is 2.5 mm.
- Advancement and deployment of the carotid stent system and embolic protection system should only be performed under fluoroscopic observation.
- When introducing the delivery system, confirm that the wire tip is free within the vessel lumen and is not directed into the vessel wall. Failure to do so may result in vessel trauma. Use the radiopaque marker on the interventional device to confirm position.
- Avoid excessive movement of the filter basket during catheter device exchanges. Excessive movement of the deployed basket may cause vessel trauma or spasm.
- One or more components of this device may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; Chemical Abstract Service (CAS) No. 7440-48-4; EC No. 231-158-0.

 Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- A component of this device may contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight: 1-methyl-2-pyrrolidone (NMP); Chemical Abstract Service (CAS) No. 872-50-4; EC No. 212-828-1. In animal tests, continuous NMP exposure has shown reproductive toxicity.

 There is no clinical study for this device in terms of reproductive toxicity in pregnant women or in men intending to father children. Effects on the developing fetus have not been studied. In the worst-case scenario, i.e. four Emboshield NAV^6 Embolic Protection devices used in a procedure and 8 hours patient exposure, the amount of NMP that may elute out is 9 times lower than NMP Permissible Daily Exposure by the ICH Q3C(R7) Guidelines, or 16 times lower than daily tolerable exposure. Occurrence of this worst-case scenario is unlikely, as the number of the devices and the procedure time are much lower than the worst-case scenario. While there is no contraindication, the risks and reproductive effects of use of this device on humans are unknown at this time. It is recommended that the need of the patient and the clinical benefit of Emboshield NAV^6 Embolic Protection System be weighed against the risk stated above.

10.0 POTENTIAL ADVERSE EVENTS

- Allergic reaction or hypersensitivity to latex, contrast agent, anesthesia, device materials (nitinol, nickel, titanium), and drug reactions to anticoagulation, or antiplatelet drugs
- Vascular access complications which may require transfusion or vessel repair, including:
 - Catheter site reactions
 - Bleeding (hemorrhage, ecchymosis, oozing, hematoma, retroperitoneal hemorrhage)
 - Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture, and laceration
 - Embolism (air, tissue, plaque, thrombotic material or device)
 - Thrombophlebitis
- Target artery complications which may require additional intervention, including:
 - Total occlusion or abrupt closure
 - Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
 - Embolism (air, tissue, plaque, thrombotic material or device)
 - Artery, stent, or filter thrombosis / occlusion thrombosis
 - Stenosis or restenosis
 - Vessel spasm
- Cardiac arrhythmias (including conduction disorders, atrial and ventricular arrhythmias)
- Cardiac ischemic conditions (including myocardial ischemia, myocardial infarction, and unstable or stable angina pectoris)
- Stroke / cerebrovascular accident (CVA) and transient ischemic attack (TIA)
- System organ failure:
 - Cardio-respiratory arrest
 - Cardiac failure
 - Cardiopulmonary failure (including pulmonary edema)
 - Renal failure / insufficiency
 - Shock
- Bleeding
- Blood cell disorders, including heparin-induced thrombocytopenia (HIT) and other coagulopathy
- Hypotension / hypertension
- Peripheral nerve injury
- Other ischemic conditions / infarct
- Target limb loss / necrosis (amputation of toe, foot, and / or leg)
- Infection – local and systemic (including post-procedural)
- Nausea and vomiting
- Chest pain
- Edema / cerebral edema and fluid overload
- Fever
- Pain, including headache
- Hyperperfusion syndrome
- Other neurologic and systemic complications
 - Seizure
- Cerebral hemorrhage
- Death

11.0 REPORTING OF ADVERSE EVENTS

For customers in the European Union, the user and / or patient should report any serious incident that has occurred in relation to the Emboshield NAV^6™ Embolic Protection System to the manufacturer and the competent authority of the European Member State in which the user and / or patient is established. For customers in other geographies, report any serious incident to the manufacturer and the regulatory authorities responsible for the jurisdiction.

12.0 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

A summary of the safety and clinical performance (SSCP) for this device is available at https://ec.europa.eu/tools/eudamed. Search for the device using the Basic UDI-DI provided in **Table 1**. This is the SSCP location after the launch of the European Database on Medical Devices (EUDAMED).

13.0 EXPECTED LIFETIME OF THE DEVICE

Emboshield NAV^6™ Embolic Protection System is intended for use during a single interventional procedure. The expected lifetime of the Emboshield NAV^6 Embolic Protection System typically ranges between 10 minutes to two hours and is dependent on the lesion and procedure complexity.

14.0 DEVICE DISPOSAL

After use, this device, its accessories and packaging should be appropriately classified for disposal (e.g., biohazard, sharps, non-hazardous waste, etc.), and carefully disposed of in compliance with facility procedures and applicable laws and regulations.

15.0 PATIENT SELECTION AND TREATMENT

Patients who are deemed, by a qualified physician, to be suitable for and within an acceptable risk profile for the indicated procedures as described in **Section 6.0 INDICATIONS** and **Section 7.0 CONTRAINDICATIONS**.

16.0 PATIENT COUNSELING AND PATIENT INFORMATION

Physicians should review warnings, precautions, and adverse events for residual risks and possible side effects or hypersensitivity risks that may be relevant when counseling patients about use of this product and should consider the following:

- Discussion of the risk / benefit issues for this particular patient
- Discussion of alteration to current lifestyle immediately following the procedure and over the long term
- Discussion of any necessary follow up
- Discuss the risks of allergic reaction or hypersensitivity to device components

17.0 CLINICIAN USE INFORMATION

17.1 Materials Required

The Emboshield NAV^6™ Embolic Protection System can only be used with the BareWire™ Filter Delivery Wire. Use of the device with any guide wire other than the BareWire Filter Delivery Wire will lead to loss of the Filtration Element during the procedure or an inability to retrieve the Filtration Element.

Confirm the compatibility of the Emboshield NAV^6 Embolic Protection System with the interventional devices before actual use.

- Minimum guiding catheter / sheath size inner diameter (I.D.) 0.070" / 1.77 mm
- 0.096" (2.44 mm) Rotating Hemostatic Valve (RHV) (optional)
- X-ACT™ or Acculink™ Carotid Stent System / Peripheral atherectomy devices (optional)
- Balloon dilatation catheter (optional)
- 1,000 u / 500 cc heparinised normal saline (sterile)
- 10-20 cc luer-lock syringe

Review **Section 2.0 HOW SUPPLIED** for a complete list of device components and ancillary items provided with the Emboshield NAV^6 Embolic Protection System.

17.2 General Preparation and Use (Aseptic Technique)

17.2.1 Sizing Information

The size of the Filtration Element for placement in the carotid arteries or the arterial vasculature of the lower extremity below the common femoral artery is determined by the diameter of the vessel at the intended deployment site. Quantitative coronary arteriography (QCA) or intravascular imaging (e.g., optical coherence tomography [OCT], intravascular ultrasound [IVUS]) is recommended.

17.2.2 Filter Loading

- **Do not move the Delivery Catheter during loading. Movement may cause misalignment of system components.**
- Fill the recessed well area with heparinised normal saline (**Figure 5**) until the Filtration Element is completely immersed (**Figure 6**).
- Fill the provided syringe with heparinised normal saline in order to flush the Delivery Catheter. Remove air from syringe. Push the syringe tip into the core of the Loading Funnel until it engages the internal diameter of the Loading Funnel (**Figure 6**). Flush until heparinised normal saline can be seen exiting at the RX exit port (**Figure 7**).
- Use the syringe with the flushing tip to remove any air bubbles from the Filtration Element. Place the white tip gently inside the Filtration Element and flush (**Figure 8**).


Figure 5


Figure 6


Figure 7


Figure 8

Load the Filtration Element into the Delivery Catheter as follows:

- Ensure that the tray is on a level surface. Hold the tray securely with one hand ensuring the well and Filtration Element are clearly visible (**Figure 9**).
- Unclip the torque device (**Figure 10**). Ensure that the torque device is securely tightened on the coating free section of the BareWire Filter Delivery Wire (**Figure 1**). The torque device can be re-tightened, if required. Carefully pull the torque device in a controlled manner to retract the BareWire Filter Delivery Wire and observe the Filtration Element being retracted through the Loading Funnel into the Delivery

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Catheter Pod. Pulling the wire too quickly or with excessive force may cause damage to the RX Delivery Catheter Pod and / or Filtration Element. The device is fully loaded when the distal marker band on the Filtration Element is at or inside the end of the delivery catheter pod and a definite stop is felt (**Figure 11**).

- **Do not continue to retract the BareWire Filter Delivery Wire against significant resistance.**
- Inspect the Delivery Catheter Pod prior to removal from the tray to ensure the Filtration Element is fully loaded (**Figure 11**).
- Remove the loaded Delivery Catheter and the BareWire Filter Delivery Wire from the tray by gripping the black handle and the Torque Device together and gently pulling them out of the tray (**Figure 12**).
- **Do not attempt to remove either the Delivery Catheter or the Filter Delivery Wire individually.**
- Inspect the Delivery Catheter to ensure the Filtration Element remains loaded properly into the pod at the end of the Delivery Catheter and that the BareWire Filter Delivery Wire is not damaged.
- Wipe the length of exposed Filter Delivery Wire from the RX port to the torque device with a sterile gauze pad soaked with heparinised normal saline.
- The Delivery System is now ready for use.


Figure 9


Figure 10


Figure 11


Figure 12

17.2.3 Using Contrast Media

During the insertion of rapid exchange catheters through guide catheters or sheaths, careful handling is required to ensure that air is not drawn into the access device.

It is therefore recommended that flushing of contrast media (or other fluids) is performed before or after insertion of the catheter, but not while the catheter is within the access device.

In the case where a contrast media injection must be performed with the catheter in place, it is essential to ensure that no air is present within the access device prior to injection.

If aspiration is to be performed prior to contrast media injection, it should be performed slowly and steadily at a rate of not more than 0.5 ml (0.5 cc) per second until the user can visually confirm that no further air is entering the aspiration syringe (**Figure 13**).

Following all instructions, slow catheter insertion, and good hemostatic valve control will help to minimize the risk of air embolization.


Figure 13

17.2.4 Preparation for Deployment

WARNINGS:

Only use the BareWire Filter Delivery Wires listed in Table 2. Use of other guide wires will lead to loss of the Filtration Element or an inability to retrieve the Filtration Element.

When using Over-the-Wire interventional equipment over 80 cm maximum length, the Filtration Element should be delivered on the BareWire Workhorse 315 cm wire ONLY.

All BareWire Filter Delivery Wires listed in Table 2 are NOT compatible with the Abbott Medical DOC™ Guide Wire Extension.

The Emboshield NAV^6 Embolic Protection System is NOT to be deployed with an access device that uses an integrated leaflet-type valve.

The Emboshield NAV^6 Embolic Protection System should be used with a hemostatic valve with minimum internal diameter of 0.096" (2.44 mm).

The Emboshield NAV^6 Embolic Protection System should be used with a guide catheter or sheath access device with an internal diameter equal to or greater than 0.070" (1.77 mm).

A device or access device that is kinked or damaged in any way should not be used. Assess all equipment before use.

Do not advance any component of the Emboshield NAV^6 Embolic Protection System against significant resistance. The cause of any resistance should be determined via fluoroscopy and remedial action taken.

Maintain a snug seal between the device and the hemostatic valve during catheter insertion. Catheter insertion should be performed slowly to minimize the risk of air entrainment. Failure to observe this may result in air being drawn into the access device through the hemostatic valve.

Perform all device manipulation using fluoroscopic visualization.

Do not torque or pull back on the RX Delivery Catheter during advancement.

Do not pull back on the BareWire Filter Delivery Wire during advancement.

Torquing the BareWire Filter Delivery Wire against resistance may cause BareWire Filter Delivery Wire damage and / or BareWire Filter Delivery Wire tip separation.

- Confirm the compatibility of the Emboshield NAV^6 Embolic Protection System with the other interventional devices intended for use in the procedure.
- The BareWire Filter Delivery Wire should be kept moistened throughout the procedure.
- A range of alternative BareWire Filter Delivery Wires is available (**Table 2**).
- The Filtration Element can only be positioned proximal to the radiopaque distal section of the BareWire Filter Delivery Wire.

Device Introduction

- Use the introducer tool provided to pass the loaded Delivery Catheter and Filter Delivery Wire through the hemostatic valve of the guide catheter or sheath. Indicator bands are positioned 90 cm and 100 cm from the tip to assist the clinician in anticipating when the catheter tip will exit the distal end of the guide catheter / introducer sheath.
- Advance the device through the guide catheter / introducer sheath to a location proximal to the lesion. Advance the Filter Delivery Wire across the lesion.

Use of an Alternative BareWire Filter Delivery Wire

- Load the device with the pre-packed BareWire Filter Delivery Wire as previously described.
- To use an alternative BareWire Filter Delivery Wire, carefully remove the torque device and the Filter Delivery Wire and dispose of the Filter Delivery Wire used to load the device (**Figure 14**).
- Remove the replacement BareWire Filter Delivery Wire from its hoop and wipe with a sterile gauze pad soaked with heparinised normal saline.
- Use an introducer tool to pass the wire through the hemostatic valve of the guide catheter or introducer sheath toward the deployment site. **(The BareWire Filter Delivery Wire should be kept moistened throughout the procedure)**.
- Advance the BareWire Filter Delivery Wire across the lesion.
- Carefully advance the Delivery Catheter over the BareWire Filter Delivery Wire (**Figure 15**). The BareWire Filter Delivery Wire will exit the Delivery Catheter at the catheter wire exit port.
- Holding the BareWire Filter Delivery Wire, advance the Delivery Catheter through the hemostatic valve toward the deployment site. Indicator bands are positioned 90 cm and 100 cm from the catheter tip to assist the clinician in anticipating when the catheter tip will exit the distal end of the Guide Catheter. Do not torque or pull back on the RX Delivery Catheter during advancement. Advance the Delivery Catheter to a location proximal to the lesion; then, continue with the Device Placement Instructions below.


Figure 14


Figure 15

Device Placement

- Position the 3 cm radiopaque portion of the Filter Delivery Wire in a suitable section of the vessel distal to the intended site of Filtration Element deployment.
- **(The Filtration Element can only be positioned proximal to the radiopaque distal section).**
- Carefully advance the Delivery Catheter over the BareWire Filter Delivery Wire to the intended filter deployment site. Do not torque or pull back on the RX Delivery Catheter during advancement.
- Once the Delivery Catheter is distal to the lesion, establish desired stent location and note where the distal edge of the stent will be placed. Using the marker bands on the Filtration Element and Delivery Catheter as guides, position the Delivery Catheter accordingly.
- Ensure that there is enough distance between the proximal tip of the Filtration Element and the most distal tip of any interventional device to be introduced over the Filter Delivery Wire to avoid contact during the procedure. Ensure that there is enough distance between the proximal end of the 0.014" (0.36 mm) Radiopaque Distal Tip of the Filter Delivery Wire and the distal tip of the Filtration Element to avoid unnecessary contact during the procedure. **Atherectomy Device / Stent Delivery System / Balloon Catheter distal tips may vary in size and radiopacity. The tip of a balloon catheter or a stent delivery system or an atherectomy device should not contact the Filtration Element (Figure 16).**
- Once the Delivery Catheter is positioned, the Filtration Element may be deployed.


Figure 16

17.2.5 Device Deployment

- Remove the red Locking Clip from the black handle of the Delivery Catheter.
- Hold the black handle steady and pull back on the white Pull Handle (**Figure 2**) to deploy the Filtration Element. Deployment is occurring when the distal marker band on the Catheter moves proximally relative to the Filtration Element marker bands.
- **WARNING: Do not attempt to move the Filtration Element after deployment and prior to the start of retrieval.**
- **WARNING: A high-pressure contrast media injection may cause Filtration Element movement in the vessel.**
- Remove the Delivery Catheter while maintaining the BareWire Filter Delivery Wire position. Silver indicator bands on the Delivery Catheter shaft help indicate the position of the BareWire Filter Delivery Wire exit port. A single silver band is located 15 cm from the exit port and twin silver bands are located 5 cm from the exit port. Open the hemostatic valve after the twin indicator bands on the hypotube shaft become visible.
- Discard the Delivery Catheter after withdrawal. Perform the required interventions over the BareWire Filter Delivery Wire.
- Maintain adequate distance between the proximal tip of the Filtration Element and the most distal tip of any interventional device.
- Monitor if excessive debris is collected in the Filtration Element such that distal profusion of dye is significantly reduced or no dye is perfusing past the filter; the Emboshield NAV^6 device may have reached its maximum capacity to contain emboli. Remove and replace the Emboshield NAV^6 device, otherwise it may be difficult to completely recover all embolic debris, and the potential for thrombus formation may increase.

Note: Filtration Element may not be within field of view.

17.2.6 Device Retrieval

WARNING: Do not attempt to use the Delivery Catheter to retrieve the Filtration Element.

- Remove the Retrieval Catheter tray from the pouch. Remove the Catheter from the tray.
- Flush the catheter lumen by injecting heparinised normal saline through the catheter's distal tip using the flushing tip provided in the tray and a standard syringe. Check that saline is exiting at the RX Port.
- Under fluoroscopy, advance the Retrieval Catheter over the BareWire Filter Delivery Wire. Indicator bands are positioned 90 cm and 100 cm from the tip to assist the clinician in anticipating when the catheter tip will exit the distal end of the guide catheter / introducer sheath.
- Advance through the stent / treated lesion until the radiopaque catheter tip contacts the proximal end of the Filtration Element. Do not torque the Retrieval Catheter during advancement.
- Grip the BareWire Filter Delivery Wire using the torque device provided in the tray, and retract the BareWire Filter Delivery Wire until the proximal radiopaque section of the wire meets the distal end of the Filtration Element.

The standard method of retrieval is as follows:

- Hold the Retrieval Catheter position steady by gripping the shaft close to the hemostatic valve.
- Pull on the tightened torque device to retract the BareWire Filter Delivery Wire until the Filtration Element is fully enclosed in the radiopaque expandable tip.
- Complete retrieval of the Filtration Element has been achieved when the distal Filtration Element marker band has been withdrawn to within 2 mm of the distal end of the radiopaque catheter tip (**Figure 17.3**).
- **Do not continue to retract the BareWire Filter Delivery Wire against significant resistance.**
- Retract the Retrieval Catheter over the BareWire Filter Delivery Wire (or remove Retrieval Catheter and BareWire Filter Delivery Wire as one unit) until the tip of the Retrieval Catheter is adjacent to the tip of the guide catheter or sheath. See **Figure 17.1**, **Figure 17.2** and **Figure 17.3** for examples of the Filtration Element prior to retraction into the Retrieval Catheter, partially retracted into the Retrieval Catheter, and fully retracted into the Retrieval Catheter.
- If significant resistance is felt, retract the guide catheter or sheath and the Retrieval Catheter together.
- Remove the Retrieval Catheter from the patient.
- If a hemostatic valve is used, ensure the valve is fully open during removal of the RX section of the Retrieval Catheter.
- **If an exceptionally large embolic load prevents complete withdrawal of the Filtration Element into the Retrieval Catheter tip, the following method of Retrieval must be employed:**
- Hold the Retrieval Catheter position steady by gripping the shaft close to the hemostatic valve.
- Pull on the tightened torque device to retract the BareWire Filter Delivery Wire and Filtration Element into the Retrieval Catheter tip.
- Do not continue to retract the BareWire Filter Delivery Wire against significant resistance.
- **WARNING: If the Filtration Element cannot be fully retracted as described above, the BareWire Filter Delivery Wire must be removed with the Retrieval Catheter.** Failure to do so may result in redeployment of the Filtration Element from the Retrieval Catheter tip. The proximal openings of the Filtration Element have been fully contained within the catheter tip once the radiopaque frame has been withdrawn into the radiopaque catheter tip (**Figure 17.2**).
- Retract the Retrieval Catheter, Filtration Element and BareWire Filter Delivery Wire as one unit until the tip of the Retrieval Catheter is adjacent to the tip of the guide catheter or sheath.
- Pull the tip of the Retrieval Catheter into and through the guide catheter / introducer sheath. If significant resistance is felt, retract the guide catheter or sheath and the Retrieval Catheter, Filtration Element and BareWire Filter Delivery Wire together.
- Remove the Retrieval Catheter from the patient. If a hemostatic valve is used, ensure the valve is fully open during removal of the RX section and tip of the Retrieval Catheter and retrieved Filtration Element.

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
Figure 17.1 **Figure 17.2** **Figure 17.3**

18.0 PRODUCT INFORMATION DISCLOSURE

Abbott has exercised reasonable care in the manufacturing of this device. Abbott excludes all warranties, whether expressed or implied, by operation of law or otherwise, including but not limited to, any implied warranties of merchantability of fitness, since handling and storage of this device, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of Abbott directly affect this device and the results obtained from its use. Abbott shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from the use of this device. Abbott neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

Reference Abbott website for patent markings: www.abbott.com/patents

™ Indicates a trademark of the Abbott group of companies.

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© 2023 Abbott. All Rights Reserved.
Abbott Medical UK Limited