VINTAGE CPM Modelling Liquid

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MIXING LIQUID

en

INDICATIONS
Mixing liquid for VINTAGE Body and Enamel porcelains to keep moisture longer than distilled water.
Higher consistency than VINTAGE Mixing Liquid for higher shape stability during the build-up.

INGREDIENTS
Distilled water, 1,3-Butanediol

NOTE ON USE
- Use VINTAGE Mixing Liquid-HC instead of distilled water or mix with distilled water to obtain the desired handling characteristics.
- If additional mixing is required to obtain desired handling, add distilled water.
- Do not use for Opaque porcelain.

CAUTION
- Do not use or keep near fire as it is flammable.
- Use in a well-ventilated place.
- Do not expose to high temperature or direct sunlight.
- Keep lid tightly closed if not in use.
- Do not use for purpose other than that described here.
- This product is for use by dental professionals only.

IMPORTANT
- Do not touch with hand. Pay particular attention to avoid contact with skin or eyes.
- In case of contact with skin, wipe off immediately with cotton wool soaked in alcohol and rinse with running water.
- In case of contact with eye, rinse immediately with plenty of water and seek medical advice.
- If skin eruption or inflammation occurs, stop using the product and seek medical advice.

STORAGE
- Keep away from heat, direct sunlight or flame. Seal the cap tightly and keep at room temperature.
- This product must be kept away from all personnel except dental professionals.

PACKAGE
VINTAGE Mixing Liquid-HC 50mL, 500mL

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VERPACKUNG
VINTAGE Mixing Liquid-HC 50mL, 500mL

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CE 0123

Location of the symbols glossary: https://www.shofu.com/en/symbols_glossary/

03.0:2020-09 72480

VINTAGE PRO

- en Instructions for use
- de Gebrauchsanleitung
- fr Mode d'emploi
- nl Gebruiksaanwijzing
- es Instrucciones de uso
- it Istruzioni per l'uso
- ru Инструкция по применению
- pl Instrukcja obsługi
- ro Instrucţiuni de utilizare
- sv Bruksanvisning
- el Οδηγίες χρήσης
- pt Instruções de uso

03.0:2020-09
72486

Manufacturer
SHOFU INC.
11 Kamitakamatsu-cho, Fukuine,
Higashiyama-ku, Kyoto 605-0983, Japan
www.shofu.com

SHOFU DENTAL CORPORATION
1225 Stone Drive, San Marcos,
California 92078-4059, USA

SHOFU DENTAL ASIA-PACIFIC PTE.LTD.
10 Science Park Road, #03-12,
The Alpha, Science Park II, Singapore 117684

EC REP
SHOFU DENTAL GmbH
An der Pönt 70, 40885 Ratingen, Germany

CE 0123

en Metalbond Porcelain

READ CAREFULLY BEFORE USE

PRODUCT DESCRIPTIONS
VINTAGE PRO is a metal bond porcelain system indicated for veneering of precious, semi-precious and non-precious alloys as well as palladium-based alloys with a CTE range of 13.5 to 14.8 x 10^{-6} K^{-1}. VINTAGE PRO is also suitable for metal frameworks made with CAD/CAM and electroforming/sintering. The shade system includes VITA* classical shades with a variety of effect shades to reproduce a wide range of shade.

INDICATIONS FOR USE
Build-up porcelain for fabrication of metal bonded ceramic restorations

PRECAUTIONS AND WARNINGS
1. If any inflammation or other allergic reactions occur on the operator while using this product, immediately discontinue use and seek medical advice.
2. Avoid contact with skin or eyes. In case of accidental contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
3. Use local dust extractor or mask etc, while grinding or polishing this product to avoid harmful influence of dust on the human body.
4. Wear protective glasses etc, while grinding or polishing this product to avoid damage to eyes.
5. Tightly close the bottle cap after each use.
6. Follow the instructions for use of each dental material, instrument or equipment to be used in conjunction with this product.
7. Do not use this product for any purposes other than specifically outlined in these instructions for use.
8. Use this product within the expiration date indicated on the package and container.
 (Example YYYY-MM-DD → Year-Month-Date of the expiration date)
9. This product is intended for use by dental professionals only.

DIRECTIONS FOR USE
Refer to the separate user's manual or detailed instructions for use.
- Use VINTAGE PRO in conjunction with alloys having coefficient of thermal expansion of 13.5 to 14.8 × $10^{-6}K^{-1}$.
- The firing conditions may vary due to the different designs and operating voltages of the porcelain furnaces. It is essential to carry out test firings before using the porcelain for actual restorations.
- Apply the first Opaque layer (BASE) in a thin layer and fire.
- The second layer of Opaque must be applied evenly in a sufficient thickness to mask the color of the underlying metal coping.
- Mix Opaque pastes sufficiently well with a plastic spatula before use.
- Opaque porcelains must be mixed with dedicated VINTAGE PRO opaque mixing liquids only. Do not mix with other VINTAGE mixing liquids.
- When Opaque porcelain is applied generously, set longer drying time at the furnace entrance and preheating time in order to burn the liquid components thoroughly.
- Select the powders based on the combination chart and build up the restoration.
- The Body build-up should be done to the final anatomical shape of the tooth in order to ensure the exact dimension for the Enamel porcelain.
- When large amount of porcelain is built up, set longer drying time at the furnace entrance in order to avoid explosion of the porcelain from the inside resulting from remaining liquid components.
- Enamel porcelain should be built up slightly oversized in consideration of firing shrinkage.
- When the area requiring an additional build-up is dry, wet the area with water by using a brush.
- Use distilled water or VINTAGE MARGIN HARDENING LIQUID for Margin porcelains.
- Use distilled water or VINTAGE CPM MODELLING LIQUID for ADD-ON porcelains; ADD-ON B and ADD-ON T. Do not use VINTAGE MIXING LIQUID or MIXING LIQUID HC for these porcelains since they affect the shade of the restoration.

STORAGE
Store at room temperature (1-30 °C/34-86 °F). Avoid direct sunlight.

CAUTION: US Federal law restricts this device to sale by or on the order of a dental professional.

Location of the symbols glossary: https://www.shofu.com/en/symbols_glossary/

* VITA is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co.KG, Bad Säckingen, Germany

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en READ CAREFULLY BEFORE USE

Mixing Liquid for VINTAGE PRO Paste Opaque

PRODUCT DESCRIPTIONS
VINTAGE PRO PASTE OPAQUE LIQUID is dedicated for VINTAGE PRO Paste Opaque.

INDICATIONS FOR USE
Mixing liquid for dilution of VINTAGE PRO Paste Opaque to obtain optimal creamy consistency

PRECAUTIONS AND WARNINGS
1. Do not handle the liquid or the mixture of this material with bare hands. This product should not come in contact with eyes or skin. In case of accidental contact with skin, immediately blot with alcohol moistened cotton ball. In case of contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
2. If any inflammation or other allergic reactions occur on the operator while using this product, immediately discontinue use and seek medical advice.
3. Avoid any source of ignition since this product is flammable.
4. Use this product in well ventilated area (several appropriate ventilations per hour are recommended).
5. Do not leave this material in areas of high temperature, such as near the heater etc. Avoid direct sunlight.
6. Tightly close the cap after each dispensing.
7. Do not use VINTAGE PRO PASTE OPAQUE LIQUID in conjunction with any materials other than VINTAGE PRO Paste Opaque.
8. Do not use this product for any purposes other than specifically outlined in these instructions for use.
9. This product is intended for use by dental professionals only.

DIRECTIONS FOR USE
1. **Wash Opaque**
 Mix VINTAGE PRO Paste Opaque with VINTAGE PRO PASTE OPAQUE LIQUID until it reaches thin consistency. Apply the mixture to the clean, dry substructure with a brush.
2. **Opaque**
 If required, mix VINTAGE PRO Paste Opaque with VINTAGE PRO PASTE OPAQUE LIQUID until it reaches creamy consistency. Apply the mixture with a brush or instrument to completely mask the surface to be veneered.

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en READ CAREFULLY BEFORE USE

Mixing Liquid for VINTAGE PRO Powder Opaque

PRODUCT DESCRIPTIONS
VINTAGE PRO POWDER OPAQUE LIQUID is dedicated for VINTAGE PRO Powder Opaque. VINTAGE PRO POWDER OPAQUE LIQUID-L (L-type liquid), provides longer working time than standard liquid (VINTAGE PRO POWDER OPAQUE LIQUID) by keeping the mixture optimally consistent.

INDICATIONS FOR USE
Mixing liquid for VINTAGE PRO Powder Opaque for obtaining creamy consistency

PRECAUTIONS AND WARNINGS
1. Do not handle the liquid or the mixture of this material with bare hands. This product should not come in contact with eyes or skin. In case of accidental contact with skin, immediately blot with alcohol moistened cotton ball. In case of contact with eyes, immediately flush the eyes with plenty of water and seek medical advice.
2. If any inflammation or other allergic reactions occur on the operator while using this product, immediately discontinue use and seek medical advice.
3. Avoid any source of ignition since this product is flammable.
4. Use this product in well ventilated areas (several appropriate ventilations per hour are recommended).
5. Do not leave these materials in areas of high temperature, such as near the heater etc. Avoid direct sunlight.
6. Tightly close the cap after each dispensing.
7. Do not use VINTAGE PRO POWDER OPAQUE LIQUID in conjunction with any materials other than VINTAGE PRO Powder Opaque.
8. Do not use this product for any purposes other than specifically outlined in these instructions for use.
9. This product is intended for use by dental professionals only.

DIRECTIONS FOR USE
1. **Wash Opaque**
 Mix VINTAGE PRO Powder Opaque with powder opaque liquids until it reaches thin consistency. Apply the mixture to the clean, dry substructure with a brush.
2. **Opaque**
 Mix VINTAGE PRO Powder Opaque with powder opaque liquids until it reaches creamy consistency. Apply the mixture with a brush or instrument to completely mask the surface to be veneered.
 ➢Use L-type liquid to extend working time.

COMPOSITION
VINTAGE PRO POWDER OPAQUE LIQUID 1,3-Butanediol
VINTAGE PRO POWDER OPAQUE LIQUID-L Glycerin, Propylene glycol, 1,3-Butanediol

STORAGE
Tightly close the cap and store at room temperature.
Keep away from direct sunlight and any source of ignition.

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