Advisor™ VL Circular Mapping Catheter, Sensor Enabled™

REF

D-AVSE-D10-V1525, D-AVSE-D20-V1525, D-AVSE-DF10-V1525, D-AVSE-DF20-V1525

Instructions for Use

EN: English 1
AR: العربية 4
BG: Български 7
CS: Čeština 11
DA: Dansk 14
DE: Deutsch 17
EL: Ελληνικά 21
ES: Español 25
ET: Eesti keel 28
FI: Suomi 31
FR : Français 34
HR: Hrvatski 37
HU: Magyar 40
IT: Italiano 43
KO: 한국어 46
LT: Lietuvių 49
LV: Latviski 52
MK: Македонски јазик 55
NL: Nederlands 58
NO: Norsk 62
PL: Polski 65
PT: Português 68
RO: Română 71
RU: Русский 74
SK: Slovenčina 78
SL: Slovenščina 81
SR: Srpski 84
SV: Svenska 87
TR: Türkçe 90
UK: Українська 93
ZF：繁體中文 97
ZH：简体中文 100

Abbott

This warranty is given by Abbott Medical. To make inquiries regarding this warranty, use the contact information for either the manufacturer or the Australian Sponsor.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott Medical warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott Medical will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott Medical to the address below and after Abbott confirms that the goods are defective.

Abbott Medical
Attn: Postmarket Surveillance
One St. Jude Medical Drive
St. Paul, MN 55117 USA
+1 651 756 2000

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott Medical's liability for the breach of the term or guarantee is limited to, at Abbott Medical's option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED.

See the Terms and Conditions of Sale for further information.

Proposition 65, a State of California voter initiative, requires the following notice:

WARNING: This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

™ Indicates a trademark of the Abbott group of companies.
‡ Indicates a third party trademark, which is property of its respective owner.

Pat. http://www.abbott.com/patents

© 2024 Abbott. All Rights Reserved.

EN: English

Advisor™ VL Circular Mapping Catheter, Sensor Enabled™

Instructions for Use

This IFU is for use only by physicians outside of the United States.

Description

The Advisor™ VL Circular Mapping Catheter, Sensor Enabled™ is a steerable, flexible, insulated electrophysiology (EP) catheter constructed of thermoplastic elastomer material and noble metal electrodes. For both bi-directional and uni-directional catheters, the shaft curvature is manipulated by the shaft deflection mechanism located on the handle at the catheter's proximal end. To adjust the curve on the uni-directional catheter, push or pull the thumb control located on the handle. To adjust the curve on the bi-directional catheter, use the actuator to deflect the catheter in either direction. The distal loop is oriented counter-clockwise as viewed from the handle. The loop is deflected by the actuator located on the handle at the catheter's proximal end. To adjust the loop diameter, rotate the loop actuator in either direction. See the label for specific loop and curve configuration. The catheter features a variable loop design that is intended to accommodate vein sizes ranging from 25 to 15 mm in diameter.

Table 1. Compatible Systems

System Device	Connect via	System Software
EnSite™ X EP System	Sensor Enabled™ Diagnostic Catheter Cable (D-ENS-AVSE-6-CBL)	v 1.1.2 or later
EnSite Precision™ System	Sensor Enabled™ Diagnostic Catheter Cable (D-AVSE-CBL22)	v 2.6.1 or later

Intended Purpose

The Advisor™ VL Circular Mapping Catheter, Sensor Enabled™ is intended for electrogram recording and cardiac stimulation during diagnostic electrophysiology studies.

Indications for Use

The Advisor™ VL Circular Mapping Catheter, Sensor Enabled™ is a steerable electrophysiology catheter with integrated sensors. The catheter is used for recording intracardiac signals and cardiac stimulation in patients with a diagnosed or suspected cardiac arrhythmia during diagnostic electrophysiology studies. The catheter can be used to map the atrial regions of the heart.

Patient Target Group

The intended patient population for the Advisor™ VL Circular Mapping Catheter, Sensor Enabled™ is patients who require recording of intracardiac signals and cardiac stimulation during diagnostic electrophysiological studies of atrial regions of the heart and who have vascular anatomy compatible with introducers that have a minimum inner diameter of 8F.

Intended User

The Advisor™ VL Circular Mapping Catheter, Sensor Enabled™ is intended for use by, or under the supervision of, physicians thoroughly trained in the techniques of percutaneous electrophysiology studies.

Clinical Benefits

The Advisor™ VL Circular Mapping Catheter, Sensor Enabled™ provides diagnostic information, measured by successful intracardiac electrogram recording and cardiac stimulation in the heart, to the physician to aid in the treatment of atrial arrhythmias and resolution of associated symptoms.

Contraindications

- The catheter is contraindicated for patients with prosthetic valves, and patients with left atrial thrombus or myxoma, or interatrial baffle or patch via transseptal approach.
- The catheter is contraindicated for the retrograde trans-aortic approach.
- The catheter is contraindicated for use in the ventricles.
- The catheter is contraindicated for use as an ablation catheter.
- The catheter is contraindicated for patients with an active systemic infection as this may increase the risk for cardiac infection.
- The catheter is contraindicated for patients unable to receive heparin or an acceptable alternative to achieve adequate anticoagulation.

Warnings

- This device should be used by, or under the supervision of, physicians thoroughly trained in the techniques of percutaneous electrophysiology studies.
- Persons with a known history of allergies to any of the materials listed below may suffer an allergic reaction to the device. Before use, counsel the patient on the materials contained in the device and discuss a thorough history of allergies. This device contains:
 - Acrylonitrile butadiene styrene (ABS)
 - Epoxy (Henkel/Hysol 2-part Medical Epoxy)
 - Pellethane
 - Platinum Iridium alloy
 - Polycarbonate (Makrolon)
 - Polyether block amide (Pebax)
- This device is intended for one time use only; do not reprocess or reuse it. Note the product "Use by" date on the package. Any attempt to resterilize and reuse this device may compromise its integrity. Adverse effects of using nonsterile components may include, but are not limited to:
 - Local or systemic infection or reaction
 - Mechanical damage
 - Inaccurate functionality
- Cardiac catheterization procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Careful consideration must therefore be given for the use of this catheter in pregnant women.
- The safety and effectiveness of the device has not been established in pregnant women or pre-pubescent children. Careful consideration must therefore be given for the use of the device in pregnant women or pre-pubescent children.
- Catheter entrapment within the heart or blood vessels or with another device is a possible complication of electrophysiology procedures. To unentangle the catheter, fully open the loop and straighten the catheter shaft, then rotate the handle clockwise.
- Vascular perforation or dissection is an inherent risk of any electrode placement. Careful catheter manipulation must be performed in order to avoid device component damage, thromboembolism, cerebrovascular accident, cardiac damage, perforation, pericardial effusion, or tamponade.
- Risks associated with electrical stimulation may include, but are not limited to, the induction of arrhythmias, such as atrial fibrillation (AF), ventricular tachycardia (VT) requiring cardioversion, and ventricular fibrillation (VF).
- Do not use force to advance or withdraw catheter when resistance is encountered.
- Do not immerse the proximal handle or cable connector in fluids; electrical performance could be affected.

Precautions

- Read directions prior to use.
- Federal law (U.S.) restricts this device to sale by or on order of a physician.
- Maintain an activated clotting time (ACT) of greater than 300 seconds at all times during use of the catheter.
- Personnel handling the electrophysiology catheter must wear gloves.
- To maintain optimal patient safety and electrode catheter integrity, do not wipe this catheter with alcohol.
- Excessive bending or kinking of the catheter may cause damage to the catheter.
- Use care to isolate any unused connector pins of the electrogram cable. This will reduce the chances of developing accidental current pathways to the heart.
- Always straighten the catheter and open the loop before insertion or withdrawal.

- Catheter advancement must be performed under fluoroscopic guidance to minimize the risk of cardiac damage, perforation, or tamponade. Compatible navigation and visualization systems may be used in conjunction with fluoroscopy.
- Do not use if the catheter appears damaged, kinked, or if there is difficulty in deflecting the distal section to achieve the desired curve. Do not use if the catheter does not hold its curve or if there is difficulty deflecting the loop. Do not use if the catheter does not hold its loop diameter.
- After use, classify this device, its accessories, and packaging for disposal appropriately (example: biohazard, sharps, non-hazardous waste, and so forth). Disposal must comply with facility procedures and applicable laws and regulations.
- Catheter materials are not compatible with magnetic resonance imaging (MRI).

Potential Adverse Events

- New or worsening of existing arrhythmia, which may include:
 - Atrial fibrillation
 - Ventricular tachycardia requiring cardioversion
 - Ventricular fibrillation
- Cardiac perforation/pericardial effusion/cardiac tamponade
- Bleeding, including access site
- Vascular access complications, which may include:
 - Femoral artery dissection
 - Dissection
 - Arteriovenous fistula
- Embolism
- Cerebrovascular accident (CVA)
- Cardiovascular injury
- Pericarditis

Directions

1. Inspect the package prior to use. Do not use if the package is open, damaged, or expired.
2. Remove the catheter from its packaging.
 - Completely remove the tray from the pouch.
 - To prevent potential damage to the loop, lift the catheter up and out of the tray.
3. Inspect the electrodes and catheter carefully for integrity and overall condition.
4. Insert the distal tip section of catheter into an 8F minimum introducer (not included):
 - Prior to insertion, deflect the catheter shaft to the straight position and open the loop.
 - Insert the catheter through the hemostasis valve.
5. Never manipulate the loop or deflectable section of the shaft while within the introducer.
6. Connect to compatible systems using the appropriate cable. Refer to the cable's instructions for use.
7. If the loop is not perpendicular to the shaft when extended from the sheath, completely retract the loop into the sheath and re-extend the catheter.
8. The catheter should be passed from a peripheral vessel to the desired endocardiac position with the aid of fluoroscopy. Compatible navigation and visualization systems may be used in conjunction with fluoroscopy.
9. To adjust the curve of the distal tip on the uni-directional catheter, push or pull the thumb control located on the handle. To adjust the curve of the distal tip on the bi-directional catheter, use the actuator to deflect the catheter in either direction. To adjust the loop diameter, rotate the loop actuator.

 NOTE: The catheter handle has an adjustable tension lock knob that allows the operator to use the actuator and deflectable section in an unlocked state or adjust the tension to where the actuator and deflectable section are locked in place. If necessary, the tension lock knob may be rotated to increase or decrease the tension.

10. Prior to withdrawal, deflect the catheter shaft to the straight position and open the loop completely.
11. After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

Connection to Other Equipment

The catheter may be connected to a commercially available EP recording system and navigation and visualization system using the connection cable. All systems must be patient isolated. For instructions regarding the use of these systems with the catheter, refer to the system's instructions for use.

Packaging and Shelf-Life

The catheter packaging is designed to prevent crushing of the product, to minimize product exposure to the atmosphere, and to provide for aseptic product transfer. It is recommended that the product remain in the unopened package until time of use. Contents are sterile if the package is unopened and undamaged. Do not resterilize. The expiration date is marked on the outside of the package. The product should be stored in a cool, dry location. The instructions for use is recyclable. Dispose of used product and packaging following standard solid biohazard waste procedures.

Incident Reporting

If, in the course of use of this device, you have reason to believe that a serious incident occurred, please report it to the manufacturer. For customers in the European Union, report the serious incident to your national authority as well as to the manufacturer.

Clinical Data Summary

A summary of the safety and clinical performance for this device is available at https://ec.europa.eu/tools/eudamed. Search for the device using the UDI-DI provided on the outer packaging of the device.








Warranty

Abbott Medical warrants that its products shall be free from defects in materials and workmanship under normal use. This warranty does not exceed the "Expiration" date stated on any product labeling. The authorized uses and approved methods of use of each of our products is set forth in the related "Instructions for Use" that accompany each product. Abbott Medical disclaims any responsibility and liability for the use of its products in a manner that has not been authorized or approved. Abbott Medical's liability under this warranty is limited to replacing its products. The foregoing warranty excludes and is in lieu of all other warranties whether expressed or implied including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Abbott Medical disclaims any liability for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this product, other than as expressly provided by specific law. Abbott Medical neither assumes nor authorizes any other person to assume for it any other or additional liability for loss, damage, or expense in connection with this product. For more details please review complete Abbott Medical warranty policy available from Abbott Medical or on the back of an Abbott Medical invoice.

Symbols

The symbols below and harmonized symbols may be found on the product or product label. For harmonized symbols, refer to the Universal Symbols Glossary at https://medical.abbott/manuals.

Symbol	Definition
MD	Medical Device
UDI	Unique Device Identification
medical.abbott/manuals	Follow instructions for use on this website
REF	Catalog number
LOT	Lot number
	Use-by date

Symbol	Definition
	Keep dry
	Keep away from sunlight
	Do not reuse
	Do not resterilize
	Do not use if package is damaged
	Sterilized using ethylene oxide
	Date of Manufacture
	Country of manufacture
	Manufacturer
	Outer diameter
	Excursion humidity limit
	Excursion temperature limit
	Authorized European representative
	Quantity
	Single sterile barrier system
	CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician.
	European conformity, affixed in accordance with the relevant provisions of European Council Regulation 2017/745 (NB 2797). Hereby, Abbott Medical declares that this device is in compliance with the relevant provisions of this regulation.
	Importer
	United Kingdom Conformity Assessment
	Responsible person in the United Kingdom
	Ukrainian mark of conformity to the technical regulations, where 101 is the assessment body.
	Number of electrodes
	Electrodes
	Spacing
	Length
	Recommended cables
	Diagnostic Catheter
	Circular Mapping Catheter
	Non-pyrogenic

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Abbott Medical UK Limited