PPL2134963 (2021-11-16)

PressureWire X Guidewire

Wireless Connection

REF

C12059, C12359

CE
2797

Instructions for Use
تعليمات الاستخدام
Указания за употреба
Návod k použití
Brugervejledning
Gebrauchsanweisung
Οδηγίες χρήσης
Instrucciones de uso
Kasutusjuhend
Käyttöohjeet
Mode d'emploi
Upute za uporabu
Használati utasítás
Istruzioni per l'uso
사용 지침
Naudojimo instrukcijos
Lietošanas pamācība
Упатство за употреба
Gebruiksaanwijzing
Bruksanvisning
Instrukcja obsługi
Instruções de Utilização
Instrucțiuni de utilizare
Инструкция по эксплуатации
Návod na použitie
Uputstvo za upotrebu
Bruksanvisning
Kullanma Talimatları
Інструкція з експлуатації
Hướng dẫn sử dụng
使用說明
使用说明

Abbott

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
PressureWire X Guidewire
Wireless Connection

REF C12059, C12359

Instructions for Use

Content

PressureWire X guidewire package includes:	**PressureWire X guidewire pouch includes:**
Sealed non-sterile outer box	Packaging coil
Instructions for use	PressureWire X guidewire
Sealed sterile pouches	Torque device
	Transmitter
	Quick Guide

1. Guidewire, refer to label for length[1]
2. Flexible length, 31 cm
3. Shaft
4. Connector
5. Radiopaque tip, 3 cm
6. Sensor element
7. Hydrophilic coating
8. PTFE coating
9. Hydrophobic insulators between contact members
10. Torque device
11. PressureWire X transmitter
12. Function indicators
13. Switch
14. Close cap
15. Open cap
16. Cap

[1]Applied part according to IEC 60601-1.

Description

The PressureWire X guidewire is a 0.014" (0.36 mm) guidewire with an integrated sensor element at the tip to enable measurements of physiological parameters. The guidewire is uniquely paired with a specific transmitter. The guidewire is available in different lengths.

Intended Use

The PressureWire guidewire is designed to fit inside a percutaneous catheter for the purpose of directing the catheter through a vessel. The signal output from the sensor is transmitted to associated equipment for analysis, calculations, and display of any physiological parameters or indices based on pressure or temperature; e.g., Fractional Flow Reserve (FFR).

Indications for Use

The PressureWire X guidewire is indicated to direct a catheter through a blood vessel and to measure physiological parameters in the heart and in the coronary and peripheral blood vessels. Physiological parameters include blood pressure. The PressureWire X guidewire can also measure blood temperature.

Contraindications

This guidewire is contraindicated for use in the cerebral vasculature.

Warnings

- No modification of this device is allowed.
- The PressureWire X guidewire is supplied sterile. Discard the guidewire if the pouch is opened or damaged, compromising the sterile barrier.
 The guidewire is designed for single use only and shall not be reused or resterilized. Adverse effects of using a non-sterile or resterilized guidewire may include, but are not limited to:
 - Local and / or systemic infection
 - Mechanical damage
 - Inaccurate readings
- Observe all guidewire movements. Whenever the guidewire is moved or torqued, the tip movement should be examined under fluoroscopy. Never push, withdraw, or torque the guidewire if it meets resistance or without observing corresponding movement of the tip; otherwise, vessel / ventricle trauma may occur.
- Torquing or excessive manipulation of the guidewire in a sharp bend, against resistance, or repeated attempts to cross a total vessel occlusion may:
 - Cause dissection or perforation of blood vessels
 - Cause vessel spasm
 - Damage and / or fracture the guidewire
- When introducing the guidewire, flush the catheter and administer anticoagulation as for a standard catheterization procedure or clotting may occur.
- Do not use the guidewire in the ventricles if the patient has a prosthetic mechanical or biological valve. It may result in damage to both the prosthesis and the guidewire, which may cause injury or death.
- Use of the PressureWire X guidewire in conjunction with interventional devices with a short rapid exchange may result in a folded or fractured guidewire.
- High frequency surgical devices must not be used on a patient at the same time as the guidewire.

Precautions

- The PressureWire X guidewire is a delicate instrument and should be handled carefully.
- Make sure that the transmitter is kept dry to ensure accurate pressure and / or temperature readings. Inaccurate readings may necessitate device replacement.
- Do not use the guidewire in conjunction with atherectomy catheters. It may damage the guidewire.
- Do not withdraw or manipulate the guidewire in a sharp-edged object. It may result in abrasion of the guidewire coating.
- Factors that may affect the accuracy of the diagnostic information include, but are not limited to:
 - Improper placement of the aortic pressure sensor.
 - Failure to achieve maximum coronary and myocardial hyperemia during FFR procedures.
 - Blood flow affected by the position of interventional devices, such as balloon catheters.
- Guidewire readings may be affected by defibrillation. Rezero the guidewire after defibrillation use.
- Do not measure pressure when the guidewire sensor element is in a sharp bend or in contact with atrial or ventricular walls. It might result in pressure artifacts.
- Do not use the PressureWire X guidewire together with another guidewire, for so-called jailed wire technique, due to difficulty in guidewire withdrawal and possible guidewire entrapment.
- Store at room temperature (15°C – 25°C) in a dry and dark place.

Adverse Events

Potential complications which may be encountered during all catheterization procedures include, but are not limited to: vessel dissection or occlusion, perforation, embolus, spasm, local and / or systemic infection, pneumothorax, congestive heart failure, myocardial infarction, hypotension, chest pain, renal insufficiency, serious arrhythmias, or death.

In addition, this device has a coating containing Polyethylene Glycol (PEG); potential allergic reactions (anaphylaxis) may occur during the interventional procedure if the patient is allergic to PEG.

Directions for Use

Use the PressureWire X guidewire in conjunction with a 6F (2 mm diameter) guiding catheter. Refer to instructions supplied with any interventional devices to be used in conjunction with the guidewire for their intended uses, contraindications, and potential complications. Note that use of the guidewire in conjunction with uncoated diagnostic catheters may reduce guidewire performance.

CAUTION: If the PressureWire X guidewire gets stuck or is damaged during the procedure, immediately disconnect the transmitter and replace the damaged guidewire, including the transmitter, with a new one. Damage may include, but is not limited to, kinks, bends, sheared polymer coating, no pressure signal, or inaccurate pressure signal. Incorrect use may result in vessel / ventricle damage, arrhythmia induction, electric shock, inaccurate pressure signals or inaccurate torque response.

CAUTION: When the PressureWire X guidewire is disconnected from the transmitter during procedure, make sure the male connector does not come into contact with conductive surfaces, thus avoiding unintentional connection with other equipment or electric shock.

Preparations

1. Open the PressureWire X guidewire packaging using aseptic technique.
2. Make sure that the packaging coil is flat on the table.
3. Fill the packaging coil with saline solution through the flush port. The saline shall cover the pressure sensor element, approximately 20 cc.

Zeroing the Guidewire

4. Activate the receiving instrument according to manufacturer's instruction.
5. Zero the PressureWire X guidewire by turning on the transmitter.
 CAUTION: Keep the guidewire stable at all times and avoid any excessive movement during zeroing.
 The transmitter's function indicator shall display a stable green light when it is successfully zeroed. Refer to Table 5 for further information.
6. Carefully remove the guidewire from the packaging coil.
7. The shapeable guidewire tip may be carefully shaped using standard tip shaping practices.
 CAUTION: Do not use a shaping tool with a sharp edge.

Diagnostic Procedure

8. Introduce the PressureWire X guidewire into the patient's blood vessels according to standard clinical practice.
 CAUTION: Avoid turning off the transmitter during the procedure. The PressureWire X guidewire is automatically zeroed when the transmitter is turned on and must be zeroed outside of the body.
9. Use the torque device to slowly steer the guidewire. Position the guidewire sensor element just outside the guiding catheter opening.
10. Verify the following to facilitate an accurate pressure reading:

- Place the aortic pressure transducer at the same height as the patient's heart.
- Pull back the insertion tool and tighten the hemostatic valve.
- Make sure that there is no remaining contrast fluid in the catheter.

11. **Verify pressures.** Verify that the pressures from the guiding catheter and the guidewire are equal.

- If the pressure difference is $\leq$ 10 mmHg, equalize the pressure signals according to compatible system's instructions. Refer to Table 6 for a list of compatible systems.
- If the pressure difference is $>$ 10 mmHg, it is recommended to:
 - Confirm the items noted in step 10.
 - If the difference is still $>$ 10 mmHg, remove the PressureWire X guidewire, re-zero, and repeat from step 8.
 - Equalize the pressure signals according to the compatible system instructions.
- If the pressure difference results in an equalization offset value, limit, or alert according to the compatible systems, refer to that equipment's instructions.

 NOTE: If equipment not identified in Table 6 is used to display pressure from the guidewire, refer to related equipment specification if additional limitations apply.

12. **Perform measurements.** Use the torque device to slowly steer the guidewire to the desired position and perform pressure measurement according to standard clinical procedures.

 NOTE: If it is difficult to reach the desired position, the guidewire may be disconnected from the transmitter for better handling.
13. If applicable, considering the results from the measurements, proceed with interventional procedures.

Interventional Procedure

14. Disconnect the transmitter.
15. Remove the torque device by loosening the screw and gently withdraw it.
16. Carefully advance the interventional device over the PressureWire X guidewire and proceed according to manufacturer's instructions.
17. Put the torque device back on the guidewire in a suitable position.
18. If it is necessary to re-connect the transmitter, carefully wipe and dry the guidewire.
19. Gently and fully insert the connector into the transmitter and secure the cap.
20. Proceed with treatment.

Removal / Re-introduction of the PressureWire X Guidewire

If it is necessary to re-introduce the PressureWire X guidewire to the patient's blood vessels during any diagnostic or interventional procedure:

1. Carefully withdraw the guidewire. Make sure that the distal flexible part of the guidewire is kept wet.
2. Carefully wipe and dry the connector.
3. Gently and fully insert the connector into the transmitter and secure the cap.
4. Repeat the procedures from step 8 and forward.

Verify Pressure Measurement

1. Carefully withdraw the PressureWire X guidewire until the sensor element is positioned just outside the catheter opening.
2. Pull back the insertion tool and tighten the hemostatic valve. Make sure that there is no remaining contrast fluid in the catheter.
3. Verify that the pressures registered by the guiding catheter and the guidewire are equal.

Finish Procedure

Carefully withdraw the PressureWire X guidewire. Turn off the transmitter.

Handling and Disposal

Handle and dispose of the PressureWire X guidewire as biohazardous material and in accordance with medical and applicable local, state and federal law and regulations. The instructions for use is recyclable.

The battery can be removed by opening the end cap using forceps, for example, in the slits under the label.

CAUTION: The transmitter is powered by a standard Lithium-ion battery. Lithium-ion batteries can contain a small concentration of 1.2-Dimethoxyethane. Do not open the battery or disassemble it. Do not expose to fire or high temperatures (> 60°C). Dispose battery according to standard local practice for Lithium-ion batteries.

Technical Specifications

Table 1. Pressure measurement

Operating pressure	-30 to +300 mmHg
Accuracy[1]	≤ 2 mmHg or ± 2%
Temperature sensitivity	0.3% of pressure change per °C
Zero drift	< 7 mmHg/hour
Operating time	≤ 3 hours

[1] The contribution to the error in average pressure from any external disturbance source will be ≤ 2 mmHg; i.e., when a disturbance source is applied in accordance with IEC 60601-1-2, the average pressure may change up to 2 mmHg. The effect of electromagnetic noise on pressure signal is less than 3 mmHg peak to peak.

Table 2. Temperature performance specification

Temperature range	34 – 42°C

Table 3. Radio signal specification

Frequency range	2.4000 – 2.4835 GHz (ISM-band)
Type	Frequency Hopping Spread Spectrum
Radiated power	1 mW peak, 70 µW average (EIRP)
Range	Up to 5 m. **NOTE**: Radio range is reduced by objects and walls, keep transmitter and receiver in line of sight wherever possible.
Delay time	< 28 ms

Table 4. Environmental conditions

Operating, ambient temperature	15°C to 35°C (59°F to 95°F) with a maximum temperature fluctuation of ± 5°C
Operating, relative air humidity	30 – 75%
Ambient pressure	425 – 850 mmHg
Transport temperature	-20°C to +50°C (-4°F to +122°F)
Transport relative humidity	25 – 90%
Storage conditions	15°C to 25°C (59°F to 77°F)

Table 5. Transmitter function indicators light

Function indicators	Explanation
Fast blinking green light	Zeroing (initializing)
Green steady light	Transmitter on, normal operation
Alternating green / yellow blinking light	Stability check failed (pulsatile pressure detected). Rezero the guidewire at stable zero pressure.
Yellow steady light	The guidewire is disconnected from transmitter. Reinsert the guidewire firmly into transmitter.
Yellow blinking light	Error. Possible causes can be that the battery level is too low for operation or other internal error. Exchange the guidewire including the transmitter.

Table 6. Compatible systems

PressureWire Receiver
ILUMIEN™ PCI Optimization System
ILUMIEN™ OPTIS™ System
OPTIS™ Integrated System
OPTIS™ Mobile System
QUANTIEN™ Measurement System
For other compatible platforms, check with your sales representative.

Compliance with Regulatory Requirements

Guidance and Manufacturer's Declaration - Electromagnetic Emissions

The PressureWire X guidewire is intended for use in the electromagnetic environment specified below. The customer or the user of the PressureWire system should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The PressureWire system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	N/A No mains connection	The PressureWire X guidewire with wireless connection is suitable for use in healthcare establishments. There may be potential difficulties ensuring electromagnetic compatibility in domestic building environments due to conducted as well as radiated disturbances.
Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A No mains connection	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

Table 7. Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The PressureWire X guidewire is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity Test	**IEC 60601 Test Level**	**Compliance Level**	**Electromagnetic Environment – Guidance**
Electrostatic discharge (ESD) EN 61000-4-2 (IEC 61000-4-2)	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material the relative humidity should be at least 30%.
Electrical fast transient / burst EN 61000-4-4 (IEC 61000-4-4)	±2 kV for power supply lines ±1 kV for input / output lines	±2 kV ±1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5 (IEC 61000-4-5)	±1 kV differential mode ±2 kV common mode	±1 kV ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variation on power supply input lines. IEC 61000-4-11	0% U_T; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% U_T; 1 cycle and 70% U_T; 25/30 cycles Single phase: at 0° 0% U_T for 5 sec @ 60 Hz (300 cycles) 0% U_T; 250/300 cycles	100% dropout in VNOM for 0.5 cycle at listed phase angles 100% dropout in VNOM for 1 cycle at 0° 30% dropout in VNOM for 25/30 cycles at 0° 100% dropout in VNOM for 5 sec 100% interrupt in VNOM for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_T is the AC mains voltage prior to application of the test level.

Table 8. Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The PressureWire X guidewire is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity Test	**IEC 60601 Test Level**	**Compliance Level**	**Electromagnetic Environment – Guidance**
			Portable and mobile RF communications equipment should be used no closer to any part of the device including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms [V1 = 3]	Recommended separation distance $d = [\frac{3.5}{V_1}]\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m [E1 = 3]	$d = [1.2]\sqrt{P}$ $d = [\frac{3.5}{E_1}]\sqrt{P}$ 80 MHz to 800 MHz $d = [1.2]\sqrt{P}$ $d = [\frac{7}{E_1}]\sqrt{P}$ 800 MHz to 2.7 GHz $d = [2.3]\sqrt{P}$

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
			where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and *d* is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,[a] should be less than the compliance level in each frequency range.[b] Interference may occur in the vicinity of equipment marked with the following symbol:
Immunity to Proximity Fields from RF wireless communications equipment IEC 60601-1-2 (Clause 8.10)	385 – 5785 MHz	9 – 28 V/m	Per IEC 60601-1-2 (Per Table 9 of Standard)

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

[a] Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

[b] Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Recommended separation distances between portable and mobile RF communications equipment and PressureWire system

The PressureWire X guidewire with wireless connection is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PressureWire X guidewire with wireless connection can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PressureWire X guidewire with wireless connection as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance to frequency of transmitter (m)		
	150 kHz to 80 MHz*	80 MHz to 800 MHz	800 MHz to 2.7 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.38	0.74
1	N/A	1.2	2.3
10	N/A	3.8	7.4
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

*Not applicable since the transmitter is not connected to mains.

NOTE: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.


NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Warranty Disclaimer

Although PressureWire X guidewire (hereafter referred to as "product") has been manufactured under carefully controlled conditions, Abbott Medical (hereafter referred to as "Abbott") has no control over the conditions under which the product is used. Abbott, therefore, disclaims all warranties, both expressed and implied, with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Abbott shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Abbott to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

Symbols

Symbol	Description
Guidewire	Guidewire
Wireless Connection	Wireless connection
	Consult instructions for use
	Caution
	Use by
	Manufacturer
	Manufacturing facility
	Date of manufacture
REF	Catalogue number
	Length
Ø	Outer diameter
LOT	Lot number
	Packaging unit
	Keep away from sunlight
	Keep dry
	Temperature limitations
EXCURSIONS	Excursion limits permitted to temperature
EXCURSIONS	Excursion limits permitted to humidity
	Defibrillation-proof Type CF equipment
	Equipment includes RF transmitter
	Do not use if package is damaged
	Do not reuse
STERILIZE	Do not resterilize

Symbol	Description
STERILE EO	Sterilized with Ethylene Oxide Gas
	Dispose of hardware in accordance with local law
CE 2797	European conformity, affixed according to the relevant provisions of MD directives 93/42/EEC and 2011/65/EU, and RE directive 2014/53/EU Annex II. Hereby, Abbott Medical declares that this device complies with the essential requirements and other relevant provisions of these directives. The full text of the European Union RE directive 2014/53/EU declaration of conformity is available at the following internet address: www.cardiovascular.abbott/euconformity.
EC REP	Authorized representative in the European Community
℞ ONLY	CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a physician.
FC	FCC identifier for the transmitter This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
IC:	IC identifier for transmitter This device complies with RSS-210 of the IC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation
RECOGNIZED COMPONENT c ETL us Intertek 5013704	cETLus valid for Canada and US. Conforms to AAMI Std. ES60601-1 & IEC Std. 60601-2-34 Certified to CSA Std. C22.2 No. 60601-1
Medical Electrical Equipment	Medical electrical equipment
	Australian Communications and Media Authority (ACMA) and New Zealand Radio Spectrum Management (RSM) Regulatory Compliance Mark (RCM)
R-NZ	New Zealand Radio Spectrum Management (RSM) compliance mark for products that have radio transmitters of conformity level A1, A2, or A3
	Japanese RF identifier for the transmitter
	KCC identifier for the receiver
Canada 310	This device complies with RSS-310 of Industry Canada. Operation is subject to the condition that this device does not cause harmful interference.
	Flush packaging coil with saline
	Connect
	Switch on transmitter

Australian Warranty. This section is provided in English only.

This warranty is given by Abbott. To make inquiries regarding this warranty, use the contact information for either the manufacturer or Australian Sponsor located on the back cover.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

The benefits to you given by this warranty are in addition to your other rights and remedies under the Australian Consumer Law.

Abbott warrants to the buyer that the goods, for a period equal to the validated shelf life of the goods (the "Warranty Period"), shall meet the product specifications established by the manufacturer when used in accordance with the manufacturer's instructions for use and shall be free from defects in materials and workmanship. Unless otherwise obligated by law or mandatory consumer guarantees, under this warranty Abbott will, at its option, replace or repair the goods at its factory, if the buyer, at the buyer's expense, returns the goods within the Warranty Period to Abbott to the address below and after Abbott confirms that the goods are defective.

Abbott Medical
Attn: Postmarket Surveillance
Return Number
Roseville Distribution Center
2305 Walnut Street
Roseville, MN 55113 USA
+1 651 756 2000

To the maximum extent permitted by law, if a mandatory term is implied by law or a mandatory consumer guarantee applies to the goods, and the goods are not of a kind ordinarily acquired for personal, domestic or household use or consumption, Abbott's liability for the breach of the term or guarantee is limited to, at Abbott's option, either replacement or repair of the goods or payment of the costs of replacing or repairing the goods.

EXCEPT FOR MANDATORY TERMS IMPLIED BY LAW, MANDATORY CONSUMER GUARANTEES, OR AS EXPRESSLY PROVIDED IN THIS WARRANTY, ABBOTT MEDICAL DISCLAIMS ANY REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED.

See the Terms and Conditions of Sale for further information.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNING: This product can expose you to chemicals including ethylene oxide, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

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© 2022 Abbott. All Rights Reserved.