Representation of the expected values on the Quo-Lab A1C Control Data Card

LEVEL 1			
Unit	**LOWER LIMIT**	**TARGET**	**UPPER LIMIT**
DCCT (%)			
IFCC (mmol/mol)			
JDS (%)			

Level 1 = normal

LEVEL 2			
Unit	**LOWER LIMIT**	**TARGET**	**UPPER LIMIT**
DCCT (%)			
IFCC (mmol/mol)			
JDS (%)			

Level 2 = abnormal

Calibration

Quality Control values are assigned using equipment calibrated using samples provided by the European Reference Laboratory, via the NGSP network. They are traceable to the IFCC reference method.

Troubleshooting

If the result of either control fails then review the Instructions For Use and repeat the control again. If the repeat test also fails then contact your local distributor or call EKF-diagnostic GmbH customer support on +49 39203 511 414 or support@ekf-diagnostic.de.

You must not use the analyzer or report any test results until the problem has been resolved.

Storage and Stability

Store unopened control kits in a refrigerator at 2 to 8°C (36 to 46°F). The expiry date is printed on the outside of the pack. The controls cannot be used after this date.

Once controls have been reconstituted they can be used for up to one month, when kept in their package in a refrigerator at 2 to 8°C (36 to 46°F) .

Protect the controls from sunlight. Always store the controls and the reconstituted controls within their original box.

Limitations of Use

For use with the Quo-Lab Analyzer and Quo-Lab A1C Test Cartridges only.

WARNING! POTENTIALLY BIOHAZARDOUS MATERIAL

Human sourced materials were used in the manufacture of this product. Each donor unit was tested for hepatitis B surface antigen (HbsAg), antibodies to hepatitis C (HCV), and antibodies to Human Immunodeficiency Viruses (HIV-1 and HIV-2), HIV-1 Antigen and Syphilis and found to be negative. **CAUTION!** Because no test method can offer complete assurance that HIV, Hepatitis B or C viruses, or other infectious agents are absent, these products must be handled as potentially infectious. Handle the Hemoglobin A1c Control with the same precautions used with patient specimens.

EKF-diagnostic GmbH
Ebendorfer Chaussee 3
39179 Barleben
Germany
T+49 39203 511 0
support@ekf-diagnostic.de
www.ekfdiagnostics.com

EKF Diagnostics for life

Quo-Lab A1C Control Kit
REF 3101-6021_VS
Instructions For Use

Read and familiarize yourself with the contents of the Instructions for Use before using the controls.

Consult the Quo-Lab A1C Control Data Card for the expected values.

Please note the modified reconstitution process of the controls.

	MEANING
	Read the instructions for use
IVD	*In-vitro*-Diagnostic Medical Device
LOT	Batch or lot number
CAL	Calibration Barcode
REF	Product Reference
CONTROL	Control
2°C 8°C	Storage Temperature
Σ	Contents sufficient for <n> Number of Tests
	Expiry date
	Manufacturer
	Caution
CE	CE mark
	Biological Hazard
	Keep away from sunlight

3101-9202-1115 EN.01.00 06/2019

Intended Use

The Quo-Lab A1C Control Kit is intended to check the general functionality and accuracy of the Quo-Lab A1C Test system.

Description of the Quo-Lab A1C Control Kit

The kit contains two controls: Level 1 (in the normal value range) and Level 2 (in the abnormal value range) controls. **Both** controls, **Level 1 and Level 2** should be used when Quality Control is being carried out.

Quality controls should be run:

- With each new lot of Test Cartridges.
- With each new shipment of Test Cartridges.
- At any time there is concern that the test result may be wrong, the Test Cartridges may not have been stored correctly, the user is unfamiliar with the test procedure or there is concern that a user may be performing the test incorrectly.

Use the controls according to the quality standards set by your organization or laboratory. Refer to the **Expected Values** on the **Quo-Lab A1C Control Data Card** to check the results obtained.

Only Quo-Lab A1C controls are suitable for use with the Quo-Lab A1C Test System. Preservatives and additives used in third party controls may interfere with the running of the A1C test and EKF-diagnostic GmbH cannot guarantee the results obtained with third party controls.

In addition to running the Quo-Lab A1C controls, it is recommended that laboratories regularly participate in an external quality control scheme.

Quo-Lab A1C Control Kit Contents

2 bottles of Level 1 Control with white lids.

2 bottles of Level 2 Control with black lids.

5 ampoules containing reconstitution solution for reconstituting the controls.

4 color coded dropper caps (2 white, 2 black).

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Quo-Lab Analyzer
Quo-Lab A1C Test Kit

Preparing the Controls

Follow these instructions to reconstitute the controls prior to use:

- Place one Level 1 control bottle and one ampoule on to a clean, dry surface.
- **Use one ampoule containing reconstitution solution for each control bottle.**
- Unscrew the cap of the control bottle. Open the ampoule by carefully breaking the plastic cap. Invert the ampoule over the control bottle. Fill the complete contents of the ampoule to the control bottle by gently squeezing the ampoule.
- When the reconstitution solution has been added, put the cap back on the control bottle. and dispose the ampoule.
- Leave the control bottle to stand for 15 minutes. At the end of this time, gently swirl the contents to ensure that they are completely reconstituted. **Avoid** vigorously shaking the bottle as this will create bubbles.
- Leave the control bottle to stand for a further 15 minutes before use.
- Remove the control bottle cap and replace it with the white dropper cap.
- Write the date prepared on the control label. The control can be used for up to one month when stored at 2 to 8°C (36 to 46°F).
- Follow the above instructions to prepare the Level 2 control using the black dropper cap.
- **The controls are now ready. Refrigerate immediately until required for use. Always store the controls within their original box.**

Using the Controls

Refer to the Quo-Lab User Manual for setting up the analyzer and the Quo-Lab A1C Test Kit Instructions For Use for preparing Test Cartridges for use.

Follow these instructions to run the controls on a Quo-Lab Analyzer:

- Place the Level 1 Control, Quo-Lab A1C Test Cartridges, the Sampling Sticks, a small plastic tray (provided in the control kit) and the Quo-Lab Analyzer on to a clean, dry, level surface. The Quo-Lab Analyzer and Test Cartridges **must** be at the same working temperature before use.
- Reconstituted controls can be used straight from the refrigerator and do not need to be at room temperature before use.
- Again swirl the bottle with the reconstituted control before use.
- Place the barcode scanner over the Level 1 Control Barcode, which is printed on the Quo-Lab A1C Control Data Card. Press the barcode scanner button until you hear a beep.
- Prepare and run an A1C test cartridge as per the Quo-Lab A1C Test Kit Instructions For Use.
- Whilst the analyzer prepares the reagent, unscrew the dropper cap, expel any control solution in the dropper into the bottle, refit the dropper cap and gently swirl the contents before use. **Avoid creating bubbles**.
- Unscrew the dropper cap, gently squeeze the bulb on the control cap, place the tip of the dropper into the control solution in the bottle and gently release the bulb on the control cap. Some of the control solution will be drawn into the glass dropper. **Avoid taking up bubbles.**
- Place a droplet of the control solution on to the small plastic tray by gently squeezing the bulb.
- Screw the dropper cap back on the control bottle.
- When the analyzer prompts the you to 'Insert Sample', touch the pointed end of the Sampling Stick to the drop of control solution on the plastic tray. **Avoid taking up bubbles**. Insert the Sampling Stick in the Test Cartridge and run the control sample, as defined in the Quo-Lab A1C Test Kit Instructions For Use.
- When the test is finished, the result will be displayed on the screen. If the result is within the acceptable range, the screen will show "Control OK". If the result is outside the range then the screen will show "Control Fail". The measured A1C value is also displayed.

The result will remain on screen until the door is opened and the used cartridge removed.

- Now run the Level 2 Control in the same way, making sure that you scan the Level 2 Control Barcode.
- Refer to the **Expected Values** on the Quo-Lab A1C Control Data Card.
- **Return the control material to the refrigerator immediately after use. Always store the controls within their original box.**
- Quality control results are stored in the analyzer memory as control samples. Control sample data files include the test lot number, control lot number, whether the control was run as a Level 1 or Level 2 control, the A1C value and if the control passed or failed. Refer to the Quo-Lab Analyzer User Manual to access results held in the memory.

Warnings and Precautions

For *in vitro* diagnostic use only.

Do not use an unopened control kit after the expiry date or if the control kit has not been stored in accordance with the recommendations.

Do not use controls beyond one month once reconstituted or if the reconstituted controls have not been stored in accordance with the recommendations.

Do not use the control kit if the vials containing the control material or dropper bottle containing the reconstitutrion solution are damaged.

Use gloves - Used Test Cartridges, sampling equipment and controls are potentially infectious.

Used Test Cartridges should be disposed of **immediately** as **clinical waste** in accordance with local / national regulations.

Protect the controls from sunlight by storing within their original box, whether unopened or in reconstituted form.

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EN - Instructions For Use

Test Cartridge Storage

Refrigerated 2 to 8°C (36 to 46°F)
Quo-Lab A1C Test Cartridges are stable until the expiry date only when stored refrigerated. The expiry date is the last day of the month stated on the cartridge, cartridge tray and kit carton.

Room Temperature Storage 15 to 25°C (59 to 77°F)
Quo-Lab A1C Test Cartridges can be stored at room temperature for up to 90 days, if still within the expiry date stated on the cartridge, cartridge tray and kit carton. Write the date of removal from the refrigerator on the cartridge tray or kit carton. Avoid exposure to direct sunlight.

Test Cartridge Preparation and Integrity Check

Allow at least 50 minutes for unopened cartridges to equilibrate to the analyzer's environmental temperature before use, even if previously stored at room temperature. **Do not** use heaters to warm the cartridges.

Check that the foil seal covering the top of the cartridge is not damaged. If the foil seal is damaged then the cartridge **must not** be used.

If problems are experienced, contact your local distributor or contact customer support on +49 39203 511 414 or support@ekf-diagnostic.de

Intended Use

The Quo-Lab A1C Test is intended for the *in vitro* quantitative determination of glycated hemoglobin (HbA1c) in whole blood obtained from finger prick or from venous blood samples collected into EDTA tubes.

The Quo-Lab A1C Test is indicated in the management and treatment of diabetes mellitus and for monitoring long term glycemic control in patients diagnosed with diabetes.

The Quo-Lab A1C Test is designed for professional use only.

Summary and Explanation of the Quo-Lab A1C Test

Human red blood cells (erythrocytes) are freely permeable to glucose present in the surrounding liquid (plasma) of the blood.

During the lifetime of the erythrocytes (normally up to 120 days) exposure to blood glucose results in binding of glucose to the hemoglobin A molecule present in erythrocytes. This is referred to as glycated hemoglobin[1].

Chronic elevated blood glucose will result in tissue and organ damage[2].

Good control of blood glucose, evident as lower HbA1c values, has been proven to result in delayed onset and slower progression of these complications[2,3].

Principle of the Test

The Quo-Lab A1C Test Cartridge contains all of the reagents required for the determination of glycated hemoglobin in whole blood.

The Quo-Lab A1C Test combines the chemical binding of boronate to glycated hemoglobin with the fluorescent quenching effect that this binding exerts on a fluorescent marker bound to the boronate molecule.

After removing the foil seal, the cartridge is placed in the Quo-Lab Analyzer.

When instructed, the user pushes a mixing bead containing the fluorescent boronate conjugate into the buffer contained within the cartridge, using the blunt end of the Sample Stick, and the analyzer automatically rehydrates it.

The user then places the blood sample, contained on the Sample Stick, gently into the top of the cartridge, breaks off and removes the handle and closes the door of the analyzer to introduce the blood sample it into the buffer.

The analyzer automatically mixes the blood sample to release the hemoglobin from the erythrocytes.

The total hemoglobin concentration is determined from the initial decrease in the fluorescent signal.

The fluorescent boronate conjugate binds to the glycated hemoglobin, which is measured by monitoring a decrease in the fluorescence of the active ingredient.

The ratio of glycated hemoglobin to total hemoglobin is determined and the result is presented in up to two user selectable units: % DCCT, mmol/mol IFCC, % JDS, eAG mg/dl or eAG mmol/l.

Linear relationships have been established between these reportable units[4,5]:

mmol/mol IFCC = (% DCCT – 2.15) x 10.929
% JDS = (0.09274 x mmol/mol IFCC) + 1.724
eAG mg/dl = (28.7 x % DCCT) - 46.7
eAG mmol/l = (1.59 x % DCCT) - 2.59

eAG values are based on a correlation study linking % DCCT to the patient's average blood glucose concentration, resulting in published formula to derive the eAG[5].

Cartridge Pack Contents

REF 3101-3051 / REF 0055

25/50 Quo-Lab A1C Test Cartridges (1/2 trays of 25 Cartridges).
1 carton containing a minimum of 25/50 Sample Sticks.
1/2 Calibration Barcodes (one on each tray of 25 Cartridges, indicated by the CAL symbol).
1 Instruction For Use leaflet.

Accessories Required (not included)

Quo-Lab Analyzer
Quo-Lab A1C Control Kit
Standard Blood Collection Equipment

Active Ingredients

Fluorescent boronate conjugate
Plastic vial containing ammonium chloride buffer, lysis agent and preservative
Sample Stick treated with EDTA and surfactant

Calibration of the Quo-Lab A1C Test

The Quo-Lab A1C System is certified by the NGSP[6].

The Quo-Lab Analyzer and A1C Test Cartridges have been calibrated using samples provided by the ERL via the NGSP network.

Results obtained using the Quo-Lab A1C System are traceable to the IFCC reference method[7].

Limitations of Use

For hemoglobin concentrations outside of the range of 6.5 g/dL to 20.4 g/dL (4.0 to 12.7 mmol/l) or for HbA1c values outside the range of 4.0 to 15.0 % DCCT (20 to 140 mmol/mol IFCC), a result will not be reported and an error message will be displayed. Consult the Quo-Lab Analyzer System User Manual for further details.

eAG values may differ significantly from a patient's glucose concentration if checked at the same time.

As with all HbA1c analyses, the results obtained are dependent upon the patient's clinical history and having red blood cells with a typical lifespan.

Warnings and Precautions

For *in vitro* diagnostic use only.

Do not use Test Cartridges after the expiry date or if the Test Cartridges have not been stored in accordance with the recommendations.

Do not use the Test Cartridge if the foil seal coving the top of the cartridge has been damaged.

The Test Cartridge contains sodium azide as a preservative. The concentration is < 0.1 %, which is below the level considered to be hazardous in normal use. In case of leakage, avoid contact with eyes and skin. In case of contact, wash with plenty of water.

Do not re-use any part of the Test Cartridge.

Use Gloves – Used Test Cartridges, sampling equipment, patient samples and controls are potentially infectious. Used Test Cartridges should be disposed of immediately as clinical waste in accordance with local / national regulations.

Operate the Quo-Lab Analyzer System according to the User Manual.

Sample Preparation and Collection

Finger Prick Blood Sample

The patient's finger must be warm, dry and clean, including being free of substances such as hand cream.

Use a single use lancet (not provided) to prick the finger.

Allow a droplet of blood to form, about the width of the Sample Stick, before sampling.

Fill the Sample Stick, avoiding air bubbles and excess sample on the outside of the tip of the Sample Stick.
Do not wipe excess blood from the Sample Stick.

Place the Sample Stick into the Test Cartridge held in the analyzer within 1 minute.

Operate the Quo-Lab Analyzer System according to the User Manual.

Venous Blood Sample

Venous samples must have been collected into EDTA tubes.

Venous samples can be stored refrigerated at 2 to 8 °C (36 to 46 °F) for up to 10 days before use and can be used without equilibration to room temperature.

Thoroughly mix the sample by inverting the tube 8 to 10 times before collecting a sample.

Place a droplet of blood, approximately 10 to 20 µl, onto a non-metallic and non-absorbent surface such as Parafilm® M[9].
Do not attempt to sample directly from the tube.

Fill the Sample Stick, avoiding air bubbles and excess sample on the outside of the tip of the Sample Stick.
Do not wipe excess blood from the Sample Stick.

Place the Sample Stick into the Test Cartridge held in the analyzer within 1 minute.

Operate the Quo-Lab Analyzer System according to the User Manual.

Expected Values

Interpret test results with careful consideration to the patient's age, ethnicity, medical history, clinical examinations and other laboratory results, see limitations of use.

Users should refer to national guidance regarding patient reference ranges, an example is as follows[8]:

The ADA recommends an HbA1c range of 4 to 5.7 % DCCT (20 to 39 mmol/mol IFCC) for non-diabetic patients.

Patients with values in the range of 5.7 to 6.0 % DCCT (39 to 42 mmol/mol IFCC) are considered to be at risk of developing diabetes and those in the range of 6.0 to 6.5 % DCCT (42 to 48 mmol/mol IFCC) are considered to be at high risk of developing diabetes.

The ADA recommends a goal of < 7 % DCCT (53 mmol/mol IFCC) for effective management of diabetes and to minimize long term complications.

Quality Control

The Quo-Lab A1C Control Kit must be used with Quo-Lab A1C Test Cartridges to ensure that the system is operating correctly.

Only Quo-Lab A1C Controls should be used with the Quo-Lab A1C Test.

Please read the IFU included with each Quo-Lab A1C Control Kit to ensure that the correct procedure is followed.

Quality control should be run:

- With each new lot of Test Cartridges.
- With each new shipment of Test Cartridges.
- At any time there is concern that the result may be wrong, the Test Cartridges may not have been stored correctly, the user is unfamiliar with the test procedure or there is a concern that a user may be performing the test incorrectly.

Use the Quo-Lab A1C Control Kit in accordance with the quality standards set by your local organization or laboratory.

If the Quo-Lab A1C controls do not give the results within the range stated in the control kit IFU, contact your local distributor or contact customer support on +49 39203 511 414 or support@ekf-diagnostic.de

Performance Characteristics

Method Comparison

Quo-Lab A1C results were compared with a laboratory HPLC system giving a linear regression of r = 0.989 (y = 0.978x + 0.278, n = 46).

Precision

Blood samples with high (12.8 %) and low (5.3 %) A1C values were tested in duplicate every morning and afternoon for twenty days on a single Quo-Lab Analyzer. The total imprecision for each level was < 3 % CV respectively.

Linearity and Test Range

High and low pooled blood samples were used to create eleven linearity standards, which were tested on the Quo-Lab A1C Test to establish the linear range. The Quo-Lab A1C Test has a linear working range of 4 to 15 % DCCT (20 to 140 mmol/mol IFCC).

If a result is outside of this range then the analyzer will display an error message. Refer to the Quo-Lab Analyzer System User Manual for further details.

Interfering Substances

No affect from bilirubin (20 mg/dl), cholesterol (500 mg/dl), creatinine (30 mg/dl), triglycerides (1,500mg/dl), uric acid (20 mg/dl), acetaminophen (20 mg/dl), ascorbic acid (3 mg/dl), dopamine (13 mg/dl), ibuprofen (40 mg/dl), salicylate (50 mg/dl), tetracycline (4mg/dl), tolazamide (100mg/dl) and tolbutamide (100 mg/dl) was found.

Hemoglobin Variants

The Quo-Lab A1C Test was found to be unaffected by labile glycated hemoglobin and the following hemoglobin variants: Hb AS, Hb AC, Hb AD, Hb AE, Hb AJ, Hb SS, Hb CC, Hb SC, and elevated fetal hemoglobin (up to 6 %).

It is possible that other substances and/or factors not listed above may interfere and cause false results.

Glossary of Terms

EDTA	Ethylenediaminetetraacetic Acid
ADA	American Diabetes Association
DCCT	Diabetes Control & Complications Trial
ERL	European Reference Laboratory
NGSP	National Glycohemoglobin Standardization Program
IFCC	International Federation of Clinical Chemistry
IFU	Instructions For Use
eAG	Estimated Average Glucose
HPLC	High Pressure Liquid Chromatography

How to Run a Quo-Lab A1C Test

Refer to images 1 to 12 on the inside of the back cover.

Step 1

Ensure that the analyzer is ready to run a test. Refer to the Quo-Lab Analyzer System User Manual for details.

Step 2

Remove a tray of cartridges from the carton, open the lid and place it next to the analyzer. Allow at least 50 minutes for unopened cartridges to equilibrate to the analyzer's environmental temperature before use, even if previously stored at room temperature.

Step 3

The analyzer prompts the user to scan the calibration barcode to start the test. The barcode is printed on the lid of the cartridge tray, indicated by the CAL symbol.

Place the barcode scanner over the calibration barcode and press the barcode scanner button, until you hear a beep. 'Insert New Cartridge' will be displayed on the analyzer screen along with the lot number scanned.

It is critical that the lot number displayed on the screen matches the lot number stated on the cartridges, otherwise incorrect results could be reported.

If the system returns to the home screen before a cartridge is inserted, rescan the barcode to restart the test sequence.

Step 4

Remove a cartridge from the tray. **Do not** handle the lower part of the cartridge containing the liquid.

The cartridge must be free from condensation before use. **Do not** wipe condensation from the cartridge.

Do not use the cartridge if the foil covering the top of the cartridge is damaged.

Step 5

Carefully remove the foil lid of the cartridge, ensuring that the reagent bead remains within the top of the cartridge. Once the foil lid has been removed, the cartridge must be used within 1 minute.

Step 6

Insert the cartridge into the analyzer. The analyzer will automatically detect the cartridge and will the display the message 'Cartridge Inserted. Running…'.

Do not remove the cartridge from the analyzer until the test is complete.

Do not close the door of the analyzer until instructed to do so in Step 11.

Step 7

When prompted to 'Insert Reagent', use the blunt end of the Sample Stick to push the reagent bead into the vial. Use one finger to hold the blue section of the top of the cartridge in place as the Sample Stick is removed.

There is a time limit of 60 seconds to insert the reagent bead once prompted. If the time runs out the test will be aborted. Remove the cartridge to reset the analyzer.

After the reagent bead has been added the analyzer will display the message 'Rehydrating Reagent'. This process takes approximately 50 seconds. Whilst the system is rehydrating the reagent, prepare the blood sample.

Step 8

For finger prick sample collection, the patient's finger must be warm, dry and clean, including being free of substances such as hand cream. Use a single use lancet (not provided) on the patient's finger to get a droplet of blood about the same width as the Sample Stick.

For venous samples previously collected into EDTA tubes, ensure that the sample is thoroughly mixed and use standard sampling equipment to place a droplet of blood, approximately 10 to 20 µl, onto a non-metallic and non-absorbent surface such as Parafilm® M[9].
Do not attempt to sample directly from the tube.

Step 9

Lightly touch the pointed tip of the Sample Stick to the top surface of the droplet of blood, as per image 9. The Sample Stick will take up the blood. The blood must completely fill the slot.

Note: Avoid taking up air bubbles and excess sample on the outside of the tip of the Sample Stick. **Do not** wipe excess blood from the Sample Stick.

Step 10

When the analyzer displays the message 'Insert Sample and Close Door', gently insert the Sample Stick, containing the blood sample, so that half of the tip is in the top of the cartridge. Bend the Sample Stick handle towards the front of the analyzer to break it off.

Do not push the Sample Stick fully into the cartridge otherwise the analyzer will report an error message. Refer to the Quo-Lab Analyzer System User Manual for further details.

Step 11

Close the analyzer door fully, using a single smooth movement, to complete the introduction of the Sample Stick into the cartridge. There is a time limit of 60 seconds to insert the Sample Stick and close the door once prompted. If the time runs out the test will be aborted. Remove the cartridge to reset the analyzer.

The used cartridge must be disposed of as clinical waste, as per local regulatory requirements. The analyzer will display a progress bar whilst the test is running.

You will be prompted to scan the Operator and / or Patient ID if these features have been activated. Refer to the Quo-Lab Analyzer System User Manual for further details.

Step 12

When the test has finished the result will be displayed on the screen, until the door is opened. Remove the cartridge and dispose of it as clinical waste, as per local regulatory requirements. The analyzer will be ready to run another test.

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1. Frantzen F, Studies on glycated haemoglobin and new analytical methods for its determination. Dr. Philos. Thesis. NTNU, Trondheim 1998. ISBN 82-7861-116-5.
2. The Diabetes Control and Complications trial Research Group, The Effect of Intensive Treatment of Diabetes on the Development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. N Engl J Med 1993; 329:977-986.
3. Sacks DB et al., Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. Clin Chem 2002; 48:436-472.
4. Hoelzel W et al., IFCC reference system for measurement of haemoglobin A1c in human blood and the national standardization schemes in the United States, Japan and Sweden: a method-comparison study. Clin Chem 2004; 50:166-174.
5. Nathan DM et al., A1c-Derived Average Glucose Study Group. Translating the A1c Assay into estimated average glucose values. Diabetes Care 2008; 31:1473-1478.
6. http://www.ngsp.org/docs/methods.pdf
7. http://www.ngsp.org/protocol.asp
8. Standards of Medical Care in Diabetes. Recommendations of the American Diabetes Association. Diabetes Care Vol. 34. Jan 2011; 1:S62-S69.
9. Parafilm® M is a trademark of Pechiney Plastics Packaging.

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