CE
IVD

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EN

Utilisation:
The diagnostic test strips PHAN® LAURA are intended for the semiquantitative analysis of the urine. The diagnostic test strips PHAN® LAURA are intended for in vitro diagnostic for professional use only.
Instruction:
For the objective evaluation of the diagnostic test strips, please, use the user manual for reader LAURA®, LAURA® Smart or LAURA® M.
Use only a freshly voided, well-mixed, uncentrifuged specimen without preservatives collected in a clean container. The urine should not be more than 4 hours old when tested.
1. Remove only as many test strips as are required, and reseal the tube immediately after use.
2. Do not touch test pads of the strips.
3. Completely immerse all reagent pads in specimen (no longer than 1-2 sec.).
4. Run edge of the strip against rim of urine container to remove excess urine. Hold the strip in horizontal position.
5. Evaluate the result using the reader LAURA®, LAURA® Smart or LAURA® M, follow enclosed instructions for that instrument.
Notes: For visual evaluation compare the tests pads to the corresponding colour scale on the label after approx. 60 sec. and for leucocytes after 120 sec. According to the different spectral sensitivity of the human eye and the optical system it is not possible to guarantee the exact correspondence between visual reading and reading by analyser.
Working reagent concentrations:
Specific gravity: poly(methylvinylether/maleic acid) 32 %; bromthymol blue 5.1 % / **Leucocytes:** indoxyl ester 0.43 %; diazonium salt 0.05 %
Nitrite: sulphanilamide 5.1 %; tetrahydrobenzo-[h]-quinoline 5.8 % / **pH:** methyl red 0.71 %; bromthymol blue 12.1 %; **Ascrobic Acid**: phosphomolybdic acid 26 %
Protein: tetrabromphenolphtalein ester 0.21 %; tetrabromphenol blue 0.35 % / **Glucose:** glucose oxidase 1.3 %; peroxidase 1.3 %; tetramethylbenzidine 21.0 %
Ketones: sodium nitroprusside 4.9 % / **Urobilinogen:** diazonium salt 2.3 % / **Bilirubin:** diazonium salt 0.75 % / **Blood:** tetramethylbenzidine 1.5 %; cumene hydroperoxide 15.2 %
Principle:
Specific gravity - The test is based on the principle of ion exchange, which runs between polyelectrolyte and ions present in the urine. Its result is a colour change of an acid-base indicator from the blue-green colour in the urine with low concentration of ions, through green and yellow-green in the urine with increased concentration of ions, to umber yellow colour. Using this test it is possible to determine the specific gravity of the urine in the range of 1.000 up to 1.030.
Leucocytes - The test is based on the enzymatic reaction. The test pad contains an indoxyl ester, that is cleaved by the granulocyte esterases. The released indoxyl reacts with a diazonium salt and violet coloration is formed. The colour intensity is proportional to the leucocytes amount in a sample of tested urine and is evaluated after 120 seconds.
Nitrites - The test is based on conversion of nitrate to nitrite by the action of certain species of bacteria contained in the urine. The colour test is based on the principle of the Griess test. Some degree of pink coloration should be interpreted as a positive nitrite test suggestive of 10^5 or more organisms in 1 ml of the urine specimen.
pH - The test is based on the double indicator principle and gives a range of colours from orange through yellow and green to blue and permits differentiation to within 0.5 pH unit in the range pH 5 to 9.
Ascorbic Acid - The test is based on the reaction of phosphomolybdic acid which is reduced by ascorbic acid to molybdenum blue. The test is not specific for ascorbic acid because the green to greyish-blue colour of the test pad is exhibited also by other strongly reducing substances present in urine, such as gentisic acid and other acetylsalicylic acid metabolites. We recommend to carry out determination of ascorbic acid in urine especially in cases, in which ascorbic acid may disturb the tests for other urine constituents, such as glucose, blood and nitrite.
Proteins - The test is based on colour change of acid-base indicator, which is caused by presence of proteins. It is particularly sensitive to albumin, but is much less sensitive to globulin, mucoprotein, haemoglobin and Bence-Jones protein.
Glucose - The test is based on the specific glucose oxidase/peroxidase reaction and is specific for D-glucose. The reagent pad does not react with other sugars, it reacts with presence of D-glucose by green to the dark green coloration.
Ketones - The test is based on the principle of Legal' s test and is more susceptible to the acetoacetic acid than to acetone. Test does not react with the ß - hydroxybutyric acid. The colour scale is calibrated for the acetoacetic acid.
Urobilinogen - The test is based on the coupling of urobilinogen with a stabilized reagent. The test is specific for urobilinogen and stercobilinogen and is not susceptible to the interfering factors known at Ehrlich's test.
Bilirubin - The test is based on the coupling of bilirubin with stabilized reagent. The reaction is not affected by pH of the urine.
Blood - The test is based on the peroxidase activity of haemoglobin which catalyzes the oxidation of the indicator due to the presence of the organic hydroperoxide contained in the diagnostic pad. The label contains two colour scales; for detection of intact erythrocytes and free haemoglobin. The test is highly sensitive to free haemoglobin and may detect its presence from concentrations corresponding approx. to 5 Ery/µl.
Compensation field - Pad, which isn't impregnated with any reagents. The compensation field is used for suppressing of the dark colour of the urine sample, because the dark colour could have the effect on the evaluation of reagent pads.
Limitations:
Specific gravity - The reaction is not affected by pH values of urine over 6.5 shift colour response towards lower values of specific gravity.
Leucocytes - In case when the urine sample is more markedly coloured (e.g. increased bilirubin content), the resulting colour could be affected by a sample coloration. The intensity of the colour reaction is increased by alkaline pH and higher urine density.
Nitrites - Before testing the patient should intake vegetable-rich meals and discontinue antibiotic therapy for 3 days prior to the test. Sensitivity of this test declines with the high specific gravity of the urine. Increased diuresis can cause the false negative results. Limited fluid intake prior testing can prevent from the excessive dilution of the urine. The test can be applied only at the fresh urine. Inaccurate results may occur at the stale urines, in which nitrite can be formed by contamination of the specimen.
Proteins - In strongly alkaline urines (pH >8) from patients on medication with quinine or quinoline containing drugs false positive reading may be obtained. False positive results may be found when the urine collection vessel contains traces of disinfectants with quarternary ammonium groups. On the other hand, in the presence of non-ionic or anionic detergents, false--negative results may occur. Do not take the colour of the dry pad into consideration.
Glucose - The reaction is independent on pH and presence of ketone bodies.
Ketones - Drugs and diagnostics on the basis of phenolphtalein or sulphophthalein may turn red to purple because of alkaline reaction of the pad.
Urobilinogen - The reaction is not affected by pH of the urine. The presence of bilirubin gives yellow colour. This colour which turns slowly to greenish-blue does not interfere with urobilinogen determination provided the reading is made 1 minute after wetting. The urine specimen should not be exposed to the direct sunlight as this promotes the oxidation of urobilinogen and thus leads to artificially low or false negative readings.
Bilirubin - The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of bilirubin and thus leads to artificially low or false negative readings. High concentrations of urobilinogen (above 100 µmol/l) interfere with the test. Also the red substances or the substances, which are turning red in low pH may interfere (e.g. Phenazopyridine).
Blood - Microbial peroxidase associated with urinary tract infection may cause a false positive reaction. The sensitivity of the test is influenced by specific gravity or by inhibitors of pharmacological origin.
All diagnostic pads do not interfere with the common concentration of the ascorbic acid. More information are stated at web site www.erbalachema.com
Reference values: The reference values are written in the Table I. The reference values is only approximate, it is recommended that each laboratory verify the reference values for their particular examined population.
The results: The results for reader LAURA® and LAURA® Smart are written in the Table I. The colour scale for visual reading is on the label.
Performance characteristics: The performance characteristics are written in the Table II. The values of analytical sensitivity were defined as the concentration of analyte, from which the 90% of results are positive. The analytical sensitivity isn't possible to use for SG and pH. The correlation is based on the comparison of the measured results on the reader LAURA® or LAURA® Smart with the quantitative method. The values Neg. and Pos. are the consensus with negative and positive results.
Please Note: Knowledge of the effects of drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is advisable to repeat the test after dis-continuing a drug. For the objective evaluation by using the analyser LAURA®, LAURA® Smart or LAURA® M the sensitivity of tests can be depended upon the variability of urines. The diagnostic test strips PHAN® LAURA can only be read by analyser LAURA®, LAURA® Smart or LAURA® M (see the scheme of the diagnostic strips above), they are not suitable for use with other instrument for analyse of urine. The analyser LAURA®, LAURA® Smart or LAURA® M use other method for evaluation of urine, than is the visual evaluation, therefore the concentration scale measured with analysers LAURA®, LAURA® Smart or LAURA® M are different from the colour scale on the label on the tube. This colour scale is determined for the orientation visual evaluation only. The semiquantitative analysis is not sufficient for the completing of diagnosis.
Storage: Keep the diagnostic test strips in tightly closed original tubes in a dry and dark place at (+2 to +30)°C. The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes in the laboratory. When stored under these conditions, test strips are stable to the expiry date given on the pack.
Waste disposal: Used strip should be treated as potentially infectious and should be liquidated in accordance with local and national regulations relating to the safe handling of such materials. Let waste recycle or put it to municipal waste.

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Utilisation:
The diagnostic test strips PHAN® AUTO are intended for the semiquantitative analysis of the urine. The diagnostic test strips PHAN® AUTO are intended for in vitro diagnostic for professional use only.
Instruction:
For the objective evaluation of the diagnostic test strips, please, use the user manual for reader LAURA XL®.
Use only a freshly voided, well-mixed, uncentrifuged specimen without preservatives collected in a clean container. The urine should not be more than 4 hours old when tested.
Notes: For visual evaluation compare the tests pads to the corresponding colour scale on the label after approx. 60 sec. and for leucocytes after 120 sec. According to the different spectral sensitivity of the human eye and the optical system it is not possible to guarantee the exact correspondence between visual reading and reading by analyser.
1. Remove only as many test strips as are required, and reseal the tube immediately after use.
2. Do not touch test pads of the strips.
3. Completely immerse all reagent pads in specimen (no longer than 1-2 sec.).
4. Run edge of the strip against rim of urine container to remove excess urine. Hold the strip in horizontal position.
Working reagent concentrations:
Specific gravity: poly(methylvinylether/maleic acid) 32 %; bromthymol blue 5.1 % / **Leucocytes:** indoxyl ester 0.43 %; diazonium salt 0.05 %
Nitrite: sulphanilamide 5.1 %; tetrahydrobenzo-[h]-quinoline 5.8 % / **pH:** methyl red 0.71 %; bromthymol blue 12.1 %; **Ascrobic Acid**: phosphomolybdic acid 26 %
Protein: tetrabromphenolphtalein ester 0.21 %; tetrabromphenol blue 0.35 % / **Glucose:** glucose oxidase 1.3 %; peroxidase 1.3 %; tetramethylbenzidine 21.0 %
Ketones: sodium nitroprusside 4.9 % / **Urobilinogen:** diazonium salt 2.3 % / **Bilirubin:** diazonium salt 0.75 % / **Blood:** tetramethylbenzidine 1.5 %; cumene hydroperoxide 15.2 %
Principle:
Specific gravity - The test is based on the principle of ion exchange, which runs between polyelectrolyte and ions present in the urine. Its result is a colour change of an acid-base indicator from the blue-green colour in the urine with low concentration of ions, through green and yellow-green in the urine with increased concentration of ions, to umber yellow colour. Using this test it is possible to determine the specific gravity of the urine in the range of 1.000 up to 1.030.
Leucocytes - The test is based on the enzymatic reaction. The test pad contains an indoxyl ester, that is cleaved by the granulocyte esterases. The released indoxyl reacts with a diazonium salt and violet coloration is formed. The colour intensity is proportional to the leucocytes amount in a sample of tested urine and is evaluated after 120 seconds.
Nitrites - The test is based on conversion of nitrate to nitrite by the action of certain species of bacteria contained in the urine. The colour test is based on the principle of the Griess test. Some degree of pink coloration should be interpreted as a positive nitrite test suggestive of 10^5 or more organisms in 1 ml of the urine specimen.
pH - The test is based on the double indicator principle and gives a range of colours from orange through yellow and green to blue and permits differentiation to within 0.5 pH unit in the range pH 5 to 9.
Ascorbic Acid - The test is based on the reaction of phosphomolybdic acid which is reduced by ascorbic acid to molybdenum blue. The test is not specific for ascorbic acid because the green to greyish-blue colour of the test pad is exhibited also by other strongly reducing substances present in urine, such as gentisic acid and other acetylsalicylic acid metabolites. We recommend to carry out determination of ascorbic acid in urine especially in cases, in which ascorbic acid may disturb the tests for other urine constituents, such as glucose, blood and nitrite.
Proteins - The test is based on colour change of acid-base indicator, which is caused by presence of proteins. It is particularly sensitive to albumin, but is much less sensitive to globulin, mucoprotein, haemoglobin and Bence-Jones protein.
Glucose - The test is based on the specific glucose oxidase/peroxidase reaction and is specific for D-glucose. The reagent pad does not react with other sugars, it reacts with presence of D-glucose by green to the dark green coloration.
Ketones - The test is based on the principle of Legal' s test and is more susceptible to the acetoacetic acid than to acetone. Test does not react with the ß - hydroxybutyric acid. The colour scale is calibrated for the acetoacetic acid.
Urobilinogen - The test is based on the coupling of urobilinogen with a stabilized reagent. The test is specific for urobilinogen and stercobilinogen and is not susceptible to the interfering factors known at Ehrlich's test.
Bilirubin - The test is based on the coupling of bilirubin with stabilized reagent. The reaction is not affected by pH of the urine.
Blood - The test is based on the peroxidase activity of haemoglobin which catalyzes the oxidation of the indicator due to the presence of the organic hydroperoxide contained in the diagnostic pad. The label contains two colour scales; for detection of intact erythrocytes and free haemoglobin. The test is highly sensitive to free haemoglobin and may detect its presence from concentrations corresponding approx. to 5 Ery/µl.
Compensation field - Pad, which isn't impregnated with any reagents. The compensation field is used for suppressing of the dark colour of the urine sample, because the dark colour could have the effect on the evaluation of reagent pads.
Limitations:
Specific gravity - The reaction is not affected by pH values of urine over 6.5 shift colour response towards lower values of specific gravity.
Leucocytes - In case when the urine sample is more markedly coloured (e.g. increased bilirubin content), the resulting colour could be affected by a sample coloration. The intensity of the colour reaction is increased by alkaline pH and higher urine density.
Nitrites - Before testing the patient should intake vegetable-rich meals and discontinue antibiotic therapy for 3 days prior to the test. Sensitivity of this test declines with the high specific gravity of the urine. Increased diuresis can cause the false negative results. Limited fluid intake prior testing can prevent from the excessive dilution of the urine. The test can be applied only at the fresh urine. Inaccurate results may occur at the stale urines, in which nitrite can be formed by contamination of the specimen.
Proteins - In strongly alkaline urines (pH >8) from patients on medication with quinine or quinoline containing drugs false positive reading may be obtained. False positive results may be found when the urine collection vessel contains traces of disinfectants with quartemary ammonium groups. On the other hand, in the presence of non-ionic or anionic detergents, false--negative results may occur. Do not take the colour of the dry pad into consideration.
Glucose - The reaction is independent on pH and presence of ketone bodies.
Ketones - Drugs and diagnostics on the basis of phenolphtalein or sulphophthalein may turn red to purple because of alkaline reaction of the pad.
Urobilinogen - The reaction is not affected by pH of the urine. The presence of bilirubin gives yellow colour. This colour which turns slowly to greenish-blue does not interfere with urobilinogen determination provided the reading is made 1 minute after wetting. The urine specimen should not be exposed to the direct sunlight as this promotes the oxidation of urobilinogen and thus leads to artificially low or false negative readings.
Bilirubin - The urine specimen should not be exposed to direct sunlight as this promotes the oxidation of bilirubin and thus leads to artificially low or false negative readings. High concentrations of urobilinogen (above 100 µmol/l) interfere with the test. Also the red substances or the substances, which are turning red in low pH may interfere (e.g. Phenazopyridine).
Blood - Microbial peroxidase associated with urinary tract infection may cause a false positive reaction. The sensitivity of the test is influenced by specific gravity or by inhibitors of pharmacological origin.
All diagnostic pads do not interfere with the common concentration of the ascorbic acid. More information are stated at web site www.erbalachema.com
Reference values: The reference values are written in the Table I. The reference values is only approximate, it is recommended that each laboratory verify the reference values for their particular examined population.
The results: The results for reader LAURA XL® are written in the Table I. The colour scale for visual reading is on the label.
Please Note: Knowledge of the effects of drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is advisable to repeat the test after dis-continuing a drug. For the objective evaluation by using the analyser LAURA XL® the sensitivity of tests can be depended upon the variability of urines. The diagnostic test strips PHAN® AUTO can only be read by analyser LAURA XL® (see the scheme of the diagnostic strips above), they are not suitable for use with other instrument for analyse of urine. The analyser LAURA XL® use other method for evaluation of urine, than is the visual evaluation, therefore the concentration scale measured with analysers LAURA XL® are different from the colour scale on the label on the tube. This colour scale is determined for the orientation visual evaluation only. The semiquantitative analysis is not sufficient for the completing of diagnosis.
Storage: Keep the diagnostic test strips in tightly closed original tubes in a dry and dark place at (+2 to +30)°C. The strips must be kept away from moisture, direct sunlight, elevated temperature and chemical fumes in the laboratory. When stored under these conditions, test strips are stable to the expiry date given on the pack.
Waste disposal: Used strip should be treated as potentially infectious and should be liquidated in accordance with local and national regulations relating to the safe handling of such materials. Let waste recycle or put it to municipal waste.

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EN

Store at (+2 to +8) °C

Specifications:
URINORM control urines are a liquid synthetic material with concentrations within the normal (URINORM N) and pathological (URINORM P) ranges.

Usage: URINORM control urines are designed for verification and confirmation of precision and accuracy of PHAN (Uro-dip 10*e*) diagnostic strips as well as LAURA, LAURA Smart, Uro-dipcheck 240*e*, Uro-dipcheck 400*e*, LAURA M urine analysers (for the visual and objective evaluation).
URINORM control urines are designed for direct usage (*ready to use*). This product is to be treated in the same manner as patient specimen and used in accordance with the instructions for PHAN diagnostic strips and LAURA, LAURA Smart, Uro-dipcheck 240*e*, Uro-dipcheck 400*e*, LAURA M analysers.

Reagents:
URINORM N – **normal (negative) level**
URINORM P – **pathological (positive) level**
URINORM control urines contain preservatives added to ensure the required stability and avoid microbial growth. **URINORM control urines do not contain human urine.**

Stability:
This product will be stable until the expiration date when stored unopened at (+2 to +8) °C; the product is to be protected against direct light and contamination.
Once opened, URINORM control urines remain stable for 30 days (blood analyte is stable for 20 days) or 10 measurements (whichever occurs first). Such used control urines are to be stored at (+2 to +8) °C – in a closed container and protected against direct light and contamination.
Mark the date of the first opening of vials on the appropriate place on the label of the vial (see cell **Opened**) containing control urines. Do not store this product in freezing appliances.

Preparation:
1. Allow the control urines to reach the room temperature (+18 to +25) °C;
2. Shake up a vial containing control urine (homogenization);
3. Open the vial carefully;
4. Immerse a PHAN - Uro-dip 10*e* diagnostic strip in the vial (1-2 seconds);
5. Remove excess urine from the strip (wipe the rim of the strip on the vial);
6. Perform a measurement (using LAURA, LAURA Smart, Uro-dipcheck 240*e*, Uro-dipcheck 400*e*, LAURA M analysers or visually);
7. Perform an evaluation (comparison with the target ranges listed in the table);
8. Close the vial containing control urine;
9. Store at (+2 to +8) °C.

Values:Values for LAURA, LAURA Smart, Uro-dipcheck 240*e*, Uro-dipcheck 400*e*, LAURA M urine analysers and for visual evaluation.

Analyt	URINORM N	URINORM P	Units
Specific gravity	1.020 - 1.030	1.000 - 1.020	
Leucocytes	negative	25 - 500 / 1+ – 3+	Leu/µl / ARB
Nitrites	negative	positive	
pH	5 - 6.5	8 - 9	
Proteins	negative	1 - 5 / 100 - 500 / 2+ – 3+	g/l / mg/dl / ARB
Glucose	negative	17 - 55 / 300 - 1000 / 3+ – 4+	mmol/l / mg/dl / ARB
Ketones	negative	1,5 - 15 / 16 - 156 / 1+ – 3+	mmol/l / mg/dl / ARB
Urobilinogen	normal	17 - 203 / 1 - 12 / 1+ – 4+	µmol/l / mg/dl / ARB
Bilirubin	negative	17 - 102 / 1 - 6 / 1+ – 3+	µmol/l / mg/dl / ARB
Blood	negative	10 - 250 / 1+ – 3+	Ery/µl / ARB

Limitations:
• This product may not be used after the expiration date specified on the package and it is not intended to be used as the standard.
• This product is not to be used after the elapse of 30 days from the day of first opening or after 10 measurements (whichever occurs first).
• If there is evidence of microbial contamination or excessive turbidity in the product, discard the vial.

Health protection:
For *in vitro* diagnostic use. To be handled solely by an authorized, trained personnel.
URINORM N control urines contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H -isothiazol-3-one (3:1).

Warning

Hazard statement:
H317 May cause an allergic skin reaction.
Precautionary statement:
P280 P280 Wear protective gloves/protective clothing/eye protection.
P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
P332 + P313 If skin irritation occurs: Get medical advice.

URINORM P control urines are not classified as hazardous.

First aid: In case of an accidental ingestion wash up the mouth and drink approximately 0.5 l of water. In case of eye contact rinse the eye quickly and thoroughly using a stream of clean water. Contaminated skin is to be washed using warm water and soap. Consult a physician in all serious cases of health damage.

Waste disposal: Residues of control urines, which are not classified as infectious material, are to be disposed of in accordance with the internal regulations and by-laws. Paper and other packaging materials are to be discarded in the form of sorted waste (paper, glass, plastic).

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EN

CE IVD

INTENDED USE
URINORM XL control urines are liquid synthetic material with concentrations in the normal (URINORM XL N) and pathological (URINORM XL P) ranges. URINORM XL control urines are designed to verify the precision and accuracy of PHAN AUTO diagnostic strips and automatic urine analyzer LAURA XL intended for *in vitro* diagnostics. This product is intended also for evaluation of sediment. For professional use only.

PRODUCT DESCRIPTION:
URINORM XL control urines are human urine-based containing constituents of human and animal origin as well as purified chemicals. URINORM XL control urines contain preservatives added to maintain product stability and avoid the microbial growth. URINORM XL urine control is a ready-to-use.

REAGENTS:
URINORM XL N– normal level 120 ml
URINORM XL P– pathological level 120 ml

STORAGE AND STABILITY:
URINORM XL control urines are stable until the expiration date unopened stored at 2–8 °C and protected from light and microbial contamination.
Once opened URINORM XL control urines are stable for 30 days at 20–25 °C and 18 months at 2–8 °C. Do not store the product in freezers.

INSTRUCTIONS:
Control urines URINORM XL should be treated the same as patient specimens.
1. Allow the control urines bottles to reach the room temperature
2. Gently invert the bottles to homogenize the sample
3. Remove the cap and pour 1–1.5 mL into the Eppendorf tube
4. Perform the QC measurement on LAURA XL urine analyzer – device automaticly evaluates the results and compare them with declared values
5. Close bottles with control urines
6. Store approprietly

DECLARED VALUES:

Table 1: Declared values for chemistry on LAURA XL analyzer:

Analyte	URINORM XL N	URINORM XL P	Units
Specific gravity	1.020–1.030	1.000–1.015	
Nitrites	negative	positive	
pH	5–6.5	7–9	
Protein	negative	0.3–5 / 30–500	g/L / mg/dL
Glucose	negative	5.5–55 / 100–1000	mmol / mg/dL
Ketones	negative	0.5–5 / 5–52	mmol / mg/dL
Urobilinogen	normal	17–203 / 1–12	µmol/L / mg/dL
Bilirubin	negative	51–102 / 3–6	µmol/L / mg/dL
Leucocytes	negative	75–500	Leu/µL
Blood	negative	50–250	Ery/µL

Table 2: Declared values for sediment on LAURA XL analyzer:

URINORM XL N			
Method	**Red blood cells**	**White blood cells**	**Crystals**
LAURA XL	0–2 hpf	0–5 hpf	none

URINORM XL P			
Method	**Red blood cells**	**White blood cells**	**Crystals**
LAURA XL	3–100 hpf	3–100 hpf	present

LIMITATIONS:
Control urines should not be used past the exspiration date on the label.
Do not use turbid control urines where there is an evidence microbial contamination.

WARNING AND PRECAUTIONS
For *in vitro* diagnostic use. To be handled solely by an authorized, trained personnel. Human source material. Treat as potencially infectious.

Hazards identification in accordance with Regulation (EC) No 1272/2008
Reagents of the kit are not classified as dangerous.

FIRST AID:
In case of an accidental ingestion wash up the mouth and drink approximately 0.5 L of water. In case of eye contact rinse the eye quickly and thoroughly using a stream of clean water. Contaminated skin is to be washed using warm water and soap. Consult a physician in all serious cases of health damage.

WASTE DISPOSAL:
Residues of control urines are to be disposed of in accordance with the internal regulations and by-laws. Paper and other packaging materials are to be discarded in the form of sorted waste (paper, glass, plastic).

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