RISHENA
瑞神安医疗
Ref: IFU-202-001.EN

Lead	VL202A VL202B VL202C
Tunneling tool	ST205

Implantable Vagus Nerve Stimulation Lead

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Table of Contents

1 Overview 1

1.1 Product Structure and Composition 1
1.2 Related products 1
1.3 Product Scope 1
1.4 Package Contents 2

1.4.1 Lead Package 2
1.4.2 Tunneling tool Package 2

2 Technical Specifications 2

2.1 Lead assembly 2
2.2 Tunneling tool 2
2.3 Product category 3
2.4 Electrical Properties 3
2.5 Transport and Storage Conditions 3
2.6 Lead Connection 3

3 Operation Instructions 3

3.1 Preoperative Preparation 3
3.2 Implant location 4
3.3 Implanting steps 4

3.3.1 Implant the Lead 4
3.3.2 Using the Tunneling tool 5
3.3.3 Placing the Lead electrodes 5
3.3.4 Connect the lead to the stimulator 9
3.3.5 Verifying impedance and stitching pockets 12

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1 Overview

1.1 Product Structure and Composition

The Implantable Vagus Nerve Stimulation Lead consists of VNS therapy Lead and the tunneling tool.

The VNS therapy lead assembly (hereinafter referred to as the lead) is composed of IS-1 connector、 silicone tube、 anchor tether、 electrodes and tie-downs. Structural diagram is shown in Figure 1.

Note that the tie-down in the diagram is a separate part used for suture fixation of the lead during the surgical operation. Each lead is packed with 3 tie-downs.


Figure 1 Diagram of VNS therapy lead assembly

The tunneling tool is composed of center shaft、 plastic cannula and insertion-tip. The structural diagram is shown in Figure 2.


Figure 2 Structural diagram of the tunneling tool

1.2 Related products

The Implantable Vagus Nerve Stimulation Lead should be used with the following products from Rishena Medical Co., Ltd..

a) Implantable Vagus Nerve Stimulator: VNS201A, VNS201B
b) Neuromodulation Programmer: VP213A/ VP213B

1.3 Product Scope

The lead combined with specific implantable vagus nerve stimulator is used for drug resistant or intractable epilepsy patients (for 6 years old and above), intended for adjunctive treatment of refractory and drug-resistant epilepsy.

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1.4 Package Contents

1.4.1 Lead Package

VNS Therapy lead	x1
Tie-downs	x3
Instruction Manual	x1
Registration form	x1
Quality Certificate	x1
Product labels	x6

1.4.2 Tunneling tool Package

Tunneling tool	x1

2 Technical Specifications

2.1 Lead assembly

a) Model: VL202; Specification: A; B; C .
b) Specifications and parameters information for the lead is presented in Table 1.

Table 1 Lead dimensions

Specification	A	B	C
Length	45cm		
Diameter	2.2mm		
Weight	2.3g		
Inner diameter of spiral electrode	1mm	2mm	3mm

Note: the size A, B, and C are used to suit different thicknesses of a patient's vagus nerve.

c) The component materials of the lead are presented in Table 2

Table 2 Lead component materials

Component	Material
IS-1 connector	316L Stainless Steel
Silicone tube	Silicone
Electrode contact	Platinum/Iridium Alloy
Suture	Polyester
Tie-down	Barium sulfate silicone
Electorates	MP35N (MP-35N alloy)

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c) The tunneling tool materials: see Table 4.

Table 4 materials of the tunneling tool components

Component	Center shaft	Threading -tip	Plastic cannula
Material	304 stainless steel	304 stainless steel	Polyurethane

2.3 Product category

Usage: Sterilized tool, single use only.
Shelf life: 3 years.

2.4 Electrical Properties

Lead resistance: see Table 5.

Table 5 Lead resistance

Connector ring 1 to terminal electrode	170Ω
Connector ring 2 to terminal electrode	170Ω

2.5 Transport and Storage Conditions

Temperature: -10°C ~ + 60°C
Atmospheric Pressure: 500hPa ~ 1060hpa
Relative Humidity: 10% ~ 85%
Keep in a dry and ventilated place.

2.6 Lead Connection

During surgical operation the lead is connected to the stimulator, also from Rishena Medical Co., Ltd, through a plug-in connector to conduct the stimulus signal to the vagus nerve. The maximum holding force of the connection is 10N.

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Please make sure that the aseptic packages of the lead and the stimulator are without any damage before opening. If damage is observed, do not use, please return the product to Rishena Medical Co., Ltd.

3.2 Implant location

Two incisions are required for the implanting, one for electrodes to be encircled onto the vagus nerve, and the other for placing the stimulator. A tunnel is created between these two incisions using the tunneling tool for the electric lead to pass through.

Surgically create a subcutaneous 'pocket' for implanting vagus nerve stimulator just below the clavicle, under the skin of the upper left chest.

Surgically open the skin to expose the left vagus nerve of the neck half-way between the clavicle and the mastoid. The spiral contact electrode tips on one end of the lead are wound around the vagus nerve, while the connector end of the lead is tunneled to the implanting 'pocket' on the chest using a subcutaneous tunneling tool (ST205), see Figure 3.


Figure 3: Placement of stimulator and lead

3.3 Implanting steps

3.3.1 Implant the Lead

1) After administering appropriate anesthesia to the patient, expose the left carotid sheath as it extends along the anterior border of the sternocleidomastoid muscle.

2) Locate and expose at least 3-cm section of the left vagus nerve, approximately half-way between the clavicle and the mastoid, where it is clear of branches (below where the superior and inferior cervical cardiac branches separate from the vagus nerve). The nerve usually lies in a posterior groove between the carotid artery and internal jugular vein.

3) Create a subcutaneous pocket in the chest below the clavicle for the stimulator.

Note: Do not expose the lead to dust or other similar particulates, because its

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silicone insulation can easily attract particulate matter.

Note: Do not soak the lead in saline or similar solution before implanting it, because this may cause the insulated portions of the connector pin to swell and become difficult to insert into the stimulator.

3.3.2 Using the Tunneling tool

The lead connector and lead body can be tunneled subcutaneously from the neck incision site to the stimulator in the chest pocket, after placement of the leads and anchor tether on the nerve and the placement of strain relief with the tie-downs.

Caution: To maximize system performance and minimize possible mechanical damage to the nerve or lead, pay careful attention to lead routing, lead stabilization, and lead placement.

Caution: Never route the lead through muscles.

Caution: Never suture the lead or lead body to muscle tissue.

Caution: Always use the tie-downs.

Caution: Do not place sutures directly on the lead body. Doing so may result in insulation damage or wire failure, causing premature failure of the lead.

To pass the tunneling tool, do the following:

1) Place the insertion-tip end of the tunneling tool through the neck incision and tunnel subcutaneously toward the chest incision, exerting force on the handle end and directing the tunneling tool as necessary.

2) After the insertion tip has passed from one incision site to the other, unscrew the tip and withdraw the shaft from the cannula, leaving the cannula extended through both incisions)

3) With the cannula in place between the two incisions, carefully insert the lead connector inside the end of the cannula at the neck incision.

4) Carefully pull the cannula, along with the lead connector, from the chest incision end until the lead connector completely exit the chest incision.

5) Remove the lead connector from the cannula, leaving the lead electrodes at the neck incision site.

6) Discard the tunneling tool after use.

3.3.3 Placing the Lead electrodes

It is very important that the surgeon implanting the VNS Therapy System be familiar with vagus nerve anatomy (see Figure 4), particularly the cardiac branches. The lead leads must not be placed on either the superior or the inferior cervical cardiac branches. Place the lead below where the superior and inferior cardiac branches separate from the vagus nerve. Stimulation of either of these two branches during the System Diagnostics (Lead Test) may cause bradycardia and/or asystole. Careful dissection laterally on the vagus nerve should aid the physician in determining proper lead placement. In most but not all patients, the main vagus nerve is the largest of the three nerves. The spiral leads and anchor tether are coiled around the nerve, beginning with the lead that is farthest from the lead bifurcation (with a blue suture) . This lead should be nearest the patient's head.

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
Figure 4 Location of the vagus nerve

Depending on the surgeon's preference, the spirals can alternately be placed by putting the anchor tether on first (distal to head), next placing the lead closest to the lead bifurcation (with yellow suture), and then placing the lead farthest from the lead bifurcation (with blue suture), see Figure 5. For correct pulse polarity, do not reverse the position of the electrodes.


Figure 5 Electrode placement position on the vagus nerve

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The spiral electrodes can be placed on the nerve as described below. As an alternative, each spiral can be placed underneath the nerve before it is spread. Some tools may be useful to separate the nerve from tissue during the procedure.

Step 1) Coil the first spiral (with blue suture) in the following manner:

a. With forceps, gently pull each end of the spiral to spread the spiral electrode, see Figure 6.


Figure 6 Spread the spiral electrode

b. Parallel to the exposed nerve, turn the spiral clockwise at a 45 degree angle to the nerve (see Figure 7).


Figure 7 Angle the spiral electrode

c. Place the turn of the spiral where the lead connects to the spiral (the section with the metal ribbon) onto the nerve (see Figure 8).

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
Figure 8 Placing the spiral electrode

d. Pass the distal suture under the nerve and back around so that it encircles the nerve (see Figure 9 and Figure 10).


Figure 9 Pass the tasteless suture under the nerve


Figure 10 Pull the sutures out from behind

e. Pass the proximal suture under the nerve and back around so that it encircles the nerve (see Figure 11).

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
Figure 11 Spiral electrode encircles the nerve

Step 2) Repeat steps a. through e. for the middle spiral electrode (with yellow suture).

Step 3) Place the anchor tether (with blue suture) around the nerve, following the same general steps as for the other two spirals.

Step 4) After all three spirals have been coiled around the nerve, verify that the two lead branch sections exit each spiral in the same direction and that the two sections are aligned parallel to each other and to the nerve. The correct placement of the two spiral lead sections and anchor tether is shown in Figure 12.


Figure 12 Placement of lead, electrodes and anchor tether

3.3.4 Connect the lead to the stimulator

Step 1) Insert the connector of the lead into the stimulator, and check the connection between the stimulator and the lead .

a) Ensure that the lead connector pin and the generator socket are clean before insert the lead connector pin into the stimulator socket.

b) Insert the torque screwdriver through the center of silicone plug into the inner hexagon groove of the fastening screw.

c) Insert the lead connector into the stimulator socket, the lead connector must be fully inserted into the bottom of the stimulator (see Figure 13).

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
Figure 13 Insert the lead connector into the stimulator

d) If the fastening screw obstructs the connector from being inserted properly, loosen the screw counterclockwise with the torque screwdriver until the connector plug can be easily and fully inserted (Figure 14).

Reminder: sufficiently loosen the fastening screw to ensure that the connector plug can be smoothly inserted fully into the stimulator socket without obstruction.

Note: Loosen 1~2 turns counterclockwise is usually enough. Avoid too many turns that may detach the fastening screw from the connector of the stimulator.


Figure 14 Loosen fastening screw with torque screwdriver

Step 2) When the lead connector is fully inserted into the stimulator socket, tighten the screw as follows:

Tighten the screw by rotate the torque screwdriver clockwise until a "click" sound is heard, then pull out the torque screwdriver (see Figure 15).

Reminder: The torque screwdriver is a disposable tool, and its performance may not be guaranteed if used multiple times.

Do not tighten the screw too hard. Excessive tightening may permanently damage the screw or the socket.

After the torque screwdriver is removed, check the silicone plug hole to ensure that the silicone plug hole has automatically closed. If the silicone seal is not fully closed, fluid infiltration may occur, and the patient may be shocked, burned, or stimulated at the stimulator implant site, it may cause intermittent or lose of stimulation altogether.

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
Figure 15 Tighten screw with torque screwdriver

Step 3) Place the stimulator in the subcutaneous pocket with the side of the RISHENA logo facing the skin. Place the stimulator carefully to ensure that the lead do not distort significantly（Figure 16）.

Reminder: to reduce the possibility of skeletal muscle stimulation, place the stimulator away from bone tissue and skeletal muscle tissue.


Figure 16 Place stimulator into pocket

Reminder: Do not wrap the lead around the front or back of the stimulator. Slack the excess lead around the outside edge of the stimulator. This can avoid increasing the depth of the subcutaneous pocket, reducing the potential risk of product replacement surgery and the possibility of lead distortion (Figure 17).


Figure 17 Wrapping lead

Step 4) Secure the stimulator in the subcutaneous pocket utilizing the suture hole of the stimulator connector.

Reminder: failure to fully immobilize the stimulator through the suture hole increases the risk of device movement or rotation, which may result in device damage, skin erosion, and accidental movement of the stimulator or lead .

Step 5) Check the functioning of the stimulator with the external programmer.

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Refer to the next section: "verifying impedance and stitching pockets."

3.3.5 Verifying impedance and stitching pockets

Step 1） Diagnose the stimulator system and check the output impedance before stitching the pocket.

Reminder: It is necessary to check before stitching the pocket. Read the resistance of stimulator with the external programmer. If the test resistance is around 3k ohm, the electrical integrity of the system is verified.

Caution:

If the external programmer finds abnormal impedance, the lead should be pulled out, replace it with the resistance rod inserted into the stimulator. Reread the resistance from the stimulator. A display of around 3k ohm indicates that the stimulator is functioning normally, and the possible problem is with the lead connection or the lead itself. Reconnect the lead and check the electrode placements. If problem persists, the lead may have to be replaced.

- Confirm that the stimulator is secured to the subcutaneous pocket through the suture hole of the connector.
- Suture and bandage the incision.

Step 2） After the incision is bandaged, re-check the resistance of the system before the patient leaves the room. This precautious step further verifies the electrical connectivity of the lead to the stimulator after all surgical procedures.

4 Safety Precautions

The warning signs and legends shown in the manual are intended to enable you to use the product safely and correctly, and to prevent injury to you and others.

Warning signs, legends and their meanings:

Warning sign	Meanings
⚠Contraindication	Risk of injury or serious injury if used incorrectly
⚠Warning	Possibility of injury or serious injury if used incorrectly
⚠Cautions	Possibility of injury or product damage if used incorrectly

Legends	Meanings
	Mandatory
	Forbidden

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⚠Contraindication	
1、**Electrotherapy:** patients implanted with the vagus nerve stimulator system shall not use short-wave electrotherapy, microwave electrotherapy or ultrasonic electrotherapy (collectively referred as electrotherapy). Because electrotherapy can heat the stimulator seal housing and damage the tissue near the implant site.	⊘
2、**Anticoagulant therapy:** patients with anticoagulant therapy may have a greater risk of postoperative complications, such as hematoma.	⊘
3、Electrocoagulation hemostasis: Electrocoagulation hemostasis may damage the lead or the stimulator, temporarily suppress the output of the stimulator or may reset the stimulator parameters (output amplitude =0 V), which requires the clinician to reprogram the stimulator. Electrocoagulation may also induce current in the lead of the stimulator system, which may cause harm or further injury to the human body. When using electrocoagulation, please follow the following warning: turn off the stimulator before electrocoagulation, do not put the lead in contact with the electrocoagulation device, only a bipolar electrocautery is recommended, if a unipolar electrocautery must be used, do not use high voltage mode. Set the electricity as low as possible, keep the current path (ground plane) as far away from the stimulator and lead as possible, and confirm whether the stimulator function is normal after electrocoagulation.	⊘
4、**Radiotherapy:** do not face the vagus nerve stimulation system to a source of radiation, such as cobalt 60 or gamma rays. If radiation therapy is needed near the vagus nerve stimulation system, shield the nerve stimulation device to prevent radiation damage.	⊘
5、**High gain output ultrasonic/lithotripsy:** it is not recommended to use high gain output ultrasonic devices such as electro-hydraulic lithotripsy for patients with implantable vagal nerve stimulation system. If lithotripsy is necessary, keep the focus at least 15 cm from the stimulator.	⊘
6、**Defibrillation/cardioversion:** when the patient is in ventricular or atrial fibrillation, the first consideration is the safety of the patient's life. External defibrillation or cardioversion may cause permanent damage to the nerve stimulation system. Use of defibrillation or cardioverter near the stimulator is not recommended. If it is necessary for external defibrillation or after cardioversion, minimize the current flowed through stimulator and conduit system as possible, the following operations: locating rod as far as possible away from stimulator, locating rod with nerve stimulation system in vertical position, with minimal power output (watts seconds) for clinical,	⊘

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dysfunction, etc.	
7、The lead can only be used in conjunction with the stimulator and programmer specified in this manual.	❗
8、Treatment-grade ionizing radiation may cause damage to the electronic components of the vagus nerve stimulator, and damage to the stimulator may not be immediately noticeable.	❗
9、Recommended precautions should be taken to prevent the harmful effects of changes in the performance of stimulator.	❗
10、Preventive measures should be taken to avoid harmful effects of certain adverse circumstances (e.g., electromagnetic interference, extreme temperature, pressure changes) on the patient.	❗
11、Medical advice should be sought before entering an environment that may adversely affect the functioning of the VNS therapy system.	❗
12、If the lead has been implanted in another patient before, it must not be reused.	⊘

⚠Cautions	
1、Before use, please check the sterile packaging of the lead. If the seal or packaging is damaged, do not use and contact the after-sale service unit.	❗
2、Medical staff shall be trained before the first operation.	❗
3、When using sharp tools near the lead, special care should be taken to avoid puncturing the insulation layer of the lead.	❗
4、Do not use salt water or other ionic liquids at the joints, as this may lead to a short circuit.	❗
5、During the transportation, storage and use of the lead, collision and vibration should be avoided as far as possible. In addition, temperature, humidity and air pressure requirements specified in this manual should be ensured.	❗

5 Troubleshooting

Fault phenomenon	Possible Cause Analysis
The stimulator detects abnormal impedance	Check whether the lead connector is fully inserted into the stimulator socket.
	Check whether the stimulator fastening screw is tightened.

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	Check whether the spiral electrodes of the lead properly encircled to the vagus nerve.
	Check whether the stimulator impedance is normal with a resistance bar, and if so, check the lead.

If problems persist, please contact service department of Rishena Medical Co., Ltd..

6 Possible Complications

In addition to the risks of surgery, implantation and usage of the vagus nerve stimulator may result in the following adverse reactions. Some adverse symptoms may require surgery.

Nausea, vomiting

Hoarse voice and uncomfortable throat

Allergies or adverse immune reactions

Change in sensitivity to a stimulus, or what some patients call discomfort (tremor or shock).

Infection

Newly occurred pain

Pain at system implants (stimulators and leads)

Hematomas, bleeding, hematomas and/or strokes

The position of the stimulator or lead may move or corrode

Possible nerve damage

Possible equipment technical problems

Instantaneous shock

If necessary, the implanted vagus stimulator system should be removed.

7 Service life

The lead has a service life of at least 10 years. A lead would require replacement if a lead discontinuity were suspected, accompanied by increased symptoms (Such as increased seizures).

Events that can shorten the useful life of the lead are as follows:

- Blunt trauma to the neck and/or any area of the body.
- Patient's excessive twisting or stretching at either the implanted lead or the pulse generator.
- Improper surgical implantation of the VNS Therapy System, such as strain relief loop is inadequate, placing sutures directly on the lead body, not using the tie-downs, etc.

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8 Product scrap treatment

If the product cannot be used or replaced due to faults, it shall be scrapped in accordance with local laws and regulations.

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MRI Safety Conditions Declaration for Implantable Vagus Nerve Stimulator (VNS)

To Whom It May Concern,

This document serves as an official declaration of the safety conditions under which a patient with an Implantable Vagus Nerve Stimulation (VNS) system may undergo a Magnetic Resonance Imaging (MRI) examination.

Implantable vagus nerve stimulator magnetic resonance environmental conditions safety information	
The Implantable vagus nerve stimulation system shall be scanned by MRI under the following conditions.	
Model and specification of Implantable Vagus Nerve Stimulator	VNS201A, VNS201B
Model and specification of Implantable Vagus Nerve Stimulation Lead	VL202A, VL202B, VL202C
Maximum number of matching electrode wires	1
System integrity	No open circuit or short circuit
Device setup	Stimulus output off
Type of magnetic resonance equipment	Cylindrical superconducting magnet placed horizontally
Multinuclear imaging	/
Static magnetic field strength	3.0T
RF frequency	128MHz±5%
Maximum spatial field gradient	≤19T/m
Maximum gradient switching rate	≤180T/m/s, per axle
RF transmitting coil	Body coil
RF receiving coil	Unlimited
Running mode of magnetic resonance equipment type	Normal operation mode
RF safety limits	Whole body SAR value should be ≤1.58W/kg
Continuous scan time limit	≤ 30 minutes
Scanning area	Whole body scan
Image artifact	The product is expected to be Implantable into the patient's chest, which may lead to image artifacts.

Address: Suite D3006-D3008, No.26, Huashan Middle Road, Xinbei District, Changzhou Jiangsu, China
Tel:+86-519-69989503 Website: www.rishena.com

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Production license No.: NO.SSYJXSCX20200033
Product registration certificate No.: NO.GXZZ20203120084
Product technical requirement No.: NO.GXZZ20203120084
Manufacturer: Rishena Medical Co., Ltd.
Production address: Suite D3006 - D3008 No.26 Huashan Middle Road Xinbei District Changzhou Jiangsu, China
Tel.: +86- 519-6998 9503
Registrant: Rishena Medical Co., Ltd.
Address: Suite D3006 - D3008 No.26 Huashan Middle Road Xinbei District Changzhou Jiangsu, China
Tel.: +86- 519-6998 9503
After sales service unit: Rishena Medical Co., Ltd.

Chinese manual Version information:Revision date：20190505 Edition C/0
IFU-202-001
This manual Version information:Revision date：20221123 Edition A/0
IFU-202-001.EN

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