Concept Medical

CONCEPT MEDICAL B.V. NETHERLANDS

MagicTouch ED

SIROLIMUS COATED BALLOON CATHETER

EN	Sirolimus Coated Balloon Catheter Instructions for Use	Page 04
EL	Καθετήρας-μπαλόνι με επικάλυψη σιρόλιμους Οδηγίες χρήσης	Page 10
DE	Sirolimus beschichteter Ballonkatheter Gebrauchsanweisung	Page 17
IT	Catetere a palloncino rivestito con Sirolimus Istruzioni per l'uso	Page 24
ES	Catéter con Balón Recubierto de Sirolimus Instrucciones de uso	Page 30
PT	Cateter de balão revestido com Sirolimus Instruções de uso	Page 37
RU	Баллонный катетер, покрытый сиролимусом Инструкции по использованию	Page 43
TR	Sirolimus Kaplı Balon Kateter Kullanım Talimatları	Page 51
SK	Balónikový katéter s povlakom Sirolimus Návod na používanie	Page 57

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ENGLISH

1.0 PRODUCT DESCRIPTION

Magic Touch - ED Sirolimus Coated Balloon Catheter is a device/drug combination product consisting of two components :1). a standard balloon catheter 2). coated with encapsulated sirolimus formulation. Generic name of Magic Touch - ED is "Sirolimus Coated Balloon Catheter" and GMDN Code is 62551.

The coating of the Magic Touch - ED utilizes Nanolute technology to deliver polymer-free sirolimus encapsulated in a Generally Recognized As Safe (GRAS) phospholipid excipient. The amphiphilic properties of the phospholipid allow the Sirolimus to remain encapsulated and protected from degradation for a time sufficient to ensure treatment of the target lesion. Upon contact of the balloon with the arterial lumen, the hydrophilic surface of the nanocarrier becomes wet, encouraging maximum drug delivery upon subsequent balloon expansion. Coating is on circumferential surface of the balloon length ensures homogeneous drug delivery.

The use of Nanolute technology provides for very low in-transit drug loss resulting in reduced shedding of particulate matter that can contribute to downstream embolization. This technology also allows for better drug diffusion and drug retention throughout the artery, as well as enhanced bioavailability due to deeper penetration of the drug into the artery tissue.

1.1 DEVICE COMPONENT DESCRIPTION

The Magic Touch - ED Sirolimus Coated Balloon Catheter is indicated for the treatment of erectile dysfunction for male patients with de novo and restenotic lesions in the internal pudendal, Common penile, dorsal, or urethral arteries that are 1.50 mm to 4.00 mm in diameter and up to 36 mm lesion length.

The balloon catheter is a monorail (Rapid Exchange) catheter with a semi-compliant inflatable balloon at the distal part and an atraumatic, tapered tip to aid in crossing the tight stenosis. The catheter is double lumen in the distal part, one lumen is used for inflation and deflation of the balloon, and the second lumen allows access to the distal part of the catheter for guide wire insertion. The maximum recommended guide wire diameter is 0.014". A luer lock fitting (Hub) allows the connection with an inflation device. The internal tubing of the balloon has two radio opaque markers (Platinum - Iridium)to provide visual reference points for balloon positioning relative to the stenosis within the vessel. The catheter working length is 140cm and is available in diameters between 1.50 mm and 4.00mm, and lengths of 10, 15, 20, 25, 30, 35 and 40 mm as shown in Table 1.


Figure 1: Sirolimus Coated Balloon Catheter

Table 1: Device Specification

Brand Name	Magic Touch - ED
Available Lengths (mm)	10, 15, 20, 25, 30, 35, 40
Available Diameters (mm)	1.50, 2.00, 2.25, 2.50, 2.75, 3.00, 3.50, 4.00
Drug Component	Sirolimus Drug (also known as Rapamycin)
Excipient	Lipoid E 80 (Egg phospholipid)
Drug Dose	1.27µg/mm^2
Delivery System Workable Length	140cm (1400mm)
Balloon Inflation Pressure	Nominal Inflation Pressure: 8 bar Rated Burst Pressure: 14 bar*
Guide wire compatibility (max.)	(max.) 0.014"
Sheath Compatibility	6F (2.00mm)
Guiding Catheter Inner Diameter	6F (1.98mm) (inner lumen ≥ 0.078")
Catheter Shaft Outer Diameter	Proximal 1.95 F (0.65 + 0.025 - 0.005 mm) / Distal 2.67F (0.89 ± 0.02 mm) / Intermediate 2.55F (0.85 ± 0.02 mm)
Catheter Configuration	Rapid Exchange (Monorail)
Radiopaque Marker Bands	90% Pt / 10% Ir alloy
Sterilization	Ethylene Oxide Gas

*Note: *Do not exceed Rated Burst Pressure (RBP)*

1.2 DRUG COMPONENT DESCRIPTION

Active Pharmaceutical Ingredient (API): Sirolimus is a macrocyclic lactone produced by Streptomyces hygroscopius.

Sirolimus appearance: White to off-white powder.

Sirolimus solubility: Freely soluble in chloroform, acetone and acetonitrile and insoluble in water and melts at approximately 183°C -185°C

Sirolimus Molecular formula: $C_{51}H_{79}NO_{13}$

Molecular weight: 914.2g/mol

CAS Registry no.: 53123-88-9

Chemical name: (3S, 6R, 7E, 9R, 10R, 12R, 14S, 15E, 17E, 19E, 21S, 23S, 26R, 27R, 34aS) - 9, 10, 12, 13, 14, 21, 22, 23, 24, 25, 26, 27, 32, 33, 34, 34a - hexadecahydro - 9, 27 - dihydroxy - 3 - [(1R) - 2 - [(1S, 3R, 4R) - 4 - hydroxy - 3 - methoxycyclohexyl] - 1 - methylethyl] - 10, 21 - dimethoxy - 6, 8, 12, 14, 20, 26 - hexamethyl - 23, 27 - epoxy - 3H - pyridol [2, 1-c] [1, 4] oxaazacyclohentriacontine - 1, 5, 11, 28, 29 (4H, 6H, 31H) - pentone.


Figure 2: Structural Formula of Sirolimus

Table 2: Product Sizes, Reference and Drug Content

Product Reference	Balloon Diameter (mm)	Balloon Length (mm)	Drug Content (±10%) (µg)
CME15010	1.50	10	60
CME20010	2.00	10	80
CME22510	2.25	10	90
CME25010	2.50	10	100
CME27510	2.75	10	110
CME30010	3.00	10	120
CME35010	3.50	10	140
CME40010	4.00	10	160
CME15015	1.50	15	90
CME20015	2.00	15	120
CME22515	2.25	15	135
CME25015	2.50	15	150
CME27515	2.75	15	164
CME30015	3.00	15	179
CME35015	3.50	15	209
CME40015	4.00	15	239
CME15020	1.50	20	120
CME20020	2.00	20	160
CME22520	2.25	20	179
CME25020	2.50	20	199
CME27520	2.75	20	219
CME30020	3.00	20	239
CME35020	3.50	20	279
CME40020	4.00	20	319
CME15025	1.50	25	150
CME20025	2.00	25	199
CME22525	2.25	25	224
CME25025	2.50	25	249
CME27525	2.75	25	274
CME30025	3.00	25	299
CME35025	3.50	25	349
CME40025	4.00	25	399
CME15030	1.50	30	179
CME20030	2.00	30	239

CME22530	2.25	30	269
CME25030	2.50	30	299
CME27530	2.75	30	329
CME30030	3.00	30	359
CME35030	3.50	30	419
CME40030	4.00	30	479
CME15035	1.50	35	209
CME20035	2.00	35	279
CME22535	2.25	35	314
CME25035	2.50	35	349
CME27535	2.75	35	384
CME30035	3.00	35	419
CME35035	3.50	35	489
CME40035	4.00	35	558
CME15040	1.50	40	239
CME20040	2.00	40	319
CME22540	2.25	40	359
CME25040	2.50	40	399
CME27540	2.75	40	439
CME30040	3.00	40	479
CME35040	3.50	40	558
CME40040	4.00	40	638

2.0 INDICATIONS

The Magic Touch - ED Sirolimus Coated Balloon Catheter is indicated for the treatment of erectile dysfunction for male patients with de novo and restenotic lesions in the internal pudendal, dorsal, or urethral arteries that are 1.50 mm to 4.00 mm in diameter and up to 36 mm lesion length.

3.0 CONTRAINDICATIONS

Use of Magic Touch - ED Sirolimus Coated Balloon Catheter is contraindicated in the following types of patients:

- Patients with a hypersensitivity to sirolimus or its structural activity related compounds.
- Patients with a known hypersensitivity to excipients with phospholipid or related origins.
- Patients whose diseased segment cannot be pre-dilated or prepared before drug coated balloon treatment.
- Lesions not suitable for interventional treatment.
- Patients who cannot receive recommended antiplatelet or anticoagulation therapy.

4.0 WARNINGS

- Contents supplied are sterilized using Ethylene Oxide (EO) process.
- Ensure that the product package is intact, and the inner package housing device has not been opened or damaged as this may indicate the sterile barrier has been broken.
- **Do not** use the product if sterile barrier is opened prior to intended use.
- Inspect the product before use for any damage during removal from the protective catheter holder and protective cap. Return to company using proper handling precautions if found to be damaged for further inspection and analysis.
- The use of product carries the risks associated with pudendal angioplasty, including vascular complications or bleeding events.
- The product shall be discarded in a proper way after use as it may have sirolimus on its surface and shall be discarded as per applicable regulations.
- Magic Touch - ED is intended for single use only.
- **Do not** re-sterilize and re-use. Reuse or re-sterilization may create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Reuse or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness and death. Concept Medical is not be responsible for any direct, incidental or consequential damages resulting from re-sterilization or reuse.
- **Do not** exceed rated burst pressure (RBP). Use of correct pressure monitoring device is recommended to prevent over pressurization. Inflation in excess of the RBP may cause the balloon to rupture.
- **Do not** use after the " Use-by Date".
- When the Magic Touch - ED is in the body, it should be manipulated only under sufficient and/or high-quality fluoroscopy. Prior to withdrawing the catheter from the lesion, the balloon must be completely deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- **Do not** use air or any gaseous medium to inflate the Magic Touch - ED. Use only the recommended inflation medium. Handle the product with caution to avoid any damage to the balloon coating or folded balloon. Avoid touching the balloon while loading it on guidewire. Always use dry Sterile gloves while loading Magic Touch - ED on Guidewire.
- Take appropriate precautions to prevent or reduce clotting when catheter is used. Flush or rinse ancillary product, if needed (except Magic Touch - ED Balloon) entering the vascular system with sterile isotonic saline or a similar solution via the guide wire access port prior to use. Consider the use of systemic heparinization.

5.0 PRECAUTIONS FOR USE

5.1 GENERAL PRECAUTIONS:

- Only physicians who have received adequate training in Interventional Cardiology, Vascular Surgeon, Interventional Radiologist and Interventional Urologist shall perform this treatment.
- The treatment shall be performed at hospitals where emergency surgical set-up is available.
- Do not expose system to organic solvents such as alcohol and detergents as it may impair with product performance in treatment.
- **Pre-dilatation of the lesion before the treatment with Magic Touch - ED is mandatory.**
- To minimize the possible introduction of air into the system, it's imperative that, before proceeding, careful attention is paid to the maintenance of tight catheter connections and through aspiration and flushing of the system.
- Never advance the Drug delivery catheter without the guide wire extending from the tip.
- Never attempt to move the guide wire when the balloon is inflated. Do not advance the catheter against significant resistance. The cause of resistance should be determined via fluoroscopy and appropriate remedial action taken.
- Magic Touch - ED is not intended for stent expansion.
- Magic Touch - ED is not intended for pre-dilation of the lesion.
- Use of Magic Touch - ED in patients and lesions outside intended use (off label)is not recommended.

 Note: The use of Magic Touch - ED for treatment of pudendal artery disease is a relatively new concept to re-open blocked or clogged pudendal (penile) artery which re-establishes blood flow. Drug release using Magic Touch - ED provides longer drug retention and bio-available and thereby prevent re-occurrence.

5.2 POST PROCEDURE PRECAUTIONS:

If a surgical procedure is anticipated that requires early discontinuation of antiplatelet therapy, the interventional surgical specialist and patient should carefully consider whether a drug coated balloon and its associated recommended antiplatelet therapy is the appropriate treatment choice.

5.3 MULTIPLE INFLATION OF DEVICE

- Multiple inflation of a Magic Touch - ED has not been evaluated in humans.
- **It is recommended that treatment using single inflation of device at least 60 second at nominal pressure.**

5.4 DRUG INTERACTION

Several drugs are known to affect metabolism of Sirolimus, and other drug interactions may be inferred from known metabolic effects. Sirolimus is known to be a substrate for both cytochrome P450 IIIA4 (CYP3A4) and P-glycoprotein as per published papers. Consideration should be given to the potential drug interaction when deciding for Magic Touch - ED Sirolimus Coated Balloon Catheter treatment in a patient who is taking a drug that could interact with Sirolimus. The effect of drug interactions on safety and effectiveness of Magic Touch - ED Sirolimus Coated Balloon Catheter has not been determined.

5.5 IMMUNE SUPPRESION POTENTIAL

Sirolimus, the active ingredient of Magic Touch - ED Sirolimus Coated Balloon Catheter, is an immunosuppressive agent that

is also available in oral formulations. The total drug concentration in blood from two inflations of Magic Touch in rabbit iliac artery, Sirolimus Coated Balloon Catheter at 24 hours was found to be 0.81ng/ml and is substantially lower than the therapeutic concentrations obtained when Sirolimus oral formulations are used as prophylaxis for renal transplant rejection (Pharmacokinetic and Histology study done at CV Path Institute, Washington DC USA Lemos et al, *Eurointervention*;9(1), 148-156, 2013).

However, when multiple Magic Touch - ED Sirolimus Coated Balloon Catheter treatment is performed in patients, systemic concentration of Sirolimus may reach immunosuppressive levels momentarily, especially in patients who also have hepatic insufficiency or who are taking drugs that inhibit CYP3A4 or P-glycoprotein. This possibility should be considered for such patients, particularly if they are also taking oral Sirolimus, other immunosuppressive drugs or are otherwise at risk of immunosuppression.

5.6 HANDLING PRECAUTIONS

- Special care must be taken not to handle or in any way disrupt the coating on balloon. This is most important while removing the catheter from the packaging, placing it over the guide wire and advancing it through large bore rotating haemostatic valve and guiding catheter hub.
- Balloon manipulation, trying to refold may damage the coating, contamination or fragmentation of coating on the balloon.

5.7 MAGIC TOUCH - ED SIROLIMUS COATED BALLOON CATHETER INFLATION PRECAUTIONS

- Magic Touch - ED is a drug delivery device to deliver drug upon appropriate inflation to the vessel.
- The vessel shall be pre-dilated with an appropriately sized balloon.
- Do not prepare or pre-inflate the balloon prior to treatment other than as directed. Use the balloon purging technique.
- Guiding catheters should have lumen sizes that are suitable to accommodate the balloon across the lesion. Wrong sizing may lead to damage of balloon coating.
- Inflating a balloon in a tight lesion may lead to dissection of the vessel distal/proximal to the lesion and may cause acute closure of the vessel leading to additional intervention (CABG, further dilatation, placement of additional stents or other intervention).
- Do not expand the balloon, if it is not properly positioned at the target lesion.
- Expansion of balloon has a potential to compromise side branch flow, hence caution needs to be taken when inflation is done.
- Inflation pressure of balloon should always be monitored. Do not exceed Rated burst pressure as listed in the compliance chart on label. Use of pressures higher than those specified on the product label may result in ruptured balloon with possible intima damage and dissection.
- Do not attempt to pull a partially inflated or inflated balloon back through the guiding catheter, as serious injury to intima layer may occur. Remove the device after complete deflation through the catheter. Use proper retrieval methods for the device as it may cause trauma at the insertion site and may lead to bleeding, haematoma or pseudo aneurysm.
- Ensure treatment device length more than or equal to lesion length with addition of 2 mm on proximal and distal end for complete coverage.

5.8 MAGIC TOUCH - ED SIROLIMUS COATED BALLOON CATHETER REMOVAL PRECAUTIONS

- Should any unusual resistance be felt at any time during either lesion access or removal of the balloon, the entire system should be removed.
- While removing the system, advance the guidewire into the pudendal/penile anatomy as far as distally possible. Tighten the rotating haemostasis valve to secure the Magic Touch - ED Sirolimus Coated Balloon Catheter with the guiding catheter; then removing the Balloon catheter system and guiding catheter as a single unit.
- Failure to follow this procedure or to apply excessive force during withdrawal can potentially result in damage to the anatomy and cause any other vascular complication.
- If it is necessary, retain the guidewire in position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.

 Note: After use, this product may be a biohazard. Handle and dispose of all such devices in accordance with accepted medical practice and applicable government regulations.

5.9 POST PROCEDURE PRECAUTIONS

- Care must be taken when the Magic Touch - ED Sirolimus Coated Balloon Catheter is being withdrawn after treatment. Post procedure treatment of Dual antiplatelet therapy (DAPT) is recommended for a period of minimum 3 months. Patients who require early discontinuation of antiplatelet therapy (e.g. secondary to active bleeding) should be monitored carefully.

6.0 DRUG INFORMATION

6.1 MECHANISM OF ACTION

The mechanism by which Magic Touch - ED Sirolimus Coated Balloon Catheter affects neo-intimal growth and arterial characteristics as observed in pre-clinical studies has not been established. Sirolimus is a potent inhibitor of T-Lymphocyte activation, Smooth Muscle Cell and Endothelial cell proliferation in response to stimulation by cytokines, growth factors and antigens. In cells, Sirolimus binds to the immunophillin, FK Binding protein-12 (FKBP-12). Sirolimus-FKBP-12 complex binds to and inhibits the activation of the Mammalian Target of Rapamycin (mTOR), a key regulatory kinase. This inhibition suppresses cytokine-driven cell proliferation, inhibiting the progress of the cell cycle from the G_1 to S phase.

6.2 PHARMACOKINETIC OF MAGIC TOUCH - ED SIROLIMUS COATED BALLOON CATHETER

Twenty-five healthy Yucatan mini swine were included. Balloon treatments with Magic Touch - SCB were performed in coronary arteries. Blood was collected from 10 animals in several time points (0 hour [h], 0.05h, 0.25h, 1h, 3h, 24h, 30 days [d] and 60d after SCB treatment) for the assessment of sirolimus concentration in whole blood. 3 animals were sacrificed at each time point (1h, 24h, 3d, 14d, 30d, and 60d) to determine sirolimus concentration in coronary arteries of treated and untreated proximal/distal segments. Downstream myocardium of SCB treated segments as well as non-target peripheral organs were also collected for the assessment. The tissues were dissected, carefully excised, and immediately frozen in liquid nitrogen for further analysis. Measurements of sirolimus blood level and tissue concentration were performed by liquid chromatography-tandem mass spectrometry (LC/MS/MS).


Figure 3A EDTA whole blood concentration


Figure 3B In-tissue sirolimus concentration

Sirolimus was quantifiable in EDTA whole blood as early as 3 minutes after SCB deployment. The maximum concentration (Cmax) in whole blood of 17.9 ± 8.6 ng/mL was recorded after 15 minutes and progressively declined thereafter and was barely detectable at 24 hours. In the treated coronary artery segments, the Cmax of Sirolimus was measured at 1 hour after SCB deployment. Two weeks after treatment, the Sirolimus concentration had dropped to less than one tenth of Cmax. The Sirolimus concentrations in the untreated proximal and distal coronary artery segments were markedly lower than those in the treated artery tissues (Figure 3B, Table 3B). In all downstream myocardial tissues as well as non-target peripheral organs, the highest Cmax of Sirolimus was measured at 1 hour after SCB deployment. At 30 days, no Sirolimus could be quantified in any of tissues (Table 3B).

Table 3A; Sirolimus concentration in EDTA whole Blood

Time (h)	0h	0.05h	0.25h	1h	3h	24h	30d	60d*
Sirolimus concentration in blood (n=7 in each)	BLOQ	11.2 ± 3.25	17.9 ± 8.59	14.7 ± 8.06	9.1 ± 5.62	0.85 ± 0.49	BLOQ	BLOQ

Values are described as mean ± standard deviation (ng/ml)
Abbreviations; BLOQ=below the lower limit of quantification (<0.1ng/ml), d=days, h=hour.
*60d time point had n=3. Other time points had n=7

Table 3B; Sirolimus concentration in vessels, downstream myocardium, and non-target peripheral organs

Time Points		1h	24h	3d	14d	30d	60d
Coronary artery (n=6)	Proximal untreated segment	192.7 ± 115.4	160.1 ± 66.2	96.3 ± 55.8	74.2 ± 84.7	30.9 ± 25.6	0.8 ± 1.9
	Treated segment	1959.4 ± 641.0	1459.9 ± 518.1	410 ± 118.6	168.1 ± 76.9	47.7 ± 38.3	10.3 ± 5.2
	Distal untreated segment	1107.1 ± 650.8	399.0 ± 212.1	113.3 ± 70.2	172.4 ± 166.6	3.2 ± 7.9	BLOQ
Downstream myocardium (n=6)		1698.1 ± 447.1	338.9 ± 349.8	25.0 ± 12.2	0.5 ± 1.3	BLOQ	BLOQ
Kidney (n=3)		56.2 ± 3.4	18.0 ± 4.4	8.3 ± 1.6	BLOQ	BLOQ	BLOQ
Liver (n=3)		61.6 ± 15.5	7.2 ± 1.5	0.9 ± 1.3	BLOQ	BLOQ	BLOQ
Lung (n=3)		284.1 ± 48.1	36.5 ± 20.5	26.6 ± 16.4	4.6 ± 3.4	BLOQ	BLOQ
Bone/rib (n=3)		3.1 ± 1.0	3.7 ± 1.1	0.9 ± 1.2	BLOQ	BLOQ	BLOQ
Spleen (n=3)		97.5 ± 19.8	12.2 ± 3.0	7.2 ± 1.6	BLOQ	BLOQ	BLOQ

Values are described as mean ± standard deviation (ng/g)
Abbreviations; BLOQ=below the lower limit of quantification (<0.25ng/ml tissue homogenate), d=days, h=hour.

6.3 DRUG INTERACTION FOLLOWING TREATMENT OF MAGIC TOUCH - ED

Drug interaction studies have not been performed with Magic Touch - ED Sirolimus due to limited systemic exposure to Sirolimus eluted from Magic Touch - ED Balloon. However, consideration should be given to the potential for both local and systemic drug interactions in the vessel wall when deciding to use Magic Touch - ED in a subject taking a drug with the known interaction with Sirolimus.

Sirolimus is metabolized by cytochrome P450 3A4 (CYP3A4) in the gut wall and liver and undergoes efflux from enterocytes of the small intestine by p-glycoprotein (P-gp). Absorption and subsequent elimination of systemically absorbed Sirolimus may be influenced by drugs that affect CYP3A4 and P-gp.

Medications that are strong inhibitors of CYP3A4 and P-gp may increase Sirolimus levels while inducers of CYP3A4 and P-gp may reduce Sirolimus metabolism in-vivo.

Sirolimus when prescribed as oral medication may interact with the drugs/foods listed below:

- P-gp Inhibitors such as but not limited to Cyclosporine, Digoxin
- Antibiotics such as but not limited to Ciprofloxacin, Ofloxacin
- Glucocorticoids
- CYP3A4 Inhibitors such as but not limited to Ketoconazole, Fluconazole
- CYP3A4 Inducers such as but not limited to Rifampin, Phenytoin, Dexamethasone
- Herbal preparations (St. John's wort- Hypericum perforatum)
- Amphotericin B such as but not limited to Abelcet, Amphocin
- Cimetidine such as but not limited to Tagamet
- Vitamins
- Vaccination- sirolimus is potent immunosuppressant hence it may affect response to vaccination making the vaccination less effective. For some period after receiving Magic Touch - ED, Sirolimus Coated Balloon Catheter use of live vaccines such as, but not limited to measles, mumps, BCG, yellow fever, varicella, Ty21a Typhoid should be avoided.
- Grapefruit/grapefruit juice: Grapefruit/grapefruit juice reduces CYP3A4 mediated metabolism of sirolimus.
- Medicine to lower your cholesterol or triglycerides
- Medicine for high blood pressure or heart problems
- Anti-seizure medicine
- Medicines used to treat stomach acid, ulcers, or other gastrointestinal problems
- Drugs that may increase sirolimus blood concentration:
 - •Antifungal agents: such as but not limited to Itraconazole, Clotrimazole
 - •Gastrointestinal prokinetic agents such as but not limited to metoclopramide
 - •Macrolide antibiotics: such as but not limited toTroleandomycin
 - •Calcium channel blockers: such as but not limited to Diltiazem, verapamil
 - •Other drugs: Danazole, HIV Protease Inhibitor (Indinavir)
- Drugs that may decrease sirolimus level:
 - •Anticonvulsants such as but not limited to Phenobarbitol
 - •Antibiotics such as but not limited to Rifabutin

6.4 CARCINOGENICITY, GENOTOXICITY AND REPRODUCTIVE TOXICOLOGY

Carcinogenicity, Genotoxicity and Reproductive Toxicology

The carcinogenicity, and reproductive toxicity of Magic Touch - ED have not been evaluated. However, carcinogenicity, and reproductive toxicity of sirolimus has been investigated in bacterial and mammalian cells in-vitro and in laboratory animals in vivo.

Carcinogenicity

Carcinogenicity studies were conducted in mice and rats. In 86-week female mouse study at sirolimus doses 30 to 120 times higher than the 2 mg daily clinical dose (adjusted for body surface area), there was a statistically significant increase in malignant lymphoma at all dose levels compared with controls.

In an another mouse study at dosages that were approximately 3 to 16 times the clinical dose (adjusted for body surface area), hepatocellular adenoma and carcinoma in males were considered sirolimus-related.

In the 104-week rat study at dosages equal to or lower than the clinical dose of 2 mg daily (adjusted for body surface area), there were no significant findings.

Reproductive Toxicity

The study of the effects of sirolimus on the different stages of the reproductive process were conducted in rats. Fertility was diminished slightly in both male and female rats following oral administration of sirolimus at doses approximately 10 times or 2 times, respectively, the clinical dose of 2 mg daily (adjusted for body surface area). Male rats showed decreased fertility and testicular atrophy (giant cells and hypospermia in the testes and epididymides). In female rats, reduced size of ovaries and uteri was observed.

Genotoxicity

The Genotoxicity test has been performed as per ISO 10993-3: 2014 for Sirolimus Coated Balloon Catheter. The results show that the product doesn`t induce any micronucleus in the bone marrow of mice, doesn`t induce chromosome aberrations in cultured human peripheral blood lymphocytes and it is non-mutagenic.

7.0 POTENTIAL ADVERSE EVENTS

Possible complications linked to the use of the balloon catheter during the procedure:

- Dissection or perforation of the pudendal artery
- Injury or rupture of the pudendal artery
- Total occlusion
- Arterial spasm
- Ventricular fibrillation
- Embolism

- Infection and or pain at the access site
- Haematoma
- Drug reaction to antiplatelet agent / anticoagulation agent / contrast medium
- Distal emboli,
- Nausea & vomiting
- Vascular complication which may requires vessel repair.
- Alopecia
- Anaemia
- Blood product transfusion
- Gastrointestinal Symptoms
- Hepatic Enzyme Changes
- Histological changes in vessel wall, including inflation, cellular damage or necrosis.
- Myalgia/Arthralgia
- Peripheral Neuropathy

Possible complications that could occur following an angioplasty procedure with balloon catheter, on short and medium term:

- Restenosis of the dilated artery
- Bleeding complications or haematoma

8.0 PUBLISHED CLINICAL DATA OF SIROLIMUS COATED BALLOON IN VASCULOGENIC ERECTILE DYSFUNCTION

A pilot clinical study was initiation (Principal Investigator Prof. Giuseppe Sangiorgi, MD FESC, FACC, FSCAI, San Gaudenzio Clinic, Italy) to evaluate the performance of Magic Touch - ED in in focal atherosclerotic lesions of the internal pudendal and dorsalis penis arteries among men with erectile dysfunction (ED) not responding to phosphodiesterase-5 inhibitors.

8.1 STUDYDESIGN

The study is observational, prospective, multicenter registry evaluating the use of Magic Touch - ED after pre-dilatation with standard balloon angioplasty for the treatment of up 2 de novo stenosis in native internal pudendal arteries.

8.2 STUDY ENDPOINTS

- Incidence of MAE (myocardial infarction, dysrhythmia, congestive heart failure, leg ischemia, renal insufficiency, bowel complication, reoperation, PTA site infection, stroke, respiratory complication, and no home discharge).
- Binary restenosis at 8 months follow-up CT-angiography.
- Delta IEF >5 from baseline to 8 mos FU

8.3 STUDYRESULTS

Table 4: Baseline characteristics

Patient's baseline characteristics, n=14	**n (%)**
Age, (Years ± SD)	59.8±10.8
BMI (kg/m^2± SD)	23.1±2.8
Active Smoker	8 (57)
Diabetes NID	5 (36)
Diabetes ID	4 (28)
Hypertension	9 (64)
Hypercholesterolemia	6 (42)
Family history for CAD	3 (21)
Coronary artery disease	4 (35)
Peripheral artery disease	3 (21)
History of prostate surgery	2 (14)
Chronic prostatitis	1 (7)
History of or current dialysis	0
Baseline IIEF-5 score	9.7±6.6
Q3: ability to achieve penetration	1.69±1.27
Q4: ability to maintain erection	1.18±1.23
Associated Venous leakage	1 (7)
Baseline medications tried before procedure	**n (%)**
Phosphodiesterase type 5 inhibitor	7 (50)
Prostaglandin	6 (42)
Intracavernosal Testosterone	3 (21)
Medications with an impact on EF	**n (%)**
β-Blockers	3 (21)
Psychotropic drugs	0
Thiazide diuretics	2 (14)

Table 5: Lesion and procedural characteristics

Target lesion (n=18)	**n (%)**
Internal iliac artery	3 (16)
Internal pudendal artery	11 (62)
Common penile artery	2 (12)
Dorsal penile artery	1 (5)
Cavernosal artery	0
Inferior gluteal artery	1 (5)
Arteries affected	**n (%)**
Left side only	5 (27)
Right side only	6 (33)
Bilateral	7 (40)
Arteries treated	**n (%)**
Left side only	9 (50)
Right side only	7 (38)
Bilateral	2 (12)
Lesion length, mm	11.9±6.6
Total lesion length, mm	17.6±13.3
RVD, mm	2.2±0.5
Diameter stenosis, %	68.0±6.5
MLD, mm	1.4±0.6
PSVb, cm/s	19.5±10.9
EDVb, cm/s	8.9±5.4
Endovascular intervention	**n (%)**
Standard balloon dilation	18 (100)
Drug-eluting stent	3 (16)
Sirolimus-coated balloon	15 (84)
No access	1 (5)
Procedural Success (stenosis <10%)	17 (94)

Treatment Success Reflected by Minimum Clinically Relevant Improvement in the International Index of Erectile Function-5 (IIEF-5) score


Figure 4: International Index of Erectile Function–5 (IIEF-5) score graph

Study results demonstrated that Sirolimus coated balloon catheter performs as intended in vasculogenic erectile dysfunction and is proven to be safe and effective treatment option.

Data source: Presented in TCT 2019. https://www.tctmd.com/slide/sirolimus-coated-balloon-treatment-vasculogenic-erectile-dysfunction-ready-prime-time***)***

9.0 INDIVUDUALIZATION OF TREATMENT

- Selection of patients for Magic Touch - ED Sirolimus Coated Balloon Catheter treatment should be carefully done looking at the risks and benefits to the patient.
- Multivariable patient conditions and outcomes modelling suggest that treatment assignment/strategy remained an independent predictor of clinical and angiographic outcomes even after adjusting for other baseline and procedural confounding variables.

10.0 PATIENT COUNSELLING INFORMATION

Physicians should consider the following in counselling patient about this product:

- Discuss the risks associated with Magic Touch - ED balloon dilation treatment
- Discuss the risks/benefits issues for this particular patient
- Discuss alteration to current lifestyle immediately following the procedure and over the long term.

11.0 HOW SUPPLIED

- STERILE: This device is sterilized with ethylene oxide gas and is non-pyrogenic. Do not use if the package is opened or damaged. Device is for one use only (single use). Do not re-sterilize.

- CONTENTS: One (1) Magic Touch - ED Sirolimus Coated Balloon Catheter (Rapid exchange delivery system).
- STORAGE: Temperature limitation: 8°C to 25°C, store at dry place. Keep away from direct sunlight.
- Disposal: Dispose of device in accordance with local regulations.

12.0 OPERATOR'S MANUAL FOR USE

Magic Touch - ED is supplied and packaged in sterile condition. The product is packed in a box with identification labels on each packing. There are two stages of packing for maintaining the integrity and sterility of device - first pack is a Tyvek pouch which is sealed and stored in a secondary laminated aluminium foil pouch for protection. It is important that integrity of packing is verified before usage of device.

12.1 ACCESS TO STERILE PACK

Open the product box pack at tear line provided at the top. This provides access to the secondary aluminium foil pouch in which the sterile package is kept. Tear open the laminated aluminium foil pouch at the cut provided on the side. The inner most Tyvek pouch holds sterile Magic Touch - ED. The Magic Touch - ED is stored in a coiled hoop. Open the Tyvek pouch and pass or drop the sterile balloon catheter into the sterile field using aseptic technique.

12.2 INSPECTION PRIOR TO USE

Before opening, carefully inspect the balloon system package, and check for damage to sterile barrier. Prior to using the device, carefully remove the system from package and inspect it for any bends, kinks or other damage. Do not use the device if any damage to the sterile package or balloon system is noted.

12.3 MATERIAL REQUIRED

Quantity	Material
N/A	Appropriate guiding catheter(s)
2-3	10-20 cc syringes
1,000 u /500 cc	Sterile Heparinized Normal Saline (HepNS)
1	0.014" (0.36 mm) diameter guidewire
1	Rotating haemostatic valve with an appropriate internal diameter (min. I.D. of 0.096" [2.4 mm])
N/A	Contrast diluted 1:1 with normal saline
1	Inflation device
1	Stopcock (3-way minimum)
1	Torque device
1	Guidewire Introducer
N/A	Appropriate anticoagulation and anti-platelet drugs

12.4 PREPARATION OF MAGIC TOUCH - ED

AVOID any manipulation of the Magic Touch - ED during preparation as this may disrupt the integrity of coating on the device.

DO NOT wipe the distal part having balloon using any media as it may damage the coating.

ANY advancement of the drug coated balloon catheter should be in deflated condition and with use of guidewire in the penile artery.

MAGIC TOUCH –ED PREPARATION

Steps

1. Prepare the inflation device with diluted contrast medium.
2. Attach the inflation device to the 3-way stopcock.
3. Attach stopcock to the balloon inflation hub.
4. Open the stopcock to balloon catheter system and leave inflation device system on neutral position.

12.5 MAGIC TOUCH - ED DELIVERY PROCEDURE

Steps

1. Prepare the vascular access site according to standard practice. Required access to the lesion has to be done using an appropriate guide wire (recommended 0.014").
2. Pre-dilate the lesion with a pre-dilatation balloon. Limit the longitudinal length of pre-dilatation by the pre-dilatation balloon to avoid creating a region of vessel injury that is outside the boundaries of the Magic Touch - ED.
3. Maintain neutral pressure on the inflation device. Open the rotating haemostatic valve as widely as possible.
4. Backload the delivery system onto the proximal portion of the guidewire while maintaining the guidewire position across the target lesion.
5. Advance the drug coated balloon over the guidewire to the target lesion. Use the radiopaque balloon markers to position the balloon across the lesion; perform check angiography shoot to confirm the position of the balloon. Close the rotating haemostatic valve at the right position of balloon catheter in lesion to lock it in place.

 Note: Should unusual resistance be felt at any time during either lesion access or removal of the drug coated balloon before inflation, the entire system should be removed as a single unit. See Precautions - 5.8 Removal precautions for removing balloon catheter system.

12.6 MAGIC TOUCH - ED DEPLOYMENT PROCEDURE

Steps

1. Before deployment, reconfirm the correct position of the drug coated balloon relative to the target lesion via the radiopaque balloon markers.
2. Attach the inflation device (only partially filled with contrast media) to a three-way stopcock and apply negative pressure to purge the balloon of air.
3. Turn the stopcock on the catheter to the off position and purge the inflation device of air. Close the side port of the stopcock.
4. Under fluoroscopic visualization, inflate the balloon to at least the nominal pressure to deliver drug from the balloon surface, but do not exceed the labelled rated burst pressure. Optimal treatment and drug delivery require the balloon to be in full contact with the artery wall, with the internal diameter matching the size of the reference vessel diameter. Higher pressures may be required for certain lesions but Rated burst pressure should not be exceeded at all times. Refer to balloon compliance chart for further details on diameters and sizes.
5. Fully cover the entire lesion and dilated balloon area (including dissections) with the Magic Touch - ED Sirolimus Coated Balloon Catheter, allowing for adequate coverage into healthy tissue proximal and distal to the lesion.
6. Keep Magic Touch - ED inflated for at least 60 seconds.
7. Deflate the balloon by pulling a vacuum with the inflation device. Make sure that the balloon is fully deflated before attempting to remove.
8. Confirm that vessel dilation is adequately done by angiographic injection through the guiding catheter.

 NOTE: It is recommended that the guidewire and/or the balloon catheter remain across the lesion until the procedure is complete. Contrast media have different viscosities and may affect the inflation/deflation time.

12.7 MAGIC TOUCH - ED REMOVAL PROCEDURE

Steps

1. Ensure that the balloon is fully deflated.
2. While maintaining the guidewire position and negative pressure on the inflation device, withdraw the drug coated balloon delivery system.

 Note: Should unusual resistance be felt at any time during either lesion access or removal of the drug coated balloon catheter before inflation, the entire system should be removed as a single unit. See Precautions – 5.8 Removal precautions for removing balloon catheter system.
3. Repeat angiography to assess the treatment area.

12.8 MAGIC TOUCH - ED COMPLIANCE CHART

Refer to labelling.

13.0 STERILISATION AND STORAGE CONDITIONS

Magic Touch - ED has been sterilized using Ethylene Oxide gas (EtO) and is non-pyrogenic. The device is for single use only. Storage temperature limitation: 8°C to 25°C. Store at dry place. Keep away from direct sunlight.

14.0 PATIENT INFORMATION

In addition to this Instructions for Use booklet, following patient specific information regarding the Magic Touch - ED is available in product box:

A Patient Treatment Card includes specific information about Magic Touch - ED (Product name, Size, Reference Number, Serial Number, LOT Number and Expiry Date) and contact details for patient and his treating physician. It is always recommended that patients treated with Magic Touch - ED keep this card in their possession and present it to their treating physician before going for any other treatment / surgical intervention.

15.0 DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

Concept Medical B.V. does not provide any expressed or implied warranty, including without limitation, any implied warranty of merchantability, fitness for a particular purpose or non-infringement, with respect to the product(s) described in this publication and Concept Medical B.V. hereby expressly disclaims any such warranties.

Descriptions or specifications in Concept Medical B.V. printed

matter, including this publication, are meant solely to generally describe the relevant product at the time of manufacture and do not constitute any express warranties.

No person has the authority to bind Concept Medical B.V. to any representation or warranty.

Under no circumstances shall Concept Medical B.V. be liable for any medical or other expenses or direct, indirect, incidental, consequential, punitive or exemplary damages with respect to the products described herein, including without limitation, in respect of any use, defect, failure or malfunction, whether based upon warranty, contract, tort or other legal theory, even if Concept Medical B.V. has been advised or should have been aware of the possibility of such damages, other than as expressly provided by specific law without possibility of exclusion.

The exclusions, disclaimers and limitations contained in this section shall be interpreted in such manner as to be effective and valid under applicable law, but in case any one or more of the provisions contained herein shall, for any reason, be held to be invalid, illegal or unenforceable in any respect, such provision shall be ineffective to the extent, but only to the extent, of such invalidity, illegality or unenforceability without invalidating the remainder of this section or any other provisions hereof.

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Žiadna osoba nemá právomoc zaviazať spoločnosť Concept Medical B.V. k akémukoľvek vyhláseniu alebo záruke.