Medacta International

Strada Regina
CH - 6874 Castel San Pietro - Switzerland
Phone +41 91 696 60 60 - Fax +41 91 696 60 66
info@medacta.ch – www.medacta.com

Ref: 75.09.039 rev. 06

CE 0476
2016

Manufactured by:
Medacta International SA
CH-6874 Castel San Pietro - Switzerland

MEDACTA SHOULDER SYSTEM

CAUTION

For a U.S. specific package insert, please go to www.medacta.com or call toll free 800-901-7836.

SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES / ΕΙΚΟΝΟΣΎΜΒΟΛΑ

	Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser / Μην επαναχρησιμοποιείτε
	Do not resterilize / Non risterilizzare / Niet resteriliseren / Bitte nicht re-sterilisieren / No resterilize / Não resterilizar / Ne pas restériliser / Να μην αποστειρωθεί ξανά
	Caution, read the accompanying documents / Attenzione, consultare i documenti di accompagnamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atenção, consultar os documentos de acompanhamento / Attention, consulter les documents d'accompagnement / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
	Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las Instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation / Οδηγίες χρήσης
	Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil / Μην εκθέτετε σε ηλιακό φως
	Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec / Φυλάσσετε σε ξηρό χώρο
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilice si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abîmé / Να μην χρησιμοποιηθεί εάν η συσκευασία είναι κατεστραμένη
	Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant / Ημερομηνία λήξης
	Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote / Numéro de lot / Αριθμός παρτίδας
	Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale / Νούμερο παραπομπής
	Sterilized with ethylene oxide / Sterilizzato con ossido di etilene / Gesteriliseerd met ethyleenoxide / Ethylenoxid-sterilisiert / Esterilizado con óxido de etileno / Esterilizado com óxido de etileno / Stérilisé à l'oxyde d'éthylène / Αποστειρωμένο με οξείδιο του αιθυλενίου
	Sterilized by irradiation / Sterilizzato mediante irradiazione / Gesteriliseerd door doorstraling / Sterilisiert durch Bestrahlung / Esterilizado por irradiación / Esterilizado por irradiação / Stérilisé par irradiation / Αποστειρωμένο με ακτινοβολία
	MR Unsafe / Non sicuro per RM / MR onveilig / MR ungeeignet / Inseguro RM / Não seguro para RM / MR non sécurisé / Μη ασφαλές στην μαγνητική τομογραφία
	Manufacturer / Produttore / Fabrikant / Hersteller / Fabricante / Fabricante / Fabricant / Κατασκευαστής
	Double sterile barrier system /Sistema a Doppia Barriera Sterile/ Dubbele steriele verpakking/ Doppeltes steriles Barrieresystem/ Sistema de doble barrera estéril/ Sistema de Dupla Barreira Estéril/ Système à barrière stérile double/ Σύστημα διπλού στείρου εμποδίου

Index

ENGLISH: MEDACTA SHOULDER SYSTEM - INSTRUCTIONS FOR USE 3
ITALIANO: MEDACTA SHOULDER SYSTEM - ISTRUZIONI PER L'USO 7
NEDERLANDS: MEDACTA SHOULDER SYSTEM - GEBRUIKSINSTRUCTIES 12
DEUTSCH: MEDACTA SHOULDER SYSTEM - GEBRAUCHSANWEISUNGEN 17
ESPAÑOL: MEDACTA SHOULDER SYSTEM - INSTRUCCIONES DE USO 22
PORTUGUÊS: MEDACTA SHOULDER SYSTEM - INSTRUÇÕES DE UTILIZAÇÃO 27
FRANÇAIS : MEDACTA SHOULDER SYSTEM - INSTRUCTIONS D'UTILISATION 32
ΕΛΛΗΝΙΚΑ: MEDACTA SHOULDER SYSTEM - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ 37

ENGLISH: MEDACTA SHOULDER SYSTEM - INSTRUCTIONS FOR USE

Important notice: the device(s) can be prescribed and implanted only by a doctor legally authorized to perform this type of surgery.

1. GENERAL

Before any surgery, the surgeon must be familiar with the sales product literature and operative technique and must carefully read these instructions for use. Patient selection is as important as implant placement or positioning. Overweight patients or unsuitable functional requirements may generate exceptional stresses and reduce the implant life. The warnings must be heeded, and the instructions for use must be strictly followed

2. PRODUCT DESCRIPTION

These instructions for use are intended for all the products described here below.

Medacta®'s product description - Anatomic

Medacta® shoulder products used in anatomic configuration are:

- Humeral diaphysis (cementless and cemented) to be used in hemi or total shoulder arthroplasty, for primary or revision surgery
- Stemless humeral metaphysis (cementless) to be used in hemi or total shoulder arthroplasty, for primary surgery
- Humeral anatomical metaphysis (cementless and cemented) to be used in hemi or total shoulder arthroplasty, for primary or revision surgery
- Humeral stem screw to be used in hemi or total shoulder arthroplasty, for primary or revision surgery
- Double eccenter to be used in hemi or total shoulder arthroplasty, for primary or revision surgery
- Stemless double eccenter to be used in hemi or total shoulder arthroplasty, for primary surgery
- Metal humeral head to be used in hemi or total shoulder arthroplasty, for primary or revision surgery
- Double eccenter/reverse metaphysis screw to be used in hemi or total shoulder arthroplasty, for primary or revision surgery
- Stemless double eccenter screw to be used in hemi or total shoulder arthroplasty, for primary surgery
- Cemented glenoid component to be used in total shoulder arthroplasty, for primary or revision surgery

Medacta®'s product description - Reverse

Medacta® shoulder products used in reverse configuration are:

- Humeral diaphysis (cementless and cemented) to be used in total shoulder arthroplasty, for primary or revision surgery
- Humeral reverse metaphysis to be used in total shoulder arthroplasty, for primary or revision surgery
- Humeral reverse liner to be used in total shoulder arthroplasty, for primary or revision surgery
- Glenoid components (baseplate, polyaxial screws, glenosphere, glenosphere screw, Glenoid Reconstruction System central screw (optional)) to be used in total shoulder arthroplasty, for primary or revision surgery
- Double eccenter/reverse metaphysis screw to be used in total shoulder arthroplasty, for primary or revision surgery.

Total Shoulder Arthroplasty - Anatomic

A shoulder prosthesis designed for use in total shoulder arthroplasty in anatomic configuration consists of: a metallic humeral diaphysis, a metallic humeral anatomical metaphysis, a humeral stem screw, a double eccenter, a metal humeral head and a cemented glenoid component.
A stemless shoulder prosthesis designed for use in total shoulder arthroplasty in anatomic configuration consists of: a metallic stemless humeral metaphysis, a double eccenter, a metal humeral head and a cemented glenoid component.

Total Shoulder Arthroplasty - Reverse

A shoulder prosthesis designed for use in total shoulder arthroplasty in reverse configuration consists of: a metallic humeral diaphysis, a metallic humeral reverse metaphysis, a polyethylene liner, a metal glenosphere, polyaxial screws and a glenoid baseplate. If using a Glenoid Reconstruction System baseplate, an optional central screw can be used to provide additional compression to the glenoid implant.

3. INTENDED USE / INDICATIONS

Indications – Anatomic

The Anatomic Medacta Shoulder System is intended for use as a hemi or total shoulder replacement. It is a single-use device for the reconstruction of the gleno-humeral joint in case of osteoarthritis, rheumatoid arthritis, traumatic arthritis, avascular necrosis as well as complex fractures, post-traumatic malunion and revision of failed arthroplasty in patients with an intact or reparable rotator cuff shoulder joint. The humeral components are for cemented and cementless use.
The Medacta Anatomic Stemless Shoulder Prosthesis is indicated for primary partial or total shoulder replacement in patients with an intact or reparable rotator cuff and a severe arthropathy. The humeral bone quality has to be sufficient in order to support the cementless biological fixation of the stemless metaphysis.

Indications - Reverse

The Reverse Shoulder Prosthesis is indicated for treatment of humeral fractures and for primary or revision total shoulder replacement in patients with a grossly rotator cuff deficient shoulder joint with severe arthropathy or a previously failed joint replacement. The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device. The glenoid baseplate is intended for cementless application with the addition of screws for fixation.

4. CONTRAINDICATIONS

Total joint replacement is contraindicated where there is:

- Local or systemic infection or sepsis
- Insufficient bone quality which may affect the stability of the implant
- Muscular, neurological, or vascular deficiencies, which compromise the affected extremity
- Any concomitant disease and dependence that might affect the implanted prosthesis
- Materials (metals, etc.) sensitivity or allergy
- Loss of ligamentous structures
- Non-functional deltoid muscle

5. WARNINGS AND PRECAUTIONS

The success of the operation depends on compliance with the operative technique supplied, and the proper use of the dedicated instruments and specially designed for that range of implants.
The trial instrumentation must be used to confirm the choice of sizes and verify the functionality of the joint. Malpositioning may reduce implant longevity and lead to early implant failure.

MRI Compatibility

Medacta implants have not been tested for heating, migration, or to determinate the specifications for conditional status in the MR (magnetic resonance) environment. The safety of the implants in the MR environment is unknown, and the scanning of patients who have the implant may result in patient injuries.

HA coated implants

HA-coated surfaces must not be implanted in contact with cement.

Cemented implants

It is essential to follow carefully the instructions for use provided by the cement manufacturer because cement handling may influence the effectiveness of implant fixation. The use of Medacta® bone cement is strongly recommended.

6. MEDACTA INTERNATIONAL IMPLANTS

Medacta® International is not responsible for the use of its implant components in combination with a component from another manufacturer (unless otherwise specified by Medacta® International in the surgical technique), therefore we advise against such a use.
The components of a shoulder prosthesis should never be reimplanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure. The operating surgeon must be aware that even a very small superficial damage, caused for instance by a sharp tool or electrocautering, can have an influence on the endurance of the device and can lead to fracture.
Always use a trial prosthesis for trial purposes. Never adapt or alter trial prostheses.
UHMWPE implants should be stored for at least three hours at 20°C (+/- 3°C) before the operation. Th e surfaces of the articulated components must be cleaned and dried before reduction. A reliable seating is only possible when both surfaces are completely intact.
Medacta® shoulder products are non-pyrogenic.

7. RISK FACTORS

The following risk factors, individually or together, may result in poor clinical outcomes:

- Strenuous physical activities (active sports, heavy physical work)
- Fretting of modular junctions
- Incorrect implant positioning
- Muscle deficiencies
- Multiple joint disabilities
- Refusal to modify postoperative physical activities
- Patient history of infections or falls
- Systemic diseases and metabolic disorders
- Local or disseminated neoplastic diseases
- Drug therapies that adversely affect bone quality, healing, or resistance to infection
- Drug use or alcoholism
- Marked osteoporosis or osteomalacia
- Patient's resistance to disease generally weakened (HIV, tumour, infections)
- Severe deformity leading to impaired anchorage or improper positioning of implants

INSTRUCTIONS FOR USE

8. PREOPERATIVE PHASE

The surgeon should verify possible patient physical limitations and mental deficiencies and should also discuss with the patient all the details of the procedure and prosthesis. The discussion should consider the limitations of the procedure and the constraints imposed by the selected implant. The factors which could limit the performance and stability of the implant, e.g. level of activity, should be set out to improve the patient's chances of avoiding complications. The necessity of following the postoperative instructions given by the surgeon should be fully understood by the patient. A stock of sterile implants of suitable sizes should be available and checked before surgery.

9. HANDLING

To avoid scratching or damaging the implants, the implants should be handled with the utmost care by qualified personnel in an environment where conditions of hygiene are controlled. The implants should be kept in their undamaged packages until needed for use. Do not use implants from opened packages, implants that are damaged, or that are beyond their expiration date.

10. SURGICAL TECHNIQUE

The surgeon should be fully familiar with the surgical technique. Supplementary information about the surgical techniques and products are available on request. Careful preoperative planning, based on X-rays, is essential. X-ray templates are available for most implants.

11. POSTOPERATIVE CARE AND FOLLOW-UP

The surgeon should caution the patient to control their level of activity and avoid excessive loads on the replaced joint and make them aware of the precautions to be taken with regards to exercise, treatments and limitations on activities, as well as avoiding exposure to magnetic fields. Periodic follow-up and X-rays are recommended to make comparisons with the immediate postoperative condition and identify implant displacement, loosening, etc. Excessive physical activity, and operated limb traumas may cause early failure of the arthroplasty through implant displacement, fracture and/or wear. If the case occurs, it is necessary to place the patient under supervision, evaluate the possible progression of the deterioration, and weigh the benefit of early revision.

ADVERSE EFFECTS AND COMPLICATIONS

12. GENERAL INFORMATIONS

Adverse effects that can occur in shoulder arthroplasty include.

- Loosening of the prosthetic components
- Prosthesis dislocation and instability
- Damage to the prosthetic implant
- Instability of the system because of inadequate soft tissue balancing
- Dissociation due to incorrect coupling of the devices
- Infection
- Local hypersensitivity
- Local pain
- Periprosthetic fractures
- Temporary or permanent nerve damage
- Fractures of the devices
- Excessive wear of UHMWPE components due to damaged articular surfaces or the presence of particles
- Additional surgery

Some adverse effects can ultimately lead to death.

General complications include:

- Venous thrombosis with/without pulmonary embolism
- Cardiovascular or pulmonary disturbances
- Haematomas
- Systemic allergic reactions
- Systemic pain

13. PACKAGING

All the implants are supplied in single-use individual packages. For components delivered sterile, the sterilization method is indicated on the label. The expiration date and package integrity must be checked to ensure that sterility of the contents has not been compromised. If the package is damaged, do not use the component. Do not resterilize.

14. INSTRUMENTS

Instruments are supplied non-sterile and must be cleaned and sterilized prior to use. Recommended cleaning, decontamination and sterilization instructions are provided on www.medacta.com

15. STORAGE

The packages must be stored in a cool, dry place, away from light.

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