Aspiration Catheter

Instruction for Use

(Please read this Instruction for Use carefully before use, paying particular attention to the various warnings and precautions)

Shanghai Heartcare Medical Technology Co., Ltd.

[Product Name]

Aspiration Catheter

[Respectfully]

- This product has been sterilized with ethylene oxide.
- This is a single-use product, do not re-sterilize or re-use.
- Do not use if the package is opened, cracked, leaked or beyond the sterilization expiration date before use.
- This Instruction for Use must be read carefully before use.
- This product should only be used by professionally trained physicians.

[Packing List]

- Aspiration catheter
- Extension tube
- Inner catheter

[Structure and composition]

Aspiration Catheter consists of Aspiration catheter, Extension tube and Inner tube. Aspiration catheter includes the accessories: Shaping Needle, Protective Sheath and Rotating Hemostasis Valve. The Inner catheter includes the accessory, Inner catheter Shaping Needle. Aspiration catheter is made of braided tube, Handle, Stainless steel Diffusion Tube and Heat shrinkable tube. The distal outer layer of Aspiration catheter is coated with hydrophilic coating to reduce the friction when the catheter moves in the blood vessel. The extension tube is composed of long high-pressure braided Tube, short high-pressure braided Tube, sliding switch, connecting tube and Luer Locking Fittings. The inner catheter consists of Handle, Diffusion Tube, and a sheath. The specific structure of Aspiration Catheter is shown in Figure 1.

Extension tube (Figure 1.2) is for connecting aspiration catheter and HeartCare Aspiration Pump and controlling the aspiration switch; Protective Sheath (Figure 1.1.3) protects the tip of aspiration catheter when delivering aspiration catheter into the puncture sheath, which can be peeled away after use. Shaping Needle (Figure 1.1.2) can shape the tip of aspiration catheter when needed. Rotating Hemostasis Valve (Figure 1.1.4) is for connecting other surgical instruments and hemostasis; Inner catheter (Fig. 1.3.1) is for guiding aspiration catheter. Inner catheter Shaping Needle(Figure 1.3.2) is for shaping the tip of Inner catheter when needed. The assembly diagram of Aspiration Catheter is shown in Figure 2.


Figure 1.1.1 Schematic diagram of the structure of the aspiration catheter

Braided Mesh Tubing; 2. Stainless Steel Diffusion Tube; 3. Heat Shrink Tubing; 4.Handle; 5.Marker; 6.Hydrophilic Coating


Figure 1.1.2 Schematic diagram of the Shaping Needle

Figure 1.1.3 Schematic diagram of the Protective Sheath


Figure 1.1.4 Rotating Hemostasis Valve

Figure 1.1 Schematic diagram of the aspiration catheter, the


Figure 1.2 Schematic diagram of the extension tube structure

7. Male taper luer locking fitting; 8. Slide switch; 9. Connecting sleeve; 10. Long high pressure Braided Tube; 11. Short high pressure Braided Tube; 12. Female taper luer locking fitting


Figure 1.3.1 Schematic diagram of the structure of the inner catheter

13.Handle;14 Diffusion Tube;15.Sheath tube;16 Marker


Figure 1.3.2 Inner Catheter Shaping Needle

Figure 1.3 Schematic diagram of the Inner catheter

Figure 1. Schematic diagram of the Aspiration Catheter


Fig. 2 Schematic diagram of the assembly of Aspiration Catheter

Table 1 aspiration catheter specification dimensions table

Model	Adapted to the specifications of the Inner Catheter	Total length L2(cm)	Effective length L1(cm)	OD1(mm)	ID1(mm)	Hydrophilic Coating Length l1(cm)	Annex
WM*ASC-8F080	N075153S	88	80	2.65	2.23	40	Rotating Hemostasis Valve Shaping Needle Protective Sheath 8F
WM*ASC-8F095		103	95				
WM*ASC-8F105		113	105				
WM*ASC-8F115		123	115				
WM*ASC-8F125		133	125				
WM*ASC-8F130		138	130				
WM*ASC-7F095	N069153S	103	95	2.40	2.10	40	Rotating Hemostasis Valve Shaping Needle Protective Sheath 8F
WM*ASC-7F105		113	105				
WM*ASC-7F115		123	115				
WM*ASC-7F125		133	125				
WM*ASC-7F132		140	132				
WM*ASC-6F115	N059153S	113	115	2.10	1.80	40	Rotating Hemostasis Valve Shaping Needle Protective Sheath 6F
WM*ASC-6F120		123	120				
WM*ASC-6F125		133	125				
WM*ASC-6F132		140	132				
WM*ASC-5F125	N049153S	133	125	1.80	1.40	40	Rotating Hemostasis Valve Shaping Needle Protective
WM*ASC-5F127		135	127				
WM*ASC-5F132		140	132				
WM*ASC-5F139		147	139				

							Sheath 6F
WM*ASC-4F153	/	161	153	1.25	0.9	40	Rotating Hemostasis Valve Shaping Needle Protective Sheath 4F

Table 2 Extension Tube Specification dimensions Table

specifications and models	Effective length L3	OD2	ID2
WM*ASC-T290	290cm	4.8mm	2.2mm

Table 3 Inner Catheter Specification dimensions Table

Model/Spec	Effective length L4(cm)	OD3(mm)	ID3（mm）	Hydrophilic Coating Length l2(cm)	Annex
N049153S	153	1.20	0.54	40	Inner Catheter Shaping Needle
N059153S	153	1.47	0.54		Inner Catheter t Shaping Needle
N069153S	153	1.70	0.69		Inner Catheter Shaping Needle
N075153S	153	1.85	1.0		Inner Catheter Shaping Needle

Table 4 Shaping Needle specification dimensions table

Name	Effective length	Diameter
Shaping Needle	L5:70mm	1.0mm
Inner Catheter Shaping Needle	L6:70mm	0.45mm

[Indications]

Aspiration Catheter

The Aspiration Catheter is indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral-M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

Extension Tube

As part of the Aspiration Catheter, extension tubing is indicated to connect the Aspiration catheters and the medical negative pressure aspiration pump of our company.

[Contraindications]

No contraindication.

[Potential Complications]

- Potential complications include, but are not limited to.
- Infection
- Hematoma at the puncture site
- Distal Embolization
- Thrombosis

- Vascular stripping
- Pseudoaneurysm formation
- Acute vascular occlusion
- Blood clot formation
- Bleeding from the puncture site
- Intracranial hemorrhage
- Arterial rupture
- Stroke and death
- Unable to remove the clot completely

[Compatibility and required accessories]

- 7F or 8F specifications are compatible with 8F or larger puncture sheaths;
- 4F is compatible with 6F guide tube;
- 5F is compatible with 7F guide tube;
- 6F is compatible with 8F guide tube;
- 5F, 6F, 7F, 8F specifications aspiration catheter is compatible with. 035 " guide wire
- 4F aspiration catheter or internal conduit compatible with. 014 " micro guide wire;
- HeartCare Aspiration Pump;

[Warning]

- This product must be used by a physician trained in appropriate interventional neurosurgical endovascular techniques and in the treatment of acute ischemic stroke.
- This is a single-use product and should not be reused. The structural integrity and/or functionality of the product may be compromised during reuse or cleaning.
- If there is significant resistance to the product during introduction into the body, the catheter should not be forcibly advanced or reversed without careful assessment of the cause of the resistance using fluoroscopy. If the cause cannot be determined, the catheter must be withdrawn. Forced advancement may result in damage to the vessel or the catheter itself.
- If the flow of fluid in the Inner Lumen is obstructed, do not clear the Inner Lumen by instillation, as this may cause catheter rupture, which could result in vascular compromise. If fluid flow in the Inner Lumen is obstructed, remove and replace the catheter.
- To minimize the risk of complications due to Hydrophilic Coating dislodgement, the following operational recommendations need to be followed.
 - Before inserting the catheter into the puncture sheath, moisten the catheter with saline; when the catheter is introduced into the puncture sheath, it is necessary to introduce the tip of aspiration catheter coaxially into the Protective Sheath to avoid scraping of the distal Hydrophilic Coating at the puncture sheath peeling area.
 - During catheter advancement, the catheter should be advanced as gently and slowly as possible, as rapid advancement may result in kinking of the catheter body.
 - If the body of the catheter breaks during use, the catheter must be replaced with a new one.
- Fluoroscopy is required to ensure that the catheter port is aligned with the thrombus and not the vessel wall to avoid damage to the vessel.

[Notes]

- The product should be stored in a cool, dry place, away from direct sunlight.
- If the package is damaged, use is prohibited.

- The product should be used within the effective date.
- When removing the product from the package, it should be carefully inspected to ensure that the catheter is intact.
- Do not expose catheter to solvents.
- Use this product in conjunction with fluoroscopy and appropriate anticoagulants.
- Twisting of the product while the catheter is Bending can result in breakage of the catheter, which in turn can result in breakage of the catheter body.
- If a device becomes lodged in the catheter, or if a serious Bending situation occurs with the catheter, all devices(aspiration catheter, guidewires, and devices)should be withdrawn together.
- To avoid the formation of blood clots, the formation of contrast crystals, or the introduction of air into the catheter, maintain a steady and continuous flow of appropriate flushing solution into the catheter.
- When suctioning, make sure that the extension tube switch is open for no longer than the minimum time required by the Thrombectomy, and do not over-suction or fail to turn off the extension tube switch when suctioning is complete.
- The total time limit for completing revascularization for patients using the Aspiration Catheter is 120 minutes.
- Aspiration Catheter applies to negative pressure range of 0~- 92 ± 7kPa;

[Instructions for use]

- To determine the vessel diameter, select the appropriate dimensions and effective length of aspiration catheter(see Table 1).
- Before withdrawing the catheter from the coil, flush the lumen of the coil with heparinized saline.
- Flushing the lumen of the catheter with heparinized saline.
- If necessary, the attached Shaping Needle can be used to change the shape of the aspiration catheter tip. The shaping position is 1-2.5cm away from the steam source, and the shaping time is about 20-30s. If the distance is too short or the time is too long, the end of the conduit head may be damaged. If damaged, please replace it before use;
- Connect the attached Rotating Hemostasis Valve to aspiration catheter Handle;
- Using tearable Protective Sheath, Aspiration catheter is slowly inserted into Protective Sheath, then Protective Sheath and aspiration catheter are slowly threaded together into the guide Inner Lumen, and after threading, Protective Sheath is removed from the puncture sheath and removed after tearing.
- Using conventional catheterization techniques under fluoroscopic guidance, the aspiration catheter is slowly advanced into the target vessel with the assistance of an appropriately sized Inner Catheter and guidewire and the port of the aspiration catheter is placed proximal to the thrombus; the aspiration catheter is then placed into the target vessel.
- After the aspiration catheter is slowly advanced to the intended target, the guidewire and the inner catheter are withdrawn from the inner lumen of the aspiration catheter; the guidewire and the inner catheter are withdrawn from the inner lumen of the aspiration catheter.
- If the tip of the inner catheter needs to be shaped, the tip is shaped under steam using Shaping Needle supplied with the inner catheter, the procedure is the same as aspiration catheter.

- Connect the extension pipe to the collection tank of HeartCare Aspiration Pump, and open HeartCare Aspiration Pump. Before use, HeartCare Aspiration Pump should run stably for at least 1 min, and ensure that the sliding switch on the extension pipe is in the position of OFF(to close).
- Connect the extension tube to Rotating Hemostasis Valve side hole(, if not use Rotating Hemostasis Valve, directly connect to the Handle of Aspiration catheter);
- The aspiration catheter is gently pushed forward so that it is gently inserted into the thrombus.
- Start aspiration, open the sliding switch on the extension tube to the position of ON (to open), and aim the tip of aspiration catheter at the thrombus for suction;
- Stop aspiration, close the sliding switch on the extension pipe to the position of OFF (and), and close the pumping pump;
- If necessary, further retrieve the clot at the discretion of the physician.
- Use a 5ml or 10ml syringe to draw about 5ml of blood from aspiration catheter to remove any thrombus that may remain in the catheter;
- Contrast agent was injected through a guide catheter to obtain angiographic images after the surgery.
- Upon completion of the procedure, the product may be potentially bio-Hazard, and the product and packaging should be disposed of in accordance with hospital or administrative regulations or local government policy.

[Sterilization]

Aspiration Catheter, extension tube and inner catheter were sterilized by ethylene oxide.

[Storage conditions]

The product should be stored in a cool, dry, clean, well-ventilated, non-corrosive gas environment.

[Conditions of transportation]

Packaging and transportation requirements in accordance with the provisions of the order contract and the number of layers of transportation stacking, should be prevented from heavy pressure, direct sunlight and rain and snow soaking.

[Date of manufacture, date of use]

See package label.

[Expiry date]

Aspiration Catheter, extension tubes are valid for 3 years under conditions consistent with storage regulations.

[Symbol Description]

Symbol	Title
	Do not re-use.
	Refer to Instruction for Use
	Keep away from sunlight
	Keep dry
	Use-by date
	Batch code
	Catalogue number
	Sterilized using ethylene oxide
	Date of manufacture
	Do not use if the package is damaged.
	Do not resterilize

[After-sales service]

Shanghai Heartcare Medical Technology Co., Ltd. (here in after referred to as "Heartcare ")is committed to providing high quality, high performance medical products for the field of neurointervention. Heartcare guarantees that our products will be free from material or manufacturing defects when received by the customer, and any remaining questions regarding the products should be directed to the company.

[Solemn declaration]

Shanghai Heartcare Medical Technology Co., Ltd. (hereinafter referred to as "Heartcare ")expressly specifies that "Aspiration Catheter" produced by our company is for single-use, and shall not be reused. Repeated use is not permitted. Heartcare will not be responsible for any damage to the product or failure of the procedure due to improper product Model selection, incorrect operation or any other man-made accidents Events. Regarding the reusability of the system, Heartcare does not recommend, indicate or imply in any way, and is not responsible for any accidents or product damage caused by reuse.

Manufacturer: Shanghai Heartcare Medical Technology Corporation Limited
Address: Building 38, No.356 Zhengbo Road, LIN-GANG Special Area, China (Shanghai) Pilot Free Trade Zone, 201422 Shanghai, China.
Tel./ Fax:: +86-021-58975055
Post Code: 201201

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Embolic Protection System

Instructions for Use

(Please read this IFU carefully before use, pay special attention to various warnings and precautions)

Shanghai HeartCare Medical Technology Co., Ltd

[Product name]

Embolic Protection System

[Warning]

- This product is sterilized using ethylene oxide.
- This product is disposable. Please do not sterilize or reuse it.
- Do not use the package if it is found to have been opened, broken, leaked or has exceeded the sterilization validity period before use.
- Read this manual carefully before use.
- The product can only be used by professionally trained doctors.

[Packing List]

- A set of embolic protection device, including a capture wire and a dual-ended catheter.
- One copy of instructions for use.
- One copy of certificate of conformance.

[Product Description]

As shown in Figure 1, the embolic protection device consists of a capture wire and a dual-ended catheter. The capture wire contains a filter (1), a guide wire (2), a distal marker coil (3), etc. and the catheter contains a catheter delivery end (4), a catheter recovery end (5), a catheter body (6), etc., in which S1 and S2 are hydrophilic coating areas.


Figure 1 Schematic Assembly Diagram of the Embolic Protection Device

As shown in Figure 2, both ends of the filter of the capture wire are equipped with a distal marker band (1-1) and a proximal marker band (1-5). The filter is woven from nickel titanium alloy wire (1-3) with shape memory, and contains a limit tube (1-2) inside, which can pass through a small catheter when being limited, and fully expand in the blood vessel to closely fit the blood vessel wall and capture broken thrombus; the mouth indicator (1-4) made of gold-plated tungsten wire is set at the opening of the filter, therefore the state of which can be clearly indicated to make it convenient for doctors' judgment and operation.


Figure 2 The structure diagram of capture wire

The wire is available in 190 cm and 190/320cm. With a score, the 190/320cm wire is 320 cm long before being broken and can be broken to 190 cm when necessary; the rod body is coated with a black polymer coating with good biocompatibility, which improves the smoothness and trafficability of the embolic protection device in the catheter and blood vessel. The distal marker coil tightly wound with platinum iridium wire fully ensures the good development of the embolic protection device in clinical use, and can effectively indicate the position of the distal end of the embolic protection device in the human body.

The catheter is used for the delivery and recovery of the capture wire. Made of different materials through welding connection, it features both flexibility and smoothness; with smaller outer diameter, larger inner lumen and hydrophilic coating at both ends, it can pass through tortuous and complex blood vessels more easily and minimize the damage to the blood vessel wall. As shown in Figure 3, the catheter delivery end (green) and recovery end (blue) are respectively equipped with soft distal white tips (4-1), (5-1) and radiopaque marker bands (4-2), (5-2), which can effectively avoid scratching blood vessels and indicate the position of the catheter in the human body. The catheter delivery end has a main wire port (4-3) and an exchange port (4-4), and the recovery end has an exchange port (5-3).


Figure 3 Schematic Diagram of the Catheter Structure

The package of this product also contains a 23 gauge blunt needle irrigator of the standard Luer taper that doctors can use for catheter rinsing and perfusion.

Table 1 Model/Specifications and Applicable Vessels of Embolic Protection Device

Model specification	Diameter of applicable vessel (mm)	Capture wire			Catheter		
		Filter diameter (mm)	Capture wire length RX (cm)	Capture wire length OTW (cm)	Outer diameter of the delivery end (mm)	Outer diameter of the recovery end (mm)	Total Length mm
SM*EPS-030-320	2.1-2.5	3.0	190	320	1.2	1.4	1800
SM*EPS-035-320	2.6-3.0	3.5	190	320	1.2	1.4	1800
SM*EPS-040-320	3.1-3.5	4.0	190	320	1.2	1.4	1800
SM*EPS-045-320	3.6-4.0	4.5	190	320	1.2	1.4	1800
SM*EPS-050-320	4.1-4.5	5.0	190	320	1.2	1.4	1800
SM*EPS-055-320	4.6-5.0	5.5	190	320	1.2	1.4	1800

SM*EPS-060-320	5.1-5.5	6.0	190	320	1.2	1.4	1800
SM*EPS-065-320	5.6-6.0	6.5	190	320	1.2	1.4	1800
SM*EPS-070-320	6.1-7.0	7.0	190	320	1.2	1.4	1800
SM*EPS-030-190	2.1-2.5	3.0	190	/	1.2	1.4	1800
SM*EPS-035-190	2.6-3.0	3.5	190	/	1.2	1.4	1800
SM*EPS-040-190	3.1-3.5	4.0	190	/	1.2	1.4	1800
SM*EPS-045-190	3.6-4.0	4.5	190	/	1.2	1.4	1800
SM*EPS-050-190	4.1-4.5	5.0	190	/	1.2	1.4	1800
SM*EPS-055-190	4.6-5.0	5.5	190	/	1.2	1.4	1800
SM*EPS-060-190	5.1-5.5	6.0	190	/	1.2	1.4	1800
SM*EPS-065-190	5.6-6.0	6.5	190	/	1.2	1.4	1800
SM*EPS-070-190	6.1-7.0	7.0	190	/	1.2	1.4	1800

[Indications]

- The embolic protection device provides distal embolization protection during vascular interventional surgery, including peripheral vascular, coronary, and carotid interventions.

[Contraindications]

- Patients with a clear history of allergy to nickel-titanium alloy materials.
- The filter shall not be deployed in a previously implanted stent.
- The product shall not be used in excessively curved blood vessels.

[Potential Complications]

Potential complications may include but are not limited to:

- Access site adverse event (for example, arteriovenous fistula, hematoma, hemorrhage, pseudoaneurysm, puncture site infection)
- Adverse reaction to antiplatelet/anticoagulation agents or contrast media
- Allergic reactions to device materials
- Amaurosis fugax (carotid indication only)
- Aneurysm
- Angina (coronary indication only)
- Arrhythmia
- Artery dissection
- Pericardial tamponade (coronary indication only)
- Death
- Device deformation, collapse, fracture, or rupture
- Device thrombosis (acute and subacute)
- Embolism of air, debris, plaque or thrombosis (evidenced by slow-flow or no-flow) from mechanical disturbance of interventional surgery, resulting in temporary ischemia or stroke
- Embolization or migration of the interventional device
- Emergency surgery

- Gastrointestinal bleeding due to anticoagulation
- Hemodynamic compromise (e.g., prolonged hypotension requiring treatment with intravenous medications)
- Infections
- Intimal flap
- Ischemia
- Myocardial infarction
- Renal failure requiring dialysis
- Repeat intervention to treatment site
- Seizure (carotid indication only)
- Significant cardiac arrhythmia requiring treatment with medications and/or transvenous pacing
- Thrombosis (acute and sub-acute)
- Vasospasm
- Vessel perforation or rupture

[Warning]

- Please avoid scratching both ends of the catheter which are covered with hydrophilic coating and keep them away from organic solvents such as alcohol.
- Do not grip the catheter at the distal white tip. Gripping the catheter's distal white tip can cause the embedded radiopaque marker band or inner lumen of the catheter to deform.
- Failure to adequately close off the primary guidewire exit port when flushing the delivery end of the catheter can result in air embolism.
- Do not undersize or oversize the filter in relation to the selected vessel diameter. Undersizing or oversizing can result in inadequate vessel wall apposition, incomplete deployment of the filter, slow-flow or no-flow, and/or debris embolization. (Please refer to Table 1 to select the capture wire of appropriate specification)
- Minimize torquing of the capture wire after it has been introduced into the patient. Torquing the capture wire is likely to result in wire whipping and tangling.
- Do not apply excessive force while manipulating the capture wire. Excessive force can lead to distal embolization of debris, vessel damage, or device damage.
- Never withdraw or move an intravascular device against any resistance until the cause of the resistance is determined. It can lead to distal embolization of debris, vessel damage, or device damage.
- Do not place the capture wire tip in a vessel that does not allow for wire movement. Placing the capture wire in a vessel that does not allow for wire movement can lead to debris embolization, vessel damage, and/or device damage.
- Minimize movement of the capture wire after Initial placement. Excessive movement can cause debris embolization, vessel damage, and/or device damage.
- Always use the catheter when repositioning or removing the capture wire. It can lead to distal embolization of debris, vessel damage, or device damage.
- Observe the capture wire frequently under fluoroscopy, and monitor the patient to verify that the filter has not become occluded with debris, resulting in slow-flow or no-flow. Recover the filter if it becomes occluded or if the flow is compromised.
- Do not exceed 60 minutes of dwell time for the capture wire. Exceeding 60 minutes of dwell time can cause occlusion, resulting in slow-flow/no-flow.
- Exercise caution when using the partial recovery method to withdraw a full filter through a deployed stent. If resistance is met, advance the filter and recovery end of the catheter distally, and draw the filter more fully into the catheter before attempting withdrawal again.
- Do not use the capture wire with mechanical rotational thrombectomy, or with atherectomy or laser devices that require use of a specific guidewire.

[Cautions]

- The product shall be stored in a cool, dry and dark place away from direct sunlight;
- Do not use the product if its package is damaged;
- Do not use the product if it is bent or damaged;
- The product shall be used prior to the shelf life;
- Exposure to an environment above 54℃ (130℉) may cause damage to the embolic protection device and accessories;
- The product must not be sterilized using high-pressure steam;
- Use this product in combination with fluoroscopy and appropriate anticoagulants;
- Please confirm that the model, specification and size of this product are suitable for specific surgery before use;
- When taking the product out of the package, check it carefully to ensure that the product is intact;
- Torquing the embolic protection device may bend or separate the embolic protection device; If the embolic protection device is seriously bent, it shall be withdrawn with other supporting devices.

[Operating Instructions]

Preparations

- Administer anticoagulation medications and monitor activated clotting time per standard institutional guidelines. Antiplatelet medications are recommended when using this device for vascular interventions.
- Performing angiography.
- Identify the location, approximately 3 cm distal to the lesion site within the vessel where the filter will be deployed. Measure the vessel diameter at this location, and select the appropriate filter. Confirm that there is a vessel length of 4 to 5 cm distal to the lesion, to allow for adequate filter deployment.
- Ensure that the introducer catheter/sheath will accommodate the embolic protection device catheter.

Table 2 Specifications of Primary Guidewire and Introducer Catheter/sheath Applicable to the Embolic Protection Device

Diameter of the filter (mm)	Diameter of the primary guidewire (in)	Introducer catheter/sheath min I.D (in)
3.0	0.014-0.018	0.066
3.5	0.014-0.018	0.066
4.0	0.014-0.018	0.066
4.5	0.014-0.018	0.066
5.0	0.014-0.018	0.066
5.5	0.014-0.018	0.066
6.0	0.014-0.018	0.066
6.5	0.014-0.018	0.066
7.0	0.014-0.018	0.066

- Using sterile technique, remove the embolic protection device from the packaging, and place in a sterile work area.
- Remove 23 gauge blunt needle irrigator and packaging hoop extension.
- Submerge the filter in heparinized saline to saturate and remove air from the filter.
- As shown in Figure 4, carefully hold the area near the distal white tip (4-1) of the catheter delivery end and the catheter body (6) to support the catheter shaft, and pull the guidewire (2) to the proximal end until the filter(1) and the distal marker coil (3) completely enter the transparent section (4-5) of the catheter delivery end.


Figure 4 Schematic Diagram of Loading the Capture Wire in Delivery End

- Flush heparinized saline through the distal white tip of the delivery end until all air is removed and fluid passes from the primary guidewire exit port. The 23 gauge blunt needle irrigator can be used, if desired.
- Gently apply pressure to the primary guidewire exit port. Flush until all air is removed and fluid passes from the exchange port of the delivery end.
 Note: If using the 320/190 cm convertible version of the capture wire as a rapid exchange wire, firmly grasp the capture wire on each side of the score of the wire, break off the proximal end of the wire and discard it.


Figure 5 Schematic Diagram of the Score of the Guidewire

Delivery and Deployment

- Using standard technique to access the vessel, advance the primary guidewire at least 5 cm beyond the distal edge of the lesion being treated.
- Ensure that the opening on the primary guidewire exit port is visible.
- Load the distal white tip of the delivery end onto the proximal end of the primary guidewire. Advance the delivery end until the primary guidewire exits at the exit port.
 Caution: Do not allow the primary guidewire to advance proximally past the primary guidewire exit port, otherwise the capture wire may be damaged.
- Gently advance the delivery end of the embolic protection device, carrying the capture wire and into the introducer catheter/sheath through the primary guidewire. Continue advancement until the radiopaque marker at the distal tip of the delivery end is at least 4 to 5 cm beyond the distal edge of the lesion.
 Caution: Clearly identify the capture wire and primary guidewire while advancing. Do not completely rotate the catheter. Rotating can cause the wires to wrap around the catheter.
- Holding the catheter stationary, withdraw the primary guidewire, leaving the delivery catheter and capture wire in place.
- Holding the catheter stationary, gently advance the capture wire until the distal radiopaque marker band on the filter aligns with the rachopaque marker on the catheter delivery end.
- Observe the filter position to ensure that the proximal radiopaque marker band is positioned at least 2 cm distal to the anticipated distal end of the interventional device that will be used to treat the lesion. Gently advance the components of the delivery end as a unit if necessary so the proximal radiopaque marker band of the filter is at least 2 cm beyond the lesion distal edge.
 Caution: Position the filter sufficiently distal to the lesion so that the interventional device can be advanced without coming within 2 cm of the proximal end of the filter.
- Holding the capture wire stationary, gently withdraw the catheter to expose and deploy the filter in the vessel. To reposition the filter (if required), advance the catheter back over the filter, adjust the position of the catheter and capture wire together, then redeploy the filter. Never attempt to move the filter when it is outside of the catheter.
 Caution: Always use the catheter when repositioning or removing the capture wire. Excessive force can lead to distal embolization of debris, vessel damage, or device damage.

- The gold radiopaque proximal mouth indicator will be visible when the filter is properly opened. If the proximal mouth indicator is not visible, attempt to redeploy the filter or to remove the filter/capture wire/primary guidewire, and the delivery catheter together. Do not reintroduce or reuse the capture wire after it has been removed from the patient.
- Confirm that the filter and proximal radiopaque marker band are in opposition to the vessel wall.
- After the filter is deployed in the desired location, carefully remove the catheter from the patient. Keep the catheter in the sterile field. Do not discard the catheter. The recovery end (blue) of the catheter is required for capture wire recovery.
- Place the catheter in sterile heparinized saline or in a sterile field, and save the catheter for recovering the capture wire. Do not bend or kink the catheter.
- The 0.014 inch capture wire can be used as the primary guidewire for other interventional devices. Always keep the capture wire stationary during use.

Recovery and Removal

- After the interventional procedure is completed, wet the distal 20 cm of the recovery end (blue) of the catheter with heparinized saline. Remove all air from the recovery end by flushing the distal tip until fluid passes from the exchange port of the delivery end.
- Load the recovery end (blue) of the catheter onto the capture wire and advance the catheter through the introducer catheter/sheath.
- Advance the catheter until the distal tip radiopaque marker aligns with the proximal radiopaque marker band on the filter.
- The recovery of filter can be divided into partial enclosure recover and full enclosure recovery.
 A-Partial enclosure recovery: Gently advance the recovery end over the filter; until the proximal portion of the filter including the mouth indicator is inside the catheter and is fully compressed against the capture wire.


Figure 6 Partial Enclosure Recovery

B-Full enclosure recovery: Gently advance the recovery end over the filter until all radiopaque markers on the filter are within the catheter. The distal marker coil can remain outside the catheter.

Figure. 6 Schematic Diagram of Full Enclosure Recovery

- Carefully remove the catheter and capture wire together. To allow the catheter to be removed without resistance, open the hemostatis valve on the introducer catheter/sheath. Use care to avoid interaction with the intervention site. Discard the embolic protection device after removal from the patient.

[Sterilization]

The embolic protection device has been sterilized using ethylene oxide.

[Storage Conditions]

- Do not store in direct sunlight.
- Keep dry.

[Transportation Conditions]

During transportation, heavy pressure, direct sunlight, rain and snow shall be prevented.

[Production Date, Service Life]

See package label.

[Shelf Life]

The shelf life is 2 years under the condition of meeting the storage regulations.

[Date of preparation/revision of the instructions for use]

Under the conditions that comply with storage regulations, the validity period is 2 years.

[Symbol Description]

	Manufacturer
	Do not re-use
	Consult instructions for use
	Keep away from sunlight
	Keep dry
	Date of Manufacture
	Use-by date
	Batch code
	Catalogue number
	Sterilized using ethylene oxide
	Do not use if package is damaged
	Do not resterilize
	The quantity of the product inside is 1 piece

[After-Sales Service]

Shanghai HeartCare Medical Technology Co., Ltd. (hereinafter referred to as "HeartCare Medical") is committed to providing high-quality and high-performance medical products for the field of neurological intervention. HeartCare Medical guarantees that there are no material or manufacturing defects when the customer receives the products produced by our company. Please ask the company directly for other questions

related to the products.

[Solemnly Declare]

Shanghai HeartCare Medical Technology Co., Ltd. (hereinafter referred to as "HeartCare Medical") clearly indicates that the "embolic protection device" produced by our company is disposable and cannot be reused. HeartCare Medical will not be responsible for any product damage or operation failure caused by improper product model/specification selection, wrong operation or any other man-made accident. Regarding the reusability of the system, HeartCare Medical does not recommend, indicate or imply in any form, and the company is not responsible for accidents or product damages caused by repeated use.

Manufacturer: Shanghai HeartCare Medical Technology Corporation Limited

Address: Building 38, No.356 Zhengbo Road, LIN-GANG Special Area, China (Shanghai) Pilot Free Trade Zone, 201422 Shanghai,China.

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Neuro Micro Guide Wire

Instruction for Use

(Please read this IFU carefully before use, and pay special attention to all kinds of warnings and precautions)

Shanghai HeartCare Medical Technology Corporation Limited

[Product Name]

Neuro Micro Guide Wire

[Warnings]

- This product is sterilized using ethylene oxide.
- This product is disposable. Please do not reuse it. Reuse, reprocessing or re-sterilization may damage the structural integrity of the device and/or cause device fault, which may accordingly lead to injury, illness or death of patients. Reuse, reprocessing or re-sterilization may also cause the risk of device contamination and/or lead to infection or cross infection of patients, including but not limited to the transmission of infectious diseases from one patient to another. Contamination of the device may lead to infection, illness or death of the patient.
- Do not use if the packaging is found to be opened, broken, leaked or has exceeded the sterilization period before use.
- Read this IFU carefully before use.
- The product can only be used by professionally trained doctors.
- After use, dispose of the product according to relevant regulations of hospitals, administrative departments and/or local governments.

[Package List]

- A set of Neuro Micro Guide Wire.
- Related accessories (including introducer, torque device and shaping mandrel).
- One copy of IFU.
- One copy of certificate of conformance.

[Product Description]

As shown in Figure 1, the Neuro Micro Guide Wire is composed of (1) Distal radiopaque coil, (2) Stainless steel coil, (3) Guide wire shaft, (4) Distal coating, (5) Proximal coating.




Figure 1 Structural Schematic Diagram of Neuro Micro Guide Wire

As shown in Figure 2, the accessories include (1) Introducer, (2) Torque device and (3) Shaping mandrel.

(1) Introducer

(2) Torque device

(3) Shaping mandrel

Figure 2 Schematic Diagram of the Accessories

The distal radiopaque coil of Neuro Micro Guide Wire is made of platinum alloy, with the length of about 3 ~ 5cm.

The guide wire can be divided into three specifications based on its proximal diameter: 0.010inch, 0.014inch, and 0.016inch. The length of the 0.010inch guide wire is 200cm; The 0.014inch and 0.016inch guide wires have two lengths, 200cm and 300cm, respectively; A 200cm length guide wire can be quickly exchanged, while a 300cm length guide wire can be quickly exchanged and coaxially exchanged.

The 010 guide wire is a stainless steel mandrel; 014 and 016 guide wires are hybrid mandrels with a nitinol mandrel attached to a stainless steel mandrel.

All specifications of guide wires are coated with polymer coating on the proximal end and hydrophilic coating on the distal end.

See Figure 1 and Table 1 for the specification and dimension of the Neuro Micro Guide Wire.


Figure 3 Schematic Diagram of Dimensioning of Neuro Micro Guide Wire

Table 1 Specification and Dimension Table of Neuro Micro Guide Wire

Model and Specification	D		L1/mm	L2/mm	L3/mm
	inch	mm			
SM*NMW-10200	0.010±0.001	0.2540±0.0254	2000±50	30±1.5	270±10
SM*NMW-14200	0.014±0.001	0.3556±0.0254	2000±50	50±2.5	370±10
SM*NMW-16200	0.016±0.001	0.4064±0.0254	2000±50	50±2.5	370±10
SM*NMW-14300	0.014±0.001	0.3556±0.0254	3000±50	50±2.5	370±10
SM*NMW-16300	0.016±0.001	0.4064±0.0254	3000±50	50±2.5	370±10

[Indication for Use]

- For use in the neurovascular system in conjunction with a diagnostic or therapeutic catheter for placement into the selected lesion site.

[Contraindications]

- Patients allergic to contrast media;
- Patient allergic to nickel-titanium alloy materials;
- Patients in whom antiplatelet and/or anticoagulation therapy is contraindicated.

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- Aneurysm rupture;
- Infection, hematoma, nerve injury and other contact site complications;
- Cerebral ischemia;
- Death;
- Embolization (catheter/device, air bubble, thrombus or debris);
- Intracranial hemorrhage;
- Pseudoaneurysm;
- Seizures;
- Stroke;
- Transient ischemic attack (TIA);
- Vasospasm;
- Vascular perforation, dissection and injury.

[Warnings]

- Improper use may lead to vascular injury. Please follow the instructions carefully. Operate the guide wire only under fluoroscopy after it enters human body.
- Do not attempt to move the guide wire without observing tip response. Do not leave the guide

wire in a prolapsed state, otherwise it may be damaged. Do not advance the wire under excessive resistance without determining the cause of the resistance under fluoroscopy, as this may cause the guide wire tip to separate, damage the catheter and/or perforate the blood vessel.

[Precautions]

- Physicians can use the product only after receiving comprehensive training in the techniques and operations related to percutaneous intervention, intravascular intervention, etc.
- The Neuro Micro Guide Wire and its accessories have been sterilized and are non-pyrogenic. The integrity of product packaging shall be checked before use to ensure no damage to the sterility of the product.
- Do not use or attempt to re-sterilize the product if the integrity of product packaging is damaged.
- The product shall be used within its shelf life.
- Do not use if the labeling is incomplete or illegible.
- Check the diameter of the matched catheter before use to ensure its compatibility with the outer diameter of the guide wire.
- The hydrophilic coating may wear out during operation due to the change in the inner diameter of the catheter tip. In case of any resistance during catheter advancing, the catheter may need to be replaced by a new one.
- Check the guide wire for any surface irregularity, bend or kinking before use. Do not use damaged and/or irregularly surfaced guide wire.
- Operate it carefully during use to reduce the possibility of rupture, bending or kinking.
- Do not withdraw or control the guide wire through a metal pin cannula in order to avoid damage to it.
- Excessively screwing down the torque device on the guide wire may lead to wear of its coating.
- Such accessories as the introducer, shaping mandrel, torque device, etc. shall not enter the human body.
- Please confirm that the model, specification, dimensions, etc. of the product are applicable to specific surgery before use.

[Instructions for Use]

- Flush the coil with heparinized saline before taking out the guide wire from the coil in order to hydrate the hydrophilic coating of the guide wire. Flush it again if the hydration effect disappears.
- Flush the catheter with heparinized saline before inserting the guide wire into the catheter in order to hydrate the hydrophilic coating of the guide wire, so that the guide wire can move smoothly in the catheter.
- Remove the guide wire from the coil slowly and gently to avoid damage to it, and check the guide wire after removing it to ensure that there is no damage.

- Gently shape the flexible tip of the guide wire in the usual way if necessary.
- Carefully introduce the guide wire into the catheter hub or hemostatic valve by the introducer.
- Use the guide wire in conjunction with the torque device if necessary. Insert the proximal end of the guide wire into the torque device during use. When the torque device is in the desired position, tighten the nut to secure the torque device. Reposition or demount the torque device by unscrewing and screwing down the nut.
- Slowly and carefully rotate the torque device connected to the guide wire, so as to rotate the tip of the guide wire and enter different blood vessels.

[Sterilization]

The guide wire is sterilized by ethylene oxide.

[Storage Conditions]

- Store at room temperature.
- Do not store in direct sunlight.
- Keep dry.

[Transport Conditions]

In transportation, it is required to prevent the product from heavy load, direct sunlight, rain and snow.

[Production Date, Service Life]

See packaging label.

[Shelf Life]

The shelf life is 3 years under the condition of meeting the storage regulations.

[Date of Preparation/Revision of the Instructions for Use]

October 25, 2025.

[Symbol Description]

1)		**Do not reuse**
2)		**Consult Instructions for Use**
3)		**Keep away from sunlight**
4)		**Keep away from rain**
5)		**Shelf Life**

6)		**Batch code**
7)		**Catalogue number**
8)		**Sterilized using ethylene oxide**
9)		**Do not use if package is damaged**
10)		**Number of Product Contained: 1**
11)		**Manufacturer**
12)		**Do not re-sterilize**
13)		**Manufacturing Date**

[Manufacturer]

Shanghai HeartCare Medical Technology Corporation Limited (hereinafter referred to as "HeartCare") is committed to providing high-quality and high-performance medical products for the field of neurological intervention. HeartCare guarantees that there is no material or manufacturing defects when the customer receives the products produced by our company. Please ask the company directly for other questions related to the products.

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(Shanghai Xinwei Medical Technology Co., Ltd.)

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Support Catheter

Instruction for Use

(Please read this Instruction for Use carefully before use, paying particular attention to the various warnings and cautions)

Shanghai HeartCare Medical Technology Co., Ltd.

[Product Name

Support Catheter

[Statements]


- This product has been sterilized with ethylene oxide.
- This is a single-use product, do not re-sterilize or re-use. Repeated sterilization or reuse can disrupt the lubricating properties of the distal end of the catheter, preventing the catheter from reaching its vascular target successfully; even severe damage to the structural integrity of the device can occur.
- Do not use if the package is opened, cracked, leaked or beyond the sterilization expiration date before use.
- This Instruction for Use must be read carefully before use.
- This product should only be used by professionally trained physicians.

[Packing List]

- One set of support catheter
- Use of one copy of Instruction for Use


[Product Description]

The support catheter consists of a catheter, introducer sheath, extension tube and Y-valve; the catheter consists of a braided mesh tube, stress diffusion tubes (inner and outer) and a handle; the outer layer of the distal end of the catheter is coated with a hydrophilic coating; the braided mesh tube is a single-lumen, three-layer tube, and the developer ring is designed at the distal end of the braided mesh tube to develop the image under X-ray. The structural composition of the support catheter is shown in Fig. 1, and the product model specifications are shown in Table 1.

1 Hub; 2 inner stress diffusion tube; 3 outer stress diffusion tube; 4 braided mesh tube;

Figure 1a) Structural Schematic Diagram of Support Catheter

5 guide sheath connectors; 6 guide sheath tube;

Figure 1b) Schematic diagram of the guiding sheath

Figure 1c) Schematic diagram of the extension tube


Figure 1d) Schematic diagram of Y valve

Fig. 1 Schematic diagram of the support catheter structure

[Model specification]

Table 1 Specification Table of Support Catheter

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[Indications]

- The product is suitable for the introduction of interventional or diagnostic devices into vessels in the neurovascular, peripheral vascular and coronary vascular systems.

[Contraindications]

- Patients who are allergic to contrast media.
- Patients for whom antiplatelet or/and anticoagulation therapy is contraindicated.

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- Acute obstruction
- Pseudoaneurysm formation
- Air embolization
- Death

- Distal embolization
- Embolus
- Hematoma or hemorrhage at puncture point
- Infection
- Intracranial hemorrhage
- Ischemia
- Neurologic deficits including stroke
- Vasospasm, thrombosis, stripping, or perforation
- Intracranial hemorrhage

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- 6F size support catheter used with 5F or less size catheter
- 7F size support catheter for use with catheters of size 5F or less
- 8F size support catheter for use with catheters 6F and below
- Use with .035“/.038” guidewire

[Warning]

- This product must be used by a physician with appropriate specialized training in interventional techniques.
- This product is a single-use product and should not be reused. After single use, the product should be disposed of according to regulations. The structural integrity and/or function of the product may be damaged during reuse or cleaning. The catheter is extremely difficult to washing completely after entering the body and may cause cross-infection if reused.
- If there is significant resistance to the product during penetration into the body, the catheter should not be forcibly advanced or reversed without careful assessment of the cause of the resistance using fluoroscopy. If the cause cannot be determined, the catheter must be withdrawn. Forced advancement may result in damage to the vessel or the catheter itself.
- If the flow of fluid in the Inner Lumen is obstructed, do not clear the Inner Lumen by instillation, as this may cause catheter rupture, resulting in vascular compromise. Remove and replace the catheter if fluid is obstructed in the Inner Lumen.
- The catheter must be threaded into the body with great care to avoid damage to blood vessels; the catheter may occlude small blood vessels; care must be taken to avoid complete interruption of blood flow.

[Notes]

- Verify that the product model number and specifications apply to the current procedure before use.
- The product should be stored in a cool, dry and dark place, away from direct sunlight.
- Prohibit use if the package is damaged.

- Prohibit use if conduit is bent or damaged.
- The product should be used before the expiration date.
- Autoclaving shall not be permitted.
- Do not adopt high-pressure steam sterilization;
- To avoid thrombus formation or contrast crystallization, appropriate flushing solutions should be continuously injected through the catheter lumen.
- Twisting the catheter may result in bending or separation of the catheter shaft; if severe bending of the catheter occurs, all instruments should be withdrawn together.

[Instructions for use]

◆ **Long Sheath Procedure**

- Select the appropriate type of support catheter product for the anatomical location.
- Slowly pull the catheter with the plate out of the bag and carefully remove the catheter, introducer sheath, extension tube and Y-valve from the plate.
- Inspect the product for bends or other damage and replace if damaged.
- Flush the lumen of the catheter and introducer sheath with heparinized saline before use.
- Rotate and secure the hemostatic valve end of the extension tube to the catheter at the handle, then insert the introducer sheath into the catheter along the hemostatic valve so that the introducer sheath is tightly snapped into place at the handle of the introducer sheath and the hemostatic valve cap, and set aside.
- Using aseptic technique, the cannula of the vascular sheath puncture needle is threaded into the vasculature, the puncture needle is pulled out by holding the puncture needle cannula still, and then the flexible end of the guidewire is inserted into the puncture needle cannula, and the guidewire is pushed to the target position.
- Hold the guidewire still and pull out the puncture needle cannula, then apply pressure to the puncture point until the catheter is inserted into the vessel.
- Holding the catheter, introducer sheath, and guidewire together, the catheter is held close to the skin; in a rotational motion, the device is moved through the skin into the vessel.
- Push the catheter to the target position under x-ray fluoroscopy and withdraw the guidewire and introducer sheath together with the catheter.
- Remove any air that may be present in the device by suctioning through the three-way valve end of the extension tube.
- Select the appropriate size guidewire for insertion into the catheter and push the guidewire to the distal end of the catheter; push the guidewire and catheter to the target vessel location.
- Introduce other instruments, along the catheter (and guidewire) into the vessel;
- Removal of the catheter after completion of the procedure requires holding the vessel at the puncture point and then slowly withdrawing the catheter.

◆ **Guiding Catheters Procedure**

- Select the appropriate type and size of the supporting catheter product according to the anatomical location.
- Check the catheter and other instruments for bending or other damage, and replace them if damaged.
- Connect the Y-valve and catheter and flush the catheter lumen with heparinized saline.
- To use, insert the guidewire into the catheter and push the guidewire until its distal end is at the distal end of the catheter.
- A puncture sheath matching the catheter's outer diameter is placed according to standard interventional practice.
- Simultaneously push the catheter and guidewire into the puncture sheath, pushing the guidewire beyond the catheter.
- The catheter and guidewire are pushed at the same time and the guidewire is utilized to find the target vessel location.
- Once the guidewire reaches the target vessel location, push the catheter to reach the desired vessel location under x-ray fluoroscopy and withdraw the guidewire.
- Other instruments, along the catheter (and guidewire) are introduced into the vessel;
- When the procedure is complete, the catheter is slowly withdrawn.

[Sterilization]

The support catheter was sterilized by ethylene oxide.

[Storage conditions

Store in a cool, dry, well-ventilated environment.

[Conditions of transportation]

Transportation should be prevented from heavy pressure, direct sunlight and rain or snow soaking.

[Date of manufacture, expiration date]


See package label.

[Expiry date]

Valid for 3 years under storage conditions.

[Symbol Description]

1）		**No secondary use**
2）		**Refer to the instructions for use**
3）		**Avoid sun exposure**

4）		**Avoid rain**
5）		**Term of use**
6）		**Batch code**
7）		**Classification Number**
8）		**Do not re-sterilize**
9）		**Sterilized by ethylene oxide**
10）		**Sterility**
11）		**Do not use if the package is damaged**
12）		**Number of products inside 1**

[After-sales service]

Shanghai HeartCare Medical Technology Co., Ltd.(hereinafter referred to as "HeartCare Medical")is committed to providing high quality, high performance medical products for the field of neurological interventions. HeartCare Medical warrants that the products manufactured by the Company will be free from material or manufacturing defects when received by the customer, and that remaining questions regarding the products should be directed to the Company. Any remaining questions regarding the product should be directed to the company.

[Solemn declaration]

Shanghai HeartCare Medical Technology Co., Ltd.(hereinafter referred to as "HeartCare Medical")clearly indicates that the "support catheter" produced by our company is for single use and should not be reused.HeartCare Medical will not be held responsible for any damage to the product, failure of the procedure, etc. Events caused by improper use of the product Model selection, operator error, or any other man-made accident. HeartCare Medical does not recommend, indicate, or imply in any way the reusability of the system, and is not responsible for accidents or damage to the product caused by reuse

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Address: Building 38, No.356 Zhengbo Road, LIN-GANG Special Area, China (Shanghai) Pilot Free Trade Zone, 201422 Shanghai, China

After-sales service unit: Shanghai Heartcare Medical Technology Corporation Limited

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(Shanghai HeartCare Medical Technology Co., Ltd.)

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