Medtronic

PulseSelect™ PSCC100

Pulsed Field Ablation Loop Catheter
Катетър с лупинг за аблация с импулсно поле
Smyčkový katétr pro ablaci pulzním polem
Løkkekateter til pulserende feltablation
Loop-Katheter für Pulsed Field Ablation
Καθετήρας βρόχου για κατάλυση παλμικού πεδίου
Catéter de espiral para ablación por campo de impulsos
Impulssvälja-ablatsiooni silmuskateeter
Pulssikenttäablaatiossa käytettävä lassokatetri
Cathéter à boucle d'ablation à champ pulsé
Kateter s petljom za ablaciju pulsirajućim poljem
Impulzusmező-ablációs hurokkatéter
Catetere circolare per ablazione a campo pulsato
Impulsinės lauko abliacijos kilpinis kateteris
Impulsu lauka ablācijas cilpveida katetrs
Ablatieluskatheter gepulseerd veld
Løkkekateter for pulsfeltablasjon
Cewnik w kształcie pętli do ablacji pulsacyjnym polem elektrycznym
Cateter de aro de ablação por campo pulsado
Cateter în buclă pentru ablație cu câmp pulsatil
Петлевой катетер для абляции импульсным полем
Slučkový katéter na abláciu pulzným poľom
Kateter z zanko, za ablacijo s pulzirajočim poljem (PFA)
Slingkateter för pulsfältsablation
Atımlı Alan Ablasyon Lup Kateteri

Technical Manual • Техническо ръководство • Technická příručka • Teknisk håndbog • Gebrauchsanweisung • Τεχνικό εγχειρίδιο • Manual técnico • Tehniline juhend • Tekninen opas • Manuel technique • Tehnički priručnik • Műszaki leírás • Manuale tecnico • Techninis vadovas • Tehniskā rokasgrāmata • Technische handleiding • Teknisk håndbok • Instrukcja techniczna • Manual técnico • Manual tehnic • Техническое руководство • Technická príručka • Tehnični priročnik • Teknisk manual • Teknik El Kitabı

CE
0123

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1 Description


The Medtronic PulseSelect Pulsed Field Ablation (PFA) loop catheter is a sterile, single use, anatomically designed, multi-electrode catheter used to detect electrical signals, pace, and ablate the pulmonary veins (*Figure 1*). The catheter is provided sterile via electron beam (e-beam) sterilization. The catheter is used to perform cardiac tissue ablation through the mechanism of irreversible electroporation. The catheter is compatible with the Medtronic PulseSelect PFA system. The catheter is designed to be used together with an introducer sheath with a minimum inside diameter of 3.33 mm (10 Fr) and a 0.81 mm (0.032 in) J-tip guide wire.

The catheter has the following 2 control features:

1. A steering lever (4, *Figure 1*) that controls bidirectional deflection of the distal shaft.
2. A slide control (3, *Figure 1*) to extend or retract the guide wire lumen and distal tip (8, *Figure 1*) and change the catheter geometry from circular to spiral to a longitudinal shape.

The slide control also allows for extension of the distal portions of the catheter to allow capture of the distal array into the capture device that is used to deliver the catheter into the sheath. The catheter is inserted into the patient's vasculature, typically through an introducer sheath via the femoral vein, then advanced into the left atrium. Once in the sheath, the catheter and guide wire are advanced through the sheath with the guide wire being guided into one of the pulmonary veins. After the catheter array exits the sheath, it is deployed into a circular shape by pulling the slide control back. The catheter array is then placed in multiple sites for cardiac tissue ablation.

Figure 1. PulseSelect PFA loop catheter

1 Spiral array	6 Connector
2 Capture device	7 Luer
3 Slide control	8 Distal tip
4 Steering lever	9 Electrode 1
5 Steering lever tension control knob	10 Electrode 9

The PulseSelect PFA loop catheter can only be used with the PulseSelect PSG100 Pulsed Field Ablation generator and the PSCIC101 PulseSelect PFA catheter interface cable. For details about the PulseSelect system and how to use it with the catheter to perform ablation procedures, see the *PulseSelect Pulsed Field Ablation System Operator's Manual*.

1.1 Contents of package

The package contains the following:

- PulseSelect PFA loop catheter
- Product documentation

Observe all symbols on the back panel of the PSG100 PulseSelect generator and carefully read all instructions before use. Observe all contraindications, warnings, and precautions noted in these instructions. Failure to read the instructions may result in less than expected outcomes.

Caution: Visually inspect the sterile packaging and catheter before use. If the sterile packaging or catheter is damaged, do not use the catheter. Contact your local Medtronic representative.

Caution: The catheter should be used only by, or under the supervision of physicians trained in left atrial ablation procedures using this catheter and the Medtronic PulseSelect PFA generator.

2 Indications for use

The PulseSelect PFA loop catheter is indicated for use in cardiac electrophysiological mapping (stimulation and recording) and for the treatment of atrial fibrillation (paroxysmal and persistent).

3 Contraindications

The PulseSelect PFA loop catheter is contraindicated for use in patients with the following conditions:

- Active systemic infections
- A known sensitivity to Heparin
- Blood clotting abnormalities
- Permanently implanted metallic objects in the left atrium

The catheter is also contraindicated in conditions where the manipulation of the catheter within the heart would be unsafe, such as intracardiac mural thrombus.

The catheter is not recommended for use in patients who cannot undergo standard anticoagulation protocol for a left-sided cardiac procedure, or who have had a recent coagulopathy or embolic event.

4 Intended purpose

The intended purpose of the PulseSelect PFA loop catheter is to detect electrical signals, ablate cardiac tissue, and perform cardiac stimulation during electrophysiology procedures.

5 Intended user

The intended user of this ablation system is a physician that performs cardiac ablation. Use of this device does not require specific qualifications or training beyond what is required per the applicable local, national, or federal regulations.

6 Intended patient population

The intended patient population for the PulseSelect PFA loop catheter includes adults (patients 18 years or older) who require or are candidates for cardiac ablation. Ablation in patients who are pregnant or breastfeeding is not recommended. For the most up-to-date medical consensus on additional considerations, refer to the applicable treatment guidelines.

7 Clinical benefit

The expected clinical benefit to the patient is the treatment of arrhythmia. This includes restoring a patient's normal heart rhythm, which has been shown to reduce symptoms and improve quality of life.

8 Summary of safety and clinical performance

The Summary of Safety and Clinical Performance (SSCP) can be found at https://ec.europa.eu/tools/eudamed. Search for the SSCP using the manufacturer and device name, and any of the following elements, as applicable: device model, reference number, catalog number, or the Basic Unique Device Identification (Basic UDI-DI) number — 0763000B00012287M.

9 Warnings and precautions

Ablation therapy hazards – To reduce the possibility of hazards associated with use of the PulseSelect PFA loop catheter:

- Use the catheter only in the recommended anatomical location.
- Maintain the catheter position during the ablation.
- If coughing occurs, reposition the catheter more proximally and review sedation management.
- Ensure electrodes are not in contact with any metal during ablations (for example the guide wire).
- Maintain substantially circular array to ensure uniform field distribution.
- If the electrode array is deployed to deliver ablation energy, avoid continuing to move the slide control forward to prevent the guide wire lumen from coming too close to the electrode array.

Air entrapment – It is recommended that the array be captured while it is submerged to help reduce the possibility of air becoming entrapped around the electrode array during capture and catheter insertion.

Anticoagulation therapy – Administer appropriate levels of anticoagulation therapy for patients undergoing left-sided and transseptal cardiac procedures to minimize the risk of thromboembolic events. It is recommended to maintain therapeutic anticoagulation levels, per guidelines for patients undergoing cardiac ablation procedures.

Biohazard disposal – Discard all used catheters and sterile accessories in accordance with hospital procedures.

Catheter integrity – Use care to avoid damage to the catheter.

- Do not bend or kink the leading end of the catheter. Doing so could cause damage to the catheter lumen and make it unusable.
- Monitor the catheter throughout the procedure. If a flash is observed in the luer, replace the catheter immediately.

Do not resterilize – Do not resterilize this device for the purpose of reuse. Resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Electrode-electrode contact – Avoid contact between electrodes. Contact between electrodes may create a short circuit.

Embolism risk – Introducing any catheter or sheath into the circulatory system entails the risk of air, gas, or thromboembolism, which can occlude vessels and lead to tissue infarction with serious consequences.

- Avoid unnecessary catheter exchanges to minimize sheath-related embolic events.
- Always advance and withdraw components slowly to minimize the vacuum created and the risk of air embolism.
- Aspirate and flush the sheath frequently to help minimize the potential for embolic events resulting from the introduction of air or clot formation within the sheath.

Fluoroscopy use during catheter placement – The use of fluoroscopy during catheter ablation procedures presents the potential for significant x-ray exposure to both patients and laboratory staff. Extensive exposure may result in acute radiation injury and increased risk for somatic and genetic effects. Only perform catheter ablation after giving adequate attention to the potential radiation exposure associated with the procedure, and taking steps to minimize this exposure. Give careful consideration before using the device in pregnant women.

For single use only – The PFA catheter is intended only to be used once for a single patient. Do not reuse, reprocess, or resterilize the PFA catheter. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Handling and placing catheters – Careful manipulation of the catheter is necessary to avoid cardiac damage, perforation, or tamponade.

- Do not use excessive force to advance, withdraw, or rotate the catheter, especially if resistance is encountered. Excessive force may lead to catheter damage and blood loss.
- Use imaging guidance during catheter advancement, manipulation, and placement.
- Vascular perforation is an inherent risk of catheter placement.
- Performing ablation with the PFA catheter array inside the sheath may result in damage to the array or the sheath and should be avoided.
- Performing steering manipulation with the PFA catheter array inside the sheath may result in damage to the catheter steering mechanism or the sheath and should be avoided.
- The PFA generator is capable of delivering significant energy. Do not touch the ablation electrodes of the PFA catheter while operating the generator.
- If the system is to be tested outside of the body, the electrode array must be immersed in saline solution in a plastic container. Never test PFA delivery in direct contact with skin.

Imaging for catheter placement – Use of imaging during catheter manipulation and placement is strongly advised. Manipulating the catheter without imaging may result in damage to cardiac and vascular structures.

Improper connection – Do not connect the PFA catheter to anything other than the PulseSelect PFA generator or use it to deliver non-PFA energy. This may cause catheter malfunction or patient injury.

Other devices, wires, or catheters – Avoid catheter entanglement with other devices, wires, or catheters, for example, intracardiac echo catheters. Failure to do so may increase the risk of entrapment of the array or damage to the array, which may affect retrieval of the device into the transseptal sheath.

Product literature – Do not operate the PFA system or connect a catheter to the generator before completely reading and understanding the product literature for the system.

Phrenic nerve injury – To reduce the potential for phrenic nerve injury, assess for proximity of the ablation catheter to the nerve using an appropriate technique such as pacing for local phrenic nerve capture or using the test pulse feature before ablation. Stop ablation immediately if phrenic nerve impairment is observed and assess for injury.

Required use environment – Cardiac ablation procedures should be performed only in a fully equipped electrophysiology laboratory.

Safety in magnetic resonance imaging (MRI) not evaluated – The PFA catheter has not been evaluated for safety and compatibility in the magnetic resonance (MR) environment. The safety of the PFA catheter in the MR environment is unknown (such as any heating, migration, or image artifact). Scanning a patient during the use of this device may result in patient injury.

Serious incident – If a serious incident related to the device occurs, immediately report the incident to Medtronic and the applicable competent authority or regulatory body.

Sheath and guide wire required – Do not attempt to advance or withdraw the catheter through the vasculature without the use of a sheath and guide wire, as it may result in damage to cardiac and vascular structures.

Sterile package inspection – Visually inspect the sterile packaging and catheter before use. If the sterile packaging or the catheter is damaged, do not use the catheter and contact your local Medtronic representative.

Storage and transport – There is no special handling, including no special storage or transport conditions, required for this device. Standard storage conditions are sufficient to safeguard the device. Store the device in the original packaging at room temperature in a dry place.

Technical manual information – If you find information in this manual that is incorrect or illegible, contact your Medtronic representative or your local competent medical authority.

Use by date – Do not use the device after the “Use by” date on the package label.

Using PFA energy near implanted devices – Implanted devices, such as pacemakers and implantable cardioverter-defibrillators (ICDs), may be adversely affected by PFA energy.

- Keep external sources of pacing and defibrillation available during ablation.
- Program pacemaker sensing parameters to asynchronous pacing to ensure that PFA energy is not sensed as an intrinsic event.
- Deactivate ICD detection during the delivery of PFA energy.
- Perform complete implantable device testing before and after ablation.
- Monitor surface and intracardiac electrograms or vital signs during PFA energy delivery to assess for device interaction. Take appropriate action if any interaction is detected.
- Refer to the appropriate implantable device technical manual for additional information.

10 Adverse events

Potential adverse events associated with cardiac catheter ablation procedures include, but are not limited to, the following conditions:

- Access site complications (such as, bruising, ecchymosis, arteriovenous fistula, hematoma, pseudoaneurysm)
- Anemia
- Arrhythmias, proarrhythmia (such as, atrial flutter, bradycardia, heart block, tachycardia)
- Bleeding, possibly requiring transfusion
- Bruising
- Cardiopulmonary arrest
- Perforation of the heart or other organs during transseptal puncture or other procedures
- Cardiac tamponade
- Catheter entrapment in cardiac structures requiring intervention
- Cerebrovascular accident [such as stroke, transient ischemic attack (TIA)]
- Chest discomfort, pain, or pressure
- Collateral damage to the conduction system or coronary vasculature
- Cough
- Death
- Embolism
- Esophageal damage (including atrial esophageal fistula)
- Hemoptysis
- Hypotension
- Hypertension
- Infections (such as, sepsis)
- Myocardial infarction or ischemia
- Nerve injury or nerve damage (for example phrenic nerve injury)
- Pericarditis or endocarditis
- Pericardial effusion
- Pneumothorax
- Pulmonary edema
- Pulmonary vein dissection
- Pulmonary vein stenosis
- Radiation injury or damage and late malignancy
- Skin laceration or puncture
- Sore throat
- Unintended complete or incomplete atrioventricular node (AV-Node) or sinus node block or damage
- Valvular insufficiency or damage

11 Instructions for use

Caution: Inspect the sterile packaging and catheter before use. If the sterile packaging or catheter is damaged, do not use the catheter and contact your local Medtronic representative.

Caution: Verify that the catheter is within its expiration date. Do not use if the product use-by date has expired.

Notes:

- Check the activated clotting times (ACT) at 10 to 15 minute intervals until therapeutic anticoagulation is achieved and then at approximately 30 minute intervals during the procedure. The anticoagulation level is recommended to be maintained at an ACT of at least 350 seconds throughout the procedure.
- The PulseSelect PFA loop catheter is an over-the-wire catheter designed to be used with a 0.81 mm (0.032 in) J-tip guide wire at all times.

11.1 Prepare the PulseSelect PFA loop catheter and the transseptal sheath

Use a standard aseptic technique when removing the product from its sterile barrier and handling the catheter. During connection of the catheter, keep the connections dry at all times. Fluid incursion may lead to system malfunction.

1. Carefully remove the catheter from the packaging and inspect the catheter before use to verify that it has not been compromised during shipping or handling.
2. Use a sterile technique to connect the catheter interface cable to the connector located on the handle of the catheter (6, *Figure 1*).
3. Pass the other end of the catheter interface cable out of the sterile field and connect it to the front panel of the PulseSelect PFA generator. (For more detailed instructions, see the *PulseSelect Operator's Manual*).
4. Flush the guide wire lumen and insert the guide wire into the catheter via the guide wire luer. The guide wire must exit the proximal end of the catheter.
5. Prepare and manage a transseptal sheath with a minimum inner diameter of 3.33 mm (10 Fr) following the directions in the manufacturer's Instructions for Use (IFU).

 Notes:

 - Aspirate and flush the sheath while the sheath is in the right atrium.
 - After gaining access to the left atrium, maintain anticoagulation per hospital protocol.


11.2 Capture the array with the capture device

1. Ensure the catheter is in a straight position with the array in a neutral position.

 Warning: The array must be captured while it is submerged in heparinized saline to wet the array and reduce the possibility of air bubbles adhering to the electrode array during capture and catheter insertion.
2. Remove the capture device from the handle and slide it up to the spiral array.
3. Verify that the guide wire extends at least 10 cm (3.9 in) from the distal tip.
4. While holding the capture device against the spiral array, slowly advance the slide control forward, as shown in *Figure 2*, while advancing the capture device.


 Note: Do not advance the capture device beyond the end of the array.

 Figure 2. Advancing the capture device

5. As the capture device is advanced forward, continue to advance the slide control forward (see *Figure 3*).


 Figure 3. Advancing the slide control while advancing the capture device

6. Hold the slide control in the forward position while capturing the tip of the spiral array (see *Figure 4*).

 Note: Ensure the slide control is fully forward before capturing the distal tip to prevent kinking.

 Figure 4. Capturing the array

11.3 Insert the catheter and deploy the array


Warning: Avoid unnecessary catheter exchanges to minimize sheath-related embolic events. Catheter exchanges follow the same steps as catheter introductions.

1. Retract the guide wire to the distal tip of the catheter.
2. With the array contained, place the capture device against the hemostasis valve and advance the catheter array into the sheath.
3. Advance the array portion of the catheter into the sheath before replacing the capture device firmly onto the handle.

 Note: It is strongly advised to perform catheter manipulation under fluoroscopic guidance or other appropriate visualization techniques.
4. Advance approximately 10-20 cm of the guide wire length into the transseptal sheath.
5. Advance the catheter through the transseptal sheath until the guide wire and array enters the left atrium.
6. Before fully deploying the catheter, verify that the guide wire is positioned at least 5 cm past the ostium in 1 of the left pulmonary veins.
7. Continue to advance the catheter until the fourth electrode on the spiral array begins to exit the sheath.
8. Slowly pull back on the slide control while continuing to advance the catheter forward to deploy the spiral array.
9. To fully deflect the catheter, make sure the steering portion of the catheter is beyond the distal tip of the sheath. This will allow for full deflection of the distal segment of the catheter, as shown in *Figure 5*.

 Note: Do not attempt to steer the catheter within the sheath.

 Figure 5. Deflecting the distal segment of the array

11.4 Pace and detect

Note: It is strongly advised to perform the remaining steps under fluoroscopic guidance or other appropriate visualization techniques.

1. Advance the guide wire into the selected pulmonary vein to support and stabilize the catheter.
2. Advance the array into position for performing diagnostic EGM recording or bipolar pacing from the electrode pairs.
3. The slide control may be advanced during pulmonary vein diagnostic recording or pacing, but it should not be advanced when performing therapeutic energy deliveries.
4. Detect pulmonary vein potentials. The generator is required to be powered on to pace or visualize electrograms when using the catheter.

11.5 Deliver therapy

To use the catheter for a pulsed field ablation procedure, follow these steps. For more detailed instructions, see the *PulseSelect Operator's Manual*.

1. Select **START PROCEDURE** on the generator touchscreen display.

 Note: The therapy setup screen will display **1500 V** and a representation of the distal end of the 9-electrode **PFA Catheter**.
2. Confirm that the PFA Catheter is **Connected**. Connection of the catheter is indicated on the PulseSelect PFA generator touchscreen display.
3. The catheter array should first be placed to isolate each pulmonary vein ostium.
4. Deliver a test pulse.

 Notes:
 - The test pulse delivers a single sub-therapeutic stimulation pulse to the electrodes. A test pulse can indicate close proximity to the phrenic nerve via diaphragmatic stimulation. A delivery tone indicates when the R-wave synchronized test pulse has been delivered.
 - Deliver the test pulse while monitoring for stimulation of the diaphragm. The respiration cycle can move the catheter into or out of contact or proximity to the phrenic nerve.
 - If the test pulse produces diaphragm stimulation, reposition the catheter array.
5. Charge the system using the charge system button on the generator, optional remote, or the left actuator pad of the optional foot switch.
6. Deliver therapy using the deliver therapy button on the generator, optional remote, or the right actuator pad of the optional foot switch. Repeat applications for additional catheter placements as necessary to isolate the pulmonary vein ostium. One application is equivalent to four pulse trains delivered.

 Note: At any point the back button on the generator or optional remote can be used to return to the therapy setup screen.
7. Complete ablation of the associated pulmonary vein carina. Following isolation of the pulmonary vein ostium and carina, the catheter array should be repositioned antrally and four pulse trains delivered at each placement until an antral lesion pattern is achieved.

8. Repeat applications for additional catheter positions as necessary.
 Note: A minimum of 8 applications per vein is recommended to achieve ostial and antral ablation. Operators should use medical judgment to achieve a wide antral ablation for pulmonary vein isolation.

11.6 Pacing and navigation between pulmonary veins

1. Advance the guide wire to provide support for the catheter during maneuvers.
 Notes:
 - Do not attempt to steer the catheter within the sheath.
 - Use the steering lever tension control knob to maintain curve configuration.
2. Retract the catheter from the antrum of the vein.
3. Advance the slide control forward as required to detect electrical signals and check for pulmonary vein isolation.
4. Advance the spiral array into the vein while rotating it in a counter clockwise rotation.
5. Retract the slide control slightly to expand the spiral array against the inside surface of the vein.
 Note: Pacing is available through all electrodes on the catheter.
6. Repeat pulmonary vein isolation for additional pulmonary veins as necessary.

11.7 Removing the catheter

1. Disconnect the catheter from the catheter interface cable.
2. Verify that the steering lever is in a neutral position and at least 5 cm (1.97 in) of guide wire is inserted into one of the left pulmonary veins.
3. Slowly advance the slide control forward 1/4 of the travel distance and slowly pull back on the catheter shaft to retract the array into the transseptal sheath.
 Note: The catheter should only be pulled back to the right atrium through the transseptal puncture while fully inside the transseptal sheath.
4. Continue to slowly advance the slide control while simultaneously retracting the catheter, until the spiral array is completely inside the transseptal sheath.
5. Slowly withdraw the catheter halfway through the sheath.
6. Place the capture device against the hemostasis valve and continue retracting the catheter until the array is located in the capture device and the guide wire has been fully removed from the hemostasis valve.

11.8 Disposal

Do not dispose of this product in the unsorted municipal waste stream. Dispose of this product according to local regulations. See http://recycling.medtronic.com for instructions on proper disposal of this product.

Caution: Discard all used catheters and sterile components in accordance with hospital procedures.

12 Specifications

Overall length	145 cm (57.09 in)
Effective length	99 cm (38.98 in)
Deflection angle	170 degrees
Catheter reach at 90 degree deflection	23 mm
Array configuration	Single circled array [25 mm (0.98 in) outside diameter]
Number of electrodes	9
Electrode spacing	3.75 mm (0.15 in)
Electrode length	3 mm (0.12 in)
Single use circuitry	Circuitry allows the generator to assess “single use”
Outside diameter at array to catheter shaft bond	3.25 mm (9.8 Fr; 0.13 in)
Outside diameter of shaft	3 mm (9 Fr; 0.12 in)
Compatibility:	
Generator	Medtronic PulseSelect PFA generator
Catheter interface cable	Medtronic PulseSelect PFA catheter interface cable
Introducer sheath minimum	3.33 mm (10 Fr; 0.13 in)
Guide wire	0.81 mm (0.032 in) J-tip
Operating environmental parameters	15 °C to 30 °C (59 °F to 86 °F) at humidity between 30% and 75%
Sterilization	Single use, sterile device, e-beam radiation

13 Medtronic limited warranty

For complete warranty information, see the accompanying limited warranty document.

14 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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© 2024 Medtronic
Medtronic and Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.
Medtronic, Inc.
Medtronic B.V.