REACH SURGICAL | ENER REACH

ENER REACH Optimus

EN **Electrosurgical Instruments** Instructions

	Reach Surgical.Inc. 120 Xinxing Road, West Zone, TEDA, 300462 Tianjin, P.R. China

Overview

Please read the following information carefully before using this device. It provides important instructions regarding proper operation, potential risks, and potential damage to the product or individuals. In the event of any abnormal situation, follow the specified instructions to prevent harm to yourself or damage to the equipment. Failure to follow these instructions may result in compromised safety, performance, warranty, or maintenance, for which the manufacturer cannot be held liable. Please read the instructions for the Transducer and Ultrasonic Shears as well before using the device.
This document is designed to assist in using this device. It is not a reference for surgical techniques.

Standard Conventions Used: Caution, WARNING, and Note Statements

Please note the following statements, categorized as CAUTION, WARNING, or NOTE, which provide essential guidance for completing tasks safely and thoroughly. These statements can be found throughout the documentation and should be read before proceeding to the next step in a procedure.
WARNING: This statement highlights an operating or maintenance procedure, practice, or condition that, if not strictly followed, could lead to personal injury or loss of life.
CAUTION: This statement alerts the user to a potentially hazardous situation that, if not avoided, may result in minor or moderate injury to the user or patient, as well as damage to the equipment or other property. It may also serve as a warning against unsafe practices. This includes the necessary precautions for the safe and effective use of the Instrument and the care required to prevent damage resulting from proper or improper use.
NOTE: This statement indicates an operating practice or condition that is essential for executing a task efficiently.

Description

The ENER REACH Electrosurgical Instrument(herein referred as ‘the device’) is used in medical operating rooms for surgical procedures for cutting, coagulating human tissue, and ligating vessels. It features two separate slots: one for ultrasonic energy and the other for radio frequency waveform output.
In radio frequency mode, the generator delivers different energy schemas depending on the connected instrument/electrosurgery device. When an advanced bipolar instrument is connected, the generator delivers RF waveforms for sealing arterial and venous vessels, lymphatics, and tissue bundles up to 7 mm in diameter. When a basic bipolar instrument is connected, constant power energy is outputted for soft tissue cutting and coagulation.
In ultrasonic energy output mode, using the TRA6 transducer in combination with Disposable Ultrasonic Shears, including CH45PD, CH36PD, CH23PD, CH14PD, the device cuts soft tissues requiring bleeding control and minimal thermal damage, and seals vessels up to a maximum diameter of 5 mm. When connected with Disposable Ultrasonic Shears, including SRB14, SRB23, SRB36, SRB45, SRE14, SRE23, SRE36, and SRE45, the device cuts soft tissues requiring bleeding control and minimal thermal damage, and seals vessels up to a maximum diameter of 7 mm. When using the USC Ultrasound Converter in combination with Disposable Ultrasonic Shears, including SFA9 and SFA17, the device cuts soft tissues requiring bleeding

control and minimal thermal damage, and seals vessels up to a maximum diameter of 5 mm.

Intended Use

The device provides radiofrequency power to drive electrosurgical handpieces that are intended to seal vessels and cut, grasp, dissect tissues.

In addition, the generator provides ultrasonic power to drive ultrasonic surgical instruments that are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired.

Indications

The device provides radiofrequency power to drive electrosurgical handpieces that are used during open surgeries or laparoscopic surgeries in general, pediatric, gynecologic, urologic, thoracic surgery to cut and seal vessels up to and including 7mm, and to cut, grasp, and dissect tissues.
In addition, the generator provides power to drive ultrasonic surgical instruments that are indicated for soft tissue incisions when bleeding control and minimal thermal injury are desired. The ultrasonic surgical instruments can be used as an adjunct to or substitute for electrosurgery, lasers, and steel Shears for cutting and/or coagulating tissue in open surgeries or laparoscopic surgeries in general, pediatric, gynecologic, urologic, thoracic, and sealing and transection of lymphatic vessels.

Intended User

The device is intended for use by healthcare professionals for surgical applications.

Intended Use Environment

The device is intended to be used in a hospital.

Intended Patient Population

This device is suitable for patients aged 3 and older who require surgical procedures involving soft tissue incisions with bleeding control and minimal thermal injury.


Clinical Benefit

The device can be used safely and effectively in ligation and division of vessels, tissue bundles, and lymphatics.

Contraindications

This device is contraindicated for bone incisions, contraceptive tubal occlusion, and coagulation procedures. It is also not intended for use in neurosurgery.

System Composition

[01] Electrosurgery handpiece	[02] OP9 Generator
[03] Ultrasonic footswitch	[04] Electrosurgical footswitch
[05] Transducer	[06] Ultrasonic Shears
[07] Electrosurgery Bipolar Converter	[08] Ultrasound Converter

[01] Electrosurgical Handpiece

The electrosurgical handpiece, connected to a generator, is operated by the surgeon during the procedure. It applies RF electrosurgical energy to tissue between its jaws to coagulate or seal. The instrument also includes a blade for tissue division. Model: OSA23, OSA37, OSA44. Refer to the instructions for detailed information. Note: The maximum operating voltage of The Osprey instrument is 250Vpk, which is compatible with the OP9 generator.

[02] OP9 Generator

The generator provides both electrosurgical energy and ultrasonic therapy energy through separate instrument connection ports. It consists of a chassis, shell, electronic board, LCD touch screen, power switch,

and interfaces for accessory connection.
Model: OP9

[03] Ultrasonic Footswitch

Used to control the on/off output of ultrasonic energy with two switches ("MIN" and "MAX").
Model: OP-FSD

[04] Electrosurgical Footswitch

Used to control the on/off output of high-frequency energy with a single switch.
Model: OP-FSS

[05] Transducer

The transducer converts electrical energy from a compatible generator into mechanical motion for the instrument blades. It is a reusable instrument with a limited service life. The transducer is non-sterile and must be sterilized according to instructions before use. Refer to the TRA6 Transducer instructions for detailed information.
Model: TRA6.

[06] Ultrasonic Shears

The Ultrasonic Shears use mechanical motion from the transducer/ Ultrasound Converter and deliver ultrasonic energy for tissue cutting or coagulation. Caution should be taken as the mechanical vibration is not detectable and could unintentionally affect non-targeted areas. The Ultrasonic Shears are sterilized with ethylene oxide. If the shelf life is exceeded or the sterilization package is damaged, the Ultrasonic Shears should not be used and should be disposed of. Refer to the instructions of the Ultrasonic Shears for detailed information.
The following Ultrasonic Shears are compatible with the OP9 generator:
PD series: CH14PD, CH23PD, CH36PD, CH45PD
SRB series: SRB14, SRB23, SRB36, SRB45
SRE series: SRE14, SRE23, SRE36, SRE45
SFA series：SFA9,SFA17

[07] Electrosurgery Bipolar Connector

Used to connect basic bipolar electrosurgical instruments to the generator for tissue coagulation.
Model: OP-BPC


Note: The rated maximum operating voltage of the basic bipolar electrosurgical instruments must be

≥250Vpk, which is compatible with the OP9 generator. Please refer to the manual of the bipolar electrosurgical instruments

[08] Electrosurgery Bipolar Connector

The Converter is used for connecting generator and open shear（SFA series shear）
Model: USC

Generator Front Panel

[01] Electrosurgical socket

Used to connect advanced bipolar electrosurgical instruments or electrosurgery bipolar Converter.

[02] Transducer socket.

Connects the Transducer or Ultrasound Converter to the Generator.


[03] LCD touch screen

Displays system information and serves as the interface for adjusting controls and settings.

[04] Power switch Button

Press to turn on the Generator; press and hold to power off.

Generator Back Panel

[05] Electrosurgical Footswitch Socket

Round socket for connecting the Electrosurgical Footswitch

[06] Footswitch of Ultrasound Surgical Equipment Socket

Round socket for connecting the Footswitch of Ultrasound Surgical Equipment.

[07] Power Socket

Connects the Power Cord to the Generator.

[08] Potential Equalization Port

If the power socket's grounding is uncertain, this port allows for connection to protected earth.

[09] Expansion Interface

Used for function expansion and maintenance.

[10] Mounting Hole

Instructions for use

Refer to the provided guidelines for the operating environment.

Unpacking

Please follow the instructions upon receiving the components below.

- Check for any visible transportation damages. If any damage is found, please contact Reach Surgical, Inc or the local agent for assistance.

Components included in the table (For detailed technical specifications and product codes, refer to Chapter 'System Technical Conditions'):

Model	Description	Component
OP9	Generator	Generator, Power cord, Instructions
TRA6	Transducer	Transducer with cable
OP-FSD	Footswitch of Ultrasound Surgical Equipment	/
OP-FSS	Electrosurgical Footswitch	/
OP-BPC	Basic Bipolar Energy Connector	/
USC	Ultrasound Converter	/

Safety Precautions

- During equipment inspection, keep the distal end of the instrument away from other apparatuses, surgical drapes, the patient, or any other objects to avoid injury.
- Implement necessary safety measures in the presence of vapors, following hospital procedures and regulations.

Ultrasonic Energy with PD/SFA Series Ultrasonic Shears

- Connect the transducer, foot switch, and PD series Ultrasonic Shear to the generator or Connect the Ultrasound Converter, foot switch, and SFA series Ultrasonic Shear to the generator
- After passing the transducer and Ultrasonic Shear test, adjust the power level (1-5) using the +/- icons on the LCD screen.
- Press the 'Min' button on the Ultrasonic Shear or Min pedal on the Foot switch to activate the ultrasonic energy at the preset power level.
- Press the 'Max' button on the Ultrasonic Shear or Min pedal on Foot switch to deliver the maximum power level of energy.
- Release the key or foot switch to stop energy output.

Ultrasonic Energy with SRB/SRE Series Ultrasonic Shears

- After passing the transducer and Ultrasonic Shear test, adjust the power level (1-5) using the +/- icons on the LCD screen.
- Press the Energy Button on the instrument or Min button on the foot pedal to activate the ultrasonic energy at the preset power level.
- Press the Energy Button with Advanced Hemostasis on the instrument to activate Advanced Hemostasis mode.
- Release the Button on the Ultrasonic Shears or the button on the Foot switch to stop energy delivery.

Electrosurgical Bipolar Energy with advanced bipolar instruments

- The LCD display shows the advanced bipolar icon and adjustable hand control icon.
- Pull the Lever on the instrument or press the footswitch pedal to activate the radio frequency.
- The generator will stop energy output when the closure is completed, indicated by the complete icon on the screen.
- If the ligation is incomplete, a warning tone will be played, and a warning icon will be displayed.

Electrosurgical Bipolar Energy with Bipolar Electrosurgical Instrument

- The LCD screen displays the power value, adjustable icons, and power level bars.
- Power settings range from 1-95 watts, with adjustable increments.
- Select the desired power output setting: Low (15 watts), Medium (30 watts), or High (60 watts).
- Press the button on the hand switch or footswitch pedal to activate the radio frequency.
- Release the button to deactivate the energy output.

Setting Recommendations

Effect Settings	Power Settings Range	Clinical application	Optimized device
Low	1-15 watts	Such as nerves/spine, hands and Facial surgery	· Devices with smaller surface areas · Microtip forceps (0.4 – 2.2 mm)
Mid	16-40 watts	Such as coagulation during head/neck, spine and anatomy	· Devices with medium surface area · Microtip forceps (1.0 – 2.2 mm) · Small flat-head laparoscopic forceps · Bipolar scissors
High	45-95 watts (5 watts increments)	Such as head/neck and plastic surgery (similar to	· Devices with large surface areas · Large flat-head laparoscopic forceps

		medium effect, but with faster effect)	

NOTE: This device is not compatible with neutral electrode.


NOTE: The Ultrasonic Shears and RF electric Shear are patient contact applied parts

Using the generator


Turn on the device

- The system is ready for operation once it's been turned on. When the generator is connected to the mains supply and the standby switch light is on, the system is ready for use.
- After pressing the Standby Button, the following image will be displayed:


- If no Transducer/ Ultrasound Converter or Shear is connected to the Generator, or if they are connected improperly, the following prompt will be indicated:


Using Ultrasonic instrument (Transducer and Ultrasonic Shear)

When transducer/ Ultrasound Converter and Ultrasonic Shear are detected, the following image will be displayed.

NOTE: if the remaining use of the Transducer is less than 10, a warning message will be prompted in screen. User should pay attention to the number of remaining uses as the Transducer needs to be replaced when the number down to zero.


Next the following icon will be displayed, user will be asked to press any button on the Shear or open shear to start shear test with open jaws.


Press the any button, and system will start transducer & shear testing.

If the test fails, a test failed result with following icon will be shown in middle of screen.


When the transducer/ Ultrasound Converter and Ultrasonic Shear test passes, the following image will be displayed:

The following working screen will be displayed when connecting Shears without Advanced Hemostasis function or single button shears.


The following working screen will be displayed when connected new SRB or SRE Shears with Advanced Hemostasis.


Output mode will be highlighted it's activated.

Using Electrosurgical instrument


When an electrosurgical instrument is detected, the following image will be displayed.

When the electrosurgical instrument is ready for use, the following image will be displayed:

The following working screen will be displayed when connecting to advanced bipolar instrument.

The following working screen will be displayed when connecting to basic bipolar instrument.


Using Electrosurgical instrument and Ultrasonic Shears

Please note that when the generator is connected to both a bipolar electrosurgical instrument and an Ultrasonic Shear, the system operates in a split-screen mode to display the current state. The device that is prioritized will be the first to activate the energy. This split-screen state is shown below:

System Settings


To access the system settings, click on the settings icon located in the upper left corner of the screen. The following system setting options are available:


Language: Clicking on the Language item will display the language selection screen:

User Configuration: Clicking on the User Configuration item will display the user configuration screen:


System Preferences: Clicking on the System Preferences item will display the system preferences screen:

System Info: Gently click on the System Info item to display the system information screen:

Diagnostics: Clicking on the Diagnostics item will display the diagnostics screen:

To exit the system settings and return to the system standby mode, press the 'Return' icon located in the top left corner of the screen.

Shutting down the system

Follow the steps below to safely shut down the system:
Press the Standby Button: Locate the Standby Button on the Generator and press it to initiate the shutdown process.
Disconnect Transducer/ Ultrasound Converter and Ultrasonic Shear: Carefully disconnect the Transducer and the Ultrasonic Shear from the system. Refer to the operating instructions for proper handling of the Transducer and Ultrasonic Shear.
Unplug the power cable: Disconnect it from supply mains
Cleaning: Clean the Generator, Foot Switch, and Transducer/ Ultrasound Converter in accordance with the specified procedures outlined in Chapter ‘Maintenance’.

Troubleshooting

The Generator incorporates various alert signals and error codes to aid in the identification and diagnosis of component faults. It is important to note that these signals and codes are designed to support clinical judgment and observation, rather than replace them.

Sound alert

Working Status	Sound Type	Possible cause and action
Self-check Status	Normal sound	/
	Abnormal sound	A hardware fault has been detected. Please contact the manufacturer for repair assistance.
Using Ultrasonic Instrument	Max Level sound	/
	Min Level sound	/
	Abnormal sound	The Ultrasonic Shear has come into contact with excessive tissue. Reduce the tissue contact with the Ultrasonic Shear. If a continuous sound persists, carefully remove any tissue that may be accumulating around the end of the Ultrasonic Shear. There is a fault detected with the Transducer and/or Ultrasonic Shear. Please refer to the manufacturer or contact support for further assistance.
	ADVANCED sound	/
	ITS sound	/
Activating Electrosurgical Instrument	Advanced bipolar warning tone	/
	Normal bipolar sound	/
	Abnormal sound	Unsuccessful Cutting: If the tissue contacted by the instrument does not meet cutting requirements, reduce the amount of tissue in contact with the instrument. If the fault tone persists, carefully clear any tissue that may be accumulating at the end of the instrument. Hardware Failure: This error indicates a

		potential failure in the circuit or connector. In the event of any fault or error, it is recommended to consult the manufacturer or seek support for appropriate guidance and resolution.

Error Codes

The Energy Platform Generator is equipped with a comprehensive fault identification system consisting of alerts and system errors. When a fault is detected, the generator emits a warning tone, displays a warning signal on the control panel, and shows a corresponding fault code on the LCD screen. Follow the steps outlined below to address the issue:

Error Codes Table

Error code	**Corresponding fault message**
Warning	Please connect Ultrasonic Shear with Transducer correctly
Warning	Make sure jaws are open while testing
Warning	Ultrasonic Shear Error Detected
Warning	The remaining time of the transducer is zero
Warning	Please activate only one button at a time
Warning	Please reduce the force applied to the Ultrasonic Shear
Warning	Ultrasonic Shear button and footswitch cannot be closed simultaneously
Warning	Button stuck, please check and continue
Warning	Adjust jaws or clamp less tissue
Warning	Remove device from tissues
System Error	System Error
System Error	Self-check Failed

If an error appears on the screen during ultrasonic testing, perform the following actions:

- Ensure that the Transducer cable is fully inserted in the correct direction.
- Check if the Ultrasonic Shear has been tightened correctly or if any tissue has accumulated around the end of the Ultrasonic Shear. Adjust the Ultrasonic Shear's tightness and carefully remove any tissue accumulation around the Ultrasonic Shear casing. (If the test is initiated prior to the operation, ensure that the Ultrasonic Shear is pointing towards the air. If Ultrasonic Shears are being used, confirm that the clamping jaw is open and not in contact with any objects.)
- If the problem persists, consider replacing the Transducer or Ultrasonic Shear.
- Proceed to the equipment's working mode.

NOTE: The Transducer may not function properly if its temperature exceeds the specified limit. In such cases, use another Transducer immediately for recovery or follow the steps below to determine the cause of the error and explore optional recovery methods:

- Allow the Transducer to cool down at room temperature for a minimum of 45 minutes. This cooling method also applies if the Transducer becomes hot after prolonged operation at high power.
- If Transducer overheating is not evident and the problem remains unresolved, contact the manufacturer's maintenance representatives for assistance.
- Apart from fuses, there are no user-serviceable parts in the Generator. For any replacement or service requirements, please get in touch with service personnel who are trained and authorized by Reach Surgical, Inc. or your local representative.

Maintenance and upgrades of the Generator should be exclusively performed by service personnel trained and authorized by Reach Surgical.

Cybersecurity precautions should be considered to prevent potential threats. The following incidents pose cybersecurity risks:

- Unauthorized access to any non-related products for the device.
- Any unauthorized network communication with the device.
- Firmware or software upgrades that have not been authorized by Reach Surgical.

In the event of any of the above incidents, please contact the sales representative of Reach Surgical, Inc. or directly reach out to Reach Surgical, Inc. at Reachquality@reachsurgical.com.

Maintenance

Cleaning and disinfecting the generator and Basic Bipolar Energy Connector

Cleaning

Clean the Generator LCD screen and Basic bipolar/ **Ultrasound Converter** in accordance with hospital procedures and regulations. Before cleaning, ensure that the main power supply of the Generator is disconnected and the Power Cord is removed from the output unit.
WARNING: Cleaning procedures must be followed carefully to avoid damaging the Generator, causing electric shock, or creating a fire hazard. Do not spill or splash liquids on or into the Generator, or immerse it in liquid.

Follow these steps for cleaning:

- Prepare a neutral pH detergent or neutral pH enzyme detergent as specified by the detergent manufacturer.
- Using a clean, soft cloth soaked with a small amount of cleaning solution, manually wipe all surfaces, including the Generator screen.
- Wipe all surfaces with a clean, soft cloth soaked in warm tap water.
- Finally, wipe all surfaces with a clean, soft cloth to ensure they are dry.

Disinfection

- If the generator becomes contaminated with blood or bodily fluids, it must be disinfected before reuse. The following chemical disinfectants have been validated for use on the generator: 70% isopropyl alcohol, 6% sodium hypochlorite, 10% hydrogen peroxide。
- Follow the manufacturer's recommendations for proper use, concentration, and contact time of the disinfectants.
- Ensure that disinfectors are configured and used according to the manufacturer's instructions.

Cleaning the Foot Switch

Clean the Foot Switch and cable after each use using the following procedure:

- Disconnect the Foot Switch from the Generator.
- Prepare a neutral pH detergent or a neutral pH enzymatic detergent according to the manufacturer’s directions
- Use a soft, clean cloth lightly moistened with the cleaning solution to manually clean the surfaces. Pay special attention to cracks and crevices.
- Wipe thoroughly using a soft, clean cloth lightly moistened with warm tap water

- Dry with a soft, clean cloth.

CAUTION: Keep the cable and Foot Switch securely connected during rinsing.

- Wipe all surfaces with a clean, soft cloth.

WARNING: Do not use an ultrasonic cleaning machine to clean the Foot Switch.

WARNING: Do not switch on the power supply of the Generator before connecting the AC power cable. Ensure all connections are dry before assembly.

Cleaning and disinfecting the transducer

Transducers are shipped non-sterile and require thorough cleaning and sterilization before each use.

Cleaning:

Thoroughly clean Transducer according to the following steps:

- Rinse the Transducer with the screw side up and clean with a soft bristle brush with purified water until there are no smudges on the surface;
- The Transducer is soaked in a pH neutral enzymatic detergent (main ingredients: protease, lipase, amylase, cellulase, pectinase and other biological enzymes, environmentally friendly surfactant, rust prevention factor and stabilizer) for a period of up to 10 minutes at an appropriate temperature 15℃~65℃. The ratio of detergent and purified water is 1:400.
- Rinse the Transducer with the screw side up with purified water for 2 minutes;
- Clean the Connecting Screw, Shear Mount Surface and Connector with an alcohol wipe.
- Soak the Transducer in 75% medical alcohol and hold and shake it for 30 times;
- Rinse the Transducer with the screw side up with purified water for 2 minutes.

Note: The use of ultrasonic cleaners is not recommended for the Transducer.

Drying:

Drying temperature: 50~70℃, drying time: 30min.

Transducer Sterilization:

Following the cleaning and drying steps above, the Transducer must be sterilized by one of the methods listed below.

Steam Sterilization (121℃)

- Sheath should be installed before sterilization. Transducer should be wrapped during sterilization. Put the Transducer into a high-temperature steam sterilization pot for sterilization, with a temperature of 121℃ and a duration of 30min.
- Drying temperature: 50~70℃, drying time: 30min.

Steam Sterilization (134℃)

- Sheath should be installed before sterilization. Transducer should be wrapped during sterilization. Put the Transducer into a high-temperature steam sterilization pot for sterilization, with a temperature of 134℃ and a duration of 10min.
- Drying temperature: 50~70℃, drying time:30min.

Safety and Functional Tests

Ensure the implementation of safety and functional tests for the Transducer, Generator, and Foot Switch in accordance with hospital procedures and regulations. For safety and function tests of other components used by multiple patients, refer to the operating instructions specific to each component.

Safety Test

Generator: Certified hospital technicians should perform a leakage current test.
Foot Switch: Inspect the pedal, cable connector, and cable for any cracks or damage. Replace any damaged components.
Other components: Check all other components as instructed in their respective operating instructions.

Functional Test

Ultrasound Mode

- Prepare the complete set of PD Ultrasonic Shear and connect the Transducer following the instructions provided in Chapter II - Installation and Operation of the Equipment.
- Verify if it is possible to enter the working state. Different Ultrasonic Shears may have different entry interfaces. Refer to Part 1, Section 3 of Chapter 2 for detailed instructions.
- Confirm the display of MIN power Level 3 and MAX power Level 5.
- Press the power increase and decrease buttons to ensure that the MIN power level can be adjusted from levels 1 to 5.
- Power on the generator and switch it to ultrasonic working mode. Verify the correct connection of the transducer and Ultrasonic Shear.
- With the jaw open, press the "MAX" button on the foot switch. The LCD screen should display the MAX power level "5," and an activation tone should sound.
- With the jaw open, press the "MIN" button on the foot switch. The LCD screen should display the MIN power level, and an activation tone should sound.

WARNING: Before activating the system, ensure that the jaw is kept away from tissues, other instruments, or any other objects to prevent injury to the user.

Bipolar Electrosurgical Mode

- Connect the advanced bipolar Electrosurgical instrument/basic bipolar instrument according to the instructions.
- Check if the system can enter the working interface. Basic bipolar mode should display the 30 power level icon, while advanced bipolar mode should display the "advanced bipolar" icon.
- Lightly touch the power increment and decrement keys under Basic bipolar mode to confirm that the power level can be adjusted between 1 and 95. Touch the Low, Medium, and High bars to switch directly

between 15, 30, and 60.

- Power off the generator and wait for 5 seconds. Then, turn on the power supply of the generator and wait for 10 seconds. Check if the Basic bipolar mode displays the 30 power level interface, and the advanced bipolar mode displays the "advanced bipolar" interface. Verify if the generator is activated according to the predetermined requirements.
- Connect the advanced and Basic bipolar Electrosurgical instruments and press the single foot switch. Check for flashing power level indications on the control panel and listen for an activation sound.

WARNING: Before activating the system, it is strictly prohibited to allow any contact of the Ultrasonic Shear with tissues, other instruments, or any other objects to prevent injury to the user.

Warnings and Precautions

System related

- Read the instructions prior to use and follow hospital guidelines for clinical practice for ultrasonic surgery, electrosurgery, gynecology, and laparoscopy.
- Minimally invasive devices may vary from manufacturer to manufacturer. If minimally invasive instruments and accessories from different manufacturers are used in a surgery at the same time, check the compatibility of instruments and accessories before the surgery and check whether the accessories inserted into the human body have a rough surface, sharp edge, or protrusion that may cause safety hazards.
- This device is intended for use by trained and licensed surgeons only. Do not use electrosurgical devices unless you have been properly trained in their use for the specific procedure that you will need to complete. Untrained use of this device can cause unintended serious injury to the patient, including bowel perforation and unconscious and irreparable tissue necrosis.
- Do not open the generator enclosure without permission to avoid possible shock hazards. Any repair and upgrade of the instrument shall be performed by a service person trained and authorized by Reach Surgical, Inc. Do not use this instrument for any purpose other than medical surgery.
- A thorough understanding of the principles and techniques involved in laser, electrosurgical, and ultrasonic procedures is essential to prevent shock and burn hazards to both the patient and medical personnel and damage to this device or other medical devices. Electrical insulation or grounding must not be compromised. Do not immerse electrosurgical devices in fluids unless the design requires it, and labeling states that they should be immersed.
- Safe and effective ultrasonic surgery and electrosurgery depend not only on the design of the equipment but also largely on many factors controlled by the operator. To improve safety and effectiveness, read, understand, and follow the instructions for use provided with the device.
- As with all energy sources (electrosurgical, laser, or ultrasound), consideration should be given to the carcinogenic and infectious risks that many tissue byproducts, such as smoke and aerosols, may present. Appropriate precautions such as safety glasses, filtration masks, and effective smoke evacuation equipment should be observed in both open and endoscopic procedures.
- After removing the device, check the tissue for hemostasis. If hemostasis is absent, appropriate methods should be used to achieve hemostasis.
- Products manufactured or distributed by companies not authorized by Reach Surgical, Inc. may not be compatible with the device. Use of such products may lead to unexpected results and may injure the user or patient.
- To reduce the risk of interference, the device and the shell be connected to an independent power circuit.
- The equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. Potential for sparking due to collision with other metallic devices. Sparks may ignite flammable gases, such as field gases.
- The must operate within the required ambient operating temperature range.
- Output socket voltage shall meet the requirements of the generator (Chapter ‘System technical conditions’). If the power supply is not connected correctly, it may damage the generator and cause electric shock or fire hazards.

- Do not use extension cords to avoid fire hazards.
- Do not turn the activation tone to an inaudible level. Activation tones can be noticed by surgical team personnel while the generator is delivering energy.
- Smoke generated during electrosurgery has been shown to be potentially injurious to patients or surgical team personnel. Suggest adequate ventilation with a surgical smoke evacuator or other means.
- There are some components in the Ultrasonic Surgical Integrated Generator that are shipped non-sterile (e.g., Transducer). Sterilize the product as required before starting system installation. For cleaning and sterilization instructions, refer to each relevant instruction.
- To avoid injury to users or patients, the Ultrasonic Shear must avoid other devices, surgical drapes, patients, or other objects before pressing the test button and during system check. Safety measures in case of aerosol (according to hospital regulations) shall be implemented in the system inspection and test method.
- Do not apply too much pressure to the jaw to avoid inhibiting the delivery of ultrasonic energy.
- To avoid injury to the user, the blade must avoid contact with tissue, other devices, or other objects before activating the system.
- If liquid is sputtered or poured on or into the generator, or the generator is spilled or poured into the liquid, it may damage the generator and cause electric shock or fire hazards.
- Sparking and heating associated with vessel closure techniques can both serve as sources of ignition. Gauze and a sponge should remain moist. Keep electrosurgical electrodes away from combustible materials and oxygen-rich (O2) environments.
- When there is significant damage to the Transducer or if any parts show signs of damage after cleaning and disinfection maintenance, discard them. Damaged parts are clearly marked to avoid misuse prior to subsequent handling.
- Disposable waste and electronic waste shall be disposed according to hospital regulations and shall not be discarded at will to avoid environmental pollution.
- Avoid using the generator close to or stacked on other equipment. If adjacent or stacked use is necessary, monitor the generator and other equipment to ensure proper operation.
- The device does not contain any operator-serviceable parts. For service, contact your Reach Surgical sales representative or service personnel.
- Check all devices connected to the system and connections prior to use. Validate that the device performs as intended. Improper connection can lead to arcing, sparking, device malfunction, or unintended surgical results.
- To avoid the risk of electric shock, this equipment must only be connected to a supply with protective earth.
- Keep as far as possible between the generator and other electronic devices (e.g., monitors). Do not cross or tie electrical wires to electronic devices. The generator may cause interference with other electronic devices.
- Use the system with caution in the presence of internal or external pacemakers or other implanted Devices(IEDs). Interference produced by electrosurgical equipment can cause a pacemaker or other device to enter an unsafe mode or permanently damage the device. Consult the device manufacturer or responsible hospital department for further information when use is planned in patients with implanted medical devices.
- Use caution if stacking instruments on top of the generator or placing the generator on top of electrical instruments. This is an unstable configuration and does not provide adequate cooling.

- If the generator fails, it may cause surgical interruption. A backup system shall be available.
- If required by local regulations, the generator should be connected to the hospital's equipotential connector using an equipotential cable.
- When the system and physiological monitoring equipment are used simultaneously on the same patient, any monitoring electrode shall be placed as far as possible from the instrument.
- In basic bipolar mode, choose the lowest possible output power for the desired effect.
- The failure of high-frequency surgical equipment may cause an unexpected increase in output power.
- The generator and basic bipolar electrosurgical connector may be invaded by water or particulate matter. In the process of use and cleaning, it is necessary to avoid the invasion of water or particulate matter.
- When the system is used in combination with an endoscope, it may increase the leakage current on the patient body. Pay attention during the procedure.
- A notice to the user and/or patient that any serious incident that has occurred in relation to the device sho uld be reported to Reach Surgical, Inc. through Reachquality@reachsurgical.com. and the competent aut hority of the Member State in which the user and/or patient is established.
- To avoid unknown risks ,do not use non-standard accessories to connect to generators.
- OP9 generator and ancillary devices should not be used when the physician believes that high-frequency bipolar or ultrasound surgery may be contrary to the best interests of the patient
- For contraindications when used with high-frequency bipolar devices or ultrasonic tool heads, please refer to the instructions for Use of each device (IFU).
-

Instrument related

During surgical procedures, it is important to be aware that the end portion of the instrument Ultrasonic Shear, Ultrasonic Shears pad, and shaft may experience elevated temperatures when ultrasound or bipolar electrocautery is applied to tissue for extended periods. To ensure safety, it is crucial to avoid any accidental contact between these instrument components and tissue, surgical drapes, surgical coats, or any unrelated areas throughout the procedure.
Furthermore, it is vital to exercise caution when positioning the bipolar instrument. Specifically, it should not be placed near or in contact with combustible materials such as gauze or surgical drapes. The bipolar instruments used in this context can generate heat during use, which could potentially lead to a fire hazard. When the bipolar electrosurgical generator is not in use, it should be properly stored in the protective sleeve or kept at a safe distance from the patient, surgical team members, and any combustible materials.
NOTE: Please refer to the respective Instructions for Use for additional warnings and precautions.

EMC Information

The product has passed the electromagnetic compatibility test, which meets the limitation requirements of IEC60601-1-2 standard for medical equipment. These restrictions provide reasonable protection against harmful interference in normal medical installations.

Product composition

Serial Number	Part Name	Model/Version No.	Remarks
1	Electrosurgical Instrument	OP9	/
2	Transducer of Ultrasound Surgical Equipment	TRA6	Compatible use
3	Ultrasonic Shear System for Single Use Ultrasonic Shear	CH14PD, CH23PD, CH36PD, CH45PD, SRB14, SRB23, SRB36, SRB45, SRE14, SRE23, SRE36, SRE45	Compatible use
4	Disposable ultrasonic high-frequency surgery unit Tissue Sealer	OSA23, OSA37, OSA44	/
5	Ultrasonic footswitch	OP-FSD	/
6	Electrosurgical footswitch	OP-FSS	/
7	Electrosurgery Bipolar Connector	OP-BPC	/
8	Ultrasound Converter	USC	/

Note：In addition to the accessories provided by our company, the use of other manufacturer's accessories may result in the degradation of the EMC performance of the ultrasonic high-frequency surgery integrated surgical system

Product cable

Serial Number	Cable Name	Length (m)	Shielded
1	Ultrasonic footswitch Cable	3	No
2	Electrosurgical footswitch Cable	3	No

3	Power cord	5	No
4	Osprey Bipolar Electrosurgical Cable	3	Yes
5	The cable of Transducer	2.9	Yes

EMC performance

This equipment may be subject to radio frequency interference caused by other medical equipment and radio communications. To prevent such interference, this product has been tested according to IEC 60601-1-2 and meets its requirements. However, Reach Surgical, Inc. does not guarantee that there will be absolutely no interference in individual installation environments.
If it is found that the device is interfered (which can be determined by turning the device on and off), the user (or maintenance personnel approved by Reach Surgical, Inc.) should try to take one or more of the following measures to solve the interference problem:
Adjust the direction or position of the device that affects it.
Increase the distance between this device and the sending device.
Use other power sources (rather than the power used to affect the equipment) to power this equipment.
Consult the supplier or service representative for other suggestions.
The manufacturer is not responsible for any interference caused by the following situations: use other interconnecting cables other than the recommended cables; alter or modify this equipment without permission. Unauthorized changes or modifications may cause the equipment lose efficacy.
All types of electronic equipment may cause electromagnetic interference to other equipment through the air or other cables connected to it. Do not use devices that can emit RF signals, such as cellular phones, radio transceivers, or radio control products, near this device, as this may cause the performance of this device to fail to meet the specified specifications. When such devices are close to this device, turn off the power of these devices. The medical personnel in charge of this equipment should instruct technicians, patients and other personnel who may be close to this equipment to fully comply with the above requirements.
To fully achieve the specified EMC performance, the user should install the product correctly according to the steps described in the manual. If there are any EMC-related problems, please contact the maintenance personnel approved by Reach Surgical, Inc.
The Transducer (with cable) and Ultrasonic Shears are defined as the applied part of the whole system.

Precautions for product installation

The equipment can be used in a hospital environment but does not include radio frequency shielding rooms around active Radio frequency surgical equipment or where magnetic resonance impact equipment is placed, because the electromagnetic disturbance intensity in these locations is high.
Separation distance and impact of fixed radio communication equipment: magnetic field strength generated by fixed transmitters, such as base stations of wireless (cellular/cordless) telephones, land mobile radio receivers, amateur radio receivers, AM and FM radio broadcasts, and TV broadcasts Generators, etc., cannot be accurately measured theoretically. To assess the electromagnetic environment generated by fixed RF

transmitters, measurement of the electromagnetic field should be considered. If the measured value of the magnetic field strength at the location of the device exceeds the corresponding radio frequency level specified in the "Anti-Interference Statement", the device should be inspected to ensure that it can operate normally. If abnormal operating conditions are found, additional measurements should be considered, such as reorienting or relocating the equipment, or using an anti-radio frequency room.

1）Use the Power Cord provided or designated by the Reach Surgical, Inc. Products equipped with a power plug should be plugged into a fixed power outlet with protective grounding. Do not use any type of adapter or connector to connect the power plug.

2）Keep this device away from other electronic devices as much as possible.

3）Follow the steps to connect the device.

General notes

(1) The specification of the cable.

The use of cables provided by the Reach Surgical, Inc will not damage the EMC performance of this product. If unspecified cables are used, the EMC performance of this equipment may be significantly reduced.

(2) Precautions for unauthorized modifications

The user shall not modify this product, otherwise the EMC performance of this product may decrease.

The modification of the product includes the following changes:

a. Cable (length, material, and wiring, etc.).

b. Equipment installation/layout.

c. Equipment configuration/components.

d. Equipment protection parts (cover opening/closing and cover fixing parts).

(3) All protective covers should be closed when operating the equipment.

This product is expected to be used in the electromagnetic environment specified below, and the purchaser and user of this product should ensure that it is used in this electromagnetic environment.

Essential Performance

Ultrasonic Basic performance: None

High-frequency Basic performance:

1. For advanced high-frequency Radio energy output, the system shall ensure that the deviation of rated power output does not exceed 20% of the standard value

2. For basic high-frequency Radio energy output, the system shall ensure that the deviation of rated power output does not exceed 20% of the standard value

Guidance and MANUFACTURER'S declaration - ELECTROMAGNETIC EMISSIONS		
The device is intended for use in the electromagnetic environment specified as follows. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance

RF emissions CISPR 11	Group1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class A	The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals
Harmonic current IEC 61000-3-2	Class A	
Voltage changes/Voltage fluctuations/ flicker IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment Specified as follows. The customer or the user of the device should assure that it is used in such an environment			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2,4,8 ,15kV air	± 8 kV contact ± 2,4,8 ,15kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, The relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line (s) ± 2 kV line (s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Ultrasound surgical Equipment requires continued operation during power mains interruptions, it is recommended that the Ultrasound surgical Equipment be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at level characteristic. Of a typical location in a typical commercial or hospital environment.

<table>
<tr><td colspan="4">Guidance and manufacturer's declaration - electromagnetic immunity</td></tr>
<tr><td colspan="4">The device aims at application under the electromagnetic environment specified as follows.
Customer or user of the device should assure that it is used under such environment</td></tr>
<tr><td>IMMUNITY
test</td><td>IEC 60601 test
level</td><td>Compliance
level</td><td>Electromagnetic environment-
guidance</td></tr>
<tr><td>Conducted RF
IEC 61000-4-6</td><td>3 Vrms
150 kHz to 80
MHz
6 V in ISM bands
between 0.15
MHz and
80 MHz</td><td>3V</td><td>Portable and mobile RF communications
Equipment should be used no closer to any part of The device , including cables, than recommended separation distance calculated from the equation for frequency of the transmitter.
Recommended separation distance
$d = [\frac{3.5}{V1}]\sqrt{P}$</td></tr>
<tr><td>Radiated RF
IEC 61000-4-3</td><td>3 V/m
80 MHz to 2.7
GHz
80 % AM at 1 kHz</td><td>3V/m</td><td>$d = [\frac{3.5}{E1}]\sqrt{P}$ 80MHz to 800MHz
$d = [\frac{7}{E1}]\sqrt{P}$ 800MHz to 2.7GHz
where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey.3 should be less than the compliance level in each frequency range.
Interference may occur in the vicinity of equipment marked with the following symbol:

The ENCLOSURE PORT of ME EQUIPMENT and ME SYSTEMS shall be tested as specified in Table 9 of IEC 60601-1-2 using the test methods specified in IEC 61000-4-3.</td></tr>
<tr><td colspan="4">Separation distances recommended between portable and mobile RF communications equipment and the devices.</td></tr>
<tr><td colspan="4">The device aims at application under an electromagnetic environment in which radiated RF disturbances are controlled. Customer or user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommendation as follows., according to the maximum output power of the communications equipment.</td></tr>
<tr><td>Rated maximum</td><td colspan="3">Separation distance according to frequency of transmitter /m</td></tr>
</table>

output power of transmitter (w)	150 kHz to 80 MHz $d = [\frac{3.5}{\mathrm{V1}}]\sqrt{\mathrm{P}}$	80 MHz to 800 MHz $d = [\frac{3.5}{\mathrm{E1}}]\sqrt{\mathrm{P}}$	800 MHz to 2.7 GHz $d = [\frac{7}{\mathrm{E1}}]\sqrt{\mathrm{P}}$
0.01	0.117	0.117	0.233
0.1	0.36999	0.36999	0.73681
1	1.17	1.17	2.33
10	3.69986	3.69986	7.36811
100	11.7	11.7	23.3

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Technical Conditions

Components required for system operation: Electrosurgical Instrument, Transducer, Ultrasonic footswitch, Electrosurgical footswitch, Ultrasonic Shear or Electrosurgery Handpieces Electrosurgery Bipolar Connector Bipolar Instrument.

Refer to the product description for this component.


Degree of protection against electric shock	Ultrasound:  Type CF Applied Part; Bipolar Electrosurgical  Defibrillation-proof Type CF Applied Part
Degree of Protection against Electric Shock	Class I
Generator to harmful infusions	Normal equipment
Degree of Protection against Harmful Ingress of Water	Footswitch: IP68
Management Category	Class IIb
Ultrasound Mode Parameter Requirements:	
Input Power	Supply voltage: 100-240V
	Supply frequency: 50Hz/60Hz Input power: 400VA
Ultrasound Mode Parameter Requirements:	
Excitation frequency	54 kHz - 57 kHz
Primary tip vibration excursion	25 μm ~ 110 μm
Ultrasonic Shear tip main acoustic output surface Product:	1.53～2.75 mm2
Secondary Crosslink Acoustic at Ultrasonic Shears Tip Output Area:	19.32～35.88 mm2
Derived Output Sound at Ultrasonic Shear Tip at Reference Tip Primary Amplitude Power	< 30W
Type of system frequency control	Independent of load, the excitation frequency is automatically adjusted continuously.
Power reserve index	Not less than 2.5
Electrosurgical Parameters Requirements:	
Working frequency	400KHz ± 5%
Maximum output voltage	⩽ 250Vpk

Advanced Bipolar mode Maximum output current	⩽ 5.5Arms
Basic bipolar mode Maximum output current	⩽ 2.2Arms
Maximum output peak-to-peak voltage	⩽ 500Vpk
Peak Factor	1.6 ± 0.4
Rated Power	Radio frequency advanced bipolar: 270 W± 20% Radio frequency Basic bipolar: 95W± 20%
Rated load	Radio frequency advanced bipolar: 30Ω Radio frequency Basic bipolar: 200Ω

Operating Environment Conditions	Temperature: 10 ° C to 30 ° C	After generator is moved from a storage environment lower than 10°C or than 40°C to a running environment, stand the generator for one hour before starting it
	Relative humidity: 35%-75%	
	Atmospheric pressure range: 800 hPa to 1060 hPa	
Transportation and Storage Conditions	Temperature: -30 ° C to + 55 ° C (generator, foot switch, adapter) Temperature: -10 ° C to + 55 ° C (Ultrasonic Shear)	
	Humidity: ⩽ 80%	
	Atmospheric pressure range: 800 hPa to 1060 hPa	
Date of manufacture	The manufacturing date can be determined by the serial number on the back panel of the generator.	
Power cord	Compliance with CCC certification requirements	
	Current rating: 10A	
Persistence Rate	Determined by transducer handpiece and Ultrasonic Shear used. For Persistence Rate information, refer to applicable Ultrasonic Shears and Transducer Handpiece Instructions or Chapter 7 – Warnings and Precautions.	
Fuse	¢5 * 20 T8AH250V	
Weight (without packaging)	Generator: nominal 8 kg	
Total Volume	OP9 generator: (length * width * height): 34 cm-34 cm-16 cm	
Disposition	Some internal components of generator, foot switch and foot switch cable contain lead. According to local Requirement and regulation for disposal. Dispose of batteries in accordance with appropriate waste disposal practices.	
AP/APG Classification	Not AP/APG equipment.	
Service life:	Service life: 7 years	
Software Release Version	V01.01	

Bipolar output waveform

Advanced bipolar output power-load curve


POWER CURVE LIMIT Max Output: 270[W]			
Load (Ω)	Lower Limit	Nominal	Upper Limit
5	83.6	104.5	125.4
10	159.84	199.8	239.76
20	216	270	324
30	216	270	324
50	204.64	255.8	306.96
100	138.72	173.4	208.08
200	78.88	98.6	118.32
500	31.6	39.5	47.4
1000	15.84	19.8	23.76

Nominal power curve of 270 watts output at rated load in accordance with the current/voltage limit of the power curve meter. The power curve represents the envelope that varies under operating condition.


Note: Advanced bipolar output power-load curve was tested in engineering mode of generator.

Basic bipolar output full power-load curve


Full Power Curve Limit [W]: Max Power 95 Watts			
Load (Ω)	Lower Limit	Nominal	Upper Limit
5	16.16	20.2	24.24
10	32.32	40.4	48.48
20	64.48	80.6	96.72
30	76	95	114
50	75.52	94.4	113.28
100	75.68	94.6	113.52
200	76	95	114
500	37.12	46.4	55.68
1000	19.76	24.7	29.64

Basic bipolar output half power-load curve

Half Power Curve Limit [W]: Max Power 50 Watts			
Load (Ω)	Lower Limit	Nominal	Upper Limit
5	16.16	20.2	24.24
10	32.32	40.4	48.48
20	40	50	60
30	40	50	60
50	39.68	49.6	59.52
100	39.84	49.8	59.76
200	39.92	49.9	59.88
500	20.48	25.6	30.72
1000	10.24	12.8	15.36

Output Power vs. Power Settings for Power in Normal Bipolar

Output Power vs. Power Set Value for Normal Bipolar [W]; Load: 200Ω			
Set Power (W)	Output Power Lower Limit (W)	Nominal Power (W)	Upper Output Power Limit (W)
5	1	5	10
10	5	10	15
15	10	15	20
20	15	20	25
30	24	30	36
45	36	45	54
60	48	60	72
75	60	75	90
95	76	95	114

Service and Warranty

Reach Surgical, Inc. warrants this product to be free from defects in material and workmanship under normal use and preventive maintenance for the warranty period specified below. Reach Surgical's obligation under this warranty is limited to repairing or replacing, at its option, any defective product or part that has been returned to Reach Surgical, Inc. or its authorized Distributor within the applicable warranty period and is found to be defective to Reach Surgical's satisfaction. This warranty does not apply to products or parts that have been:





- Adversely affected due to use with unauthorized devices manufactured or distributed by parties not authorized by Reach Surgical, Inc.
- Repaired or altered outside Reach Surgical's factory, if it affects the stability or reliability of the device as determined by Reach Surgical.
- Subjected to improper use, negligence, or accident.
- Used in a manner inconsistent with the design, use parameters, instructions, and guidelines for the product or with industry-accepted functional, operational, or environmental standards for similar products.



Warranty Periods

Basic Bipolar Energy Converter (OP-BPC)/ Ultrasound Converter(USC): 1 year for components and labor.
Generator (OP9): 1 year for components and labor.
Foot Switch/Power cord: 1 year for components and labor.

This warranty is the exclusive remedy for the original purchaser and replaces all other warranties, express or implied, including warranties of merchantability and fitness for a particular purpose. Reach Surgical, Inc. shall not be liable for any special, incidental, or consequential damages, including damages resulting from loss of use, profits, business, or goodwill, except as expressly provided by applicable law.
Reach Surgical, Inc. does not authorize any person to assume any additional liability in connection with the sale or use of its products. There are no warranties that extend beyond the terms stated herein.
Reach Surgical, Inc. reserves the right to make changes to its products without incurring any obligation to retroactively apply those changes to previously sold or built products.

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
EN/English

Overview

Please read the following information carefully before using this instrument. It provides important instructions regarding proper operation, potential risks, and potential damage to the product or individuals. In the event of any abnormal situation, follow the specified instructions to prevent harm to yourself or damage to the equipment. Failure to follow these instructions may result in compromised safety, performance, warranty, or maintenance, for which the manufacturer cannot be held liable.
This document is designed to assist in using this instrument. It is not a reference for surgical techniques.

Nomenclature (Illustration 01)

[01] Jaws	[02] Shaft
[03] Rotation Knob	[04] Handle
[05] Lever	[06] Cutting Trigger
[07] Cable	[08] Connector

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Description

The ENER REACH Enerseal Electrosurgery Handpieces (hereinafter referred to as the **instrument**) are designed for use with the ENER REACH Optimus Electrosurgical Instruments OP9 generator, which includes vessel sealing capability. For more detailed information, refer to the generator's instructions or contact the relevant department at Reach Surgical, Inc.
These instructions assume that the operator possesses knowledge and proficiency in correctly setting up and

operating the associated Reach Surgical, Inc. generator. For setup instructions and additional warnings and precautions, please consult the user's guide of the generator.
The instrument functions by utilizing radiofrequency (RF) electrosurgical energy to create a seal on vascular structures such as vessels and lymphatics, as well as tissue bundles positioned between the instrument's jaws. The surgeon can activate a blade within the instrument to divide tissue as needed. Shorter shaft lengths are typically used for open procedures, while longer shaft lengths are typically employed for laparoscopic procedures.

Instrument Product Codes

Product Code	Shaft Length (cm)	Diameter (mm)
OSA23	23	Φ5
OSA37	37	Φ5
OSA44	44	Φ5

Intended Use

The instrument is used in medical institutions for cutting and coagulation of tissues during surgery.

Indications

The instrument is used in medical institutions for cutting and coagulation of tissues during open surgeries or laparoscopic surgeries in general, pediatric, gynecologic, urologic, thoracic surgery. It can be used to close arteriovenous vessels and lymphatic vessels with a diameter of no more than 7mm and tissue bundles suitable for the size of instrument jaws.

Contraindications

The instrument is contraindicated for bone incisions, contraceptive tubal occlusion, and coagulation procedures. It is also not intended for use in neurosurgery.

Intended User

The instrument is intended for use by healthcare professionals for surgical applications.

Intended Use Environment

The instrument is intended to be used in a hospital.

Intended Patient Population and medical condition to be treated

Patients aged 3 and older who need surgery in which soft tissue incisions with bleeding control and minimal thermal injury are required.

Clinical Benefits

The instrument can be used safely and effectively in ligation and division of vessels, tissue bundles, and lymphatics.

Instructions for use

Unpacking and connection

- Remove instrument from tray by firmly pulling on the Handle. Do not pull on the instrument's Jaws or Cable.
- Insert the Connector into the Electrosurgical socket on the generator. Follow the instructions in the generator user's guide to complete the setup procedure.

Tissue Manipulation and Dissection

- The instrument can be used to manipulate and dissect tissue with the jaws either open or closed.

Rotating the jaws

- Spin the Rotation Knob to rotate the jaws.

CAUTION: Do not turn the Rotation Knob when the Lever is fully pulled as this may damage the instrument

Closing and opening the jaws

- To close the jaws, place the tissue between the jaws and pull the Lever. The first click indicates the jaws are fully closed.
- To open the jaws, release the Lever.

CAUTION: Do not turn the Rotation Knob when the Lever is fully pulled as this may damage the instrument.
CAUTION: Continue pulling the Lever after the first click will activate energy and start sealing.

Sealing and cutting with Hand Activation

- After the jaws are fully closed and the first click heard, pull the Lever further until a second click heard, this will activate energy delivery, hold until the seal cycle is complete.
- Open the jaws to release tissue or use the Cutting Trigger to cut the tissue.

NOTE: To use this mode, ensure hand-activation is enabled on the generator before use.
NOTE: The user should inspect the seal before cutting the vessel or tissue.
NOTE: To apply a double seal, clamp overlap the edge of the existing seal. The second seal should be distal to the first seal to increase seal margin.
NOTE: A tone with multiple pulses indicates that the seal cycle was not completed. Refer to the Troubleshooting section possible causes and corrective actions. Do not cut tissue until you have verified that there is an adequate seal.
CAUTION: Pull the Lever firmly until the sealing is complete. The Lever does not latch onto the activation position.
CAUTION: A continuous tone sounds to indicate the activation of RF energy. When the activation cycle is

complete, there will be a two-pulsed Seal-Cycle-Complete tone sounds to indicate RF output stops.

Sealing with Footswitch Activation

- After the jaws are fully closed and the first click heard.
- Pull the Lever further until a second click is heard, while keep pulling the Lever, press and hold the footswitch pedal to activate energy delivery until the seal cycle is complete.
- Release from the footswitch, and release the Lever to open jaws, or use the Cutting Trigger to cut the tissue.

NOTE: To use this mode, ensure hand-activation is disabled on the generator before use.
NOTE: The user should inspect the seal before cutting the vessel or tissue.
NOTE: To apply a double seal, clamp overlap the edge of the existing seal. The second seal should be distal to the first seal to increase seal margin.
NOTE: A tone with multiple pulses indicates that the seal cycle was not completed. Refer to the Troubleshooting section possible causes and corrective actions. Do not cut tissue until you have verified that there is an adequate seal.
CAUTION: Pull the Lever firmly until the sealing is complete. The Lever does not latch onto the activation position.
CAUTION: A continuous tone sounds to indicate the activation of RF energy. When the activation cycle is complete, there will be a two-pulsed Seal-Cycle-Complete tone sounds to indicate RF output stops.
WARNING: Activating energy delivery with a footswitch with jaws not fully closed may result in improper sealing and increase thermal spread to tissue outside the application site. Make sure the jaws are fully closed during energy delivery.

Cutting Tissue

- To cut the tissue, first pull the Lever, make sure the first click heard and the jaws are fully closed.
- Pull the Cutting Trigger activate the cutting mechanism.
- Release the Cutting Trigger to retract the blade.
- Release the Lever to open jaws.

WARNING: Energy-based devices, such as electrosurgical pencils or ultrasonic scalpels which generates thermal spread should not be used to transect tissues sealed by this instrument.
WARNING: Make sure the seal cycle is completed before cutting.
CAUTION: Do not use the cutting mechanism if clips, staples, or other metal objects are located between the jaws, as this could damage the blade.
CAUTION: Using excessive force to pull the Lever can result in inadvertently activate of energy delivery.
CAUTION: If the cutting trigger does not automatically retract to its original position, push the cutting trigger to reset it.

Cleaning the instrument

- The jaws should be kept clean for optimum performance. Build-up of eschar may impact the seal and/or cutting effectiveness. Wipe jaw surfaces and edges with a wet gauze pad and remove any residue tissue from blade track and jaw hinge area.

WARNING: Inspect the jaws prior to cleaning, make sure that the blade is not deployed.

WARNING: Do not pull the Lever or the Cutting Trigger while cleaning the jaws, as this may cause injury to the user.

CAUTION: Do not pull the Lever with gauze between the jaws as this may damage the instrument.

CAUTION: Do not clean the instrument scratch pad or other abrasives.

Post procedure

- Dispose of the Instrument and its accessories into appropriate containers. Do not reuse Instrument.

Warnings and Precautions

- The instrument is sterile, and single use only. Discard it after use. Do not sterilize or reuse. Attempts to clean sterilize these devices without appropriate regulatory authorization may result in bio-incompatibility, infection, or product failure risks to the patient.
- This instrument is intended for use ONLY with the compatible generators from Reach Surgical. Use of this instrument with other generators may not result in the desired tissue effect, may result in injury to the patient or the user, or may cause damage to the instrument.
- Do not use this instrument unless properly trained to use it in the specific procedure being performed. Use of this equipment without such training may result in serious unintended patient injury.
- Exercise caution when using the system in the presence of internal or external pacemakers or other implanted devices. Interference produced by electrosurgical equipment can cause a pacemaker or other device to enter an unsafe mode or sustain permanent damage. Consult the device manufacturer or the responsible hospital department for further information when planning to use the instrument in patients with implanted medical devices.
- When using this instrument with an energized endoscope, the leakage current from the instrument and the endoscope are additive. The patient may be exposed to unexpected levels of leakage current if this instrument is used with an energized endoscope that is not a type CF applied part.
- Before inserting the instrument through the cannula, inspect its outer surfaces to ensure there are no rough or sharp edges that could damage tissue.
- Contact between an active instrument electrode and any metal objects (hemostats, staples, clips, retractors, etc.) may increase current flow and result in unintended surgical effects or insufficient energy deposition.
- The safe and effective use of RF energy depends on various factors solely under the control of the operator. It is important to read, understand, and follow the operating instructions provided.
- This Instrument is supplied sterile, EO sterilized. If sterile package is damaged, DO NOT use.
- Verify the instrument model is correct before the procedure.
- Electric Shock Hazard - Do not connect a wet instrument to the generator.
- Position instrument cords to prevent contact with the patient or other cords. Avoid wrapping cords around metal objects, as this can cause currents that may lead to shock, fire, or injury to the patient or surgical team.
- Carefully inspect all connections between the instrument and generator before use. Improper connections can result in arcing, sparks, accessory malfunction, or unintended surgical effects.
- Inspect the instrument and cords for breaks, cracks, nicks, or other damage before use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team, or cause damage to the instrument. Do not use a damaged instrument.

- Avoid using the instrument in the presence of flammable anesthetics or oxidizing gases, such as nitrous oxide (N20) and oxygen, or in close proximity to volatile solvents (e.g., ether or alcohol) to prevent the risk of explosion.
- Due to concerns regarding the potential carcinogenic and infectious nature of electrosurgical by-products (such as tissue smoke plume and aerosols), it is advisable to use protective eye wear, filtration masks, and effective smoke-evacuation equipment during both open and minimally invasive procedures.
- Connect adaptors and accessories to the electrosurgical unit only when the unit is turned off or in standby mode. Failure to adhere to this instruction may result in injury or electrical shock to the patient or operating personnel.
- Avoid placing fingers between the Lever and the Handle, or between the jaws. It may result in injury to the user.
- Avoid placing cable between the Lever and the Handle, or between the jaws.
- Use caution when handling the instrument between uses to prevent accidental activation. Do not place the instrument on the patient or on the drape when not in use.
- **Fire Hazard**: Avoid placing the instrument near or in contact with flammable materials such as gauze, surgical drapes, or flammable gases. Activating or hot instruments may cause a fire. When not in use, store the instruments in a clean, dry, highly visible area that is not in contact with the patient. Accidental contact with the patient may result in burns.
- For minimally invasive procedures, do not use hybrid trocars that have both metal and plastic components. Capacitive coupling of RF current may cause unintended burns.
- For minimally invasive procedures, verify the combability of the trocar size before use to easy insertion and extraction of the instrument. Carefully insert and withdraw the instrument through the cannula to avoid damage to the device and/or injury to the patient.
- Close the jaws before insertion/extraction through the access point.
- Do not use the instrument to ligate vessels larger than **7 mm** in diameter.
- Do not bend the shaft. If the instrument shaft is visibly bent, do not use it. Discard it, and replace it with a new instrument.
- Do not place the vessel and/or tissue near the jaw hinge. Always place the vessel and/or tissue at the center part of the jaws.
- Conductive fluids (e.g., blood or saline) in direct contact with or near the instrument may carry electrical current or heat, which can cause unintended burns to the patient. Aspirate fluid near the instrument jaws before activating the instrument.
- Keep the external surface of the instrument jaws away from adjacent tissue while activating the electrosurgical instrument system to avoid unintended injury.
- During a seal cycle, energy is applied to the tissue between the jaws. This generates heat and can cause water to convert into steam, which may result in unintended injury to tissues near the jaws. Exercise caution during surgical procedures in confined spaces.
- Release tension on the tissue when sealing and cutting to ensure proper function.
- Use caution when grasping, manipulating, sealing, and dividing large tissues.
- Do not attempt to seal or cut over clips or staples, as this may compromise seals. Contact between an active electrode and any metal objects may result in alternate site burns or incomplete seals.
- The surfaces of the jaws may remain hot enough to cause burns after the RF current is deactivated.

- Avoid using the instrument in the presence of flammable anesthetics or oxidizing gases, such as nitrous oxide (N20) and oxygen, or near volatile solvents (such as ether or alcohol) to prevent the risk of explosion.
- Due to concerns regarding the potential carcinogenic and infectious nature of electrosurgical by-products, such as tissue smoke plume and aerosols, it is recommended to use protective eyewear, filtration masks, and effective smoke-evacuation equipment in both open and minimally invasive procedures.
- Do not activate the instrument while the instrument jaws are in contact with or near other instruments, including metal cannulas, as localized burns to the patient or physician may occur.
- Exercise caution during surgical procedures involving patients with certain types of vascular pathology (such as atherosclerosis or aneurysmal vessels). For optimal outcomes, apply the seal to unaffected vasculature.
- Do not close the jaws if there's excessive tissue between the jaws.

Environmental Conditions for Transport and Storage

Temperature: -10°C ~ 55°C
Relative Humidity: ≤ 80 %
Air Pressure: 800 hPa ~ 1060 hPa

Expiration Date

The Instrument is sterilized by Ethylene Oxide. The expiration date is labeled on the package. The period of validity is 5 years from sterilization date. Do not use this Instrument beyond its expiration date.

How Supplied

This Instrument is supplied sterile for single patient use. Discard after use.

Troubleshooting

The following is a list of troubleshooting suggestions for situations encountered when using the instrument with compatible Reach Surgical generators. For details on specific situations, refer to generator instructions for reference.

Alerts

When an alert occurs, energy delivery stops. After the alert condition has been corrected, energy delivery will be immediately available. When an alert condition occurs, energy delivery stops, the generator produces a sequence of pulsed tones, and an alert will be displayed on the generator. Do Not Cut the Vessel The user should inspect the seal site and instrument before proceeding.

Troubleshooting steps


- Stop activation by releasing the Lever or Footswitch pedal.
- Open the jaws and inspect if the seal is complete.
- Follow the suggested corrective actions on the generator screen, the generator quick reference card, or in the generator user's guide.




- If possible, reposition the instrument and regrasp tissue in a location that overlaps the previous seal, then reactivate the seal.




Possible reasons and solutions

Inadequate tissue between the jaws: Open the jaws and grasp more tissues. Close the jaws and reactivate the seal.
Excessive tissue between the jaws: Open the jaws, reduce the amount of tissue grasped. Close the jaws and reactivate the seal.
Activating on a metal object: Avoid grasping objects, such as staples, clips, or encapsulated sutures in the jaws of the instrument.
Dirty jaws: Use a wet gauze pad to clean surfaces and edges of instrument jaws.
Excess Fluids in the Surgical Field: Reduce or remove excess fluids from around the instrument jaws.
Activation switch released before seal complete tone: The footswitch pedal or Lever was released before the seal cycle was complete.
A notice to the user and/or patient that any serious incident that has occurred in relation to the device should b e reported to Reach Surgical, Inc. through Reachquality@reachsurgical.com. and the competent authority of t he Member State in which the user and/or patient is established.

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
REACH SURGICAL | iREACH Magnum

IREACH MAGNUM PLUS Staplers


	Reach Surgical, Inc. 120 Xinxing Road,West Zone,TEDA,300462 Tianjin,P.R. China

Illustrations


01


02


03


04


05

06


07


08


09

10

11

Before you use this Stapler, please read the following contents carefully.

This document is designed to assist in using this product. It is not a reference for surgical techniques.

This Stapler is designed, inspected and manufactured for single procedure only. Reuse or reprocessing may result in product failure or cause patient injury. Do not reuse, reprocess or resterilize this Stapler.

Device Description

The Endoscopic Motorized Cutting Staplers (hereinafter referred to as the Stapler) are sterile, single patient use instruments that simultaneously cut and staple tissue. There are six staggered rows of staples, three on either side of the cut line. The Shaft can be rotated freely in both directions and an articulation mechanism enables the distal portion of the Shaft to pivot to facilitate lateral access to the operative site.

The Staplers are packaged with a primary Battery Pack that must be installed prior to use. There are specific requirements for disposing of the Battery Pack. Refer to the Battery Pack Disposal section.

The Staplers are packaged without Reloads and must be loaded prior to use. A Staple Retaining Cap on the Reload protects the loading unit staples points during shipping and transportation. The Stapler’s lock-out feature is designed to prevent a used or improperly installed Reload from being refired or a Stapler from being fired without a Reload.

Caution:

Do not load the Stapler more than 16 times.

The Stapler can be fired for a maximum of 16 times.

Use of staple line reinforcement material may reduce the maximum number of firings.

Intended Use

This instrument is intended for transection, resection of tissues and/or creation of anastomoses.

Indications

This instrument is intended to be used with the Reloads for transection, resection, and/or

creation of anastomoses. This instrument has applications in open and minimally invasive surgeries including thoracic, abdominal, gynecological and urological surgeries. It is used for transection and resection of lungs, bronchial tissue, intestines, stomach, urethra, kidney, uteru.

Intended User

This instrument is used for healthcare professionals who use this instrument for surgical purposes.

Intended Use Environment

This instrument is intended to be used in a hospital.

Intended patient population:

Patients aged 3 and older requiring resection and reconstruction of organs and tissues in th e thoracic and abdominal cavities.

Clinical benefit:

•Shorter operative time;

•Less intraoperative blood loss;

•Reduced postoperative complications such as anastomotic leak.

Contraindications

1. Do not use the Stapler on edematous tissue, tissue with a muscle layer that is too thick, or tissue with bad healing ability.

2. Do not use the Stapler on tissue with suspected cancer cell residual.

3. Do not use the Stapler on the aorta.

4. Do not use the Stapler on ischemic or necrotic tissue.

5. Do not use the Stapler on major vessels without making provision for proximal and distal control.

6. Do not use on tissue that cannot be properly compressed to the closed staple height or easily compressed to less than the closed staple height.

7.The Stapler is not intended for use when surgical stapling is contraindicated.

Compatibility Information

The Stapler is designed to be compatible with Battery Pack, Curved Tip Accessory and Reloads referred in this document. Third party accessories are not compatible.

When the Stapler is used for minimally invasive surgery, a trocar is needed.

Note: The Stapler is designed to be passed through a 12mm diameter trocar.

The Stapler is designed to be used with IREACH MAGNUM PLUS Reloads and ENDO REACH REC Reloads (hereinafter referred to as the Reloads) manufactured by Reach Surgical, Inc. Product Code and tissue thickness reference are listed as follows:

Product Codes and Specifications

Table 1 Product Codes and specifications

Product Code	**Staple Line length (mm)**	**Shaft Length (mm)**	**Trocar Compatibility** (mm)	**Battery included in the package**
IMS60L	60	365	12	YES
IMS60M	60	265	12	YES
IMS60S	60	215	12	YES
IMS45L	45	365	12	YES
IMS45M	45	265	12	YES
IMS45S	45	215	12	YES
IMS60L-0	60	365	12	NO
IMS60M-0	60	265	12	NO
IMS60S-0	60	215	12	NO

Product Code	**Staple Line length (mm)**	**Shaft Length (mm)**	**Trocar Compatibility** (mm)	**Battery included in the package**
IMS45L-0	45	365	12	NO
IMS45M-0	45	265	12	NO
IMS45S-0	45	215	12	NO

Table 2 Reload Product Codes

Product Code	Tissue Thickness (mm)	Staple Line Length (mm)	Color	Open Staple Height (mm)	Closed Staple Height (mm)
IMST45TAN	Vascular/Thin	45	Tan	2.0/2.5/3.0	0.75 /1.0/1.25
IMST45PUL	Medium/Thick	45	Purple	3.0/3.5/4.0	1.25/1.5/1.75
IMST45BLK	Extra Thick	45	Black	4.0/4.4/4.4	1.75/2.0/2.25
IMST60TAN	Vascular/Thin	60	Tan	2.0/2.5/3.0	0.75 /1.0/1.25
IMST60PUL	Medium/Thick	60	Purple	3.0/3.5/4.0	1.25/1.5/1.75
IMST60BLK	Extra Thick	60	Black	4.0/4.4/4.4	1.75/2.0/2.25
REC45GRA	Extra Thin	45	Gray	2.0	0.75
REC45WHT	Thin tissue	45	White	2.5	1.0
REC45BLU	Regular	45	Blue	3.5	1.5
REC45GLD	Regular/thick	45	Gold	3.8	1.75

REC45GRN	Thick tissue	45	Green	4.1	2.0
REC45BLK	Extra thick	45	Black	4.4	2.2
REC60GRA	Extra Thin	60	Gray	2.0	0.75
REC60WHT	Thin	60	White	2.5	1.0
REC60BLU	Regular	60	Blue	3.5	1.5
REC60GLD	Regular/thick	60	Gold	3.8	1.75
REC60GRN	Thick	60	Green	4.1	2.0
REC60BLK	Extra thick	60	Black	4.4	2.2

Table 3 Accessory

Trade Name	**Product Code**
IREACH Powered Battery	PB-A
IREACH Curved Tip Accessory	IMB1

Instructions for Use

Nomenclature (Illustration 01)

[A] Anvil Jaw	[B] Reload Jaw
[C] Staple Mark	[D] Cut Mark
[E] Reload Alignment Slot	[F] Knife Blade Indicator
[G] Articulation Joint	[H] Shaft
[I] Articulation Lever	[J] Rotation Knob
[K] Knife Reverse Switch	[L] Safety Release Button
[M] Anvil Release Button	[N] Manual Override

[O] Battery Pack	[P] Closing Trigger
[Q] Firing Trigger	[R] Battery Indicator
[S] Staple Retaining Cap	[T] Reload
[U] Alignment Tab	[V] Blade
[W] Proximal Mark	[X] Curved Tip Accessory

Prior to use

1. Verify the Stapler and Reload codes to be used as well as their compatibility, warnings and precautions.
2. Verify the integrity of the package for all the Stapler and accessories. Do not use if the package is damaged.
3. Verify the integrity of the antimicrobial package for all Stapler and accessories prior to using the Stapler. All Stapler and accessories are sold sterile.

Preparing the Stapler for Use

4. Using sterile technique, remove the Stapler, Battery Pack, and Reload from their respective packages. To avoid damage, do not flip the Stapler, Battery Pack, or Reload into the sterile field.
5. Install the Battery Pack. The Battery Pack must be installed prior to use. It can be inserted in only one orientation. (Illustration 05) Ensure that the Battery Pack is fully inserted into the Stapler. An audible click can be heard when the Battery Pack is fully inserted. The Battery Pack should not have any visible gap when the Battery Pack is fully inserted.
6. Ensure that the jaws are open and the Articulation Joint is straight, not articulated (Illustration 01).
7. Ensure that the Staple Retaining Cap on the Reload is in place; do not use the Reload if the Staple Retaining Cap is missing.

Caution: Tissue thickness should be carefully evaluated prior to using the Stapler. Refer to the codes and specifications of Single Use Loading Units for Endoscopic Linear Cutting Staplers for proper Reload selection.

Loading

8. Load the Reload by sliding it against the top of the jaw until the Reload Alignment Tab snaps in the Reload Alignment Slot. Remove the Staple Retaining Cap and discard. (Illustration 02) The Stapler is now loaded and ready for use.

Caution: Make sure that the knife does not touch the bottom of the jaw as this prevents the knife from being moved or it can become deformed when loading the Reload. Otherwise, the Stapler may be locked out and cannot be fired.

Caution: After removing the Staple Retaining Cap, observe the surface of the loaded Reload. The Reload must be replaced with another Reload if any colored drivers are visible as it may indicate absence of staples.

Using the Stapler

9. Close the jaws of the Stapler by squeezing the Closing Trigger until it locks in place (Illustration 06). An audible click indicates that the Closing Trigger and the jaws are locked. When the jaws of the stapler are closed, the Firing Trigger will be exposed.

Caution: Do not press the Safety Release Button at this time. Otherwise the Stapler may be partially or completely fired and will need to be reloaded before using on tissue.

10. Visually inspect the Stapler to ensure proper Reload seating. Introduce the Stapler into the body cavity through a trocar of the appropriate size or through an incision. When using a trocar, the jaws must be visible past the trocar sleeve before opening the jaws.

Caution: For insertion and removal, the jaws of the Stapler must be closed and straight, in-line with the Shaft of the Stapler. Failure to have the Stapler jaws in the straight position will result in difficult insertion or withdrawal of the Stapler and may result in damage to the Stapler or trocar.

Caution: When placing the Stapler through the trocar or incision, avoid inadvertently pressing the safety release button and pulling the Firing Trigger. Otherwise the Stapler may be partially or completely fired and will need to be reloaded before using on tissue. If the Stapler is partially fired, remove the Stapler and push the Knife Reverse Switch to return the knife to home position. Remove the Reload and insert a new one. If the Stapler is fired completely, the knife will return to home position automatically. Remove the Stapler and insert a new Reload (see Reloading the Stapler).

Caution: At any time, if the Knife Reverse Switch does not return the knife to its initial position and the jaws will not open:

A). First, ensure the Battery Pack is securely installed and the stapler has power; then, try the Knife Reverse Switch again.

B). If the knife still does not return, use the Manual Override. After the Manual Override system is used, the Stapler is disabled and cannot be used for any subsequent firings. To use the Manual Override, remove the access panel labeled “Manual Override” on the top of the Stapler handle. The reset lever and the forced return lever will be exposed. Pull the Reset Lever (marked ①) perpendicular to the stapler handle, move the Return Lever (marked ②) forward and backward until it can no longer be moved. The knife will now be in the initial position. This can be verified by viewing the position of the Knife Blade Indicator on the underside of the Reload Jaw. Discard the Stapler.

11. Once in the cavity, open the jaws by pulling the Anvil Release Button (Illustration 07).
12. If needed, rotate the jaws by pushing on the fins of the Rotation Knob with the index finger using a downward or upward pressure (Illustration 04). The Stapler Shaft will be rotated freely in either direction.
13. To articulate the jaws inside the body cavity, pivot the Articulation Lever clockwise or counter-clockwise to articulate the jaws, ensuring that it stays within the field of view. (the jaws must be open in order to articulate the Stapler). (Illustration 04)

Caution: The Stapler can only achieve a maximum articulation angle of 65°. When the maximum angle is reached, the force will increase, indicating that the maximum angle has been reached. Avoid applying excessive force to the Articulation Lever as damage to the Stapler may occur.

14. Position the Stapler around the tissue to be stapled.

Caution: Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the Reload, particularly in the base of the jaws, may result in an incomplete staple line. The Staple Mark on Reload Jaw designate the ends of the staple line. The line on the Reload Jaw that reads “cut” references the cut line on the Stapler.

Caution: When positioning the Stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the Stapler jaws. Firing over an obstruction may result in an incomplete cutting action, improperly formed staples, and/or inability to open the jaws.

15. After positioning the jaws, close the jaws by squeezing the Closing Trigger until it locks (Illustration 06). An audible click indicates that the jaws are fully closed. When the jaws are fully closed, the Safety Release Button is ready to be activated.

Caution: Ensure tissue has not been packed (extended) to the Proximal Mark on the Stapler. Tissue forced into the Stapler proximal to the Proximal Mark may be transected without staples. When firing across thick tissue, holding the jaws in place for 15 seconds after closing and prior to firing may result in better compression and staple formation.

Caution: If the trigger is difficult to lock:

A) Ensure that the proper Reload selection has been made (Refer to the Reload Selection Chart).

B) Reposition the Stapler and reduce the amount of tissue to be clamped.

C) When clamping across thick tissue, holding the jaws in place for 15 seconds after closing and prior to firing may result in better compression and staple formation.

Caution: If the clamping mechanism becomes inoperative and the jaws do not clamp on tissue, do not fire the Stapler. Remove and do not continue to use the Stapler.

16. Press the Safety Release Button on either side of the Stapler.

17. Fire the Stapler by pulling the Firing Trigger, so that the motor will be activated audibly

(Illustration 08). Continue to depress the trigger until the motor stops (audible feedback), in order to reach the end of the tissue transection.

Caution: Since the motor may stop if it stalls, it is important to do a visual check to ensure that the Knife Blade Indicator, on the underside of the Reload Jaw, has reached the end of the tissue transection.

Caution: The firing process may come to stop and the audible tone will be heard if there's excessive tissue between jaws. In such situation, use the Knife Reverse Switch to return the Blade to its original position and reload. It is advised to clamp less tissue and select proper Reload according to the thickness of the tissue.

Caution: The motor will slow down when device encounters thick tissues.

18. To complete the firing sequence, release the Firing Trigger and the motor will return the knife to home position where the motor will stop. In this position, the Stapler is locked out until the jaws are opened and re-closed.

Caution: If it is necessary to interrupt the firing sequence or it is interrupted inadvertently by releasing the trigger during the firing sequence, pull the Firing Trigger again to continue. The status of the transection can be determined by observing the Knife Blade Indicator on the underside of the Reload Jaw throughout the firing action. This can be done as often as necessary until the knife reaches its end where releasing the trigger automatically returns the knife to home position.

Caution: If the Stapler locks out, the motor will stop. Release the Firing Trigger and slide the Knife Reverse Switch forward to return the knife to the home position. While keeping the jaws closed, remove the Stapler from the body cavity, open the jaws and reload in order to continue. To open the jaws, press the Anvil Release Switch on either side of the Stapler while pulling on the Closing Trigger. Slowly release the Closing Trigger while pressure is still applied on the Anvil Release Switch. Follow the instructions for Reloading the Stapler.

Caution: At any time, if the Knife Reverse Switch does not return the knife to home position and the jaws will not open:

a. First, ensure the Battery Pack is securely installed and the Stapler has power. Then, try the Knife Reverse Switch again.

b. If the knife still does not return, use the Manual Override.

Caution: After the Manual Override system is used, the Stapler is disabled and cannot be used for any subsequent firings. To use the Manual Override, remove the access panel labeled “Manual Override” on the top of the Stapler handle. The Reset Lever and the forced Return Lever will be exposed. Pull the Reset Lever (marked ①) perpendicular to the Stapler handle, Move the Return Lever (marked ②) forward and backward until it can no longer be moved. The knife will now be in the home position. This can be verified by viewing the position of the Knife Blade Indicator on the underside of the Reload Jaw. Discard the Stapler.

Caution: Incomplete firing may result in malformed staples, incomplete cut line, bleeding, and/or difficult removal of the Stapler.

Caution: Crossing of staple lines may shorten the life of the Stapler.

Caution: If the firing mechanism becomes inoperative, do not continue to use the Stapler.

19. To open the jaws, squeeze the Closing Trigger, then simultaneously press the Anvil Release Switch on either side of the Stapler. (Illustration 07).

Caution: If the jaws do not automatically open after the Anvil Release Switch is pressed, first ensure that the knife is in the home position. The position of the knife can be determined by observing the Knife Blade Indicator under the Reload Jaw. If the Knife Blade Indicator is not in the home position or the position of the knife cannot be determined, slide the Knife Return Switch to activate the motor and return the knife to home position. Try to open the jaws again using the Anvil Release Switch. If the jaws do not open at this point, then gently push the Closing Trigger (1) upward (away from the handle) until both firing and Closing Triggers return to their original positions.

20**.** Gently pull the Stapler away from the transected tissue and ensure it is released from the jaws.

Caution: Examine the staple lines for hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.

21. Before removing the articulating Stapler, move the jaws away from any obstruction inside the body cavity while keeping the jaws open and within the field of view, and pivot the

Articulation Lever of the Stapler clockwise or counter-clockwise to return the jaws to the straight position manually.

Caution: For insertion and removal of articulating Stapler, the jaws must be straight, parallel to the Shaft of the Stapler. Failure to have the jaws in the straight position will result in difficult insertion or withdrawal of the articulating Stapler and may result in damage to the Stapler.

22. To remove the Stapler from the cavity, keep the jaws close and completely withdraw the Stapler from the body cavity or the trocar.

23. Push the Anvil Release Button to open the Stapler jaws (Illustration 07).

24. Push upward (toward the Anvil) to unsnap the spent Reload from the Reload Jaw. Discard the used Reload. (Illustration 03)

Caution: Before loading a new Reload, hold the Stapler in a vertical position, with Anvil and Reload Jaw completely submerged in sterile solution. Swish vigorously and then wipe the inside and outside surfaces of the Anvil and Reload Jaw to clean any unused staples from the Stapler. Do not use the stapler until it has been visually inspected to confirm there are no staples on the Anvil or Reload Jaw.

Caution: Before unloading the Reload, make sure the jaws and the Shaft are in a straight line. Otherwise, the Stapler may be damaged.

25. When needed, insert the Curved Tip Accessory into the Anvil Jaw (Illustration 10). After use, hold and remove it from the Anvil Jaw. (Illustration 11)

26.Reload and use the Stapler by repeating steps 6-24. The Stapler can be fired up to a total of 16 times.

Instrument and Reload Disposal

Once the instrument has been used, depending on your local regulations, it can be disposed of according to hospital procedure after removing the Battery Pack. If the instrument requires decontamination prior to disposal, follow the hospital protocol and local regulation regarding stapler disposal. The used instrument and Reload also share the same disposal method and are processed as per the Battery Pack’s cleaning and disinfection instructions below.

Battery Pack Disposal

The Battery Pack will automatically drain itself if it remained in the instrument.

The Battery Pack should be removed from the instrument prior to disposal. The Battery Pack and the instrument should be disposed separately or handled according to local regulations. If the Battery Pack requires cleaning and disinfection prior to disposal, follow the hospital protocol or the Battery Cleaning and Disinfection instructions below.

Prior to installation into the Instrument

If Battery Pack needs to be disposed of prior to installation into Stapler (e.g. product is beyond expiration date indicated on the package, Battery Pack is dropped), first install the Battery Pack into the instrument, then remove after the Battery Pack is depleted.

After Use

The Battery Pack must be removed from the instrument prior to disposal.

If Battery Pack requires decontamination prior to disposal, follow the hospital protocol or the Battery Pack Cleaning and Disinfection instructions below.

To remove the Battery Pack, squeeze the release tabs and pull the Battery Pack straight back.

Note: It is not necessary to disassemble the Battery Pack.

Caution: The disposable Battery Pack is not allowed to be charged or disassembled after use.

Battery Pack Cleaning and Disinfection

Warning: Do not use hospital autoclaves to sterilize or disinfect Battery Pack.

· Manual Cleaning

Step 1 Remove the Battery Pack from the Stapler before cleaning.

Note: Battery Pack should not be submerged in water or cleaning solutions.

Step 2 Clean the Battery Pack surfaces with a neutral pH detergent or neutral pH enzymatic detergent, prepared according to the manufacturer’s instructions.

Step 3 Use soft bristle brush to manually clean the Battery Pack with the cleaning solution.

Step 4 Ensure areas containing crevices are scrubbed thoroughly.

Step 5 Wipe off detergent thoroughly with lukewarm tap water.

Step 6 Perform visual inspection to determine if debris is removed.

Step 7 Repeat cleaning as necessary to obtain a visually clean Battery Pack.

· Chemical Disinfection

Disinfectants should be prepared and used according to the manufacturer's recommendations. It is recommended that the chemical disinfectant be wiped off with tap water.

- 70% Isopropyl alcohol
- 10% Bleach (sodium hypochlorite solution)

Standard Conventions Used

The Use of Caution, Warning, and Note Statements

Information relative to the operation of this Stapler will be supplied in the form of a Caution, Warning, or Note statement. These statements are found throughout the document. These statements should be read before continuing to the next step in a procedure.

Warning: A Warning statement indicates an operating, practice, or condition that, if not strictly observed, could result in personal injury or loss of property.

Caution: A Caution statement indicates an operating, practice, or condition that, if not strictly observed, could result in damage to or destruction of the Stapler.

Note: A Note statement indicates an operating, practice, or condition that is necessary to execute a task efficiently.

Warnings and Precautions

• Minimally invasive procedures should be performed only by persons having adequate training and familiarity with minimally invasive techniques. Consult medical literature relative to techniques, complications, and hazards prior to performance of any minimally invasive procedure.

• Minimally invasive Staplers may vary in diameter from manufacturer to manufacturer. When minimally invasive Staplers and accessories from different manufacturers are used together in a procedure, verify compatibility prior to initiation of the procedure.

• When using other technologies in the procedure, observe the precautions suggested by the

original equipment manufacturer to avoid the hazards associated with their use.

• Prior to using the Stapler, check the endoscope or endoscopic accessories inserted into the human body for rough surfaces, sharp edges or protrusions that may cause safety hazards.

• Failure to properly follow the instructions may lead to serious surgical consequences, such as leakage or disruption.

• Do not load the stapler more than 16 times for a maximum of 16 firings per Stapler.

• The Stapler can only be matched with the Battery Pack and Reload manufactured by Reach Surgical.

• Tissue thickness should be carefully evaluated prior to using the Stapler. Refer to the Table 2 in this manual for proper Reload selection.

• After removing the Staple Shipping Wedge, observe the Reload surface of the new Reload. The Reload must be replaced with another Reload if any colored drivers are visible. (If colored drivers are visible, the Reload may not contain staples.)

• For insertion and removal of Stapler, the jaw must be straight, parallel to the Shaft of the Stapler. Failure to have the jaws in the straight position will result in difficult insertion or withdrawal of the Stapler and may result in damage to the Stapler.

• When placing the Stapler through the trocar or incision, avoid inadvertently pressing the Safety Release Button. the Stapler may be partially or completely fired and will need to be reloaded before using on tissue.

• Do not attempt to articulate the jaws when clamping the tissues.

• The Stapler can only achieve a maximum articulation angle of 65°. When the maximum angle is reached, the force will increase indicating the maximum angle has been reached. Avoid applying excessive force to the Articulation Lever, otherwise it may cause damage to the Stapler.

• Ensure that the tissue lies flat and is positioned properly between the jaws. Any “bunching” of tissue along the Reload, particularly in the crotch of the jaws, may result in an incomplete staple line.

• When positioning the Stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the jaws. Firing over an obstruction may result in

incomplete cutting action, improperly formed staples, and/or inability to open the jaws.

• Ensure tissue has not extended proximal to the proximal black line on the Reload. Tissue forced into the Reload proximal to the black line may be transected without staples.

• If the trigger is difficult to be squeezed, reposition the Stapler and take a smaller amount of tissue. Ensure that the proper Reload selection has been made.

• If the trigger becomes inoperative and the jaws do not clamp on tissue, do not fire the Stapler. Remove and do not continue to use the Stapler.

• Attempting to force the trigger to complete the firing stroke with too much tissue between the jaws, or with dense/thick tissue between the jaws, may result in motor stall and the knife will stop.

• Since the motor may stop if it stalls, it is important to do a visual check to ensure that the I beam on the underside of the Reload Jaw, has reached the distal end of the tissue transection.

• If the Stapler locks out, the motor will stop. Push the Knife Reverse Switch forward to return the I beam to the home position. In this position, the Stapler should be removed, opened, and reloaded in order to continue.

• If the motor does not run when pushing the Knife Reverse Switch, the Manual Override panel on top of Handle should be opened to access Manual Override which will be used to backward Knife blade to its home position. Then the Stapler can not be used any more.

• Incomplete firing may result in malformed staples, incomplete cut line, bleeding, and/or difficult removal of the Stapler.

• If the firing mechanism becomes inoperative, do not continue to use the Stapler.

• If the jaws do not open after pushing the trigger forward, first ensure that the knife is in the home position. The position of the knife can be determined by observing the Knife Blade Indicator under the Reload Jaw. If the Knife Blade Indicator is not in the home position or the position of the knife cannot be determined, push the Knife Reverse Switch to return the knife to its home position.

• Examine the staple lines for pneumostasis/hemostasis and proper staple closure. Minor bleeding can be controlled with manual sutures or other appropriate techniques.

• Prior to loading a new Reload, hold the Stapler in a vertical position, with Anvil and Reload

Jaw completely submerged in sterile solution. Swish vigorously and then wipe the inside and outside surfaces of the Anvil and Reload Jaw to clean any unused staples from the Stapler. Do not use the Stapler until it has been visually inspected to confirm there are no staples on the Anvil or Reload Jaw.

• Before removing the Stapler, be sure that tissue is cleared from the jaws and then close the jaws.

• When selecting the Reload, careful consideration should be given to existing pathologic conditions as well as any pre-surgical treatment, such as radiotherapy, that the patient may have undergone. Certain conditions or preoperative treatments may cause change in tissue thickness that would exceed the indicated range of tissue thickness for the standard choice of Reload.

• Do not use hospital autoclaves to sterilize or disinfect Battery Pack and Stapler.

• Use of any other type of battery other than the battery supplied with the Stapler may result in increased EMISSIONS or decreased IMMUNITY of the Stapler.

• Avoid use of the Stapler adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, monitor the Stapler and the other equipment to assure normal operation.

• Do not modify this device and Battery Pack without authorization from the manufacturer.

• Device or Battery Pack which comes into contact with bodily fluids may require special disposal handling to prevent biological contamination.

• Use of accessories other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this Stapler and result in improper operation.

• This Stapler and Battery Pack are packaged and sterilized for single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the Stapler and/or lead to device failure that in turn may result in patient injury, illness or death. Also, reprocessing or resterilization of single use Stapler may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the Stapler may lead to injury, illness, or death of the patient.

• The Stapler and the Reload should not be used on tissues such as liver or spleen, as the

compression of such tissues may cause damage when the jaws are closed.

• If the hemostasis of the staple line cannot be clearly observed, this Stapler should not be used.

• If Stapler and auxiliary devices from different manufacturers are used in one operation, the compatibility of Stapler with devices from different manufacturers must be checked, and the insulation as well as grounding must be checked.

• User should not try to load the Reload while squeezing the trigger.

• The purchaser or user of the Stapler should use the Stapler under a specified electromagnetic environment. Follow the directions of the following tables when using the device. Otherwise, it may cause that the Stapler cannot work properly.

•The Stapler cannot be operated under oxygen enriched environment.

• A notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to Reach Surgical, Inc. through Reachquality@reachsurgical.com. and the competent authority of the Member State in which the user and/or patient is established;

Side effects

Potential complications related to the use of Stapler, Reload and Battery Pack include: hemorrhage, tissue injury, introduction of non-sterile surface or pathogen transfer, inflammatory or accidental tissue reaction, electrical shock, property damage or environmental damage. In addition, incomplete suture, inability to cut or device damage may cause accidental injury, prolongation of operation time or change of operation method.

Technical Parameters

Nominal Voltage: Battery Pack 12V.

The Endoscopic Motorized Cutting Staplers have a power rating of 40 W.

The Stapler is not resistant to fluids ingress and is classified per IEC 60601-1 as IPX0.

The Stapler needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document. Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

Essential Performance

1. The Stapler and the Reload should be able to be assembled and disassembled smoothly. The movement of moving and rotating parts should be flexible, without any stucking or loosening. The connection between the Stapler and the Reload should be reliable, without any loosening.

2. The safety mechanism of the Stapler should be opened and closed smoothly; it should be safe and reliable. Trigger return spring should have enough elasticity, the trigger should be able to quickly reset when releasing the trigger.

3. The firing and articulating operation should be safe and reliable. The electronic control system and driving motors should be safe and reliable, without any harmful noise, electricity, radiation to the user as well as nearby devices.

4. The Manual Override of the Stapler should be safe and reliable, without any stucking or loosening. Once the Battery Pack could not output enough power to Stapler, the Stapler should be able to finish retraction smoothly by Manual Override back to its original position, without any stuck.

Storage Requirements

Temperature: 5℃～35℃

Relative Humidity: 0%-70%

Air Pressure: 500 hPa～1060 hPa

Transport Requirements

Temperature: -10℃～54℃

Relative Humidity: 0%-70%

Air Pressure: 500 hPa～1060 hPa

Operating Environment Requirements

Temperature: 10℃～40℃

Relative Humidity: 35%-75%

Air Pressure: 800 hPa～1060 hPa

Software information

The Stapler is controlled by an embedded software program with the version of software V01. The software program is intended to detect the current. If the current is too high, the software will cut off the circuit to avoid motor damage.

Expiration Date

The Stapler, Curved Tip Accessory and Battery Pack are sterilized by ethylene oxide and the expiration date is labeled on the package of the Stapler. The Stapler should not be used after its expiration date.

Guidance and manufacturer's declaration for EMC

WARNING: Use of the Stapler adjacent to or with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used not closer than 30 cm (12 inches) to any part of the Endoscopic Motorized Cutting Staplers, including cables specified by the MANUFACTURER. Otherwise, the degradation of the performance of this equipment could result.

NOTE: The EMISSIONS characteristics of the Stapler make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), the stapler might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

• The Stapler must be used in a specified electromagnetic environment. Follow the directions of the following tables when using the device.

Guidance and manufacturer's declaration – electromagnetic emission		
The Endoscopic Motorized Cutting Staplers are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:		
Emission test	Compliance	Electromagnetic environment -

		guidance
RF emission CISPR 11	Group 1	The Endoscopic Motorized Cutting Staplers use RF energy only for its internal function. Therefore, its RF emissions are low and there is little possibility of producing interference to nearby electronic equipment.
RF emission CISPR 11	Class A	The Endoscopic Motorized Cutting Staplers are suitable for use in professional healthcare facilities.
Harmonic distortion IEC 61000-3-2	N/A	
Voltage fluctuations and flicker IEC 61000-3-3	N/A	


Guidance and manufacturer's declaration – Electromagnetic immunity			
The Endoscopic Motorized Cutting Staplers are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	The floor should be wood, concrete or ceramic. If the floor is covered with synthetic material, the relative humidity should be at least 30%.



Electrical fast transients / bursts IEC 61000-4-4	±2 kV 100 kHz repetition frequency	N/A	Battery powered and no signal line >3m
Surges IEC 61000-4-5	±1 kV line-to-line ±2 kV line-to-ground	N/A	Battery powered and no signal line >30m or going out to outdoor
Voltage dips IEC 61000-4-11 Voltage interruptions IEC 61000-4-11	0% U_T; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315° 70% U_T; 25 cycles at 0° 0% U_T; 250 cycle	N/A	Battery powered
Rated power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz	30 A/m 50Hz	The power frequency magnetic field should have the characteristics for use in a typical place in a typical commercial or hospital environment.
Note: U_T refers to the AC voltage of the power supply before the test voltage is applied.			


Guidance and manufacturer’s declaration –Electromagnetic immunity			
The Endoscopic Motorized Cutting Staplers are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted disturbances induced by RF fields IEC 61000-4-6 Radiated RF EM fields IEC 61000-4-3	3V 0.15 MHz～80 MHz 6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1kHz 3 V/m 80 MHz～2.7 GHz 80% AM at 1 kHz	N/A 3 V/m 80 MHz ～ 2.7 GHz 80% AM at 1 kHz	

Guidance and manufacturer’s declaration –Electromagnetic immunity						
The Endoscopic Motorized Cutting Staplers are intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in this electromagnetic environment:						
Immunity to RF wireless communications equipment (IEC 61000-4-3)						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380—390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430—470	GMRS 460, FRS 460	FM ±5kHz deviation 1kHz sine	2	0.3	28
710	704—787	LTE Band 13, 17	Pulse modulation 217Hz	0.2	0.3	9
745						
780						
810	800—960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
870						
930						
1720	1700—1990	GSM 1800; CDMA 1900;	Pulse modulation	2	0.3	28
1845						

1970		GSM 1900; DECT; LTE Band 1, 3, 4, 25, UMTS	217Hz			
2450	2400—2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217Hz	2	0.3	28
5240	5100—5800	WLAN 802.11 a/n	Pulse modulation 217Hz	0.2	0.3	9
5500						
5785						

	EN Sterilization batch
	EN Peel Here
	EN HDPE recyclable
	EN Recyclable
	EN Refer to instruction manual
	EN Do not use if package is damaged.
	EN Do not resterilize
	EN Manufacturer
	EN Date of manufacture
	EN Serial number

	EN Batch code
	EN Use-by date
	EN Fragile, handle with care
	EN Keep dry
	EN Keep away from sunlight
	EN Up
	EN Do not re-use
	EN Caution
	EN Catalogue number

	EN Storage temperature limit
	EN Storage humidity limitation
	EN Single sterile barrier system
	EN Country of manufacture
	EN Medical device
	EN Unique device identifier
	EN Sterilized by Ethylene Oxide.