Medtronic

Vectris™ SureScan™ MRI 1x8 Subcompact	977A160 977A175 977A190
Vectris™ SureScan™ MRI 1x8 Compact	977A260 977A275 977A290

Lead kit

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Souprava elektrody

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! USA Rx only

CE 0123

Explanation of symbols on product or package labeling

Refer to the appropriate product for symbols that apply.

CE 0123 Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union acts.

Importer

STERILE	EO

Sterilized using ethylene oxide

Medtronic and the Medtronic logo are trademarks of Medtronic. ™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

Table of contents

Description 7
Neurostimulation systems with SureScan MRI Technology 7
Intended purpose 7
Package contents 7
Device specifications 8
Implanted components and MRI scans 10
Implant criteria for full-body MRI scan eligibility 10
When changing components 11
Identifying implanted Vectris SureScan MRI leads 11
Instructions for use 13
Preparing for surgery 13
Placing a percutaneous lead 13
Testing stimulation intraoperatively 15
Anchoring the lead and creating strain-relief loops 16
Removing an anchor from the lead 21
Creating a subcutaneous pocket for the neurostimulator 23
Tunneling the lead for direct connection to a neurostimulator 24
Appendix: Trial procedures 26
Tunneling the percutaneous extension for the trial 26
Connecting the percutaneous extension to the lead and closing the incision for the trial 27
Preparing components for the trial 28
Removing the percutaneous extension after the trial 29

Information available for the system:

The information for prescribers manual provides information about contraindications, warnings, precautions, adverse events, sterilization, patient selection, individualization of treatment, and component disposal. For customers in Japan, the appropriate package insert provides information about safety, contraindications, warnings, precautions, and adverse events.

The indications sheet provides information about indications and related information. For customers in Japan, the appropriate package insert provides information about indications.

The system eligibility and battery longevity manual describes programming considerations and provides battery longevity information to aid in the appropriate neurostimulator selection.

MRI guidelines provide information about any MRI conditions and MRI-specific contraindications, warnings, and precautions for MRI scans with the neurostimulation system.

Product manuals, such as programming guides, recharging guides, and implant manuals provide device descriptions, package contents, device specifications, product-specific warnings and precautions, and instructions for use.

! USA **The clinical summary provides information about the clinical study results for the neurostimulation system.**

Description

The Medtronic Vectris™ SureScan™ MRI 1x8 Subcompact Lead Kit Models 977A160, 977A175, and 977A190 and the Vectris™ SureScan™ MRI 1x8 Compact Lead Kit Models 977A260, 977A275, and 977A290 are part of a neurostimulation system for pain therapy.

Neurostimulation systems with SureScan MRI Technology

When a Medtronic neurostimulation system with SureScan MRI Technology is implanted as directed (see "Implant criteria for full-body MRI scan eligibility" on page 10), a patient's full body may be eligible for MRI scans under specific conditions, ie, any part of the patient's anatomy can be scanned when specific conditions are met.

For the MRI conditions and MRI-specific warnings and precautions for conducting an MRI scan, refer to the *MRI guidelines for Medtronic neurostimulation systems for chronic pain instructions for use* manual.

Intended purpose

The lead is an implanted component of a neurostimulator system intended to conduct electrical stimulation from a neurostimulator to the intended target, as part of a neurostimulation system for pain therapy.

The accessories are intended to facilitate the secure placement of a lead.

Package contents

- Lead
- Injex™ Bumpy Anchor with anchor dispenser tool
- Guide wire
- Needles (epidural, modified Tuohy):
 - curved tip, 14 gauge, 8.9 cm
 - straight tip, 14 gauge, 8.9 cm
- Stylets:
 - bent, 0.30 mm diameter (green stylet handle and blue cap; inserted in the lead)
 - straight, 0.30 mm diameter (green stylet handle and green cap)
- Tunneling tools: tunneling rod, tunneling tip (narrow), tube (passing straw)
- Product literature
- ! USA Warranty card
- ! USA Registration form

Device specifications

The lead has electrodes on the distal end; the proximal (connector) end fits into an 8-conductor connector. A stylet has been inserted into the proximal end of the lead to aid in positioning the lead.


Figure 1. *Electrode, lead contact, and stylet (shown with bent stylet).*

Table 1. Device specifications for leads: Models 977A1 and 977A2 series[a,b]

Description	Models 977A160, 977A175, 977A190	Models 977A260, 977A275, 977A290
Connector	Octapolar, in-line	Octapolar, in-line
Conductor resistance[c,d]	1.6 Ω per cm	1.6 Ω per cm
Length	60, 75, 90 cm	60, 75, 90 cm
Diameter (lead body)	1.3 mm	1.3 mm
Distal (electrode) end		
Number of electrodes	8	8
Electrode shape	Cylindrical	Cylindrical
Electrode length	3.0 mm	3.0 mm
Electrode spacing	1.5 mm	4.0 mm
Electrode distance	34.5 mm	52.0 mm
Distal tip distance	Not applicable	Not applicable
Radiopaque tip marker	No	No
Proximal (connector) end		
Lead contact length	1.5 mm	1.5 mm
'0' contact length	3.0 mm	3.0 mm
Lead contact spacing	2.8 mm	2.8 mm
Lead contact distance	22.5 mm	22.5 mm
Stylet handle length	29.0 mm	29.0 mm
Stylet wire diameter	0.30 mm	0.30 mm
Expected lifetime	Five (5) years	Five (5) years

[a] All measurements are approximate.

[b] The numbers that follow the base model number in a series denote the lead length in centimeters.

[c] Electrical resistance of this device only.

[d] Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude.

Table 2. Material of components in the Models 977A1 and 977A2 series lead package

	Component	Materials[a,b,c]
Materials and substances to which the patient can be exposed	Lead	Polyurethane, platinum-iridium, epoxy
	Anchor	Silicone rubber, barium sulfate
	Anchor dispenser tool[d]	Polytetrafluoroethylene, stainless steel[e]
	Guide wire[d]	Stainless steel[e]
	Needles[d] (curved, straight)	Stainless steel,[e] nickel-plated brass,[e] epoxy
	Tunneling rod[d]	Stainless steel[e]
	Tunneling tip[d]	Stainless steel[e]
	Tube (passing straw)	Polytetrafluoroethylene

[a] The materials and substances to which the patient can be exposed are ordered from the most to least amount.

[b] Tested for category 1A and category 1B carcinogenic, mutagenic, or toxic for reproduction (CMR) substances, and endocrine disrupting chemicals (EDC).

[c] Does not contain natural rubber latex.

[d] This device may contain a substance: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0 defined as CMR 1B, in a concentration above 0.1% weight by weight.

[e] Contains nickel.

Implanted components and MRI scans

Implant criteria for full-body MRI scan eligibility

Caution: To allow a patient full-body MRI scan eligibility under specific conditions, implant a Medtronic neurostimulation system with SureScan MRI Technology as follows:

- Use only SureScan MRI neurostimulation system components (eg, leads and neurostimulators).

 Note: Extension model numbers in the 3708 series (eg, 37081, 37082, 37083) are not full-body MRI scan eligible.
- Implant the neurostimulator in the buttocks, abdomen, or flank (ie, the lateral and posterior region between the ribs and pelvis).
- Place the lead tip(s) in the spinal epidural space.

- Explant any previously abandoned pain leads or extensions that may be in the patient (ie, leads or extensions, or portions of, that are not connected to a neurostimulator).

 Note: Confirm MRI compatibility of any other implanted medical devices. Other implanted medical devices may limit or restrict MRI scans.
- Enter all component model number and implant location information using the clinician programmer.

If the above implant criteria are not met, the patient will not have a neurostimulation system with full-body MRI scan eligibility. MRI scan eligibility will be restricted.

For the MRI conditions and MRI-specific warnings and precautions for conducting an MRI scan, refer to the *MRI guidelines for Medtronic neurostimulation systems for chronic pain instructions for use* manual. MRI scans performed under different conditions can result in patient injury or damage to the implantable device.

When changing components

⚠ **Warning**: Before explanting and replacing an existing neurostimulator, print a report from the existing neurostimulator that shows implanted and any abandoned component information. When adding, changing, or removing neurostimulators, leads, extensions, and accessories, always program up-to-date component model numbers, implant locations, and any abandoned component information to re-establish MRI-scan eligibility.

If this information is not updated or is entered incorrectly, MRI scan-type eligibility data will be inaccurate, and the patient is at risk for one of the following:

- The patient is allowed to have an MRI scan inappropriate for the implanted components, which could cause tissue heating, resulting in tissue damage or serious patient injury.
- The patient is unnecessarily restricted from having an MRI scan.

Identifying implanted Vectris SureScan MRI leads

To identify an implanted Vectris SureScan MRI lead, use x-ray imaging and look for a transition of the lead body diameter size near the lead electrodes. The portion of the lead body with the braiding has a wider radiographic diameter than the lead body at the lead-electrode (distal) end (Figure 2).

Note: X-ray identification of the lead should only be used to confirm the accuracy of the component and implant location information entered into the neurostimulator. Do not use

x-ray identification to conclude that the entire neurostimulation system is full-body MRI scan eligible.


Figure 2. *A Vectris SureScan MRI lead (with braiding) compared with a lead without the braiding.*

Instructions for use

Implanting physicians should be experienced in epidural-access procedures and should be thoroughly familiar with all product labeling.

Cautions:

- Do not bend, kink, or stretch the lead or extension, which may damage the component.
- Do not use any instrument to handle the lead. The force may compress the lead, resulting in the inability to insert or withdraw the stylet.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

Cautions:

- Do not bend or kink the stylet, which may damage the component.
- When reinserting a stylet, use care because too much pressure on the stylet could damage the lead, resulting in intermittent or loss of stimulation.

Caution: Do not abandon pain leads or any portions of leads or extensions in the patient when changing components or permanently explanting the neurostimulator. Abandoned components may prevent the patient from being allowed MRI scans in the future due to concerns of lead electrode heating that can result in tissue damage.

Preparing for surgery

Before opening the lead package, verify the model number, use-by date, lead length, and connector type.

Note: Select a lead length that is long enough for strain-relief loops and direct connection to a SureScan MRI neurostimulator.

Warning: Medtronic has sterilized the contents of the package according to the process indicated on the package label. Visually inspect the sterile packaging for damage that may invalidate device sterility before opening, and inspect the components before use. Do not use a component if the following circumstances have occurred:

- The storage package or sterile seal has been pierced or altered and component sterility cannot be guaranteed and infection may occur.
- The component shows signs of damage and the component may not function properly.
- The use-by date has expired and component sterility cannot be guaranteed and infection may occur.

Placing a percutaneous lead

Warning: As with any spinal procedure, the risk of serious injury to the patient (eg, hemorrhage, hematoma, or paralysis) increases as location of the selected needle

insertion site progresses up the vertebral column—from a lower risk at a lumbar location to a higher risk at a cervical location. Select a vertebral location that provides the widest and easiest access to the epidural space during needle insertion to reduce the risk of serious patient injury resulting from direct trauma to the spinal cord.

⚠ **Warning**: Do not make custom bends to the distal (electrode) end of the lead with or without the stylet in the lead. Bending the lead at the distal end can break the conductor wires and render the respective electrode inoperable, which may require additional surgery to restore therapy. If the patient receives an MRI scan, broken conductor wires increase the risk of tissue heating, which could cause tissue damage or serious patient injury.

⚠ **Caution**: Do not place the lead electrodes in locations outside of the epidural spine. If other lead electrode locations are used (eg, neck or head), MRI scan eligibility will be restricted. The patient will not be full-body eligible for MRI scans.

⚠ **Caution:** To reduce the risk of lead damage that may result in intermittent stimulation or loss of stimulation, requiring additional surgery to replace the lead:

- use only the modified Tuohy needle supplied in the lead kit or compatible accessory kits (3550-14, 3550-28, 3550-43).
- use a shallow needle-insertion angle (45° or less) when inserting or withdrawing the needle into or out of the epidural space.

1. Make an incision at the needle-entry site to the depth of the subcutaneous fascia.
2. Using a paramedial approach and under fluoroscopy, insert the needle (included in the kit) into the epidural space at the appropriate angle until you encounter resistance from the ligamentum flavum.

 ⚠ **Caution**: Use a paramedial approach for needle placement so that lead placement is not on the midline. Over time, ligament or spinous process movement may cause damage to a lead implanted on the midline, resulting in loss of stimulation and requiring additional surgery to restore therapy.

3. Confirm needle location under fluoroscopy.
4. After rotating the needle so that the beveled edge faces cephalad, remove the needle stylet.
5. Advance the needle and confirm entry into the epidural space (eg, using the loss-of-resistance technique with air or sterile [United States Pharmacopeia—USP] water).

 ⚠ **Caution:** Do not use contrast media or a saline flush. Contrast media may obscure the field of view and a saline flush may increase the difficulty of lead placement.

6. For a second lead, repeat steps 1 to 5 noting these recommendations:
 - Implant the second lead parallel to the first lead and approximately 1 to 3 mm lateral of the physiological midline.

- Introduce the second lead one vertebral space below the first lead to help prevent nicking or cutting the first lead and to allow sufficient space for suturing both lead anchors.
- Stagger the lead tips or place them several vertebral spaces apart, depending on the position that produces the most effective paresthesia.

7. After inserting the guide wire through the needle, advance the guide wire no farther than 1 to 3 cm past the needle tip. Next, remove the guide wire from the needle.

 Note: If the guide wire track deviates from the intended pathway, steering and manipulating the lead will be more difficult.

8. Using fluoroscopy, slowly insert the lead through the needle and advance the lead to the initial target placement site. A stylet may need to be reinserted.

 Notes:

 - The bent stylet is keyed so that the flat notch of the stylet handle faces the same direction as the distal bent tip.
 - Ensure the lead electrodes are completely out of the distal end of the needle when changing a stylet. If the lead electrodes still partially remain in the needle, the force needed to insert the stylet may increase.
 - The force to insert and withdraw the stylet increases as the needle angle increases.
 - Before advancing the lead, ensure the stylet is fully inserted.
 - If resistance is felt while steering (with the stylet inserted), rotate the lead and stylet or pull back slightly to ensure the lead is not stuck. Use short, firm movements to advance the lead.

9. After verifying the lead electrode position under fluoroscopy (anterior-posterior and lateral views), compare that location with the location that is expected to have the highest probability of paresthesia coverage.

 Notes:

 - To increase lead stability, insert enough lead length to extend at least three vertebral bodies into the epidural space.
 - Position the lead so that the center bipole electrodes (eg, electrodes 3 and 4 on an octapolar lead) will be active. If lead migration occurs, effective stimulation may be regained through programming adjustments instead of surgical repositioning.

Testing stimulation intraoperatively

1. Refer to the instructions associated with the external devices (for example, external neurostimulator or trialing cable) used in intraoperative test stimulation.
2. Identify optimal stimulation parameters, beginning at a pulse width of 210 to 240 µs and a rate of 30 Hz.

 Note: Ensure that the patient can provide immediate feedback.

3. Increase the amplitude (V) or intensity (mA) beginning from the lowest setting while asking the patient closed-ended questions to identify the perception threshold (the setting at which the patient first perceives paresthesia), the discomfort threshold (the setting at which paresthesia is beyond the patient's tolerance), and the paresthesia coverage.

 Note: If good paresthesia coverage is not attained, change electrode settings before repositioning the lead to confirm the direction of lead movement.

4. If two leads were placed, repeat steps 2 – 3 for the second lead; optimize paresthesia coverage using both leads.
5. In the patient's chart, document the lead position that provided appropriate stimulation coverage (ie, record the settings and patient responses and include a fluoroscopic image of the final lead position).

Anchoring the lead and creating strain-relief loops

Note: The Injex anchor in this lead kit has been specifically designed for securing these leads.

△ **Caution**: Use only Injex anchors to secure the SureScan MRI lead in this kit. Other anchors are not recommended as they may damage the lead (eg, Titan or Twist-Lock anchors) or may not provide sufficient anchoring strength (eg, EZ, two-wing, three-wing anchors), resulting in lead migration. Lead damage or lead migration may require additional surgery to restore therapy.

1. Prepare the anchor site by making a 5 to 7 cm (2.0 to 2.8 in) longitudinal incision around the needle shaft, dissecting down to the supraspinous ligament, and establishing hemostasis.
2. Disconnect the stylet from the lead. Partially expose but do not completely withdraw the stylet.
3. While maintaining lead position, use minimal force to remove the needle and stylet.

 △ **Caution:** Use minimal traction to remove the needle because quick or sudden removal may dislodge the lead.

4. Thread 2-0 nonabsorbable suture material (such as silk or some types of braided polyester mesh) through the ligament or fascia in preparation for tying the anchor to secure the lead to the anchor site tissue in step 11.

 △ **Cautions:**

 - Do not use polypropylene suture material on silicone components. Polypropylene may damage the component, resulting in component failure.
 - Do not overtighten ligatures on the anchor or connector boot. Ligatures that are too tight may damage the component.

5. Locate the anchor dispenser tool with a preloaded anchor (Figure 3).


Figure 3. *Anchor dispenser tool with the anchor.*

6. Thread the connector end of the lead into the metal shaft of the anchor dispenser tool until the lead exits past the handle and through the slot (Figure 4).


Figure 4. *Threading the lead through the metal shaft.*

7. With one hand, grasp a portion of the lead that exited the slot, and with the other hand, gently advance the anchor dispenser tool until the anchor is at the anchoring site. Use care to maintain lead position.

Cautions:

- The anchoring site should be as close as possible to where the lead emerges from the fascia. Anchors not secured to the ligament or fascia in the anchor site may result in lead migration.
- Do not use medical adhesives with the anchors in this kit. Medical adhesives may not provide sufficient anchoring strength, resulting in lead migration.

Notes:

- If desired, embed the distal tip of the anchor into the anchor site tissue.
- Lubricants (eg, sterile water) are not required for anchor placement using the anchor dispenser tool.

8. Use fluoroscopy to verify lead location.
9. While maintaining lead position with one hand, position the fingers and thumb of the other hand on the anchor dispenser tool similar to holding a syringe (Figure 5).

 To release the anchor, use the fingers to pull the handle of the anchor dispenser tool towards the thumb. The anchor will slide off the metal shaft onto the lead.

 Note: If the lead or anchor is not dry—after releasing the anchor, maintain lead position for at least 15 seconds to allow for anchor compression and adequate retention on the lead.


Figure 5. *Releasing the anchor.*

10. Grasp the lead near the anchor and use your other hand to gently slide the anchor dispenser tool off from the end of the lead.

Note: If the anchor needs to be removed, use the anchor removal tool. Refer to "Removing an anchor from the lead" on page 21.

11. Tie two or more ligatures to secure the anchor to the ligament or fascia using the suture material prepared in step 4.

⚠ **Warning**: Do not tie ligatures directly around the lead body. Ligatures around the lead body can damage the lead body or conductor wires, resulting in a loss of therapy. During an MRI scan, the patient could be at risk of tissue heating, resulting in tissue damage or serious patient injury.

Cautions:

- Do not use polypropylene suture material on silicone components. Polypropylene may damage the component, resulting in component failure.
- Do not overtighten ligatures on the anchor or connector boot. Ligatures that are too tight may damage the component.

12. While maintaining lead position, create a strain-relief loop with the lead body near the secured anchor (Figure 6).


Figure 6. *Strain-relief loop made with the lead body.*

Cautions:

- Do not pull the lead or extension taut. Allow enough slack in the lead or extension to accommodate patient movement. Pulling the lead or extension taut may result in a short or open circuit or migration of implanted components.
- Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 7). A circular loop decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.


Figure 7. *Coiling excess lead.*

13. Ensure that the lead has not moved during the anchoring process. If the lead has moved, reestablish paresthesia coverage using program settings or minor lead repositioning as required.

 Notes:

 - For lead repositioning, refer to "Placing a percutaneous lead" on page 13.
 - If lead placement is for a trial and not for system implant, go to "Appendix: Trial procedures" on page 26 and continue from there.

Removing an anchor from the lead

An anchor removal tool is provided in the 97791 Injex Bumpy Anchor Accessory Kit.

If an anchor must be removed from the lead, proceed as follows:

1. Cut any suture material that is securing the anchor.
2. If necessary, expose the embedded anchor tip.
3. Locate the anchor removal tool (Figure 8).


Figure 8. *Anchor removal tool.*

4. Position the anchor removal tool near the anchor with the slitting blade facing in the direction of the lead electrodes and the bottom of the track aligned with the lead body.
5. Starting at the wide end of the track (opposite the slitting blade), slide a finger along the underneath side of the lead to press the lead into the track (Figure 9).


Figure 9. *Anchor removal tool placement on the lead.*

6. While holding the lead with one hand, use the other hand to push the anchor removal tool about one-half to two-thirds of the way through the anchor. Take care not to disrupt lead placement.
 - Keep the top of the anchor removal tool handle parallel with the anchor axis when slitting the anchor wall.

 Caution: Do not hold the anchor removal tool at a sharp angle and then use excessive force to push the tool through the anchor. Using a sharp angle and excessive force at the same time could result in dislodgement of the lead and/or damage to the lead or anchor removal tool. Dislodgement of or damage to the lead may require additional surgery to restore therapy.

 - Do not pull back on the lead to create counter tension. Pulling on the lead may dislodge it.
7. Push the anchor removal tool through the remainder of the anchor, while using the forceps to grasp a slit segment of the anchor, if needed (Figure 10).


Figure 10. *Using a grasping tool during anchor removal.*

8. Gently remove the anchor removal tool.
9. Remove any remaining suture material that was securing the anchor.
10. Gently remove the anchor from the lead.
11. Inspect the lead for any damage near the anchor site.
12. Discard the anchor and anchor removal tool after the anchor has been removed from the lead.
13. Use fluoroscopy to verify lead position.

Creating a subcutaneous pocket for the neurostimulator

Notes:

- Refer to the neurostimulator implant manual for the proper subcutaneous pocket depth.
- The neurostimulator should be placed on the opposite side of the body from another active implanted device and should be placed preferably on the right side of the body to allow for future placement of cardiac devices on the patient's left side.
- If a patient has a condition in the neurostimulator implant area that may require MRI scans, place the neurostimulator in another appropriate implant location to avoid image distortion from the neurostimulator at the anatomy of interest during an MRI scan.

1. Choose a neurostimulator implant location in the buttocks, abdomen, or flank (Figure 11) so that the patient will be full-body eligible for MRI scans.


Figure 11. *Neurostimulator implant locations for full-body MRI scan eligibility*

2. Make a small incision that will allow the neurostimulator to fit tightly in the subcutaneous pocket.
3. Use blunt dissection to create a subcutaneous pocket just large enough that the neurostimulator and excess lead or extension coils fit tightly.

 Caution: To prevent device inversion, do not make the neurostimulator pocket any larger than necessary to fit the neurostimulator and excess lead or extension. Device inversion may result in component damage, lead dislodgement, skin erosion, or stimulation at the implant site, requiring repeat surgery to restore therapy.

Tunneling the lead for direct connection to a neurostimulator

1. After simulating the tunneling route, mark the patient's skin at the neurostimulator pocket.

 Cautions:

 - When routing the lead, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
 - When multiple leads are implanted, route the leads so the area between them is minimized (Figure 12). If the leads are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).


Figure 12. *Routing for multiple leads.*

2. After assembling the tunneling tools, begin at the neurostimulator pocket and tunnel subcutaneously to the lead incision site.

 Caution: Proceed slowly when the tunneling tool approaches the lead site. If excess force is used, the patient could experience additional trauma when resistance to tunneling suddenly ceases.

 Notes:

 - Deep tunneling is not desirable.
 - Avoid the lower thoracic ribs.
 - If the tunneling tool does not extend to the lead incision site, make an intermediate incision.

3. Unscrew the tunneling tip; then, remove the tunneling tool, leaving the passing straw in place in the tunnel.
4. Gently insert the proximal (connector) end of the lead through the passing straw to the neurostimulator pocket.
5. After carefully removing the passing straw, verify that the location of the lead has not changed.
6. Refer to the neurostimulator implant manual for instructions on connecting the lead to the neurostimulator.

Appendix: Trial procedures

If these leads were implanted for a trial:

- A 1x8 extension kit is required.
- The extension connects to the external devices (for example, external neurostimulator or trialing cable) for the trial period to ensure that lead sterility is not compromised.
- The lead length chosen for the trial should be sufficient for direct connection to the neurostimulator at system implant.
- Record the lead models and tip locations for use at system implant to ensure that MRI scan eligibility is accurate.
- Perform the procedures in "Preparing for surgery" on page 13 through "Anchoring the lead and creating strain-relief loops" on page 16.

Note: Use the procedure "Removing an anchor from the lead" on page 21 only if needed.

⚠ **Warning**: Physicians should not prescribe MRI for patients undergoing a trial or test stimulation period, or who have any neurostimulation system components that are not fully implanted. Explant all trial components if an MRI scan is required. MRI has not been tested on trial components and may cause heating of the lead electrodes, resulting in tissue damage or serious patient injury.

⚠ **Warning**: Do not tie ligatures directly around the lead body. Ligatures around the lead body can damage the lead body or conductor wires, resulting in a loss of therapy. During an MRI scan, the patient could be at risk of tissue heating, resulting in tissue damage or serious patient injury.

Tunneling the percutaneous extension for the trial

1. Use blunt dissection to create a small, subcutaneous pocket for the lead-extension connection.

 Note: Ensure that the pocket is large enough to accommodate the lead-extension connection and excess, coiled lead length.

2. After simulating a tunneling route, mark the patient's skin at the lead-extension connection pocket and at the percutaneous extension exit site (exit site). The exit site should avoid the intended neurostimulator pocket site and be at least 10 cm (4.0 in) lateral to the lead-extension connection pocket.
3. After assembling the tunneling tools, make an incision at the exit site.
4. Begin at the exit site and tunnel subcutaneously to the pocket.

 Notes:

 - Deep tunneling is not desirable.
 - Avoid the lower thoracic ribs.

5. Unscrew the tunneling tip; then, remove the tunneling tool, leaving the passing straw in place in the tunnel.
6. Gently insert the proximal end of the extension (ie, the end that will be externalized) through the passing straw from the pocket to the exit site.

Connecting the percutaneous extension to the lead and closing the incision for the trial

Caution: Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site, intermittent stimulation, or loss of stimulation.

1. Slip the narrow end of the connector boot onto the lead body with the open end of the connector boot facing outward, exposing the lead contacts.

 Notes:

 - Irrigate the connector boot with a nonionic antibiotic solution. Dry the interior of the boot.
 - A clear and white radiopaque boot are provided in the extension kit for distinguishing the leads in two-lead implants.

2. Wipe the lead body and extension setscrew connector junction with sterile gauze. If necessary, use sterile (USP) water or a nonionic antibiotic solution. Dry all connections.
3. Insert the lead connector fully into the extension setscrew connector junction.

 Notes:

 - Each lead contact must be aligned under each extension connector contact.
 - During insertion, some resistance is typical.

4. Use the torque wrench from the extension kit to tighten the setscrew, which completes the electrical circuit with the lead contacts.

 Cautions:

 - To avoid overtightening, do not use a hex wrench to tighten an extension setscrew. Overtightening an extension setscrew may damage the lead contacts and cause an open or short circuit, resulting in intermittent or loss of stimulation.
 - Discard the torque wrench after making all connections. Reusing a torque wrench may result in undertightening or overtightening and subsequently, intermittent or loss of stimulation.

5. Carefully pull the proximal end of the extension (ie, the end that will be externalized) through the tunnel to the exit site.
6. Carefully remove the passing straw from the tunnel.
7. Connect the proximal (externalized) end of the percutaneous extension to the external device(s) used for checking impedances.

 Refer to the instructions associated with the external device(s) for connection information.

8. Confirm the integrity of the lead-extension connection by checking the impedance of all electrode pairs with the clinician programmer, then disconnect the external device(s).
9. After sliding the connector boot over the lead-extension connection, tie a ligature between the rings on the wide end of the connector boot.
10. Insert the coiled lead body in the lead-extension connection pocket, under the connection, leaving as much slack as possible in the lead body between the anchor and the lead-extension connection.

Cautions:

- Do not pull the lead or extension taut. Allow enough slack in the lead or extension to accommodate patient movement. Pulling the lead or extension taut may result in a short or open circuit or migration of implanted components.
- Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 13). A circular loop decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.


Figure 13. *Coiling excess lead.*

11. Close the lead incision, coiling the percutaneous extension under a dressing with the connector contacts exiting the bandage.

Preparing components for the trial

The patient's response to the system and the efficacy of the system should be evaluated in both a clinical and a home setting. Different electrode configurations should be evaluated at various parameter settings during the trial period.

1. Place a gauze bandage on the skin where the extension and external device(s) used for the trial will be placed on the patient.
2. Use the settings from "Testing stimulation intraoperatively" on page 15 to program the external neurostimulator for the trial.

3. Tape the extension and the external device(s) together to secure the connection.
4. Tape the entire assembly to the patient's skin, allowing for strain relief.
5. Proceed with the trial.

Removing the percutaneous extension after the trial

1. Refer to the instructions associated with the external device(s) to disconnect the device(s) from the percutaneous extension(s).
2. Using minimal traction, pull the percutaneous extension until approximately 1 to 2 cm (0.4 to 0.8 in) of the extension body is exposed from under the skin.

 Caution: When severing the percutaneous extension, use minimal traction on the extension because strong force may dislodge the lead.

3. Use sterile scissors to cut the exposed segment of the extension and allow the implanted part of the extension to retract under the skin.
4. Repeat steps 2 and 3 for the second extension, if applicable.
5. Discard the severed part of the percutaneous extension.
6. To remove the remaining segment of the percutaneous extension, find and expose the lead-extension connection, then make an incision large enough to grasp the lead body while maintaining lead position.

 Notes:

 - Blunt dissection may be needed to expose the lead-extension connection.
 - Use care to avoid damaging the coiled lead body in the pocket site.

7. While maintaining lead position, carefully remove the lead-extension connection from the incision.
8. Disconnect the lead from the percutaneous extension:
 a. Remove the ligature, then slip the connector boot off the connection.

 Caution: Do not use sharp instruments near the lead. Nicking or cutting the insulation may cause loss of stimulation or component failure requiring surgical replacement.

 b. Use the torque wrench from the INS kit to loosen the setscrew by turning the wrench counterclockwise.
 c. Gently separate the extension setscrew connector from the lead.

 Caution: If resistance is felt while removing the lead from the extension connector, first, loosen (but do not remove) the setscrews to ensure the lead contacts are not engaged. Next, inspect the lead connector for damage (eg, flattening of the lead contacts, stretching of the lead), which may result in intermittent or loss of stimulation.

 d. Remove and discard the connector boot; if necessary, carefully slit the boot.

9. Remove and discard the remaining segment of the percutaneous extension.
10. For a second percutaneous extension, repeat steps 6 to 9.
11. If proceeding with system implant, refer to "Creating a subcutaneous pocket for the neurostimulator" on page 23.

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
Use by

Lot number

Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union acts.

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Importer

Single sterile barrier system with protective packaging inside

Medtronic and the Medtronic logo are trademarks of Medtronic. ™* Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

Table of contents

Description 9

Additional contraindication 9

Intended purpose 9

Intended users 9

Package contents 9

Device specifications 10

Instructions for use 13

Preparing for surgery 13
Placing a percutaneous lead 13
Testing stimulation intraoperatively 15
Removing the needle and closing 16
Preparing components for the trial 16
Removing components after the trial 17

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Description

The Medtronic Vectris™ 1x8 Subcompact Trial Lead Kit Model 977D160 and Vectris™ 1x8 Compact Trial Lead Kit Model 977D260 are part of a neurostimulation system for pain therapy.

Additional contraindication

The Medtronic Vectris 1x8 Subcompact Model 977D160 and Vectris 1x8 Compact Model 977D260 Trial Leads are contraindicated for long-term implantation. These leads **MUST BE REMOVED** within thirty (30) days of implant.

Intended purpose

The trial lead is intended to conduct electrical stimulation from an external neurostimulator to the intended target as part of a neurostimulation system used for trial stimulation for pain therapy.

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Intended users

The primary users of trial lead products are implanting physicians who are experienced in spinal procedures. Implanting physicians include neurosurgeons, anesthesiologists, and orthopedic spine surgeons. Managing clinicians should be experienced in the diagnosis and treatment of chronic intractable pain and should be familiar with using the neurostimulation system. Managing clinicians include anesthesiologists, nurses, and pain management physicians. Medtronic field personnel provide therapy and product training to clinicians and patients.

The patient ultimately receives the therapy and feels the paresthesia and pain relief. The patient also interacts with their patient programmer and the recharger to manage their therapy. A percentage of patients have a caregiver that adjusts their therapy.

Package contents

- Lead
- Guide wire
- Needles (epidural, modified Tuohy):
 - curved tip, 14 gauge, 8.9 cm
 - straight tip, 14 gauge, 8.9 cm
- Stylets:
 - bent, 0.30 mm diameter (green stylet handle and blue cap; inserted in the lead)
 - straight, 0.30 mm diameter (green stylet handle and green cap)
- Product literature

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Table 1. Device specifications for leads: Models 977D160 and 977D260[a]

Description	Models 977D160	Models 977D260
Connector	Octapolar, in-line	Octapolar, in-line
Conductor resistance[b,c]	1.6 Ω per cm	1.6 Ω per cm
Length	60 cm	60 cm
Diameter (lead body)	1.3 mm	1.3 mm
Distal (electrode) end		
Number of electrodes	8	8
Electrode shape	Cylindrical	Cylindrical
Electrode length	3.0 mm	3.0 mm
Electrode spacing	1.5 mm	4.0 mm
Electrode distance	34.5 mm	52.0 mm
Distal tip distance	Not applicable	Not applicable
Radiopaque tip marker	No	No
Proximal (connector) end		
Lead contact length	1.5 mm	1.5 mm
'0' contact length	3.0 mm	3.0 mm
Lead contact spacing	2.8 mm	2.8 mm
Lead contact distance	22.5 mm	22.5 mm
Stylet handle length	29.0 mm	29.0 mm
Stylet wire diameter	0.30 mm	0.30 mm
Expected lifetime	The trial leads have a minimum expected lifetime that exceeds the 30-day trial period.	The trial leads have a minimum expected lifetime that exceeds the 30-day trial period.

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[c] Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude.

Table 2. Material of components in the Models 977D160 and 977D260 lead package

	Component	Materials[a,b,c]
Materials and substances to which the patient can be exposed	Lead	Polyurethane, platinum-iridium, epoxy
	Guide wire[d]	Stainless steel[e]
	Needles[d] (curved, straight)	Stainless steel[e], nickel-plated brass[e], epoxy

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[b] Does not contain natural rubber latex.

[c] Tested for category 1A and category 1B carcinogenic, mutagenic, or toxic for reproduction (CMR) substances, and endocrine disrupting chemicals (EDC).

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[e] Contains nickel. Nickel is used in mixtures of metals, such as stainless steel, and is a known allergen. Discuss allergies or other intolerances related to the materials with the patient before the procedure.

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Placing a percutaneous lead

Warning: As with any spinal procedure, the risk of serious injury to the patient (eg, hemorrhage, hematoma, or paralysis) increases as location of the selected needle insertion site progresses up the vertebral column—from a lower risk at a lumbar location to a higher risk at a cervical location. Select a vertebral location that provides the widest and easiest access to the epidural space during needle insertion to reduce the risk of serious patient injury resulting from direct trauma to the spinal cord.

Caution: To reduce the risk of lead damage that may result in intermittent stimulation or loss of stimulation, requiring additional surgery to replace the lead:

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 Caution: Use a paramedial approach for needle placement so that lead placement is not on the midline. Over time, ligament or spinous process movement may cause damage to a lead implanted on the midline, resulting in loss of stimulation and requiring additional surgery to restore therapy.

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 Caution: Do not use contrast media or a saline flush. Contrast media may obscure the field of view and a saline flush may increase the difficulty of lead placement.

6. For a second lead, repeat steps 1 to 5 noting these recommendations:
 - Implant the second lead parallel to the first lead and approximately 1 to 3 mm lateral of the physiological midline.
 - Introduce the second lead one vertebral space below the first lead to help prevent nicking or cutting the first lead.
 - Stagger the lead tips or place them several vertebral spaces apart, depending on the position that produces the most effective paresthesia.
7. After inserting the guide wire through the needle, advance the guide wire no farther than 1 to 3 cm past the needle tip. Next, remove the guide wire from the needle.

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 Notes:

 - The bent stylet is keyed so that the flat notch of the stylet handle faces the same direction as the distal bent tip.

- Ensure the lead electrodes are completely out of the distal end of the needle when changing a stylet. If the lead electrodes still partially remain in the needle, the force needed to insert the stylet may increase.
- The force to insert and withdraw the stylet increases as the needle angle increases.
- Before advancing the lead, ensure the stylet is fully inserted.
- If resistance is felt while steering (with the stylet inserted), rotate the lead and stylet or pull back slightly to ensure the lead is not stuck. Use short, firm movements to advance the lead.

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Removing the needle and closing

After the lead is in place, the physician can remove the needle, close, and use medical tape to secure the lead to the skin.

1. While holding the lead in place, disconnect the stylet from the lead. Partially expose, but do not completely withdraw the stylet.
2. While maintaining lead position, use minimal force to remove the needle and stylet.

 Caution: Use minimal traction to remove the needle because quick or sudden removal may dislodge the lead.

3. If two leads were placed, repeat steps 1 – 2 for the second lead.
4. Connect the external device(s) used for the trial period.

 Refer to the instructions associated with the external device(s) for connection information.
5. Verify stimulation parameters to ensure that the lead has not moved. If the lead has moved, reposition it.
6. Close the lead exit site using standard suture technique.
7. Apply a dressing over the exit site of the lead.

 Note: An antimicrobial dressing may be used at the lead exit site.
8. Place a sterile gauze pad between the skin and the lead. Place the lead over the sterile pad lateral to midline.
9. Use medical tape to secure the lead to the skin.

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1. Place a gauze bandage on the skin where the lead and external device(s) used for the trial will be placed on the patient.
2. Use the settings from "Testing stimulation intraoperatively" on page 15 to program the external neurostimulator for the trial period.
3. Tape the lead and the external device(s) together to secure the connection.
4. Tape the entire assembly to the patient's skin, allowing for strain relief.
5. Proceed with the trial.

Warning: Physicians should not prescribe MRI for patients undergoing trial stimulation or who have any neurostimulation system components that are not fully implanted. Explant all trial stimulation components if an MRI scan is required. MRI has not been tested on trial stimulation components and may cause heating of the lead electrodes, resulting in tissue damage or serious patient injury.

Removing components after the trial

1. Verify that the external neurostimulator is turned off.
2. Refer to the instructions associated with the external device(s) used for the trial for information on removing.
3. Remove all tape from the length of the lead body.
4. Open the original lead incision.
5. Use your thumb and index finger to grasp the lead as close as possible to the lead exit site.
6. Carefully retract the lead from the epidural space.
7. If the trial was successful, proceed with the chosen therapy.
8. If the trial was unsuccessful, clean the exit site(s) and close.

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CE 0123

M954166A046 Rev B

Printing instructions: *Refer to doc# A88898 for Printing Instructions*

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Temperature limitation

Lot number

Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Union acts.

Authorized Representative in the European Community

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Table of contents

Description 9
Neurostimulation systems with SureScan MRI Technology 9

Intended purpose 9

Package contents 9

Device specifications 9
Specify SureScan MRI 5-6-5 Surgical lead 10
Specify SureScan MRI 2x8 Surgical lead 13

Implanted components and MRI scans 15
Implant criteria for full-body MRI scan eligibility 15
When changing components 16
Identifying implanted Specify SureScan MRI leads 16

Instructions for use 18
Preparing for surgery 19
Placing a surgical lead 19
Testing stimulation intraoperatively 20
Implant procedure: connecting the lead directly to a neurostimulator 20
Creating strain-relief loops 21
Anchoring the lead (optional) 23
Creating a subcutaneous pocket for the neurostimulator 23
Tunneling the lead for direct connection to a neurostimulator 24

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Description

The Medtronic Specify™ SureScan™ MRI Lead Kits, Models 977C165, 977C190, 977C265, and 977C290, are part of a neurostimulation system for pain therapy.

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When a Medtronic neurostimulation system with SureScan MRI Technology is implanted as directed (see "Implant criteria for full-body MRI scan eligibility" on page 15), a patient's full body may be eligible for MRI scans under specific conditions, ie, any part of the patient's anatomy can be scanned when specific conditions are met.

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Package contents

- Lead
- Anchors:
 - two-wing (2)
 - three-wing (2)
 - EZ (2)
- Passing elevator
- Lead blank
- Tunneling tools: tunneling rod, tunneling tip (narrow), tube (passing straw)
- Torque wrench
- Product literature
- ! USA Warranty card
- ! USA Registration form

Device specifications

The lead has 16 electrodes on the distal end; the two proximal ends fit into 8-conductor connectors (ie, a Medtronic neurostimulator). One lead body has a white marker band to indicate electrodes 0-7.

Specify SureScan MRI 5-6-5 Surgical lead


Figure 1. *Specify SureScan MRI 5-6-5 Models 977C165 and 977C190 Surgical Lead (electrode side shown).*

Table 1. Device specifications[a] for leads: Models 977C165 and 977C190

Description	Models 977C165, 977C190
Connector	Octapolar, in-line
Shape	Contoured
Conductor resistance[b,c]	<77 Ω (65 cm)
	<107 Ω (90 cm)
Length (includes lead paddle)	65, 90 cm
Diameter	1.3 mm
Distal (electrode) end	
Number of electrodes	16
Electrode shape	Rectangular
Electrode size (width x length)	1.5 mm x 4.0 mm
Electrode stimulating area	6.0 mm²
Electrode spacing (edge to edge)	
in-line spacing	4.5 mm
row spacing	1.0 mm
Lead paddle length	64.2 mm
Lead paddle width	10 mm
Lead paddle thickness	2.0 mm
Proximal (connector) end	
Lead contact length	1.5 mm
Lead contact spacing (center to center)	2.8 mm
Lead contact distance	22.5 mm
Expected lifetime	Five (5) years

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[b] Electrical resistance of this device only.

[c] Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude or intensity.

Table 2. Material of components in the Models 977C165 and 977C190 lead package

	Component	Materials[a,b,c]
Materials and substances to which the patient can be exposed	Lead	Polyurethane, silicone rubber, platinum-iridium alloy, silicone rubber adhesive
	Anchors (Two-wing, three-wing, EZ)	Silicone rubber
	Passing elevator	Acetal resin, silicone rubber, barium sulfate
	Lead blank	Silicone rubber, barium sulfate
	Tunneling rod[d]	Stainless steel[e]
	Tunneling tip[d]	Stainless steel[e]
	Tube (passing straw)	Polytetrafluoroethylene

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[e] Contains nickel.

Specify SureScan MRI 2x8 Surgical lead


Figure 2. *Specify SureScan MRI 2x8 Models 977C265 and 977C290 Surgical Lead (electrode side shown).*

Table 3. Device specifications[a] for leads: Models 977C265 and 977C290

Description	Models 977C265, 977C290
Connector	Octapolar, in-line
Shape	Contoured
Conductor resistance[b,c]	<77 Ω (65 cm)
	<107 Ω (90 cm)
Length (includes lead paddle)	65, 90 cm
Diameter	1.3 mm
Distal (electrode) end	
Number of electrodes	16
Electrode shape	Rectangular
Electrode size (width x length)	1.5 mm x 4.0 mm
Electrode stimulating area	6.0 mm²
Electrode spacing (edge to edge)	
in-line spacing	1.0 mm
row spacing	1.0 mm
Lead paddle length	56.4 mm
Lead paddle width	7.6 mm
Lead paddle thickness	1.9 mm
Proximal (connector) end	
Lead contact length	1.5 mm
Lead contact spacing (center to center)	2.8 mm
Lead contact distance	22.5 mm
Expected lifetime	Five (5) years

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[b] Electrical resistance of this device only.

[c] Resistance is proportional to length: long lengths have higher resistance that may limit the amplitude or intensity.

Table 4. Material of components in the Models 977C265 and 977C290 lead package

	Component	Materials[a,b,c]
Materials and substances to which the patient can be exposed	Lead	Polyurethane, silicone rubber, platinum-iridium alloy, silicone rubber adhesive
	Anchors (Two-wing, three-wing, EZ)	Silicone rubber
	Passing elevator	Acetal resin, silicone rubber, barium sulfate
	Lead blank	Silicone rubber, barium sulfate
	Tunneling rod[d]	Stainless steel[e]
	Tunneling tip[d]	Stainless steel[e]
	Tube (passing straw)	Polytetrafluoroethylene

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⚠ **Caution**: To allow a patient full-body MRI scan eligibility under specific conditions, implant a Medtronic neurostimulation system with SureScan MRI Technology as follows:

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Identifying implanted Specify SureScan MRI leads

To identify an implanted Specify SureScan MRI lead, use x-ray imaging and look for a transition of the lead body diameter size where the braiding ends inside the paddle (Figure 3).

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
Figure 3. *A Specify SureScan MRI lead (with braiding) compared with a lead without the braiding (5-6-5 lead shown).*

Instructions for use

Implanting physicians should be experienced in open spinal surgical procedures and should be thoroughly familiar with all product labeling.

Note: Implant this lead after the patient completes a test stimulation period to determine the effectiveness of stimulation in managing pain. The Specify SureScan MRI 5-6-5 and 2x8 surgical leads are not intended for test stimulation outside of the operating room.

⚠ **Warning:** Evaluate the geometry and suitability of the spinal canal (by using an appropriate imaging method) to ensure that the spinal canal can safely accommodate insertion and permanent placement of the surgical lead paddle and body. The following factors should be considered during this evaluation:

- Adequate space within the spinal canal
- Pre-existing conditions or surgical procedures that may have altered the patient's epidural anatomy
- Potential for postimplant tissue changes around the lead
- Changes in epidural space associated with typical spinal movement

Inadequate space to accommodate the lead may result in localized compression of the spinal cord, increasing the risk of spinal cord injury.

Note: Surgical planning considerations should include insertion of the lead at a shallow angle. A steep angle may cause compression of the spinal cord and possible neurological deficits.

⚠ **Warning:** Safety and effectiveness of this lead has not been established for use with more than one neurostimulator. If this lead is connected to more than one neurostimulator, damage to nerve tissue may result.

△ **Caution:** Implanting the surgical lead paddle in anatomical locations that have a high degree of mobility (such as the cervical spine) is not recommended. The stresses associated with these locations may reduce the functional survival time of the lead, requiring reprogramming and/or early replacement of the lead to restore effective therapy.

△ **Cautions:**

- Do not bend, kink, or stretch the lead, which may damage the component.
- Use only smooth forceps on the lead. Do not use sharp-edged instruments (eg, hemostat), which may nick or cut the insulation.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

△ **Caution**: Do not abandon pain leads or any portions of leads or extensions in the patient when changing components or permanently explanting the neurostimulator. Abandoned components may prevent the patient from being allowed MRI scans in the future due to concerns of lead electrode heating that can result in tissue damage.

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Placing a surgical lead

Caution: Do not place the lead electrodes in locations outside of the spinal epidural space. If other lead electrode locations are used (eg, non-epidural neck or head), MRI scan eligibility will be restricted. The patient will not be full-body eligible for MRI scans.

Consider the following suggestions when performing the laminotomy:

Note: A laminotomy is preferred over a laminectomy because surgical leads are designed to fit in the epidural space under the vertebral laminae.

- Perform a laminotomy 1-2 vertebral spaces inferior to the desired final location of the surgical lead paddle using standard surgical technique.
- Use standard surgical technique to minimize bleeding and maintain hemostasis to reduce the risk of epidural hematoma.
- Make the laminotomy width consistent with the lead paddle's width.

 Caution: Ensure that the laminotomy is wide enough to accommodate the surgical lead paddle. Failure to do so may result in forcing the lead into the epidural space, which may damage the lead.

Note: The use of the passing elevator and the use of the lead blank are optional.

Warning: Introduce the passing elevator, the lead blank, and the surgical lead at a shallow angle. Using an angle that is too steep could cause contusion to the spinal cord. If resistance is encountered, do not force the passing elevator, lead blank, or surgical lead. Forced insertion into the epidural space could cause spinal cord injury, which could result in neurological deficits.

1. The passing elevator and the lead blank may be used prior to lead insertion.

 Optional use of the passing elevator:

 a. Slowly and with minimal force, introduce the passing elevator into the epidural space along the midline, guiding the entry so that the approach angle is as shallow as possible.
 b. Remove the passing elevator.

 Optional use of the lead blank:

 a. Slowly and with minimal force, introduce the radiopaque lead blank into the epidural space to ensure proper size and location of the lead site.
 b. Remove and discard the lead blank.

 ⚠ **Caution:** Make bony edges smooth, as needed, to prevent damage to the lead.

2. Position the surgical lead:
 a. Using smooth forceps to handle the lead paddle, carefully position the lead paddle in the epidural space.

 Note: Make sure the stimulating electrodes face the dura mater.
 b. Advance the lead paddle cephalad until the entire paddle is in the epidural space.
 c. For bilateral pain, place the lead close to the midline.
3. Verify the lead position under fluoroscopy (AP and lateral views).

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1. Refer to the instructions associated with the external device(s) used in intraoperative test stimulation.
2. Program stimulation parameters, beginning at a pulse width of 210 to 240 µs and a rate of 30 Hz.

 Note: Ensure that the patient can provide immediate feedback.
3. Increase the amplitude (V) or intensity (mA) beginning from the lowest setting while asking the patient closed-ended questions to identify the perception threshold (the setting at which the patient first perceives paresthesia), the discomfort threshold (the setting at which paresthesia is beyond the patient's tolerance), and the paresthesia coverage.

 Note: If good paresthesia coverage is not attained, change electrode settings before repositioning the lead to confirm the direction of lead movement.
4. In the patient's chart, document the lead position that provided appropriate stimulation coverage (ie, record the settings and patient responses and include a fluoroscopic image of the final lead position).

Implant procedure: connecting the lead directly to a neurostimulator

Following intraoperative test stimulation, use this procedure to connect the surgical lead directly to a Medtronic neurostimulator.

This procedure includes the following sections:

- Creating strain-relief loops
- Anchoring the lead (optional)
- Creating a subcutaneous pocket for the neurostimulator
- Tunneling the lead for direct connection to a neurostimulator

Creating strain-relief loops

Caution: Use strain-relief loops as directed when implanting surgical leads. Implanting surgical leads without strain-relief loops may reduce the functional survival time of the lead, requiring reprogramming and/or early replacement of the lead to restore effective therapy.

1. While maintaining lead position, create a strain-relief loop with each of the lead bodies, by coiling the lead bodies near the lead paddle (Figure 4).


Figure 4. *Strain-relief loops made with the lead bodies (shows optional use of anchors).*

Cautions:

- Do not pull the lead taut. Allow enough slack in the lead to accommodate patient movement. Pulling the lead taut may result in a short or open circuit or migration of implanted components.
- Coil the excess lead body into a circular loop greater than 2 cm (0.8 in) in diameter. Do not use a U-shaped loop or bend (Figure 5). A circular loop

decreases the possibility of electromagnetic interference and its effects and prevents kinking or damaging the lead body.


Figure 5. *Coiling excess lead.*

2. Insert the strain-relief loops into the lead pocket (Figure 6).


Figure 6. *Insert strain-relief loops into lead pocket (shows optional use of anchors).*

3. If desired, use the external neurostimulator to verify that the location of the lead and the pattern of stimulation have not changed.

Anchoring the lead (optional)

The use of anchors is at the clinician's discretion. Proper use of anchors requires the use of strain-relief loops in order to reduce the stresses placed on the lead. The use of non-Medtronic anchors is not recommended and may limit or restrict MRI scans.

△ **Caution:** Do not tie ligatures around the lead body or place sutures through the lead paddle, which may damage the insulation.

Notes:

- Each lead body should be anchored separately.
- The anchor can be lubricated with sterile water.
- The EZ anchor can be trimmed to a shorter length.
- Use 2-0 nonabsorbable suture (such as silk or some types of braided polyester mesh) to secure the anchor.

1. While maintaining lead position, slide the anchor onto the proximal end of each lead body.
2. Leave slack in the lead body and create a strain relief loop between the lead paddle and the anchor.
3. Secure an anchor to each lead body by tying ligatures around the grooves in the middle of the anchor.

 △ **Cautions:**

 - Do not use polypropylene suture material on silicone components. Polypropylene may damage the component, resulting in component failure.
 - Do not overtighten ligatures on the anchor. Ligatures that are too tight may damage the component.

4. Secure the anchors in place by suturing the anchors to the supraspinous ligament or deep fascia.
 - **For the EZ anchor:** Make a ligature, pass the suture around the anchor and make another ligature. Use a minimum of 4 sutures. Suture the anchor to the supraspinous ligament or deep fascia.
 - **For the winged anchor:** Suture the anchor to the supraspinous ligament or deep fascia after securing the anchor to the lead.

 Note: Laboratory testing has shown that injecting silicone medical adhesive between the anchor and lead after tying the ligatures increases the anchor-to-lead holding force.
5. If desired, use the external neurostimulator to verify that the location of the lead and the pattern of stimulation have not changed.

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1. Choose a neurostimulator implant location in the buttocks, abdomen, or flank (Figure 7) so that the patient will be full-body eligible for MRI scans.


Figure 7. *Neurostimulator implant locations for full-body MRI scan eligibility*

2. Make a small incision that will allow the neurostimulator to fit tightly in the subcutaneous pocket.
3. Use blunt dissection to create a subcutaneous pocket just large enough that the neurostimulator and excess lead coils fit tightly.

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Note: If a longer tunneling tool is desired, a Medtronic tunneling tool accessory kit may be used.

1. After simulating the tunneling route, mark the patient’s skin for the neurostimulator pocket.

 Cautions:

- When routing the lead, avoid sharp bends or kinks, which may break the wires. Broken wires may create an open circuit, resulting in loss of stimulation or component failure and requiring surgical replacement.
- When routing two lead bodies, assure that the area between them is minimized (Figure 8). If the lead bodies are routed in a loop and the patient is exposed to some sources of electromagnetic interference (eg, theft detectors), the patient may perceive a momentary increase in stimulation, which some patients have described as uncomfortable stimulation (jolting or shocking sensation).


Figure 8. *Routing for two lead bodies.*

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Caution: Proceed slowly when the tunneling tool approaches the lead site. If excess force is used, the patient could experience additional trauma when resistance to tunneling suddenly ceases.

Notes:

- Deep tunneling is not desirable.
- Avoid the lower thoracic ribs.
- If the tunneling tool does not extend to the lead incision site, make an intermediate incision.

3. Unscrew the tunneling tip; then, remove the tunneling tool, leaving the passing straw in place in the tunnel.
4. Gently insert the proximal (connector) end of both lead bodies through the passing straw to the neurostimulator pocket.
5. After carefully removing the passing straw, verify that the location of the lead has not changed.
6. Refer to the neurostimulator implant manual for instructions on connecting the lead to the neurostimulator.

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All Rights Reserved
Medtronic B.V.
Medtronic International Ltd.