IFU, en, ru: 445 mm x 200 mm (pa: 435 mm x 190 mm)

KD-*SPIKE*

Spike for multiple withdrawals or infusions of fluids

en – Instructions for use

Intended use:
KD-*SPIKE* is a fluid transfer set, intended to establish a needle-free access to vials/containers to enable sterile withdrawal and/or transfer of fluids intended for infusion, administration or introduction into the body.

Indications:
KD-*SPIKE* is provided in sterile condition for short term use for fluid transfer by multiple needle-free withdrawals or infusions of fluids in a single vial/container and for channelling medication that is intended for parenteral administration in combination with other medical devices.

Intended users:
Qualified healthcare professionals.

Description of the product:
KD-*SPIKE* is available in the following variants:
Green color (for general drugs):
- REF 775714 with snap fit cap, without fluid filter, with 0.1 µm air filter
- REF 775776 with swabable valve, without fluid filter, with 0.1 µm air filter
- REF 775837 with snap fit cap, with integrated 2-way valve, without fluid filter, with 0.1 µm air filter*
- REF 775899 with snap fit cap, with integrated 2-way valve, without fluid filter, with 0.45 µm air filter*


Blue color (for general drugs):
- REF 775721 with snap fit cap, with 5 µm fluid filter, with 0.1 µm air filter
- REF 775783 with swabable valve, with 5 µm fluid filter, with 0.1 µm air filter
- REF 775905 with snap fit cap, with integrated 2-way valve, with 5 µm fluid filter, with 0.1 µm air filter*
- REF 775844 with snap fit cap, with integrated 2-way valve, with 5 µm fluid filter, with 0.45 µm air filter*

Red color (for chemotherapeutic drugs):
- REF 775738 with snap fit cap, with 5 µm fluid filter, with 0.1 µm air filter
- REF 775790 with swabable valve, with 5 µm fluid filter, with 0.1 µm air filter
- REF 775851 with snap fit cap, with integrated 2-way valve, with 5 µm fluid filter, with 0.2 µm air filter*

* Validated with syringes with a nominal capacity up to 20 ml
* Integrated 2-way Valve prevents free flow when the system is turned upside down

Materials used:
Acrylonitrile butadiene styrene (ABS), Polypropyine (PP), Polytetrafluorethylene (PTFE) membrane (for Air filter), Polyethylene Terephthalate (PET) (for fluid filter), ETT copolymer (for spikes for chemotherapeutic drugs), colorant, Polycarbonate (PC), silicone (for swabable valve).

Figure No. 1 **Figure No. 2** **Figure No. 3** **Figure No. 4**

Contraindications:
Not to be used for patients with known hypersensitivity to any of the materials used.

Instructions for use:
1. Swab the vial/container with preferred solution and allow it to dry.
2. Open the individual package and take out the device.
3. Remove the protective cover from the spike.
4. Hold the body of the spike and insert the spike into the center of the vial/container ensuring that the vial/container is on a flat, solid surface (see Figure no 1 with (A) for variants with snap-fit cap and (B) for variants with swabable valve).
5. Ensure the spike is entering vertically and not tilted or at an angle.
6. Once the spike is inserted into the vial/container, proceed as follows.
 A. Variants with snap-fit cap: Open the snap-fit cap.
 B. Variants with swabable valve: Swab the surface of the valve with a sterile alcohol wipe and let it dry.
7. Firmly connect a syringe to the free end of the spike by using a clockwise twisting motion until the syringe enters the spike (see figure no 2 with (A) for variants with snap-fit cap and (B) for variants with swabable valve).
8. Either invert the vial/container and withdraw the desired amount of fluid (see figure no 3 with (A) for variants with snap-fit cap and (B) for variants with swabable valve) or leave the vial/container upright on a flat, solid surface and inject the desired amount of fluid.
9. If used for withdrawal, turn the vial/container after fluid withdrawal has been completed so that it stands upright on a flat, solid surface.
10. Disconnect the syringe from the spike by using an anticlockwise twisting motion (see figure no 4 with (A) for variants with snap-fit cap and (B) for variants with swabable valve).
11. Close the snap-fit cap (in variants with snap-fit cap).

⚠ Cautions:
- Use of this product is restricted to qualified healthcare professionals.

Read instructions before use.
Consult the instructions for use of the medication intended to be used with KD-*SPIKE*.
KD-*SPIKE* has not been tested on resistance against saline containing more than 0.9 % benzyl alcohol. Higher concentrations of benzyl alcohol and other solvents (e.g. polyethylene glycol) as well as different inorganic, surface-active substances might lead to stress and tension cracks in plastics. For these substances only a single aspiration with the spike is recommended. Visually inspect and check the expiry date as well as the product and package integrity prior to use. Improper transport and handling may cause structural and/or functional damage to device or packaging.
Do not use if the individual package is already opened or damaged.
Do not use after the expiry date.
Use the product immediately after opening the individual package.
Do not insert the spike at an angle or in a tilted fashion.
Do not hold the vent while inserting the spike into the vial/container.
Do not inject fluid while the vial/container is inverted.
Do not hold the spike when the vial is inverted without syringe.
To prevent valve damage do not use needles or blunt cannula to access the swabable valve.
For single vial use only. Re-use and cleaning of the product creates a potential risk for patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
Discard after single vial use. Do not recycle, reprocess or re-sterilize.

Storage conditions:
Store in a cool and dry place. Do not expose to heat or direct sunlight.

Disposal:
Discard the product after single vial use in accordance with pertinent national and institutional guidelines.

Notice:
Any serious incident that has occurred in relation to the device must be reported to KD Medical GmbH Hospital Products and to the competent authority of the EU Member State (or state with the same rules (Regulation (EU) 2017/745 on medical devices)) in which the user and/or patient is established.

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