INSTRUCTIONS FOR USE

Disposable Circular Staplers

BEFORE USING THE PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY. IMPORTANT!

This booklet is designed to assist in using this product. It is not a reference to surgical techniques.

This device was designed, tested and manufactured for single use only. Reusing or reprocessing this device may lead to its failure and subsequent patient injury. Reprocessing and/or re-sterilizing this device may create the risk of contamination and patient infection. Do not reuse, reprocess, or re-sterilize this device.

Device description

The Disposable Circular Stapler simultaneously places two staggered, circular rows of staples in the target tissue and resects the tissue inside with a circular blade. The stapler has good anastomosis capability. Compressive stress on the anastomosed tissue should not be less than 3.6 kPa. The staple line should have no leakage or splitting. The safety lock can be engaged and disengaged flexibly, maintaining safe use. The spring of the handle has good elasticity and the handle can return to its original position quickly upon release.


Indications

The Disposable Circular Stapler is intended for end-to-end, end-to-side, and side-to-side anastomosis during digestive tract reconstruction for the esophagus, stomach,and intestines.

Contraindications

- Severe tissue edema;
- Muscle layer of the organ wall is too thick or too thin;
- Lumen is too narrow.

Illustration and Nomenclature

1. Anvil
2. Protector
3. Cartridge
4. Shaft
5. Indicator
6. Fixed Handle
7. Adjusting Knob
8. Safety
9. Active Handle

Device models

DH14, DH17, DH19, DH21, DH24, DH26, DH29, DH32, DH34

Instructions For Use

1. Take the sterile instrument out of the package. Turn the adjusting knob counterclockwise until the protector is visible (see Fig. 1).

 Pull out the anvil (see Fig. 2) and remove the protector (see Fig. 3).

Fig. 1


Fig. 2

Fig. 3

2. Purse-string suture the target tissue as usual but with the stapler anvil inside. Firmly pull the suture string (see Fig.4) and cut off excess tissue, leaving 2-3 mm of tissue beyond the sutures. The thickness of the tissue layers between the anvil and stapler should be between 1.0 mm and 2.5 mm, or the device may not perform successfully.


Fig. 4

3. Assemble the anvil onto the central rod until an audible “click” is heard (see Fig.5). Turn the adjusting knob clockwise to clamp the device on the target tissue. Ensure that the red index in the indicator is within the green range (see Fig. 6).


Fig. 5


Fig. 6

4. After double-checking that the index is in the right position, release the safety (see Fig. 7). Keep the stapler parallel to the tissue segments and squeeze the handles uniformly (see Fig. 8) until the active handle comes into contact with the fixed handle and can be squeezed no further. A cutting “crunch” sound will be heard. Release the handles to complete the procedure, then restore the safety.

Fig. 7



Fig. 8

5. Turn the adjusting knob 2 to 3 rounds counterclockwise to open the anvil (see Fig. 9). Slowly and carefully rotate the device 90° back and forth while pulling, keeping the distal end of the stapler as still as possible, to remove the device from the body (see Fig. 10).

2 to 3 rounds

Fig. 9


Fig. 10

6. Check the resected tissue donuts to ensure the completed staple lines. In the case of a staple line defect, reinforcement is necessary.

Precautions and Warnings: (Please read all information carefully before surgery)

1. Check the sterilized package carefully. Do not use if the package is damaged.
2. Check the product shelf life carefully. Do not use if expired.
3. Endoscopic procedures should be performed only by surgeons with adequate training and familiarity with endoscopic techniques. Before any endoscopic procedures, consult the medical literature relative to techniques, complications, and hazards. Individual consumers shall not use it by themselves.
4. Preoperative radiotherapy may result in changes to tissue. These changes may, for example, cause the tissue thickness to exceed the indicated range for the selected staple size, and cause anastomotic leakage. Careful consideration should be given to any pre-surgical treatment the patient may have undergone and to selection of staple size.
5. Thoroughly clean the tissue segments to be anastomosed before making the purse-string suture.
6. Uneven purse-string suturing will cause bleeding, which must be treated on a case by case basis.
7. When assembling the anvil onto the central rod, ensure that the safety is engaged to avoid accidental firing. Otherwise, the blade may be accidentally exposed and the staples prematurely deployed (partially or fully).
8. The tissue segments to be anastomosed should be positioned evenly between the anvil and cartridge.
9. Do not over-tighten the adjusting knob. Keep the tissue segments in proper orientation.
10. Each firing must be fully completed. Do not fire the device partially, or it may lead to staple malformation, incomplete staple line, bleeding and leakage, or difficulty when removing the device.
11. Caution: After firing, release the active handle and re-engage the safety immediately to avoid an accidental second cut.
12. Once the anastomosis is completed, please check the resected tissue to ensure that there is no leakage or bad hemostasis. If a staple line defect is found, please place a suture on the defective

part as remedy.

13. Dispose of the used devices in accordance with local medical and biological waste disposal requirements.

Packing

The devices are packed in a sterile blister box and paper box.

Environmental conditions for transport and storage

Store the products in a ventilated, room-temperature, dry, and non-corrosive space with humidity under 80%.

 Caution

 Sterilized using ethylene oxide

 Do not use if package is damaged and consult instructions for use

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Add: Building 1, the 3rd floor of Building 3, the 2nd, 3rd, 5th floor of Building 5 & Room 7-3 of Building 7, No.66, Hehuan Road, Zhonglou Economic Development Area, Changzhou, 213000 Jiangsu, P.R. China

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Instructions for Use

CKH

Disposable Linear Staplers and Reloads

Indication

Disposable Linear Stapler and Reload is automatically placing a double-staggered row of staples to suture the tissue. It is simply used to substitute the manual suturing and widely applied for the efficiency and convenience, shortened time of the surgery it may provide with.

Schematic View and Nomenclature


Fig. 1

I.Jaws 2. Anvil 3. Cartridge 4. Cartridge Housing 5. Retaining Pin

6. Active Handle 7. Fixed Handle 8. Retaining Pin Thumb Button 9. Release Button

Product Specifications

Instrument	Cartridge	Tissue Thickness	Color of Cartridge	Height of Staple
CKH303	CH303	1.5mm	blue	3.8mm
CKH304	CH304	2.0mm	green	4.5mm
CKH453	CH453	1.5mm	blue	3.8mm
CKH454	CH454	2.0mm	green	4.5mm
CKH603	CH603	1.5mm	blue	3.8mm
CKH604	CH604	2.0mm	green	4.5mm
CKH753	CH753	1.5mm	blue	3.8mm
CKH754	CH754	2.0mm	green	4.5mm
CKH903	CH903	1.5mm	blue	3.8mm
CKH904	CH904	2.0mm	green	4.5mm

1.Take out the sterilized instrument from the package, Place the intended tissue between the cartridge and the anvil (see Fig.2).


Fig. 2

2.Grasp the handles (see Fig.3) to move backward the active handle to about half of the movable range and release the active handle and keep it.

Adjust the tissue between the cartridge and the anvil to make the tissue adequate and no overlapping.

Further press the active handle to the end (see Fig.4) until the click sound can be heard. The intended tissue is fraped in between the cartridge and the anvil.

Release the active handle to the original place automatically and clamp the tissue tightly between the anvil and cartridge.


Fig. 3


Fig. 4

3. Grasp the handles again, and press the active handle to the end (see Fig.5). Then fire the instrument and suture the tissue. The active handle is locked and complete the firing.


Fig. 5

4.Cut the excessive tissue by scalpel. Press

release button (see Fig.6), and the active handle will be automatically opened. withdraw the retaining pin thumb button and remove the instrument.


Fig. 6

Reloading

1. If one surgical procedure need suturing several times (see Fig 7), move backward the retaining pin thumb button to the end and take out the used cartridge and change the new one (see Fig. 7).


Fig. 7

2. Repeat use steps No.1-4 for another suturing

Warning and Precautions

1. When closing the instrument completely, the thickness of the sutured tissue shall be 1.0-2.5mm.
2. After removing the instrument, examine the staple line for hemostasis and proper staple closure. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.
3. It is single-use instrument and can be only used for one surgical procedure.
4. Make sure not to fall the excessive staples into the body if the width of tissue is smaller than the width of instrument.

Caution

STERILE EO — Sterilized using ethylene oxide

Do not use if package is damaged and consult instructions for use

Do not re-use

Jiangsu Brightness Medical Devices Co., Ltd.
ADD: Building 1, the 3rd floor of Building 3, the 2nd, 3rd, 5th floor of Building 5 & Room 7-3 of Building 7, No.66, Hehuan Road, Zhonglou Economic Development Area, Changzhou, 213000 Jiangsu, P.R. China
Tel: +86-519-83908835 Fax: +86-519-83908831

EU REP **Lepu Medical (Europe) Cooperatief U.A.**
Add: Abe Lenstra Boulevard 36, 8448 JB, Heerenveen, The Netherlands
Tel: +31-515 573399 Fax: +31-515 760020

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
Instructions of Disposable
Linear Staplers and Reloads

【Product Name】Disposable Linear Staplers and Reloads

【Components】

1. Anvil 2.Cartridge 3.Cartridge Housing 4.Shaft 5.Fixed Handle
6. Firing Handle 7.Closure Handle 8.Release Button 9.Staple

Form of Disposable Linear Staplers and Reloads

【Product Specifications】

Instrument	Cartridge	Tissue Thickness	Color of Cartridge	Height of Staple
SAF302	SF302	1.0mm	white	2.5mm
SAF303	SF303	1.5mm	blue	3.8mm
SAF304	SF304	2.0mm	green	4.5mm
SAF452	SF452	1.0mm	white	2.5mm
SAF453	SF453	1.5mm	blue	3.8mm
SAF454	SF454	2.0mm	green	4.5mm
SAF602	SF602	1.0mm	white	2.5mm
SAF603	SF603	1.5mm	blue	3.8mm
SAF604	SF604	2.0mm	green	4.5mm
SAF752	SF752	1.0mm	white	2.5mm
SAF753	SF753	1.5mm	blue	3.8mm
SAF754	SF754	2.0mm	green	4.5mm
SAF902	SF902	1.0mm	white	2.5mm
SAF903	SF903	1.5mm	blue	3.8mm
SAF904	SF904	2.0mm	green	4.5mm

【**Main Technical Index**】Disposable Linear Stapler and Reload is constituted by anvil, cartridge, shaft, closure handle, firing handle, fixed handle, release button, staples. Material: Anvil, Cartridge housing, shaft is made by 06Cr19Ni10, cartridge housing, closure handle, firing handle, fixed handle, release button are made by ABS, staples are made by TA1G or TA2G.

Characters: Cartridge and shaft should match each other and disassembled smoothly, and all moveable parts should be pushed smoothly without stuck and loose. Compressive stress should be endured≧3.6Kpa and no leakage and split. Staples heads should be sharp and acuity. it should be installed into cartridge easy, and the staples should not fall off during swinging the instrument. The spring should be elastic enough.

【**Application**】Used in the surgery of esophagus, stomach and digestive tract


【**Usage**】During the incision surgery or partly cutting surgery of esophagus, stomach and digestive tract, disposable Linear Stapler and Reload can be used as a substitute of handwork. See its instructions as below: ①Take it out form the sterile package, then remove the protector, place the target tissue between anvil and cartridge (as picture 1), clasp the closure handle to pull it back until it reaches its half range, then release the closure handle (not restoration) and adjust the place of target tissue which is clamped by anvil and cartridge to prevent overlap (as picture 2).

② Continue to clasp the closure handle again after adjust the tissue to make it reaches the end when there is a click sound which means the target tissue is clamped by anvil and cartridge well. Then release the closure handle but it is fixed now (as picture 3).


③ Clasp firing handle until it reaches the end, the staples in the cartridge will be fired and suture the tissue. Then release the firing handle, it will restore to its original place (as picture 4).

④ Cutting the superfluous tissue by knife (as picture 5).


⑤ Push the release button, the closure handle will be restored automatically, then stapler will be opened and exited.


Picture 1


Picture 2

Picture 3


Picture 4

Picture 5

P. S.: Methods of cartridge replacement when times of suture are necessary in a same surgery: Take out the cartridge directly and replace a new one, then it can be reused again (as picture 6).

Picture 6

【Attentions】

1. Surgeons who operate must be qualified.
2. There is a safety ability of it, the staples will be fired without cartridge implant.
3. Tissue's thickness should be less than 1.0mm-2.5mm when the instrument is totally clamped.
4. Please check the cutting point tissue, if damage or bleeding is found, you should do some repairing suture or electrocoagulation.
5. Please pay attention that the surplus staples should not be left inside patient's body when tissue's width is narrow than the cartridge's.
6. Disposable Linear Stapler and Reload is single use only and sterilized by EO, the package is sterile. and there is "Sterile" marked in package, and there are marks of STERILE EO , sterilization date and warnings of no usage of broken-packed product.
7. Check the package before use and never use it if the package is broken.
8. Please destroy it after use according to regulations.

【**Packing**】 Inner package is a sterile and seal blister box ,and middle package is a paper box, there is one piece of product in every paper box together with quality certification and use instruction, the out package is a corrugated paper carton.

Instructions of marks in the package: STERILE EO Sterilized using ethylene oxide

Do not re -use

【**Storage**】

Put the products in ventilated, drying and non-corrosive room with the humidity under 80%.

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ADD: Building 1, the 3rd floor of Building 3, the 2nd, 3rd, 5th floor of Building 5 & Room 7-3 of Building 7, No.66, Hehuan Road, Zhonglou Economic Development Area, Changzhou, 213000 Jiangsu, P.R. China

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Intended Purpose

This device is intended for use in abdominal, thoracic, and gynaecologic surgery for resection, transection, and creation of anastomosis.

Description

This is a sterile, disposable, single-patient use device. It places two triple rows of staggered staples and simultaneously cuts tissue between the two rows, completing the tissue anastomosis. Staple size and staple line length are determined by the selection of the single-use reloads. The device may be used in open or laparoscopic surgery.

Any model of the JQ series stapler may be used with any model of JQ series and JJII series reloads. JQ series staplers may be reloaded and fired up to 25 times in a single procedure.

The JQ series reloads offer staples of uniform height, with options of 2.0 mm, 2.5 mm, 3.5 mm, 4.0 mm, or 4.8 mm.

The JJII series reloads offer staples with progressive heights. The triple rows at either side of the cut-line are symmetrical, with the row closest to the cut-line having the shortest staples and the farthest row having the tallest staples. The staple height increases by 0.5mm per row. The staple height of the closest row may be selected from the options of 2.0mm, 3.0mm, or 4.0mm, depending on tissue thickness. Progressive height staples enable better blood circulation in the anastomosis area.

The short model of JQ staplers is designed for open procedures only. The instructions, warnings, and precautions referring to trocars or endoscopic procedures are not applicable to short model staplers.

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Indications


The device is indicated for abdominal, thoracic, and gynaecologic surgeries that use stapled resection, transection, or creation of anastomosis. Examples of procedures include: gastrectomy, appendectomy, pulmonary wedge resection, etc.

Contraindications


- Do not use the device where surgical stapling is not appropriate;
- Do not use the short shaft JQ-B & JQIIB model staplers for endoscopic procedures; they are designed for open procedures only;
- Do not use the device where adequacy of haemostasis cannot be visually verified after application;
- Do not use the device on the aorta and major blood vessels;
- Do not use the device on ischemic, oedematous, friable, or necrotic tissue, as this could result in staple line leak or bleeding;
- Do not use the reloads on tissue whose thickness cannot be comfortably compressed to less than the maximum value, or can be easily compressed to less than the minimum value of the corresponding "tissue thickness range" from the table in "Reload Variants Specifications". Failure to respect the appropriate tissue thickness may lead to tissue damage or anastomosis leakage;

Stapler Schematic Illustration

JQ model stapler:

JQII model stapler:


1. **Pin:** connection/disconnection between the stapler and the reload;
2. **Internal rod:** supports internal mechanical parts;
3. **Shaft:** penetration through trocar sleeve or incision;
4. **Rotation collar:** to rotate the shaft freely 360° in both directions;
5. **Safety button:** safety lock to prevent unintentional firing, to be pressed just before firing;
6. **Black return knob:** to open the reload jaws and return the knife to home position;
7. **Fixed handle:** to hold the stapler with one hand;
8. **Movable handle:** for JQ model, to close or fire the reload; for JQII model, to open/close or fire the reload;
9. **Articulating lever:** to pivot the distal portion of the articulating model reload left or right, up to 45° on each side;
10. **Unload button:** to unload the reload.


JQII model offers a different handling option. A closed reload may be opened by pushing the moveable handle forward, enabling the surgeon to single-handedly reposition the device on a different site. In comparison, the JQ model stapler requires the surgeon to pull the return knob with another hand.

Reload Schematic Illustration


JQZ model straight reload:


JJIIZ model straight reload:

JQD model articulating reload:

JJIID model articulating reload:

1. Cartridge: with two triple rows of staples; **2. Anvil:** forms the staples into "B" shapes; **3. Cartridge holder:** houses the cartridge; **4. Knife position indicator:** indicates the knife position; **5. Shaft protector:** protects mechanical parts.

Stapler Variants Specifications

Staplers models and variants are coded according to the following rules:

(JQ/JQII)–(A/B/L)

- Stapler type: A=Standard, B=Short, L=Long
- Series and model code

Stapler model & variant	Shaft length (mm, ±10)	Total stapler length (mm, ±10)	Specifications
JQ-A	155	370	JQ stapler with standard-length shaft
JQ-B	105	320	JQ stapler with short shaft
JQ-L	255	470	JQ stapler with long shaft
JQIIA	155	378	JQII stapler with standard-length shaft
JQIIB	105	328	JQII stapler with short shaft
JQIIL	255	478	JQII stapler with long shaft

Reload Variants Specifications

JQ series reload models and variants are coded according to the following rules:

JQ(Z/D)–(30/45/60)–(2.0/2.5/3.5/4.0/4.8)

- Open staple height (mm)
- Staple line length (mm)
- Reload model: Z=straight, D=articulating

JJII series reload models and variants are coded according to the following rules:

JJII(Z/D)–(30/45/60)–(1/2/3)

Code	Open staple height (mm)
1	2.0, 2.5, 3.0
2	3.0, 3.5, 4.0
3	4.0, 4.5, 5.0

- Staple line length (mm)
- Reload model: Z=straight, D=articulating

Reload model: JQZ, straight reload						
Reload variant (code)	Cartridge colour	Open staple height (mm,±0.2mm)	Tissue thickness range (mm)	Staple line length (mm,±2.0mm)	Reload length (mm±5.0mm)	Number of staples
JQZ-30-2.0	Gray	2.0	0.75 – 1.00	30.0	170	48
JQZ-30-2.5	White	2.5	1.00 – 1.50			
JQZ-30-3.5	Blue	3.5	1.50 – 2.00			
JQZ-30-4.0	Gold	4.0	1.70 – 2.20			
JQZ-30-4.8	Green	4.8	2.00 – 2.50		178	
JQZ-45-2.0	Gray	2.0	0.75 – 1.00	43.0	184	66
JQZ-45-2.5	White	2.5	1.00 – 1.50			
JQZ-45-3.5	Blue	3.5	1.50 – 2.00			
JQZ-45-4.0	Gold	4.0	1.70 – 2.20			
JQZ-45-4.8	Green	4.8	2.00 – 2.50		193	
JQZ-60-2.0	Gray	2.0	0.75 – 1.00	60.0	199	90
JQZ-60-2.5	White	2.5	1.00 – 1.50			
JQZ-60-3.5	Blue	3.5	1.50 – 2.00			
JQZ-60-4.0	Gold	4.0	1.70 – 2.20			
JQZ-60-4.8	Green	4.8	2.00 – 2.50		209	

Reload model: JQD, articulating reload						
Reload variant (code)	Cartridge colour	Open staple height (mm,±0.2mm)	Tissue thickness range (mm)	Staple line length (mm,±2.0mm)	Reload length (mm±5.0mm)	Number of staples
JQD-30-2.0	Gray	2.0	0.75 – 1.00	30.0	205	48
JQD-30-2.5	White	2.5	1.00 – 1.50			
JQD-30-3.5	Blue	3.5	1.50 – 2.00			
JQD-30-4.0	Gold	4.0	1.70 – 2.20			
JQD-30-4.8	Green	4.8	2.00 – 2.50		213	
JQD-45-2.0	Gray	2.0	0.75 – 1.00	43.0	217	66
JQD-45-2.5	White	2.5	1.00 – 1.50			
JQD-45-3.5	Blue	3.5	1.50 – 2.00			
JQD-45-4.0	Gold	4.0	1.70 – 2.20			
JQD-45-4.8	Green	4.8	2.00 – 2.50		225	
JQD-60-2.0	Gray	2.0	0.75 – 1.00	60.0	233	90
JQD-60-2.5	White	2.5	1.00 – 1.50			
JQD-60-3.5	Blue	3.5	1.50 – 2.00			
JQD-60-4.0	Gold	4.0	1.70 – 2.20			
JQD-60-4.8	Green	4.8	2.00 – 2.50		241	

Reload model: JJIIZ, straight reload								
Reload variant (code)	Cartridge colour	Open staple heights (mm, ±0.2mm 1st row: nearest to cutline, 3rd row: farthest to cutline)			Tissue thickness range (mm)	Staple line length (mm,±2.0mm)	Reload length (mm±5.0mm)	Number of staples
		1st row	2nd row	3rd row				
JJIIZ301	Tan	2.0	2.5	3.0	0.88 – 1.50	30.0	202	48
JJIIZ302	Purple	3.0	3.5	4.0	1.50 – 2.25			
JJIIZ303	Black	4.0	4.5	5.0	2.25 – 3.00			
JJIIZ451	Tan	2.0	2.5	3.0	0.88 – 1.50	43.0	217	66
JJIIZ452	Purple	3.0	3.5	4.0	1.50 – 2.25			
JJIIZ453	Black	4.0	4.5	5.0	2.25 – 3.00			
JJIIZ601	Tan	2.0	2.5	3.0	0.88 – 1.50	60.0	232	90
JJIIZ602	Purple	3.0	3.5	4.0	1.50 – 2.25			
JJIIZ603	Black	4.0	4.5	5.0	2.25 – 3.00			

Reload model: JJIID articulating reload								
Reload variant (code)	Cartridge colour	Open staple heights (mm, ±0.2mm 1st row: nearest to cutline, 3rd row: farthest to cutline)			Tissue thickness range (mm)	Staple line length (mm,±2.0mm)	Reload length (mm±5.0mm)	Number of staples
		1st row	2nd row	3rd row				
JJIID301	Tan	2.0	2.5	3.0	0.88 – 1.50	30.0	202	48
JJIID302	Purple	3.0	3.5	4.0	1.50 – 2.25			
JJIID303	Black	4.0	4.5	5.0	2.25 – 3.00			
JJIID451	Tan	2.0	2.5	3.0	0.88 – 1.50	43.0	217	66
JJIID452	Purple	3.0	3.5	4.0	1.50 – 2.25			
JJIID453	Black	4.0	4.5	5.0	2.25 – 3.00			
JJIID601	Tan	2.0	2.5	3.0	0.88 – 1.50	60.0	232	90
JJIID602	Purple	3.0	3.5	4.0	1.50 – 2.25			
JJIID603	Black	4.0	4.5	5.0	2.25 – 3.00			

Implanted Materials

The staples are made of titanium or titanium alloy Ti-3Al-2.5V.

Warnings and Precautions

- The device is intended to be used by surgeons who have adequate training and familiarity with open and laparoscopic techniques.
- Prior to any procedure, surgeons must read and carefully follow this IFU, consult the medical literature relative to techniques, complications and hazards, be familiar with the structures of and around the intended stapling site. Do not use the device where adequacy of haemostasis cannot be visually verified after application.
- Consider other options if the patient's tissue is oedematous (swollen with fluid), friable (tissue that readily tears, fragments, or bleeds when gently palpated or manipulated), necrotic (dead tissue), or if it has received any preoperative chemotherapy or radiation, as the staples may be less likely to securely close tissue. Use in these circumstances may result in staple line leak or bleeding, prolonged surgery time, or inflammatory reaction.
- The device is packaged and sterilized by Ethylene Oxide (EO) for use in a single procedure only. Check the integrity of the sterilized package before use. Discard if the package is not sealed. Use of an unsealed device may lead to pathogen transfer, contamination of the patient, infection, or inflammation.
- Check the date of validity before use, and discard if the device is expired. Use of an expired device may lead to pathogen transfer, contamination of the patient, infection, or inflammation.
- Reprocessing, re-sterilization, or reuse of a disposable device may damage the device's structural integrity, pose a risk of contamination / cross-infection to the patient, result in injury, illness or even death of patient.
- Carefully evaluate the tissue thickness. Have a range of staple sizes available and select the appropriate size reload for the tissue type and thickness. Overly thick or thin tissue may result in unacceptable staple formation, which may lead to staple line leak, bleeding, or prolonged surgery time.
- Do not attempt to remove the protective cap until the reload is properly loaded onto the stapler.
- Before loading or unloading, the jaws must always be set in open position. Additionally, for articulating JQD or JJIID reloads, the jaws must be parallel to the shaft (neutral position). Attempting to load/unload a reload with closed or articulated jaws may damage the load/unload mechanism and make multi-reloading no longer possible.
- Before introducing the stapler into the trocar or incision, always perform a close-open cycle of the jaws to confirm proper loading and clamping mechanism.
- Prior to introducing or removing the stapler from the trocar or incision, the jaws must always be in closed position. Additionally, articulating models must be in neutral position. Failure to set the jaws in closed neutral position will result in difficult insertion or withdrawal of the device.
- For articulating reloads, always use the stapler's articulating lever to regulate the articulation angle. Do not exert direct force on the jaws, as this may damage the articulating mechanism and make unloading difficult or even impossible without damaging the stapler.
- The anvil must be completely past the trocar (if a trocar is used) prior to opening the reload within the body cavity.
- Placement of tissue proximal to the tissue stops (on the reload) may result in stapler malfunction. Any tissue extending beyond the black "cut" line on the reload will not be transected, all tissue placed within the "cut" line will be transected.
- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the jaws, which could lead to staple line leak, bleeding, or prolonged surgery time.
- Ensure that the tissue lies flat and is evenly distributed inside the stapling range prior to closing the jaws. "Bunching" tissue may result in staple malformation, poor haemostasis and/or leakage.
- Always completely fire the reload. Partial firing leads to incomplete cut and/or staple malformation, which may result in poor haemostasis and/or leakage.
- The device is equipped with a safety interlock to prevent an empty or partially-fired reload from being fired a second time. Do not attempt to override the safety interlock, as it will cause device malfunction.
- If the clamping or firing mechanism fails to work (excessive resisting force), do not continue the operation. Remove and carefully check the stapler-reload structural integrity following the IFU. The use of excessive force may result in device damage, injury, or tissue damage.
- After firing and removal of the device, always inspect the staple line and the surrounding site for air sealing / haemostasis of the anastomosis line and to see whether the staples are closed properly. Minor bleeding or leakage may be treated by electrocautery, manual suture or other appropriate methods.
- In a single procedure, any JQ series stapler may be reloaded and fired up to 25 times with any JQ or JJII reload.
- Dispose of all opened devices whether used or unused. Devices which come into contact with bodily fluids require special disposal handling to prevent biological contamination, relevant local regulations shall be respected.
- The implantable staples are MR conditional, see the MRI Safety Information section for more information.

Instructions for Use

Prior to use

1. Carefully evaluate the length and thickness of the tissue to be excised, then select one or several reloads with the appropriate staple sizes from the table in "Reload Variants Specifications".
2. Check the integrity of the sterile barrier of the packaging, and the device expiration date. Do not use any device with damaged sterile barrier or device past expiration date.
3. Open the blister box of the stapler, ensure that the black return knob on the stapler is pulled back completely, and that the articulating lever is in neutral position.
4. Open the blister box of the reload, check that the jaws are in open position, and the protective cap remains wedged. Discard the reload if the cap is not in place.

Loading the reload

5. Align the "LOAD" indicator on the reload with the "LOAD" indicator on the stapler shaft, push the reload into the stapler's pin and twist clockwise 45° relative to the stapler, so that the reload locks into place and an audible click is heard. Do not squeeze the movable handle during the loading, or prior to removing the protective cap.
6. Remove the reload protective cap, check for the presence of all staples.
7. Perform a close-open cycle to confirm proper loading: squeeze the movable handle once to close the jaws, then pull back the black return knob to open the jaws. For JQII model staplers, an additional way to open the jaws is to push the movable handle forward.

During use

8. Close the jaws, ensure the articulating lever's neutral position, and insert the device into an appropriately sized trocar. Note: the green JQ reload (4.8mm) and black JJII reload (4.0/4.5/5.0mm) must be inserted through a 15mm trocar sleeve (or larger with a converter), while other reloads are designed for introduction through a 12mm trocar sleeve (or larger with a converter).
9. Once inside the body cavity, ensure that the anvil is completely visible. Open the jaws by pulling back the black return knob (do not squeeze the movable handle at same time). If necessary, use the shaft's 360° rotation collar and (for articulating reload) the articulating lever to articulate the reload up to 45° in either direction. Position the device across the target tissue, ensure that the tissue to be transected lies flat and evenly, with no obstructions within the jaws. Close the jaws by fully squeezing the movable handle. It is recommended to wait at least 15s before firing. The jaws of the device may be repositioned on the tissue prior to firing by pulling back the black return knob, allowing the jaws to open.
10. Push the green safety button, squeeze the movable handle fully and repeatedly until the knife position indicator reaches the distal end to completely fire the reload. The total number of squeezes depends on the reload length (30mm, 45mm or 60mm). Failure to completely fire the reload will lead to incomplete cut and/or staple malformation, which may result in poor haemostasis and/or leakage.
11. Once the reload has been completely fired, it is recommended to keep the jaws closed and wait at least 15s. Open the jaws by pulling the black return knob back completely. Gently remove the device from the tissue, close the jaws, ensure the articulating lever's neutral position, and withdraw the device out of the body cavity.
12. Carefully inspect the staple line for proper staple formation and haemostasis. Minor bleeding or leakage may be treated by electrocautery, manual suture or other appropriate methods.

Unloading the reload

13. Ensure that the jaws are open and in neutral position. Press the "UNLOAD" button back towards the stapler, twist the reload counter clockwise 45° and remove the reload from the shaft of the stapler.
14. Additional reloads (up to 25 times in total) may be then loaded and fired in a single procedure. A safety interlock prevents an empty reload from being fired a second time. Do not attempt to override the safety interlock.

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- Static magnetic field (B0) no more than 3.0T;
- Maximum allowable spatial gradient magnetic field (∇B) is 30T/m (3,000 gauss/cm);
- Magnetic induced torque is less than worst-case gravitational torque;
- In 3.0T or 1.5T MR system, whole-body average Specific Absorption Rate (SAR) no more than 2W/kg, RF exposure no more than 15mins;
- Artifact information, sample causes artifacts up to 3.8mm in 3.0T MR system and 3.2mm in 1.5T MR system.

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Refer to package label.

Storage and Handling Conditions

Store at room temperature, avoid prolonged exposure to elevated temperatures. Do not expose to temperatures above 30°C or below 10°C, do not expose to humidity above 70% or below 35%. Handle with care, do not over-stack or over-palletize.

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Symbol Description

	Date of manufacture		Manufacturer
REF	Catalogue number	EU REP	Authorized representative
LOT	Batch code	STERILE EO	Sterilized using ethylene oxide
UDI	Unique device identifier		Single sterile barrier system with protective packaging outside
MD	Medical device		Single sterile barrier system
	Use-by date		Do not re-sterilize
	Caution		Do not re-use
	Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Upper limit of temperature		Keep away from sunlight
	Fragile, handle with care		Keep dry
	Humidity limitation	MR	MR conditional

Jiangsu Brightness Medical Devices Co., Ltd.
Building 1, the 3rd floor of Building 3, the 2nd, 3rd, 5th floor of Building 5 & Room 7-3 of Building 7, No.66, Hehuan Road, Zhonglou Economic Development Area, Changzhou, 213000 Jiangsu, P.R. China

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Disposable Linear Cutter Staplers and Reloads
Instruction for Use

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

【Main Technical Information】

Material: Staple: TA1,TA2, Anvil: 06Cr19Ni10 stainless steel, Cutter: 6Cr13 stainless steel, Handle, Releasing Button: ABS.

Capability: The quality disposable Linear cutter is able to provide good performance of cutting and suturing, the sutured tissue is allowed under pressure of more than 3.6Kpa without leaking and tearing appearing. Clear scale mark should be presented on the instrument, and its retractable handle can open and close flexibly. Its safety lock can secure from mis-firing, ensure to provide safe performance. The products is sterilized and can be used directly after unpacking.

Instrument Description

1 Releasing Button 2 Retractable handle 3 Trigger 4 Cartridge Rack

5.Cartridge 6 Scalpel 7 Anvil 8.Anvil Handle 9.Staple

Remarks: 1pcs cartridge is along with each Disposable Linear Cutter as standard kit when delivery.

【Application】 Applied for gastrointestinal and pulmonary surgery.

【Instruction for Use】

1.After Press the Trigger to open the retractable handle, separate the cutting cartridge and anvil, and then remove the protector from the cartridge gently.

2. Place tissue between cutting assembly and anvil, then press the retractable handle to the locked the part(Note: The" Click" sound will be heard after the retractable handle press completely in place).

3. Carry the disposable linear cutter with one hand while press&push forward the releasing button gently to fire position by thumb and index finger of another hand.(The frictional force is not so obvious when the cutting haven't acted). And then keep push toward the releasing button completely to end position. (The frictional force increase when push the trigger that means transection and suturing acting to the tissue). According to the operation habits, push forward one side releasing button from initial fire position to select being fired by left hand or right hand.

4. Don't lift and pull the instrument immediately after withdraw the releasing button to initial position. The instrument should only be taken out after separate the cutting assembly and the anvil.

5. Replace the cartridge: withdraw the releasing button completely and press the trigger to open the retractable handle. And the replace the new cartridge to continue the surgery operation.

【Warning and Precautions】

①. The user should only be qualified surgeon.

②. The disposable linear cutter is only for single use, adopt ethylene oxide sterilization, sterilized and sealed package when delivery,marked the warnings of STERILE EO ， ， date of sterilization, and damaged package is prohibited from using on the package box.

③. According to the surgery need,1 to 3pcs cutting cartridges is required to one surgery.

④. Check the seal integrity of package before use, damaged package and over the period of valid shall be prohibited from using.

⑤. If the instrument is found rusted before use, it shall be prohibited strictly from using.

⑥. Please destroy after use according to regulation.

【Storage】 Instrument should be placed in ventilated, dry and non-corrosive indoor environment, and prevent it from being crushed, worn and impacted.

Jiangsu Brightness Medical Devices Co., Ltd

Address: Building 1, the 3rd floor of Building 3, the 2nd, 3rd, 5th floor of Building 5 & Room 7-3 of Building 7, No.66, Hehuan Road, Zhonglou Economic Development Area, Changzhou, 213000 Jiangsu, P.R. China

Tel: +86-519-83908835 Fax: +86-519-83908831

EU REP **Lepu Medical (Europe) Cooperatief U.A.**

Address: Abe Lenstra Boulevard 36, 8448 JB, Heerenveen, The Netherlands

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Instructions for Use
Disposable Linear Cutter Staplers and Reloads

【Product Name】 Disposable Linear Cutter Staplers and Reloads

【CE certificate No】 HD 60086113 0001


【Product Specification】

Instrument Code No.	Cartridge Code No.	Tissue Thickness	Number of staples	Cartridge Color
GH553	GK553	1.5mm	56	Blue
GH554	GK554	2.0 mm	56	Green
GH753	GK753	1.5mm	76	Blue
GH754	GK754	2.0 mm	76	Green
GH1003	GKI003	1.5 mm	100	Blue
GH1004	GKI004	2.0 mm	100	Green


Application

Applied for gastrointestinal and pulmonary surgery. Description

Instrument Description

Exploded View

1.Cartridge 2.Cartridge Housing 3.Anvil

4.Protector 5.Scale 6.Handle 7.Button

Instructions for Use

Select instrument with blue or green cartridge according to tissue thickness.Take out the sterilize instrument from the package, pull out active handle to divide the instrument into two parts.

remove the protector from the cartridge, place tissue between anvil and cartridge and close active handle. Push button forward, the instrument will suture the tissue and cut it in between. Release the Stapleline push-button to the original position, pull out active handle to divide the instrument into two. Do no retire the instrument. Take out the sutured and cut tissue, examine the staple line for hemostasis and proper staple closure. If hemostasis is not present, appropriate techniques should be used to achieve hemostasis.

If there is a need to change new cartridge, release push-button to the original position and pull out active handle. Take out the used cartridge by Lilting up the end of cartridge and reload the new one. remove the protector from the cartridge.

The instrument is now ready for use.

Warning and Precautions

1. The instrument can be used for one patient several times, dispose it after using. Do not resterilize.
2. Cartridge is only single use.
3. The use of the instrument should only be done by well experienced surgeon.
4. Make sure to remove the protector from the cartridge before firing instrument.

Contraindication

1. Don’t use the instrument with blue cartridge on any tissue that requires excessive force to compress to 1.5mm or on any tissue that compresses easily to below 1.5mm.
2. Don’t use the instrument with green cartridge on any tissue that requires excessive force to compress to 2.0 mm or on any tissue that compresses easily to below 2.0mm.
3. Don’t use on ischemic or necrotic tissue.
4. Don’t use on pulmonary vessel.
5. Don’t use on uterine vessels.

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Jiangsu Brightness Medical Devices Co., Ltd.

BL DISPOSABLE ENDOCUTTER & DISPOSABLE ENDOCUTTER RELOAD

CE 0344

Instructions for Use

PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE USE. FAILURE TO PROPERLY FOLLOW THE INSTRUCTIONS MAY LEAD TO SERIOUS SURGICAL CONSEQUENCES, SUCH AS STAPLE MALFORMATION, LEAKAGE, OR DISRUPTION OF OPERATION.


This document aims to assist surgeons in using this product; it is not a reference to surgical procedure.

This device is for single-use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reprocessing and/or re-sterilization of this device may create the risk of contamination and patient infection. Do not reuse, reprocess or re-sterilize this device.


Description

This is a sterile, disposable single-patient use device intended for use in abdominal, thoracic, and gynaecologic surgery for resection, transection, and creation of anastomosis. It places two triple rows of staggered staples and simultaneously cuts tissue between the two rows, completing the tissue anastomosis. Staple size and staple line length are determined by the selection of the single-use reloads. The device may be used in open or laparoscopic surgery.

BLR reloads offer staples of uniform height or of progressive heights, as illustrated below.

1) Cartridge; 2) Cartridge holder;
3) Staple pusher; 4) Staples of uniform height

1) Cartridge; 2) Cartridge holder;
3) Staple pusher; 4) Staples of progressive heights

Uniform staple height reloads are available in 2.0mm, 2.5mm, 3.5mm, 3.8mm, 4.1mm, or 4.3mm, while progressive staple height reloads are available in 2.5/3.0/3.5mm and 3.8/4.1/4.3mm. Progressive staple height reloads enable better blood circulation in the anastomosis.

45mm staplers may only be used with 45mm reloads; 60mm staplers may only be used with 60mm reloads.

The BL stapler may be reloaded and fired up to 12 times in a single procedure.

Intended Users and Target Patient Group

The device is intended to be used by surgeons who have adequate training and familiarity with open and laparoscopic techniques. The target patient population consists of adults undergoing abdominal, thoracic, or gynaecologic surgery with stapled resection, transection, or creation of anastomosis.

Indications


The device is indicated for abdominal, thoracic, and gynaecologic surgeries that use stapled resection, transection, or creation of anastomosis. It simultaneously places two triple rows of staggered staples and cuts tissue between the two rows, completing the tissue anastomosis. Examples of procedures include: gastrectomy, appendectomy, pulmonary wedge resection, etc.

Contraindications


- Do not use the device where surgical stapling is not appropriate;
- Do not use the device where adequacy of haemostasis cannot be visually verified after application;
- Do not use the device on the aorta and major blood vessels;
- Do not use the device on ischemic, oedematous, friable, or necrotic tissue, as this could result in staple line leak or bleeding;
- Do not use the reloads on tissue whose thickness cannot be comfortably compressed to less than the maximum value, or can be easily compressed to less than the minimum value of the corresponding “tissue thickness range” from the table in “Reload Variants Specifications”. Failure to respect the appropriate tissue thickness may lead to tissue damage or anastomosis leakage;

Stapler Schematic Illustration

BLE model stapler:

BLC model stapler:


1. **Anvil**: to form the staples into “B” shape;
2. **Reload holder**: to house the reload;
3. **Shaft**: penetration through trocar sleeve or incision;
4. **Rotation collar**: to rotate the shaft freely 360° in both directions;
5. **Articulating fins**: apply lateral force to articulate the jaws to the desired angle, release to maintain them in the desired angle, and pull to reset to neutral position;
6. **Manual knife reverse switch**: allows to pull back the knife when squeezing the firing handle;
7. **Knife position window**: to show the knife position in a character wheel;
8. **Knife direction indication**: to indicate the knife direction (forward/backward);
9. **Release button**: to release the anvil and open the jaws;
10. **Fixed handle**: to hold the stapler;
11. **Closing handle**: to clamp the jaws on the tissue;
12. **Firing handle**: to fire the reload, or (when the knife reverse switch is activated) to pull back the knife.

The reload holder may be protected to facilitate the stapler's movement in the body cavity.

Back-view of a basic/bare reload holder structure

Back-view of a protected reload holder structure

Stapler Variants Specifications

BL series stapler models and variants are coded according to the following rules:


BL – (E/C) – (B/P) – (S/M/L) – (45/60)

- Jaw size: 45 = use with 45mm reload, 60 = use with 60mm reload
- Shaft length: S = short, M = medium, L = long
- Reload holder structure: B = basic/bare, P = protected
- Anvil tip shape: E = normal end, C = curved end

The BLC model curved anvil tip provides ease of manipulation and better access underneath tissue/vessels.


Stapler model & variant (code)	Shaft length (mm, ±15mm)	Total stapler length (mm, ±20mm)	Specifications
BLEBS45	190	455	Normal anvil tip, basic reload holder, 45mm, short shaft
BLEBS60	190	468	Normal anvil tip, basic reload holder, 60mm, short shaft
BLEBM45	250	515	Normal anvil tip, basic reload holder, 45mm, medium shaft
BLEBM60	250	528	Normal anvil tip, basic reload holder, 60mm, medium shaft
BLEBL45	350	615	Normal anvil tip, basic reload holder, 45mm, long shaft
BLEBL60	350	628	Normal anvil tip, basic reload holder, 60mm, long shaft
BLEPS45	190	455	Normal anvil tip, protected reload holder, 45mm, short shaft
BLEPS60	190	468	Normal anvil tip, protected reload holder, 60mm, short shaft
BLEPM45	250	515	Normal anvil tip, protected reload holder, 45mm, medium shaft
BLEPM60	250	528	Normal anvil tip, protected reload holder, 60mm, medium shaft
BLEPL45	350	615	Normal anvil tip, protected reload holder, 45mm, long shaft
BLEPL60	350	628	Normal anvil tip, protected reload holder, 60mm, long shaft
BLCBS45	190	455	Curved anvil tip, basic reload holder, 45mm, short shaft
BLCBS60	190	468	Curved anvil tip, basic reload holder, 60mm, short shaft
BLCBM45	250	515	Curved anvil tip, basic reload holder, 45mm, medium shaft
BLCBM60	250	528	Curved anvil tip, basic reload holder, 60mm, medium shaft
BLCBL45	350	615	Curved anvil tip, basic reload holder, 45mm, long shaft
BLCBL60	350	628	Curved anvil tip, basic reload holder, 60mm, long shaft
BLCPS45	190	455	Curved anvil tip, protected reload holder, 45mm, short shaft
BLCPS60	190	468	Curved anvil tip, protected reload holder, 60mm, short shaft
BLCPM45	250	515	Curved anvil tip, protected reload holder, 45mm, medium shaft
BLCPM60	250	528	Curved anvil tip, protected reload holder, 60mm, medium shaft
BLCPL45	350	615	Curved anvil tip, protected reload holder, 45mm, long shaft
BLCPL60	350	628	Curved anvil tip, protected reload holder, 60mm, long shaft

Reload Schematic Illustration

Reload Variants Specifications

BLR series reload models and variants are coded according to the following rules:

Reload model: BLRS, smooth surface reload								
Reload variant (code)	Open staple heights (mm, ±0.2mm, 1st row: closest to cutline, 3rd row: farthest to cutline)			Cartridge colour	Tissue thickness range (mm)	Staple line length (mm, ±2.0mm)	Reload length (mm, ±5.0mm)	Number of staples
	1st row	2nd row	3rd row					
BLRSA45	2.0			Gray	0.75 – 1.0	49.3	76	70
BLRSW45	2.5			White	1.0 – 2.0	49.3	76	70
BLRSB45	3.5			Blue	1.5 – 2.4	49.3	76	70
BLRSY45	3.8			Gold	1.8 – 3.0	49.3	76	70
BLRSG45	4.1			Green	2.0 – 3.3	49.3	76	70
BLRSH45	4.3			Black	2.3 – 4.0	49.3	76	70
BLRSA60	2.0			Gray	0.75 – 1.0	61.3	88	88
BLRSW60	2.5			White	1.0 – 2.0	61.3	88	88
BLRSB60	3.5			Blue	1.5 – 2.4	61.3	88	88
BLRSY60	3.8			Gold	1.8 – 3.0	61.3	88	88
BLRSG60	4.1			Green	2.0 – 3.3	61.3	88	88
BLRSH60	4.3			Black	2.3 – 4.0	61.3	88	88
BLRSL45	2.5	3.0	3.5	Tan	1.0 – 2.4	49.3	76	70
BLRSM45	3.8	4.1	4.3	Purple	1.8 – 4.0	49.3	76	70
BLRSL60	2.5	3.0	3.5	Tan	1.0 – 2.4	61.3	88	88
BLRSM60	3.8	4.1	4.3	Purple	1.8 – 4.0	61.3	88	88
Reload model: BLRP, gripping surface reload								
Reload variant (code)	Open staple heights (mm, ±0.2mm 1st row: closest to cutline, 3rd row: farthest to cutline)			Cartridge colour	Tissue thickness range (mm)	Staple line length (mm, ±2.0mm)	Reload length (mm, ±5.0mm)	Number of staples
	1st row	2nd row	3rd row					
BLRPA45	2.0			Gray	0.75 – 1.0	49.3	76	70
BLRPW45	2.5			White	1.0 – 2.0	49.3	76	70
BLRPB45	3.5			Blue	1.5 – 2.4	49.3	76	70
BLRPY45	3.8			Gold	1.8 – 3.0	49.3	76	70
BLRPG45	4.1			Green	2.0 – 3.3	49.3	76	70
BLRPH45	4.3			Black	2.3 – 4.0	49.3	76	70
BLRPA60	2.0			Gray	0.75 – 1.0	61.3	88	88
BLRPW60	2.5			White	1.0 – 2.0	61.3	88	88
BLRPB60	3.5			Blue	1.5 – 2.4	61.3	88	88
BLRPY60	3.8			Gold	1.8 – 3.0	61.3	88	88
BLRPG60	4.1			Green	2.0 – 3.3	61.3	88	88
BLRPH60	4.3			Black	2.3 – 4.0	61.3	88	88
BLRPL45	2.5	3.0	3.5	Tan	1.0 – 2.4	49.3	76	70
BLRPM45	3.8	4.1	4.3	Purple	1.8 – 4.0	49.3	76	70
BLRPL60	2.5	3.0	3.5	Tan	1.0 – 2.4	61.3	88	88
BLRPM60	3.8	4.1	4.3	Purple	1.8 – 4.0	61.3	88	88

Implanted Materials

The staples are made of titanium alloy Ti-3Al-2.5V.

Warnings and Precautions

- The device is intended to be used by surgeons who have adequate training and familiarity with open and laparoscopic techniques;
- Prior to any procedure, surgeons must read and carefully follow this IFU, consult the medical literature relative to techniques, complications and hazards, be familiar with the structures of and around the intended stapling site. Do not use the device where adequacy of haemostasis cannot be visually verified after application.
- Consider other options if the patient's tissue is oedematous (swollen with fluid), friable (tissue that readily tears, fragments, or bleeds when gently palpated or manipulated), necrotic (dead tissue), or if it has received any preoperative chemotherapy or radiation, as the staples may be less likely to securely close tissue. Use in these circumstances may result in staple line leak or bleeding, prolonged surgery time, or inflammatory reaction.
- The device is packaged and sterilized by Ethylene Oxide (EO) for use in a single procedure only. Check the integrity of the sterilized package before use. Discard if the package is not sealed. Use of an unsealed device may lead to pathogen transfer, contamination of the patient, infection, or inflammation.
- Check the date of validity before use, and discard if the device is expired. Use of an expired device may lead to pathogen transfer, contamination of the patient, infection, or inflammation.
- Reprocessing, re-sterilization, or reuse of a disposable device may damage the device’s structural integrity, pose a risk of contamination / cross-infection to the patient, result in injury, illness or even death of patient.
- Carefully evaluate the tissue thickness. Have a range of staple sizes available and select the appropriate size reload for the tissue type and thickness. Overly thick or thin tissue may result in unacceptable staple formation, which may lead to staple line leak, bleeding, or prolonged surgery time.
- Do not attempt to remove the protective cap until the reload is properly loaded onto the stapler.
- Before loading or unloading, the jaws must always be open and parallel to the shaft (neutral position). Attempting to load/unload a reload with closed or articulated jaws may damage the load/unload mechanism, making multi-reloading no longer possible;
- Before introducing the stapler into the trocar or incision, always perform a close-open cycle of the jaws to confirm proper loading and clamping mechanism;
- Prior to introducing or removing the stapler from the trocar or incision, the jaws must always be in closed neutral position. Failure to set the jaws in closed neutral position will result in difficult insertion or withdrawal of the device;
- The anvil must be completely past the trocar (if a trocar is used) prior to opening the reload within the body cavity;
- Placement of tissue proximal to the tissue stops (on the reload) may result in stapler malfunction. Any tissue extending beyond the black “cut” line on the reload will not be transected, all tissue placed within the “cut” line will be transected;
- When positioning the stapler on the application site, ensure that no obstructions such as clips, stents, guide wires, etc. are within the jaws. Firing over an obstruction may result in incomplete cutting action, improperly formed staples, and/or inability to open the jaws, which could lead to staple line leak, bleeding, or prolonged surgery time.
- Ensure that the tissue lies flat and is evenly distributed inside the stapling range prior to closing the jaws. “Bunching” tissue may result in staple malformation, poor haemostasis and/or leakage;
- Always completely fire the reload. Partial firing leads to incomplete cut and/or staple malformation, which may result in poor haemostasis and/or leakage;
- The device is equipped with a safety interlock to prevent an empty or partially-fired reload from being fired a second time. Do not attempt to override the safety interlock, as it will cause device malfunction;
- If the clamping or firing mechanism fails to work (excessive resisting force), do not continue the operation, remove and carefully check the stapler-reload structural integrity following the IFU. The use of excessive force may result in device damage, injury, or tissue damage;

- After firing and removal of the device, always inspect the staple line and the surrounding site for air sealing / haemostasis of the anastomosis line and to see whether the staples are closed properly. Minor bleeding or leakage may be treated by electrocautery, manual suture or other appropriate methods;
- In a single procedure, a BL series stapler may be reloaded and fired up to 12 times with BLR reloads;
- Dispose of all opened devices whether used or unused. Devices which come into contact with bodily fluids require special disposal handling to prevent biological contamination, relevant local regulations shall be respected.
- The implantable staples are MR conditional, see the ***MRI Safety Information*** section for more information.

Instructions for Use

Prior to use

1. Carefully evaluate the length and thickness of the tissue to be excised, then select one or several reloads with the appropriate staple sizes from the table in “Reload Variants Specifications”. Keep in mind that 45mm staplers may only be used with 45mm reloads, while a 60mm stapler may only be used with 60mm reloads.
2. Check the integrity of the sterile barrier of the packaging, and the device expiration date. Do not use any device with damaged sterile barrier or device past expiration date.
3. Open the blister box of the stapler; ensure that the jaws are open.
4. Get familiar with the articulating action by pulling on the articulating fins and simultaneously applying a lateral pressure on the jaws with the other hand or against a supporting surface. The jaws pivot to a maximum angle of 45° in either direction, at which point the pivoting force increases to signal the end of the articulating arc. Maintain the lateral pressure while releasing the articulating fins to lock the jaws at one of the predetermined 15°/30°/45°angles, then release the lateral pressure. To return to straight (neutral) position, just pull on the articulating fins to release the lock, the angle of the jaws will automatically reset.
5. Open the blister box of reload, look for the presence of a protective cap. Discard the reload if the cap is not in place.

Loading the reload

6. With the jaws open, insert the reload by sliding it against the bottom of the jaw until the reload alignment tab stops in the alignment slot and snaps securely in place.
7. Remove and discard the protective cap. Examine the reload surface for the presence of all staples. Discard and replace the reload if any staples are missing.
8. Perform a close-open cycle to confirm proper loading. To close the jaws, squeeze the closing handle until it locks in place and an audible click is heard. To open the jaws, squeeze the closing handle, press and maintain the anvil release button, then slowly release the closing handle.

During use

9. Close the jaws, ensure that they are in neutral position, and insert the device into a 12mm trocar or incision. Avoid inadvertently pulling the firing handle. If the firing handle is unintentionally pulled, the device will be locked and it will no longer be possible to open the jaws by pressing the anvil release button. To recover from this state, remove the device and push down the red manual knife reverse switch. The knife direction indicator will display a backward arrow, indicating that the knife is in backward mode. Squeeze the firing handle completely, the knife position indicator will display zero, indicating that the knife is back to home (starting) position. Press the anvil release button to open the jaws and replace the reload.
10. Once inside the body cavity, ensure that the anvil is completely visible. Open the jaws. Use the shaft rotation collar to rotate the jaws 360° freely, and the articulating fins to articulate the jaws up to 45° in either direction. A body structure, organ, or another device may be used as a supporting surface to apply a lateral pressure on the jaws and to bend the jaws to one of the predefined 15°/30°/45° angles. When using body structures or organs as supporting surface, particular attention should be placed to the visual cues and tactile feedback received from the device. Do not apply frontal pressure or excessive lateral pressure to the tissue as tissue damage or tissue trauma may occur.
11. Position the device around the tissue to be stapled and ensure that the tissue lies flat and evenly inside the cutting range, with no obstructions within the jaws. Close the jaws by squeezing the closing handle completely. It is recommended to wait at least 15s before firing. If the closing handle is difficult to lock, reposition the device and take a smaller amount of tissue.
12. Squeeze the firing handle repeatedly to fully fire the reload: three complete strokes to cut/staple the full length of tissue within the cutting/stapling range, and 1 complete stroke to send the knife back to home position. For each stroke, squeeze fully and release the firing handle with continuous and smooth movement. The knife position indicator moves with each stroke. After the third stroke, the knife direction indicator will display a backward arrow,

indicating that the knife is in backward mode. At the end of the fourth stroke, the knife position indicator will display 0, indicating that the knife is back to home position.

13. Once the reload has been completely fired, it is recommended to keep the jaws closed and wait at least 15s. Open the jaws, gently remove the device from the tissue, close the jaws, pull on the articulating fins to reset the jaws to neutral position, and withdraw the device out of the body cavity through the trocar sleeve.

14. Carefully inspect the staple line for proper staple formation and haemostasis. Minor bleeding or leakage may be treated by electrocautery, manual suture or other appropriate methods.

Unloading the reload

15. Ensure that the jaws are open and in neutral position. Push the used reload up toward the anvil to unsnap the reload from the jaw, discard the used reload.

16. Additional reloads (up to 12 times in total) may be then loaded and fired in a single procedure. A safety interlock prevents an empty reload from being fired a second time. Do not attempt to override the safety interlock.

Clinical Benefits

The device offers the possibility of minimally invasive surgery and its associated lower amounts of intra-operative blood loss and smaller incisions. Smaller incisions result in less incisional infections, shorter recovery time, less post-operative pain, and smaller scars.

Undesirable Side Effects/Residual Risks

Undesirable side effects and risks associated with the device are the same as those related to stapling surgery and other suture methods: potential leakage, infection, bleeding, haematoma, stricture, abscess, ulcers, herniation, reoperation, conversion to open surgery, death, reflux, pneumonia, pain, allergy, pancreatitis, ileus, delayed gastric emptying, perforation, vomiting, effusion, fistula, empyema, recurrence, chylothorax, stroke, acute respiratory distress syndrome, conversion, readmission, gastrointestinal dysfunction. Users must communicate these risks to patients.

MRI Safety Information

The staples are MR conditional 

For whole-body MR examinations: a person implanted with the staples may be safely scanned anywhere in the body at 3.0T under the following conditions. Failure to follow these conditions may result in injury.

- Static magnetic field (B0) no more than 3.0T;
- Maximum allowable spatial gradient magnetic field (▽B) is 30T/m (3,000 gauss/cm);
- Magnetic induced torque is less than worst-case gravitational torque;
- In 3.0T or 1.5T MR system, whole-body average Specific Absorption Rate (SAR) no more than 2W/kg, RF exposure no more than 15mins;
- Artifact information, sample causes artifacts up to 3.8mm in 3.0T MR system and 3.2mm in 1.5T MR system.

Packaging

The device is packaged in a sterile blister box and paper box.

Duration of Validity

The device is valid for 5 years.

Production Date

Refer to package label.

Storage and Handling Conditions

Store at room temperature, avoid prolonged exposure to elevated temperatures. Do not expose to temperatures above 30°C or below 10ºC, do not expose to humidity over 70% or below 35%. . Handle with care, do not over-stack or over-palletize.

Sterilization Method

Ethylene oxide sterilization



Incident Reporting

Any serious incidents related to the device should be reported to the manufacturer and the competent authority of the country in which the user and/or patient is established.

Summary of Safety and Clinical Performance (SSCP)

The SSCP is available in the European database on medical devices (Eudamed). Use the Basic UDI-DI 69492989LS0001W5 to find the device on the website: https://ec.europa.eu/tools/eudamed.

Symbol Description

	Date of manufacture		Manufacturer
	Catalogue number		Authorized representative
	Batch code		Sterilized using ethylene oxide
	Unique device identifier		Single sterile barrier system with protective packaging outside
	Medical device		Single sterile barrier system
	Use-by date		Do not re-sterilize
	Caution		Do not re-use
	Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Upper limit of temperature		Keep away from sunlight
	Fragile, handle with care		Keep dry

	Humidity limitation		MR conditional

Jiangsu Brightness Medical Devices Co., Ltd.

Building 1, the 3rd floor of Building 3, the 2nd, 3rd, 5th floor of Building 5 & Room 7-3 of Building 7, No.66, Hehuan Road, Zhonglou Economic Development Area, Changzhou, 213000 Jiangsu, P.R. China

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