Jiangsu Brightness Medical Devices Co., Ltd.

SKIN STAPLERS AND REMOVERS

CE 0197

Instructions for Use

PLEASE READ THE FOLLOWING INSTRUCTIONS CAREFULLY BEFORE USE. FAILURE TO PROPERLY FOLLOW THE INSTRUCTIONS MAY LEAD TO SERIOUS SURGICAL CONSEQUENCES, SUCH AS STAPLE MALFORMATION, LEAKAGE, OR DISRUPTION OF OPERATION.

This document aims to assist surgeons in using this product; it is not a reference to surgical procedure.

This device is for single-use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reprocessing and/or re-sterilization of this device may create the risk of contamination and patient infection. Do not reuse, reprocess or re-sterilize this device.

Intended purpose

This device is intended for use for skin suturing of wounds and surgical incisions.

Description

This is a sterile, disposable, single-patient use device. It places one single staple on the applied skin surface, completing the skin suturing. Staple size and number is determined by the selection of the stapler type. The device may be used on superficial skin.

PF(X) series staplers may fire up to 55 times in a single procedure according to the specification chosen. They offer staples of uniform width and height, with options of width 9.8 mm and leg length 2.7mm, or width 12.9mm and leg length 3.4mm.

Intended Users and Target Patient Group

The device is intended to be used by surgeons who have adequate training and familiarity with skin techniques. The target patient population consists of adults requiring skin suturing of wounds and surgical incisions.

Indications


The device is indicated for skin suturing of wounds and surgical incisions.

Contraindications

- Do not use the device where surgical stapling is not appropriate;
- Do not use the device when it is not possible to maintain at least a 5 mm distance from the stapled skin to underlying bones, vessels, or internal organs;
- If the tissue's tension is too high, make an approximating-suture before closing skin with skin stapler.

Stapler Variants Specifications

Staplers models and variants are coded according to the following rules:

For example, PFW35 means PF stapler with 35 W size staples.

Stapler specifications basic dimensions (unit: mm)

Type	Specification	Size	Total length (±5)	Staples					Remover
				Diameter (±0.1)	Before formation		After formation	Staple quantity (0~+3)	
					Width(±0.2)	Height(±0.1)	Width x Height (±0.5)		
PF	15, 25, 35, 45, 55	R	140	0.6	9.8	2.7	5.5 x 3.5	15, 25, 35, 45, 55	PHC-R、PHU-R
		W		0.6	12.9	3.4	7.0 x 4.5		PHC-W、PHU-W
PFX	15, 25, 35, 45, 55	R	168	0.6	12	2.9	6.0 x4.0	15, 25, 35, 45, 55	PHC-R、PHU-R
		W		0.6	13.8	3.6	7.0 x 4.5		PHC-W、PHU-W

【Main technical specification】

Material: The activator consists of stainless steel 06Crl9Ni10 according to the GB/T1220-2007, the staples consist of stainless steel 022Cr17Ni12Mo2 according to GB/T 4240-2019. The handle consists of ABS according GB/T12672-2009. The staple remover consists of ABS and stainless steel 06Cr19Ni10.

Performance characteristics: Non-toxic, harmless and sterile, non-caustic, good biocompatibility, good metal stability, quick accomplishment of sutures, good site of wound edges, applicable for different kinds of wounds and cuts.

Easy removal of staples by using staple remover.

【Application】

used to suture wound and surgical cuts.

【Instructions for use】

Instruction of Skin Stapler

1.Hold skin stapler steady, keep the staples vertical to the skin, center the activator lightly over the wound edge, press skin stapler on skin.

2.Compress active handle completely until the staple formation is completed. (Do not release the handle before the staple has been formed.)

3.Release active handle until the handle is back in home position Place skin stapler to the next position of closure.
4.Repeat actions as mentioned above.

Instruction of Staple Remover

1.Place remover's two hooks in-between staples and tissue
2.Close the staple remover and remove staples from tissue.

【Warning and Precautions】

1.If the tissue's tension is too high, make an approximating-suture before closing skin with skin stapler.
2.If the distance between the stapled skin and underlying vital structure is less than 5 mm, the skin stapler should not be used!
3.The skin stapler is for skin closure only.
4.Please check the package's sealing before use. Do not use the instrument if the package is damaged.
5.Use skin stapler in valid time of sterilization.
6.For single use only!
7.Non-resterillzable product.

【Package】

Skin stapler is packed in a blister box and sterilized with EO. One piece per package. The outer package is a paper box and labeled with

■ Sterilized using ethylene oxide

■ Do not re-use

There is also a "forbidden to use if package is damaged" warning on package.
Valid time of sterilization is five years.
The skin stapler remover is disposable. It is packed sterile.

【Storage】

The depository should be dry, cool and well ventilated. Prevent instrument from damage.

Jiangsu Brightness Medical Devices Co., Ltd.
ADD:Building 1, the 3rd floor of Building 3, the 2nd, 3rd, 5th floor of Building 5 & Room 7-3 of Building 7, No.66, Hehuan Road, Zhonglou Economic Development Area, Changzhou, 213000 Jiangsu, P.R. China
Tel: +86-519-83908835 Fax: +86-519-83908831

EU REP **Lepu Medical (Europe) Cooperatief U.A.**
Add: Abe Lenstra Boulevard 36, 8448 JB, Heerenveen, The Netherlands
Tel: +31-515 573399 Fax: +31-515 760020
Date of issue: 06.05.2025 Revision no: A/2

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