Regenera Activa AMT


hbw srl
human brain wave

EUROPEAN PATENTED RIGENERA TECHNOLOGY

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RIGENERACONS (PRODUCT CODE: 7945AS)

INTENDED PURPOSE

Rigeneracons device is intended for the mechanical disaggregation and filtering of small fragments of soft tissues such as dermis and hypodermis, in order to obtain a suspension of autologous micrografts.

USAGE NOTE

The device should be used for the disaggregation of n. 2-3 punches of soft tissues such as dermis, hypodermis collected using a 2.5 mm biopsy punch. The device must be used in combination with the Rigenera machine N4SA EVO that permit its rotation activating the tissue disaggregation.

INDICATIONS

- Rigeneracons device must be used in a medical environment by professionals informed about user instructions.
- Rigeneracons device is used to generate autologous micrografts that in turn are able to promote tissue repair.

CONTRAINDICATIONS

- The devices are sterile and disposable. Do not reuse to avoid the risk of infections.
- No other contraindications are reported.

STORAGE AND USE OF ADIPECONS DEVICE

- Storage temperature: Keep in a dry place and away from heat sources.
- Transportation temperature: Keep in a dry place and away from heat sources.
- Handling temperature: Keep in a dry place and away from heat sources.
- How to maintaint the product: keep the device in the package until to use.
- How to dispose: according to legislation in force.

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