IntraNova® Intra-Aortic Balloon Pump

User Manual

Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.

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Manufacturer Information

Registrant / Manufacturer

Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.

Registrant/Production Location

1401, Building D1, Nanshan i Park
No.1001, Xueyuan Avenue, Taoyuan Street
Nanshan District
Shenzhen, Guangdong, China

Tel: +86 755 23251894
WWW.LTCARDIO.COM

Device information and Support

Product Name: Intra-Aortic Balloon Pump
Mode: IAP-50
After-sales Service Provider:
Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.
After-Sales Service Hotline: 400-803-5503

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QUICK START GUIDE ---------- 4

1. DEVICE OVERVIEW ---------- 5

1.1. General Description ---------- 5
1.2. Product Configuration List ---------- 5
1.3. Product Functions ---------- 5
1.4. Safety Features ---------- 5
1.5. Intended Use ---------- 6
1.6. Indications ---------- 6
1.7. Contraindications ---------- 6
1.8. Operator Requirements ---------- 6
1.9. Compatible Equipment ---------- 7
1.10. Transportation and Storage Conditions ---------- 7
1.11. Operating Environment ---------- 7
1.12. Terms, Abbreviations, and Definitions ---------- 8

2. WARNINGS, CAUTIONS AND SAFETY INFORMATION ---------- 10

2.1. WARNINGS ---------- 10
2.1.1. Routine Use and Maintenance ---------- 10
2.1.2. Clinical Use ---------- 12
2.2. PRECAUTIONS ---------- 17
2.3. Safety Information ---------- 21
2.3.1. Measures to Prevent User Misoperation ---------- 21
2.3.2. Abnormal conditions ---------- 21

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5. INSTALLATION INFORMATION ---------- 131

6. USER MAINTENANCE ---------- 132

6.1. Replacing the Helium Cylinder ---------- 132
6.2. Refilling the Helium Cylinder ---------- 133
6.3. Battery Maintenance ---------- 133
6.4. Backup Batteries and Filter Maintenance ---------- 134
6.5. Removing the Main Unit from the Cart ---------- 135
6.6. Installing the Main Unit on the Cart ---------- 135
6.7. Fuse Maintenance ---------- 135
6.8. User Cleaning and Maintenance ---------- 135
6.8.1. Cleaning the System and Accessories ---------- 135
6.8.2. Cleaning the Display Screen ---------- 135
6.8.3. Cleaning and Disinfection of Cables ---------- 136
6.9. Preventive Maintenance Plan ---------- 136
6.9.1. Preventive Maintenance Plan A ---------- 136
6.9.2. Preventive Maintenance Plan B ---------- 137
6.10. Software Maintenance ---------- 138
6.11. Warranty ---------- 138

7. EXTERNAL MONITORING DEVICE INTERFACES ---------- 140

7.1. ECG Requirements from External Monitoring Devices ---------- 140
7.2. Arterial Pressure Requirements from External Monitoring Devices ---------- 140
7.3. ECG Input Description from Monitoring Devices ---------- 140
7.4. Arterial Pressure Input Description from Monitoring Devices ---------- 140
7.5. ECG Output Description to Monitoring Devices ---------- 140
7.6. Description of Low-Level Pressure Output to Monitoring Devices ---------- 141

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8.12. Invasive Blood Pressure Sensor Cable ---------- 147
8.13. Display Cable ---------- 147
8.14. Power Cord ---------- 147
8.15. Electromagnetic Compatibility ---------- 148
8.16. Equipotential Bonding Terminal ---------- 151
8.17. Specifications of the Fuse Device ---------- 152
8.18. Main Safety Characteristics of the Device ---------- 152
8.19. Software Release Version ---------- 152
8.20. Standards Compliance ---------- 152
8.21. Essential performance ---------- 153

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Warning

MR Unsafe

Indicates approximate lithium battery charge level

Indicates approximate remaining helium level

Alarm paused

High-priority alarm

Medium-priority alarm

Low-priority alarm


Indicates a detected physiological trigger event

Indicates an internally generated trigger event by the IABP

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

Humidity limitation

Equipotentiality

Do not dispose of this product as unsorted municipal waste. Dispose of in accordance with local regulations

Defibrillation-proof type CF applied part


Product information

Date of manufacture

AA/LR6 battery

Serial number

Product number

Batch code

PN

Material code

Patient monitoring equipment

Optical fiber connector

Orientation for optical fiber connector connection


USB port

Non-ionizing electromagnetic radiation

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

WARNING

Do not operate this device until you have thoroughly read this manual and all related warnings and precautions. If you are already familiar with this product, the Quick Start Guide may be used. If you are not familiar with this product, continue reading the remainder of this manual.

Quick Start Guide

1) Turn on the IABP power inlet switch. After the status display illuminates, press the IABP power button "⏻".
2) Ensure that the Helium Cylinder is open and confirm that the helium level is sufficient.
3) Establish the ECG and pressure signal connections.
4) If a direct pressure source from a sensor is used, zero the sensor.
5) Vent the Disposable Blood Pressure Sensor or the fiber-optic sensor to atmosphere.
6) Press and hold the "Zero" button.
7) Close the sensor to atmospheric pressure.
8) Confirm that the operating mode is set to Automatic.
9) Connect the IAB catheter and the helium extension tubing to the pneumatic module.
10) Press the "Start" button to initiate pumping. In response, the IABP will automatically inflate the balloon and initiate counterpulsation. If necessary, use the "Deflation" control to fine-tune the balloon deflation timing.
11) Verification of Augmentation Pressure Alarm Settings：Approximately 3 minutes after augmentation assistance is initiated, confirm that the augmentation pressure alarm is set to approximately 10 mmHg below the patient's augmentation pressure. If necessary, adjust the alarm value using the Up "▲" or Down "▼" buttons in the "Augmentation Pressure Alarm" threshold settings menu.

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1. Device Overview

1.1. General Description

Lifetech Cardio IntraNova® Intra-Aortic Balloon Pump (hereinafter referred to as the IABP) consists of a Main Console, a Display, a Cart, ECG Lead Wires, and an Invasive Blood Pressure Sensor Cable.

IntraNova® IABP is a cardiac assist device used in conjunction with Intra-Aortic Balloon Catheters. The Intra-Aortic Balloon Catheter (hereinafter referred to as the IAB catheter) is inserted into the arterial system. The balloon inflates during cardiac diastole and deflates during cardiac systole, thereby increasing coronary artery perfusion and reducing left ventricular workload to provide support to the left ventricle.

1.2. Product Configuration List

1) Main Console ×1
2) Display×1
3) Cart×1
4) ECG Lead Wires×1
5) Invasive Blood Pressure Sensor Cable×1
6) Instructions For Use×1
7) Others: Display Cable ×1, Power Cord ×1, Helium Cylinder ×1, and Luer Caps×5

The fully assembled IABP, including the Main Console, Display, and Cart, is shown in Figure 1.


Figure 1

1.3. Product Functions

The functional features of the Lifetech Cardio IntraNova® IABP are summarized as follows:

- Automatic operating mode
- Semi-automatic operating mode
- Augmentation function
- Report output

1.4. Safety Features

Prolonged inflation time protection: During augmentation assistance, if the IABP detects that the inflation time exceeds 1500 ms, the inflation protection procedure is activated, triggering is suspended, and the balloon is automatically deflated.

Power-loss deflation: when the IABP loses power, the balloon automatically deflates.

R-wave deflation: if an early R wave is detected while the balloon is inflated, the IABP immediately initiates deflation.

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Screen lock: in any mode, pressing the screen lock button locks the touchscreen. After the IABP is powered on and a augmentation cycle has been initiated, the touchscreen automatically locks after 2 minutes of inactivity by default.

1.5. Intended Use

This product is intended for use in conjunction with specific models of Intra-Aortic Balloon Catheters manufactured by the company and is clinically indicated for intra-aortic balloon counterpulsation therapy in adult patients.

This product is intended to be operated by qualified physicians in medical institutions.

1.6. Indications

Specific indications for the Intra-Aortic Balloon Pump include, but are not limited to, the following conditions：

1) Acute coronary syndrome
2) Cardiac and non-cardiac surgery
3) Complications of heart failure

1.7. Contraindications

Absolute contraindications:

1) Hemodynamically significant or severe aortic valve insufficiency
2) Aortic aneurysm or aortic wall disease (including dissection)

Relative contraindications:

1) Atherosclerosis or severe peripheral vascular disease (PVD)
2) End-stage disease or sepsis
3) Severe coagulopathy

WARNING

The Lifetech Cardio IABP is not a substitute for patient monitoring equipment.

1.8. Operator Requirements

The primary intended users of the Intra-Aortic Balloon Pump are physicians and nurses in cardiology, critical care medicine, emergency medicine, and anesthesiology who have received training in the theory, techniques, and clinical aspects of counterpulsation therapy. Users shall not operate the IABP before reading the operating instructions, including warnings and precautions. Users shall read and understand the language descriptions of the system control methods and features. Users must have normal vision and be able to clearly observe the waveforms displayed on

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the device screen.

Lifetech Cardio recommends that healthcare personnel receive appropriate training before operating the IABP. Training may be provided onsite by local sales and/or clinical representatives. For additional information, education, and training, please visit the Lifetech Cardio website.

Users may operate the IABP in a seated or standing position. The display angle may also be adjusted to accommodate different operator positions for viewing and operating the device.

1.9. Compatible Equipment

Lifetech Cardio Intra-Aortic Balloon Catheter:

Model		Balloon Membrane Volume
Standard Series	Fiber-Optic Series	
IAB-01-30	IAB-02-30	30cc
IAB-01-40	IAB-02-40	40cc
IAB-01-50	IAB-02-50	50cc

Disposable Blood Pressure Sensor: use sensors that comply with YY0781-2010 and its recommended performance parameters, and that meet the following requirements: (1) sensor excitation of 5 V DC; (2) sensitivity of 5 μV/V/mmHg; and (3) an RJ11 connector or compatibility with the sensor interface of the Lifetech Cardio IAP-50 invasive blood pressure sensor cable via a standard adapter.

Disposable ECG Electrodes: use silver/silver chloride disposable ECG electrodes that are approved for market use through registration.

Invasive Blood Pressure Monitor

ECG Monitor

1.10. Transportation and Storage Conditions

Temperature: −20°C to 60°C;

Humidity: 5% to 95% RH;

Altitude: ≤4000 m.

1.11. Operating Environment

Operating temperature: 5°C to 40°C;

Operating humidity: 5% to 95% RH;

Operating altitude: ≤4000 m.

This device is intended for extracorporeal use and does not come into direct contact

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with human tissue. Applicable medical environments include operating rooms, emergency rooms, and hospital wards.

1.12. Terms, Abbreviations, and Definitions

—The abbreviation “IABP” shall be considered equivalent to the Lifetech Cardio “Intra-Aortic Balloon Pump.”

—The abbreviation “IAB catheter” shall be considered equivalent to the Lifetech Cardio “Intra-Aortic Balloon Catheter.”

—The term “medical device” shall be considered equivalent to ME equipment.

—The term “user” shall be equivalent to the operator.

—The term “safety” shall be equivalent to basic safety and essential performance.

—The term “user interface” shall be equivalent to the operator–device interface.

—The term “A pacemaker trigger” refers to atrial pacemaker triggering.

—The term “V/A-V pacemaker trigger” refers to ventricular/atrial–ventricular pacemaker triggering.

—The term “pressure trigger” refers to blood pressure triggering.

—The term “reference line” refers to the pressure value corresponding to the current reference line position displayed.

—The term “pneumatic test” refers to verification of the integrity of the pneumatic module.

—The term “augmentation ratio” refers to the ratio of counterpulsation cycles to heart rate cycles, such as 1:1 or 1:2.

—The term “trigger source” refers to a signal capable of triggering counterpulsation, such as ECG signals or blood pressure signals.

—The term “inflation start time” refers to the time interval from the trigger event to the start of balloon inflation during counterpulsation therapy.

—The term “deflation start time” refers to the time interval from the trigger event to the start of balloon deflation during counterpulsation therapy.

—The term “signal source” refers to the origin of the trigger signal. For the ECG channel, signal sources include leads I, II, III, aVR, aVL, aVF, V, and external. For the

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blood pressure channel, signal sources include fiber-optic, fluid-filled (hydraulic), and external sources.

—The term “inflation timing step” refers to the timing step between the trigger event and the start of inflation.

—The term “deflation timing step” refers to the timing step between the trigger event and the start of deflation.

—The abbreviation “AP” stands for Arterial Pressure.

—The abbreviation “BPW” stands for Balloon Pressure Waveform.

—The abbreviation “Aug” stands for Augmentation.

—The abbreviation “Device Info” stands for Device Information

—The abbreviation “Ext” stands for External

—The abbreviation “Semi” stands for Semi-automatic

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2. Warnings, Cautions and Safety Information

2.1. WARNINGS

2.1.1. Routine Use and Maintenance

Device operation: Do not operate this device until you have read this manual and all related warnings and precautions. If you are already familiar with this product, refer to the Quick Operation Guide. If you are not familiar with this product, continue reading the remainder of this manual.

Device operation: Do not operate this device until you have read this manual and all related warnings and cautions. If you are already familiar with this product, refer to the Quick Operation Guide. If you are not familiar with this product, continue reading the remainder of this manual.

Packaging and transport: Compressed gas (Helium Cylinders) and lithium batteries are classified as dangerous goods. Packaging must be identified, packed, labeled, classified, marked, and documented in accordance with applicable regulations, and appropriate primary or secondary packaging must be used for dangerous goods transportation.

Risk of electric shock: There are no user-serviceable parts inside this device (except fuses). Do not open the device enclosure. Maintenance must be performed only by service personnel trained and certified by Lifetech Cardio.

Fuses: Always disconnect the external power supply before replacing fuses.

Use of non-approved components: Use of components or cables not supplied or approved by Lifetech Cardio may increase electromagnetic emissions, reduce electromagnetic immunity, or cause malfunction of the IABP. It may also delay recovery after defibrillator discharge. Use only power cords with 3C certification or supplied by Lifetech Cardio.

Pneumatic module: Wear surgical gloves when replacing the pneumatic module to avoid contact with residual condensate or other bodily fluids. Dispose of used safety disks in accordance with hospital medical waste regulations.

Device servicing: Do not modify the device. Except for fuses, there are no user-serviceable parts. Do not disassemble or repair the IABP yourself. Any disassembly or repair not performed by Lifetech Cardio may result in malfunction or compromise safety.

Liquid ingress protection: The operator may wipe the device with a cloth dampened

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with clean water. Do not immerse the device in any liquid. Liquid ingress may impair device function and safety.

Fire and explosion hazard: This device is not intended for use in environments containing flammable or explosive gases, such as anesthetic gases, oxygen, or nitrous oxide, nor near flammable or explosive materials.

Monitoring: The Lifetech Cardio IABP does not replace patient monitoring equipment.

Not MRI compatible: Magnetic resonance imaging (MRI) is a medical imaging technique that uses magnetic fields to generate views of internal body structures. Do not perform MRI scans on patients who are being treated with the Lifetech Cardio IABP. Do not use the IABP near an MRI scanner or during MRI scanning. MRI scanning may induce tachycardia and may result in serious injury, or cause system malfunction or damage.

Installation and placement: Do not place the Lifetech Cardio IABP adjacent to or stacked with other equipment. If adjacent or stacked use is unavoidable, observe and verify that the device operates normally in the intended configuration.

When placing the main unit on the cart, ensure that the locking latch is engaged to secure the main unit.

Do not place the device in a location where the disconnection device is difficult to access or operate.

Device location: When positioning the IABP, ensure that the AC plug remains accessible. If an AC power issue occurs, disconnect the IABP from the AC power input by unplugging it from the wall outlet.

Grounding: To reduce the risk of electric shock, this device must be connected to a power supply with protective grounding.

Power disconnection: Do not place the device in a location where the power cord cannot be easily unplugged, to ensure that power can be disconnected promptly in an emergency.

Alarms: Do not perform alarm testing while the IABP is being used for patient treatment.

Do not set the alarm volume too low, otherwise alarms may not be audible above the background noise in the intended use environment of the IABP.

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Waste disposal: Waste materials, residues, and the ME device and its accessories must be disposed of properly in accordance with local regulations at the end of their intended service life.

2.1.2. Clinical Use

Physiological signals: Operating the IABP below the minimum amplitude or value of the patient's physiological signals may result in incorrect outcomes.

Unexpected shutdown: In the event of a sudden unexpected shutdown, disconnect the IAB catheter from the IABP and ensure that the IAB catheter is deflated.

Leak test: Do not perform a pneumatic module leak test while the patient's IAB catheter is connected to the IABP.

Internal trigger: Do not use internal triggering when the IABP is being used for patient treatment.

Pressure triggering: If necessary, after each manual adjustment of the pressure trigger threshold, reassess and adjust the inflation and deflation timing selections.

If the manual pressure trigger threshold option is used, adjust the threshold when the systolic pulse amplitude changes continuously. Changes in patient condition or initial zeroing of a Fiber-Optic IAB catheter may cause such amplitude changes. Regardless of the reason for adjusting the trigger threshold, inflation and deflation timing must be reassessed.

Augmentation pressure alarm: The augmentation pressure alarm automatically set at startup supports IAB alarms (leak and IAB catheter alarms) at higher heart rates. Therefore, do not manually disable this alarm.

When weaning by reducing augmentation pressure, the reduction must not lower the IAB status indicator below 50%.

Preventive maintenance: Do not perform preventive maintenance while the IABP is connected to a patient.

Screen cleaning: Do not clean the screen while the IABP is connected to a patient.

Components: To ensure proper system performance and defibrillator protection, use only cables and components approved by Lifetech Cardio, including lead wires, IAB catheters, system software, and pneumatic modules. Use only Lifetech Cardio-supplied ECG Lead Wires and Invasive Blood Pressure Sensor Cables.

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Balloon leakage: If balloon membrane leakage is suspected, do not inflate the IAB catheter using a syringe or any other method.

Continued assistance with a leaking IAB catheter may result in formation of large blood clots within the balloon, potentially causing difficulty in balloon removal and requiring surgical intervention.

Balloon perforation: Balloon perforation may indicate vascular conditions that cause balloon abrasion or perforation, which may also damage a replacement balloon.

System standby: To avoid thrombus formation, avoid leaving the balloon inactive (neither inflating nor deflating) inside the patient for more than 30 minutes.

Monitoring: Do not leave the patient unattended during IABP therapy.

Interference: External bedside monitors providing ECG signals to the IABP must have electrosurgical interference suppression capability.

Placement environment: Under no circumstances should patients with an IAB catheter or the IABP be placed in a hyperbaric chamber.

Patient heart rate: Maintain the patient's heart rate below 210 bpm to avoid system mis-triggering due to ECG waveform changes.

Displayed heart rate values may fluctuate significantly during arrhythmias.

Cable routing: Route cables neatly and keep them away from the patient's neck to prevent strangulation. Do not place cables on floors or walkways to reduce tripping hazards.

Hose routing and the placement of the helium inflator must be safe. Keep walkways clear to reduce the risk of injury.

Power cords must be routed safely.

Heart rate and timing selection: At certain heart rates and timing selections, catheter alarms may be interrupted. For further details, refer to this manual.

When pressure is used as the trigger source, balloon deflation must be continuously adjusted to ensure deflation is completed during cardiac systole. Late deflation may reduce systolic pulse pressure and delay detection of the pressure signal. To achieve consistent and reliable pressure triggering, the system relies on clear and timely systolic pulses. If balloon deflation conflicts with systolic ejection when pressure is used as the trigger source, triggering may be delayed or missed, potentially resulting in loss of synchronization.

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External monitoring equipment: Output signals from external monitoring equipment must meet specified requirements. Refer to the External Monitoring Equipment Connection section (Chapter 7).

Battery: The internal lithium battery and backup battery may only be replaced or removed by service personnel trained and certified by Lifetech Cardio. Improper replacement may result in device malfunction.

Use only Lifetech Cardio lithium batteries.

Lithium batteries and backup batteries pose risks of fire, explosion, and severe burns. Do not disassemble, crush, heat above 60°C (140°F), or incinerate.

Insulation: The ECG lead wires and invasive blood pressure channels used by this device employ standard insulation. Ensure that the conductive parts of ECG electrodes do not contact other conductive parts, including ground. When connected to a patient, do not connect any non-insulated accessories to the IABP or to the ECG or invasive blood pressure inputs. Ensure that the total enclosure leakage current of all connected units does not exceed 300 μA. If any covers are removed from the device, do not simultaneously touch the patient and any part of the electrical equipment.

The power cord and interface cables (i.e., non-patient cables) may share the same ground. Therefore, unplugging the power cord may not fully isolate the IABP if non-patient interface cables are connected.

Defibrillator use: Exercise extreme caution when a defibrillator is used on a patient. During defibrillation, do not touch the patient, the table, or any part of the IABP.

Risk of electric shock: When connecting the IABP to other medical or non-medical electrical equipment, securely connect the equipotential bonding connector to a safety ground to minimize the risk of excessive leakage current and electric shock.

Do not place sockets used with the IABP or its cables on the floor. Connect IABP cables only to sockets designated for IABP use. Do not cover the sockets. Do not connect other equipment to sockets used for the IABP, as system leakage current may increase.

The operator must not touch the device and the patient simultaneously to prevent excessive leakage current.

Patients with pacemakers: In patients with pacemakers, pacing pulses may be counted by the heart rate meter during asystole or arrhythmia. Do not rely solely on heart rate alarms. Closely monitor patients with pacemakers. For information on pacing pulse suppression, refer to Section 8.1 of this manual.

Trigger source selection: Pressure triggering is not recommended in the presence

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of persistent irregular heart rates or tachyarrhythmias. When pressure triggering is used, promptly and adequately adjust deflation to ensure completion before cardiac systole, and continuously monitor the patient.

If possible, use ECG or pressure triggering during CPR. This helps synchronize counterpulsation with the rate and rhythm of chest compressions. In automatic mode, ECG (R-wave) or pressure signals are automatically selected as the trigger source based on relative signal quality. If neither ECG nor pressure signals provide reliable triggering, the device's internal trigger frequency may be used by selecting semi-automatic mode and setting the trigger source to internal.

Do not remain in internal trigger mode when the patient has cardiac ejection.

Helium Cylinder replacement: Helium Cylinders may only be replaced by personnel familiar with high-pressure gas cylinder handling.

Alarm sound: When multiple IABP systems are in use nearby, confirm the corresponding visual indicators to correctly identify the source of alarm sounds.

In high-noise transport environments, users should rely primarily on visual alarm messages.

Instruments: Do not reuse single-use devices.

Electromagnetic interference: Electromagnetic fields may affect device performance; therefore, nearby equipment must comply with EMC requirements. Mobile phones, X-ray, and MRI equipment are potential sources of interference.

Use of the IABP near high-frequency electrosurgical, defibrillation, or shortwave therapy equipment may result in electromagnetic interference that affects normal operation and poses a risk to the patient.

Except for transducers and cables sold by the manufacturer as internal component spare parts, use of unauthorized accessories, sensors, or cables may increase emissions or reduce immunity of the device or system.

Alarm sound pressure level adjustment: Adjust the alarm sound pressure level with caution. Alarm sound levels below ambient noise may impair recognition of alarm conditions.

Augmentation pressure alarm: Setting the augmentation pressure alarm limit too low or disabling it will render the alarm ineffective.

IAB catheter: The Lifetech Cardio IABP may only be used with IAB catheter models specified by Lifetech Cardio.

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Blood pressure sensors: Disposable Blood Pressure Sensors used with the IABP must comply with YY0781-2010 and its recommended specifications. Refer to Section 1.9 for detailed requirements.

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2.2. PRECAUTIONS

Power on: When repeatedly switching the device on and off, wait at least 10 seconds before restarting.

After completion of the power-on self-test, the device enters power-saving mode if idle for more than 5 minutes.

Installation: Before use, mount the Main Console and Display on the Cart and ensure that the Cart locks secure both components.

Display cable: The cable connecting the Main Console and the Display is not detachable.

Removal: When removing the Main Console and Display while not in use, place them on a level surface to prevent damage due to movement or falling. Handle the Main Console and Display with care to avoid damage caused by impact, pushing, pulling, sitting, or stepping on them.

Electrostatic discharge: To eliminate the effects of static electricity on the device and the patient, the operator should discharge any static electricity before use, for example by touching a large conductive grounded surface.

Mains-powered equipment: If equipment connected to mains power (220 V) is used, ensure that it is properly grounded.

RF communication equipment: Keep handheld RF communication equipment (including antenna cables and external antennas) as far as possible from the cables of the Lifetech Cardio IABP, typically no less than 30 cm, otherwise device malfunction or damage may occur.

RF immunity: As the Lifetech Cardio IABP is only immune to specified RF signals, operation near RF transmitters may cause abnormal device operation.

Inspection: Operators shall routinely check the IAB catheter helium extension tubing for condensation. Excessive condensation may affect system performance and indicates that pneumatic module maintenance is required.

Do not operate the device if ventilation openings or speaker vents are blocked.

Maintenance: To maintain specified performance, this product requires regular preventive maintenance. Maintenance includes periodic cleaning of electronic components using cooling airflow.

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Lithium battery: When powered on, the lithium battery charges only when the AC plug is connected to a live AC outlet and the battery charging icon is displayed. When powered off, charging occurs only when the power button flashes green. The battery may remain charging even when the system is not in use.

If long-term storage is required, or if the ambient temperature exceeds the operating range, refer to Battery Maintenance in the User Maintenance section (Section 6.3).

After storage at low temperatures, allow the Lifetech Cardio IABP to remain at room temperature for at least 30 minutes before operating on lithium battery power.

Before shipping the Lifetech Cardio IABP via commercial transport, the lithium battery must be disconnected.

Lithium batteries of the Lifetech Cardio IABP must be transported by carriers qualified in dangerous goods transportation. Packages must be properly packed and labeled.

The lithium battery and backup battery must be properly maintained and inspected regularly. Refer to Battery Maintenance in User Maintenance (Section 6.3).

If a mains power interruption occurs while sufficient internal battery capacity is available, the device automatically switches to internal battery power and continues normal operation without changes to operating mode, user settings, or stored data.

When operating on internal battery power, monitor the battery charge level to prevent insufficient power.

Helium: Use medical-grade helium only.

The Helium Cylinder should be closed when not in use for an extended period.

If the Helium Cylinder is correctly installed and no leakage is present, the cylinder valve may remain open when the system is not in use.

When refilling helium for the IABP, follow hospital biohazard procedures to avoid contact with pathogens.

The internal Helium Cylinder contains helium at pressures up to 15 MPa. Before commercial transport, the cylinder must be emptied and removed from the device.

Information display area: If more than two alarms or messages are present, additional items will be collapsed. Use the Up "▲" or Down "▼" buttons to view them.

Diagnostic function: The device has no intended diagnostic function.

Power supply: Prior to emergency use, if powered by an AC inverter, qualified maintenance personnel must verify proper compatibility. When operating correctly, the icon and indicator area should not display "Operating on lithium battery."

When AC power operation is required, ensure the plug is connected to a live AC outlet and that "Operating on lithium battery" is not displayed.

When positioning the IABP, ensure the AC plug remains accessible. If an AC power

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issue occurs, unplug the device from the wall outlet to disconnect AC power.

Leads: Conductive parts of applied electrodes and associated connectors (including neutral electrodes) must not contact other conductive parts, including ground.

Do not touch the exposed end of the fiber-optic balloon catheter cable or allow it to contact other surfaces, as this may damage or contaminate the sensor interface.

Pull handle: The pull handle is intended only for moving the device; do not use it to lift the device.

Timing selection: Do not set timing based on highlighted inflation interval markers. Set timing by evaluating the arterial pressure waveform.

Display: The displayed ECG signal is automatically scaled to provide optimal screen presentation. Due to automatic scaling, low-amplitude ECG signals may appear with normal amplitude on the display. The ECG amplitude of interest shall be interpreted according to the indicated scale.

The displayed AP signal is automatically scaled and offset to provide optimal screen presentation. Due to automatic scaling, low-amplitude AP signals may appear with normal amplitude on the display. Arterial pressure parameters shall be interpreted according to the indicated scale.

If abnormal display behavior occurs, check whether the cable connections between the main unit and the display are loose or disconnected, and confirm that all cables are properly connected.

IAB catheter: To ensure reliable operation of the automatic inflation system and correct IAB catheter inflation pressure, it is critical not to alter the total internal volume of the balloon and extracorporeal tubing together with the helium extension line. Use of tubing with a different length or inner diameter from that specified for the Lifetech Cardio intra-aortic balloon catheter and accessories will alter the inflation pressure and may result in automatic inflation failure. Such conditions shall be avoided.

When connecting the IAB catheter to the device, ensure that the Luer connector is fully seated and securely tightened. During movement, take care to maintain secure connection between the IAB catheter and the device.

Transport: When removing the pump main unit from the trolley or returning it to the trolley, ensure that the trolley wheels are locked.

The IV pole shall not be used as a handle or lifting point. During transport or movement of the system, only the trolley handle shall be used.

The pull handle shall not be used to lift the system. Only the designated lifting points and handles shall be used.

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During transport, secure the system to prevent impact or mechanical damage.

Operator position: When an alarm occurs, the operator shall be positioned facing the display to clearly view visual alarm indicators and obtain alarm information for appropriate action.

Display: The touchscreen of the display is a pressure-sensitive surface. Do not allow contact with sharp objects or strong chemical substances.

Environmental considerations: Environments exceeding normal atmospheric pressure may damage the Lifetech Cardio intra-aortic balloon pump. To ensure proper operation, observe the following precautions during operation, transport, and storage.

- The device shall be used in clinical settings only by trained professionals.
- Prior to operation, the operator shall properly discharge static electricity.
- Do not subject the device to physical damage such as mechanical impact.
- Do not place the device in environments outside the operating and storage temperature and humidity ranges specified in this manual (see Device Specifications).
- Do not place the Lifetech Cardio intra-aortic balloon pump within the patient's reach.
- Do not immerse the Lifetech Cardio intra-aortic balloon pump or any of its components in water.
- Do not sterilize the device using gamma radiation, steam, or similar methods.
- If there is a sudden change in ambient temperature or humidity, allow the Lifetech Cardio intra-aortic balloon pump to reach ambient conditions before use.
- Do not place liquids on top of the device. Ensure that saline containers and tubing are not suspended directly above the IABP. If spillage occurs, wipe it off immediately and service the device to prevent hazards.
- Cables shall be replaced promptly if their performance changes due to aging or environmental conditions.
- To prevent condensation, when moving the IABP from a cold area to a warm area, allow it to warm up. It is recommended to leave the device in the warmer area for at least 1 hour before use.

Electromagnetic compatibility information: Refer to **8.15 Electromagnetic Compatibility**. The device shall be used only in the electromagnetic environment specified in that section.

Potential adverse environmental factors include, but are not limited to, those listed above. In principle, the device shall be operated with priority given to patient health and safety, minimizing environmental impact on the device.

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Timing: Improper timing settings may result in suboptimal therapeutic effect or cause patient injury.

Trigger source: Improper trigger source settings may result in suboptimal therapeutic effect or cause patient injury.

2.3. Safety Information

2.3.1. Measures to Prevent User Misoperation

Screen lock:

To prevent unintended contact, the display screen can be locked manually. In addition, the system includes an automatic lock function: after initial power-on, once the IABP has initiated an assist cycle, the display screen will automatically lock after 2 minutes of inactivity. To unlock the screen, press the corresponding button.

When the screen is locked, all currently displayed items will appear “grayed out” but will not be.

2.3.2. Abnormal conditions

In the event of user interface errors, application logic errors, or system availability failures, continued operation of the device cannot be guaranteed. It is recommended to open the Settings menu, select Shutdown, and then select Confirm to power off the IABP. After waiting 10 seconds, press and release the main unit power button “ ”to restart the device. If the above shutdown procedure cannot be performed, press and hold the main unit power button “ ”for 5 seconds to force shutdown, wait 10 seconds, and then restart the device.

When no waveform is selected for report output or when shutting down the software, a secondary confirmation dialog will be displayed.

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3. User Interface

The user interface of the Lifetech Cardio IntraNova® IABP is described in detail below.

3.1. Display and Operation

This section of the operating instructions describes the display interface features and control methods. For additional operating instructions, refer to Chapter 4, Device Operation.

PRECAUTION: The touchscreen of the Display is a pressure-sensitive surface. Do not allow contact with sharp objects or strong chemical substances.

The display layout of the Lifetech Cardio IABP is shown below.

Figure 3-1: Example of Display Layout

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

Depending on the current balloon pressure waveform settings (located in Settings → Display, see Section 3.1.24.1), the IABP can display two waveforms (ECG and AP) or three waveforms (ECG, AP, and BPW). On the display, the waveforms are arranged from top to bottom as ECG, AP, and BPW. From left to right, each waveform is drawn by continuously replacing old data with new data.

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3.1.1. Alarm Sound Pause / Cancel Alarm Sound Pause

1 Alarm Icon

Figure 3-2: Example of Alarm Icon, Information Display Area, and the Icons/Indicators Area

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

Pressing the alarm icon temporarily silences the active alarm sound for 60 seconds, and the icon changes to the alarm sound pause icon . This control does not override the alarm condition. If the alarm condition is not resolved within 60 seconds, the alarm sound will be reactivated. While the audio is muted, the alarm message remains displayed with its active visual attributes. If a new alarm occurs during the alarm sound pause period, the alarm sound will be immediately reactivated. In addition, pressing the alarm sound pause icon again during the 60-second pause period will immediately restore the alarm sound.

3.1.2. Information Display Area

2 Information Display Area

When there are no alarms or messages, the center area displays time information. When alarms or messages are present, the center area displays alarm or message information, and the time information is shown on the left side.


Figure 3-2 ： Two Examples of the Information Display Area

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

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If the total number of alarms and messages in the information display area exceeds two, alarms or messages with lower priority will be collapsed. Press the down arrow in the information display area to view the collapsed alarms or messages. When expanded, up to four alarms or messages can be displayed. Additional alarms or messages can be viewed by pressing the down arrow to scroll through the list.

Alarm and message priorities are classified as high-priority alarms, medium-priority alarms, low-priority alarms, and messages. All items are displayed according to their predefined priority order, with higher-priority alarms shown first and lower-priority items shown afterward. Once an alarm condition is corrected, the corresponding alarm or message will be cleared. Selecting an alarm opens its help information.

PRECAUTION:

If the total number of alarms and messages in the information display area exceeds two, alarms or messages with lower priority will be collapsed. Press the Up "▲" or Down "▼" buttons in the information display area to view the collapsed alarms or messages

3.1.3. Icons and Indicators Area

3 Alarm Disable Icon

The alarm disable icon has two states.

When the gas leakage alarm or the front-end helium circuit leakage alarm is disabled, the alarm disable icon displays a crossed solid line.

When the catheter alarm is temporarily paused, the alarm disable icon displays a crossed dashed line.

4 Helium Indicator Icon

The helium indicator icon indicates the approximate remaining helium level. When the Helium Cylinder is full, the indicator is shown as a fully filled green cylinder. As helium is consumed, the fill level decreases accordingly. When the helium level becomes low, the fill color changes from green to yellow (30%) and then to red (10%), and continues decreasing until the cylinder is empty.

PRECAUTION:

When the helium cylinder level becomes too low, a "low helium" message will be displayed; the helium indicator will turn red and replacement will be required

5 Battery Status Icon

The battery status icon indicates the current status of the lithium battery. When

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AC power is connected and the battery is being charged, a lightning symbol appears on the battery icon. When AC power is connected and the battery is fully charged, the lightning symbol is not displayed. When the device is powered by the lithium battery, the lightning symbol disappears. As the battery charge level decreases, the fill level inside the battery icon gradually decreases accordingly. Similarly, as the battery charge level drops, the battery icon color transitions from green to yellow (30%) and then to red (10%).

The numeric battery charge level and the current power source (AC or DC) are displayed to the right of the lithium battery icon.

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3.1.4. ECG Waveform

Figure 3-3: Example of ECG Waveform

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

6 ECG Lead

The ECG lead field displays the selected (active) ECG signal source and its status.

In automatic operating mode, the system automatically selects an ECG signal source (I, II, III, or External). This field displays the currently automatically selected (active) ECG signal source, for example: Lead: Auto-II.

In semi-automatic operating mode, the user can select the ECG signal source (I, II, III, aVR, aVL, aVF, V, or External). This field displays the ECG signal source currently selected by the user based on the connected electrodes, for example: Lead: Semi-auto-II.

- Regardless of the ECG signal source, ECG signals must be acquired via the RA, LA, and LL electrodes. If a single lead becomes disconnected, a lead-off message is displayed in this field, for example: –LA Lead Fault. This message indicates the currently selected ECG signal source and the lead that requires replacement to restore ECG signal acquisition. If multiple leads are disconnected, the system displays a general lead-off or lead fault message without specifying an individual lead. If the cable is disconnected, –No Cable is displayed.

 Lead-off current level: < 100 nA.
- When the External source is selected and no cable is connected, the message –External – No Cable is displayed.。

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7 Filter

The filter field displays the current filter setting.

8 ECG Signal Source Selection

- In automatic operating mode, the system automatically selects an ECG signal source (I, II, III, or External), which is displayed between the left and right arrow buttons. No other selection is permitted. This field displays the currently automatically selected (active) ECG signal source, for example: Lead: Auto-II.
- In semi-automatic operating mode, the user can use the “<<<” and “>>>” buttons to scroll left or right to select an available ECG signal source (I, II, III, aVR, aVL, aVF, V, or External). Selection must be based on the connected electrodes. The current setting is displayed between the left and right arrow keys. This field also displays the ECG signal source currently selected by the user, for example: Semi-auto-II.

9 Pacemaker Spike Pulses

When pacemaker triggering is used (see Section 4.1.6), the ECG waveform area displays pacemaker spike pulses as gray vertical lines extending upward from the ECG waveform.

10 ECG Waveform

The ECG waveform is drawn from left to right. The oldest data are erased and replaced with new data. To facilitate identification of the most recent data, a vertical erase line is displayed. The newest data appear to the left of the line, and the oldest data appear to the right.

11 ECG Scale

The ECG scale is located immediately to the left of the displayed ECG waveform. This scale can be used to determine ECG amplitude. The indicator bar has a height of 1 cm, and the annotation indicates the current scaling factor of the displayed ECG waveform.

PRECAUTION:

The displayed ECG signal is automatically scaled to provide optimal screen presentation. Due to automatic scaling, low-amplitude ECG signals may appear with normal amplitude on the display. The ECG amplitude of interest shall be interpreted according to the indicated scale.

12 ECG Markers

ECG markers are based on the length and position beneath the ECG waveform.

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These horizontal lines indicate the calculated inflation (assist) timing intervals. Display of these markers can be enabled or disabled in Settings → Display → ECG Markers (see Section 3.1.24.1).

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3.1.5. Trigger Signal and Heart Rate Parameters

Figure 3-4: Example of Trigger Signal and Heart Rate Parameters

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

13 Trigger Indicator

- When ECG, pressure, V/A-V pacemaker, or A pacemaker is selected as the trigger source, this indicator is displayed as a heart-shaped icon and flashes once (on/off) for each detected trigger event.
- When internal triggering is selected, this indicator is displayed as a diamond-shaped icon and flashes with each internal trigger event.

14 Heart Rate

The heart rate field displays the average heart rate over the most recent four cardiac cycles in beats per minute (bpm). When the heart rate is below 15 bpm, the IABP displays zero (0). When no signal is present or the signal is invalid, "---" is displayed. For technical specifications, refer to the heart rate meter in Section 8.4.

WARNING: Displayed heart rate values may fluctuate significantly during arrhythmias.

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3.1.6. Arterial Pressure Waveform

Figure 3-5: Example of Arterial Pressure Wavefor

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

15 AP Source Status

- This field displays the currently selected pressure signal source and its status (if applicable).
- In semi-automatic operating mode, if the pressure source selected by the operator is detected as disconnected, the IABP displays “Signal Source Disconnected.”
- In automatic operating mode, if no pressure source is available, the IABP displays “No Signal Source.”
- If the Disposable Blood Pressure Sensor has not been zeroed, “Not Zeroed” is displayed. If zeroing is successful, “Zeroed” is displayed. If zeroing fails, “Zeroing Failed” is displayed.
- If fiber-optic sensor identification fails, “Identification Failed” is displayed. If the sensor is identified but not zeroed, “Connected – Not Zeroed” is displayed. If zeroing fails, “Connected – Zeroing Failed” is displayed. If zeroing is successful, “Connected – Zeroed” is displayed.

16 AP Source Selection

This field allows selection of the AP Source, as described below.

- In semi-automatic operating mode, the operator can use the “<<<” and “>>>” buttons to select the AP Source (Fiber, Fluid-filled, or External). The current selection is displayed between the left and right arrow buttons (the example above shows External pressure selected). The operator can

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manually change the AP Source.

- In automatic operating mode, the device automatically selects the AP Source. Although user selection is permitted, the system automatically chooses the most appropriate pressure source. The priority order is: Fiber > Fluid-filled (Zeroed) > External > Fluid-filled (Not Zeroed) > Fluid-filled (Not Connected). When the AP Source is set to “Fiber,” the fiber-optic sensor is used as the direct pressure signal source and has priority over a directly connected Disposable Blood Pressure Sensor. The fiber-optic balloon sensor connector must be disconnected to enable input from the blood pressure sensor.
- When the AP Source is set to External (in the signal source menu), this field displays “External,” indicating that the signal is provided by an external AP Source (such as a patient monitor).

17 Arterial Pressure Waveform Scale

The arterial pressure waveform scale is located to the left of the displayed arterial waveform. The scale is a vertical line with five equally spaced tick marks. Numeric annotations indicate the maximum, mean, and minimum arterial pressure values visible within the waveform area. The measurement unit is mmHg, with a display range of –30 to 300 mmHg. This scale facilitates measurement of arterial pressure amplitude. If the pressure has not been zeroed, the arterial pressure waveform scale is not displayed.

18 Arterial Pressure Waveform

The arterial pressure waveform is drawn from left to right. The oldest data are erased and replaced with new data. To facilitate identification of the most recent data, a vertical erase line is displayed. The newest data appear to the left of the line, and the oldest data appear to the right.

In standby mode, the assist timing intervals are continuously calculated and the corresponding intervals are highlighted on the arterial pressure waveform.

During assist operation, highlighted inflation timing intervals are not displayed. If it is necessary to temporarily view the highlighted inflation timing intervals during assist, press the “Show Inflation Interval” button located in Settings → Display (see Section 3.1.24.1).

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Pressing the “Show Inflation Interval” button highlights the inflation timing interval for 15 seconds. Pressing the button again refreshes the display.


Figure 3-6: Example of Highlighted Inflation Timing Interval

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

PRECAUTION:

Do not set timing based on highlighted inflation interval markers. Set timing by evaluating the arterial pressure waveform.

When pressure is selected as the trigger source, thin horizontal lines are overlaid on the arterial pressure waveform. On the rising edge of systole, these arterial pressure trigger event markers provide a visual indication of the current arterial pressure trigger threshold. These markers do not indicate the exact timing of the trigger event and are displayed with a slight time delay.

19 Reference Line with Numeric Pressure Readout

On the arterial pressure waveform, the user can control a horizontal cursor to measure arterial pressure at any point. If the pressure indication is valid, a numeric value corresponding to the pressure at the reference line position (in mmHg) is displayed to the right of the line. When the sensor is not zeroed, the numeric readout of the reference line is not displayed. The user can move the reference line into the balloon pressure waveform area. When the reference line is located within the balloon pressure waveform area, the numeric value is not displayed. The display and position of the reference line can be controlled in the reference line control area. When the reference line is enabled for the first time, it is positioned at the mean pressure value and can be moved using the up or down buttons (see Reference Line Control Area).

In this area, the Up and Down buttons are used to move the reference line to the desired measurement point on the waveform. If the reference line is at the highest position of the arterial pressure waveform area, the Up button is grayed out, indicating that the reference line cannot be moved higher. Similarly, if the reference line is at the lowest position of the balloon pressure waveform area, the Down button is grayed out, indicating that the reference line cannot be moved lower. Pressing the Up or Down buttons moves the reference line accordingly. Pressing and holding the Up or Down button accelerates the movement of the

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reference line, allowing rapid traversal across the waveform area.

3.1.7. Arterial Pressure–Related Parameters

20 Arterial Pressure–Related Parameters

This area displays systolic pressure, diastolic pressure, and mean pressure in mmHg (see Figure 3-6). For directly connected arterial pressure values, the valid display range is –30 to 300 mmHg (inclusive). If the arterial pressure exceeds 300 mmHg, “300” is displayed. If the arterial pressure is invalid or not zeroed, three dashes “---” are displayed in place of the numeric value. During assist operation, if the IAB ratio is set to 1:2 or 1:3, the unassisted systolic pressure and unassisted diastolic pressure are also displayed.

3.1.8. Balloon Pressure Waveform

Figure 3-7: Example of Balloon Pressure Waveform

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

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21 Standby Time

Standby time (in minutes) is displayed only when the IABP is in standby mode. When the IABP enters standby mode for any reason, the elapsed standby time is displayed and cleared when assist is restarted. If the standby time is less than 1 minute, a “<1 minute” message is displayed until the standby duration reaches 1 minute or longer.

If the IABP remains in standby mode for 10 minutes, a “Standby Time Too Long” message remains displayed in the information display area. If the IAB catheter remains in the deflated state for an extended period, standby time alert tones provide a reminder of the risk of thrombus formation. A 30-second double-beep tone is emitted at 10, 15, and 20 minutes. After 20 minutes, the 30-second double-beep tone sounds every 2 minutes until assist is initiated.

22 Balloon Status Indicator

The balloon status indicator depicts inflation and deflation of the intra-aortic balloon during assist. It does not provide quantitative measurement results, nor is it intended for such calculations.

23 Balloon Pressure Waveform

The balloon pressure waveform is drawn from left to right. The oldest data are erased and replaced with new data. To facilitate identification of the most recent data, a vertical erase line is displayed. The newest data appear to the left of the line, and the oldest data appear to the right.

3.1.9. Augmentation Pressure Information Area

Figure 3-8: Example of Augmentation Pressure Information Area

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

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24 “Aug” Value

The Aug value indicates the peak augmentation pressure in mmHg. The valid display range is 0 to 300 mmHg (inclusive). If the arterial pressure exceeds 300 mmHg, “300” is displayed. The augmentation pressure value is not displayed under the following conditions:

- The IABP is in standby mode
- When the fiber-optic AP Source is selected and the fiber sensor is not connected or not zeroed
- When the fluid-filled AP Source is selected and the Disposable Blood Pressure Sensor is not connected or not zeroed
- When the external AP Source is selected and the external cable connector is not inserted.

25 “Aug Alarm Off” Status

This field displays the status of the augmentation pressure alarm. If the augmentation pressure alarm is disabled, the alarm disabled icon is displayed.。”

3.1.10. “Zero” Button

Figure 3-9: Example of Sidebar Menu / Buttons

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

26 “Zero” button

The “Zero” button applies only to the two types of direct pressure signal sources: Fluid-filled and Fiber-optic.

When the signal source is Fluid-filled, press and hold the “Zero” button to perform zeroing.

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When the signal source is Fiber-optic, confirm that the “Fiber Zeroing Method” setting is correctly configured before performing zeroing.

- To perform in vivo zeroing, select “In Vivo Zeroing” in Settings → System Options → Tools → Fiber Zeroing Method, then press and hold the “Zero” button to execute (see Section 3.1.24.4, System Options Menu).
- To perform in vitro zeroing, select “In Vitro Zeroing” in Settings → System Options → Tools → Fiber Zeroing Method, then press and hold the “Zero” button to execute (see Section 3.1.24.4, System Options Menu).

PRECAUTION: Switching the fiber zeroing method is not permitted during augmentation therapy. In this case, the corresponding buttons in the Fiber Zeroing Method menu are grayed out and unavailable.

PRECAUTION: The factory default fiber zeroing method of the IABP is In Vivo Zeroing. Once updated, the selected zeroing method is saved and will not change after device shutdown or restart.

When the system is using an external pressure signal source, the Zero button is grayed out and unavailable.。

3.1.11. “Start” / “Pause” Button

Figure 3-10: Example of “Start” / “Pause” Button

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

27 “Start” / “Pause” Button

If the IABP is not providing assist, pressing this button starts assist (if the current settings allow), and the button displays “Start” to indicate that assist is active. If the IABP is already providing assist, pressing this button pauses assist, and the button displays “Pause” to indicate that assist has been paused.

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When the device is used for the first time after power-on, pressing the **Start** button automatically initiates the IABP, accelerating assist activation. The device automatically performs the following steps:

- Automatic deflation and inflation of the IAB.
- Initiation of assist.

If the IABP is providing assist and pressure is selected as the trigger source, pressing and holding this button sends a request to resynchronize arterial pressure. Timing is immediately resynchronized. During this process, temporary interruption of heart rate–based assist may occur.

3.1.12. “Operating Mode” Menu

Figure 3-11 Example of Operating Mode Menu

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

28 “Mode” Button

The IABP provides two operating modes: Automatic (for basic operation) and Semi-Automatic (for advanced operation). The current operating mode is indicated by a gray background on the corresponding mode button. The default operating mode is Automatic.

“AUTO” Button: This button sets the IABP to Automatic operating mode.

“Semi” Button: This button sets the IABP to Semi-Automatic operating mode. In Semi-Automatic mode, five trigger sources can be selected: ECG, Pressure, Pacer A, Pacer V/AV, and Internal. For details, see Section 4.1.3.2, Semi-Automatic Operating Mode.

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3.1.13. Trigger Signal Source Selection

Figure 3-12 Example of Trigger Source Selection Menu

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

29 "Trigger" Button

After selecting Semi-Automatic mode, the following trigger sources can be selected: ECG, Pressure, Pacer A, Pacer V/AV, and Internal. Regardless of operating mode, the active trigger source is displayed in the trigger source selection field.

3.1.14. Reference Line Control Area

30 Up "▲" or Down "▼" buttons

The Up or Down buttons move the reference line. Pressing and holding a directional button accelerates the movement of the reference line, allowing rapid traversal across the waveform area.

PRECAUTION

The reference line can also be moved into the balloon pressure waveform area, in which case the numeric pressure value to the right of the reference line is not displayed.

PRECAUTION

If arterial pressure is invalid or not zeroed, the numeric pressure value to the right of the reference line is not displayed.

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3.1.15. Screen Lock

Figure 3-13: Example of Bottom Menu / Buttons (Example 1)

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

31 "🔒"Lock and"🔓"Unlock buttons

To prevent unintended touches, pressing this button"🔒" locks the display screen. The lock can be activated manually. In addition, the system includes an automatic screen lock function.

When the screen is locked, all displayed items appear grayed out but remain active and unaffected. Touch operations are disabled. Press the button"🔓" again to unlock the screen.

Upon initial power-on, after the IABP completes one assist cycle, the display screen automatically locks after 2 minutes of inactivity. Press the button to unlock the screen.

When high-, medium-, or low-priority alarms occur, or when a standby time pop-up appears, the display screen automatically unlocks and remains unlocked for the duration of the alarm condition.

The automatic screen lock time can be adjusted in Settings → System Options → Tools.

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3.1.16. Timing Selection Menu

32 **“Timing” button**

Pressing the “**Timing**” button displays the inflation timing selection indicator and the deflation timing selection indicator.


Figure 3-14 Example of Timing Selection Menu

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

1) Inflation Timing Selection Indicator

This graphical indicator provides visual feedback for the selected time between the trigger event and the start of inflation. The inflation timing selection range is –4.0 to 4.0 (Semi-Automatic mode).

In Semi-Automatic mode, the inflation timing can be adjusted using the left and right directional buttons of the inflation timing selection indicator.

In Automatic mode, the word “Auto” appears in the center of the indicator, indicating that the inflation timing is set automatically and does not require operator intervention.

The inflation timing selection control adjusts the duration between the trigger event and the start of inflation. Press and hold the left “◀” or right “▶” directional button to move inflation to the desired timing point, then release the control.

- Press the “Reset” button below the inflation timing selection indicator to quickly reset the inflation timing to 0.
- The left“◀”/right “▶” directional buttons and the Reset button are adjustable only in Semi-Automatic mode. In Automatic mode, these buttons are disabled because inflation timing selection is set automatically.
- Movement of the inflation timing selection indicator depends on the user’s adjustment of the inflation timing.

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2) Deflation Timing Selection Indicator

The deflation timing selection indicator adjusts the timing at which balloon deflation begins. Press and hold the left "◀" or right "▶" directional button to move deflation to the desired timing point, then release the control. This graphical indicator provides visual feedback for the selected deflation timing. The deflation timing selection range is –4.0 to 4.0 (Semi-Automatic and Automatic modes).

- Press the "Reset" button below the deflation timing selection indicator to quickly reset the deflation timing to 0.
- In Semi-Automatic mode, the deflation timing can be adjusted using the left and right directional buttons of the deflation timing selection indicator.
- Deflation timing selection is adjustable in all operating modes except when pressure is selected as the trigger source in Automatic mode. In this case, the word "Auto" appears in the center of the indicator, indicating that deflation timing is set automatically and operator intervention is not permitted.
- When the ECG trigger source is used (Automatic mode), the user may be allowed to modify the automatically determined deflation point, although this is not required. Unlike inflation timing selection, deflation timing selection can vary significantly. This option allows the user to set deflation timing based on experience.
- The deflation timing selection control also allows the user to change the deflation timing selection mode. Moving the deflation timing selection indicator to the far right (end) position manually enables R-wave Deflation mode. The label "R-wave Deflation" will appear on the deflation timing selection indicator, and an R-wave deflation message will appear in the balloon pressure waveform status area. R-wave Deflation mode delays IAB deflation until the next valid R-wave trigger event occurs.
- When R-wave Tracking is enabled (see Section 3.1.24.4, R-wave Tracking), the IABP performs Automatic R-wave Deflation whenever an unpredictable rhythm pattern (e.g., atrial fibrillation) is detected. The deflation timing selection indicator automatically moves to the far right (end) position and displays "Auto R-wave Deflation." During Auto R-wave Deflation, when a predictable rhythm is restored, the deflation timing selection returns to its previous value. The IABP then continues monitoring for unpredictable rhythms.
- If a cardiac cycle is detected earlier than expected, deflation will begin immediately, overriding the preset deflation point.

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Press the“Assist Ratio” button to open the “Assist Ratio” setting menu. The current ratio is indicated by a blue background on the corresponding button. The default assist ratio is 1:1.

Select the IABP assist ratio using the controls in the Assist Ratio menu. Available selections include:

- 1:1 – Counterpulsation on every assisted beat.
- 1:2 – Counterpulsation on every two assisted beats.
- 1:3 – Counterpulsation on every three assisted beats.

3.1.18. “Augmentation Pressure Level” Menu

34 “Augmentation” button

The Augmentation Pressure Level menu includes the current pressure level display and Up/Down buttons, allowing adjustment under any operating mode and any trigger source. The Up and Down buttons increase or decrease the augmentation pressure level in a range of 1 to 10 units, with an increment of 1 unit per step, for a total of 10 levels. When the level is at the maximum, the Up button is grayed out, indicating it cannot be increased further. Similarly, when the level is at the minimum (Off), the Down button is grayed out, indicating it cannot be decreased further.

3.1.19. “Pressure Threshold” Menu

35 “Thresholds” button

In Automatic mode, the word “Auto” is displayed under “Pressure Threshold,” and the user is not allowed to adjust the pressure threshold.

In Semi-Automatic mode, the pressure threshold can be increased or decreased using the Up and Down buttons within a range of 7 to 30 mmHg, in 1 mmHg increments. The current setting is displayed between the Up “▲” and Down “▼” buttons. After adjustment, the label “Auto” changes to “Manual.”

3.1.20. “Augmentation Pressure Alarm” Threshold Setting Menu

36 “Aug Alarm” button

After assist is started, the augmentation pressure alarm is automatically set after 3 minutes (zeroing is required if “Hydraulic” or “Fiber” pressure signal source is used). The alarm limit is automatically determined by the program and is set to 10 mmHg below the patient’s augmentation pressure. This function also allows the user to adjust the augmentation pressure alarm limit (the pressure that triggers the alarm, in mmHg).

The Up “▲” and Down “▼” buttons increase or decrease the augmentation pressure alarm limit by 2 mmHg per step, within a range of 60 to 200 mmHg. The current setting is displayed between the Up and Down buttons. Press and hold the Up or Down button to set the alarm limit; holding the button accelerates the adjustment to quickly reach the desired value. If the augmentation pressure

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alarm is Off, pressing the Up or Down button will automatically set the alarm to 10 mmHg below the patient's augmentation pressure.

To disable the augmentation pressure alarm, adjust the alarm limit beyond the 60–200 mmHg range to turn the alarm Off (displayed as "OFF"). The augmentation pressure alarm icon will indicate the disabled status (see Section 3.1.9, Augmentation Pressure Information Area).

About **3 minutes after assist is initiated** and valid pressure values are available, verify that the augmentation pressure alarm has been automatically set to approximately **10 mmHg below** the patient's augmentation pressure.

WARNING:

Setting the augmentation pressure alarm limit too low or disabling it will render the alarm ineffective.

3.1.21. "IAB Fill" Button

Figure 3-15 Example of Bottom Menu / Buttons (Example 2)

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

37 "IAB Fill" Button

Press and hold this button for 2 seconds to start the inflation procedure. If augmentation is in progress, augmentation will be interrupted. The automatic inflation procedure evacuates the gas in the tubing (IAB catheter and helium extension line) and replaces it with pure helium.

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After inflation is completed, a “Inflation Complete” message is issued and remains displayed until augmentation is resumed.

Because the IABP automatically evacuates and replaces the gas in the tubing every 2 hours, operator intervention is usually not required. Therefore, this button is not used frequently. Typically, it is used to address alarms when corrective actions are prompted by the help screen.

3.1.22. “Freeze” Button

38 **“Freeze” Button**

Press this button to **freeze** all waveforms displayed on the screen. Press again to **unfreeze** all waveforms. This button is disabled when the print function is enabled. If the user is adjusting timing selection or presses the waveform print button, the display will automatically unfreeze.

3.1.23. “” Report Output Button

39 **“”Button**

When no storage device is inserted and the device is not connected to a printer, the report output button “” is grayed out and unavailable. When the system detects a storage device or printer, the report output button becomes available. Press the button “” to initiate report output according to Settings → Report Output → Output Method (see Section 4.7), to generate a PDF to the storage device or print the document.

3.1.24. “Settings” Menu

40 **“Settings” Button**

Press the “Settings” button to open the Settings menu, where the user can access the Display, Sound, Report Output, System Options, Alarm Log, and About submenus. By default, the Display menu is always shown first. Regardless of the selected submenu, the Power Off option is always located in the lower left corner.


Figure 3-16：Example of Display Menu

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3.1.24.1. “Display” Menu

The Display menu allows selection of scan speed, balloon pressure waveform, ECG markers, and inflation interval, as described below:

- **Scan Speed**: Sets the scan speed for displayed waveforms and the printer. Options: 25 mm/s and 50 mm/s. Default: 25 mm/s.
- **Balloon Pressure Waveform**: Enables or disables balloon waveform display. Options: On and Off. Default: On. For balloon pressure waveform information, see Section 3.1.8.
- **ECG Markers**: Enables or disables ECG marker display. Options: On and Off. Default: On. For ECG marker information, see Section 3.1.4.
- **Inflation Interval**: The inflation interval is the marked and highlighted segment on the arterial pressure waveform and is used to view the approximate balloon inflation period. During assist, the highlighted inflation interval is not displayed. To temporarily view it during assist, press “Show Inflation Interval” to display or refresh the inflation interval highlight for 15 seconds.

3.1.24.2. “Sound” Menu

Allows adjustment of alarm volume, adjustment of QRS beep volume or turning off the QRS beep, and enabling/disabling the standby reminder tone.

- **Alarm Volume**: Manually sets the alarm volume with 8 levels (minimum 1, maximum 8). Default volume level is 6. Each adjustment plays a sound sample to allow the operator to evaluate the volume.
- **QRS Beep Volume**: Sets the QRS trigger beep volume. Options: Off, Low, Medium, and High. Default: Off.
- **Standby Reminder Tone:** Enables or disables the reminder tone for excessive standby time. Options: On and Off. Default: On.

WARNING:

Adjust the alarm sound pressure level with caution. Alarm sound levels below ambient noise may impair recognition of alarm conditions.

3.1.24.3. “Report Output” Menu

The Report Output menu allows selection of waveforms, scheduled output, and output method, as described below:

- Waveforms: Defines the vertical layout of printed waveforms. Multiple selections are supported: one waveform, two waveforms, or all three waveforms may be selected. Default: all three waveforms. For printer format details, see Section 4.7.
- Scheduled Output: Sets the interval for automatic waveform printing. Options: Off, 15 min, 30 min, 1 h, 2 h, and 4 h. Default:

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Off.

- Output Method: Selects output to a storage device for saving or to a connected printer for printing. Default: output to the storage device.

3.1.24.4. “System Options” Menu

This menu allows selection of R-wave tracking, catheter alarms, pacer detection settings, fiber manual zero adjustment, and date/time settings, as described below:

- R-wave Tracking: Enables or disables R-wave tracking to track unpredictable rhythm patterns. Options: On and Off. Default: On.
- Catheter Alarms: Enables or disables catheter alarms. Options: Catheter Alarms On, Leak Alarm Off, and Pause Catheter Alarms. Default: Catheter Alarms On. The Pause Catheter Alarms setting lasts only 60 minutes; if needed, the operator must re-enter and re-set it. To enable Pause Catheter Alarms, press and hold the pause button for 2 seconds. Pausing catheter alarms suspends alarms for IAB circuit air intake, IAB circuit leakage, and excessive balloon pressure.
- Pacer Detection: Under any operating mode and regardless of the selected trigger source, the pacer detection level can be adjusted using the left/right directional buttons. There are 4 pacer detection levels, displayed numerically in the area between the left/right buttons.
- Fiber Manual Zero Adjustment: The left/right directional buttons can be used to adjust the current fiber baseline pressure, shifting the overall fiber pressure upward or downward.
- Date/Time: Sets the current time and date. Time is displayed in 24-hour format with two-digit hours and two-digit minutes. The date is displayed as an 8-digit number (4-digit year, 2-digit month, 2-digit day). The Up “▲”/Down “▼” directional buttons adjust all fields. After setting the time and date, press Confirm to save and apply the changes.

3.1.24.4.1 “Tools” Submenu of “System Options”

Press “Tools” to enter the Tools menu, where you can select the fiber zeroing method, the helium circuit front-end leak alarm, and the auto screen-lock time, as described below:

- Fiber Zeroing Method: Select “In-vivo Zeroing” or “In-vitro Zeroing” depending on the scenario. If the balloon has not been implanted in the patient and is fully exposed to ambient air, select “In-vitro Zeroing.” If the balloon has been advanced to the target vessel position in the patient,

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select “In-vivo Zeroing.”

PRECAUTION: Switching the fiber zeroing method is not allowed during counterpulsation therapy. In this case, the corresponding button in the “Fiber Zeroing Method” menu will be grayed out and unavailable.

PRECAUTION: The factory default fiber zeroing method is In-vivo Zeroing. Once updated, the zeroing method is saved and will not change due to power-off or restart of the IABP.

- Helium Circuit Front-end Leak Alarm: The helium circuit front-end leak alarm can be turned On/Off. When the alarm “Helium circuit front-end leak detected” has not occurred, this alarm cannot be turned off. It can only be turned off after the alarm “Helium circuit front-end leak detected” occurs. Each time the device is powered on, the helium circuit front-end leak alarm is set to On by default.

 WARNING:

 If an IABP has a helium front-end leak, continuous helium supply cannot be guaranteed. For long-term counterpulsation therapy, address the helium leak promptly or replace with another available Lifetech Cardio IABP.

 WARNING:

 After the helium pathway front-end leakage issue has been resolved, promptly set the “Helium pathway front-end leakage” alarm to “On”. The IABP will re-detect whether leakage is present at the front end of the helium pathway.

- Auto Screen-lock Time: The auto screen-lock time can be adjusted within a range of 0–120 min. Setting it to 0 disables auto screen-lock; 0 is displayed as “OFF.”

3.1.24.4.2 “Device Information” Submenu of “System Options”

- Device Information summarizes the device total operating time, device total assist time, operating time of the current fiber (module), cycle count of the current safety disk, battery cycle count, and battery capacity.

3.1.24.5. “Alarm Log” Menu

Lifetech Cardio maintains a non-editable alarm log here, containing high-, medium-, and low-priority alarms that occurred since the current power-on, with storage of up to 300 entries. If the power-off duration exceeds 3 minutes, the alarm log prior to shutdown will be cleared. The log can be

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accessed in this menu. The alarm log is sorted from most recent to oldest by the actual alarm occurrence time.

3.1.24.6. "About" Menu

- The About menu displays the product model, software version, company, release date, copyright, and other information.
- Pneumatic Test: A pneumatic test can also be performed in the About screen. Press "Pneumatic Test" to enter the pneumatic test interface. On the test screen, press "Start Test." Follow the prompt: when "Please use a Luer plug to block the IAB catheter port" appears, tightly insert the Luer cap into the IAB catheter port. Then tap "Blocked." When the test is completed, the leak test result will be displayed. If the system passes, a green pass message will appear; otherwise, failure information will be displayed. See Section 4.2.

 If there is any doubt about the integrity of the pneumatic module, please contact Lifetech Cardio service.
- Power-line Filter: Press this button to enable/disable the mains-frequency (power-line) filter on the ECG channel.

3.1.24.7. "Power Off" Button

Regardless of which submenu is selected, the "Power Off" button is always located in the lower-left corner of the Settings menu. Press "Power Off" and a dialog "Please confirm whether to power off" will pop up. Press "Confirm" to shut down the device.

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3.2. Panels

A schematic of the front panel interfaces is shown in Figure Figure 3-17。


Figure 3-17 A schematic of the front panel interfaces

1 Status Display

The status display shows the lithium battery level, remaining helium volume, and operating time.

2 ECG Input Connector

The ECG input connector is a 12-pin waterproof connector for the accessory ECG cable.

3 External ECG Input Connector

The external ECG input connector is used to connect an external ECG signal cable. See Chapter 7 for external monitoring device interfaces.

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4 Blood Pressure Input Connector

The blood pressure input connector is used to connect an invasive blood pressure sensor cable.

5 External Blood Pressure Input Connector

The external blood pressure input connector is used to connect an external pressure signal cable. See Chapter 7 for external monitoring device interfaces.

6 ECG Output Connector

The ECG output connector is used to connect an external ECG signal cable. See Chapter 7 for external monitoring device interfaces.

7 Blood Pressure Output Connector

The blood pressure output connector is used to connect an external pressure signal cable. See Chapter 7 for external monitoring device interfaces

8 Fiber-Optic Sensor Connector

This dedicated connector is used to detect fiber-optic sensor insertion and to read fiber-optic parameter information.

9 Main Console Power Button“ ”

The button “ ” includes a power indicator symbol.

When the device is powered off, pressing the button once starts the system, accompanied by a beep indicating system startup.

When powered on, a short press “ ” opens the Power Off menu.

When powered on, pressing and holding the button “ ” for 5 seconds forces the device to shut down.

The green indicator light is on when operating on AC power.

The green indicator light flashes when charging the lithium battery while the device is powered off.

The orange indicator light is on when operating on lithium battery power.

The orange indicator light flashes when the lithium battery is at a low charge level.

When the device is powered off and the lithium battery is fully charged, the indicator light is off.

10 IAB Catheter Connector

The IAB catheter connector is used to connect the IAB catheter extension tube. It is a standard female Luer connector, white in color.

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
Figure 3-18 A schematic of the rear panel interfaces

11 Main Console Video Interface

The Main Console Video Interface is a Type-C connector used to connect the display.

12 USB interface

The USB interface is a USB-A female connector used for report and log output (see Section 4.7).

13 Equipotential Bonding Terminal

The equipotential bonding terminal is a circular M6 socket used for device grounding or equipotential bonding with other hospital equipment.

14 Power Inlet

The power inlet is used to connect the power cord.

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15 Power Inlet Switch

The power inlet switch includes an indicator light and symbol. When connected to a valid AC power source and switched on, the indicator remains lit; it turns off when switched off.

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3.3. Helium Cylinder Compartment

Sealing Gasket: This gasket is used to replace the old gasket when changing the Helium Cylinder, improving airtightness at the cylinder connection.

Helium Cylinder Knob: Used to open or close the Helium Cylinder. When the cylinder is closed, helium cannot be supplied to the device.

T-handle: Used to secure or remove the Helium Cylinder. When securing the cylinder, ensure the T-handle is fully tightened to prevent unnecessary helium leakage. When removing the cylinder, ensure the cylinder knob has been rotated to the closed position.

Helium Cylinder Compartment Door Release: Used to open the Helium Cylinder compartment door. Pressing this switch opens the compartment door.

3.4. Display

The Display mainly consists of a touch screen, control buttons, and a video interface.


Figure 3-19 Display schematic

1) The display is equipped with a locking mechanism that allows it to be secured to the main unit or cart. When locked, it cannot be detached from the cart or main unit. When unlocked, it can be easily removed.
2) Buttons: One display power button, one brightness “+” button, and one brightness “–” button.
 “ ⏻ ” Power button: Turns the display on or off
 “+” button: Increases the display screen brightness.
 “–” button: Decreases the display screen brightness.
3) Brightness adjustment range: 0–100.

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4) Display video interface: Type-C female connector used to connect the display cable.

3.5. IV Pole

The IABP system includes an IV pole, which can be used to conveniently place saline containers and infusion tubing when required.

PRECAUTION:
Do not place liquids on top of the device. Ensure that saline containers and tubing are not suspended directly above the IABP. If spillage occurs, wipe it off immediately and service the device to prevent hazards.

PRECAUTION:
The IV pole shall not be used as a handle or lifting point. During transport or movement of the system, only the trolley handle shall be used.

3.6. System Default Settings

If the system is powered off for less than 3 minutes, user settings and system options are automatically saved in internal memory. If the power-off time exceeds 3 minutes, the system automatically restores the default settings. The default settings are listed below.

User Settings:

User Settings	Default
Operating mode	Automatic
Trigger signal	ECG
ECG lead	II
Internal trigger rate	80
AP Source:	Fiber-optic
Pressure threshold	Default
Pacemaker detection	2
Augmentation ratio	1:1
Augmentation pressure level	10
Inflation timing selection	Automatic
Deflation timing selection	Automatic
Augmentation pressure alarm	OFF(For more information, see Section 3.1.20)

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Reports

User Settings	Default
Waveforms	ECG & AP & BPW
Scheduled reports	OFF
Output format	PDF output

System options

User Settings	Default
R-wave tracking	Enabled
Catheter alarm	Enabled
Leak alarm	Enabled
Date and time	Refer to the most recent setting

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4. Device Operation

Before operating the Lifetech Cardio IntraNova® IABP system, users must be familiar with its controls and functions (see Chapter 3, User Interface) and have a comprehensive understanding of IABP therapy.

4.1. Operating Instructions

4.1.1. Power-On Procedure

1) Connect one end of the AC power cord to the device power inlet and the other end to an AC outlet. Ensure the power cord is not under tension when connected.
2) Plug the power cord into a properly grounded, compatible AC outlet. Do not use adapters that defeat the grounding conductor. If auxiliary equipment is used with the IABP system, ensure proper grounding of that equipment as well.

 WARNING: To reduce the risk of electric shock, this device must be connected to a supply mains with protective earth.
3) Turn on the power inlet switch; the green indicator illuminates steadily. After values appear on the status display, press the main unit power button. A beep indicates that the system has entered the startup sequence.
4) At this time, the IABP system begins subsystem self-tests. Verify that a self-test completion message appears in the information display area.

 If any electrical or pneumatic test fails, an alarm will be displayed on the screen or status display indicating the failed item. Power-cycle the system (power off for at least 10 seconds). If the message reappears, record the alarm name and contact Lifetech Cardio service.

4.1.2. Opening the Helium Cylinder Valve

Check the Helium Cylinder valve. If it is not open, open the valve and confirm the helium pressure. The Helium Cylinder is located at the rear of the unit (marked with a Helium Cylinder icon).

1) Press the Helium Cylinder compartment door release to open the side panel.
2) Turn the Helium Cylinder knob counterclockwise to open the cylinder.
3) Ensure that the helium level indicated on the display or status screen is greater than 1%.

 PRECAUTION: If the Helium Cylinder is properly installed with no leakage and the system will be used shortly, the cylinder valve may remain open. For long-term non-use, close the Helium Cylinder valve.
4) Ensure that no “low helium” message is displayed in the information area. If such a message appears, confirm the cylinder is open, check the helium indicator and pressure gauge, and inspect the cylinder clamp for leakage. Replace the cylinder if necessary. See Section 6.1 for details.

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4.1.3. Operating Modes

Selection of the operating mode determines the IABP functions and user interface layout. For example, in Automatic mode many controls are not required, whereas in Semi-automatic mode these controls are enabled.

Controls	Automatic mode (default)	Automatic mode (User selected)
Trigger signal and lead source	Automatic	User selected
IAB Timing selection	Automatic	User sets initial timing and periodically evaluates and adjusts; Trigger or lead loss
Trigger or lead loss	Automatic detection of a new signal source	IABP alarms and user intervention required

The “Operating Mode” button selects Automatic or Semi-automatic control for patient therapy. Selecting a mode cancels the previous selection. Switching from Automatic to Semi-automatic interrupts assistance and places the system in Standby; press Start to resume. Switching from Semi-automatic to Automatic automatically resumes assistance.

4.1.3.1. Automatic Operating Mode

Press the Automatic button to set the IABP system to Automatic mode. In this mode, operation is fully automated.

- The IABP automatically selects the optimal ECG lead, pressure channel, and trigger source (ECG or pressure). Inflation and deflation timing are automatically set without user intervention. When ECG is selected as the trigger, deflation timing may be fine-tuned using the deflation timing control.
- In Automatic mode, the operator cannot select or deselect the trigger source. If the active trigger is interrupted, the system automatically selects the next best trigger source and retimes if necessary. If a higher-priority trigger becomes available, it is automatically selected and retimed. ECG is the preferred trigger source.
- The IABP software automatically tracks changes in patient heart rate or rhythm and adjusts inflation and deflation accordingly. When an unpredictable rhythm with a valid ECG trigger is detected, the system automatically selects Automatic R-wave deflation, moving the deflation timing indicator to the far right and displaying the corresponding message.

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- When switching from Automatic to Semi-automatic mode, assistance is automatically interrupted. Press Start to resume assistance.

WARNING:

Whenever possible, use ECG or pressure triggering during CPR to synchronize assistance with chest compression rate and rhythm. In Automatic mode, ECG or arterial pressure is automatically selected based on signal quality. If neither provides reliable triggering, use the internal rate trigger by selecting Semi-automatic mode, choosing Internal trigger, and adjusting the rate with the left/right arrow buttons.

4.1.3.2. Semi-automatic Operating Mode

Press the Semi-automatic button to set the IABP system to Semi-automatic mode.

- In this mode, the operator selects the appropriate ECG lead and trigger source, then sets inflation and deflation timing. The software tracks heart rate or rhythm changes and adjusts timing accordingly.
- When the patient rhythm is unpredictable and a valid ECG trigger is available, R-wave deflation is automatically selected and the Automatic R-wave deflation status is displayed.
- In this mode, loss of the trigger source generates a No Trigger Source alarm and audible alert. Trigger source selection is not automatic and must be performed manually.
- Switching from Automatic to Semi-automatic interrupts assistance to emphasize the need for operator intervention. When changing trigger sources, timing must be re-evaluated. Press Start to resume assistance.
- Similarly, if a new trigger source is selected without setting timing, the system immediately interrupts assistance and enters Standby. Press Start to resume assistance.

4.1.4. ECG Acquisition

To ensure optimal triggering, high-quality ECG signals must be used. ECG signals may be obtained directly from electrodes or from compatible external monitoring equipment. See Chapter 7 for external interface requirements.

Direct ECG acquisition requires ECG lead wires and electrodes. To minimize noise, only ECG lead wires provided by Lifetech Cardio may be used.

Electrode type and application technique are the most important factors affecting signal quality. High-quality electrodes are recommended to ensure optimal baseline stability, rapid recovery after defibrillation, and minimal motion artifact.

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To obtain high-quality ECG signals with minimal noise and baseline drift, the following electrode application techniques are recommended:

- Shave the patient's skin or select hair-free electrode sites.
- Properly prepare the electrode sites according to the electrode manufacturer's instructions.

PRECAUTION: This may include thoroughly cleaning the skin with alcohol to remove oils, then drying with a coarse towel, or using an electrode skin-prep agent.

- Use NMPA-approved, commercially available disposable **s**ilver/silver chloride (Ag/AgCl) ECG electrodes, as they can provide high-quality electrode contact immediately upon application.
- If solid-gel (hydrogel) electrodes must be used, apply them to the patient as early as possible before therapy to improve electrode contact.
- Electrodes may dry out when exposed to air for extended periods. Do not remove electrodes from the package too far in advance of use.
- Do not use electrodes that are expired, have damaged packaging, are torn, or have been previously opened.
- For patients with dry skin, poor perfusion, or diabetes, consider using special low-impedance electrodes or electrodes with a larger contact area whenever possible.

4.1.4.1. Use of ECG Electrodes

WARNING: Use only ECG lead wires provided by Lifetech Cardio. Use of any other lead wires may result in abnormal system function.

- When obtaining ECG signals directly from skin electrodes:

a. Ensure the ECG lead cable is securely inserted into the ECG input connector. Connect each ECG lead wire to the corresponding electrode. The table below shows the number of ECG electrodes and leads. To provide trigger lead selection options, five electrodes are recommended.

(#)Electrodes used (AHA standard)	(#) Electrodes used (IEC standard)	(#) Electrodes used (IEC standard)
RA, LA, LL, RL, V	R, L, F, N, C	I, II, III, aVR, aVL, aVF, V

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b. Place the electrodes at the appropriate sites on the patient


Figure 4-1 Electrode placement (AHA standard)


Figure4-2 Electrode placement (IEC standard)

c. Insert the ECG lead cable into the ECG input connector on the panel (Section 3.2).

d. Check the ECG waveform on the display and confirm that the ECG lead field (Item 6, Section 3.1.4) shows the desired lead. Also verify trigger consistency.

- In Automatic mode, the IABP system automatically selects a lead (I, II, III, or External) that provides reliable triggering. If another lead is required:

a. Set the IABP system to Semi-automatic mode, and set the trigger source to ECG in the Trigger menu.

b. In the ECG Signal Source menu, select the required ECG lead based on the connected electrodes (I, II, III, aVR, aVL, aVF, V, or External).

4.1.4.2. ECG Output from External Bedside Monitors

1) When using the high-level ECG output of an external monitor, connect

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the interface cable to the External ECG Input connector (Item 3, Section 3.2).

2) Set the ECG signal source to External in the ECG Signal Source menu (Item 8, Section 3.1.4).
3) Confirm in the ECG lead field (Item 6, Section 3.1.4) that External is being read, and verify that the displayed ECG waveform quality is acceptable. Also check trigger consistency.

4.1.4.3. ECG Troubleshooting

Symptom	Possible Cause	Recommended Corrective Action
ECG noise	Electrode/electrode lead fault.	Check electrode contact; replace electrodes.
ECG baseline noise	Device not configured for the correct mains frequency.	Contact Lifetech Cardio service personnel.
Intermittent ECG	Electrode/electrode lead/patient cable fault.	Check electrode contact; replace electrodes. Check or replace the patient cable.
Lead fault message on display	ECG patient lead fault.	Check electrode contact; replace electrodes. Check or replace the patient cable.
Motion artifact	Electrode fault.	Check electrode contact; replace electrodes.
ECG signal too weak	Poor electrode placement or quality.	Try changing the lead configuration (Lead II typically provides the largest R-wave).
ESU interference	Poor electrode adhesion, improper ECG cable routing, or incorrect ECG lead wires.	Verify correct electrode placement; refer to electrode placement guidelines. In the OR, use OR lead wires equipped with ESU noise filtering.
	Interface to the monitor lacks ESU suppression.	Use direct-contact ECG electrodes and OR lead wires.

4.1.5. Arterial Pressure Acquisition

Arterial pressure signals can be obtained directly or indirectly from the high-level output of a compatible external monitor. See Chapter 7 for connection requirements for external monitor outputs.

Direct arterial acquisition sources may be a fiber-optic IAB catheter or a disposable arterial pressure transducer placed in an external artery (e.g., IAB lumen/radial artery).

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PRECAUTION: When connecting the IAB catheter to the device, ensure the Luer fitting is fully seated and secured. During movement, take care to maintain the connection between the IAB catheter and the device.

4.1.5.1. Direct Arterial Acquisition

Connection requirements between the IAB catheter and the IABP vary depending on the balloon catheter selected.

When using a Lifetech Cardio fiber-optic IAB catheter, insert the catheter connector directly into the IABP system fiber sensor input to transmit arterial signals via an optical connection. Alternatively, a pressure extension line, disposable pressure transducer, and invasive pressure cable may be used to connect to the arterial pressure input; however, in this case the catheter connector must not be connected to the fiber sensor input.

When using a standard (non-fiber) IAB catheter, connect the arterial pressure signal to the IABP system via a pressure extension line, disposable pressure transducer, and invasive pressure cable to the arterial pressure input connector.

PRECAUTION: Before connecting the IAB catheter, read the catheter Instructions for Use. The following steps provide general guidance for connecting a Lifetech Cardio IAB catheter to the IABP system.

1) Fiber-Optic IAB Catheter

PRECAUTION:

Do not touch the exposed end of the fiber-optic IAB catheter cable or allow it to contact any surface. This may damage or contaminate the sensor mating surface.

a. Hold the fiber sensor connector with one hand, ensuring the side with the triangle marking faces upward.

b. Remove the protective cap from the fiber-optic balloon catheter connector. Do not touch the exposed end of the optical sensor cable.

c. With the other hand, remove the protective cover from the fiber sensor port (Item 8, Section 3.2).

d. Insert the fiber-optic IAB catheter connector into the IABP Fiber Sensor Port until it locks into place.

PRECAUTION: Ensure that the triangle marking on the fiber-optic IAB catheter connector is aligned with the panel.

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e. Connect the male Luer fitting of the fiber-optic IAB catheter accessory to the female Luer fitting of the Lifetech Cardio IAB catheter extension line. Then connect the male Luer fitting to the IAB Catheter Port on the IABP system (Item 10, Section 3.2).

PRECAUTION: To ensure reliable operation of the automatic inflation system and correct balloon inflation pressure, the total volume of the balloon membrane, extracorporeal tubing, and catheter extension line must not be altered. Using tubing with a different length or inner diameter from that supplied with the Lifetech Cardio IAB catheter and accessories will change the balloon inflation pressure and may result in automatic inflation failure or assist malfunction. Such situations must be avoided.

PRECAUTION: The pneumatic module supports all Lifetech Cardio adult intra-aortic balloon catheters and accessories.

f. The IABP response differs depending on the current operating mode, as described below:

Automatic Operationg Mode

- If the fiber-optic IAB catheter is the only valid arterial source connected to the IABP system, the pressure signal source in the AP Source menu (Item 16, Section 3.1.6) is automatically set to Fiber. The AP Source status field displays “Fiber” (Item 15, Section 3.1.6), and the fiber sensor is used as the preferred AP Source.
- If another valid arterial source is connected before the fiber-optic IAB catheter is connected, once the IABP system tests and confirms the fiber-optic catheter, the fiber source is automatically given priority as the preferred AP Source.

Semi-automatic Operationg Mode

- Ensure that the AP Source is set to Fiber in the Pressure Signal Source menu (Item 16, Section 3.1.6).
- Once the AP Source is set to Fiber, arterial signals from the fiber-optic IAB catheter are automatically detected and used as the AP Source by the IABP system.
- The AP Source status field displays “Fiber” (Item 15, Section 3.1.6).

g. Zeroing of the fiber-optic pressure sensor is described below:

Initially, the AP Source status field (Item 15, Section 3.1.6) displays Not Zeroed, indicating that sensor zeroing is required.

- Ex-vivo zeroing: Set the fiber zeroing mode to Ex-vivo Zeroing in the System Options submenu, and ensure the fiber-optic balloon is fully exposed to air. Press and hold the Zero button (Item 26, Section 3.1.10) for at least 2 seconds. The automatic zeroing procedure is performed. Upon suc-

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- In-vivo zeroing: Set the fiber zeroing mode to In-vivo Zeroing in the System Options submenu after the fiber-optic balloon has been advanced to the target vessel in the patient. Press and hold the Zero button (Item 26, Section 3.1.10) for at least 2 seconds. The automatic zeroing procedure is performed. Upon successful zeroing, the AP Source status field displays Zeroed.

Patient blood pressure values are valid only after correct zeroing of the fiber sensor. If the sensor is not zeroed or zeroing fails, all pressure values in the arterial parameter area will display "---".

When the AP Source is Fiber, the system automatically corrects the fiber pressure values at appropriate times to ensure measurement accuracy.

PRECAUTION: If fiber-optic IAB catheter identification fails while a disposable pressure transducer is connected and selected as the arterial source, a "Fiber Sensor Fault" error will be displayed.

2) Disposable Pressure Transducer

PRECAUTION: Before setting up and flushing the disposable pressure transducer, read the Instructions for Use supplied with the IAB catheter. The following steps provide general guidance for acquiring arterial signals using a disposable pressure transducer.

a. Connect one end of the invasive pressure transducer cable to the disposable pressure transducer, then connect the other end to the Blood Pressure Input on the IABP system (Item 4, Section 3.2). For details on invasive pressure cables, see Section 8.12.

b. Determine the measurement site: the inner lumen of the intra-aortic balloon or another artery (e.g., radial artery). Patients receiving counterpulsation therapy should have arterial pressure monitored.

c. Connect the arterial line (including the transducer) to the arterial access site, then connect the corresponding pressure cable to the transducer and to the IAB catheter.

d. The IABP system response differs according to the current operating mode, as described below:

PRECAUTION: If a fiber-optic IAB catheter is connected to the IABP system, the fiber sensor has priority over the disposable pressure transducer and becomes the direct arterial source. To use a disposable pressure transducer, the fiber-optic IAB catheter connector must be disconnected.

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Automatic Operationg Mode

- If the disposable pressure transducer is the only valid arterial source connected to the IABP system, the AP Source in the Pressure Signal Source menu (Item 16, Section 3.1.6) is automatically set to Hydraulic, and hydraulic pressure is used as the preferred AP Source.
- AP Source priority is: Fiber > Hydraulic (zeroed) > External > Hydraulic (not zeroed). If an external arterial source is connected before the disposable transducer, the external source remains preferred until the disposable transducer is zeroed. During zeroing, the AP Source must be set to Hydraulic.

Semi-automatic Operationg Mode

- Ensure the pressure signal source is set to the desired source: Hydraulic, Fiber, or External (Item 16, Section 3.1.6).
- Once the AP Source is set to Hydraulic, arterial signals from the disposable pressure transducer are automatically detected and used as the AP Source by the IABP.
- The label “AP” is displayed in the AP Source field above the arterial pressure waveform (Item 15, Section 3.1.6).

e. Zeroing of the disposable pressure transducer is described below:

Initially, the AP Source field on the arterial waveform (Item 15, Section 3.1.6) displays Not Zeroed, indicating that sensor zeroing is required.

- Open the transducer to atmospheric pressure.
- Press and hold the Zero button for at least 2 seconds to initiate the automatic zeroing procedure. (The zeroing function can compensate the disposable pressure transducer with an offset range up to ±150 mmHg.) After successful zeroing, all numeric pressure values on the right side of the display show zero (±2 mmHg), and the AP Source status field displays Zeroed.

 Patient blood pressure parameters are valid only after correct zeroing of the disposable pressure transducer. If zeroing fails, all values in the arterial parameter area display “---”. If the transducer offset exceeds ±150 mmHg (as may occur with sensor failure), “---” is also displayed.

PRECAUTION: If the disposable pressure transducer offset exceeds ±150 mmHg, automatic zeroing is not possible. The Not Zeroed message indicates a fault condition. In this case, the disposable pressure transducer may be incompatible or defective and must not be used.

To prevent zeroing at inappropriate times, the pressure zero button is disabled under the following conditions:

- Pulsations detected on the arterial waveform

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- Disposable pressure transducer not connected
- AP Source set to External

f. Close the transducer to atmosphere and verify that the arterial waveform is displayed and that systolic, diastolic, and mean pressures appear in the arterial parameter area.

4.1.5.2. Pressure Output from External Monitors

This brief procedure assumes that the monitoring line has been established and the external monitor has already been zeroed.

1) Connect the interface cable from the external monitor to the External Blood Pressure Input on the IABP system (Section 3.2).
2) The IABP system response differs depending on the current operating mode, as described below.

Automatic Operationg Mode

- If the external bedside monitor is the only valid arterial source connected to the IABP system, the AP Source in the Signal Source menu (Item 15, Section 3.1.6) is automatically set to External, and the external bedside monitor is used as the preferred AP Source.
- If a valid fiber or hydraulic (zeroed) AP Source is connected to the IABP system before the external monitor is connected, the fiber or hydraulic (zeroed) source remains the preferred AP Source. If the external source is required, disconnect the fiber or hydraulic (zeroed) source. The user may also manually switch to the external source in the Signal Source menu.

Semi-automatic Operationg Mode

- Set the AP Source to External in the Signal Source menu (Item 16, Section 3.1.6).
- Once the AP Source is set to External, arterial signals from the external bedside monitor are automatically detected and used as the AP Source by the IABP system.

3) Verify that the displayed pressure waveform quality is acceptable.

4.1.5.3. Pressure Source Troubleshooting

1) Fiber-Optic IAB catherter

Symptom	Possible Cause	Recommended Corrective Action
Fiber cable connected but no signal displayed	Loose connection between the fiber cable and the fiber sensor port	Reinsert the sensor connector into the fiber sensor port and ensure it is fully seated.

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Pressure readings from the fiber IAB differ from those measured by the disposable pressure transducer.	Different pressure sources are measuring pressure at different anatomical locations.	None required.
	Environmental conditions may have changed since the last calibration.	If applicable, re-zero other pressure sources. While the IABP system is providing assist, press the Zero button to recalibrate the fiber pressure (the fiber zeroing mode must be set to In-vivo Zeroing).
Fiber cable connected but no signal displayed	The fiber pressure sensor has not been calibrated, preventing accurate display.	If the fiber sensor cable was just inserted, the pump will automatically calibrate after assist begins. If initial zeroing fails, press the Zero button while the IABP system is providing assist to recalibrate the fiber pressure (fiber zeroing mode must be In-vivo Zeroing).

2) Disposable Pressure Transducer

Symptom	Possible Cause	Recommended Corrective Action
The waveform from the disposable pressure transducer is not available to the pump.	The disposable pressure transducer input is not activated while the fiber IAB sensor cable is connected.	Remove the fiber sensor cable from the pump.
Continuous low-amplitude, high-frequency noise appears on the arterial	The pressure monitoring tubing may be in contact with the	Reroute the pressure tubing to avoid direct contact with

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waveform, affecting automatic operation.	Main Console, causing resonance with the pump compressor.	the main unit. Consider mounting the transducer and flush bag on a separate IV pole.
Excessive attenuation of the arterial waveform or low pulse amplitude. Possible cause: Excessive air bubbles or excessive bends in the pressure tubing.	Excessive air bubbles or excessive bends in the pressure tubing.	Flush the pressure line to remove air bubbles. Use tubing no longer than 8 feet (5 feet preferred) and maintain 300 mmHg pressure in the flush bag.

4.1.6. Trigger Source Selection

In Automatic Opertating Mode, the trigger source is selected automatically. In Semi-automatic mode, the trigger source must be selected in the Trigger Signal menu. If the trigger source is lost or invalid, a No Trigger alarm message and audible alert are generated. For details, see Section 3.1.13 Trigger Source Selection.

PRECAUTION: In Semi-automatic mode, switching from one trigger source to another automatically interrupts assist. This design alerts the operator to reassess timing selection and, if necessary, adjust inflation and deflation timing before resuming assist. Press the Start button to resume assist.

PRECAUTION: Trigger Source: Inappropriate trigger source selection may result in suboptimal therapy or patient injury.。

4.1.6.1. ECG Triggering

When ECG is used as the trigger signal, the IABP system triggers on the patient's R wave. Each detected trigger event causes the heart icon in the trigger heart rate parameter area to flash on and off (Item 13, Section 3.1.5). The ECG signal may originate from patient electrodes or an external patient monitor.

- When ECG is used as the trigger source and pacemaker pulses are detected, these pulses are ignored and the IABP system continues to trigger on the patient's R wave. Detected pacemaker pulses are enhanced and displayed as pacemaker spikes on the ECG waveform (Item 9, Section 3.1.4). If necessary, pacemaker detection sensitivity

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may be increased or decreased by adjusting the pacemaker detection level in the Threshold menu (see Section 3.1.24.4).

- In Automatic mode, the system automatically selects the optimal ECG signal source (I, II, III, or External). If these trigger sources become invalid, pressure may be automatically selected as the trigger source. Electrosurgical noise, motion artifact, or electrode loss may cause automatic trigger source switching.
- In Semi-automatic mode, use the "<<<" and ">>>" buttons to scroll left or right to change the ECG signal source (I, II, III, aVR, aVL, aVF, V, or External) in the Signal Source menu (Item 8, Section 3.1.4).
- In all cases, the ECG signal source field (Item 6, Section 3.1.4) displays the currently active ECG lead.

4.1.6.2. Pressure Triggering

When pressure is used as the trigger source, the IABP system triggers on the patient's arterial blood pressure.

PRECAUTION:

When the ECG signal is of good quality, this is generally a secondary (less preferred) trigger source. Whenever possible, use ECG triggering.

The IABP system triggers on the upstroke of systole on the patient's arterial waveform. The signal source can be Fiber, Hydraulic, or External.

- When Hydraulic is used as the pressure trigger source, assist may be started with the transducer not zeroed to enable rapid initiation. However, zeroing is required to display numeric arterial parameters.
- In Automatic mode, the system selects the best available arterial signal source. Fiber is preferred. If fiber is not connected, Hydraulic is preferred. If neither is available, External arterial input is selected automatically. If all valid signal sources are lost, a No Trigger alarm and audible alert are generated. Electrosurgical noise, motion artifact, or electrode loss may also cause automatic trigger source switching.
- In Automatic mode, the signal source can be selected manually. If the operator-selected pressure source becomes unavailable or its triggering is lost, the IABP system automatically selects another pressure source, if available.
- In Semi-automatic mode, use the Left and Right buttons to select the pressure signal source for pressure triggering (Fiber/Hydraulic or External) (Item 16, Section 3.1.6).

1) Pressure Threshold

- In Automatic Operating mode, the trigger source is selected automatically. In Semi-automatic mode, the trigger source must be selected in the Trigger Signal menu. If the trigger source is lost or invalid, a No Trigger alarm message and audible alert are generated. For

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details, see Section 3.1.13 Trigger Source Selection.

- In certain clinical settings or environments, fixed sensitivity may be advantageous. In Semi-automatic mode, when pressure is used as the trigger source, the Pressure Threshold Up "▲" or Down "▼" buttons allow increasing or decreasing the pressure threshold from 7 to 30 mmHg in 1 mmHg increments. The current setting is displayed between the Up and Down buttons (see Item 35, Section 3.1.19).

WARNING:

If the manual pressure trigger threshold option is used, adjust the threshold whenever the systolic pulse amplitude changes continuously. Changes in patient condition or initial zeroing of the fiber IAB catheter by the IABP system can alter the amplitude. After any threshold adjustment, inflation and deflation timing must be re-evaluated.

2) Using Pressure Triggering

The IABP system can trigger on the rising timing of systole and ignore arterial pressure increases caused by balloon counterpulsation. This protective mechanism helps prevent false triggers, but it may mask premature pulses associated with a sudden increase in heart rate. The system can detect and adapt to such increases. However, very large changes (e.g., doubling of heart rate) may still cause failure to detect every systole. If the displayed heart rate is half of the patient's actual heart rate, the trigger mechanism must be corrected. Although rare, if this occurs, press and hold the Start button to correct it. This action interrupts heart rate assist and restores correct trigger detection. If no action is taken, correct triggering will be restored at the next 60-second automatic synchronization check.

(1) Automatic Synchronization Check During Pressure Triggering

An important feature of the IABP system is the 60-second periodic synchronization check during pressure triggering. When the IAB ratio is 1:1, the pump pauses for one cardiac cycle per minute. The patient's natural cardiac cycle and interval parameters, unaffected by counterpulsation, are displayed for evaluation. This periodic check ensures correct operation and timing interval determination, and the intervals may be corrected if necessary. During assist, the user may press and hold the Start button at any time to initiate an automatic synchronization check. If rare loss-of-synchronization events are observed, this method can immediately re-synchronize triggering and timing selection. When the check is initiated by pressing and holding the Start button, the timer for the next check is automatically reset to 60 seconds.

(2) Arrhythmias During Pressure Triggering

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The IABP system can automatically adapt to persistent random arrhythmias, such as atrial fibrillation. Such rhythms have no predictable pattern and may cause early systolic ejection that can unavoidably conflict with balloon deflation, reducing the stability and consistency of systolic pulse detection. The system automatically detects such rhythms and minimizes the likelihood of conflicts.

In Semi-automatic Operating Mode, the information display area shows “Irregular Pressure Triggering”. This message indicates that the system will automatically compensate by advancing deflation to avoid interfering with systolic ejection. Therefore, the user must not adjust deflation timing control. Once the patient returns to normal rhythm and the system automatically returns to standard timing selection, adjusting deflation timing may impair triggering performance.

WARNING:

When pressure is used as the trigger source, balloon deflation must be continuously optimized so that deflation is completed before systole. Late deflation reduces systolic pulse pressure and delays detection. To achieve consistent and reliable pressure triggering, the system relies on a distinct and timely systolic pulse. If balloon deflation conflicts with systolic ejection while using pressure triggering, triggering may be delayed or lost, potentially resulting in loss of synchronization.

WARNING:

Pressure triggering is not recommended in persistent irregular heart rates or tachyarrhythmias. When pressure triggering is used, promptly and sufficiently adjust deflation so that deflation is completed before systole, and continuously monitor the patient. If “Irregular Pressure Triggering” is displayed, do not attempt to adjust deflation control; the system will automatically compensate by advancing deflation to avoid interfering with systolic ejection.

4.1.6.3. V/A-V Pacemaker Triggering

When V/A-V pacemaker triggering is used as the trigger source, the IABP system automatically determines whether a corresponding ventricular or A-V sequential pacemaker is present. In either case, the system triggers on the pacemaker’s ventricular beat. This ignores the patient’s QRS and arterial pulse. If applicable, the ECG (Section 3.1.4) displays ventricular or A-V information accordingly. V/A-V pacemaker triggering is available only in Semi-automatic mode.

- When V/A-V pacemaker triggering is used, the patient must have 100% ventricular pacing or 100% A-V sequential pacing that is detectable (i.e., no demand pacing). The ECG signal may come from

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patient skin electrodes or an external ECG source.

- When the patient is paced and the resulting QRS response is too weak to serve as a reliable trigger source (e.g., in the operating room), V/A-V pacemaker triggering is commonly selected.

- V Pacemaker Identification Criteria:
- If the ventricular pacing interval is fixed and the heart rate is below 180 bpm, the system identifies it as a V pacemaker.
- A-V Pacemaker Identification Criteria:

 If the atrial-ventricular interval is between 50–250 ms, the system identifies it as an A-V pacemaker. However, to correctly identify the pacemaker at higher pacing rates, the atrial-ventricular interval must shorten.

 When a ventricular/A-V sequential pacemaker is used as the trigger source, ensure that the IABP system can reliably detect the pacemaker spikes. If a ventricular pacemaker is used, verify that an enhanced ventricular pacing pulse is visible in every cardiac cycle.

 If an A-V sequential pacemaker (A-V pacemaker) is used, verify that two enhanced pacing pulses are visible in every cardiac cycle.

 If atrial or ventricular pacing pulses are not detected, or if the detected signals are intermittent and may affect heart rate calculation, consider adjusting the pacemaker detection level (Section 3.1.24.4).

 If 100% paced rhythm is present (i.e., no demand pacing), the system selectively triggers on the ventricular pulse of any pacemaker type. If the criteria for any pacemaker type are not met, a trigger alarm is generated and a message to check pacemaker timing selection is displayed.

PRECAUTION:

V/A-V pacemaker triggering does not function in the presence of ESU interference. Assist is temporarily interrupted and automatically resumes when interference decreases. When V/A-V pacemaker triggering is used, detection of ESU interference generates a trigger interference alarm message.

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WARNING:

For patients with pacemakers, the heart rate meter may count pacing pulses during asystole or arrhythmias. Do not rely solely on heart rate meter alarms. Closely monitor patients with pacemakers. For the system's capability to suppress pacing pulses, see Section 8.1 of this manual.

4.1.6.4. A Pacemaker Triggering

When A pacemaker (atrial pacemaker) triggering is used, the patient's QRS is actually used as the trigger event. The IABP system triggers on the patient's R wave and suppresses atrial pacemaker artifact. The ECG signal may originate from patient electrodes or an external ECG source. This mode is recommended only when atrial pacemaker trailing artifact interferes with R-wave detection. Fixed-rate or demand atrial pacemakers may be used in this mode. A pacemaker triggering is available only in Semi-automatic mode.

4.1.6.5. Internal Triggering

Internal triggering may be used when no mechanical cardiac cycle is present (e.g., during cardiopulmonary bypass or cardiac arrest). When Internal is used as the trigger source, the IABP triggers at the fixed rate set by the internal trigger rate (see Section 3.1.13). The Left and Right buttons for internal trigger rate control increase or decrease the rate from 40 to 120 bpm in 5 bpm increments. The current setting is displayed between the Left and Right buttons. Internal triggering is available only in Semi-automatic mode.

- This is a typical mode used in the operating room during coronary artery bypass grafting (CABG). It is used when the patient is on cardiopulmonary bypass and no QRS waveform is present. During full cardiopulmonary bypass, this assist state is recommended to reduce the risk of thrombus formation on the balloon membrane.
- When the internal trigger source is used, the system continues to monitor R-wave activity via the ECG patient cable. If a usable QRS is present, internal triggering must not be used. In this mode, if a valid R wave is detected, the IAB will automatically deflate (to prevent asynchronous assist), and an alarm message indicating that ECG was detected while using internal triggering will be displayed, accompanied by an audible alert. If reliable R-wave activity returns, the system should switch back to ECG triggering to achieve correct timing selection.

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Do not remain in Internal Trigger mode when the patient has cardiac ejection.

PRECAUTION:

The internal trigger rate can be changed only when the trigger source is set to Internal.

4.1.6.6. Trigger Troubleshooting

Symptom	Possible Cause	Recommended Corrective Action
No triggering occurs	In Automatic Operating Mode, valid pressure and ECG trigger sources are not present or are lost.	Check all pressure lines and verify that the photoelectric (fiber-optic) cable connections are intact. If a disposable pressure transducer is in use, confirm that it is not open to atmosphere. Connect or reapply electrodes and check that all ECG cables/electrode connections are secure.
	The pulse pressure is unsuitable for pressure triggering	If a disposable pressure transducer is being used, flush the fluid line to attempt to restore pressure.
Abnormal triggering	The paced ECG acquired from an external monitor has insufficient bandwidth.	Refer to Chapter 7 on external monitor connection.
	Failure to correctly recognize pacemaker pulses.	Adjust pacemaker detection according to Section 3.1.24.4.
	During ECG triggering, the atrial pacemaker trailing artifact is too large.	Select A pacemaker triggering.
	Demand pacing in V/A-V pacemaker mode.	Select ECG or pressure triggering.
	Pacemaker spike overlaps the R wave.	Consider switching to pressure triggering mode (see Section 3.1.13).
Under pressure triggering, the system triggers every other cardiac cycle.	Pressure triggering requires resynchronization.	Press and hold the Start button to resynchronize.

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In Semi-automatic mode, the system skips assist cycles during pressure triggering.	Deflation timing is set too late.	If necessary, readjust deflation timing selection and press and hold the Start button to resynchronize.

4.1.7. Selecting the IAB Assist Ratio

When evaluating timing selection or reducing assist, it is often necessary to change the IAB assist ratio (i.e., to suspend IAB treatment for some beats). Selecting the IAB assist ratio determines the heart rate ratio for IAB assistance. Options are: **1:1** (assist every beat), **1:2** (assist every two beats), or **1:3** (assist every three beats). To select the desired IAB assist ratio, press the corresponding button in the IAB assist ratio menu (see Section 3.1.17).

The selected IAB assist ratio determines how arterial values are calculated and displayed. In 1:1 mode, all beats are assisted and each pressure reading displays an assisted value.

When 1:2 or 1:3 is selected, assisted and unassisted systolic and diastolic pressures are displayed numerically.

PRECAUTION:

Improper timing settings may fail to achieve the expected therapeutic effect or may cause patient injury.

4.1.8. Timing Selection

Timing selection means selecting the inflation and deflation points on the arterial waveform.

In Automatic Operating Mode, the IABP automatically initializes and continuously adjusts timing selection. In Semi-automatic Operating Mode, the user must initialize and periodically check IAB timing selection.

In Standby mode, the inflation time interval is highlighted on the arterial waveform. The highlighted inflation interval is an approximate indication of balloon inflation time.

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
Figure 4-3：Example of a highlighted inflation interval starting at the dicrotic notch

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

4.1.8.1. Inflation

- Occurs at the dicrotic notch.
- Displayed as a “V” shape.
- Ideally, augmented pressure exceeds systolic pressure.

4.1.8.2. Deflation

- Occurs immediately before the next systolic event.
- Results in reduced assisted end-diastolic pressure.
- Results in reduced assisted systolic pressure.


Figure 4-4：Example waveform showing an IABP frequency of 1:2

The examples provided are for illustrative purposes only and do not represent actual clinical conditions.

4.1.9. Starting Assist

1) Press the Start button and confirm that automatic inflation is in progress. The IAB Fill button is highlighted to indicate inflation, and the balloon catheter circuit will be evacuated and filled with an appropriate amount of helium.

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- Assist begins after automatic inflation completes successfully. The information display area shows an “Inflation Complete” message, the IAB Fill button is no longer highlighted, and the Start button is highlighted to indicate assist is active.
- In Automatic mode, the pump sets augmented pressure to the maximum and starts running, and timing selection is determined automatically.
- Balloon timing selection may be fine-tuned to achieve optimal counterpulsation.

PRECAUTION:

After power-on self-test is completed, if the system remains idle for more than 5 minutes, it enters power-saving mode.

PRECAUTION:

Press the Up/Down buttons in the augmented pressure level menu (Section 3.1.18) to increase or decrease augmented pressure step by step.

2) The augmented-pressure alarm limit is displayed between the Up and Down buttons in the augmented-pressure alarm control area (see Section 3.1.20). About 3 minutes after starting assist, verify that the alarm setting is approximately 10 mmHg lower than the patient’s augmented pressure. If a disposable pressure transducer or fiber sensor is used, zeroing must be performed in advance for automatic alarm setting.
3) If adjustment of the alarm limit is required, open the Augmented Pressure Alarm threshold setting menu. Then press the Up “ ▲ ” or Down “ ▼” button to adjust the alarm limit.

PRECAUTION:

When used correctly, the augmented-pressure alarm can serve as an important backup for internal monitoring alarms. Due to the dynamic nature of the gas delivery system, leak- and catheter-related alarm operations must not be performed when the patient is critically ill (see Section 4.4 Alarms and Messages for details). By setting the augmented-pressure alarm, the system can monitor the level of assist and warn the user if augmented pressure decreases.

4.1.10. IAB Fill

To evacuate and reinflate the balloon catheter circuit with helium every 2 hours, the IABP automatically pauses assist. Assist automatically resumes when the automatic inflation cycle is completed. By pressing and holding the “IAB Fill” button, the operator can initiate IAB inflation at any time. This action resets the 2-hour automatic inflation timer. After the 2-hour automatic inflation timeout (when in Standby mode), automatic inflation will be executed one minute after returning to Assist mode.

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During air transport, automatic inflation is performed automatically if the local atmospheric pressure decreases or increases by 25 or 50 mmHg, respectively. Such pressure changes may occur when altitude changes by approximately more than 300 m; in these cases, automatic inflation is triggered to ensure the balloon pressure is appropriate for local conditions.

PRECAUTION:

After assist is started, the device initiates automatic **"**IAB Fill**"** every 2 hours (automatic evacuation and replacement of the gas in the catheter).

To maintain optimal augmented pressure, the IABP system replenishes helium lost due to diffusion one hour after each automatic inflation cycle.

PRECAUTION:

If the automatic inflation program cannot properly evacuate and inflate the pneumatic module, an Automatic Inflation Failure message will be displayed and an alarm tone will sound. Tap the alarm to view corrective actions (for help screen information, see Section 4.3).

4.1.11. Weaning the Patient from IABP Support

Weaning can be achieved by gradually and continuously reducing the IAB frequency or the IAB volume displacement (augmented pressure), or by using both methods.

During weaning, it is recommended to frequently assess hemodynamic parameters and the patient's condition.

WARNING:

When weaning by reducing augmented pressure, do not reduce the augmented pressure to a level below Level 5.

4.1.12. Safety-Related Operations

Procedures for safely terminating operation of the IABP system

1) Confirm that the device is not providing assist.
2) Open the Settings menu, tap "Power Off", then tap "Confirm" to shut down; or press and hold the main unit power button " ⏻ " for 5 s to shut down.
3) Turn off the power inlet switch so that it is in the OFF position.

To disconnect from mains power:

Turn off the power inlet switch, or unplug the power cord.

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4.2. Pneumatic Module Leak Test

This test is used to verify the integrity of the pneumatic module (see Section 3.1.24.6). It is recommended to perform this test before and after each use. This test requires operator interaction. In particular, during the test, the operator will be instructed to occlude the pneumatic module's IAB catheter port. The IABP system will prompt the user to perform these actions at the appropriate time.

WARNING:

Do not perform the pneumatic module leak test when the IAB catheter inside the patient is connected to the pump.

In the Settings menu, open the About submenu and select Pneumatic Test. On the test screen, press Start Test and follow the prompts. When "Please block the IAB catheter port with a Luer plug" appears, firmly insert a Luer cap into the IAB catheter port. Then tap Blocked. When the test is complete, the leak-test result will be displayed. If the system passes, a green Test Passed message will appear; otherwise, failure information will be displayed.

If you have any concerns about the integrity of the pneumatic module, please contact Lifetech Cardio service/support.

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4.3. Help Screen

The Help screen is used to view alarm/message descriptions and alarm configuration. It is context-sensitive and provides help based on what is shown on the display. For example, help for an alarm/message is available only when an alarm condition or message exists. When one or more alarms/messages appear, the user may tap an alarm to enter its Help screen, which provides guidance for the current alarm/message.

Browsing the Help screen:

1) The information display area shows all current alarms/messages in a predefined priority order, with higher-priority items displayed above lower-priority items. The priority order is: high-priority alarms, medium-priority alarms, low-priority alarms, and messages.
2) If there is only one alarm or message, a single item is displayed in the information area. If there are multiple alarms or messages, multiple items are displayed. Tap a specific alarm/message item to show the corresponding Help screen.
3) Useful help information can be viewed on the Help screen.
4) Tap a blank area of the screen to close the Help screen.

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4.4. Alarm and Message Examples

The IABP uses audible and visual alarm signals to require immediate operator response, prompt operator response, or draw the operator's attention. Audible alarm signals (sound sequences) are primarily used to attract attention. These signals are carefully designed to convey the urgency of response and identify the pump. The sound can effectively cover a relatively large area around the operator. The inherent alarm latency (time from alarm trigger to alarm display) is within 200 ms.

When an alarm occurs, the background of the information display area flashes in a corresponding color as an indicator. High-priority alarms cause a red background to flash at 1.7–2.3 Hz; medium-priority alarms cause a yellow background to flash at 0.4–0.6 Hz; low-priority alarms cause a yellow background to remain steadily lit.

Press the alarm icon " " in the upper-left of the screen, or touch the flashing alarm lamp in the upper-right of the display with your palm, to temporarily silence the active alarm sound for 60 seconds, and the icon will change to the Alarm Sound Paused icon "". This control does not override the alarm. If the alarm condition is not resolved within 60 seconds, the alarm sound will be re-enabled. While audio is silenced, the alarm message remains displayed with its active visual indications. If a new alarm occurs during the 60-second silence period, the alarm sound is immediately re-enabled. In addition, if the "" alarm sound pause icon is pressed during the 60-second silence period, sound will resume immediately (see Section 3.1.1).

In this alarm system, pressure-related alarm states are generated based on time and pressure.

When two or more alarm states of the same priority generate alarm signals (including but not limited to internal ordering and factors affecting alarm generation), there will be only one alarm sound for that priority level.

This alarm system does not change the preassigned priority or relative priority order (e.g., escalation or de-escalation) of a specific alarm state.

Because the device assists the heart, certain alarm delays are related to the cardiac cycle; alarm delay may vary as heart rate changes.

Adjusting the volume changes the volume of existing alarm tones.

When alarms of multiple priority levels exist, only the highest-priority alarm sound is emitted to avoid nuisance.

If multiple alarms are displayed, high-priority alarms appear at the top, medium-priority above low-priority alarms and messages, and low-priority alarms above messages. If multiple alarms share the same priority, they are ordered by occurrence time; the more recent the alarm, the higher it appears.

When multiple alarms of the same priority occur simultaneously, the audio/visual

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signals are the same as for a single alarm at that priority level. When alarms of different priorities occur simultaneously, the highest-priority alarm is displayed.

Alarm performance is the same in all use scenarios.

If the system is powered off for less than 3 minutes, alarm settings are automatically saved in internal memory. If the power-off time exceeds 3 minutes, default alarm settings are automatically restored (see Section 3.6).

When an alarm is turned off or paused, a prompt indicator remains displayed and cannot be dismissed, but no alarm sound is emitted.

PRECAUTION:

If the total number of alarms/messages in the information display area exceeds two, the later items will be collapsed. Press the Down button in the information display area to view collapsed alarms/messages. When expanded, up to four alarms/messages can be shown; press the Down button again to switch to the next page for additional items.

PRECAUTION:

When an alarm occurs, the operator should be positioned facing the display so the visual alarm signals can be seen clearly, enabling appropriate actions based on the alarm information.

WARNING:

Do not set the alarm volume too low; otherwise, alarms may not be heard above background noise where the IABP is used.

Visual alarm signals include icons and text. Displayed information further emphasizes urgency (color and icon shape), while specific text messages identify the alarm. At a distance of 1 meter from the screen, the text message can be clearly observed from the operator's position. In addition, context-related help information may be displayed to supplement the alarm, aiding understanding and guiding resolution.

For ease of operator attention and understanding, alarms and messages are classified as high-priority, medium-priority, low-priority alarms, and messages. They are displayed according to the priority of the condition indicated.

High-priority alarms are displayed first at the top of the information display area. These alarms require immediate operator response. High-priority alarms emit a high-priority alarm tone and, in most cases, interrupt counterpulsation. If multiple alarm conditions occur, all current alarms/messages are displayed in a predefined

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priority order, with high-priority above low-priority. The Help screen area is arranged to reflect the priority order.

When a medium-priority alarm is displayed, IAB assist is not interrupted depending on the device status, and a medium-priority alarm tone is emitted. Medium-priority alarms require immediate operator response.

When a low-priority alarm is displayed, assist is not interrupted and a low-priority alarm tone is emitted. Low-priority alarms require the operator's attention.

Messages usually call for attention to a condition but do not require immediate action and may persist for a longer time.

WARNING:

The augmented-pressure alarm that is automatically set at power-on supports IAB alarms (leak and balloon catheter alarms) at higher heart rates. Therefore, do not manually disable this alarm.

WARNING:

If multiple pumps are in use nearby, confirm the corresponding visual indications to ensure the source of an alarm sound is correctly identified.

4.4.1. High-Priority Alarms

High-priority alarms indicate conditions requiring immediate operator response. In most high-priority alarms, assist will be interrupted. A red alarm icon with three exclamation marks appears in the information display area, indicating a high-priority alarm. All high-priority alarms share the same audible tone: a sequence of five notes that repeats in a loop.

Figure 4-5: High-Priority Alarms Icon

Insufficient Negative Pressure in the Vacuum Reservoir During Inflation
A component fault may exist in the inflation system.
1) Re-inflate the IAB. 2) If the alarm message persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

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Operating modes: All

Trigger sources: All

Cause: Insufficient negative pressure in the vacuum reservoir while evacuating gas from the helium circuit.

System response: Deflation / IAB deflation.

Alarm reset: Press “Start” button, or press and hold “IAB Fill” button for 2 seconds to retry automatic inflation.

Insufficient Drive Pressure During Inflation
A component fault may exist in the inflation system.
1) Re-inflate the IAB. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Insufficient drive pressure while evacuating gas from the helium circuit.

System response: Deflation / IAB deflation

Alarm reset: Press “Start” button, or press and hold “IAB Fill” button for 2 seconds to retry automatic inflation.

Helium Evacuation Timeout During Inflation
IAB catheter is loose.
1) Confirm the IAB extension tubing size is appropriate and all connections between the IAB catheter and IABP are secure, with tubing unobstructed. 2) Re-inflate the IAB. 3) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.
A component fault may exist in the inflation system.
1) Re-inflate the IAB. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

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Cause: Helium evacuation timeout during inflation.

System response: Deflation / IAB deflation.

Alarm reset: Press “Start” button, or press and hold “IAB Fill” button for 2 seconds to retry automatic inflation.

Inflation Timeout
IAB catheter is loose.
1) Confirm the IAB extension tubing size is appropriate and all connections between the IAB catheter and IABP are secure, with tubing unobstructed. 2) Re-inflate the IAB. 3) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.
The Helium Cylinder is not providing sufficient helium.
1) Confirm there is sufficient helium in the cylinder and the cylinder valve is open. 2) Re-inflate the IAB. 3) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.
A component fault may exist in the inflation system.
1) Re-inflate the IAB. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

Alarm attributes: Other alarm; Interlock alarm.

Operating modes: All

Trigger sources: All

Cause: Inflation timeout.

System response: Deflation / IAB deflation.

Alarm reset: Press “Start” button, or press and hold “IAB Fill” button for 2 seconds to retry automatic inflation.

Isolation Valve 1 Failure During Inflation.
A component fault may exist in the inflation system.
1) Re-inflate the IAB. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

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Operating modes: All

Trigger sources: All

Cause: Isolation Valve 1 failure during inflation.

System response: Deflation / IAB deflation.

Alarm reset: Press “Start” button, or press and hold “IAB Fill” button for 2 seconds to retry automatic inflation.

Inflation Completion Timeout
IAB catheter is loose.
1) Confirm the IAB extension tubing size is appropriate and all connections between the IAB catheter and IABP are secure, with tubing unobstructed. 2) Re-inflate the IAB. 3) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.
The Helium Cylinder is not providing sufficient helium.
1) Confirm there is sufficient helium in the cylinder and the cylinder valve is open. 2) Re-inflate the IAB. 3) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.
A component fault may exist in the inflation system.
1) Re-inflate the IAB. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Inflation completion timeout.

System response: Deflation / IAB deflation.

Alarm reset: Press “Start” button, or press and hold “IAB Fill” button for 2 seconds to retry automatic inflation.

Vacuum Valve not Actuating
Vacuum valve actuation failure during assist mode.

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1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: The system detects that after the positive-pressure valve has been opened for a period of time, the vacuum valve is not energized.

System response: The system is disabled and the pneumatic module vents to atmosphere / IAB deflation. Physiological signals remain displayed. Press “Start” button, or press and hold “IAB Fill” button to re-enter self-test; assist can resume after a successful self-test.

Alarm reset: Attempt to resolve by pressing “Start”.

IAB Disconnection
Disconnection of the IAB catheter or IAB extension tubing detected during assist mode.
1) Reconnect the IAB catheter and IAB extension tubing. 2) Press “Start” button to re-inflate the IAB catheter and resume augmentation.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: The IAB catheter or IAB extension tubing became disconnected during assist mode.

System response: Deflation / IAB deflation.

Alarm reset: After reconnecting the IAB catheter, press “Start” button, or press and hold “IAB Fill” button for 2 seconds to initiate automatic inflation.

IAB Circuit Leak
Gas leakage or rapid loss detected in the IAB circuit during assist mode.
1) Check the tubing for signs of blood. If blood is detected, immediately stop augmentation and notify the physician; remove the IAB catheter per the manufacturer’s instructions. 2) If no blood is detected, confirm all tubing connections are secure. Press “IAB Fill” button for 2 seconds to perform automatic inflation, then press “Start” button to resume augmentation.

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3) If the patient has fever or tachycardia, consider initiating automatic inflation before the routine 2-hour cycle to increase inflation frequency. 4) If the alarm persists and no blood is present in the IAB circuit, consider entering “Settings → System Options” and setting “Leak Alarm to Off”.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Cumulative transfer gas loss exceeds the nominal dynamic limit of 5 cc/hour. This alarm is activated only when the IAB Fill duration is ≥ 80 ms and the deflation duration is ≥ 250 ms.

System response: Deflation / IAB deflation.

Alarm reset: Press and hold “IAB Fill” button for 2 seconds; the message and audible alert will clear.

Air Ingress into IAB Circuit.
Air entry into the IAB circuit detected during assist mode.
1) Verify all connections are leak-free. 2) Press “Start” button to perform automatic inflation and resume augmentation. 3) If the alarm persists, contact Lifetech Cardio service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Cumulative transfer gas relative to the most recent automatic inflation volume exceeds 5 cc. This alarm is activated when the IAB inflation duration is ≥ 80 ms and the deflation duration is ≥ 250 ms.

System response: Deflation / IAB deflation.

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ECG Detected
A valid ECG detected while in internal trigger mode.
1) Select a patient signal trigger source (e.g., ECG or pressure) via the Trigger Signal menu. 2) Verify timing. Press “Start” button to resume augmentation.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: Semi-automatic

Trigger source: Internal

Cause: ECG activity detected continuously for 4–6 seconds while in internal trigger mode.

System response: Assist continues, but the device immediately deflates upon detection of an R-wave.

Alarm reset: Automatically resets when ECG activity ceases, or manually exit internal trigger mode.

Balloon Pressure Too High
Restriction of the IAB catheter or tubing detected during assist mode.
1) Check the IAB catheter tubing, extracorporeal tubing, and IAB extension tubing for restrictions and remove any obstruction. 2) Press “Start” button to resume augmentation.
IAB membrane not fully deployed.
1) Aspirate to ensure no blood reflux into the extracorporeal tubing. 2) Using a syringe, manually inflate and deflate 30 cc via the IAB connector. 3) Press “Start” button to re-inflate the IAB and resume augmentation.
Part of the IAB remains within the sheath after insertion.
1) Check the markings on the IAB catheter to confirm the balloon is fully advanced beyond the sheath. If not, reposition the sheath and catheter per the manufacturer’s instructions. 2) Press “Start” button to resume augmentation.
Safety diaphragm leakage (cross-leak).
1) Replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

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Trigger sources: All

Cause: Refer to the possible causes listed above.

System response: Deflation / IAB deflation.

Alarm reset: Press “Start” button to resume assist and reset.

Drive Pressure Too High
Possible component fault in the inflation system detected during assist mode.
1) Press “Start” button to attempt to resume augmentation. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Reservoir pressure regulation failure.

System response: Standby / IAB deflation

Alarm reset: Press “Start” button to resume assist and reset.

Drive Pressure Insufficient
Possible component fault in the inflation system detected during assist mode.
1) Press “Start” button to attempt to resume augmentation. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact Lifetech Cardio service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Solenoid valve not actuating or reservoir pressure regulation failure.

System response: Standby / IAB deflation

Alarm reset: Press “Start” button to resume assist and reset.

Abnormal Negative Pressure in Vacuum Reservoir
Possible component fault in the inflation system detected during assist mode.
1) Press “Start” button to attempt to resume augmentation. 2) If the alarm persists, replace the device and contact Lifetech Cardio service personnel.

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Operating modes: All

Trigger sources: All

Cause: Abnormal negative pressure in the vacuum reservoir.

System response: Awaiting augmentation / IAB deflation.

Alarm reset: Automatically resets when negative pressure recovers.

No Trigger in Automatic Mode
ECG and pressure signals are absent or insufficient for triggering.
1) Reconnect or reset ECG electrodes and check all ECG cables/leads. 2) Check pressure channels and fiber-optic IAB sensor cable connections; ensure disposable transducers are not open to air. 3) If pressure signal is weak, aspirate and flush the pressure line.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: Automatic

Trigger sources: ECG or pressure

Cause: ECG or pressure trigger unavailable, lost, or poor quality for a defined duration in automatic mode.

System response: Awaiting augmentation / IAB deflation.

Alarm reset: Automatically resets when trigger is restored.

No trigger in Semi-Automatic Mode
An electrode is disconnected or improperly positioned.
1) Reconnect or reposition ECG electrodes.
Insufficient signal acquisition.
1) Select another ECG lead or external ECG source via the ECG Signal Source selector.
Selected trigger source is inappropriate.
1) Select another trigger source in the Semi-automatic menu. Press “Start” button to resume augmentation.
Pacemaker pulses not detected in V/A-V pacemaker trigger mode.
1) Adjust pacemaker detection level to enhance pacemaker pulse recognition.

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Trigger sources: R-wave–based (ECG or atrial pacemaker) or pacemaker-pulse–based (ventricular / AV sequential pacemaker).

Cause: ECG or pacemaker trigger unavailable, lost, or poor quality for a defined duration during assist.

System response: Awaiting augmentation / IAB deflation.

Alarm reset: Automatically resets when trigger is restored.

No Pressure Trigger
Pressure waveform absent or flat.
1) Confirm the desired AP Source is selected: Direct or External. Use the Pressure Signal Source selector to choose the appropriate source. 2) Check all pressure channels and verify the fiber-optic IAB catheter sensor cable is properly connected. If using a disposable pressure transducer, ensure it is not open to air. 3) If using a disposable transducer, aspirate and flush the pressure line to clear the lumen. If unsuccessful, occlude the lumen, and use another AP Source (e.g., radial artery) with a disposable transducer and zero it, or connect the external monitor pressure output to the IABP via cable.
The patient’s pulse pressure is insufficient for pressure triggering.
1) If using a disposable pressure transducer, aspirate and flush the pressure line to improve the pressure signal. 2) If possible, select another trigger source in Semi-automatic mode or switch to Automatic mode; press “Start” button to resume augmentation. 3) If pressure triggering is required, lower the trigger threshold using the down button for pressure threshold.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: Semi-automatic

Trigger sources: Pressure

Cause: In pressure trigger mode, a valid trigger is unavailable, lost, or of poor quality for a defined duration (not applicable to automatic mode).

System response: Awaiting augmentation / IAB deflation.

Alarm reset: Automatically resets when a valid pressure trigger is detected or another trigger source is selected.

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Trigger Interference
ESU (electrosurgical unit) interference detected in V/A-V pacemaker trigger mode.
1) Augmentation will resume automatically after interference clears. 2) If the condition persists, select Pressure trigger in the Semi-automatic menu. 3) Check timing settings and press “Start” button to resume augmentation.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: Semi-automatic

Trigger sources: Ventricular / AV sequential pacemaker

Cause: Sustained ESU noise detected when pacemaker trigger source is selected.

System response: Awaiting augmentation / IAB deflation.

Alarm reset: Automatically resets when ESU interference stops, or by manually canceling the pacemaker trigger source.

Check Pacemaker Timing
The patient is not being paced 100%.
1) When the pacemaker operates in demand pacing mode, always use ECG triggering. 2) Select ECG trigger in the Semi-automatic menu. 3) Press “Start” button to resume augmentation.
The system cannot recognize the atrial-ventricular sequential pacing interval.
1) The typical A-V interval (about 170 ms) may shorten when the pacing rate exceeds 130–135 bpm.
Ventricular pacemaker rate is set too high.
1) Confirm the ventricular pacing rate does not exceed 180 bpm.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: Semi-automatic

Trigger sources: Ventricular / AV sequential pacemaker

Cause: Beat-to-beat pacing interval variation exceeds 6.25% or ventricular pacing rate exceeds 180 bpm for a defined duration.

System response: Awaiting augmentation / IAB deflation.

Alarm reset: Automatically resets when acceptable pacemaker timing conditions are met, or by manually canceling the pacemaker trigger source.

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Augmentation Pressure Below the Set Limit.
The patient’s hemodynamic status has changed.
1) Assess the patient’s hemodynamics and attempt optimization. 2) If necessary, press the down button for the augmentation pressure alarm to lower the alarm limit to 8–10 mmHg below the patient’s augmentation pressure.
Augmentation pressure alarm limit is set too high.
1) Press the down button to lower the alarm limit to 8–10 mmHg below the patient’s augmentation pressure.
“Augmentation pressure level” is set too low.
1) Press the up button in the Augmentation Pressure Level menu to increase IAB augmentation pressure until it reaches the maximum setting.
Balloon catheter malposition.
1) Confirm the IAB catheter position and reposition if necessary.
Timing error present.
1) Assess whether inflation timing is too late or deflation timing is too early, and correct as needed.
Balloon leak.
1) Check for blood in the tubing. If blood is present, stop augmentation immediately and notify the physician; remove the IAB catheter per the catheter manufacturer’s instructions.

Alarm attributes: Physiological alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Augmentation pressure decreases below the alarm limit for a defined duration.

System response: Not affected.

Alarm reset: Automatically resets when augmentation pressure rises above the alarm limit.

System Temperature Too High
Over-temperature condition detected.

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1) Check whether there is a high-temperature source around the device and remove it promptly if present. 2) Open Settings, tap Power Off, then tap Confirm to shut down the IABP. 3) Wait 10 seconds. 4) Press and release the main power button to power on the IABP. 5) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Over-temperature detected in the processor module.

System response: The system is disabled and the pneumatic module vents to atmosphere/ IAB deflation

Alarm reset: Try resolving by power-cycling the device.

Host–Controller Communication Interrupted
Internal communication fault.
1) If the alarm persists for 2 minutes, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Internal electronic communication failure

System response: Deflation and system disabled / IAB deflation. Physiological signals remain visible. Press “Start” button or long-press “IAB Fill” button to re-enter self-test; assistance resumes after a successful self-test.

Alarm reset: Refer to the measures above.

Atmospheric Pressure Out of Range
The system detects atmospheric pressure beyond the allowable range.
1) Confirm the ambient pressure is within the normal range. 2) Open Settings, tap Power Off, then Confirm to shut down the IABP. 3) Wait 10 seconds. 4) Press and release the main power button to turn on the IABP. 5) If the alarm persists, replace with another available Lifetech Cardio

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IABP and contact service personnel.

Alarm attributes: Other alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Atmospheric pressure exceeds the allowable range.

System response: Deflation and system disabled / IAB deflation.

Alarm reset: Try resolving by power-cycling the device.

Suspected Blood in the Helium Pathway
Suspected IAB catheter leak.
1) Check the tubing for blood. If blood is present, stop augmentation immediately and notify the physician; remove the IAB catheter per the manufacturer's instructions. 2) If no blood is found, press "Start" button to enter self-test, then re-inflate the IAB catheter and resume augmentation. 3) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: IAB catheter leak causing blood to return into the system.

System response: Deflation and system disabled / IAB deflation.

Alarm reset: Try resolving by power-cycling the device.

4.4.2. Medium-Priority Alarms

Medium-priority alarms require immediate operator attention. They do not interrupt augmentation but may indicate corrective actions are needed. A yellow alarm icon with two exclamation marks appears in the message area. All medium-priority alarms share a common audio tone.

Figure 4-6: Medium-Priority Alarm Icon

Insufficient Drive-Side Negative Pressure During Pause State
Possible component fault while assistance is paused.
1) As needed, press "Start" button to attempt resuming augmentation, or long-press "IAB Fill" button to remain/re-enter pause state.

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2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Drive-side negative pressure cannot be maintained while the pump is paused.

System response: Deflation / IAB deflation.

Alarm reset: Press “Start” button or long-press “IAB Fill” button for 2 s to retry automatic inflation.

Abnormal Positive-Pressure Reservoir During Pause State
Possible component fault while assistance is paused.
1) As needed, press “Start” button to attempt resuming augmentation, or long-press “IAB Fill” button to remain/re-enter pause state. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Positive-pressure reservoir regulation failure while the pump is paused.

System response: Deflation / IAB deflation.

Alarm reset: Press “Start” button or long-press “IAB Fill” button for 2 s to retry automatic inflation.

Abnormal Vacuum Reservoir Pressure During Pause State
Possible component fault while assistance is paused.
1) As needed, press “Start” button to attempt resuming augmentation or long-press “IAB Fill” button to remain/re-enter pause state. 2) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Vacuum reservoir pressure regulation failure while the pump is paused.

System response: Deflation / IAB deflation.

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Alarm reset: Press “Start” button or long-press “IAB Fill” button for 2 s to retry automatic inflation.

Catheter Alarm Paused–Augmentation Pressure Alarm Off
When catheter alarms are paused, the augmentation pressure alarm is set to Off.
1) Set the augmentation pressure alarm to 12the patient’s augmentation pressure.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Catheter alarms are paused when the augmentation pressure alarm is turned off.

System response: Not affected.

Alarm reset: After the user re-enables the augmentation pressure alarm.

No Available AP Source
In automatic mode, no direct or external arterial AP Source detected.
1) If using a fiber-optic IAB catheter, ensure the sensor cable is connected; the system will auto-detect within 20 s. 2) If using a disposable pressure transducer, ensure the invasive pressure cable connects the transducer and the IABP. If the alarm persists, consider replacing the cable. 3) If still unavailable, connect the external monitor’s arterial pressure output to the IABP via cable.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: Automatic

Trigger sources: ECG or Pressure

Cause: ECG signal is normal, but no direct or external pressure source detected.

System response: Augmentation continues with conservative timing selection.

Alarm reset: Automatically resets when a pressure source becomes available.

No ECG signal in Auto mode
ECG not detected or ECG quality is poor in automatic mode.
1) Ensure electrodes are properly connected and optimize placement

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to improve ECG; replace electrodes if necessary. 2) Check all cables/leads for proper connection.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: Automatic

Trigger sources: ECG

Cause: Poor ECG signal quality detected in automatic mode.

System response: Augmentation continues with conservative timing selection.

Alarm reset: Automatically resets when good ECG or pressure signals are detected.

Pressure Sensor SP2 Abnormal
Internal electronic component issue.
1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel for internal repair.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Internal electronic component issue.

System response: Not affected.

Alarm reset: Refer to the corrective actions above.

Pressure Sensor SP3 Abnormal
Internal electronic component issue.
1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel for internal repair.

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Cause: Internal electronic component issue.

System response: Not affected.

Alarm reset: Refer to the corrective actions above.

Pressure Sensor SP4 Abnormal
Internal electronic component issue.
1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel for internal repair.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Internal electronic component issue.

System response: Not affected.

Alarm reset: Refer to the corrective actions above.

Pressure Sensor SP5 Abnormal
Internal electronic component issue.
1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel for internal repair.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Internal electronic component issue.

System response: Not affected.

Alarm reset: Refer to the measures above.

Power Management Module Communication Failure
Internal component communication fault.

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1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel for internal repair.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Internal component communication failure within the augmentation pump.

System response: Not affected.

Alarm reset: Refer to the measures above.

Fiber-Optic Sensor Failure.
Communication failure or recognition error between the fiber-optic sensor signal and the IABP.
1) Firmly grasp the fiber-optic sensor connector; unplug it from the input port and reinsert it. 2) Check the fiber-optic cable for visible kinks and straighten them if present. 3) If the issue persists, disconnect the fiber-optic sensor connector. 4) Provide an alternative AP Source (e.g., radial artery) by connecting a disposable pressure transducer and zeroing it, or connect the external monitor’s arterial pressure output to the IABP via cable. 5) If an external arterial pressure is available and the IABP is in semi-automatic mode, press the Pressure Signal Source toggle button to manually switch the pressure input to External. Confirm triggering and timing settings are appropriate. 6) If no alternative AP Source is available, consider replacing the IAB catheter.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Communication failure or recognition error between the fiber-optic sensor signal and the IABP.

System response: With ECG triggering, augmentation continues, but fiber-optic pressure waveforms and readings are not provided.

Alarm reset: When the fiber-optic sensor connector is removed.

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System Temperature High
High system temperature detected.
1) Check for high ambient heat sources around the device and eliminate them promptly if present. 2) Open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 3) Wait 10 seconds. 4) Press and release the main power button to turn on the IABP. 5) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: High system temperature detected.

System response: Not affected. If the temperature continues to rise, the system will exit assistance.

Alarm reset: The message clears when the temperature returns to the allowable range.

Critically Low Lithium Battery Level
Lithium battery level detected below 10%.
1) Connect the system to AC power.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Low lithium battery level detected.

System response: Not affected.

Alarm reset: Automatically clears the message and silences the tone when AC power is restored.

Failed to Save Log (File Write Failure)
File write failure due to an internal system fault.
1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP.

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4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Internal system fault causing file write failure.

System response: Not affected.

Alarm reset: Refer to the measures above.

Insufficient Available Helium
Helium Cylinder is closed.
1) Open the Helium Cylinder.
Remaining helium supports fewer than 2% of the expected automatic inflation cycles.
1) Replace the Helium Cylinder.

Alarm attributes: Other alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: In hospital configuration use, helium supply falls below the preset reserve determined by cylinder pressure.

System response: Not affected.

Alarm reset: Automatically resets when helium supply is restored.

Gas Pump Control Abnormal
A component fault may exist in the system.
1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel to service internal components.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: The gas pump did not reduce speed or stop as commanded.

System response: Not affected.

Alarm reset: Refer to the measures above.

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Leak Detected at the Front End of the Helium Pathway
The system detects a leak at the interface between the high-pressure Helium Cylinder and the bracket.
1) Replace the cylinder gasket and reinstall the Helium Cylinder. Ensure there is no leak at the cylinder-to-bracket connection. 2) Consider going to Settings → System Options → Tools and setting the Helium Front-End Leak alarm to Off. 3) Once the leak is resolved, promptly set the Helium Front-End Leak alarm to On; the IABP will recheck for leaks. 4) If the leak persists, replace with another available Lifetech Cardio IABP at an appropriate time and contact service personnel to service internal components. WARNING: Lifetech Cardio notice: If an IABP device has a helium front-end leak, continuous helium supply cannot be guaranteed. For long-term counterpulsation therapy, address the helium leak promptly or replace with another available Lifetech Cardio IABP.
Internal helium sensor failure.
1) Replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: See the possible causes above.

System response: Not affected.

Alarm reset: Refer to the measures above.

Lithium Battery not Charging
In AC mode, the lithium battery is not charging.
1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel to service internal components.

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Cause: Lithium battery not charging in AC mode.

System response: Not affected.

Alarm reset: Refer to the measures above.

Helium Pathway Humidity too High
Possible leakage from the balloon catheter.
1) Check the catheter tubing for blood. If blood is present, stop augmentation immediately and notify the physician; remove the IAB catheter per the manufacturer's instructions. 2) Confirm the IAB extension tubing size is appropriate and all connections are secure; ensure the tubing is unobstructed and not restricted.
A component fault may exist in the system.
1) After this use, open "Settings", tap "Power Off", then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel to repair internal components.

Alarm attributes: Other alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause：Helium pathway humidity exceeds the allowable limit.

System response: Not affected.

Alarm reset: Refer to the measures above.

Helium Sensor Abnormal
A component fault may exist in the system.
1) After this use, open "Settings", tap "Power Off", then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel to repair internal components.

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Cause: Helium sensor abnormal.

System response: Not affected.

Alarm reset: Refer to the measures above.

Humidity Sensor Abnormal
A component fault may exist in the system.
1) After this use, open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel to repair internal components.

Alarm attributes: Other alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Humidity sensor abnormal.

System response: Not affected.

Alarm reset: Refer to the measures above.

PM Hardware Failure
A critical component of the IABP has failed.
1) open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Non-interlock alarm.

Operating modes: All.

Trigger sources: All.

Cause: One or more electrical components failed detection during system operation.

System response: Not affected.

Alarm reset: Try resolving the issue by repeatedly powering off and on.

LC Hardware Failure
A critical component of the IABP has failed.
1) Open “Settings”, tap “Power Off”, then confirm to shut down the IABP.

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2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: One or more electrical components failed detection during system operation.

System response: Not affected.

Alarm reset: Try resolving the issue by repeatedly powering off and on.

Disposable Blood Pressure Transducer Failure
Detected open circuit or short circuit in the disposable or invasive blood pressure transducer cable.
1) Replace with a usable disposable blood pressure transducer or invasive pressure cable.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Disposable blood pressure transducer failure detected.

System response: Not affected.

Alarm reset: Automatically resets when the transducer or cable returns to normal or is disconnected.

4.4.3. Low-Priority Alarms

Low-priority alarms indicate conditions requiring operator attention. A dark blue icon with one exclamation mark indicates a low-priority alarm. All low-priority alarms share a common audio tone.

Figure 4-7： Low-Priority Alarm Icon

Software Version Mismatch
Internal software version mismatch detected at startup.

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1) Replace with another available Lifetech Cardio IABP and contact service personnel to update the software.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Software version mismatch detected at startup.

System response: Physiological signals can be displayed, but assistance cannot be provided.

Alarm reset: Refer to the measures above.

Sensor SP2 Communication Failure
Internal sensor communication fault detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Sensor SP2 communication fault detected at startup.

System response: Physiological signals can be displayed; pressing “Start” button or long-pressing “IAB Fill” button re-enters self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Sensor SP3 Communication Failure
Internal sensor communication fault detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm.

Operating modes: All.

Trigger sources: All.

Cause: Sensor SP3 communication fault detected at startup.

System response: Physiological signals can be displayed; pressing “Start” button or long-pressing “IAB Fill” button re-enters self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

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Sensor SP4 Communication Failure
Internal sensor communication fault detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Sensor SP4 communication fault detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Sensor SP5 Communication Failure.
Internal sensor communication fault detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm.

Operating modes: All.

Trigger sources: All.

Cause: Sensor SP5 communication fault detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

On-Board Pressure Sensor Failure
On-board pressure sensor communication fault detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

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Cause: On-board pressure sensor communication fault detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Humidity Sensor Failure
Humidity sensor communication fault detected at startup.
1) Open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm.

Operating modes: All.

Trigger sources: All.

Cause: Humidity sensor communication fault detected at startup.

System response: Not affected.

Alarm reset: Refer to the measures above.

Sensor SP2 Calibration Failure
Internal sensor calibration failure detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Sensor SP2 calibration failure detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Sensor SP3 Calibration Failure
Internal sensor calibration failure detected at startup.

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1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Sensor SP3 calibration failure detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Sensor SP4 Calibration Failure.
Internal sensor calibration failure detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Sensor SP4 calibration failure detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Sensor SP5 Calibration Failure
Internal sensor calibration failure detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Sensor SP5 calibration failure detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available

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after passing self-test.

Alarm reset: Refer to the measures above.

Pump Not Operating
A component fault may exist in the system.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Pump not operating properly detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Positive-Pressure Reservoir Pressure Abnormal During Startup Self-Test
A component fault may be detected in the system at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm.

Operating modes: All.

Trigger sources: All.

Cause: Positive-pressure reservoir pressure abnormal detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Vacuum Reservoir Pressure Abnormal During Startup Self-Test
A component fault may be detected in the system at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

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Operating modes: All

Trigger sources: All

Cause: Vacuum reservoir pressure abnormal detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Pump Speed Regulation Failure
Pump speed regulation failure detected at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Pump speed regulation failure detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

High-Flow Valve Leakage Detected
A component fault may be detected in the system at startup.
1) Open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: High-flow valve leakage detected at startup.

System response: Not affected.

Alarm reset: Refer to the measures above.

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Vacuum Proportional Valve Regulation Failure
A component fault may be detected in the system at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm.

Operating modes: All.

Trigger sources: All.

Cause: Vacuum proportional valve regulation failure detected at startup.

System response: Physiological signals can be displayed. Press “Start” key or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Vacuum Solenoid Valve not Actuating
A component fault may be detected in the system at startup.
1) Open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel to service internal components.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Vacuum solenoid valve not actuating detected at startup.

System response: Not affected.

Alarm reset: Refer to the measures above.

Leak Detected When the Drive Side is Under Negative Pressure
A component fault may be detected in the system at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

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Trigger sources: All

Cause: Leak detected at startup when the drive side is under negative pressure.

System response: Physiological signals can be displayed. Press “Start” key or long-press "IAB Fill" button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Positive-Pressure Proportional Valve Regulation Failure
A component fault may be detected in the system at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Positive-pressure proportional valve regulation failure detected at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Positive-Pressure Solenoid Valve not Actuating
A component fault may be detected in the system at startup.
1) Open “Settings”, tap “Power Off”, then confirm to shut down the IABP. 2) Wait 10 seconds. 3) Press and release the main power button to turn on the IABP. 4) If the alarm persists, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Technical alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Positive-pressure solenoid valve not actuating detected at startup.

System response: Not affected.

Alarm reset: Refer to the measures above.

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Leak Detected When the Drive side is Under Positive Pressure
A component fault may be detected in the system at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Leak detected at startup when the drive side is under positive pressure.

System response: Physiological signals can be displayed. Press “Start” button or long-press “IAB Fill” button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Self-Test Timeout
Self-test timeout at startup.
1) Press “Start” button to enter self-test. 2) If the alarm persists after self-test, replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm.

Operating modes: All.

Trigger sources: All.

Cause: Self-test timeout at startup.

System response: Physiological signals can be displayed. Press “Start” button or long-press "IAB Fill" button to re-enter self-test. Assistance is available after passing self-test.

Alarm reset: Refer to the measures above.

Device Serial Number Error
Software detected an error in the device serial number.
1) Replace with another available Lifetech Cardio IABP and contact service personnel.

Alarm attributes: Other alarm; Interlock alarm.

Operating modes: All

Trigger sources: All

Cause: Device serial number error detected by software at startup.

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System response: Physiological signals can be displayed, but assistance is not available.

Alarm reset: Refer to the measures above.

Irregular Pressure Triggering
The patient’s rhythm varies excessively, making it impossible to accurately predict the next trigger event.
1) The system automatically compensates by early deflation to avoid interfering with systolic ejection. Therefore, do not adjust deflation timing. 2) If the ECG is also unstable and arrhythmia persists, consider ECG triggering for more reliable triggering.
The rhythm is regular, but the IAB deflation timing is set too late, interfering with systolic pressure detection.
1) Confirm the ECG rhythm is regular. If so, set deflation timing earlier to improve pressure-trigger consistency.

Alarm attributes: Other alarm; Non-interlock alarm

Operating mode: Semi-automatic

Trigger source: Pressure

Cause: Pressure triggering abnormal due to: 1) arrhythmia; 2) deflation timing too late, interfering with pressure pulse detection.

System response: The system automatically deflates based on a shorter beat-to-beat interval. Early deflation helps avoid interference at higher heart rates.

Alarm reset: Automatically resets when trigger intervals return to normal.

Low Lithium Battery Level
Lithium battery level detected below 20%.
1) Connect the system to AC power.

Alarm attributes: Technical alarm; Non-interlock alarm

Operating modes: All

Trigger sources: All

Cause: Lithium battery level is low.

System response: Not affected.

Alarm reset: Refer to the measures above.

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The system detected a report export failure.
1) Ensure the removable storage device inserted into the IABP is secure and fully inserted, and that it has more than 5 MB of free space. 2) Or confirm the device is connected to a printer and the printer is powered on. 3) Press the “Report” button to export the report again. 4) If the alarm persists, contact Lifetech Cardio service personnel.

Alarm attributes: Other alarm; Interlock alarm

Operating modes: All

Trigger sources: All

Cause: Failed to export the report to removable storage, or the printer is not connected / malfunctioning.

System response: Not affected.

Alarm reset: Refer to the measures above.

4.4.4. Messages

Messages provide system information to the operator. Messages are text-based. A message icon will appear in the information display area.

Figure 4-8 Message Icon

Messages provided by the IABP include, but are not limited to, the following:

Balloon-End Pressure Abnormal in Pause State
Balloon-end leakage detected in Pause state.
1) Confirm the IAB extension tubing connected to the IAB catheter and IABP is of the correct size, all connections are secure, and the tubing is unobstructed. 2) If augmentation is required, press “Start” button to resume augmentation.

Operating modes: All

Trigger sources: All

Cause: Balloon-end leakage detected while assistance is paused.

System response: Deflation / IAB deflation.

Message reset: The message will disappear after either long-pressing "IAB Fill" button for 2 seconds or pressing “Start” button.

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Standby Time too Long
The IABP has remained in Standby for too long.
1) Check whether pump operation can be resumed. 2) Press “Start” button to resume augmentation. WARNING: To avoid thrombus formation, avoid leaving the IAB catheter inactive (neither inflation nor deflation) for more than 30 minutes.

Operating modes: All

Trigger sources: All

Cause: The IABP system has been in Standby for at least 10 minutes.

System response: Not affected.

Message reset: Press “Start” button to clear the message.

Self-Test not Completed, Command not Allowed
Self-test is not complete; commands are temporarily unavailable.
1) Wait for the IABP system to complete self-test.

Operating modes: All

Trigger sources: All

Cause: A user command was received while the pump is in self-test.

System response: Not affected.

Message reset: Clear the message after self-test is completed.

Self-Test in Progress
The device has entered self-test.
1) Wait for the IABP system to complete self-test.

Operating modes: All

Trigger sources: All

Cause: The IABP system is in self-test.

System response: Not affected.

Message reset: Clear the message after self-test is completed.

Self-Test Completed
The pre-startup self-test has been completed.
No relevant help information.

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Trigger sources: All

Cause: The IABP system self-test completed.

System response: Not affected.

Message reset: Clear the message after starting assistance.

IAB Circuit Leak
Gas leakage or rapid loss in the IAB circuit detected during assistance.
1) Check for any signs of blood in the tubing. If blood is present, stop counterpulsation immediately, notify the physician, and remove the IAB catheter per the catheter manufacturer's IFU. 2) If no blood is detected, confirm all tubing connections are secure. Long-press "IAB Fill" button for 2 seconds to perform automatic inflation, then press "Start" button to resume counterpulsation. 3) If the patient has fever or tachycardia, consider initiating automatic inflation earlier than the routine 2-hour interval to increase inflation frequency.

Operating modes: All

Trigger sources: All

Cause: Leak alarm is disabled while gas leakage or rapid loss exists in the IAB circuit.

System response: Not affected because the alarm has been overridden.

Message reset: The message and tone will disappear after holding "IAB Fill" button for 2 seconds.

Catheter Alarms Paused
The user is detected to have paused catheter alarms.
1)Turn catheter alarms back on when needed. WARNING: To prevent irreversible harm due to catheter alarms being disabled, the IABP will automatically re-enable catheter alarms 60 minutes after they are turned off.

Operating modes: All

Trigger sources: All

Cause: This message occurs when catheter alarms are set to "Pause".

System response: Not affected.

Message reset: Refer to the measures above.

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The leak alarm is set to Off.
1) If permitted, set the Leak alarm to On.

Operating modes: All

Trigger sources: All

Cause: The user set the leak alarm to Off.

System response: Not affected.

Message reset: Clear the message after the user sets the leak alarm to On.

Front-End Helium Line Leak Alarm Set to Off
The front-end helium line leak alarm is set to Off.
1) If permitted, set the Front-end Helium Line Leak alarm to On.

Operating modes: All

Trigger sources: All

Cause: The user set the Front-end Helium Line Leak alarm to Off.

System response: Not affected.

Message reset: Clear the message after the user sets the Front-end Helium Line Leak alarm to On.

Inflation Completed
IAB inflation has been completed.
No relevant help information.

Operating modes: All

Trigger sources: All

Cause: The system successfully completed inflation.

System response: Not affected.

Message reset: Clear the message after starting assistance.

Operating on Lithium Battery Power
The system is running on lithium battery power. If possible, switch to AC power. If AC power is properly connected but this message still appears, a power fault may exist. Replace with another available Lifetech Cardio IABP and contact service personnel.
No relevant help information.

Operating modes: All

Trigger sources: All

Cause: The IABP system is operating on lithium battery power.

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System response: Not affected.

Message reset: Clear the message after connecting AC power.

Alarm Tone Volume is Low, Please Increase the Volume
The alarm tone volume is set too low.
No relevant help information.

Operating modes: All

Trigger sources: All

Cause: The alarm tone volume is set too low.

System response: Not affected.

Message reset: Increase the volume, then clear the message.

No Trigger in Auto Mode
ECG and arterial pressure signals are absent or insufficient for triggering.
1) Reconnect or reset the ECG electrodes and check that all ECG cables/leads are properly connected. 2) Check the corresponding pressure channel and ensure the fiber-optic IAB catheter sensor cable is properly connected. If using a disposable blood pressure transducer, confirm it is not open to air. 3) If the pressure signal is too weak, attempt to improve signal quality by aspirating and flushing the pressure line.

Operating mode: Auto

Trigger source: ECG or Pressure

Cause: No trigger detected while the IABP is in Standby.

System response: No action; message displayed with no alarm tone.

Reset: Automatically resets when a valid trigger is detected or when “Start” is pressed.

No Trigger in Semi-Auto Mode
An electrode is disconnected or incorrectly positioned.
1) Reconnect or reposition the ECG electrodes.
Signal acquisition is inadequate.
1) Use the ECG Signal Source switch button to select another ECG lead or an external ECG signal.
The selected trigger source is inappropriate.

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1) Select another trigger source in the Semi-automatic menu. Press “Start” button to resume augmentation.
No pacer pulse detected in Pacer V/AV trigger mode.
1) Adjust the pacer detection level appropriately to enhance pacer pulse recognition.

Operating mode: Semi-automatic

Trigger source: R-wave source trigger (ECG & Pacer A) or pacer-pulse source trigger (Pacer V/AV).

Cause: No trigger detected while the IABP is in Standby.

System response: No action; message displayed with no alarm tone.

Reset: Automatically resets when a valid trigger is detected or when “Start” is pressed.

No Pressure Trigger
Pressure waveform is absent or flat.
1) Confirm the required AP Source (Direct or External) has been selected. Use the Pressure Signal Source switch button to select the appropriate AP Source. 2) Check all pressure channels and ensure the fiber-optic IAB catheter sensor cable is properly connected. If using a disposable blood pressure transducer, confirm it is not open to air. 3) If using a disposable blood pressure transducer, you may aspirate and flush the pressure line to clean the lumen. If unsuccessful, seal the lumen. Then either connect a disposable blood pressure transducer to another AP Source (e.g., radial artery) and zero it, or connect the external monitor’s arterial pressure output signal to the IABP via cable.
The patient’s pulse pressure is insufficient for pressure triggering.
1) If using a disposable blood pressure transducer, try aspirating and flushing the pressure line to improve the pressure signal. 2) If possible, select another trigger source in the Semi-automatic menu or switch to “Auto” mode. Press “Start” button to resume augmentation. 3) If pressure triggering is required, lower the trigger threshold using the “Threshold” Down button.

Operating mode: Semi-automatic

Trigger source: Pressure.

Cause: No trigger detected while the IABP is in Standby.

System response: No action; message displayed with no alarm tone.

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Reset: Automatically resets when a valid trigger is detected or when “Start” is pressed.

4.4.5. Alarm Test

This section is intended to allow the operator to verify that the alarm system is functioning. Testing may only be performed when patient support is not required:

Disconnect the balloon helium extension tube from the device balloon port so the port is open to air, then press “Start” button. Within 20 s, a high-priority visual and audible alarm “Helium purge timeout during inflation” will occur. Pressing “Alarm Silence” will mute the alarm for 60 s. The user can verify the alarm system is functioning.

WARNING:

Do not use the alarm test while the IABP is being used for patient treatment.

4.4.6. Sudden Power-Off Alarm

If the IABP powers off unexpectedly, it will emit a sudden power-off alarm tone (e.g., the power plug is unplugged and the lithium battery fails to supply power). The alarm repeats for about 90 s or until the IABP restarts, whichever comes first. No visual indicator is provided. The alarm audio is powered by its dedicated backup battery.

4.5. Alarm Log

The alarm log is used to record information related to alarm faults.

The alarm log records fault information of the alarm system. The log can:

1) Record and identify the occurrence of high-, medium-, and low-priority alarm states, including date and time.
2) Not include when the alarm system lost power; power-off time is not recorded in the log.
3) After a complete power loss (mains and/or internal power) lasting 3 minutes, the log contents are cleared; logging restarts upon the next power-on.
4) When the log is full, the alarm system discards the oldest data.
5) Record the start time of each alarm state via the log.
6) Record technical alarm states for service and maintenance. This log cannot be reset or edited by the operator.

4.6. Clinical Considerations

4.6.1. Condensation

During augmentation, fine mist or small droplets may occasionally be observed

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in the IAB catheter helium extension tubing. This mist is condensation. During each inflation/deflation cycle, the IABP uses its condensation-clearing system to remove water vapor; no operator intervention is required.

If larger droplets appear in the helium extension tubing, remove them by suction. Ensure the extension tubing is disconnected from both the balloon and the pneumatic module.

If the patient is transferred from another IABP system that does not have condensation-clearing capability, ensure all condensation droplets are removed from the helium extension tubing before connecting to this IABP.

Operators should routinely check the IAB catheter helium extension tubing for condensation. Excessive accumulation can affect system performance. Excessive condensation indicates that the pneumatic module requires maintenance.

4.6.2. Operating Procedures

ECG: If reliable ECG acquisition is altered or impeded, several methods may resolve the issue: reposition or replace ECG electrodes, and verify correct patient cable connections. In Semi-auto mode, select an alternative lead configuration. For proper electrode use, follow the techniques recommended in Chapter 4 (ECG acquisition).

Pressure: When using an arterial pressure monitoring line, maintain continuous flushing per standard hospital procedures. When monitoring via the IAB lumen, follow the specific maintenance instructions provided with the IAB. Proper flushing to keep the line patent and setting the stopcock to the correct position can prevent most pressure waveform damping/tailing issues.

Atrial fibrillation: Ensure R-wave tracking is enabled. No special pump considerations are required in Auto mode. In Semi-auto mode, select ECG as the trigger source. By automatically holding IAB inflation until the next R-wave is detected, the system can track irregular rhythms in the most stable manner. Pressure triggering is not recommended in atrial fibrillation.

Premature beats: If a premature R-wave is detected, the IABP automatically deflates and assists the premature beat. In Semi-auto mode, if the premature beat amplitude is small, select ECG triggering to maximize reliable triggering by minimizing amplitude differences between normal QRS complexes and ectopic beats. No special adjustments are required.

Ventricular fibrillation with cardiac arrest: During defibrillation, the IABP has protec-

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tive measures and is fully isolated from the patient and defibrillator electrodes. However, operators must keep a distance from the pump during defibrillation. If the IABP is disconnected from ground (e.g., AC outlet), it is especially important to comply with the above requirements when operating the IABP.

Asystole or prolonged cardiac arrest: If possible, use ECG or pressure triggering during CPR. This helps synchronize assistance with chest compression rate and rhythm.

Trigger source selection: In Auto mode, ECG (R-wave) or pressure is automatically selected as the trigger source based on relative signal quality. If neither ECG nor pressure can provide reliable triggering to start Auto mode, the IABP can be triggered by its internal clock. Select Semi-auto mode and set the trigger source to Internal. Default internal rate is 80 bpm, adjustable from 40 to 120 bpm in the Trigger Signal menu.

Extended standby and thrombus formation on the balloon membrane: The standby time reminder tone is enabled by default at power-on. If the IABP is in Standby, standby time is displayed on the balloon pressure waveform. After 10 minutes, a “standby too long” alarm will be displayed. Check whether assistance can be resumed; if yes, press “Start” to resume.

Hazardous patient leakage current: Connecting multiple medical devices to the patient may cause cumulative leakage currents that exceed limits allowed by safety standards.

4.6.3. Use in an Electrosurgical Environment

The IABP incorporates electrosurgical interference suppression technology to minimize ESU noise and reduce disruption to system operation. The system can suppress ESU noise but cannot eliminate it completely. During ESU use, sparking may occur in tissue, generating noise that extends into ECG and pressure signal frequency ranges. Because the system must pass these frequencies, some ESU noise may interfere with ECG signals, especially at higher ESU power settings. When using specified third-party medical sensors and cables, pressure signals are typically unaffected by ESU noise. Refer to “1.9 Compatible Devices/Accessories” for sensors meeting Lifetech Cardio requirements. The IABP is designed to suppress ESU noise and, when recommended cables and sensors are used, to protect the patient from burns at the ECG and pressure system input connections.

Limit the noise energy power. The interference level is directly related to ESU power settings; reduce power as much as possible while achieving the in-

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tended effect. With ESU noise present, successful ECG triggering largely depends on proper patient preparation and ESU use.

Following the listed guidelines can minimize electrical energy from ESU that enters the IABP via ECG input, typically enabling stable ECG triggering. If the clinician prefers pressure triggering in the OR, select “Semi-auto” mode to achieve this, then select “Pressure” in the Trigger Signal menu. Return to Auto mode as appropriate.

Under ESU interference, the following alarms will be suspended until the interference stops:

- Drive pressure too high.
- IAB connection disconnected.
- Balloon pressure too high.
- Air ingress into the IAB circuit.
- IAB circuit leak.

Once ESU interference stops, suspended alarms should automatically resume. If ESU interference has stopped but any of the above alarm conditions remain, the current alarm will be displayed.

When using the IABP in an electrosurgical environment, the following techniques are recommended to minimize ESU interference:

- Keep ECG and invasive BP transducer cables as perpendicular as possible to ESU cables.
- Keep ECG electrodes as far as possible from the surgical site.
- Place ECG electrodes approximately equidistant from the surgical site to minimize potential differences between electrodes.
- Place all ECG electrodes on the same plane (anterior or posterior) to minimize potential differences.
- Keep all fluid pressure tubing away from ESU cables as much as possible.
- Place the ESU return (grounding) pad directly under the surgical site.
- Use as little electrosurgical equipment as possible.

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4.6.4. Use During Cardiopulmonary Bypass

During full cardiopulmonary bypass, select Semi-auto mode and set “Internal” as the trigger source so the IABP can inflate and deflate the IAB. The internal rate can be adjusted in the Trigger Signal menu. Inflation and deflation timing controls should be set to their default mid positions.

When “Internal” is the trigger source, if ECG activity is present, the IAB will deflate immediately to avoid asynchronous counterpulsation and systolic interference.

When weaning the patient from bypass, the IABP is used to support cardiac function. If ECG stabilizes, it can be used as the trigger source. Observe the blinking heart icon in the trigger rate parameter to confirm stable ECG triggering; the icon should blink with each detected trigger event.

If the patient’s ECG amplitude is insufficient for triggering, using pressure as the trigger source is another feasible option. The IABP will automatically adjust the pressure threshold down to a minimum of 7 mmHg. In Semi-auto mode, the user may also manually control the pressure threshold level.

As bypass flow decreases and the heart begins generating cardiac output, check whether timing selection is correct. Because patient condition may alter the electromechanical relationship, frequent reassessment is necessary.

PRECAUTION:

During or after surgery, if the patient’s core temperature rises by more than 6°C within 2 hours, while providing assistance press and hold the “Zero” button for 2 seconds to re-zero the fiber-optic IAB catheter (select in vivo zeroing in Settings).

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4.7. Reports

4.7.1. Report Output Operation

Press the “Report” button to output a report according to current output settings. If “PDF” is selected, a PDF will be generated and saved to the user’s removable storage device. If “Print” is selected, the connected printer will print the current report (see Section 3.1.24.3).

PRECAUTION:

For print output, connect a printer to the USB port (see Section 3.2). For a better printing experience, EPSON L-series printers are recommended. The connected printer should comply with GB 4943.1 (recommended).

When output is set to “Print”, pressing “” starts report generation. If the IABP detects no printer connected, a pop-up will prompt that no printer is connected and allow the user to reselect Print, PDF output, or cancel. If printing fails during printing, a pop-up will prompt the user to retry printing, save to external storage (PDF output), or cancel printing.

When output is set to “PDF”, pressing “Report” starts report generation. If the IABP detects no external storage inserted, a pop-up will prompt that no external storage is inserted and allow the user to choose Print, PDF output, or cancel. If insufficient space is detected on the external storage during saving, a pop-up will prompt insufficient storage space and allow the user to replace the storage device, choose print output, or cancel output.

Whether the user selects Print or PDF output, data from the first failed PDF output attempt is not lost. The subsequently output report uses the data from the initial PDF output attempt.

If necessary, printing can be automatically initiated upon alarm events or at user-selected fixed intervals. Print initiation, format, and length depend on Settings → Report (see Section 3.1.24.3). Print speed (25 or 50 mm/s) is controlled by scan speed settings (see Section 3.1.24.1, Settings → Display).

4.7.2. Report Format

Waveforms can be exported as a patient waveform PDF for copying, or printed. For print format, see Settings → Report. Three-trace, two-trace, or single-trace formats are available. Supported waveform selections include: ECG, AP, BPW, i.e., ECG+AP+BPW, ECG+AP, ECG+BPW, AP+BPW, ECG only, AP only, and BPW only. See report output options in Section 3.1.24.3 (Settings → Report Output).

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For three- or two-waveform settings, two or three waveforms (from ECG, AP, BPW) can be recorded simultaneously for the duration specified in the report output tab.

PDF filename format: waveform type_YYYYMMDD-HHMMSS.pdf

Example: ECG_AP_BPW_20231114-102731.pdf

4.7.3. Annotations and Examples

Upon completion of recording, corresponding annotations will be printed at the beginning of the report. The annotations describe patient heart rate and filter information. Relevant system settings are also included, such as IAB assist ratio, operating mode and trigger source selection, ECG lead, and sweep speed.

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5. Installation Information

After the Intra-Aortic Balloon Pump is delivered to the purchasing hospital or other buyer, Lifetech Cardio will dispatch qualified engineers or designated representatives for on-site installation.

Before installation, ensure the following conditions are met:

Power supply: AC 220 V, 50 Hz, with protective earth grounding.

Installation environment: Ensure the device is placed on a stable, level surface with sufficient space for operation and heat dissipation.

Ensure an adequate helium supply is available for installation.

Installation service must be performed by personnel trained by our company.

It is recommended that medical staff involved in installation acceptance be physicians or technical professionals who have received relevant system training. The installation engineer shall provide on-site training to key operators.

Prior to installation, verify the device and accessories against the product configuration list, confirm completeness, and check whether the intra-aortic balloon pump was damaged during transportation.

WARNING:

This device must be connected to a power supply with protective earth grounding.

Read the Instructions for Use thoroughly and receive operational training before use.

Unauthorized modification of the device is prohibited.

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6. User Maintenance

Lifetech Cardio may, upon request, provide circuit diagrams, component lists, drawings/annotations, calibration procedures, or other information that helps maintenance personnel repair ME equipment parts that are designated by the manufacturer as repairable by authorized maintenance personnel.

6.1. Replacing the Helium Cylinder

The IABP is equipped with an internal Helium Cylinder. The internal Helium Cylinder provides fresh helium for inflating the IAB catheter. The internal Helium Cylinder can support continuous use for 3 months.

[WARNING]

Replacement of the Helium Cylinder shall only be performed by personnel familiar with high-pressure gas cylinder operation.

PRECAUTION:

Only medical helium may be used.

The Helium Cylinder must be replaced when the pressure drops below the preset level. If a “low helium” message is displayed during operation, or the helium pressure gauge indicates the supply has entered the red zone, the cylinder needs to be replaced. The IABP therapy does not need to be interrupted when performing the above operation. However, the cylinder must be replaced as soon as possible to avoid automatic inflation failure and delayed counterpulsation.

Replacing the Helium Cylinder:

1) Open the side panel of the Helium Cylinder compartment.
2) Turn the Helium Cylinder knob clockwise to ensure the cylinder is fully closed.
3) Slowly loosen the T-handle of the cylinder bracket counterclockwise (some helium may escape) until the cylinder can be removed. Support the bottom of the cylinder during loosening to prevent it from falling.
4) Remove the Helium Cylinder.
5) Install a usable Helium Cylinder, aligning the protruding locating point on the cylinder bracket with the locating hole on the cylinder.
6) Fully tighten the T-handle of the bracket clockwise.
7) Install the cylinder knob onto the new Helium Cylinder, then turn the knob counterclockwise to slowly open the cylinder.
8) Close the Helium Cylinder panel.
9) Observe the helium indication on the status display and verify that the helium pressure is within acceptable limits.

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When disposing of used Helium Cylinders, comply with current local regulations and applicable recycling requirements.

6.2. Refilling the Helium Cylinder

When refilling the cylinder, Lifetech Cardio Helium Cylinders with standard post valves must be used. Depending on local conditions, additional local requirements may apply to Helium Cylinder refilling.

6.3. Battery Maintenance

If AC power cannot be supplied to the IABP (intentionally or due to power failure), the system will automatically switch to lithium battery power. When AC power is available, the IABP will automatically use AC power. Meanwhile, if the lithium battery level is low, it will begin charging.

In the information display area, the lithium battery icon indicates the current charging status of the lithium battery (Section 3.1.3).

PRECAUTION:

If a mains power interruption occurs while internal lithium battery power is available, the device will automatically switch to internal lithium battery supply to continue normal operation, and the operating mode, all operator settings, and stored data will not change.

PRECAUTION:

When the device is powered on, the lithium battery charges only when the AC plug is connected to a live AC outlet and the charging icon is shown on the display. When the device is powered off, the lithium battery charges only when the power button "⏻"flashes green. The system can remain in a charging state even when not in use.

Using lithium battery power

- When operating on lithium battery power, the "operating on lithium battery power" message is displayed in the information display area (Section 3.1.2).
- As the lithium battery level decreases, the fill level within the battery icon gradually decreases to reflect the condition. Similarly, the battery icon color transitions from green to yellow and then to red as the battery level drops.
- When the lithium battery level is low, the medium-priority alarm message "Lithium battery critically low" will remain displayed in the information display area

Battery operating time may decrease due to factors such as age, storage temperature, and discharge cycles. If recharging is not performed promptly as recommended and the battery repeatedly undergoes deep discharge cycles, permanent damage may occur.

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PRECAUTION:

When AC power is required for operation, ensure the power plug is connected to a live AC outlet and the lithium battery icon indicates charging. Also ensure that the message “Operating on lithium battery power” is not shown in the information display area.

To achieve optimal lithium battery performance and expected service life, follow the guidelines below:

When the device is not in use, keep the lithium battery fully charged by plugging the device into an AC outlet. If the device is stored for an extended period (2 months or longer) without AC power to maintain the battery level, or if the storage temperature exceeds the maximum operating temperature, disconnect the lithium battery from the main unit. Because the battery self-discharges, fully recharge the disconnected battery at least once every 6 months.

High temperature is extremely harmful to lithium battery life. Do not operate the device if the ambient temperature exceeds the maximum operating temperature.

If a low battery warning appears, recharge within 24 hours to prevent battery damage.

If the lithium battery leaks, deforms, is damaged, or shows abnormal conditions, stop using it.

Replace the lithium battery when total capacity drops below 80% of the rated capacity or when the cycle count exceeds 1500. The battery cycle count and capacity can be viewed in Settings → System → Device Info.

6.4. Backup Batteries and Filter Maintenance

Backup batteries: The main unit has two AA 1.5 V alkaline batteries built in to power the sudden shutdown alarm when an unexpected shutdown occurs (during operation, mains power is interrupted and the lithium battery cannot supply power normally). When stored and used under specified conditions, the expected service life is 4 years.

Backup battery maintenance: Replace the backup batteries after 4 years of storage and use under specified conditions. Stop using them if leakage, deformation, damage, or abnormal conditions occur.

Replacement of backup alkaline batteries: Open the battery cover in the Helium Cylinder compartment door to replace the backup batteries.

Filter replacement maintenance: After opening the battery cover in the Helium Cylinder compartment door, the air-side filter of the vacuum tank can be replaced. Remove the small panel at the IAB catheter connector to replace the corresponding balloon-end filter. Remove the decorative panel of the main unit to replace the filter

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cotton inside the pump-side filter of the vacuum tank.

6.5. Removing the Main Unit from the Cart

Disconnect the video port and video cable from the main unit. From beneath the main unit, release the latch on the cart. Hold the handle and slowly pull the main unit outward to about halfway. Hold the handles on the top and front of the main unit and remove the console from the cart.

PRECAUTION:

When removing the main unit from the cart or installing it onto the cart, ensure the cart wheels are locked.

6.6. Installing the Main Unit on the Cart

Disconnect the cable between the display and the main unit and place the display on the cart. Release the latch at the bottom of the cart. With one hand holding the main unit handle and the other supporting the bottom, push the main unit along the cart rails until the cart pins fully engage into the main unit. Then lift the latch to secure the main unit. Route the display cable through the hole in the cart tabletop, fix the display on the disk damper of the tabletop, and connect the display to the main unit with the display cable.

[WARNING]

When installing the main unit on the cart, ensure the latch is lifted to securely lock the main unit in place.

6.7. Fuse Maintenance

Fuse specification: T5AH/250V. Each device uses two fuses.

6.8. User Cleaning and Maintenance

6.8.1. Cleaning the System and Accessories

Clean the device using a damp sponge and mild soapy water or an ammonia-containing cleaner. Do not use organic solvents or corrosive cleaners. Prevent cleaning agents from entering vents or connector openings. Patient-contact parts (such as ECG lead wires and the disposable blood pressure transducer) shall be kept clean and disinfected. Always follow standard hospital procedures for cleaning and disinfection control. If applicable, contact the disposable blood pressure transducer manufacturer for recommendations on disinfection procedures.

6.8.2. Cleaning the Display Screen

Clean the display screen carefully to prevent scratches. Dust and dirt particles

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can be blown off or wiped with a soft cloth. Use a screen cleaning spray or a lint-free microfiber cloth moistened with solution to remove fingerprints and smudges. Do not use alcohol or solvents containing chlorinated hydrocarbons. Do not apply cleaning agents directly onto the screen.

6.8.3. Cleaning and Disinfection of Cables

1) Disconnect the cable from the device before cleaning the cable.
2) Replace the cable if it is damaged or shows signs of aging.
3) The cable surface may be wiped with a soft cloth dampened with clean water; if necessary, wipe with 75% alcohol, then air-dry naturally or wipe with a clean, dry cloth.
4) Wipe gently during cleaning to avoid damaging the cable.
5) Do not immerse the cable in liquid and do not wet the connector pins.
6) Do not use oils, corrosive liquids, or other solutions to clean or disinfect the cable.
7) Clean the cable before disinfection.
8) After cleaning and disinfection, inspect the cable to confirm it is normal and intact before reuse.

6.9. Preventive Maintenance Plan

6.9.1. Preventive Maintenance Plan A

Preventive Maintenance Plan A shall be performed by clinical users or Lifetech Cardio trained and certified technical personnel.

Required actions	Before and after each use	Monthly
1. Clean the device and accessories as required; inspect patient cables*, pneumatic module Luer connectors for cracks, and check the power cord for damage;	√	
2. Perform the pneumatic test (see Section 3.1.24.6);	√	
3. Check the battery status icon (see Section 3.1.3).	√	
4. Check the operation of automatic inflation and the helium supply;		√
5. Check for alarms and messages related to low available helium, lithium battery, and sensors;		√

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6. Connect a disposable blood pressure transducer to the invasive blood pressure transducer cable and then to the device, apply pressure using a syringe and zero the transducer, and verify that the device displays a blood pressure waveform and any related alarms.		√

Patient-contact components (such as ECG lead wires and disposable blood pressure transducers) shall be kept clean and disinfected.

Always follow standard hospital operating procedures for cleaning and disinfection.

To verify the alarm system, refer to the instructions in Section 4.4.5 Alarm Test.

6.9.2. Preventive Maintenance Plan B

Preventive Maintenance Plan B shall be performed by Lifetech Cardio trained and certified technical personnel.

Required actions	Actions every 12 months or 2,500 hours	Actions at 6,000,000 inflation/deflation cycles or 4 years *	Actions every 4 years:	Actions every 5,000 hours
1. Perform the visual inspection checklist.	√			
2. Remove the rear cover of the main unit and clean internal dust using a vacuum cleaner.	√			
3 Replace the safety diaphragm in the pneumatic module		√		
4 Replace the dehumidification components in the pneumatic module			√	
5 Replace the filter cotton on the vacuum tank pump side, the air-side filter, and the balloon-end filter.				√
6 Service the pump assembly.				√

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7 Replace the alkaline batteries supplying power for the unexpected shutdown alarm.			√	
8 Verify that the pump cooling fan and power supply fan are operating properly	√			
9 Replace the lithium battery* when total capacity drops below 80% of nominal capacity or when cycle count exceeds 1,500*	√			
10 Perform fiber optic testing and clean fiber optic components if necessary.	√			

* Whichever comes first

6.10. Software Maintenance

When new software versions are released, provide corresponding software upgrade services to users, as applicable

6.11. Warranty

From the date of delivery, Lifetech Cardio warrants that the product is free from defects in workmanship and materials for 12 months, except exclusions below:

(1) Disposable products are warranted to be free from defects in workmanship and materials for one year from purchase or first use (whichever comes first).

(2) Filters and tubing are consumables and are not covered by the warranty.

The safety diaphragm originally supplied with the IABP is also warranted for 12 months from delivery.

The recommended preventive maintenance is the user's responsibility. The associated parts and labor are not covered under warranty.

Lifetech Cardio shall not be liable for any direct or indirect incidental, special, or consequential loss, damage, or expense arising from the use of its products. Under

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these warranty terms and the purchaser's exclusive remedies specified herein, Lifetech Cardio's warranty obligation is limited, at its discretion, to repair at the factory or service/replacement provided through an authorized Lifetech Cardio distributor, for products and services that are found to have defects in materials or workmanship under normal use and for which Lifetech Cardio is responsible.

No agent, employee, or representative of Lifetech Cardio is authorized to make any assertion, representation, or warranty on behalf of Lifetech Cardio regarding its products, and any such assertion, representation, or warranty shall not be enforceable.

This warranty expressly replaces any and all other express or implied warranties, including warranties of merchantability or fitness for a particular purpose, and any other implied warranties assumed by the seller.

Damage to products or components caused by misuse, negligence, failure to perform preventive maintenance, accidents, use of accessories and IAB catheters not supplied by Lifetech Cardio, or any modifications made by the customer is not covered by this warranty. Lifetech Cardio makes no warranty for accessories sold and warranted by other companies.

Under the terms of this warranty, upon Lifetech Cardio's authorization, the device or any accessory claimed to be defective must be returned to Lifetech Cardio with freight prepaid. Lifetech Cardio shall not be responsible for any loss or damage during transportation.

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7. External Monitoring Device Interfaces

The IABP can acquire patient ECG and arterial pressure signals from compatible external monitoring devices, and can also provide ECG output signals and arterial pressure output signals to bedside patient monitors; the external monitoring device signal connection ports are located on the IABP panel (see Section 3.2 Panel Description).

Warning: Any external bedside monitor providing ECG signals to the IABP in the operating room must have electrosurgical interference suppression capability.

7.1. ECG Requirements from External Monitoring Devices

Scale factor *:1 V/mV ±10%

Pacemaker enhancement: Enabled (if supported by the monitoring device) to enhance pacemaker pulse detection and blanking.

* Required for correct triggering of the IABP

PRECAUTION: Due to bandwidth limitations of external monitoring device outputs, direct patient leads are preferred to optimize pacemaker detection and blanking

PRECAUTION: Proper pacemaker detection may require the external monitoring device to reinsert pacing events into its ECG signal.

7.2. Arterial Pressure Requirements from External Monitoring Devices

Scale factor*:1 V/100 mmHg ±4%

* Required for correct triggering and accurate pressure measurement by the IABP

7.3. ECG Input Description from Monitoring Devices

Interface specification: 6.35 mm (1/4 inch) audio jack with a green indicator symbol, used to connect the external ECG signal cable (see the external ECG input interface in Section 3.2).

7.4. Arterial Pressure Input Description from Monitoring Devices

Interface specification: 6.35 mm (1/4 inch) audio jack with a red indicator symbol, used to connect the external pressure signal cable (see the external blood pressure input interface in Section 3.2).

7.5. ECG Output Description to Monitoring Devices

Interface specification: 6.35 mm (1/4 inch) audio jack with a green indicator symbol (see ECG output interface in Section 3.2);

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Transmission delay (QRS complex delay): maximum 25 ms;

Output amplitude/scale factor (at 10 Hz): output-to-input signal ratio of 1 V/mV ±20%;

Pacemaker enhancement: ECG output includes pacemaker enhancement, with enhanced pacing pulses summed with the ECG signal.

7.6. Description of Low-Level Pressure Output to Monitoring Devices

Interface specification: 6.35 mm (1/4 inch) audio jack with a red indicator symbol, used to connect the external pressure signal cable (see the blood pressure output interface in Section 3.2).

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8. Technical Specifications

8.1. Triggering

Trigger Sources	Parameters	Requirements
QRS wave triggering	ECG/Pacer A (R-wave detection mode)	Minimum threshold: 80 ±20 μV
	Trigger blanking (R-wave detection mode)	ECG: 44 ms Pacer A: 100 ms Tolerance ±20 ms
	Tall T-wave suppression (R-wave detection mode)	Suppresses all T waves with a suppression duration of 180 ms, Q-T interval ≤350 ms, and amplitude ≤120% of the QRS amplitude
	Pacemaker pulse suppression (R-wave detection mode)	Suppresses all pacing pulses with amplitudes from ±2.0 mV to ±700 mV and pulse widths from 0.1 ms to 2.0 ms, without trailing or overshoot.
Pressure triggering	Pressure threshold	Automatica threshold adjustment: Automatically adjusted to 3/8 of the pulse pressure difference between systolic pressure and end-diastolic pressure (minimum 7 mmHg); Manual adjustment range :7 mmHg to 30 mmHg; accuracy ±3 mmHg for both modes; Default: Automatic; Adjustment increment: 1 mmHg.
Pacemaker triggering	Pacer V/AV	If the ventricular pacing interval is fixed and heart rate is below 180 bpm, it is identified as a ventricle pacemaker; If the atrial-to-ventricular interval is 50–250 ms and is less than or equal to two-thirds of the ventricular-to-atrial interval, it is identified as an AV sequential pacemaker.
Internal triggering	Internal rate	Range: 40 bpm to 120 bpm Default: 80 bpm Increment: 5 bpm.

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8.2. ECG Channel

8.2.1. ECG Signal Input

Parameters	Requirements
ECG leads	Automatic: I, II, III, external Semi-Automatic: I, II, III, aVR, aVL, aVF, V or external default: II
QRS detection	Shall respond to QRS signals with an amplitude of at least 0.25 mV to 5 mV and a duration (width) of 70 ms to 120 ms.
Input dynamic range and differential offset voltage*:	DC offset voltage minimum ±300 mV Linear input minimum ±5.0 mV
Input impedance	≥2.5 MΩ
Defibrillation protection*	After withstanding a defibrillation voltage, the device shall return to the previous operating mode and resume normal operation within 5 seconds, without loss of any operator settings or stored data, and shall continue to perform its intended functions.
Lead fault detection	Detects faults when any energized electrode lead is open-circuited
Input noise	Less than 30 μV p-p
Common-mode rejection	Greater than 80 dB
Signal reproduction accuracy	For input signals within ±5 mV and with a slew rate up to 125 mV/s, the reconstructed output deviation shall not exceed ±20% of the nominal output value or ±100 μV, whichever is greater.
Frequency response*	Bandwidth at least 0.67 Hz to 40 Hz.

Note: Items marked with * comply with the requirements of GB 9706.227-2021.

8.2.2. ECG Waveform Display

Parameters	Requirements
Automatic waveform amplitude scaling	Applicable up to ±5 mV (minimum) without clipping
Resolution (display)	1 bpm
Sweep speed	25 mm/s or 50 mm/s (default 25 mm/s)
Gain indicator*	Provides a gain indicator to represent the amplitude of a 1 mV input voltage. When a (1±0.01) mV input signal is applied on an appropriate lead, the displayed output amplitude variation shall be within ±10%
Pacemaker pulse display capability*	Capable of displaying ECG signals in the presence of pacing pulses (pulse amplitude ±2 mV to ±700 mV, pulse width 0.5 ms to 2.0 ms)
Channel height*	Permanent display vertical height at least 20 mm; Non-permanent display at least 30 mm;
Baseline recovery*	After application of a 1 V peak-to-peak, 50/60 Hz overload voltage for at least 1 s, the device returns to normal operation within 3 s.

Note: Items marked with * comply with the requirements of GB 9706.227-2021.

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8.2.3. ECG Signal Output

Parameters	Requirements
Bandwidth	0.5 to 135 Hz
Transmission delay (QRS complex delay)	Maximum 25 ms (at 10 Hz)
Output-to-input signal ratio (scale factor)	1 V/mV, ±20%

8.3. Pressure Channel

8.3.1. Pressure Signal Input

Signal Source	Parameters	Requirements
Hydraulic sensor input	Pressure range	−30 to 300 mmHg
	Accuracy	±4 mmHg or 4% (whichever is greater)
	Frequency response	DC to 42 Hz (attenuation relative to 1 Hz not exceeding 3 dB)
	Input sensitivity	5.0 μV/V/mmHg (nominal)
	Sensor excitation	+5 V DC ±5% (maximum supply current 25 mA)
Fiber-optic sensor input	Pressure range	−30 to 300 mmHg
	Accuracy	±4 mmHg or 4% (whichever is greater)
	Frequency response	DC to 26 Hz (attenuation relative to 1 Hz not exceeding 3 dB).
External pressure input	Pressure range	0 to 300 mmHg
	Accuracy	±4 mmHg or 4% (whichever is greater)

8.3.2. Pressure Waveform Display

Parameters	Requirements
Automatic scaling	Arterial pressure waveforms are automatically scaled and displayed in the waveform window; within the scaled vertical range, numerical pressure scale values are displayed at the lowest, middle, and highest grid lines
Reference lines	−30 to 300 mmHg

8.3.3. Digital Pressure Display

Parameters	Requirements
Resolution	1mmHg
Pressure range	-30~300mmHg
Accuracy	±4 mmHg or 4% (whichever is greater)
Update rate	2 seconds (nominal)

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Hydraulic sensor	Zeroing shall be possible for offsets not exceeding ±150 mmHg, with waveform display offset within ±2 pixels (±0.4 mm) relative to the zero baseline
Fiber-optic pressure sensor	Zero adjustable, with waveform display offset within ±12 pixels (±2.2 mm) relative to the zero baseline
Automatic zeroing time	< 3 seconds

8.3.5. Pressure Output Signal

Parameters	Requirements
Technical parameters	0~300mmHg
Accuracy	±4% or ±4 mmHg (whichever is greater)
Scale factor	10.0 mV/mmHg ±4%
Frequency response	Hydraulic sensor: DC to 42 Hz ±15% (attenuation relative to 1 Hz not exceeding 3 dB); Fiber-optic sensor: DC to 26 Hz ±15% (attenuation relative to 1 Hz not exceeding 3 dB)

8.4. Heart Rate Meter

A heart rate indicator is displayed on the monitor screen, showing the current average heart rate in numerical form, in beats per minute (bpm). When no signal is present or the signal is invalid, "---" is displayed.

Range: 15 to 220 bpm;

Resolution (display): 1 bpm;

Accuracy: ±3% or ±3 bpm, whichever is greater;

Step response time: ≤10 seconds.

8.5. Augmentation Function

Parameters	Requirements
Operating modes	Automatic mode: trigger sources include ECG trigger, pressure trigger; Semi-automatic mode: trigger sources include ECG trigger, pressure trigger, Pacer A trigger, Pacer V/AV trigger, internal trigger
Augmentation ratio	Adjustable 1:1, 1:2, 1:3
Augmentation pressure level	Adjustable from level 1 to 10, default level 10
Inflation timing adjustment	Range −4.0 to 4.0, increment 0.1
Deflation timing adjustment	Range −4.0 to 4.0, increment 0.1
Electrosurgical interference suppression (ESIS)	After exposure to the electromagnetic field of high-frequency surgical equipment, the IABP can return to its original operating mode within 10 seconds without loss of stored data.

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8.6. Alarm Sound Level

Alarm sound pressure level: adjustable from greater than 35 dBA to less than 85 dBA.

8.7. Display

Display specification: 15.6-inch touch screen display.

Resolution: 1920 × 1080.

Viewing angle: 89° horizontal / 89° vertical.

8.8. Report

Report size: A4 format.

Recording speed: selectable 25 mm/s or 50 mm/s.

8.9. Power Supply

8.9.1. Mains Power

Input voltage: 220 V ±10%.

Input frequency: 50 Hz ±2 Hz.

Power consumption: 400 VA.

8.9.2. Internal Battery

Provides continuous power to the IABP when mains power is disconnected, suitable for emergency scenarios such as patient transport between wards or operating rooms, or sudden mains power failure.

Battery type: Lithium iron phosphate ($LiFePO_4$).

Operating voltage: Typical 25.6 V, minimum 20 V, maximum 29.2 V.

Rated capacity: 12 Ah.

Charging time: 6 hours to reach at least 90% charge.

Operating time: at 120 bpm and ambient temperature of 22 ±5 °C, a new battery operates for a minimum of 180 minutes.

The internal battery is used when mains power is disconnected, for emergency use such as patient transfer between wards/operating rooms or sudden mains power failure.

Note: when operating on internal battery power, battery charge level shall be monitored to prevent insufficient power.

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Working pressure	15 MPa
Standard	DOT-3AL
Durability	When fully charged with helium, supports continuous operation for 3 months (inflation once every 2 hours, continuous operation 24 hours per day)

8.11. ECG Lead Wires

Integrated 5-lead cable, reusable

Total length 3020 mm;

Snap-type electrode connection;

AHA lead color coding;

Right arm (white, RA);

Left arm (black, LA);

Left leg (red, LL);

Chest (brown, V);

Right leg (green, RL).

IEC lead color coding:

Right arm (red, R),

Left arm (yellow, L),

Left leg (green, F),

Chest (white, C),

Right leg (black, N).

Reusable, non-sterilizable.

8.12. Invasive Blood Pressure Sensor Cable

Total length: 3250 mm.

Connector specification: RJ11 connector at the sensor end.

Reusable, non-sterilizable.

8.13. Display Cable

Total length: 1 m.

8.14. Power Cord

Plug type: three-flat-pin straight plug.

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8.15. Electromagnetic Compatibility

Lifetech Cardio IABP complies with YY 9706.102-2021, Chapter 202 of YY 9706.234-2021, and Chapter 202 of GB 9706.227-2021.

The device shall be used in accordance with specific EMC precautions.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
Lifetech Cardio IABP (model IAP-50) is intended for use in the electromagnetic environment specified below. The user should ensure it is used in such an environment:		
Emission test	Compliance	Electromagnetic environment - guide
RF emissions GB 4824	Group 1	The device uses RF energy only for internal functions. Its RF emissions are low and unlikely to interfere with nearby electronic equipment.
RF emissions GB 4824	Class A	The device is suitable for use in all establishments except those di-rectly connected to the public low-voltage power supply network.
Harmonic emissions GB 17625.1	Not applicable	
Voltage fluctuations/flicker emissions GB 17625.2	Not applicable	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
Lifetech Cardio IABP (model IAP-50) is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Assessment
Electrostatic discharge (GB/T 17626.2)	±6 kV contact discharge and ±8 kV air discharge	±6 kV contact discharge and ±8 kV air discharge	Floors should be wood, con-crete, or ce-ramic tile. If covered with synthetic mate-rial, relative humidity should be at least 30%
Electrical fast transient/burst (GB/T 17626.4)	±2 kV for power supply lines and ±1 kV for input/output lines	±2 kV for power supply lines and ±1 kV for input/output lines	Mains power quality shall be that of a commercial or hospital environment.

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Surge (GB/T 17626.5)	±1 kV line-to-line and ±2 kV line-to-earth	±2 kV line-to-earth ±1 kV line-to-line and	Mains power quality shall be that of a commercial or hospi-tal environment
Voltage dips, short interruptions, and voltage variations on power supply input lines (GB/T 17626.11)	<5%U_T (>95% dip in U_T) for 0.5 cycle 40%U_T （60% dip in U_T）for 5 cycle 70%U_T （30% dip in U_T）for 25 cycle <5%U_T (>95% dip in U_T) for 5 s	<5%U_T (>95% dip in U_T) for 0.5 cycle 40%U_T （60% dip in U_T）for 5 cycle 70%U_T （30% dip in U_T）for 25 cycle <5%U_T (>95% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Model IAP-50 image intensifier requires continued operation during power mains interruptions, it is rec-ommended that the Model IAP-50 Image Intensifier be pow-ered from an uninter-ruptible power supply or a battery.
Power frequency magnetic field (GB/T 17626.8)	3 A/m	3 A/m	Power fre-quency mag-netic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF dis-turbances (GB/T 17626.6-2017) Radiated RF dis-turbances (GB/T 17626.3)	3 V (rms) 150 kHz to 80 MHz (excluding ISM bands) 10 V/m (rms) 150 kHz to 80 MHz (for ISM bands) 10 V/m 80MHz to 2.5 GHz	3 V (rms) 10 V/m rms 10V/m	Portable and mobile RF communica-tions equipment should be used no closer to any part of the 6108 Tem-porary Pacemaker, in-cluding ca-bles, than the recommended separation dis-tance calculat-ed from the applicable for-mula. Recommended sepa-ration dis-tance: d = [1.17]√Pd = [1.20]√Pd = [1.20]√P (80 MHz – 800 MHz)d = [2.30]√P (800 MHz – 2.5 GHz)

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			Where: P = max rated output power of the transmitter in watts (W)d = recommended separation dis-tance in meters (m) The field strengths of fixed RF trans-mitters should be determined by an electro-mag-netic site survey and should be less than the appli-cable com-pli-ance level in each frequency range. Interference may oc-cur in the vicinity of equipment marked with the following sym-bol. 

Note 1: At 80 MHz and 800 MHz, the separation distance formula for the higher frequency range applies.

Note 2: These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

[a] The ISM frequency bands between 150 kHz and 80 MHz refer to: 6.765–6.795 MHz, 13.553–13.567 MHz, 26.957–27.283 MHz, and 40.66–40.70 MHz

[b] For ISM bands from 150 kHz to 80 MHz and frequencies from 80 MHz to 2.5 GHz, compli-ance levels are designed to reduce the likelihood of interference from port-able/mobile RF communications equipment accidentally brought into the patient area. A 10/3 factor is applied to the separation distance in these cases.

[c] Fixed transmitters, such as base stations for radio (cellular/cordless), land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, may not be pre-dictable. A site survey should be considered to assess the electromagnetic environment.

If the measured field strength at the location exceeds the applicable RF compliance level, the device should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as reorient-ing or relocating the device..

[d] Across the 150 kHz–80 MHz frequency range, field strength should be less than

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3 V/m.

Recommended separation distances between portable/mobile RF communi-cations equipment and the Lifetech Cardio Intra-Aortic Balloon Pump				
Lifetech Cardio IABP (model IAP-50) is intended for use in an RF-controlled electromag-netic environ-ment. The user can help prevent electromagnetic interference by main-taining the recommended minimum distance between the device and any RF transmit-ter ac-cording to its output power				
Maximum Rated Output Power of Transmit-ter (W)	Corresponding Separation Distances by Frequency (m)			
	150 kHz～80 MHz (excluding ISM bands) d=[1.17]$\sqrt{P}$	150 kHz～80 MHz (excluding ISM bands) d=[1.20] $\sqrt{P}$	80MHz～ 800MHz d=[1.20] $\sqrt{P}$	800MHz～2.5 GHz d=[2.30] $\sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.37	0.38	0.38	0.73
1	1.17	1.20	1.20	2.30
10	3.70	3.79	3.79	7.27
100	11.70	12.00	12.00	23.00
For transmitters not listed above, the recommended separation distance d (in me-ters) can be calculated using the appropriate formula, where P is the transmitter's maximum rated output power in watts (W). Note 1: At 80 MHz and 800 MHz, use the formula for the higher frequency band. Note 2: The ISM bands from 150 kHz to 80 MHz include: 6.765–6.795 MHz, 13.553–13.567 MHz, 26.957–27.283 MHz, and 40.66–40.70 MHz. Note 3: A 10/3 factor is used in the calculation of separation distance in ISM bands (150 kHz–80 MHz) and 80 MHz–2.5 GHz, to reduce the risk of interference in patient areas. Note 4: These guidelines may not apply in all situations. Electromagnetic propaga-tion is affected by absorption and reflection from structures, objects, and people.				

8.16. Equipotential Bonding Terminal

During operation, leakage currents to earth generated by the device flow into the protective earth conductor of the power supply installation, resulting in different potentials at different points along the protective earth conductor. For example, protective earth points at different outlets within the same room may have different potentials. The potential difference depends on the current flowing between the two points and the impedance of the protective earth conductor between them. If more than one device is used in the room, a potential difference may exist between the grounded enclosures of these devices, which may pose a safety risk to personnel in

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the vicinity. If these devices are equipped with equipotential bonding terminals (typically connected to the protective earth of the device), connecting each device's equipotential bonding terminal to the equipotential bonding busbar of the power supply installation via equipotential bonding conductors can reduce the impedance between enclosures and thereby reduce the potential difference. This is the purpose of equipotential bonding.

8.17. Specifications of the Fuse Device

The IABP contains two fuses, rated at 5×20 mm / T5AH / 250 V.

8.18. Main Safety Characteristics of the Device

Classification by electric shock protection type: Class I equipment when connected to mains power; internally powered equipment when operated using the internal battery.

Classification by degree of protection against electric shock: Type CF applied parts.

Classification by degree of protection against liquid ingress: IPX0.

Classification by safety when used with flammable anesthetic gases mixed with air or with flammable an-esthetic gases mixed with oxygen or nitrous oxide: Not applicable

Classification by mode of operation: continuous operation.

Rated voltage and frequency: 220 V, 50 Hz.

Input power: 400 VA.

Application section with protection against defibrillation discharge: CF type application

Signal output or input section: Yes

Non-permanently installed equipment.

Description: when used in combination with the Intra-Aortic Balloon Catheter, the ECG Lead Wires, Invasive Blood Pressure Sensor Cable, and IAB catheter constitute the applied parts of the device.

8.19. Software Release Version

The software release version of the Lifetech Cardio Intra-Aortic Balloon Pump is version 1.

8.20. Standards Compliance

GB 9706.1-2020 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance

YY 9706.102-2021 Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic

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compatibility – Requirements and tests.

YY 9706.234-2021[1] Medical electrical equipment – Part 2-34: Particular requirements for the basic safety and essential performance of invasive blood pressure monitoring equipment

YY 9706.108-2021 Medical electrical equipment – Part 1-8: General requirements for basic safety and essential performance – Collateral standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems.

[1] Testing was performed using a disposable pressure transducer manufactured by Leister (model DPT-248).。

8.21. Essential performance

1）Defibrillation protection

2）Interruption of the power supply / SUPPLY MAINS to ME EQUIPMENT

3）Protection against depletion of battery

4）Essential performance of ME equipment as specified in 201.12.101 of GB 9706.227-2021*

5）Electrosurgical interference

6）Accuracy of pressure measurements

7）Detection of TRANSDUCER and TRANSDUCER cable fault

8）Detection of disconnected catheter

9）Counterpulsation based on input signals and user settings, or activation of a no-trigger alarm

* Within the essential performance requirements of ME equipment in 201.12.101 of GB 9706.227-2021, the following clauses are not applicable to this product: “201.12.1.101.5 Multi-channel crosstalk, 201.12.1.101.6 Gain control and stability, 201.12.1.101.14 Synchronized pulse for cardioversion, and Aspect ratio in 201.12.1.101.16”.

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